




### Low pH Stability with HALO® LPH-C18

294



#### TEST CONDITIONS:

**Column:** HALO 90 Å LPH-C18, 2.7 µm 2.1 x 50 mm

**Part Number:** 92822-416

**Mobile Phase A:** Water, 1% TFA (pH: 1)

**Mobile Phase B:** Acetonitrile

Gradient:	Time	%B
	0.0	20
	7.50	20
	7.51	5
	45.00	5
	47.00	100
	51.00	100
	51.01	20
	60.00	20

**Flow Rate:** 0.5 mL/min

**Pressure:** 108 bar

**Temperature:** 60 °C

**Detection:** UV 254 nm, PDA

**Injection Volume:** 0.4 µL

**Sample Solvent:** 25/75 ACN/ Water


**Data Rate:** 100 Hz

**Response Time:** 0.025 sec.

**Flow Cell:** 1 µl

**LC System:** Shimadzu Nexera X2

A separation of parabens is performed on a HALO 90 Å LPH-C18 column under low pH (pH 1) and high temperature conditions compared to a standard C18 SPP column. Due to the sterically protected ligand, the LPH-C18 column can withstand these conditions and maintain stable retention times while other columns show decreased retention over time indicating a loss of stationary phase.



#### AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17• D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de

