


# ACE<sup>®</sup>

UHPLC and HPLC Columns

## Food and Beverage Applications Guide



**ACE<sup>®</sup>**  
UHPLC and HPLC Columns



## Ultra-Inert Base Deactivated UHPLC / HPLC Columns

### Contents

Application Index	1
Analyte Index	2 - 4
Application Notes	5 - 90

### Send us your application and receive a free ACE column

Send us your application on an ACE column and help extend our applications database. Your proven method will enable your chromatography colleagues to benefit and if we select your application for publication we'll send you a **FREE** ACE analytical column of your choice.

To submit your application e-mail us at: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# ACE<sup>®</sup> Food and Beverage Applications: Application Index

Application	Pages	Application	Pages
Additives and intense sweeteners	5	Organic acids 2	47
Argicultural pesticides	6	Organic acids 3	48
Amino acids in peas	7	Organic acids 4	49
Amino acids and biogenic amines in wine and beer	8	Organophosphorus flame retardants in water by LC-MS/MS	50
Aminoglycosides in eggs	9	Organophosphorus isomer flame retardants in water	51
Annatto	10	Paraben preservatives	52
Anthocyanins from sambucus nigra (elderberry)	11	Perfluoro acids by LC-MS/MS	53
Appetite suppressants by LC-MS	12	Perfluoroalkyl substances by ion pairing LC-MS/MS	54
Arsenolipids from edible seaweed	13	Perfluorinated compounds in water by LC-MS/MS	55
Artificial colours (water soluble)	14	Pesticides (250 analytes) by LC-MS/MS	56
Artificial food colouring	15	Pesticides (47 analytes) by LC-MS/MS	60
Artificial sweeteners global method	16	Pesticides in water	61
Artificial sweeteners (stevia glycosides)	17	Phenolic compounds in ground water and landfill leachates	62
Avenacins	18	Phenols and phenoxy acid herbicides	63
Benfluralin (herbicide)	19	Phenols in purple coneflower	64
Caffeine and metabolites	20	Phytoestrogens in hop extracts	65
Chloramphenicol in milk	21	Preservatives 1	66
Chocolate analysis	22	Preservatives 2	67
Coffee metabolite profiling	23	Senosides in traditional Chinese medicine	68
Cyclodextrin encapsulated flavour compounds in beer	24	Spice analysis for illegal dyes	69
Defensins (human) in saliva matrix	25	Sulphonamides	70
Ethanol extract from seed cover (Acacia Farnesiana)	26	Tocopherols	71
Flavone and dibucaine	27	Trifluralin (Herbicide)	72
Flavonoids	28	Tryptic Digest of BSA	73
Garlic analysis 1	29	Vanillins 1	74
Garlic analysis 2	30	Vanillins 2	75
Green tea metabolite profiling	31	Veterinary steroids by LC-MS/MS	76
Herbicide impurity profile	32	Vitamin D2/D3	77
Hippuric acid	33	Vitamin D (25-hydroxy) by LC-MS/MS	78
Insulins	34	Vitamins -fat soluble	79
Isoflavones	35	Vitamin in fruit juice by fast LC-MS	80
Isoflavones in red clover and soy extract	36	Vitamins –water soluble (gradient 1)	81
Liquorice extracts	37	Vitamins –water soluble (gradient 2)	82
Malachite green	38	Vitamins –water soluble (isocratic 1)	83
Melamine	39	Vitamins –water soluble (isocratic 2)	84
Microcystins from blue/green algae in drinking water	40	Vitamins –Water soluble	85
Milk proteins	41	Vitamins –Water soluble	86
Mycotoxins by LC-MS/MS	42	Water soluble vitamins in green vegetables by LC-MS/MS	87
Mycotoxins/aflatoxins from peppers	43	Water soluble vitamins/polar molecules	88
Nitrofurans metabolites by LC-MS/MS	44	Whey proteins from whole milk	89
Ochratoxin A	45	Wine acids	90
Organic acids 1	46		

# ACE® Food and Beverage Application: Analyte Index

Analyte	Pages	Analyte	Pages	Analyte	Pages
Abamectin	56	n-Butyl paraben	52	Demeton S-methylsulfone	6, 60
Acetamidiprid	6, 57, 60	Cadaverine	8	Demeton S-methylsulfoxide	6, 60
Acephate	6, 57, 60	Cadusafos	57	Demeton S-sulfone	57
Acesulfame K	5, 16, 22	Caffeic acid	49	Deoxyalliin	29
Acetyl-daizdin	36	Caffeine	5, 12, 20, 22	Desethyltriazine	61
Acetyl-genistin	36	Caftaric acid	64	N-Desmethylsibutramine	12
Acetyl-glycitin	36	Carbaryl	6, 57, 60	Dialifos	57
Aflatoxin G2	42, 43	Carbendazim	6, 57, 60	Diazinon	57
Aflatoxin G1	42, 43	Carbofuran	6, 57, 60	Diazinon OA	57
Aflatoxin B2	42, 43	Carboxin	57	Dibucaine	27
Aflatoxin B1	42, 43	Carfentrazone ethyl	57	Dichlormid	57
α-Alanine	8	Carmoisine	14	2,4-Dichlorophenol (2,4-DCP)	24, 63
β-Alanine	8	α-Caseins	41	2,6-Dichlorophenol	63
Agmatine	8	β-Caseins	41	2,4-Dichlorophenoxyacetic acid (2,4-D)	63
Ajoene	30	Catechin	22	4(2,4-Dichlorophenoxy)butyric acid (2,4-DB)	63
Alachlor	61	Chloramphenicol	21	Dichlorprop (2,4-DP)	63
Aldicarb	6, 57, 60	Chlorantranilprole	57	Dichlorvos	57
Aldicarb sulfone	6, 57, 60	Chlordiazepoxide	12	Dicrotophos	6, 57, 60
Aldicarb sulfoxide	6, 57, 60	Chlorfenvinphos	57	Didesmethylsibutramine	12
Alitame	16	Chlorogenic acid	22, 64	Dienestrol	76
Allethrin	57	Chlorimuron ethyl	57	Diethofencarb	57
Alliin	29, 30	4-Chloro-2-methylphenoxyacetic acid (MCPA)	63	Diethylpropion	12
Alliin	29	4-Chloro-2-methylphenoxybutyric acid (MCPB)	63	Diethylstilbestrol	76
Allura Red	14	2-Chlorophenol	63	Diflubenzuron	57
Amaranth	14, 15	4-Chlorophenol	63	Dihydrostreptomycin	9
Ametoctradin	57	Chlorpyrifos	57	Dimethenamid	57
L-2-Aminoadipic acid	8	Chlorpyrifos OA	57	Dimethoate	6, 57, 60
p-Aminobenzoic acid	81, 82, 85, 86, 88	Cholecalciferol (vitamin D3)	77, 79	Dimethomorph	6, 57, 60
Ammonium ion	8	Cichoric acid	64	3,4-Dimethoxycinnamic acid	47, 48
Arginine	7, 8, 40	Cinnamic acid	47, 48	Dinotefuran	6, 57, 60
L-Ascorbic acid (vitamin C)	46, 80, 81, 82, 83, 84, 87, 88, 90	Citric acid	80, 90	Dioxacarb	57
Asparagine	7, 8	Clethodim	57	Dioxathion	57
Aspartame	5, 16	Clofentezine	6, 57, 60	Disulfoton sulfone	57
Aspartic acid	8	Cloransulam methyl	57	Disulfoton sulfoxide	57
Atrazine	57, 61	Clothianidin	6, 57, 60	Diuron	57
Avenacin A-1	18	CMPP	63	DMST (dimethylsulftoluidide)	57
Avenacin A-2	18	Coumaphos	57	Dodine	57
Avenacin B-1	18	p-Coumaric acid	49	Dulcin	16
Avenacin B-2	18	o-Cresol	62, 63	Dulcoside A	17
Azinphos ethyl	57	Crotoxyphos	57	Emamectin	57
Azinphos methyl	57	Crufomate	57	Ephedrine	12
Azinphos methyl OA	57	Cyanazine	57	Epicatechin	22
Azoxystrobin	57	Cyanidin-3,5-diglucoside	61	Ergocalciferol (vitamin D2)	77, 79
Bendiocarb	57	Cyanidin-3-glucoside	11	Erythrosine	14
Benfluralin	19	Cyanidin-3-sambubioside	11	Ethiofencarb	57
Benomyl	6, 60	Cyanidin-3-sambubioside-5-glucoside	11	Ethiofencarb sulfone	57
Benoxacor	57	Cyanocobalamin (vitamin B12)	80, 81, 82, 85, 86	Ethiofencarb sulfoxide	57
Benzoic acid	46, 67	Cyantranilprole	57	Ethion	57
Benzyl paraben	52	Cyazofamid	57	Ethion monoxon	57
Bifenazate	57	Cyclamate	5, 16	Ethiprole	57
Biochanin A	36	Cyflufenamid	57	Ethofumesate	57
d-Biotin (vitamin B7)	80, 81, 82, 85, 86	Cyfluthrin	6, 60	Ethoprop	57
Bis(2-ethylhexyl) phosphate (BEHP)	50	Cymoxanil	57	Ethyl paraben	22, 52, 66, 67
Bis(1,3-dichloro-2-propyl) phosphate (BDPC)	50	Cyphenothrin	57	4-Ethylphenol	24
Bitertanol	57	Cyproconazoles	6, 60	Ethyl vanillin	74, 75
9'-cis-Bixin	10	Cyprosulfamide	57	Etofenprox	57
9'-trans-Bixin	10	Cyromazine	57	Etoxazole	57
α-Boldenone	76	Cysteine	8	Eugenol	74, 75
β-Boldenone	76	Daidzein	35, 36	Famoxadone	57
Boscalid	57	Daidzin	35, 36	Fenamidon	57
2-Bromophenol	24	DEF	57	Fenamiphos	57
Bupirimate	57	Defensin HNP-1	25	Fenamiphos sulfone	57
Buprofezin	57	Defensin HNP-2	25	Fenamiphos sulfoxide	57
Butter Yellow	69	Defensin HNP-3	25	Fenazaquin	57
i-Butyl paraben	52	Deisopropyltriazine	61	Fenbuconazole	6, 60
				Fenfluramine	12

# ACE® Food and Beverage Application: Analyte Index

Analyte	Pages	Analyte	Pages	Analyte	Pages
Fenhexamid	57	Imazosulfuron	58	Microcystin YR	40
Fenobucarb	57	Imidacloprid	6, 58, 60	Monocrotophos	6, 58, 60
Fenoxaprop ethyl	57	Imiprothrin	58	Monolinuron	58
Fenoxycarb	57	Indaziflam	58	Myclobutanil	58
Fenpropimorph	57	Indoxacarb	6, 58, 60	Neohesperidin dihydrochalcone	5, 16
Fenpyroximate	58	Ipconazole	58	Neotame	16
Fensulfthion	58	Iprovalicarb	58	Nicosulfuron	58
Fenuron	58	Isoamylamine	8	Nicotinamide	80, 87, 88
Ferulic acid	47, 48, 49	Isofenphos	58	Nicotine	6, 60
Flavone	27	Isoleucine	8	Nicotinic acid	80, 81, 82, 83, 84
Flonicamid	58	Isonicotinamide	88	Nitrofurazone	44
Flubendiamide	6, 60	Isoproc carb	58	Norbixin	10
Fluazifop P butyl	58	2-Isopropylphenol	24	Norflurazon	58
Fludioxonil	58	Isoproturon	58	Norflurazon desmethyl	58
Flufenoxuron	58	Isoxanthohumol	65	α-Nortestosterone	76
Flufenpyr ethyl	58	Kaempferol	28	β-Nortestosterone	76
Flumetsulam	58	Kanamycin	9	Ochratoxin A	42, 45
Flumiclorac pentyl	58	Kresoxim methyl	58	Omethoate	6, 58, 60
Fluometuron	58	α-Lactalbumin	89	Orange II	69
Fluopicolide	58	Lactic acid	90	Orlistat	12
Fluopyram	58	Lactofen	58	Ornithine	8
Fluoxastrobin	58	β-Lactoglobulin A	41, 89	Orotic acid	88
Fluoxetine	12	β-Lactoglobulin B	41, 89	Oxalic acid	90
Fluridone	58	Lenacil	58	Oxamyl	6, 58, 60
Flusilazole	58	Leptophos OA	58	Oxamyl oxime	58
Fluthiacet methyl	58	Leucine	8, 40	Oxydemeton methyl	58
Flutolanil	58	Leucomalachite green	38	Oxydemeton methyl sulfone	58
Flutriafol	58	Linuron	58	Pantothenic acid (vitamin B5)	80, 81, 82, 85, 86, 87
Fluxapyroxad	58	Lorcaserin	12	Para Red	69
Folic acid (vitamin B9)	81, 82, 85, 86, 87	Lysine	7, 8	Parathion methyl OA	58
Folpet	6, 60	Malachite green	38	Parathion OA	58
Forchlorfenuron	58	Malathion	58	Paraxanthine	20
Formetanate	6, 58, 60	Malathion OA	58	Pencycuron	6, 58, 60
Formononetin	36	Malic acid	90	Pendimethalin	61
Fosthiazate	58	Mandipropamid	6, 58, 60	Penflufen	58
Furaldalone	44	MCPA (4-Chloro-2-methylphenoxyacetic acid)	63	Penthiopyrad	58
Furazolidone	44	MCPB (4-Chloro-2-methylphenoxybutyric acid)	63	Perfluorobutanoic acid (PFBA)	54, 55
GABA	8	Mefenpyr diethyl	58	Perfluorobutylsulfonic acid (PFBS)	53, 54
Genistein	35, 36	Melamine	39	Perfluorodecanoic acid (PFDA)	53, 54, 55
Genistin	35, 36	Menadione (vitamin K3)	79	Perfluoroheptanoic acid (PFHpA)	53, 54
Gentamicin	9	Mepanipyrim	58	Perfluorohexanoic acid (PFHxA)	53, 54, 55
Glutamic acid	8	Mesotrione	58	Perfluorohexylsulfonic acid (PFHxS)	53, 54, 55
Glutamine	8	Metaflumizone	58	Perfluorononanoic acid (PFNA)	53, 54, 55
Glycine	8	Metalaxyl	58	Perfluorooctanesulfonic acid (PFOS)	53, 54, 55
Glycitein	36	Metaldehyde	58	Perfluorooctanoic acid (PFOA)	53, 54, 55
Glycitin	36	Metanil Yellow	69	Perfluoropentanoic acid (PFPeA)	54
Guaiaicol	22, 74, 75	Metconazole	58	Phenol	62, 63
Heptafluorobutyric acid	53	Methamidophos	6, 58, 60	Phenothrin	58
Hexaconazole	6, 58, 60	Methidathion	58	Phentermine	12
Hexaflumuron	6, 60	Methiocarb	58	Phenthoate	58
Hexestrol	76	Methiocarb sulfone	58	Phenylalanine	8
Hexythiazox	58	Methiocarb sulfoxide	58	Phenylethylamine	8
Hippuric acid	33, 87	Methionine	8	Phenolphthalein	12
Histamine	8	Methomyl	6, 58, 60	Phorate OA	58
Histidine	7, 8	4-Methoxycinnamic acid	47, 48	Phorate OA sulfone	58
HT-2-toxin	42	Methoxyfenozide	58	Phorate OA sulfoxide	58
4-Hydroxybenzaldehyde	74, 75	2-Methoxy-4-vinylphenol	24	Phorate sulfone	58
4-Hydroxybenzoic acid	22, 66, 67, 74, 75	Methyl paraben	22, 52, 66, 67	Phorate sulfoxide	58
3-Hydroxycarbofuran	57	2-Methylphenol	24	Phosalone	58
Hydroxyproline	8	Methyltestosterone	76	Phosmet	58
Hydroxystanozolol	76	Metolcarb	58	Phosmet OA	58
5-Hydroxythiabenzazole	58	Metribuzin	58	Phosphamidon	58
25-Hydroxy vitamin D2	78	Mevinphos	58	Phoxim	58
25-Hydroxy vitamin D3	78	Microcystin RR	40	Phthalic acid	66, 67
Hygromycin-B	9	Microcystin LR	40	Picoxystrobin	58
Hypoxanthine	8			Pirimicarb	59

# ACE® Food and Beverage Application: Analyte Index

Analyte	Pages	Analyte	Pages	Analyte	Pages
Pirimicarb desmethyl	59	Spermidine	8	γ-Tocopherol	71
Pirimiphos methyl	59	Spinetoram	59	δ-Tocopherol	71
Ponceau 4R	14, 15	Spinosad A	59	Topramezone	59
Prallethrin	59	Spinosad D	59	Toluidine Red	69
Prochloraz	6, 59, 60	Spiroclorfen	59	T-2-toxin	42
Profoxydim	59	Spiromesifen	59	Tri-m-tolyl	51
Proline	8	Spiromesifen alcohol	59	Tri-o-tolyl	51
Promecarb	59	Spirotetramat	59	Tri-p-tolyl	51
Prometryn	61	Spiroxamine	59	Triadimefon	59
Propamocarb	59	Steviolbioside	17	Triadimenol	59
Propaquizafop	59	Stevioside	17	Triazophos	59
Propargite	6, 59, 60	Streptomycin	9	Trenbolone	76
Propazine	61	Sucralose	5, 16	Triethyl phosphate	50
Propetamphos	59	Sudan I	69	Tribenuron methyl	59
Propoxur	59	Sudan II	69	Trichlorfon	59
n-Propyl paraben	22, 52, 66, 67	Sudan III	69	Tricyclazole	59
Prosulfuron	59	Sudan IV	69	Trifloxystrobin	59
Putrescine	8	Sudan Black	69	Triflumizole	59
Pymetrozine	59	Sudan Orange G	69	Trifluralin	72
Pyraclostrobin	59	Sudan Red 7B	69	Triforine	6, 59, 60
Pyraflufen ethyl	59	Sudan Red B	59	Triisopropyl phosphate (TIPP)	50, 51
Pyrazophos	59	Sudan Red G	69	2,4,6-Trimethylphenol	24, 62
Pyridaben	59	Sulfachloropyridazine	70	Trimethyl phosphate (TMP)	50
Pyridalyl	59	Sulfadiazine	70	Tri-n-butyl phosphate (TBP)	50
Pyridoxal (vitamin B6)	81, 82, 83, 84, 88	Sulfadimethoxine	70	Tri-n-propyl phosphate (TPPrP)	50, 51
Pyridoxamine (vitamin B6)	81, 82, 83, 84, 88	Sulfamerazine	70	Triphenyl phosphate (TPP)	50
Pyridoxine (vitamin B6)	80, 81, 82, 83, 84, 85, 86, 87, 88	Sulfamethazine	70	Tris(2-butoxyethyl) phosphate (TBEP)	50
		Sulfamethoxazole	70	Tris(2-chloroethyl) phosphate (TCEP)	50
Pyrimethanil	59	Sulfamonomethoxine	70	Tris((2R)-1-chloro-2-propyl) phosphate (TCPP)	50
Pyriproxyfen	59	Sulfamoxole	70	Tris(1,3-dichloro-2-propyl) phosphate (TDCPP)	50
Quercetin	22, 28	Sulfapyridine	70		
Quercitrin	28	Sulfoxaflor	59	Tryptamine	8
Quinalphos	59	Sulprofos	59	Tryptophan	8
Quinine	5	Sunset Yellow	14, 15	Tyramine	8
Quinoxifen	59	Taleranol	76	Tyrosine	8, 40
Quizalofop ethyl	59	Tartaric acid	90	Valine	8
Rebaudioside A	17	Tartrazine	15	Vanillic acid	74, 75
Rebaudioside B	17	TCMTB (Busan)	59	o-Vanillin	75
Rebaudioside C	17	Tebufenozide	59	Vanillin	22, 74, 75
Rebaudioside D	17	Tebufenpyrad	59	Vitamin A	79
Rebaudioside F	17	Tebuthiuron	59	Vitamin A acetate	79
Red 2G	14	Tepaloxymdim	59	Vitamin B1 (thiamine)	80, 81, 82, 83, 84, 87, 88
Resmethrin	59	Terbufos sulfone	59	Vitamin B2 (riboflavin)	80, 81, 82, 85, 86, 87
Retinol (vitamin A)	79	Terbufos sulfoxide	59	Vitamin B5 (pantothenic acid)	80, 81, 82, 85, 86, 87
Rhodamine B	69	Terbutylazine	59, 61	Vitamin B6	80, 81, 82, 83, 84, 85, 86, 87, 88
Riboflavin (vitamin B2)	80, 81, 82, 85, 86, 87	Terbutryn	61	Vitamin B7 (d-biotin)	80, 81, 82, 85, 86
Rimonabant	12	Tetrachlorvinphos	59	Vitamin B9 (folic acid)	81, 82, 85, 86, 87
Rimsulfuron	59	Tetramethrin	59	Vitamin B12 (cyanocobalamin)	80, 81, 82, 85, 86
Rotenone	59	Theobromine	5, 20, 22	Vitamin C (ascorbic acid)	46, 80, 81, 82, 83, 84, 87, 88, 90
Rubusoside	17	Theophylline	5, 20, 22	Vitamin D2 (ergocalciferol)	77, 79
Rutin	28, 36	Thiabendazole	6, 59, 60	Vitamin D3 (cholecalciferol)	77, 79
Saccharin	5, 16, 22	Thiacloprid	6, 59, 60	Vitamin E (tocopherol)	71, 79
Saffluenacil	59	Thiamethoxam	6, 59, 60	Vitamin E acetate	79
Sebuthylazine	61	Thiamine (vitamin B1)	80, 81, 82, 83, 84, 87, 88	Vitamin K1	79
Sedaxane	59	Thifensulfuron methyl	59	Vitamin K3	79
Senenoside A	68	Thiobencarb	59	Xanthine	20
Senenoside B	68	Thioctic acid	81	Xanthohumol	65
Serine	8	Thiodicarb	6, 59, 60	α-Zearalenol	76
Sertraline	12	Thionazin	59	β-Zearalenol	76
Sethoxydim	59	Thiophanate methyl	6, 60	Zearalenone	76
Sibutramine	12	Threonine	8	Zeranol	76
Simazine	59, 61	Tobramycin	9		
Sinapic acid	47, 48	α-Tocopherol	71, 79		
Sorbic acid	22, 66, 67	β-Tocopherol	71		
Spectinomycin	9				

# Separation of Additives and Intense Sweeteners


Application #AN2950

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.0 mm  
Part Number: ACE-121-2504  
Mobile Phase: A: H<sub>2</sub>O  
B: MeCN  
C: 1% TFA in H<sub>2</sub>O

Time (mins)	%A	%B	%C
0	88	2	10
25	50	40	10
30	30	60	10
35	88	2	10

Flow Rate: 1 mL/min  
Temperature: 30 °C  
Detection: ELSD



Reproduced with permission of Durham County Council, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)




## Conditions


Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 50 x 2.1 mm  
Part Number: CORE-25A-0502U  
Mobile Phase: A: 0.1% formic acid + 5 mM ammonium formate in H<sub>2</sub>O/MeOH (90:10 v/v)  
B: 0.1% formic acid + 5 mM ammonium formate in H<sub>2</sub>O/MeOH (10:90 v/v)

Time (mins)	%B
0.00	0
1.00	0
15.00	100
18.00	100
18.05	0
20.00	0


Flow Rate: 0.4 mL/min  
Injection: 20 µL  
Temperature: 40 °C  
Detection: Agilent 6420 Triple Quadrupole MS, +ve mode ESI, Dynamic MRM




1. Carbaryl  
(*m/z* 202.10 → 145.10)




2. Dimethomorphs  
(*m/z* 388.10 → 301.10)



3. Cyproconazoles  
(*m/z* 292.10 → 70.00)



4. Fenbuconazole  
(*m/z* 337.10 → 70.00)



Also analysed under same conditions:

Acephate, Acetamiprid, Aldicarb, Aldicarb sulfone, Aldicarb sulfoxide, Benomyl, Carbendazim, Carbofuran, Clofentezine, Clothianidin, Cyfluthrin, Demeton S-methylsulfone, Demeton S-methylsulfoxide, Dicrotophos, Dimethoate, Dinotefuran, DMA, DMPF, Flubendiamide, Folpet, Formetanate, Hexaconazole, Hexaflumuron, Imidacloprid, Indoxacarb, Mandipropamid, Methamidophos, Methomyl, Monocrotophos, Nicotine, Omethoate, Oxamyl, Pencycuron, Prochloraz, Propargite, Thiabendazole, Thiacloprid, Thiamethoxam, Thiodicarb, Thiophanate methyl and Triflorine

Reproduced with permission of Kent Scientific Services  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.





# Free Amino Acids in Peas (*Pisum sativum*) by HPLC-HRAM-MS


Application #AN2660

## Conditions

Column: ACE 3 AQ  
Dimensions: 150 x 3.0 mm  
Part Number: ACE-116-1503  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeCN

Time (mins)	%B
0	0
10	100

Flow Rate: 0.4 mL/min  
Injection: 5 µL  
Temperature: 30 °C  
Detection: Exactive Orbitrap high resolution MS  
ESI positive ion mode  
Capillary temperature: 350 °C



Reproduced with permission of Fera Science Ltd, York, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

# Amino Acids and Biogenic Amines in Wine and Beer

Application #AN2800


## Conditions

Column: ACE 5 C18-HL  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-321-2546  
Mobile Phase: A: 25 mM acetate buffer pH 5.8 in H<sub>2</sub>O  
B: MeCN/MeOH (80:20 v/v)

Flow Rate: 0.8 mL/min  
Injection: 20 µL  
Temperature: 16 °C  
Detection: UV, 269 nm, 280 nm and 300 nm

Time (mins)	%B
0.0	45
20.0	60
30.5	17
33.5	17
65.0	40
73.0	72
78.0	82
82.0	100
85.0	100

1. Aspartic acid
2. Glutamic acid
3. Asparagine
4. Serine
5. Hydroxyproline
6. Glutamine
7. Histidine
8. Glycine
9. Threonine
10. β-Alanine
11. Arginine
12. α-Alanine
13. GABA
14. Proline
15. Histamine
16. Tyrosine
17. Ammonium ion
18. Agmatine
19. Valine
20. Methionine
21. Cysteine
22. Isoleucine
23. Tryptophan
24. Leucine
25. Phenylalanine
26. Ornithine
27. Lysine
28. Spermidine
29. Tyramine
30. Putrescine
31. Tryptamine
32. Cadaverine
33. Phenylethylamine
34. Isoamylamine
- I.S. L-2-Aminoadipic acid



Reproduced with permission of Instituto de la Vid y el Vino de Castilla-La Mancha, Spain  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Aminoglycosides in Eggs


Application #AN1920

## Conditions

Column: ACE Excel 2 C18-PFP  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1010-1002U  
Mobile Phase: A: 20 mM HFBA in H<sub>2</sub>O/MeCN (98:2 v/v)  
B: 20 mM HFBA in MeCN/H<sub>2</sub>O (98:2 v/v)

Time (mins)	%B	Curve
0.0	5.0	-
2.0	15.0	6
4.5	19.0	6
5.5	19.5	8
6.0	22.0	6
7.0	35.0	6
9.0	48.0	8
9.5	5.0	6

Flow Rate: 0.4 mL/min  
Temperature: 40 °C  
Detection: Positive ESI MRM (transitions as shown)  
Sample: Extraction at low pH, clean up with WCX SPE cartridge  
Egg sample spiked at 100 µg/kg (CCα)



Key  
GNT Gentamicin  
TBM Tobramycin  
KNM Kanamycin  
DHS Dihydrostreptomycin  
STR Streptomycin  
HGB Higromycin-B  
SPT Spectinomycin

Reproduced with permission of Veterinary Drug Residues Department, Centro Nacional de Alimentación (CNA, AECOSAN), Madrid, Spain

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.


# Detection of Annatto

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2840

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: MeCN/0.16% acetic acid in H<sub>2</sub>O (70:30 v/v)  
Flow Rate: 1.2 mL/min  
Temperature: Ambient  
Detection: UV-Vis, 478 nm



Reproduced with permission of Lincolne Sutton & Wood, Norwich, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Anthocyanins from *Sambucus Nigra* (Elderberry)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


## Application #AN2750

### Conditions


Column: ACE UltraCore 5 SuperC18  
Dimensions: 150 x 4.6 mm  
Part Number: CORE-5A-1546U  
Mobile Phase: A: 5% formic acid in H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0	5
35	10
55	65
65	65


Flow Rate: 1 mL/min  
Temperature: 40 °C  
Detection: UV-Vis, 525 nm




1. Cyanidin-3-sambubioside-5-glucoside




2. Cyanidin-3,5-diglucoside



3. Cyanidin-3-sambubioside



4. Cyanidin-3-glucoside



Reproduced with permission of University of Vila Real, Portugal  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Appetite Suppressants by LC-MS


Application #AN1960

## Conditions

Column: ACE Excel 2 SuperC18  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1011-1002U  
Mobile Phase: A: 10 mM ammonium acetate pH 9.35 with ammonium hydroxide  
B: 10 mM ammonium acetate pH 9.35/MeCN (10:90 v/v)

Time (mins)	%B
0.0	11.11
1.0	11.11
21.0	100
23.0	100

Flow Rate: 0.5 mL/min  
Injection: 2 µL  
Temperature: 25 °C  
Detection: MS



1. Caffeine
2. Ephedrine
3. Phentermine
4. Phenolphthalein
5. Chlordiazepoxide
6. Lorcaserin
7. Fenfluramine
8. Fluoxetine
9. Diethylpropion
10. Sertraline
11. Didesmethylsibutramine
12. Rimonabant
13. N-Desmethylsibutramine
14. Sibutramine
15. Orlistat

Reproduced with permission of US FDA  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Arsenolipids from Edible Seaweed by LC-ICP-MS and LC-ESI-MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN1970

## Conditions

Column: ACE 3 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-111-1546  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeOH

Time (mins)	%B
0	0
20	100
45	100

Flow Rate: 1 mL/min  
Injection: 100 µL  
Temperature: 45 °C  
Detection: Split ratio: 75% ESI-MS: 25% ICP-MS  
Thermo Scientific Element 2 ICP-MS  
Mode: Organic mode  
Medium resolution  
Thermo Scientific Orbitrap Discovery  
Positive ESI mode  
Spray voltage: 4.5 kV  
Capillary temperature: 320 °C  
Capillary voltage: 42 V

Arsenic-containing hydrocarbon:


*m/z* 389 [M + H]<sup>+</sup> for C<sub>21</sub>H<sub>46</sub>AsO

Arsenic-containing phospholipids:

*m/z* 959 [M + H]<sup>+</sup> for C<sub>45</sub>H<sub>89</sub>AsO<sub>14</sub>P (C16:0/C16:0)

*m/z* 987 [M + H]<sup>+</sup> for C<sub>47</sub>H<sub>93</sub>AsO<sub>14</sub>P (C18:0/C16:0)

*m/z* 1015 [M + H]<sup>+</sup> for C<sub>49</sub>H<sub>97</sub>AsO<sub>14</sub>P (C20:0/C16:0)



Reproduced with permission of Trace Element Speciation Laboratory, University of Aberdeen, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Water Soluble Artificial Colours

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN3010

## Conditions

Column: ACE 3 C18  
Dimensions: 100 x 4.6 mm  
Part Number: ACE-111-1046  
Mobile Phase: A: 3 mM tetrabutylammonium bromide and 5 mM KH<sub>2</sub>PO<sub>4</sub> in H<sub>2</sub>O  
B: 5 mM tetrabutylammonium bromide in MeOH

Time (mins)	%B
0	45
20	70
30	45
40	45

Flow Rate: 0.8 mL/min  
Injection: 10 µL  
Temperature: Ambient  
Detection: UV-Vis, 420 nm, 520 nm and 600 nm



1. Amaranth
2. Sunset Yellow
3. Allura Red
4. Red 2G
5. Ponceau 4R
6. Carmoisine
7. Erythrosine

Reproduced with permission of Lincolne Sutton & Wood, Norwich, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Separation of Artificial Food Colours

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2960

## Conditions

Column: ACE 3 C18  
Dimensions: 100 x 4.6 mm  
Part Number: ACE-111-1046  
Mobile Phase: A: 3.1 mM tetrabutylammonium bromide and 5 mM KH<sub>2</sub>PO<sub>4</sub> in H<sub>2</sub>O  
B: 5 mM KH<sub>2</sub>PO<sub>4</sub> in MeOH

Time (mins)	%B
0	45
12	60
25	45

Flow Rate: 0.8 mL/min  
Injection: 10 µL  
Temperature: Ambient  
Detection: UV-Vis, 480 nm



1. Tartrazine
2. Amaranth
3. Sunset Yellow
4. Ponceau 4R

Reproduced with permission of Lincolne Sutton & Wood, Norwich, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Artificial Sweeteners Global Method


Application #AN1980

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: H<sub>2</sub>O  
B: MeCN  
C: 0.1% TFA

Time (mins)	%A	%B	%C
0	88	2	10
25	50	40	10
30	30	60	10
35	88	2	10

Flow Rate: 1 mL/min  
Injection: 50 µL  
Temperature: 30 °C  
Detection: Corona CAD



Data courtesy of Durham County Council Scientific Services, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Artificial Sweeteners (Stevia Glycosides)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN1020

## Conditions

Column: ACE Excel 2 SuperC18  
Dimensions: 150 x 2.1 mm  
Part Number: EXL-1011-1202U  
Mobile Phase: A: 10 mM sodium dihydrogen phosphate pH 2.8 in H<sub>2</sub>O  
B: 10 mM sodium dihydrogen phosphate pH 2.8 in H<sub>2</sub>O/MeCN (80:20 v/v)


Time (mins)	%B
0	39.5
4	48.0

Flow Rate: 0.6 mL/min  
Injection: 1 µL  
Temperature: 50 °C  
Detection: UV, 200 nm



Stevioside  
*m/z* = 804.87

Rebaudioside  
*m/z* = 967.01



1. Rebaudioside D
2. Rebaudioside A
3. Stevioside
4. Rebaudioside F
5. Rebaudioside C
6. Dulcoside A
7. Rubusoside
8. Rebaudioside B
9. Steviolbioside

Reproduced with permission of Shimadzu Europe GmbH. Analysed on Shimadzu Nexera X2.  
<http://images2.advanstar.com/PixelMags/lcgc-eu/pdf/2013-10-sp.pdf>

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns


# Detection of Avenacins

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2740

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: H<sub>2</sub>O/MeOH (30:70 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 225 nm  
Sample: Partially purified extract from oat root



Reproduced with permission of Sainsbury Laboratory, John Innes Centre, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Detection of the Herbicide Benfluralin

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2880

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: H<sub>2</sub>O/MeOH (15:85 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 254 nm



1. Benfluralin



# Caffeine and Metabolites

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2010



## Conditions

Column: ACE 5 SuperC18  
Dimensions: 150 x 4.6 mm  
Part Number: EXL-1211-1546U  
Mobile Phase: A: 20 mM ammonium acetate pH 7.0 in H<sub>2</sub>O  
B: 20 mM ammonium acetate pH 7.0 in MeCN/H<sub>2</sub>O (90:10 v/v)


Time (mins)	%B
0	2
45	15
48	15
49	2

Post time 10 minutes


Flow Rate: 1 mL/min  
Injection: 1 µL  
Temperature: 60 °C  
Detection: UV, 273 nm




1. Xanthine




2. Theobromine



3. Paraxanthine



4. Theophylline



5. Caffeine



# Chloramphenicol in Milk by LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2030



## Conditions

Column: ACE 3 C18  
Dimensions: 50 x 2.1 mm  
Part Number: ACE-111-0502  
Mobile Phase: A: H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0.00	10
0.05	10
2.50	95
3.00	95
3.10	10
4.50	10


Flow Rate: 0.5 mL/min  
Injection: 10 µL  
Detection: Bruker EVOQ Elite triple quad MS  
VIP heated-ESI temperature: 400 °C  
Cone gas temperature: 350 °C  
Spray voltage: -4500 V

TIC of 3 MRMs of 0.02 ppb chloramphenicol spiked in milk matrix



Chloramphenicol

MRM chromatograms of 0.05 ppb chloramphenicol in milk



Reproduced with permission of Bruker UK Ltd  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Chocolate Analysis


Application #AN2040

## Conditions

Column: ACE Excel 2 C18-Amide  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1012-1002U  
Mobile Phase: A: 10 mM ammonium formate pH 2.8 in H<sub>2</sub>O  
B: 10 mM ammonium formate pH 2.8 in MeCN/H<sub>2</sub>O (90:10 v/v)

Time (mins)	%B
0.0	5
1.5	85

Flow Rate: 1.2 mL/min  
Temperature: 42 °C  
Detection: UV, 254 nm



Reproduced with permission of Shimadzu Europa GmbH. Analysed on Shimadzu Nexera X2.  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.





# Metabolite Profiling of Coffee by LC-MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2590

## Conditions


Column: ACE Excel 1.7 C18-Amide  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1712-1002U  
Mobile Phase: A: 0.01% formic acid in H<sub>2</sub>O  
B: 0.01% formic acid in MeCN

Time (mins)	%B
0.0	3
2.5	10
8.0	100
8.5	3
10.0	3

Flow Rate: 0.5 mL/min

Detection: Exactive accurate mass MS system  
ESI in negative ion mode

Sample: Analytes between *m/z* 70-800 monitored  
Metabolites from coffee extracted into cold water by vortexing for 20 mins.  
Samples filtered prior to injection onto column and modular Accela LC system



Reproduced with permission of School of Pharmacy, University of Nottingham, UK

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)




**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Cyclodextrin-Encapsulated Flavour Compounds in Beer


Application #AN2300

## Conditions

Column: ACE 3 C18  
Dimensions: 150 x 4.0 mm  
Part Number: ACE-111-1504  
Mobile Phase: 0.1% phosphoric acid in MeOH/H<sub>2</sub>O (53:47 v/v)  
Flow Rate: 0.5 mL/min  
Injection: 20 µL  
Temperature: 35 °C  
Detection: UV, 225 nm



1. 2-Methylphenol
2. 2-Bromophenol
3. 2-Methoxy-4-vinylphenol
4. 4-Ethylphenol
5. 2,4-Dichlorophenol
6. 2,4,6-Trimethylphenol
7. 2-Isopropylphenol



Reproduced with permission of FlavorActiV Ltd, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Human Defensins in a Saliva Matrix

Application #AN1270

## Conditions


Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 50 x 3.0 mm  
Part Number: CORE-25A-1503U  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeCN

Time (mins)	%B
0	2
2	2
17	50
19	95
20	95

Flow Rate: 0.6 mL/min  
Sample Preparation: SPE on C18  
Detection: Synapt G1 QToF +ESI MS  
Sampling cone voltage: 40 V  
Source temperature: 150 °C  
Capillary voltages: 4.8 kV  
Extraction cone voltages: 41 kV  
Desolvation temperature: 500 °C  
Acquisition: 100-2000 m/z



## Defensin Human Neutrophil Peptides

1. HNP-1  
(30 amino acid residues)
2. HNP-2  
(29 amino acid residues)
3. HNP-3  
(30 amino acid residues)



Extracted ion current chromatograms  
(sum of multiply protonated ions  $[M+3H]^{3+}$ ,  $[M+4H]^{4+}$  and  $[M+5H]^{5+}$ )

## Mass spectra



Reproduced with permission of School of Science, University of Greenwich, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Ethanol Extract from Seed Cover (*Acacia Farnesiana*)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2900

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: MeOH  
B: H<sub>2</sub>O

Time (mins)	%B
0.0	85
2.5	85
60.0	50
62.5	50
70.0	85

Flow Rate: 2 mL/min  
Temperature: Ambient  
Detection: UV, 230 nm



Reproduced with permission of School of Pharmacy, University of Sunderland, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Detection of Flavone and Dibucaine

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns



Application #AN2850

## Conditions


Column: ACE 3 C18  
Dimensions: 30 x 4.6 mm  
Part Number: ACE-111-0346  
Mobile Phase: A: 6.5 mM ammonium acetate in H<sub>2</sub>O  
B: MeCN  
C: MeOH

Time (mins)	%A	%B	%C
0.0	80	10	10
5.2	0	50	50
5.6	0	0	100

Flow Rate: 2 mL/min  
Temperature: 60 °C  
Detection: UV, 200-450 nm



1. Flavone



2. Dibucaine

Reproduced with permission of Johnson & Johnson Pharmaceutical Research & Early Development, a division of Janssen-Cilag, S.A, Toledo, Spain

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)





**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Analysis of Flavonoids


Application #AN2810

## Conditions


Column: ACE 5 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-121-1546  
Mobile Phase: MeCN/0.1% formic acid in H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 1 mL/min  
Injection: 1 µL  
Temperature: Ambient  
Detection: UV, 254 nm




1. Rutin



2. Quercetin



3. Quercitrin



4. Kaempferol





# Analysis of Garlic (I)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2820

## Conditions


Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: H<sub>2</sub>O/MeOH (50:50 v/v)  
Flow Rate: 1 mL/min  
Injection: 20 µL  
Temperature: 30 °C  
Detection: UV, 210 nm



1. Allicin



2. Alliin



3. Deoxyalliin

Reproduced with permission of Neem Biotech, Cardiff, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Analysis of Garlic (II)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns



Application #AN2830

## Conditions


Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: H<sub>2</sub>O  
B: MeCN

Time (mins)	%B
0	40
20	100
25	100

Flow Rate: 1 mL/min  
Injection: 20 µL  
Temperature: 30 °C  
Detection: UV, 254 nm



1. Allicin



2. Ajoene

Reproduced with permission of Neem Biotech, Cardiff, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Metabolite Profiling of Green Tea by LC-MS

Application #AN2580

## Conditions

Column: ACE Excel 1.7 C18-Amide  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1712-1002U  
Mobile Phase: A: 0.01% formic acid in H<sub>2</sub>O  
B: 0.01% formic acid in MeCN


Time (mins)	%B
0.0	3
2.5	10
8.0	100
8.5	3
10.0	3

Flow Rate: 0.5 mL/min

Detection: Exactive accurate mass MS system  
ESI in negative ion mode

Analytes between *m/z* 70-800 monitored

Sample: Metabolites from green tea extracted into cold water by vortexing for 20 mins. Samples filtered prior to injection onto column and modular Accela LC system



Reproduced with permission of School of Pharmacy, University of Nottingham, UK

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



# Impurity Profile of a Herbicide

Application #AN2130

## Conditions


Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 150 x 4.6 mm  
Part Number: CORE-25A-1546U  
Mobile Phase: A: MeCN/H<sub>2</sub>O/TFA (5:95:0.05 v/v/v)  
B: MeCN/TFA (99.9:0.05 v/v)

Time (mins)	%B
0	10
3	10
35	100
55	100
56	10
60	10

Flow Rate: 0.6 mL/min  
Injection: 10 µL  
Temperature: 25 °C  
Detection: UV, 240 nm

Technical Grade Herbicide

Active component  
↓





# Separation of Hippuric Acid

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2760

## Conditions

Column: ACE 5 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-121-1546  
Mobile Phase: 10 mM KH<sub>2</sub>PO<sub>4</sub> pH 3.5 in H<sub>2</sub>O/MeCN (15:85 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 254 nm



1. Hippuric acid



2. 2-Methylhippuric acid




### Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: 0.1% TFA in H<sub>2</sub>O/MeCN (71:29 v/v)  
B: 0.1% TFA in H<sub>2</sub>O/MeCN (68:32 v/v)

Time (mins)	%B
0	10
16	90

Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 215 nm



1. Bovine insulin
2. Human insulin
3. Porcine insulin




# Separation of Isoflavones

Application #AN2970

## Conditions

Column: ACE 5 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-121-1546  
Mobile Phase: MeCN/0.1% formic acid in H<sub>2</sub>O (35:65 v/v)  
Flow Rate: 1 mL/min  
Injection: 1 µL  
Temperature: Ambient  
Detection: UV, 254 nm



1. Daidzin
2. Genistin
3. Daidzein
4. Genistein



# Isoflavones in Red Clover and Soy Extract

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN1130

## Conditions

Column: ACE 3 C18-AR  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-119-1502  
Mobile Phase: A: Acetic acid in H<sub>2</sub>O pH 2.8  
B: 0.6% Acetic acid in MeCN

Time (min)	%B
0	15
7	15
27	75

Flow Rate: 0.35 mL/min  
Injection: 3 µL  
Temperature: 25 °C  
Detection: UV, 254 nm



1. Daidzin
2. Glycitin
3. Rutin (Int. Standard)
4. Genistin
5. Acetyl-Daidzin
6. Acetyl-Glycitin
7. Daidzein
8. Glycitein
9. Acetyl-Genistin
10. Genistein
11. Formononetin
12. Biochanin A

K. Weinfurter et al. Forsch. Komplementmed. 21 (Suppl.1): 45 (2014)  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Fingerprinting of Liquorice Extracts

Application #AN2090


## Conditions

Column: ACE 3 C18-PFP  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-1110-1502  
Mobile Phase: A: Ammonium acetate in H<sub>2</sub>O pH 4  
B: MeOH


Time (mins)	%B
0	10
1	10
11	15
55	90
60	100

Flow Rate: 0.4 mL/min  
Injection: 2 µL  
Temperature: 40 °C  
Detection: UV, 260, 280, 330 and 350 nm  
Sample: Plant material ground to a fine powder in pestle and mortar. Powdered material extracted into methanol by ultrasonication for 30 minutes, followed by centrifugal filtration.

Methanolic liquorice leaf extract at different wavelengths



Comparison of methanolic extracts at 260 nm



Reproduced with permission of Unilever, Colworth Science Park, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.





# Detection of Malachite Green

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2860

## Conditions

Column: ACE 5 C18  
Dimensions: 150 x 3.0 mm  
Part Number: ACE-121-1503  
Mobile Phase: 10 mM oxalic acid pH 2.9 in H<sub>2</sub>O/MeCN (80:20 v/v)  
Flow Rate: 0.4 mL/min  
Temperature: Ambient  
Detection: UV-Vis, 618 nm



1. Malachite green



2. Leucomalachite green

Reproduced with permission of Central Science Laboratory, York, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)





# Detection of Melamine using Ion-Pairing Reagent

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2510

## Conditions

Column:	ACE 5 C8
Dimensions:	150 x 4.6 mm
Part Number:	ACE-122-1546
Mobile Phase:	5 mM heptafluorobutyric acid:MeCN (95:5 v/v)
Flow Rate:	1 mL/min
Injection:	5 µL
Temperature:	Ambient
Detection:	UV, 240 nm



1. Melamine



# Microcystins from Blue/Green Algae in Drinking Water


Application #AN1190

## Conditions


Column: ACE Excel 2 C18  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-101-1002U  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: MeCN

Time (mins)	%B
0.0	30
1.0	30
7.0	95
7.1	30
10.0	30

Flow Rate: 0.4 mL/min  
Injection: 50 µL  
Temperature: 40 °C  
Sample: 0.05 ppb  
Detection: Bruker EVOQ Elite triple quad MS  
VIP heated-ESI temperature: 350 °C  
Cone gas temperature: 200 °C  
Spray voltage: 4500 V (+)  
Collision gas: argon 1.5 mTorr



1. Microcystin RR (MW 1038)  
(*m/z* 520 → 135)
2. Microcystin LR (MW 995)  
(*m/z* 498 → 135)
3. Microcystin YR (MW 1045)  
*m/z* 523 → 135)



Variants	R	L
Microcystin-LR	Leucine	Arginine
Microcystin-RR	Arginine	Arginine
Microcystin-YR	Tyrosine	Arginine

Reproduced with permission of Bruker UK Ltd  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

# Separation of Milk Proteins


Application #AN1540

## Conditions

Column: ACE 5 C18-300  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-221-1502  
Mobile Phase: A: 0.01% TFA in H<sub>2</sub>O  
B: 0.01% TFA in MeCN

Time (mins)	%B
0.0	33
5.0	33
9.0	35
18.0	37
22.0	40
27.5	41
28.0	41
43.0	43

Flow Rate: 0.2 mL/min  
Temperature: 45 °C  
Detection: UV, 214 nm



Reproduced with permission of The Chemical Analysis Facility, University of Reading, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)




## Conditions

Column: ACE Excel 2 C18-AR  
 Dimensions: 50 x 2.1 mm  
 Part Number: EXL-109-0502U  
 Mobile Phase: A: 1 mM ammonium acetate, 0.5% acetic acid in H<sub>2</sub>O  
 B: 1 mM ammonium acetate, 0.5% acetic acid in 95% MeOH

Time (mins)	%B
0.0	40
1.0	40
2.4	60
6.8	87

Flow Rate: 0.6 mL/min  
 Injection: 2 µL  
 Temperature: 40 °C  
 Detection: AB SCIEX triple quad 5500  
 Positive ESI mode  
 Source temperature: 500 °C  
 IonSpray voltage: 5500 V





1. Aflatoxin G2  
(*m/z* 331.1 → 313.1)
2. Aflatoxin G1  
(*m/z* 329.0 → 243.1)
3. Aflatoxin B2  
(*m/z* 315.1 → 287.0)
4. Aflatoxin B1  
(*m/z* 313.1 → 285.0)
5. HT-2-toxin  
(*m/z* 442.2 → 263.1)
6. T-2-toxin  
(*m/z* 484.2 → 305.1)
7. Ochratoxin A  
(*m/z* 404.1 → 239.0)

Reproduced with permission of Biotage GB Ltd  
 © 2015 ACT Ltd. All rights reserved. All products are available worldwide.




## Conditions


Column: ACE 3 C18-PFP  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-1110-1546  
Mobile Phase: H<sub>2</sub>O/MeOH (60:40 v/v)  
Flow Rate: 1 mL/min  
Injection: 100 µL  
Temperature: 45 °C  
Detection: Fluorescence, Ex λ 362 nm, Em λ 425 nm




1. Aflatoxin G2



2. Aflatoxin G1



3. Aflatoxin B2




4. Aflatoxin B1



## Conditions

Column:	ACE 3 C18
Dimensions:	50 x 2.1 mm
Part Number:	ACE-111-0502
Mobile Phase:	MeOH/0.5 mM ammonium acetate in H <sub>2</sub> O (50:50 v/v)
Flow Rate:	0.2 mL/min
Injection:	20 µL
Temperature:	Ambient
Detection:	ESI MS/MS (+ve mode)
Sample:	Metabolites derivatised with 2-nitrobenzaldehyde to form nitrophenyl derivatives, prior to LC-MS analysis




# Detection of Ochratoxin A

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2870

## Conditions

Column: ACE 5 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-121-1546  
Mobile Phase: MeCN/H<sub>2</sub>O/Acetic acid (51:47:2 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: Fluorescence –  $\lambda_{\text{ex}}$  333 nm,  $\lambda_{\text{em}}$  443 nm



Reproduced with permission of R-Biopharm Rhone Ltd, Glasgow, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns



# Separation of Organic acids

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2780

## Conditions


Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 5.7 in H<sub>2</sub>O/MeOH (70:30 v/v)  
Flow Rate: 1 mL/min  
Temperature: 22 °C  
Detection: UV, 220 nm




1. L-Ascorbic acid




2. Maleic acid



3. Acetylsalicylic acid



4. Benzoic acid



5. Salicylic acid

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

  
**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns





# Separation of Organic Acids

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN1570

## Conditions


Column: ACE 3 C18-Amide  
Dimensions: 250 x 2.1 mm  
Part Number: EXL-1112-2502U  
Mobile Phase: 20 mM H<sub>3</sub>PO<sub>4</sub> in MeOH/H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 0.21 mL/min  
Injection: 5 µL  
Temperature: 20 °C  
Detection: UV, 210 nm




1. Sinapic acid




2. Ferulic acid



3. 3,4-Dimethoxycinnamic acid



4. Cinnamic acid



5. 4-Methoxycinnamic acid



# Fast Separation of Organic Acids

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2200

## Conditions


Column: ACE 3 C18-Amide  
ACE Excel 1.7 C18-Amide  
Dimensions: 250 x 2.1 mm  
50 x 3.0 mm  
Part Number: 250 x 2.1 mm (EXL-1112-2502U), 50 x 3 mm (EXL-1712-0503U)  
Mobile Phase: 20 mM H<sub>3</sub>PO<sub>4</sub> in MeOH/H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 0.21 mL/min (250 x 2.1 mm)  
0.8 mL/min (50 x 3.0 mm)  
Injection: 5 µL (250 x 2.1 mm)  
2 µL (50 x 3.0 mm)  
Temperature: 20 °C  
Detection: UV, 210 nm




1. Sinapic acid




2. Ferulic acid




3. 3,4-Dimethoxycinnamic acid



4. Cinnamic acid



5. 4-Methoxycinnamic acid



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Separation of Phenolic Acids

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN3030

## Conditions


Column: ACE 5 C18  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-121-1546  
Mobile Phase: MeCN/0.1% formic acid in H<sub>2</sub>O (20:80 v/v)  
Flow Rate: 1 mL/min  
Injection: 1 µL  
Temperature: Ambient  
Detection: UV, 254 nm



1. Caffeic acid



2. p-Coumaric acid



3. Ferulic acid



# Organophosphorus Flame Retardants in Water by LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN1240


## Conditions

Column: ACE 3 C18  
Dimensions: 100 x 2.1 mm  
Part Number: ACE-111-1002  
Mobile Phase: A: 0.05 mM ammonium formate + 0.005% formic acid in H<sub>2</sub>O  
B: MeOH/MeCN (95:5 v/v)

Time (mins)	%B
0.1	50
12.0	90
13.0	100
15.0	100
15.1	50
20.0	50

Flow Rate: 0.25 mL/min  
Injection: 80 µL  
Temperature: 25 °C  
Detection: MS/MS

	Q1 Mass	Q3 Mass
1. Trimethyl phosphate (TMP)	141	109
2. Triethyl phosphate (TEP)	183	127
3. Triisopropyl phosphate (TiPP)	225	99
4. Tri-n-propyl phosphate (TPrP)	225	99
5. Tri-n-butyl phosphate (TBP)	267	211
6. Tris(2-chloroethyl) phosphate (TCEP)	285	223
7. Tris((2R)-1-chloro-2-propyl) phosphate (TCPP)	327	99
8. Tris(1,3-dichloro-2-propyl) phosphate (TDCPP)	431	99
9. Bis(1,3-dichloro-2-propyl) phosphate (BDCP)	321	99
10. Triphenyl phosphate (TPP)	327	215
11. Tris(2-butoxyethyl) phosphate (TBEP)	399	299
12. Bis(2-ethylhexyl) phosphate (BEHP)	323	99



Reproduced with permission of Ontario Ministry of the Environment and Climate Change, Canada  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns


For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Isomeric Organophosphorus Flame Retardants in Water

## Conditions



Column:	ACE 3 C18 ACE 3 C18-PFP
Dimensions:	100 x 2.1 mm
Part Number:	ACE-111-1002, ACE-1110-1002
Mobile Phase:	A: 0.05 mM ammonium formate + 0.005% formic acid in H <sub>2</sub> O B: MeOH/MeCN (95:5 v/v)
	<b>Time (mins)</b> <b>%B</b>
	0.1                    50
	12.0                  90
	13.0                  100
	15.0                  100
	15.1                  50
	20.0                  50
Flow Rate:	0.25 mL/min
Injection:	80 µL
Temperature:	25 °C
Detection:	MS/MS

TiPP and TPrP





ACE C18

TOTP, TPTP and TMTP





1. Triiso-propyl phosphate (TiPP)  
(*m/z* 225 → 99)
2. Tri-n-propyl phosphate (TPrP)  
(*m/z* 225 → 99)
3. Tri-o-tolyl phosphate (TOTP)  
(*m/z* 369 → 91)
4. Tri-p-tolyl phosphate (TPTP)  
(*m/z* 369 → 91)
5. Tri-m-tolyl phosphate (TMTP)  
(*m/z* 369 → 91)

ACE C18-PFP




## Conditions


Column: ACE 3 Phenyl  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-115-1502  
Mobile Phase: 25 mM ammonium acetate pH 6.8 in H<sub>2</sub>O/MeOH (50:50 v/v)  
Flow Rate: 0.2 mL/min  
Injection: 2 µL  
Temperature: 40 °C  
Detection: UV, 240 nm




1. Methyl paraben




2. Ethyl paraben




3. n-Propyl paraben



4. i-Butyl paraben



5. n-Butyl paraben



6. Benzyl paraben



# Perfluoro Acids by LC-MS/MS


Application #AN1280

## Conditions

Column: ACE Excel 2 C18  
Dimensions: 50 x 2.1 mm  
Part Number: EXL-101-0502U  
Mobile Phase: A: 2 mM ammonium acetate, 0.1% acetic acid/MeCN (95:5 v/v)  
B: 2 mM ammonium acetate, 0.1% acetic acid/MeCN (5:95 v/v)

Time (mins)	%B
0.0	25
0.5	25
5.5	95
7.5	95
8.0	25
10.0	25

Flow Rate: 0.5 mL/min  
Injection: 20  $\mu$ L  
Temperature: 40  $^{\circ}$ C  
Detection: AB SCIEX triple quad 5500  
Negative ESI MRM  
Source temperature: 450  $^{\circ}$ C  
IonSpray voltage: -2400 V



1. Heptafluorobutyric acid ( $m/z$  212.9  $\rightarrow$  168.9)
2. Perfluorohexanoic acid ( $m/z$  313  $\rightarrow$  268.9)
3. Perfluorobutylsulfonic acid ( $m/z$  299  $\rightarrow$  79.9)
4. Perfluoroheptanoic acid ( $m/z$  363  $\rightarrow$  319)
5. Perfluorooctanoic acid ( $m/z$  413  $\rightarrow$  368.9)
6. Perfluorohexylsulfonic acid ( $m/z$  399  $\rightarrow$  80)
7. Perfluorononanoic acid ( $m/z$  463  $\rightarrow$  419)
8. Perfluorodecanoic acid ( $m/z$  513  $\rightarrow$  469)
9. Perfluorooctanesulfonic acid ( $m/z$  499  $\rightarrow$  80)

Reproduced with permission of Biotage GB Ltd, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Perfluoroalkyl Substances by Ion-Pairing LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2560


## Conditions

Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 50 x 2.1 mm  
Part Number: CORE-25A-0502U  
Mobile Phase: A: 2 mM ammonium acetate + 5 mM 1-methylpiperidine in MeOH/H<sub>2</sub>O (5:95 v/v)  
B: 2 mM ammonium acetate + 5 mM 1-methylpiperidine in MeOH/H<sub>2</sub>O (95:5 v/v)

Time (mins)	%B
0.0	10
0.3	10
1.0	20
1.5	50
5.0	80
10.0	80
13.0	100
16.0	100

Flow Rate: 0.3 mL/min  
Injection: 5 µL  
Temperature: 35 °C  
Detection: Agilent 6430 triple quad MS  
ESI in negative ion mode  
Capillary voltage: 3000 V  
Nebulizer pressure: 50 psi

1. PFBA  
(*m/z* 213 → 169)
2. PFPeA  
(*m/z* 263 → 219)  
(*m/z* 263 → 175)
3. PFBS  
(*m/z* 299 → 99)  
(*m/z* 299 → 80)
4. PFHxA  
(*m/z* 313 → 269)  
(*m/z* 313 → 119)
5. PFHpA  
(*m/z* 363 → 319)  
(*m/z* 363 → 169)
6. PFHxS  
(*m/z* 399 → 99)  
(*m/z* 399 → 80)
7. PFOA  
(*m/z* 413 → 369)  
(*m/z* 413 → 169)
8. PFNA  
(*m/z* 463 → 419)  
(*m/z* 463 → 169)
9. PFOS  
(*m/z* 499 → 99)  
(*m/z* 499 → 80)
10. FOSA  
(*m/z* 498 → 498)  
(*m/z* 498 → 78)
11. PFDA  
(*m/z* 513 → 469)  
(*m/z* 513 → 269)



Reproduced with permission of Department of Analytical Chemistry, Faculty of Science and Technology, University of the Basque Country, Spain

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Perfluorinated Compounds in Water by LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2260

## Conditions

Column: ACE Excel 1.7 C18  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-171-1002U  
Mobile Phase: A: 2 mM ammonium acetate, 0.1% formic acid in H<sub>2</sub>O/MeCN (90:10 v/v)  
B: 2 mM ammonium acetate, 0.1% formic acid in H<sub>2</sub>O/MeCN (10:90 v/v)

Time (mins)	%B
0.0	25
0.5	25
3.5	70
4.0	100
5.5	100
6.0	25
9.0	25

Flow Rate: 0.5 mL/min  
Injection: 10 µL  
Temperature: 40 °C  
Detection: AB SCIEX triple quad 5500  
Negative ESI MRM  
Source temperature: 450 °C  
IonSpray voltage: -2400 V



1. Perfluorohexanoic acid  
(*m/z* 313.0 → 268.9)
2. Perfluorobutanesulfonic acid  
(*m/z* 299.0 → 79.9)
3. Perfluorooctanoic acid  
(*m/z* 413.0 → 368.9)
4. Perfluorohexanesulfonic acid  
(*m/z* 399.0 → 80.0)
5. Perfluorononanoic acid  
(*m/z* 463.0 → 419.0)
6. Perfluorodecanoic acid  
(*m/z* 513.0 → 469.0)
7. Perfluorooctanesulfonic acid  
(*m/z* 499.0 → 80.0)

Reproduced with permission of Biotage GB Ltd, UK. For extraction conditions see Biotage Application Note AN839

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# 250 Pesticide Screen using LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN3060


Page 1 of 4

## Conditions

Column: ACE Excel 2 C18  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-101-1002U  
Mobile Phase: A: 10 mM ammonium formate + 0.05% formic acid in H<sub>2</sub>O  
B: 10 mM ammonium formate + 0.05% formic acid in MeOH

Time (mins)	%B
0.00	2
0.25	30
10.00	100
12.00	100
12.50	2
14.50	2

Flow Rate: 0.5 mL/min  
Temperature: 50 °C  
Detection: TSQ Quantiva triple quad MS  
Positive mode HESI  
Spray voltage: 3500 V  
Ion transfer tube temperature: 350 °C  
Vaporizer temperature: 300 °C



Reproduced with permission of Centre for Analytical Chemistry, California Department of Food and Agriculture, USA

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# 250 Pesticide Screen using LC-MS/MS

Application #AN3060

Page 2 of 4

Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z	Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z
3-OH Carbofuran	2.25	[M+H] <sup>+</sup>	238.1	181.2	163.1	Cyrosulfamide	3.30	[M+H] <sup>+</sup>	375.1	135.1	254.1
5-OH Thiabendazole	1.66	[M+H] <sup>+</sup>	218.0	147.2	191.1	Cyromazine	1.15	[M+H] <sup>+</sup>	167.1	125.2	68.2
Abamectin	9.45	[M+NH <sub>4</sub> ] <sup>+</sup>	890.5	305.3	567.5	DEF	9.20	[M+H] <sup>+</sup>	315.1	169.0	113.0
Acephate	1.26	[M+H] <sup>+</sup>	184.0	143.1	125.1	Demeton-S sulfone	2.55	[M+H] <sup>+</sup>	291.1	235.1	263.1
Acetamiprid	2.24	[M+H] <sup>+</sup>	223.1	126.1	90.1	Dialifos	7.46	[M+H] <sup>+</sup>	394.0	208.1	181.0
Aldicarb	2.95	[M+NH <sub>4</sub> ] <sup>+</sup>	208.1	116.1	89.0	Diazinon	7.12	[M+H] <sup>+</sup>	305.1	169.1	153.2
Aldicarb sulfone	1.44	[M+NH <sub>4</sub> ] <sup>+</sup>	240.1	148.0	86.0	Diazinon OA	5.32	[M+H] <sup>+</sup>	289.1	153.2	233.1
Aldicarb sulfoxide	1.37	[M+NH <sub>4</sub> ] <sup>+</sup>	224.1	132.0	89.1	Dichlormid	3.85	[M+H] <sup>+</sup>	208.0	140.0	81.2
Allethrin	8.33	[M+H] <sup>+</sup>	303.2	135.1	123.1	Dichlorvos	3.63	[M+H] <sup>+</sup>	221.0	109.1	127.0
Ametoctradin	7.64	[M+H] <sup>+</sup>	276.2	149.1	176.2	Dicrotophos	1.87	[M+H] <sup>+</sup>	238.1	112.2	193.1
Atrazine	4.64	[M+H] <sup>+</sup>	216.1	174.0	104.0	Diethofencarb	5.53	[M+H] <sup>+</sup>	268.2	124.1	180.2
Azinphos ethyl	6.30	[M+H] <sup>+</sup>	346.0	132.1	233.0	Diflubenzuron	6.66	[M+H] <sup>+</sup>	311.0	158.0	141.0
Azinphos methyl	5.14	[M+H] <sup>+</sup>	318.0	132.0	124.9	Dimethenamid	5.70	[M+H] <sup>+</sup>	276.1	244.1	168.2
Azinphos methyl OA	2.98	[M+H] <sup>+</sup>	302.0	132.2	160.0	Dimethoate	2.23	[M+H] <sup>+</sup>	230.1	199.0	125.0
Azoxystrobin	5.59	[M+H] <sup>+</sup>	404.1	372.1	344.1	Dimethomorph	5.76, 6.07	[M+H] <sup>+</sup>	388.1	301.0	165.1
Bendiocarb	3.72	[M+H] <sup>+</sup>	224.1	167.1	109.1	Dinotefuran	1.36	[M+H] <sup>+</sup>	203.1	129.1	114.2
Benoxacor	5.23	[M+H] <sup>+</sup>	260.1	134.1	120.1	Dioxacarb	2.26	[M+H] <sup>+</sup>	224.1	123.1	167.1
Bifenazate	6.27	[M+H] <sup>+</sup>	301.1	198.0	170.1	Dioxathion	8.10	[M-C <sub>4</sub> H <sub>10</sub> O <sub>2</sub> PS <sub>2</sub> ] <sup>+</sup>	271.1	97.0	125.0
Bitertanol	7.41	[M+H] <sup>+</sup>	338.2	269.3	99.1	Disulfoton sulfone	4.59	[M+H] <sup>+</sup>	307.0	261.1	125.0
Boscalid	5.85	[M+H] <sup>+</sup>	343.0	307.0	140.0	Disulfoton sulfoxide	4.49	[M+H] <sup>+</sup>	291.0	185.1	213.1
Bupirimate	6.68	[M+H] <sup>+</sup>	317.2	210.2	237.3	Diuron	4.82	[M+H] <sup>+</sup>	233.0	72.1	160.0
Buprofezin	8.24	[M+H] <sup>+</sup>	306.1	201.1	106.1	DMST	3.90	[M+H] <sup>+</sup>	215.1	106.1	151.0
Cadusafos	7.58	[M+H] <sup>+</sup>	271.1	159.0	131.0	Dodine	7.56	[M+H] <sup>+</sup>	228.3	186.3	60.1
Carbaryl	4.07	[M+NH <sub>4</sub> ] <sup>+</sup>	219.1	145.1	127.0	Emamectin	8.57	[M+H] <sup>+</sup>	886.5	158.1	126.1
Carbendazim	2.10	[M+H] <sup>+</sup>	192.1	160.1	132.1	Ethiofencarb	4.27	[M+H] <sup>+</sup>	226.1	107.1	169.1
Carbofuran	3.77	[M+H] <sup>+</sup>	222.1	165.2	123.2	Ethiofencarb sulfone	1.90	[M+NH <sub>4</sub> ] <sup>+</sup>	275.1	107.1	201.1
Carboxin	3.97	[M+H] <sup>+</sup>	236.1	143.0	93.0	Ethiofencarb sulfoxide	1.98	[M+H] <sup>+</sup>	242.1	107.1	185.0
Carfentrazone ethyl	6.88	[M+H] <sup>+</sup>	412.0	346.1	366.0	Ethion	8.31	[M+H] <sup>+</sup>	385.0	199.1	143.0
Chlorantraniliprole	5.24	[M+H] <sup>+</sup>	484.0	286.0	194.0	Ethion monoxon	6.73	[M+H] <sup>+</sup>	369.0	199.0	143.0
Chlorfenvinphos	7.21	[M+H] <sup>+</sup>	359.0	170.0	99.1	Ethiprole	5.77	[M+NH <sub>4</sub> ] <sup>+</sup>	413.9	351.0	255.0
Chlorimuron ethyl	5.73	[M+H] <sup>+</sup>	415.1	186.0	83.0	Ethofumesate	5.55	[M+H] <sup>+</sup>	287.1	121.1	241.1
Chlorpyrifos	8.47	[M+H] <sup>+</sup>	349.9	198.0	97.0	Ethoprop	6.46	[M+H] <sup>+</sup>	243.1	173.0	131.0
Chlorpyrifos OA	6.65	[M+H] <sup>+</sup>	334.0	278.0	197.9	Etofenprox	9.75	[M+NH <sub>4</sub> ] <sup>+</sup>	394.2	177.2	107.1
Clethodim	7.71	[M+H] <sup>+</sup>	360.3	164.1	136.1	Etoxazole	8.73	[M+H] <sup>+</sup>	360.2	141.0	304.2
Clofentezine	7.38	[M+H] <sup>+</sup>	303.0	138.1	102.0	Famoxadone	7.24	[M+NH <sub>4</sub> ] <sup>+</sup>	392.2	331.1	238.0
Cloransulam methyl	4.13	[M+H] <sup>+</sup>	430.0	398.1	370.0	Fenamidone	5.76	[M+H] <sup>+</sup>	312.1	236.1	92.2
Clothianidin	1.99	[M+H] <sup>+</sup>	250.0	169.1	132.0	Fenamiphos	6.71	[M+H] <sup>+</sup>	304.1	217.1	202.0
Coumaphos	7.07	[M+H] <sup>+</sup>	363.0	227.1	307.1	Fenamiphos sulfone	4.10	[M+H] <sup>+</sup>	336.1	266.1	188.1
Crotoxyphos	5.86	[M+NH <sub>4</sub> ] <sup>+</sup>	332.1	127.1	193.1	Fenamiphos sulfoxide	3.96	[M+H] <sup>+</sup>	320.1	233.1	171.1
Crufomate	6.77	[M+H] <sup>+</sup>	292.1	236.1	108.1	Fenazaquin	9.21	[M+H] <sup>+</sup>	307.2	161.2	57.2
Cyantraniliprole	4.33	[M+2+H] <sup>+</sup>	475.0	286.0	444.1	Fenhexamid	6.39	[M+H] <sup>+</sup>	302.1	178.0	97.2
Cyazofamid	6.52	[M+H] <sup>+</sup>	325.1	108.1	261.2	Fenobucarb	5.49	[M+H] <sup>+</sup>	208.1	95.0	152.0
Cyflufenamid	7.42	[M+H] <sup>+</sup>	413.1	295.1	203.0	Fenoxaprop ethyl	8.04	[M+H] <sup>+</sup>	362.1	288.1	91.1
Cymoxanil	2.48	[M+H] <sup>+</sup>	199.1	128.1	111.1	Fenoxycarb	6.80	[M+H] <sup>+</sup>	302.1	88.1	116.1
Cyphenothrin	9.27	[M+NH <sub>4</sub> ] <sup>+</sup>	393.2	151.2	123.2	Fenpropimorph	6.42	[M+H] <sup>+</sup>	304.3	147.2	119.1



Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z	Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z
Fenprothion	8.90	[M+H] <sup>+</sup>	422.2	366.1	214.2	Mepanipyrim	6.21	[M+H] <sup>+</sup>	224.1	106.2	77.1
Fensulfothion	4.89	[M+H] <sup>+</sup>	309.0	235.0	281.1	Mesotrione	2.01	[M+H] <sup>+</sup>	340.1	228.1	104.1
Fenuron	2.17	[M+H] <sup>+</sup>	165.1	72.1	77.1	Metaflumizone	8.30	[M+H] <sup>+</sup>	507.1	178.0	287.1
Flonicamid	1.66	[M+H] <sup>+</sup>	230.1	203.0	98.0	Metalaxyl	4.91	[M+H] <sup>+</sup>	280.1	220.1	192.1
Fluazifop P butyl	8.12	[M+H] <sup>+</sup>	384.1	282.2	328.2	Metaldehyde	2.02	[M+NH <sub>4</sub> ] <sup>+</sup>	194.1	62.2	45.3
Fludioxonil	5.76	[M+NH <sub>4</sub> ] <sup>+</sup>	266.1	158.1	131.0	Metconazole	7.32	[M+H] <sup>+</sup>	320.2	70.1	125.0
Flufenoxuron	8.79	[M+H] <sup>+</sup>	489.0	158.1	141.1	Methamidophos	1.16	[M+H] <sup>+</sup>	142.0	94.2	125.1
Flufenpyr ethyl	6.72	[M+H] <sup>+</sup>	409.1	335.0	307.0	Methidathion	4.97	[M+NH <sub>4</sub> ] <sup>+</sup>	320.0	145.1	85.1
Flumetsulam	2.03	[M+H] <sup>+</sup>	326.1	129.1	109.0	Methiocarb	5.64	[M+H] <sup>+</sup>	226.1	169.2	121.1
Flumiclorac pentyl	8.13	[M+NH <sub>4</sub> ] <sup>+</sup>	441.1	308.1	354.1	Methiocarb sulfone	2.35	[M+NH <sub>4</sub> ] <sup>+</sup>	275.0	122.1	201.1
Fluometuron	4.31	[M+H] <sup>+</sup>	233.1	72.2	46.3	Methiocarb sulfoxide	2.10	[M+H] <sup>+</sup>	242.1	185.1	122.1
Fluopicolide	6.00	[M+H] <sup>+</sup>	383.0	173.0	145.0	Methomyl	1.61	[M+H] <sup>+</sup>	163.1	106.1	88.1
Fluopyram	6.33	[M+H] <sup>+</sup>	397.1	173.0	208.0	Methoxyfenozide	6.04	[M+H] <sup>+</sup>	369.2	149.1	313.1
Fluoxastrobin	6.40	[M+H] <sup>+</sup>	459.1	427.2	188.1	Metolcarb	3.28	[M+H] <sup>+</sup>	166.1	109.1	94.1
Fluridone	5.32	[M+H] <sup>+</sup>	330.1	309.1	290.0	Metribuzin	3.59	[M+H] <sup>+</sup>	215.1	187.1	131.1
Flusilazole	6.77	[M+H] <sup>+</sup>	316.1	247.2	165.1	Mevinphos	2.70	[M+NH <sub>4</sub> ] <sup>+</sup>	242.1	193.1	127.1
Fluthiacet methyl	6.88	[M+H] <sup>+</sup>	404.0	344.0	273.9	Monocrotophos	1.71	[M+H] <sup>+</sup>	224.1	193.0	127.0
Flutolanil	5.95	[M+H] <sup>+</sup>	324.1	262.0	282.0	Monolinuron	4.16	[M+H] <sup>+</sup>	215.1	126.1	148.1
Flutriafol	4.74	[M+H] <sup>+</sup>	302.1	70.1	123.1	Myclobutanil	6.15	[M+H] <sup>+</sup>	289.1	125.0	70.1
Fluxapyroxad	6.02	[M+H] <sup>+</sup>	382.1	342.1	314.1	Nicosulfuron	3.45	[M+H] <sup>+</sup>	411.1	182.0	213.0
Forchlorfenuron	4.78	[M+H] <sup>+</sup>	248.1	129.1	93.1	Norflurazon	4.98	[M+H] <sup>+</sup>	304.0	160.0	140.0
Formetanate HCl	1.26	[M+H] <sup>+</sup>	222.0	165.1	120.0	Norflurazon desmethyl	4.43	[M+H] <sup>+</sup>	290.0	179.0	140.0
Fosthiazate	4.40	[M+H] <sup>+</sup>	284.1	104.1	228.1	Omethoate	1.33	[M+H] <sup>+</sup>	214.0	183.0	125.0
Hexaconazole	7.29	[M+H] <sup>+</sup>	314.1	158.9	70.0	Oxamyl	1.48	[M+NH <sub>4</sub> ] <sup>+</sup>	237.1	72.0	90.0
Hexythiazox	8.51	[M+H] <sup>+</sup>	353.1	228.0	168.0	Oxamyl oxime	1.34	[M+H] <sup>+</sup>	163.1	72.1	90.1
Imazalil	5.14	[M+H] <sup>+</sup>	297.1	159.1	255.1	Oxydemeton methyl	1.57	[M+H] <sup>+</sup>	247.0	169.1	109.1
Imazosulfuron	5.28	[M+H] <sup>+</sup>	413.0	153.0	156.1	Oxydemeton methyl sulfone	1.62	[M+H] <sup>+</sup>	263.0	169.0	109.0
Imidacloprid	1.96	[M+H] <sup>+</sup>	256.1	209.1	175.0	Parathion methyl OA	3.10	[M+H] <sup>+</sup>	248.0	202.0	109.1
Imiprothrin	6.34	[M+H] <sup>+</sup>	319.2	151.1	123.1	Parathion OA	4.61	[M+H] <sup>+</sup>	276.1	220.1	248.1
Indaziflam	6.58	[M+H] <sup>+</sup>	302.2	158.1	145.1	Pencycuron	7.50	[M+H] <sup>+</sup>	329.1	125.1	89.1
Indoxacarb	7.75	[M+H] <sup>+</sup>	528.1	249.0	150.1	Penflufen	6.95	[M+H] <sup>+</sup>	318.2	234.1	141.0
Ipconazole	7.81	[M+H] <sup>+</sup>	334.2	70.1	125.0	Penthiopyrad	7.05	[M+H] <sup>+</sup>	360.1	177.1	276.1
Iprovalicarb	6.31	[M+H] <sup>+</sup>	321.2	119.1	186.2	Phenothrin	9.56	[M+H] <sup>+</sup>	351.2	183.1	168.0
Isofenphos	7.39	[M+H] <sup>+</sup>	346.1	217.0	245.1	Phenthoate	6.81	[M+H] <sup>+</sup>	321.0	247.1	79.1
Isoprocarb	4.67	[M+H] <sup>+</sup>	194.1	95.1	152.2	Phorate OA	5.10	[M+H] <sup>+</sup>	245.0	75.2	47.2
Isoproturon	4.79	[M+H] <sup>+</sup>	207.2	72.2	165.2	Phorate OA Sulfone	2.51	[M+H] <sup>+</sup>	277.0	155.0	127.0
Kresoxim methyl	6.90	[M+H] <sup>+</sup>	314.1	267.2	222.1	Phorate OA Sulfoxide	2.31	[M+H] <sup>+</sup>	261.0	153.0	81.0
Lactofen	8.22	[M+NH <sub>4</sub> ] <sup>+</sup>	479.1	344.1	223.0	Phorate Sulfone	4.61	[M+H] <sup>+</sup>	293.0	114.9	171.0
Lenacil	4.67	[M+H] <sup>+</sup>	235.1	153.1	136.1	Phorate Sulfoxide	4.49	[M+H] <sup>+</sup>	277.0	170.9	199.0
Leptophos OA	7.75	[M+2+H] <sup>+</sup>	396.9	155.1	364.9	Phosalone	7.35	[M+H] <sup>+</sup>	368.0	182.0	111.1
Linuron	5.46	[M+H] <sup>+</sup>	249.0	182.1	160.1	Phosmet	5.21	[M+H] <sup>+</sup>	318.0	160.1	133.1
Malathion	5.92	[M+H] <sup>+</sup>	331.0	127.1	285.1	Phosmet OA	3.12	[M+H] <sup>+</sup>	302.0	160.0	133.0
Malathion OA	3.89	[M+H] <sup>+</sup>	315.1	127.1	99.0	Phosphamidon	3.43	[M+H] <sup>+</sup>	300.1	127.1	174.1
Mandipropamid	5.94	[M+H] <sup>+</sup>	412.1	328.2	356.2	Phoxim	7.25	[M+H] <sup>+</sup>	299.1	77.2	129.1
Mefenpyr diethyl	7.26	[M+H] <sup>+</sup>	373.1	327.1	160.0	Picoxystrobin	6.79	[M+H] <sup>+</sup>	368.1	145.0	115.0



# 250 Pesticide Screen using LC-MS/MS

Application #AN3060

Page 4 of 4

Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z	Analyte	R <sub>t</sub>	Adduct	Precursor Ion m/z	Quant Ion m/z	Conf Ion m/z
Pirimicarb	4.24	[M+H] <sup>+</sup>	239.2	182.1	72.0	Spiromesifen	8.66	[M+NH <sub>4</sub> ] <sup>+</sup>	388.1	273.1	187.0
Pirimicarb Desmethyl	2.71	[M+H] <sup>+</sup>	225.1	168.2	72.1	Spiromesifen Alcohol	5.01	[M+H] <sup>+</sup>	273.2	187.1	179.1
Pirimiphos Methyl	7.34	[M+H] <sup>+</sup>	306.1	164.2	108.1	Spirotetramat	6.38	[M+H] <sup>+</sup>	374.2	302.3	216.2
Prallethrin	7.69	[M+H] <sup>+</sup>	301.2	133.0	151.2	Spiroxamine	5.95	[M+H] <sup>+</sup>	298.3	144.2	100.2
Prochloraz	7.39	[M+H] <sup>+</sup>	376.0	308.1	70.1	Sulfoxaflor	2.39	[M+NH <sub>4</sub> ] <sup>+</sup>	295.2	174.1	154.1
Profoxydim	7.71, 9.00	[M+H] <sup>+</sup>	466.2	280.0	180.0	Sulprofos	8.56	[M+H] <sup>+</sup>	323.0	219.1	139.1
Promecarb	5.88	[M+H] <sup>+</sup>	208.1	109.0	151.1	TCMTB	5.48	[M+H] <sup>+</sup>	239.0	180.0	136.0
Propamocarb	1.41	[M+H] <sup>+</sup>	189.1	102.0	144.0	Tebufenozide	6.78	[M+H] <sup>+</sup>	353.2	133.0	104.8
Propaquizafop	8.21	[M+H] <sup>+</sup>	444.1	299.2	371.2	Tebufenpyrad	8.19	[M+H] <sup>+</sup>	334.2	117.1	145.1
Propargite	8.74	[M+NH <sub>4</sub> ] <sup>+</sup>	368.2	231.2	175.1	Tebuthiuron	3.89	[M+H] <sup>+</sup>	229.1	172.0	116.0
Propetamphos	6.13	[M+H] <sup>+</sup>	282.1	138.1	156.1	Tepraloxydim	4.10, 6.19	[M+H] <sup>+</sup>	342.2	250.1	166.1
Propoxur (S)	3.69	[M+H] <sup>+</sup>	210.1	168.2	111.1	Terbufos Sulfone	5.46	[M+H] <sup>+</sup>	321.0	115.0	143.0
Prosulfuron	5.29	[M+H] <sup>+</sup>	420.1	167.1	141.1	Terbufos Sulfoxide	5.49	[M+H] <sup>+</sup>	305.1	97.0	187.0
Pymetrozine	1.44	[M+H] <sup>+</sup>	218.1	105.1	78.1	Terbutylazine	5.71	[M+H] <sup>+</sup>	230.1	174.1	104.1
Pyraclostrobin	7.30	[M+H] <sup>+</sup>	388.1	163.1	194.1	Tetrachlorvinphos	6.86	[M+2+H] <sup>+</sup>	366.9	127.1	206.0
Pyraflufen Ethyl	7.13	[M+H] <sup>+</sup>	413.0	339.0	253.1	Tetramethrin	7.91, 8.10	[M+H] <sup>+</sup>	332.2	164.1	135.1
Pyrazophos	7.31	[M+H] <sup>+</sup>	374.1	222.2	194.1	Thiabendazole	2.48	[M+H] <sup>+</sup>	202.0	175.0	131.1
Pyridaben	9.22	[M+H] <sup>+</sup>	365.1	309.0	147.1	Thiacloprid	2.55	[M+H] <sup>+</sup>	253.0	126.1	99.1
Pyridalyl	10.21	[M+2+H] <sup>+</sup>	492.0	110.9	164.0	Thiamethoxam	1.65	[M+H] <sup>+</sup>	292.0	211.1	181.1
Pyrimethanil	5.45	[M+H] <sup>+</sup>	200.1	107.1	168.1	Thifensulfuron Methyl	3.28	[M+H] <sup>+</sup>	388.0	167.1	205.0
Pyriproxyfen	8.39	[M+H] <sup>+</sup>	322.1	96.0	227.1	Thiobencarb	7.46	[M+H] <sup>+</sup>	258.1	125.0	89.0
Quinalphos	6.78	[M+H] <sup>+</sup>	299.1	163.1	147.1	Thiodicarb	4.34	[M+H] <sup>+</sup>	355.1	163.2	88.1
Quinoxifen	8.50	[M+H] <sup>+</sup>	308.0	197.1	214.1	Thionazin	4.74	[M+H] <sup>+</sup>	249.1	193.1	97.0
Quizalofop Ethyl	8.01	[M+H] <sup>+</sup>	373.1	299.2	255.1	Topramezone	1.63	[M+H] <sup>+</sup>	364.1	334.1	125.1
Resmethrin	9.40	[M+H] <sup>+</sup>	339.2	128.1	171.1	Triadimefon	6.07	[M+H] <sup>+</sup>	294.1	197.0	225.0
Rimsulfuron	3.94	[M+H] <sup>+</sup>	432.1	182.1	139.0	Triadimenol	6.25	[M+H] <sup>+</sup>	296.1	70.2	99.0
Rotenone	6.71	[M+H] <sup>+</sup>	395.2	213.2	192.1	Triazophos	6.19	[M+H] <sup>+</sup>	314.1	162.1	119.1
Saflufenacil	5.32	[M+H] <sup>+</sup>	501.1	349.1	198.0	Tribenuron Methyl	4.59	[M+H] <sup>+</sup>	396.1	155.1	181.1
Sedaxane	6.20, 6.54	[M+H] <sup>+</sup>	332.2	159.0	139.0	Trichlorfon	2.26	[M+H] <sup>+</sup>	256.9	109.0	221.0
Sethoxydim	8.03	[M+H] <sup>+</sup>	328.2	178.0	220.1	Tricyclazole	2.80	[M+H] <sup>+</sup>	190.0	163.1	136.1
Simazine	3.66	[M+H] <sup>+</sup>	202.1	104.1	132.1	Trifloxystrobin	7.78	[M+H] <sup>+</sup>	409.1	186.2	206.2
Spinetoram	8.14	[M+H] <sup>+</sup>	748.5	142.1	203.1	Triflumizole	7.87	[M+H] <sup>+</sup>	346.1	278.0	73.0
Spinosad A	7.69	[M+H] <sup>+</sup>	732.5	142.1	98.0	Triforine	5.23	[M+2+H] <sup>+</sup>	434.9	213.0	98.2
Spinosad D	8.10	[M+H] <sup>+</sup>	746.5	142.1	98.0	Zoxamide	7.09	[M+H] <sup>+</sup>	336.0	187.0	159.0
Spirodiclofen	8.91	[M+H] <sup>+</sup>	411.1	313.1	71.1						




## Conditions


Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 50 x 2.1 mm  
Part Number: CORE-25A-0502U  
Mobile Phase: A: 0.1% formic acid + 5 mM ammonium formate in H<sub>2</sub>O/MeOH (90:10 v/v)  
B: 0.1% formic acid + 5 mM ammonium formate in H<sub>2</sub>O/MeOH (10:90 v/v)

Time (mins)	%B
0.00	0
1.00	0
15.00	100
18.00	100
18.05	0
20.00	0


Flow Rate: 0.4 mL/min  
Injection: 20 µL  
Temperature: 40 °C  
Detection: Agilent 6420 Triple Quadrupole MS, +ve mode ESI, Dynamic MRM




1. Carbaryl  
(*m/z* 202.10 → 145.10)




2. Dimethomorphs  
(*m/z* 388.10 → 301.10)



3. Cyproconazoles  
(*m/z* 292.10 → 70.00)



4. Fenbuconazole  
(*m/z* 337.10 → 70.00)



Also analysed under same conditions:

Acephate, Acetamiprid, Aldicarb, Aldicarb sulfone, Aldicarb sulfoxide, Benomyl, Carbendazim, Carbofuran, Clofentezine, Clothianidin, Cyfluthrin, Demeton S-methylsulfone, Demeton S-methylsulfoxide, Dicrotophos, Dimethoate, Dinotefuran, DMA, DMPF, Flubendiamide, Folpet, Formetanate, Hexaconazole, Hexaflumuron, Imidacloprid, Indoxacarb, Mandipropamid, Methamidophos, Methomyl, Monocrotophos, Nicotine, Omethoate, Oxamyl, Pencycuron, Prochloraz, Propargite, Thiabendazole, Thiacloprid, Thiamethoxam, Thiodicarb, Thiophanate methyl and Triflorine

Reproduced with permission of Kent Scientific Services  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.




## Conditions

Column: ACE 3 C18  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-111-1502  
Mobile Phase: A: 0.1 M ammonium acetate in H<sub>2</sub>O  
B: MeCN

Time (mins)	%B
0	10
40	80
47	90
49	10

Flow Rate: 0.3 mL/min  
Injection: 25 µL  
Temperature: 40 °C  
Detection: UV, 220 nm (Pendimethalin at 245 nm)  
Sample: 0.05 µg/L standards in MeCN/H<sub>2</sub>O 10:90 v/v



1. Deisopropylatrazine
2. Desethylatrazine
3. Simazine
4. Cyanazine
5. Atrazine
6. Internal standard
7. Sebuthylazine
8. Propazine
9. Terbutylazine
10. Prometryn
11. Terbutryn
12. Alachlor
13. Pendimethalin

Reproduced with permission of Amt der Tiroler Landesregierung, Chemisch-technische Umweltschutzanstalt, Innsbruck, Austria

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.




# Detection of 13 Phenolic Compounds in Ground Water & Landfill Leachates

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN3070

## Conditions

Column: ACE Excel 3 C18-Amide  
Dimensions: 150 x 4.6 mm  
Part Number: EXL-1112-1546U  
Mobile Phase: 0.1% formic acid v/v in H<sub>2</sub>O/MeCN (65:35 v/v)  
Flow Rate: 1 mL/min  
Injection: 10 µL  
Temperature: 30 °C  
Detection: UV, 274 nm



Reproduced with permission of Environmental Scientific Group Ltd, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Separation of Phenol and Phenoxy Acid Herbicides


Application #AN2290

## Conditions

Column: ACE 3 C18-PFP  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-1110-1546  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0.0	10
20.0	100

Flow Rate: 1 mL/min  
Injection: 10 µL  
Temperature: 35 °C  
Detection: UV, 280 nm



1. Phenol
2. o-Cresol
3. 2-Chlorophenol
4. 4-Chlorophenol
5. 2,6-Dichlorophenol
6. 6-CP
7. 2,4-D
8. MCPA
9. PCOC
10. 2,4-DCP
11. 2,4-DP
12. CMPP
13. 2,4-DB
14. MCPB

Reproduced with permission of Nufarm UK Ltd  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



# Phenols in Purple Coneflower (*Echinacea Purpurea*)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2920

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: 0.1% H<sub>3</sub>PO<sub>4</sub> in H<sub>2</sub>O  
B: MeCN

Time (mins)	%B
0	10
13	22
14	40

Flow Rate: 1.5 mL/min  
Injection: 10 µL  
Temperature: 35 °C  
Detection: UV, 330 nm



1. Caftaric acid
2. Chlorogenic acid
3. Cichoric acid

Reproduced with permission of Bionorica Extracts SL, Spain  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# LC-MS/MS of Phytoestrogens from Hop Extract


Application #AN1160

## Conditions


Column: ACE 3 C18-AR  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-119-1546  
Mobile Phase: A: 1% formic acid in MeCN  
B: 1% formic acid in MeOH  
C: 1% formic acid in H<sub>2</sub>O  
D: MeOH

Time (mins)	%A	%B	%C	%D
0	56	0	44	0
8	51	5	44	0
10	51	5	44	0
17	95	5	0	0
22	95	0	0	5


Flow Rate: 0.6 mL/min  
Detection: TSQ-Quantum triple quad ESI  
Spray voltage: -4500 V  
Precursor ion: 355.4 [M+H]<sup>+</sup>  
MRM transition ions: 179 and 299  
Collision energy: 28 and 16 V



1. Isoxanthohumol  
LOQ 0.07 µg/mL



2. Xanthohumol  
LOQ 0.01 µg/mL





# Separation of Preservatives

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2230

## Conditions


Column: ACE Excel 1.7 C18  
Dimensions: 50 x 3.0 mm  
Part Number: EXL-171-0503U  
Mobile Phase: 20 mM potassium phosphate pH 2.5 in MeCN/H<sub>2</sub>O (30:70 v/v)  
Flow Rate: 0.43 mL/min  
Injection: 0.7 µL  
Temperature: 20 °C  
Detection: UV, 230 nm




1. Phthalic acid




2. 4-Hydroxybenzoic acid




3. Sorbic acid



4. Methyl paraben



5. Ethyl paraben



6. Propyl paraben

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)


# Separation of Seven Preservatives

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN3040

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: MeCN/50 mM KH<sub>2</sub>PO<sub>4</sub> pH 4.4 in H<sub>2</sub>O (40:60 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 230 nm



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

  
**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Sennosides in Traditional Chinese Medicine

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns



Application #AN1390

## Conditions


Column: ACE 3 C18-PFP  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-1110-1546  
Mobile Phase: A: 0.75% acetic acid in H<sub>2</sub>O  
B: MeCN/MeOH (90:10 v/v)

Time (mins)	%B
0	9
23	28
40	28

Flow Rate: 0.6 mL/min  
Temperature: 35 °C  
Detection: UV, 271 nm  
Sample: Herbal tea bag containing Folium Sennae, Peppermint, Folium Mori, Folium Nelumbinis, Glycyrrhiza Uralensis and Lelang Grass Rhizome



Sennoside A



Sennoside B

Reproduced with permission of Vercotech Inc., Taiwan  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Illegal Dyes in Spices


Application #AN2910

## Conditions

Column: ACE 3 C18  
Dimensions: 100 x 4.6 mm  
Part Number: ACE-111-1046  
Mobile Phase: A: H<sub>2</sub>O  
B: MeOH  
C: 0.06 M Tetrabutylammonium bromide and 0.5 M KH<sub>2</sub>PO<sub>4</sub> in H<sub>2</sub>O pH 2.55

Time (mins)	%A	%B	%C
0	45	50	5
45	3	92	5
65	3	92	5
66	45	50	5
75	45	50	5

Flow Rate: 1 mL/min  
Injection: 10 µL  
Temperature: Ambient  
Detection: UV-Vis, 420 nm, 520 nm and 600 nm



1. Rhodamine B
2. Orange II
3. Metanil Yellow
4. Butter Yellow
5. Para Red
6. Sudan Orange G
7. Toluidine Red
8. Sudan I
9. Sudan Red G
10. Sudan II
11. Sudan Black
12. Sudan III
13. Sudan Red 7B
14. Sudan Red B
15. Sudan IV

# Separation of Sulfonamides


Application #AN1610

## Conditions

Column: ACE Excel 3 C18-PFP  
Dimensions: 150 x 4.6 mm  
Part Number: EXL-1110-1546U  
Mobile Phase: A: H<sub>2</sub>O  
B: MeCN  
C: 10% formic acid

Time (mins)	%A	%B	%C
0	84	15	1
12	74	25	1
14	59	40	1
16	84	15	1
18	84	15	1

Flow Rate: 1 mL/min  
Detection: UV, 268 nm



1. Sulfadiazine
2. Sulfapyridine
3. Sulfamerazine
4. Sulfamoxole
5. Sulfamethazine
6. Sulfamonomethoxine
7. Sulfachloropyridazine
8. Sulfamethoxazole
9. Sulfadimethoxine

Reproduced with permission of National Research Institute of Agriculture and Food Technology  
INIA, Spain

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.






# Separation of Tocopherols

Application #AN2790

## Conditions

Column: ACE 5 SIL  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-127-2546  
Mobile Phase: Hexane/IPA (98:2 v/v)  
Flow Rate: 1 mL/min  
Injection: 1  $\mu$ L  
Temperature: Ambient  
Detection: UV-Vis, 450 nm



1.  $\gamma$ -Tocopherol
2.  $\alpha$ -Tocopherol
3.  $\beta$ -Tocopherol
4.  $\beta$ -Tocopherol
5.  $\delta$ -Tocopherol





# Detection of the Herbicide Trifluralin

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2890

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: H<sub>2</sub>O/MeOH (15:85 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 254 nm



1. Trifluralin

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# BSA Tryptic Digest Profiling

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2000

## Conditions

Column: ACE 5 C18-300  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-221-1546  
Mobile Phase: A: 1% TFA in H<sub>2</sub>O  
B: 1% TFA in MeCN/H<sub>2</sub>O (1:1 v/v)

Time (mins)	%B
0	4
5	4
25	20
45	20
75	40
95	65
115	70
120	4

Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 214 nm



Reproduced with permission of Department of Food Science, University of Reading, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Separation of Vanillins

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns



Application #AN1620

## Conditions


Column: ACE 3 C18-Amide  
Dimensions: 150 x 4.6 mm  
Part Number: EXL-1112-1546U  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeCN

Time (mins)	%B
0.0	30
10.0	55
10.5	55
15.0	30
Post time 5 minutes	


Flow Rate: 1 mL/min  
Injection: 5 µL  
Temperature: 40 °C  
Detection: UV, 260 nm




1. Vanillic acid




2. 4-Hydroxybenzoic acid




3. Vanillin




4. 4-Hydroxybenzaldehyde



5. Guaiacol



6. Ethyl Vanillin



7. Eugenol

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

  
**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Fast Separation of Vanillin Compounds

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2240

## Conditions


Column: ACE Excel 1.7 C18-Amide  
Dimensions: 50 x 3.0 mm  
Part Number: EXL-1712-0503U  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeCN

Flow Rate: 1.3 mL/min  
Injection: 1 µL  
Temperature: 45 °C  
Detection: UV, 260 nm


Time (mins)	%B
0.00	25
1.32	75
1.49	75
1.60	25




1. Vanillic acid




2. 4-Hydroxybenzoic acid




3. Vanillin




4. 4-Hydroxybenzaldehyde




5. Guaiacol




6. o-Vanillin



7. Ethyl vanillin



8. Eugenol



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns



## Application #AN1830

### Conditions


Column: ACE UltraCore 2.5 SuperC18  
Dimensions: 100 x 2.1 mm  
Part Number: CORE-25A-1002U  
Mobile Phase: A: 0.01 mM ammonium fluoride + 0.001% formic acid  
B: MeCN

Time (mins)	%B
0.0	25
0.5	25
7.0	35
7.5	35
10.5	60
12.5	90

Flow Rate: 0.5 mL/min  
Temperature: 45 °C  
Detection: Positive or negative ESI  
MRM data




1.  $\alpha$ - and  $\beta$ -Zearalenol  
( $m/z$  319.17  $\rightarrow$  275.12)




2. Diethylstilbestrol-d8  
( $m/z$  275.23  $\rightarrow$  245.09)

Also analysed in -ESI:  
Talaranol and zeranone-d4  
Talaranol and zeranone  
Zearalenone  
Hexestrol  
Diethylstilbestrol  
Dienestrol



3.  $\alpha$ - and  $\beta$ -Boldenone  
( $m/z$  287.17  $\rightarrow$  121.12)





4.  $\alpha$ - and  $\beta$ -Nortestosterone  
( $m/z$  275.23  $\rightarrow$  109.09)

Also analysed in +ESI:  
Hydroxystanozolol  
Hydroxystanozolol-d3  
Methyltestosterone  
Methyltestosterone-d3  
 $\beta$ -Nortestosterone-d3  
 $\beta$ -Trenbolone  
 $\alpha$ -Trenbolone




## Conditions

Column: ACE Excel 2 C18-Amide  
Dimensions: 50 x 3.0 mm  
Part Number: EXL-1012-0503U  
Mobile Phase: 100% MeCN  
Flow Rate: 0.43 mL/min  
Injection: 2 µL  
Temperature: 20 °C  
Detection: UV, 265 nm



1. Ergocalciferol (D2)



2. Cholecalciferol (D3)

# 25-Hydroxy Vitamin D in Serum by LC-MS/MS

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns







Application #AN2390

## Conditions


Column: ACE Excel 2 C18-PFP  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1010-1002U  
Mobile Phase: A: 2 mM ammonium acetate, 0.1% formic acid in H<sub>2</sub>O  
B: 0.1% formic acid in MeOH

Time (mins)	%B
0.0	75
3.0	100
4.0	100

Flow Rate: 0.4 mL/min  
Injection: 15 µL  
Temperature: 40 °C  
Detection: Quattro Premier XE triple quad MS  
MRM positive ESI mode  
Desolvation temperature: 450 °C  
Ion source temperature: 150 °C



- 1a. 25-OH Vitamin D2  
(*m/z* 395.5 → 269.5)
- 1b. 25-OH Vitamin D2  
(*m/z* 395.5 → 119.2)
- 2. d6-25-OH Vitamin D3 (IS)  
(*m/z* 389.6 → 263.5)



- 3a. 25-OH Vitamin D3  
(*m/z* 383.5 → 257.5)
- 3b. 25-OH Vitamin D3  
(*m/z* 383.5 → 107.2)

Reproduced with permission of Biotage GB Ltd  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)




# Separation of Fat Soluble Vitamins

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


## Conditions

Column: ACE 3 C18-Amide  
Dimensions: 150 x 4.6 mm  
Part Number: EXL-1112-1546U  
Mobile Phase: MeOH/MeCN (90:10 v/v)  
Flow Rate: 1 mL/min  
Temperature: 20 °C  
Detection: UV, 280 nm


## Application #AN2420




1. Menadione (Vitamin K3)




2. Retinol (Vitamin A)




3. Vitamin A acetate




4. Ergocalciferol (Vitamin D2)




5. Cholecalciferol (Vitamin D3)




6. Vitamin E acetate



7.  $\alpha$ -Tocopherol (Vitamin E)



8. Vitamin K1



# Fast LC-MS Analysis of Multivitamin Fruit Juice

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2610


## Conditions

Column: ACE Excel 3 C18-PFP  
Dimensions: 100 x 2.1 mm  
Part Number: EXL-1110-1002U  
Mobile Phase: A: 15 mM formic acid, adjusted to pH 3.8 with ammonia solution  
B: MeOH

Time (mins)	%B
0.00	1
1.00	1
3.00	8
3.10	25
6.00	50
6.50	50
6.51	1
9.00	1


Flow Rate: 0.4 mL/min  
Temperature: 30 °C  
Detection: LCMS-8040 triple quad MS  
ESI positive mode (ESI negative for ascorbic and citric acid)  
DL temperature: 250 °C  
Heart block temperature: 400 °C

Standard



1. Thiamine (Vitamin B1)  
(*m/z* 266.10 → 122.15)
2. Pyridoxine (Vitamin B6)  
(*m/z* 170.20 → 152.15)
3. Nicotinic acid (Vitamin B3)  
(*m/z* 124.00 → 78.00)
4. Pantothenic acid (Vitamin B5)  
(*m/z* 220.30 → 90.05)
5. Cyanocobalamin (Vitamin B12)  
(*m/z* 678.50 → 147.05)
6. Riboflavin (Vitamin B2)  
(*m/z* 377.20 → 243.10)
7. Biotin (Vitamin B7)  
(*m/z* 245.10 → 227.05)
8. Nicotinamide (Vitamin B3)  
(*m/z* 123.20 → 80.05)
9. Ascorbic acid (Vitamin C)  
(*m/z* 175.10 → 114.80)
10. Citric acid  
(*m/z* 191.10 → 87.15)

Juice sample



Reproduced with permission of Shimadzu Europe GmbH  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Separation of 14 Water Soluble Vitamins

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2940

## Conditions

Column: ACE 5 C8  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-122-2546  
Mobile Phase: A: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 2.5 in H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0.0	0
3.0	0
16.5	45
19.5	80

Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 205 nm



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

# Separation of 13 Water Soluble Vitamins using a Gradient

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN2930

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: A: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 3.0 in H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0	3
5	3
45	45
50	80

Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 205 nm



1. Pyridoxamine
2. Thiamine (Vitamin B1)
3. L-Ascorbic acid (Vitamin C)
4. Nicotinic acid
5. Pyridoxal
6. Impurity
7. Pyridoxine
8. p-Aminobenzoic acid
9. Pantothenic acid (Vitamin B5)
10. Folic acid (Vitamin B9)
11. Cyanocobalamin (Vitamin B12)
12. d-Biotin (Vitamin B7)
13. Riboflavin (Vitamin B2)




# Separation of Seven Water Soluble Vitamins using Isocratic Conditions (I)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2990

## Conditions

Column: ACE 5 C18  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-121-2546  
Mobile Phase: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 3.0 in H<sub>2</sub>O/MeOH (97:3 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 205 nm



1. Pyridoxamine
2. Thiamine (Vitamin B1)
3. L-Ascorbic acid (Vitamin C)
4. Nicotinic acid
5. Pyridoxal
6. Impurity
7. Pyridoxine




# Separation of Seven Water Soluble Vitamins using Isocratic Conditions (II)

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN2980

## Conditions

Column: ACE 5 C8  
Dimensions: 250 x 4.6 mm  
Part Number: ACE-122-2546  
Mobile Phase: 50 mM KH<sub>2</sub>PO<sub>4</sub> pH 2.5 in H<sub>2</sub>O/MeOH (97:3 v/v)  
Flow Rate: 1 mL/min  
Temperature: Ambient  
Detection: UV, 205 nm



© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)



## Application #AN1870

### Conditions


Column: ACE 3 C18-AR  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-119-1546  
Mobile Phase: A: 20 mM potassium phosphate pH 2.83 in H<sub>2</sub>O  
B: 20 mM potassium phosphate pH 2.83 in MeOH/H<sub>2</sub>O (50:50 v/v)

Flow Rate: 1.5 mL/min  
Injection: 1 µL  
Temperature: 40 °C  
Detection: UV, 205 nm


Time (mins)	%B
0	20
15	70




1. Pyridoxine  
(Vitamin B6)




2. Pantothenic acid  
(Vitamin B5)




3. p-Aminobenzoic acid




4. Folic acid  
(Vitamin B9/ Vitamin M)



5. D-Biotin  
(Vitamin B7/ Vitamin H)



6. Cyanocobalamin  
(Vitamin B12)



7. Riboflavin  
(Vitamin B2)



## Application #AN1880

### Conditions


Column: ACE Ultracore 2.5 SuperPhenylHexyl  
Dimensions: 50 x 2.1 mm  
Part Number: CORE-25B-0502U  
Mobile Phase: A: 20 mM KH<sub>2</sub>PO<sub>4</sub>, pH 2.7  
B: 20 mM KH<sub>2</sub>PO<sub>4</sub>, pH 2.7 in MeOH/H<sub>2</sub>O (50:50 v/v)

Time (mins)	%B
0.00	20
1.50	60
3.00	70
3.75	70
4.50	20
Post time 4.5 minutes	


Flow Rate: 0.4 mL/min  
Injection: 1 µL  
Temperature: 40 °C  
Detection: UV, 205 and 254 nm




1. Pyridoxine




2. Pantothenic acid




3. p-Aminobenzoic acid




4. Cyanocobalamin



5. d-Biotin



6. Folic acid



7. Riboflavin




# Water Soluble Vitamins in Green Vegetables by LC-MS/MS

## Conditions


Column: ACE 3 C18  
Dimensions: 100 x 2.1 mm  
Part Number: ACE-111-1002  
Mobile Phase: A: 10 mM ammonium acetate in H<sub>2</sub>O, pH 4.5  
B: 0.1% acetic acid in MeOH  
C: 0.3% acetic acid in MeOH

Time (mins)	%A	%B	%C
0	90	10	0
3	90	10	0
4	50	0	50
7	50	0	50
10	0	100	0


Flow Rate: 0.2 mL/min  
Injection: 10 µL  
Temperature: 20 °C  
Detection: TSQ triple quad MS; SRM mode  
-ESI for vitamin C  
+ESI for B vitamins




## Application #AN1860




1. Ascorbic acid




2. Thiamine




3. Pyridoxine




4. Nicotinamide




5. Pantothenic acid



6. Hippuric acid (IS)



7. Folic acid



8. Riboflavin

Reproduced with permission of Foodomics Laboratory, CIAL-CSIC, Madrid, Spain  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.


# Water Soluble Vitamins and Polar Molecules

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN1850

## Conditions

Column: ACE 3 C18-AR  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-119-1546  
Mobile Phase: 0.1 % phosphoric acid in H<sub>2</sub>O/MeOH (96.5:3.5 v/v)  
Flow Rate: 1 mL/min  
Injection: 2 µL  
Temperature: 22 °C  
Detection: UV, 260 nm



1. Pyridoxamine (Vitamin B6)
2. Thiamine (Vitamin B1)
3. Isonicotinamide
4. Nicotinamide
5. L-Ascorbic acid (Vitamin C)
6. Orotic Acid
7. Hypoxanthine
8. Pyridoxal (Vitamin B6)
9. Pyridoxine (Vitamin B6)
10. p-Aminobenzoic acid

© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)

# Whey Proteins from Whole Milk

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns


Application #AN3000

## Conditions

Column: ACE 3 C4-300  
Dimensions: 150 x 2.1 mm  
Part Number: ACE-213-1502  
Mobile Phase: A: 0.5% formic acid in H<sub>2</sub>O  
B: 0.5% formic acid in MeCN

Time (mins)	%B
0	35
16	43
17	80
20	80
21	35
31	35

Flow Rate: 0.4 mL/min  
Injection: 10 µL  
Temperature: 40 °C  
Detection: ESI-MS (+ve)



1.  $\alpha$ -Lactalbumin
2.  $\beta$ -Lactoglobulin B
3.  $\beta$ -Lactoglobulin A

Reproduced with permission of University of Vienna, Austria  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.



For more information contact your local ACE distributor or visit  
[www.ace-hplc.com](http://www.ace-hplc.com) or email: [info@ace-hplc.com](mailto:info@ace-hplc.com)




**ACE**<sup>®</sup>  
UHPLC & HPLC  
Columns

## Conditions


Column: ACE Excel 3 C18-Amide  
Dimensions: 250 x 2.1 mm  
Part Number: EXL-1112-2502U  
Mobile Phase: 40 mM ammonium phosphate pH 2.5 in H<sub>2</sub>O  
Flow Rate: 0.21 mL/min  
Injection: 5 µL  
Temperature: 25 °C  
Detection: UV, 214 nm




1. Oxalic acid




2. Tartaric acid




3. Malic acid



4. Lactic acid



5. Ascorbic acid



6. Citric acid



For further details please contact:




**Hichrom Limited**

Tel: +44 (0)118 930 3660

Fax: +44 (0)118 932 3484

[technical@hichrom.co.uk](mailto:technical@hichrom.co.uk)

[www.hichrom.co.uk](http://www.hichrom.co.uk)



**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: [info@mz-at.de](mailto:info@mz-at.de), [www.mz-at.de](http://www.mz-at.de)