



Determination of aromatic hydrocarbons in diesel fuel

The aromatic hydrocarbons content of diesel fuel is a factor that can affect exhaust emissions and fuel combustion characteristics, as measured by cetane number. This application shows fast and robust method for the determination of monoaromatic, di-aromatic, and polyaromatic hydrocarbon contents in diesel fuels and petroleum distillates boiling in the range from 150 to 400 °C.



Calibration standard solution on ARION® column




Test solution on ARION® column



Determination of aromatic hydrocarbons in diesel fuel

Column	ARION® NH ₂ , 5 µm
Dimensions	150 mm × 4.6 mm
Part number	ARI-5736-LK46
Mobile phase	Heptane Isocratic elution
Flow rate	0.8 mL/min
Temperature	30 °C
Injection volume	10 µL
Detection	RID
Analytes	1. Non-aromatic hydrocarbons (NAHs) 2. Mono-aromatic hydrocarbons (MAHs) 3. Di-aromatic hydrocarbons (DAHs) 4. Polycyclic aromatic hydrocarbons (PAHs)



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de

This application was developed
in co-operation with ORLEN UniCRE a.s.

