

How to Choose a Column For Your Sample

In choosing the best column for your application, there are several factors to consider:

- Resolution of peaks of interest
- Analysis time
- Selectivity (elution order of peaks)
- Durability

In general, the **Coregel 87H3** is a good starting point when considering a Transgenomic organic acids analysis column since it combines good resolution of many common organic acids, with high durability. However, the many combinations of polymer cross-linkage, particle size and column sizes offered by Transgenomic allow customers to choose a column with features most important to maximize their separation needs.


The chart below is a general guideline for choosing your Transgenomic column:

Resolution (Highest to Lowest)	Analysis Time (Faster to Slower)	Selectivity (Higher to Lower)	Durability (More to Less)
ICSep ION300	ICSep USP L-17	ICSep ION300	ICSep Coregel 107H
ICSep Coregel 64H	ICSep Coregel 87H1	ICSep Coregel 64H	ICSep Coregel 87H3
ICSep WA-1	ICSep ION310	ICSep ORH801	ICSep Coregel 87H1
ICSep Coregel 87H3	ICSep ARH601	ICSep WA-1	ICSep WA-1
ICSep ORH801	ICSep ORH801	ICSep Coregel 87H3	ICSep ORH801
ICSep Coregel 107H	ICSep Coregel 107H	ICSep Coregel 107H	ICSep ION310
ICSep ARH601	ICSep Coregel 87H3	ICSep ARH601	ICSep USP L-17
ICSep ION310	ICSep WA-1	ICSep ION310	ICSep ARH601
ICSep Coregel 87H1	ICSep Coregel 64H	ICSep Coregel 87H1	ICSep Coregel 64H
ICSep USP L-17	ICSep ION300	ICSep USP L-17	ICSep ION300

Another valuable tool for selecting a column is to use retention charts of many common organic acids, sugars, and alcohols, and the innumerable possible combinations of the compounds in mixtures, it is impossible to have a totally comprehensive retention chart to conclusively guide a chemist through all applications. Since we stock most of the common organic acids, sugars and alcohols in our laboratory, we frequently will generate a sample chromatogram of your particular sample mixture, or will actually run your sample to ensure that our columns will separate your sample. Please contact Transgenomic for any questions regarding column selection.

Compound	ION300	ORH801	Coregel 64H	Coregel 107H	Coregel 87H3
1 Malic	12.0	7.5	11.5	10.2	10.0
2 Malonic	11.1	7.4	10.8	9.8	9.8
3 cis-Aconitic	7.5	7.4	11.2	9.9	9.6
4 Adipic	21.8	13.7	21.9	17.5	16.5
5 Formic	17.1	11.1	16.5	15.0	14.6
6 Maleic	8.1	5.2	7.7	7.6	7.8
7 Ascorbic	7.8	5.0	7.5	7.4	10.5
8 Butyric	—	18.5	—	24.1	—
9 Glycolic	15.4	9.9	14.8	13.4	12.8
10 Glycoxalic	11.5	7.4	11.1	10.1	9.8
12 Citric	9.5	6.3	9.1	8.4	8.3
13 Tartaric	10.0	6.4	9.6	8.9	8.7
14 Nicotinic	—	—	—	—	26.4
15 Propionic	20.1	15.2	21.7	19.3	18.6
16 Succinic	9.2	9.8	8.8	12.7	8.7
18 Oxalic	6.6	4.5	6.5	6.7	6.9
19 Sorbic	18.7	12.0	17.9	16.0	15.3
23 Acrylic	23.3	—	22.5	19.9	19.4
26 Isobutyric	—	—	—	21.8	21.4
27 Lactic	16.0	10.3	15.5	13.6	13.1
28 Shikimic	15.5	9.7	14.7	12.5	12.1
29 Fumaric	16.4	10.3	15.7	14.0	13.1
30 Glutaric	18.2	11.3	17.2	14.7	14.0
31 Pyruvic	9.3	6.0	9.0	8.6	8.7
32 Acetic	—	12.4	18.4	16.6	—
33 Proponal	28.7	28.9	—	—	—
36 Quinic	12.6	7.9	12.1	10.7	10.4

Test Conditions: 2.5 mM H₂SO₄, 50 °C, Flow Rate 0.5 mL/minute, 100 ppm



MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

AUTHORIZED DISTRIBUTOR

Transgenomic Offers More Column Options

Transgenomic has a wider variety of polymeric HPLC columns than our competitors. The vast majority of columns on the market for ion-exclusion analysis of organic acids utilize an 8% cross-linked polystyrene-divinylbenzene co-polymer. Although the 8% cross-linked polymeric columns provide an excellent separation of many typical samples containing organic acids, sugars and alcohols, Transgenomic has determined that by varying the cross-linkage and particle sizes of the polymers, certain separations are enhanced tremendously. Although ion-exclusion of organic acids is the main separation mechanism for samples, the size exclusion limit of the polymers is controlled by the percent cross-linking. The additional size exclusion mechanism for Transgenomic columns provides unique selectivity for some components. For maximum efficiency, Transgenomic also provides columns with smaller particle sizes and tighter size distributions than the typical industry standard of 9 micron particles. By providing additional choices, Transgenomic allows you to maximize separations in terms of speed, selectivity and resolution of the sample components.


The chart below illustrates the variety of columns available from Transgenomic compared to similar columns offered by other companies:

Transgenomic	6% xI	6.4% xI	7% xI	8% xI	10% xI	Particle Size (uM)
ICSep ARH601	x					6.5
ICSep Coregel 107H					x	8
ICSep Coregel 64H		x				10
ICSep Coregel 87H1				x		9
ICSep Coregel 87H3				x		9
ICSep ION300	x					7
ICSep ION310				x		8
ICSep ORH801			x			9
ICSep USP L-17				x		8
ICSep WA-1				x		9
Bio-Rad						
Fast Analysis Column				x		9
Fermentation Monitoring				x		9
Aminex HPX-87H				x		9
Phenomenex						
Rezex RFQ-Fast Acid				x		8
Rezex RHM-Monosaccharide				x		8
Rezex ROA-Organic Acid				x		8
Shodex						
KC811				N/A		6
Sugar SH1011				N/A		6
Sugar SH1821				N/A		6
Supelco						
SUPELCOGEL™ HPN 59304-U				N/A		9
SUPELCOGEL™ 6-10H				N/A		9
SUPELCOGEL™ HPN 59346				N/A		9


New Applications

The following are new methods developed using our **ICSep** columns based on customer application requests.


Analysis of 1,3 – propanediol in the presence of other organic acids and methanol on the **Coregel 87H3** column




Separation of formaldehyde in the presence of methanol and formic acid on the **Coregel 87H3** column. UV detection at 210 nm.



Analysis of thiodiglycol **Coregel 87H3** column



Analysis of four sugars in the presence of two organic acids with the **ION300** column



Separation of galactose and the aminosugar glucosamine with the **Coregel 87H3** column

