

## Concise Separations Columns for Amino Acid Analysis

- Rugged polymeric substrate, stable in pH range of 0 to 14
- High efficiency and resolution
- Reproducibility lot-to-lot and column-to-column
- Available for both physiological samples (Li<sup>+</sup> format) and protein hydrolysate samples (Na<sup>+</sup> form)
- Post column derivitization detection

Ion-exchange chromatography is a popular technique for the analysis of amino acids because both retention times and quantification are highly reproducible regardless of the sample matrix. This unique matrix insensitivity is important when comparing results from different patients or batches of protein hydrolysate.

Amino acids are zwitterions; at low pH, they are positively-charged and are bound to the resin by their attraction to the negatively-charged ion-exchange sites. Almost all the contaminants, i.e. matrix, are eluted at the void. The amino acids are then selectively eluted by increasing the pH and salt concentration with different buffers. With few exceptions, the order of elution follows the isoelectric point of the amino acids, i.e. acidic amino acids first, then neutral and basic. Because the separation and the ensuing post-column reaction of amino acids are devoid of contaminants, amino acid analyses via ion-exchange chromatography are highly reproducible.


**CONCISE SEPARATIONS AMINO ACID ANALYSIS COLUMN COMPARISON CHART**

PHASE	CROSS-LINKAGE	IONIC FORM	PARTICLE SIZE (µm)	KEY SAMPLES	COMMENTS
AMINOSep AA-911	8	Sodium	9	Designed for complicated samples from protein hydrolysates. Increased polymer bed yields better resolving power	Popular column—Higher capacity than the Beckman columns
AMINOSep AA-511	10	Sodium	5	Designed for faster analysis than the AA-911 but still gives high resolution	Most popular AA column
Lithium Amino Acid (6300 & 7300 systems)	10	Lithium	6	Designed for use with the Beckman Coulter® 6300 and 7300 Amino Acid Analyzers using either the Beckman or Pickering Lithium buffer systems	Ideal for Physiological amino acid analysis
Sodium Amino Acid (6300 & 7300 systems)	10	Sodium	5	Designed for use with the Beckman Coulter® 6300 and 7300 Amino Acid Analyzers using either the Beckman Coulter or Pickering Sodium buffer systems	Ideally suited for routine hydrolysate analysis, Extremely rugged polymer
Sodium Amino Acid (for System Gold)	10	Sodium	5	Designed for use with the Beckman Coulter® System Gold Amino Acid Analyzer	Ideal for the separation of hydrolysate amino acids.

DESCRIPTION	SIZE	AMINOSep AA-911	AMINOSep AA-511	LITHIUM AMINO ACID (6300/7300)	SODIUM AMINO ACID (6300/7300)	SODIUM AMINO ACID (System Gold)
Column	4.6 x 250 mm	AAA-99-8553	—	—	—	—
Column (Waters)	4.6 x 250 mm	AAA-99-8553W	—	—	—	—
Column	4.6 x 150 mm	—	AAA-99-7554	—	—	—
Column	4.6 x 120 mm	—	AAA-99-6554	—	—	—
Column	4.0 x 200 mm	—	—	—	—	AAA-99-6310
Column	4.0 x 120 mm	—	—	—	AAA-99-6312	—
Column	4.0 x 100 mm	—	—	AAA-99-6311	—	—
Guard Kit (1 holder, 2/pk cartridges)		AAA-99-2353	AAA-99-2354	AAA-99-2311	AAA-99-2312	AAA-99-2312
Guard Cartridges (2/pk)		AAA-99-1353	AAA-99-1354	AAA-99-1311	AAA-99-1312	AAA-99-1312


## AMINOSep AA-511 Sodium Column

0.5 ml/min with Sodium Citrate A, B, C at 82 °C, Fluorescence Detector




## AMINOSep AA-511 High Speed Sodium Column

0.5 ml/min with Sodium Citrate A, B, C at 60 °C, Fluorescence Detector




## Lithium Amino Acid Column

20 mL/hr with LiA, LiB, LiC at 33-60-77°C, Fluorescence Detector




## AMINOSep AA-911 Sodium Column

0.5 ml/min with Sodium Citrate A, B, C at 82 °C, Fluorescence Detector



## Sodium Amino Acid Column (System Gold)

26 ml/hr with NaE, NaF, NaD at 50-65-77°C, Fluorescence Detector



## Sodium Amino Acid Column (6300/7300 systems)

26 ml/hr with NaE, NaF, NaD at 50-65-77°C, Fluorescence Detector

