


# $\mu$ SPEed<sup>®</sup> Cartridges




## ...more than micro-SPE

Disposable micro separation cartridges for ultra clean  $\mu$ SPE, fractionation and HPLC separations

# μSPEed Cartridges (with valve)

μSPEed cartridges offer a unique opportunity to revolutionise micro SPE and micro fractionation. Using a one-way check valve, the sample is aspirated through the cartridge and analytes focused on to the top of a sorbent bed. Using an analytical syringe, relatively high pressures can be generated meaning sorbents of  $\leq 3\mu\text{m}$  can be used.

The μSPEed cartridges perform like a short HPLC column where separation occurs not just during digital extraction but in the sample clean-up. This results in cleaner sample extractions, compound fractionation, high concentration factors and even LC isocratic and stepped elution for Direct MS Infusion directly from the cartridge. *(Patented)*



## FEATURES

- Achieves narrow, low-volume, high concentration elution bands
- Higher concentration factors with narrow, low-volume elution bands enable analyte elution in just  $\mu\text{L}$  of solvent. Concentration factors of 1000:1 can be achieved in minutes
- Improved reproducibility with no sample pass through
- Does not have the capacity limitations of SPME
- Reduces or eliminates solvent blowdown – cost efficient
- Small particle media bed for cleaner extracts and better analyte separation capability
- Can be configured for "analytical" separation, allowing direct connection to an analytical detector (eg. Mass Spectrometer)
- Full samples can be transferred to and from the syringe/fitting, reducing the risk of cross-contamination



## TYPICAL OPERATION SEQUENCE

STEP 1: [ACTIVATION] Aspirate and Dispense Conditioning Solvent (min 50 $\mu\text{L}$ )


STEP 2: [CONDITIONING] Aspirate and Dispense (100 $\mu\text{L}$ ) Conditioning Solvent (min 50 $\mu\text{L}$ )

STEP 3: [SAMPLE LOAD] Aspirate Sample to Trap Analyte(s)

STEP 4: [SAMPLE DISPENSE] Dispense Sample and needle Wash to Trap Analyte(s)

STEP 5: [WASH] Aspirate and Dispense Wash (typically same as Conditioning) Solvent

STEP 6: [ELUTION] Aspirate and Dispense Elution Solvent(s)



## Multiple Use

μSPEed cartridges can be used multiple times depending on sample matrix and SOP requirements. To reuse, aspirate and dispense 100 $\mu\text{L}$  organic conditioning solvent multiple times prior to beginning the next sequence.

## AUTOMATED HIGH PRESSURE SYRINGE CONNECTION



































The cartridges feature a high pressure, low dead volume connection for fast attachment and disconnection of the cartridge from the syringe and **automated operation**.

- Eprep's automated high-pressure connector allows a seal up to 1500psi. Full sample volume can be delivered from the syringe due to the ultra-low dead volume.
- The connector is ideal for robotic automation applications. The simple push-pull fittings of male and female components connect and disconnect without twisting or turning the parts. **There is no need for complex automation apparatus, processes or programming.** During a sequenced workflow operation, multiple fittings can be easily picked up and dropped off at the required stations.



High Pressure Connection



# μSPEed vs Other SPE Methods

	μSPEed	SPE	MEPS/SPEmx	SPME
<b>Sorbent Size</b>	 ≤3μm	 40-60μm	 40-60μm	 Coated Fibre
<b>Typical Sample Volume</b>	 10μL-10mL	 2-3mL	 50μL-5mL	 2mL
<b>Time/Speed</b>	 Very Fast <sup>1</sup>	 Very Slow	 Fast	 Slow
<b>Extraction Efficiency</b>	 Excellent <sup>2</sup>	 Poor	 Very Poor	 Very Poor
<b>Concentration Factor</b>	 Very High(x2-x1000)	 Very Low(x1-x10)	 High	 Low
<b>Price</b>	 <\$1 (multi use)	 \$3-\$5	 <\$1 (multi use)	 \$2
<b>Total Solvent</b>	 200μL	 10mL	 500μL	n/a
<b>Evaporation Step</b>	 Never Required	 Required	 Rarely Required	n/a
<b>Automation</b>	 Yes	 Difficult	 Limited	 Yes


<sup>1</sup> μSPEed is faster than conventional SPE. Because smaller particle size sorbents are used in μSPEed, it performs like a high resolution HPLC column with elution in a narrow precise band. Resolution of targeted compounds in conventional SPE is very limited due to large particle size and inefficient packing of these sorbents.

<sup>2</sup> Smaller particle size sorbent gives greater efficiency of extraction, with elution conditions becoming far less critical. Also the kinetics of the extraction is less critical with the greater efficiency. Where conventional large particle SPE requires precise conditions, μSPEed is far less method critical leading to greater reproducibility.

## μSPEed Application Examples




### PFAS by Automated $\mu$ SPEd/PFAS




5  $\mu$ g/L of PFBA, PFBS, PFHxA, PFHxS, PFOA and PFOS that has undergone a ten-times pre-concentration

### Recovery of Drug Mix from Serum and Blood




Selective Ion Monitoring at 0.02ppm after  $\mu$ SPEd extraction from serum and blood

### $\mu$ SPEd Plasma Opiate Standard Separation




Direct Cartridge to MS of Opiate Standards

### $\mu$ SPEd Clean-up Blood Metanephrines




Extraction of Metanephrines in plasma using  $\mu$ SPEd

### On-column Derivatisation of Fatty Acids




Stacked chromatograms of  $\mu$ SPEd derivatised FAMES in olive oil

### $\mu$ SPEd extraction of polyphenols in teas




$\mu$ SPEd traction of 8 selected phenolics compounds from tea at 0.1-1 $\mu$ g/mL

### $\mu$ SPEd Recovery of 4-nitrotoluene



$\mu$ SPEd extraction of 4-nitrotoluene at 101-110% recoveries and five point calibration with a linear regression of 0.9994

### 2min Trypsin Digest of Bovine Serum Albumin (BSA)



Overlaid chromatograms at 215 nm of undigested BSA (blue), and the resulting 2 minute BSA digestion on trypsin  $\mu$ SPEd-Cxyl

## μSPEed Cartridges (valve)

### Ordering Information

Part No	Code	Description
μSPEed Cartridges		
<b>Silica Based</b>		
01-10110	μSPEed, C18RPS-3μm/120Å (Pkt 10)	3μm/ 120Å ODS spherical silica packing with high acidic resistance suitable for general organic compound applications.
01-10115	μSPEed, Silica-3μm/120Å (Pkt 10)	3μm/120Å spherical bare silica packing. High purity silica for normal and hiliic applications
<b>Speciality Silica Based</b>		
01-10118	μSPEed, PFAS-3μm/120Å (Pkt 10)	3μm/120Å PFAS spherical silica packing. PFAS specific applications
01-10185	μSPEed, Cxyl-3μm (Pkt 10)	3μm Carboxyl spherical inert silica packing. Customisable chemistry applications
<b>Polymer Based</b>		
01-10150	μSPEed, PS/DVB -3μm/ 300Å (Pkt 10)	3μm/ 300Å spherical, crosslinked polystyrene divinyl benzene
01-10151	μSPEed, PS/DVB RP-3μm/ 300Å (Pkt 10)	3μm/ 300Å Phenyl (RP) spherical, crosslinked polystyrene divinyl benzene
01-10155	μSPEed, PS/DVB SAX-3μm/ NP (Pkt 10)	3μm/Non-Porous SAX spherical, crosslinked polystyrene divinyl benzene
01-10156	μSPEed, PS/DVB SCX-3μm/ NP (Pkt 10)	3μm/Non Porous SCX spherical, crosslinked polystyrene divinyl benzene

## SPEmx Cartridges (without valve)

SPEmx cartridges are packed with 40-60 micron sorbents and do not include the one-way check valve. Sample is loaded from the bottom of the sorbent bed making them ideal for “dirty” matrices such as blood, plasma and saliva. Loading mechanism makes the separation process digital as per conventional MEPS/SPE.



### Ordering Information

Part No	Code	Description
SPEmx Cartridges		
01-10205	SPEmx, C4-Silica/SPE (Pkt 10)	40-60μm C4 spherical silica
01-10209	SPEmx, C18-Silica/SPE (Pkt 10)	40-60μm C18 spherical silica


## Detachable Needles for μSPEed Syringes

Needles fitted with a high-pressure connection hub can be used with μSPEed syringe for automated operation on ePrep. Needles are used in a Workflow for syringe filling, priming and washing prior to cartridge or filter connection.




### Ordering Information

Part No	Code	Description
Filters		
01-10990	Ndle-Hub (Pkt 25)	Needle μSPEed Hub 50 26g Domed (Pkt 25)



[www.eprep.com.au/uspeed](http://www.eprep.com.au/uspeed)

Eprep Pty Ltd  
14/35 Dunlop Rd  
Mulgrave VIC 3170  
AUSTRALIA  
(T) +61-(0)3-9574 3600  
(E) [info@eprep.com.au](mailto:info@eprep.com.au)



**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20

e-mail: [info@mz-at.de](mailto:info@mz-at.de)

[www.mz-at.de](http://www.mz-at.de)

ePrep® is a patented registered trademark of EPREP Pty Ltd.  
μSPEed® is a patented registered trademark of EPREP Pty Ltd.  
μXact<sup>3</sup>® is a patented registered trademark of EPREP Pty Ltd.  
Other trademarks are trademarks of their respective owners

*Specifications are subject to change without notice*

Publication No.94-12002-02 Mar 2018  
© Copyright 2018, Eprep Pty Ltd