

Applications

AUTHORIZED DISTRIBUTOR


MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de

ProteoSil Application

Peptide and Protein



Analyte:

| | |
|---------------------------------|----------|
| 1. Angiotensin II | 100 mg/L |
| 2. Conalbumin | 100 mg/L |
| 3. α -Chymotrypsinogen A | 300 mg/L |
| 4. Ovalbumin | 500 mg/L |

Conditions


System : HPLC
Column : ProteoSil 200-C18 UHLC PEEK (1.9 μ m, 50 x 2.1 mm I.D)
Eluent : A) 0.1% HCOOH in H₂O
 B) 0.1% HCOOH in CH₃CN

| Time (min) | B% |
|------------|----|
| 0 | 30 |
| 1.5 | 80 |
| 1.51 | 30 |
| 4.5 | 30 |

Flow Rate : 0.6 mL/min
Col. Temp. : 40 °C
Detection : UV 280 nm
Injection Vol. : 10 μ L
Sample : Protein



Peptide



Analyte:

| | |
|--------------------------|---------|
| 1. Gly-Tyr | 50 mg/L |
| 2. Val-Tyr-Val | 50 mg/L |
| 3. Angiotensin II | 50 mg/L |
| 4. Methionine enkephalin | 50 mg/L |
| 5. Leucine enkephalin | 50 mg/L |

Conditions


System : HPLC
Column : ProteoSil 200-C18 (1.9 μ m, 100 x 2.1 mm I.D)
Eluent : A) 0.1% HCOOH in H₂O
 B) 0.1% HCOOH in CH₃CN

| Time(min) | B% |
|-----------|----|
| 0 | 5 |
| 0.5 | 30 |
| 3.0 | 40 |
| 3.5 | 40 |
| 3.51 | 5 |
| 10.0 | 5 |

Flow Rate : 0.3 mL/min
Col. Temp. : 40 °C
Detection : UV 280 nm
Injection Vol. : 5 μ L
Sample : Protein

Applications

Analysis of Peptides and Proteins




Conditions

Column : ProteoSil 200-C8 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) 0.1% TFA in CH₃CN
 B) 0.1% TFA in H₂O
 A/B = 20/80 – 20 min – 55/s45, v/v
 Flow Rate : 1.5 mL/min
 Col. Temp. : 40 °C
 Detection : UV 220 nm
 Injection Vol. : 5 μ L

Analyte

1. Ribonuclease A (0.2 mg/mL)
 2. Insulin (0.2 mg/mL)
 3. Lysozyme (0.2 mg/mL)
 4. BSA (0.2 mg/mL)

Analysis of Protein (ProteoSil 300-C8)



Conditions


Column : ProteoSil 300-C8 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) CH₃CN/0.05 % TFA in H₂O = 80/20, v/v
 B) CH₃CN/0.05 % TFA in H₂O = 10/90, v/v
 A/B = 0/100 – 30 min – 100/0, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm
 Injection Vol. : 20 μ L

Analyte

1. DL-Phenylalanine (1.01 mg/mL)
 2. Cytochrome C (0.11 mg/mL)
 3. Lysozyme (0.07 mg/mL)
 4. BSA (0.21 mg/mL)
 5. α -Chymotrypsinogen A (0.08 mg/mL)
 6. Ovalbumin (0.30 mg/mL)

Applications

Analysis of Protein (Proteosil 300-SEC)




Conditions

Column : Proteosil 300-SEC (5 μ m, 250 x 4.6 mm I.D.)
 Eluent : 0.1 M Na₂HPO₄ (pH 6.9, NaH₂PO₄)
 Flow Rate : 0.5 mL/min
 Col. Temp. : 30 °C
 Detection : UV 220 nm
 Injection Vol. : 20 μ L

Analyte

1. Thyrogloblin (0.25 mg/mL)
 2. BSA (0.25 mg/mL)
 3. Insulin Chain A (0.25 mg/mL)
 4. Oxytocin (0.21 mg/mL)

Analysis of Peptides and Proteins



Conditions


Column : ProteoSil 300-C18 (5 μ m, 150 x 4.6 mm I.D.)
 Eluent : A) 0.05% TFA in (CH₃CN/H₂O = 90/10, v/v)
 B) 0.05 % H₂O
 A/B = 20/80 – 40 min – 40/60, v/v
 Flow Rate : 1.0 mL/min
 Col. Temp. : 30 °C
 Detection : UV 280 nm
 Injection Vol. : 20 μ L

Analyte

1. Oxytocin (0.05 mg/mL)
 2. Methionine Enkephalin (0.11 mg/mL)
 3. Leucine Enkephalin (0.11 mg/mL)
 4. Angiotensin II (0.05 mg/mL)
 5. Angiotensin I (0.16 mg/mL)
 6. Insulin (0.05 mg/mL)
 7. Insulin Chain B (0.10 mg/mL)

Applications

Analysis of Peptides and proteins (Proteosil 300-C18)




Conditions

Column : Proteosil 300-C18 (5 μ m, 150 x 4.6 mm I.D.)
Eluent : A) CH₃CN/0.05 % TFA = 80/20, v/v
B) CH₃CN/0.05 % TFA = 10/90, v/v
A/B = 0/100 - 30 min - 100/0 - 10 min - 100/0, v/v
Flow Rate : 1.0 mL/min
Col. Temp. : 30 °C
Detection : UV 280 nm

Analyte

1. DL-Phenylalanine (FW 165)
2. Cytochrome C (FW 13,000)
3. Lysozyme (FW 14,000)
4. BSA (FW 66,000)
5. α -Chymotrypsinogen A (FW 25,600)
6. Ovalbumin (FW 45,000)

Analysis of Peptides and proteins (Proteosil 300-C8)



Conditions


Column : Proteosil 300-C8 (5 μ m, 150 x 4.6 mm I.D.)
Eluent : A) 0.1 % TFA in (CH₃CN/0.1 % TFA = 90/10, v/v)
B) 0.1 % TFA
A/B = 20/80 - 25 min - 60/40 - 5 min - 60/40, v/v
Flow Rate : 1.0 mL/min
Col. Temp. : 30 °C
Detection : UV 280 nm

Analyte

1. Ribonuclease A (FW 13,700)
2. Insulin (FW 6,000)
3. Cytochrome C (FW 13,000)
4. Lysozyme (FW 14,000)
5. BSA (FW 66,000)
6. STI (FW 21,000)
7. Ovalbumin (FW 45,000)

Applications

Analysis of BSA Digests




Conditions

Column : Proteosil 200-C18
(1.9 μm , 150 x 2.1 mm I.D.)
Eluent : A) 0.1% TFA in CH_3CN
B) 0.1% TFA in H_2O
A/B = 10/90 - 30 min - 50/50 - 0.1 min - 90/10
- 5 min - 90/10 - 0.1 min - 10/90 - 15 min
Flow Rate : 0.2 mL/min
Col. Temp. : 40 °C
Detection : UV 210 nm
Injection Vol. : 10 μL

Analyte

Tryptic Digest of BSA (0.5 mg/mL)

Analysis of Peptides and proteins (Proteosil 300-C4)



Conditions


Column : Proteosil 300-C4 (5 μm , 150 x 4.6 mm I.D.)
Eluent : A) 0.2 % HCOOH in ($\text{CH}_3\text{CN}/\text{H}_2\text{O}$ = 90/10, v/v)
B) 0.2 % HCOOH
A/B = 20/80 - 20 min - 80/20, v/v
Flow Rate : 1.0 mL/min
Col. Temp. : 30 °C
Detection : UV 280 nm

Analyte

1. Neurotensin (FW 1673)
2. Leucin Enkephalin (FW 556)
3. Cytochrome C (FW 12,000)
4. Insulin (FW 6,000)
5. BSA (FW 66,000)
6. Myoglobin (FW 17,000)
7. Creatine amidinohydrolase (FW 43,000)
8. Ovalbumin (FW 45,000)
9. Creatinine amidohydrolase (FW 170,000)

Applications

Analysis of Oligonucleotides



Conditions


Column : ProteoSil 200-C18
(3 μ m, 100 x 3.0 mm I.D.)
Eluent : A) 5 mM TEAA in H₂O (pH 6.5)/CH₃CN
= 80/20, v/v
B) 5 mM TEAA in H₂O (pH 6.5)
A/B = 5/95 – (60 min) – 50/50, v/v
Flow Rate : 0.8 mL/min
Col. Temp. : 40 °C
Detection : UV 260 nm
Injection Vol. : 1 μ L
Sample : Standard

Analyte

1. CATGACGTTCTGATGCT
(18 mer, M.W. 5465.61)
2. CCATGACGTTCTGATGCT
(19 mer, M.W. 5754.79)
3. TCCATGACGTTCTGATGCT
(20 mer, M.W. 6058.99)

Applications

Analysis of IgG




Conditions

Column Cat. No.: 5020-10818
Column : MonoSelect RP-mAb (20 x 2.1 mm I.D.)
A) 0.1%TFA in CH₃CN
Eluent : B) 0.1%TFA in H₂O
A/B = 5/95 – 5 min – 90/10 – 0.1 min – 5/95 – 3
Flow Rate : 0.3 mL/min
Col. Temp. : 80 °C
Detection : UV 210 nm
Injection Vol. : 5 µL
Sample : Standard

Analyte

1. IgG 0.1 mg/mL

Analysis of Abraxane



Conditions

Column : MonoSelect nPEC (50 x 3.0 mm I.D.)
Column Cat. No.: 5020-10816
Eluent : A) CH₃OH
B) 50 mM Na₂SO₄ in 20 mM Acetate buffer (pH 4.6)
A/B = 30/70, v/v
Flow Rate : 1.0 mL/min
Col. Temp. : Room temperature (approx. 25 °C)
Detection : UV 260 nm
Injection Vol. : 10 µL

Analyte

1. Abraxane
(Nanoparticle albumin-bound
paclitaxel)
2. Paclitaxel

* Unknown peak