

New!



HPLC, LC/MS Columns

# InertSustain® Cyano

The Most Reliable and Reproducible Cyano Column

## Physical Properties

- Silica : Newly Developed ES Silica Gel
- Particle Size : 3  $\mu\text{m}$ , 5  $\mu\text{m}$
- Surface Area : 350  $\text{m}^2/\text{g}$
- Pore Size : 100  $\text{\AA}$  (10 nm)
- Pore Volume : 0.85 mL/g
- Bonded Phase : Cyanopropyl Groups
- End-capping : Yes
- Carbon Loading : 8 %
- pH Range : 2~7.5
- USP Code : L10



# InertSustain® Cyano

## The Most Reliable and Reproducible Cyano Column

In general, the stability and reproducibility of the Cyano phase available in the market are poor. Many batch-to-batch or lot-to-lot reproducibility issues are occurring at many laboratories.

The InertSustain Cyano columns were developed to resolve these problems and are designed using the most modern LC column technology available providing them to be extremely inert, stable and reproducible.


The InertSustain Cyano columns are highly recommended for all pharmacopeia methods requiring a Cyano phase to be used. (Ex: USP L10)

## Benefits

- Endlessly reproducible from column-to-column and batch-to-batch
- Highly recommended for all pharmacopeia methods requiring a Cyano phase to be used (Ex: USP L10)
- Highly inert packing material results in less tailing of peaks for virtually any type of analytes
- A new selectivity option for method development due to the multiple retention mechanisms
- Originally shipped in reversed-phase solvents and is ready to use for reversed-phase methods
- Can be used for both reversed-phase separations as well as normal-phase separations

## Comparison of Batch-To-Batch Reproducibility


As proven below, InertSustain Cyano provide exceptional reproducibility from batch-to-batch even with those challenging strong basic compounds such as Dextromethorphan or Berberine.




## Highly Inert Packing Material

As shown below, InertSustain Cyano columns provide symmetric peaks for strong bases and chelating compounds, delivering highly stable chromatograms for qualitative and quantitative analysis.

**Figure 1. Comparison of Peak Shapes using Strong Basic Compounds**




**Figure 2. Comparison of Peak Shape using Strong Chelating Compound**



# InertSustain<sup>®</sup> Cyano


## Highly Stable Chemistry for Normal-Phase Separations

Although the InertSustain Cyano columns are originally shipped in reversed-phase solvents, it can also be used for normal-phase separations by properly equilibrating the column with ethanol or 2-propanol prior to the analysis.



## A New Selectivity Option


The InertSustain Cyano column uses multiple retention mechanisms which can lead to achieving the desired separation where a C18 or Phenyl column failed to separate. The Cyano column provides different separation pattern and decreased retention for hydrophobic compounds when comparing to a C18 or Phenyl column.



Applications

**Analysis of Nortriptyline Hydrochloride Capsules  
[USP Method]**

N : 7988  
Tailing factor : 1.10



USP Column : 5 µm, 150 x 4.6 mm I.D. (L10)

System suitability requirements:

Efficiency (N) : > 500

Tailing factor : < 3.0


Sample Conc. : 0.38 mg/mL (in Methanol)

Mobile Phase : ACN : CH<sub>3</sub>OH : 12 mM Potassium phosphate (pH 6.7)  
= 40 : 43 : 17

Flow Rate : 2.5 mL


Detection : UV 239 nm

Injection : 5 µL



**Analysis of Sertraline Hydrochloride  
[USP Method]**

Tailing factor : 1.24



USP Column : 5 µm, 150 x 4.6 mm I.D. (L10)

System suitability requirements:

Tailing factor : < 2.0

Sample Conc. : 0.050 mg/mL (in Mobile Phase)


Mobile Phase : CH<sub>3</sub>OH : 0.1% (v/v) Phosphoric acid  
= 1 : 1

Flow Rate : 1.5 mL

Detection : UV 210 nm

Colum Temp. : 30 °C

Injection : 10 µL




# InertSustain<sup>®</sup> Cyano

## Applications


### Analysis of Tetracaine Hydrochloride Ophthalmic Solution [USP Method]

N : 9679  
Tailing factor : 1.16



USP Column : 5  $\mu$ m, 150 x 4.6 mm I.D. (L10)  
System suitability requirements:  
Efficiency (N) : > 500  
Tailing factor : < 2.0

Sample Conc. : 0.1 mg/mL (in Water)  
Mobile Phase : ACN : 10 mM Ammonium phosphate (pH 3.0)  
= 30 : 70  
Flow Rate : 2.0 mL  
Detection : UV 280 nm  
Injection : 10  $\mu$ L



## Ordering Information

### InertSustain<sup>®</sup> Cyano Analytical Columns

| HP Series<br>Particle Size: 3 $\mu$ m<br>Max. Operating Pressure:<br>50 MPa (500 Bar) | Length / I.D. (mm) | 2.1        | 3.0        | 4.6        |
|---|--------------------|------------|------------|------------|
|   | 30                 | 5020-89459 | 5020-89465 | 5020-89471 |
| 50  | 5020-89460         | 5020-89466 | 5020-89472 |            |
| 75  | 5020-89461         | 5020-89467 | 5020-89473 |            |
| 100   | 5020-89462         | 5020-89468 | 5020-89474 |            |
| 150   | 5020-89463         | 5020-89469 | 5020-89475 |            |
| 250   | 5020-89464         | 5020-89470 | 5020-89476 |            |

\* End-fittings are 1/16" Waters-compatible.

| Particle Size: 3 $\mu$ m | Length / I.D. (mm) | 2.1        | 3.0        | 4.0        | 4.6        |
|--------------------------|--------------------|------------|------------|------------|------------|
|                          | 30                 | 5020-89374 | 5020-89381 | 5020-89388 | 5020-89395 |
| 50                       | 5020-89375         | 5020-89382 | 5020-89389 | 5020-89396 |            |
| 75                       | 5020-89376         | 5020-89383 | 5020-89390 | 5020-89397 |            |
| 100                      | 5020-89377         | 5020-89384 | 5020-89391 | 5020-89398 |            |
| 150                      | 5020-89378         | 5020-89385 | 5020-89392 | 5020-89399 |            |
| 250                      | 5020-89379         | 5020-89386 | 5020-89393 | 5020-89400 |            |

| Particle Size: 5 $\mu$ m | Length / I.D. (mm) | 2.1        | 3.0        | 4.0        | 4.6        |
|--------------------------|--------------------|------------|------------|------------|------------|
|                          | 30                 | 5020-89251 | 5020-89258 | 5020-89265 | 5020-89272 |
| 50                       | 5020-89252         | 5020-89259 | 5020-89266 | 5020-89273 |            |
| 75                       | 5020-89253         | 5020-89260 | 5020-89267 | 5020-89274 |            |
| 100                      | 5020-89254         | 5020-89261 | 5020-89268 | 5020-89275 |            |
| 150                      | 5020-89255         | 5020-89262 | 5020-89269 | 5020-89276 |            |
| 250                      | 5020-89256         | 5020-89263 | 5020-89270 | 5020-89277 |            |

\* End-fittings are 1/16" Waters-compatible.

\* Max. Operating Pressure: 20 MPa (200 Bar)

## Ordering Information

### Cartridge Guard Column E

| I.D. of the Analytical Column Applicable (mm) | Length (mm) | I.D. (mm) | Replacement Cartridge E Guard Column |                  | Cartridge E Holder / Cartridge Set       |            |
|---|-------------|-----------|--------------------------------------|------------------|--|------------|
|   |             |           | (2 EA.)                              |                  | (2 Cartridge E Guard Columns & 1 Holder) |            |
|   |             |           | Particle Size                        |                  | Particle Size                            |            |
|   |             |           | 3 µm                                 | 5 µm             | 3 µm                                     | 5 µm       |
| 1.0   | 10          | 1.0       | 5020-89449                           | 5020-89355       | 5020-89450                               | 5020-89356 |
| 1.5,2.1                                       |             | 1.5       | 5020-89451                           | 5020-89357       | 5020-89452                               | 5020-89358 |
| 2.1,3.0                                       |             | 3.0       | 5020-89447                           | 5020-89353       | 5020-89448                               | 5020-89354 |
| 4.0,4.6                                       |             | 4.0       | 5020-89445                           | 5020-89351       | 5020-89446                               | 5020-89352 |
| 2.1,3.0                                       | 20          | 3.0       | 5020-89455                           | 5020-89361       | 5020-89456                               | 5020-89362 |
| 4.0,4.6                                       |             | 4.0       | 5020-89453                           | 5020-89359       | 5020-89454                               | 5020-89360 |
| Holder for Cartridge Guard Column E           |             |           |                                      | For 10 mm Length |  | 5020-08500 |
|   |             |           |                                      | For 20 mm Length |  | 5020-08550 |

\* End-fittings are 1/16" Waters-compatible.

\* Max. Operating Pressure: 20 MPa (200 Bar)



Cartridge Guard Column E

## Worldwide Ordering Information

### GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan

Phone: +81-3-5323-6620

Fax: +81-3-5323-6621

Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)

Web: [www.glsciences.com](http://www.glsciences.com)

### GL Sciences B.V.

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands

Phone: +31 (0)40 254 95 31

Email: [info@glsciences.eu](mailto:info@glsciences.eu)

Web: [www.glsciences.eu](http://www.glsciences.eu)

### GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503

Phone: 310-265-4424

Fax: 310-265-4425


Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)

Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

### International Distributors

Visit our Website at

[www.glsciences.com/distributors](http://www.glsciences.com/distributors)



#### AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17• D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: [info@mz-at.de](mailto:info@mz-at.de), [www.mz-at.de](http://www.mz-at.de)

The GL Sciences name, the GL Sciences logo and the following registered trademark or trademark are the property of GL Sciences Inc.

InertSustain  
Inertsil

All other trademarks or service marks are the property of their respective owners.

The specification and the column type are subject to change without notice due to continual improvements.