

Analytical HPLC Column Introduction

Grace® Key Column Families

Single-Source Solution for Discovery to Recovery Applications

HPLC is commonly used in a wide range of applications, including drug discovery and purification for the pharmaceutical and biotechnology industries, environmental analysis, forensics, petrochemical analysis, food, cosmetics, and vitamins. The combined premier product lines of Alltech®, Davisil®, Flexit™, Grom™, Jones Chromatography™, Modcol®, and Vydac® create a single-source solution for HPLC columns and accessories from discovery to recovery.

Our column families include reversed-phase, normal-phase, HILIC, ion-exchange, ion-exclusion, size-exclusion, and affinity stationary phases for small- and large-molecule separations, and our column formats maximize their performance. Our product range includes standard and custom columns for analytical separations, preparative phases for scale up, bulk packings that customers can pack in their own columns, and accessories to maintain separation ruggedness and quality.

To help select the appropriate column for your application, we describe key column families and highlight unique phases within these families. Whether the most important factor is analysis speed, column bleed, pH stability, resolution, adjustable selectivity, or analyte molecular weight, Grace offers a column to suit your application.




more info

For information about Grace® analytical hardware formats, see page 31.

more info

For information about prep columns, see pages 148–171.




MZ-Analysentechnik GmbH
 Barcelona-Allee 17 • D-55129 Mainz
 Tel +49 6131 880 96-0
 Fax +49 6131 880 96-20
 e-mail: info@mz-at.de
 www.mz-at.de

AUTHORIZED DISTRIBUTOR

www.discoverysciences.com

Column Families Suitable for Small Molecules (<2000 Molecular Weight)

VisionHT™ Ultra High-Pressure Columns



12,000psig pressure-rated columns, with sub 2µm media for high-efficiency, high-speed separations

VisionHT™ columns are optimized for high throughput and ultra high-pressure LC applications. Mechanically stable 1.5µm media and ultra-low volume hardware delivers new separation benefits with excellent stability. A variety of phases available.

Differentiated Phases: C18, C10-B, C18-P, C18-HL, HILIC, Silica

Specifications: Spherical silica, monomerically bonded, proprietary endcapping, 100, 120Å pore size

Formats: Ultra High-Pressure

Alltech® Alltima™ HP Columns



Premium quality, exceptionally stable columns, without phase bleed

Alltech® Alltima™ HP phases are made from high purity silica to eliminate peak tailing. Proprietary bonding eliminates the problem of column bleed for MS and ELSD detection. Full range of phases with pH stability from 1 to 10.

Differentiated Phases: C18 for classic reversed-phase separations, EPS for extended polar selectivity, C18 HiLoad for extra reversed-phase retention, C18 Amide with low bleed, and HILIC.

Specifications: High purity spherical silica, monomerically bonded, endcapped, 100, 120, 190Å pore size

Formats: Microbore, Expedite™, Rocket™, Solvent-Reducer, Analytical, Prep

Alltech® Prevail™ Columns



Reversed-phase columns for use with 100% organic to 100% aqueous mobile phases

Prevail™ phases are designed for wettability with aqueous and organic mobile phases, and are especially useful for broad gradients. Retain highly polar analytes with aqueous mobile phases as well as hydrophobic analytes in organic mobile phases.

Differentiated Phases: Carbohydrate ES, Organic Acid


Specifications: Spherical silica, monomerically bonded, endcapped, 110Å pore size

Formats: Microbore, Expedite™, Rocket™, Solvent-Reducer, Analytical, Prep

Analytical HPLC Column Introduction

Column Families Suitable for Small Molecules
(<2000 Molecular Weight) (continued)

Alltech® Platinum™ Columns



Alternate retention and selectivity for polar and nonpolar analytes available in sub 2µm particle sizes

Alltech® Platinum™ phases have a controlled polar silica surface which provides enhanced polar selectivity (EPS) in addition to reversed-phase retention. C18 and EPS C18 columns are a complimentary duo to separate polar and nonpolar analytes.

Differentiated Phases: C18, EPS C18

Specifications: Spherical silica, monomerically bonded, proprietary endcapping, 100Å pore size

Formats: Microbore, Expedite™, Rocket™, Solvent-Reducer, Analytical, Prep

GraceSmart™ Columns



Classic reversed-phase and premium performance at an exceptional value

GraceSmart™ HPLC phases use high purity silica. Which translates into symmetrical peaks for acids and bases, and predictable reversed-phase selectivity. Whether routine analysis or new method development, use individually tested GraceSmart™ columns to get premium performance at exceptional value.

Differentiated Phases: C18

Specifications: Spherical silica, monomerically bonded, endcapped, 120Å pore size

Format: Analytical

GraceAlpha™ Columns



A New Silica Generation

New silica from Grace with high porosity surface and dense core ideal for scale-up applications. While high porosity surface increases mass transfer results in increased column efficiencies and loading capacity, dense core and highly spherical shape yields mechanically robust particle.

Differentiated Phases: C18, C8, and Silica

Specifications: High purity spherical silica with high porosity surface and dense core

Formats: Analytical, Prep

Column Families Suitable for Large Molecules
(>2000 Molecular Weight)

Vydac® MS Columns



Next generation Vydac® media from the leader in peptide and protein separations

Vydac® MS columns provide unique selectivity and exceptional protein recovery. They require less TFA in the mobile phase for good peak shape, which increases microbore sensitivity of peptides and proteins by reducing “quenching”. Vydac® MS columns are also applicable to hydrophobic proteins and peptide mapping.

Differentiated Phases: C18, C8, C4

Specifications: Spherical silica, polymerically and monomerically bonded, endcapped, 300Å pore size

Formats: Capillary, Microbore, Expedite™, Rocket™, Solvent-Reducer, Analytical, Prep

Vydac® ProZap™ Columns



High-speed, high-efficiency columns for fast protein analysis

Short 10 and 20mm ProZap™ columns with 1.5µm particles provide high-speed protein and life science separations. Sharp, efficient peaks maximize method sensitivity.

Differentiated Phases: C18

Specifications: Sub 2µm spherical silica, proprietary bonding and endcapping, 500Å pore size

Format: Expedite™

Vydac® Everest® Columns



For peptide separations, peptide mapping and proteomics applications

High-capacity Everest® C18 columns are recommended as the first choice in peptide separations. Everest® offer high capacity to resolve complex samples such as protein digests.

Differentiated Phases: C18

Specifications: Spherical silica, monomerically bonded, endcapped, 300Å pore size

Formats: Capillary, Microbore, Expedite™, Analytical

HPLC Column Selection

A Comparison of Reversed-Phase Columns

Based on the widely accepted work and data of Drs. Lloyd Snyder and John Dolan^{1,2} Grace has developed this column selection tool for choosing reversed phase HPLC columns based on peak capacities and column selectivity of polar and nonpolar compounds. Typically, chromatographers choose HPLC columns by comparing physical characteristics, such as surface area and carbon load. Often, this does not provide enough information about selectivity or capacity for adequate column selection. This chart provides a reliable means of choosing HPLC columns based on acidic, basic, and hydrophobic character.

The Snyder/Dolan column test procedure has been described in a series of publications. Based on retention data for a series of standard mixtures and the same separation conditions (50% acetonitrile/buffer; pH 2.8 and 7.0; 35°C), reversed-phase columns are characterized by five column-selectivity parameters: hydrophobicity (H), steric interaction (S*), hydrogen-bond acidity (A), basicity (B), and relative silanol ionization or cation-exchange capacity (C). Here we have chosen to graphically highlight data for H, A, and C, with C results at pH 7.0. Hydrophobicity (H) is often the primary analyte interaction with reversed phase columns and indicates overall capacity. Secondary interactions are often polar interactions with basic analytes. The degree of unprotonated base interaction (A) and protonated base interaction (C) with the packing material is measured and represented here.

Directions for Using the Column Chart

The chart lists the columns in descending order of hydrophobic capacity (H). To find similar HPLC columns to test as back-up columns, follow these steps.

- 1) Find the column you are currently using and note neighboring columns which have similar (H) capacity factors.
- 2) Compare the values for interaction of polar compounds (A and C).

If there is more than one choice for a back-up column, then compare your actual sample to the test probes. If your sample is nonpolar, then place more emphasis on hydrophobic values. If your sample is basic (polar), then pay special attention to A and C and determine if your sample will be protonated (A) or unprotonated (C) and place greater emphasis on one of these values.

Key to Chart

- Hydrophobic Indicator
- ▲ Hydrogen bonding Indicator—pH 2.8 (Protonated under acidic conditions)
- Cation Exchange Indicator—pH 7.0 (Unprotonated under neutral conditions)

| Manufacturer | Column | Selectivity Parameters | |
|-----------------|--------------------------------|------------------------|-------|
| ZirChrom | ZirChrom®-PBD C18 | H, A, C | 1.284 |
| YMC | J'Sphere® H80 C18 | H, A, C | 1.132 |
| Restek | Allure® C18 | H, A, C | 1.116 |
| Phenomenex | Ultrasorb® ODS (30) | H, A, C | 1.114 |
| YMC | YMC® Pack Pro C18 RS | H, A, C | 1.114 |
| Grace (Alltech) | Adsorbosphere™ UHS C18 | H, A, C | 1.103 |
| Thermo/Hypersil | Hypersil® BetamaxNeutral C18 | H, A, C | 1.098 |
| Agilent | Zorbax Extend C18 | H, A, C | 1.098 |
| Agilent | Zorbax C18 | H, A, C | 1.089 |
| Beckman | Ultrasphere® ODS | H, A, C | 1.085 |
| Grace (Alltech) | Alltima™ HP C18 High Load | H, A, C | 1.080 |
| Agilent | Zorbax Rx-18 | H, A, C | 1.077 |
| Agilent | Zorbax Eclipse XDB-C18 | H, A, C | 1.077 |
| Supelco | Ascentis® C18 | H, A, C | 1.077 |
| Macherey Nagel | Nucleodur® C18 Gravity | H, A, C | 1.056 |
| Grace (Grom) | Grom™ Sapphire 110 C18 | H, A, C | 1.055 |
| Restek | Restek® Ultra C18 | H, A, C | 1.055 |
| Varian | OmniSpher™ 5 C18 | H, A, C | 1.055 |
| Grace (Vydac) | Denali® 120 C18 | H, A, C | 1.052 |
| Waters | Symmetry® C18 | H, A, C | 1.052 |
| Akzo Nobel | Kromasil® 100-5C18 | H, A, C | 1.051 |
| Waters | Nova-Pak® C18 | H, A, C | 1.049 |
| Thermo/Hypersil | Hypersil® 100 C18 | H, A, C | 1.048 |
| MacMod/ACT | ACE® 5 C18-HL | H, A, C | 1.045 |
| ZirChrom | ZirChrom®-EZ C18 | H, A, C | 1.040 |
| Grace (Grom) | Grom™ Sil 120 ODS-5 ST | H, A, C | 1.035 |
| Dionex | Acclaim® 120 C18 | H, A, C | 1.032 |
| Waters | Sunfire™ C18 | H, A, C | 1.031 |
| Agilent | Zorbax Eclipse Plus C18 | H, A, C | 1.030 |
| Merck | Superspher® 100 RP-18e | H, A, C | 1.030 |
| Shiseido | CAPCELL™ C18 AG120 | H, A, C | 1.030 |
| Grace (Grom) | Grom™ Sil 120 ODS-3 CP | H, A, C | 1.029 |
| Waters | Delta-Pak™ C18 100A | H, A, C | 1.028 |
| Macherey Nagel | Nucleodur® Isis | H, A, C | 1.023 |
| Phenomenex | Prodigy™ ODS (3) | H, A, C | 1.023 |
| Phenomenex | Synergi™ Hydro-RP C18 | H, A, C | 1.022 |
| Phenomenex | Luna™ C18 | H, A, C | 1.018 |
| Supelco | Supelcosil™ LC-18 | H, A, C | 1.018 |
| YMC | YMC® Pro C18 | H, A, C | 1.015 |
| Phenomenex | Onyx™ Monolithic C18 | H, A, C | 1.012 |
| Bischoff | ProntoSIL™ SphenBOND 80-5-ODS2 | H, A, C | 1.010 |
| Grace (Jones) | Apex™ II C18 | H, A, C | 1.008 |
| Shiseido | CAPCELL™ C18 UG120 | H, A, C | 1.007 |
| GL Sciences | Inertsil® ODS-2 | H, A, C | 1.007 |

| Manufacturer | Column | Selectivity Parameters | |
|-----------------|-----------------------------|------------------------|-------|
| Merck | LiChrospher® 100 RP-18 | H, A, C | 1.006 |
| Bischoff | ProntoSIL™ 120-5-C18 H | H, A, C | 1.005 |
| Shiseido | CAPCELL™ C18 M G | H, A, C | 1.005 |
| Grace (Jones) | Genesis® 120 C18 | H, A, C | 1.005 |
| Bischoff | EU Reference Column C18 | H, A, C | 1.004 |
| Grace (Alltech) | Allsphere™ ODS2 | H, A, C | 1.004 |
| Merck | Purospher® STAR RP18e | H, A, C | 1.003 |
| Merck | Chromolith® RP18e | H, A, C | 1.003 |
| Phenomenex | Luna® C18(2) | H, A, C | 1.002 |
| Varian | Pursuit® C18 | H, A, C | 1.001 |
| MacMod/ACT | ACE® 5 C18 | H, A, C | 1.000 |
| Tosoh | TSKge® Super-ODS | H, A, C | 0.998 |
| Agilent | Zorbax StableBond 80A C18 | H, A, C | 0.996 |
| Phenomenex | Prodigy™ ODS(2) | H, A, C | 0.995 |
| Thermo/Hypersil | Hypersil® BDS C18 | H, A, C | 0.993 |
| Grace (Alltech) | Alltima™ C18 | H, A, C | 0.993 |
| Grace (Vydac) | Vydac® Everest® C18 | H, A, C | 0.993 |
| Thermo/Hypersil | Hypersil® Beta Basic-18 | H, A, C | 0.993 |
| GL Sciences | Inertsil® ODS-3 | H, A, C | 0.990 |
| Grace (Alltech) | Adsorbosphere™ C18 | H, A, C | 0.989 |
| Phenomenex | Synergi™ Max-RP C18 | H, A, C | 0.989 |
| Shiseido | CAPCELL™ C18 SG120 | H, A, C | 0.987 |
| Grace (Jones) | Apex™ I C18 | H, A, C | 0.985 |
| Thermo/Hypersil | Hypersil® ODS-2 | H, A, C | 0.985 |
| Grace (Alltech) | Alltima™ HP C18 | H, A, C | 0.985 |
| Waters | Xterra® MS C18 | H, A, C | 0.984 |
| Waters | Symmetry® 300 C18 | H, A, C | 0.984 |
| Supelco | Discovery C18 | H, A, C | 0.984 |
| Supelco | Supelcosil™ LC-18-DB | H, A, C | 0.979 |
| Waters | Spherisorb® S5 ODSB | H, A, C | 0.975 |
| Thermo/Hypersil | Hypersil® Bio Basic-18 | H, A, C | 0.974 |
| Thermo/Hypersil | Hypersil® ODS | H, A, C | 0.974 |
| Bischoff | ProntoSIL™ 120-5-C18-AQ | H, A, C | 0.974 |
| Grace (Jones) | Genesis® 300 C18 C18 | H, A, C | 0.974 |
| Bischoff | Prontosil™ 200-5-C18 AQ | H, A, C | 0.974 |
| Agilent | Zorbax C8 | H, A, C | 0.974 |
| Tosoh | TSK ge® ODS-80Ts | H, A, C | 0.971 |
| Waters | Resolve C18 | H, A, C | 0.968 |
| Phenomenex | Gemin® C18 110A | H, A, C | 0.967 |
| Grace (Alltech) | Econosil™ C18 | H, A, C | 0.966 |
| Phenomenex | Aqua® C18 | H, A, C | 0.966 |
| YMC | YMC® ODS-AQ C18 | H, A, C | 0.965 |
| Waters | Spherisorb® ODS-2 | H, A, C | 0.962 |
| Macherey Nagel | Nucleosil® 100-5-C18 HD C18 | H, A, C | 0.961 |

| | | | |
|-----------------|----------------------------------|-----|-------|
| Grace (Jones) | Genesis® 120 AQ C18 | ⊕ ⊕ | 0.960 |
| Macherey Nagel | Nucleodur® Pyramid | ⊕ ⊕ | 0.958 |
| ThermoHypersil | Hypersil® Elite C18 | ⊕ ⊕ | 0.958 |
| Dionex | Acclaim® 300 C18 | ⊕ ⊕ | 0.957 |
| Bischoff | ProntoSIL™ 300-5-C18 H | ⊕ ⊕ | 0.956 |
| Waters | Delta-Pak™ C18 300A | ⊕ ⊕ | 0.955 |
| Bischoff | ProntoSIL™ HyperSORB 120 ODS | ⊕ ⊕ | 0.951 |
| ThermoHypersil | Hypersil® PAH C18 | ⊕ ⊕ | 0.949 |
| Bischoff | ProntoSIL™ 120-5-C18 Agplus | ⊕ ⊕ | 0.947 |
| Phenomenex | Jupiter® 300 C18 | ⊕ ⊕ | 0.945 |
| Varian | Polaris® C18-Ether | ⊕ ⊕ | 0.943 |
| Waters | Atlantis™ T3 C18 | ⊕ ⊕ | 0.941 |
| Tosoh | TSKgel® 80Ts QA | ⊕ ⊕ | 0.940 |
| Grace (Alltech) | Alltima™ C18-WP | ⊕ ⊕ | 0.938 |
| Waters | YMC® Hydrosphere C18 | ⊕ ⊕ | 0.937 |
| Grace (Alltech) | Brava™ BDS C18 | ⊕ ⊕ | 0.938 |
| Bischoff | ProntoSIL™ 60-5 C8 SH | ⊕ ⊕ | 0.929 |
| Varian | Polaris® C18-A | ⊕ ⊕ | 0.928 |
| YMC | J'Sphere® M80 C18 | ⊕ ⊕ | 0.926 |
| Agilent | Zorbax Eclipse XDB-C8 | ⊕ ⊕ | 0.919 |
| Waters | Atlantis® oC18 b | ⊕ ⊕ | 0.917 |
| Supelco | Ascentis® Express C8 | ⊕ ⊕ | 0.915 |
| Thermo | Hypersil® GOLD aQ | ⊕ ⊕ | 0.915 |
| Phenomenex | Selectosil™ C18 | ⊕ ⊕ | 0.911 |
| Merck | LiChrosorb® RP-18 | ⊕ ⊕ | 0.909 |
| Grace (Vydac) | Vydac® 218TP C18 | ⊕ ⊕ | 0.909 |
| Waters | Acquity UPLC® BEH Shield RP18 EP | ⊕ ⊕ | 0.907 |
| Macherey Nagel | Nucleosil® C18 | ⊕ ⊕ | 0.906 |
| Agilent | Zorbax StableBond 300A C18 | ⊕ ⊕ | 0.905 |
| Grace (Alltech) | Prospere® C18 300 | ⊕ ⊕ | 0.903 |
| Grace (Vydac) | Vydac® 201TP C18 | ⊕ ⊕ | 0.901 |
| Supelco | Ascentis® C-8 | ⊕ ⊕ | 0.899 |
| Waters | Nova-Pak® C8 | ⊕ ⊕ | 0.899 |
| Waters | Symmetry® C8 | ⊕ ⊕ | 0.893 |
| YMC | YMC® Pro C8 | ⊕ ⊕ | 0.890 |
| Agilent | Zorbax Eclipse Plus C8 | ⊕ ⊕ | 0.889 |
| Phenomenex | Luna™ C8(2) | ⊕ ⊕ | 0.889 |
| Grace (Alltech) | Prevail™ C18 | ⊕ ⊕ | 0.888 |
| Grace (Alltech) | Prospere™ 100 C18 | ⊕ ⊕ | 0.883 |
| Grace (Alltech) | Alltima™ AQ EP | ⊕ ⊕ | 0.882 |
| ThermoHypersil | Hypersil® GOLD C18 | ⊕ ⊕ | 0.881 |
| Phenomenex | Synergi™ Fusion-RP EP | ⊕ ⊕ | 0.879 |
| Phenomenex | Luna™ C8 | ⊕ ⊕ | 0.875 |
| Grace (Grom) | Grom™ Sil 120 Octyl-6 MB C8 | ⊕ ⊕ | 0.872 |
| Grace (Jones) | Apex™ I C8 | ⊕ ⊕ | 0.869 |
| Shiseido | CAPCELL™ C18 A Q | ⊕ ⊕ | 0.867 |
| Macherey Nagel | Nucleosil® 100-5-C8 HD | ⊕ ⊕ | 0.865 |
| Akzo Nobel | Kromasil® 100-5C8 | ⊕ ⊕ | 0.864 |
| Grace (Jones) | Genesis® 120 EC C8 | ⊕ ⊕ | 0.863 |
| Grace (Alltech) | Prevail™ Amide EP | ⊕ ⊕ | 0.862 |
| Macherey Nagel | Nucleosil® ODS | ⊕ ⊕ | 0.860 |
| Restek | Ultra AQ C18 | ⊕ ⊕ | 0.857 |
| Waters | Sunfire™ C8 | ⊕ ⊕ | 0.856 |
| Waters | Acquity UPLC® BEH C8 | ⊕ ⊕ | 0.855 |
| Shiseido | CAPCELL™ PAK C8 UG120 | ⊕ ⊕ | 0.854 |
| Waters | Symmetry® Shield C18 | ⊕ ⊕ | 0.850 |
| Grace (Alltech) | Alphabond™ C18 | ⊕ ⊕ | 0.845 |
| Supelco | Ascentis® RP-Amide | ⊕ ⊕ | 0.843 |
| Merck | Purospher® RP-18 | ⊕ ⊕ | 0.841 |
| Supelco | Discovery BIO Wide pore C8 | ⊕ ⊕ | 0.839 |
| Grace (Grom) | Grom™ Sapphire 110 C8 | ⊕ ⊕ | 0.835 |
| ThermoHypersil | Hypersil® Beta Basic-8 | ⊕ ⊕ | 0.834 |
| Grace (Alltech) | Alltima™ HP C8 | ⊕ ⊕ | 0.834 |
| ThermoHypersil | Hypurity® C8 | ⊕ ⊕ | 0.833 |
| Supelco | Discovery C8 | ⊕ ⊕ | 0.832 |
| GL Sciences | Inertsil® C8-3 C8 | ⊕ ⊕ | 0.830 |
| MacMod/ACT | ACE® 5 C8 | ⊕ ⊕ | 0.830 |
| Grace (Jones) | Genesis® 120 C8 | ⊕ ⊕ | 0.829 |
| Thermo | Hypersil® GOLD C8 | ⊕ ⊕ | 0.825 |
| Phenomenex | Oryx™ Monolithic C8 | ⊕ ⊕ | 0.824 |
| Tosoh | TSKgel® Super-Octyl | ⊕ ⊕ | 0.824 |
| Grace (Alltech) | Prevail™ Select C18 | ⊕ ⊕ | 0.822 |
| ThermoHypersil | Hypersil® Bio Basic-8 | ⊕ ⊕ | 0.821 |
| YMC | YMC® Basic C18 | ⊕ ⊕ | 0.821 |
| Grace (Alltech) | Econosphere™ C18 | ⊕ ⊕ | 0.818 |
| Tosoh | TSKgel® Octyl-80Ts | ⊕ ⊕ | 0.814 |
| Whatman | Partisil™ ODS(3) | ⊕ ⊕ | 0.810 |

| | | | |
|-----------------|----------------------------------|-----|-------|
| Macherey Nagel | Nucleodur® Sphinx RP | ⊕ ⊕ | 0.805 |
| ThermoHypersil | Aquasil™ C18 | ⊕ ⊕ | 0.805 |
| MacMod/ACT | ACE® AQ EP | ⊕ ⊕ | 0.804 |
| Waters | Xterra® MS C8 | ⊕ ⊕ | 0.803 |
| Phenomenex | Luna™ C5 | ⊕ ⊕ | 0.800 |
| Waters | MicroBondapak C18 | ⊕ ⊕ | 0.798 |
| Agilent | Zorbax StableBond 80A C8 | ⊕ ⊕ | 0.795 |
| Agilent | Zorbax Rx-C8 | ⊕ ⊕ | 0.792 |
| Grace (Alltech) | Platinum™ C18 | ⊕ ⊕ | 0.786 |
| Grace | VisionHT™ C18 | ⊕ ⊕ | 0.786 |
| Phenomenex | Luna™ Phenyl-Hexyl | ⊕ ⊕ | 0.782 |
| Grace (Alltech) | Alltima™ C18-LL | ⊕ ⊕ | 0.780 |
| Bischoff | ProntoSIL™ 120-5-C18 ace-EPS | ⊕ ⊕ | 0.772 |
| Grace (Vydac) | Vydac® 218MS C18 | ⊕ ⊕ | 0.770 |
| Waters | Acquity UPLC® BEH phenyl | ⊕ ⊕ | 0.764 |
| Waters | Spherisorb® C8 | ⊕ ⊕ | 0.763 |
| YMC | J'Sphere® L80 C18 | ⊕ ⊕ | 0.762 |
| Bischoff | ProntoSIL™ 300-55-C18 ace-EPS | ⊕ ⊕ | 0.762 |
| Dionex | Acclaim® Organic Acid C18 | ⊕ ⊕ | 0.761 |
| Waters | Xterra® C18 RP | ⊕ ⊕ | 0.757 |
| Grace (Alltech) | Alltima™ C8 | ⊕ ⊕ | 0.756 |
| Whatman | Partisil™ C8 | ⊕ ⊕ | 0.749 |
| Merck | LiChrospher® 60 RP-Select B C18 | ⊕ ⊕ | 0.747 |
| Bischoff | ProntoSIL™ 120-5 C8 SH | ⊕ ⊕ | 0.739 |
| Grace (Alltech) | Allsphere™ ODS1 | ⊕ ⊕ | 0.733 |
| Waters | Symmetry® Shield C8 | ⊕ ⊕ | 0.730 |
| ThermoHypersil | Hypurity® C4 | ⊕ ⊕ | 0.713 |
| MacMod/ACT | ACE® 5 C4-300 | ⊕ ⊕ | 0.710 |
| Varian | Polaris® C8-Ether | ⊕ ⊕ | 0.705 |
| Bischoff | ProntoSIL™ 60-5-Phenyl | ⊕ ⊕ | 0.705 |
| Waters | Nova-Pak® Phenyl | ⊕ ⊕ | 0.704 |
| Macherey Nagel | Nucleosil® 100-5-C18 Nautilus | ⊕ ⊕ | 0.702 |
| Agilent | Zorbax StableBond 300A C8 | ⊕ ⊕ | 0.701 |
| Bischoff | ProntoSIL™ SpheriBOND 80-5-ODS1 | ⊕ ⊕ | 0.700 |
| ThermoHypersil | Fluophase® RP F | ⊕ ⊕ | 0.698 |
| Phenomenex | Jupiter® 300 C4 | ⊕ ⊕ | 0.698 |
| Grace (Alltech) | Prospere® 300 C4 | ⊕ ⊕ | 0.689 |
| Bischoff | Prontosil™ 60-5-C4 | ⊕ ⊕ | 0.686 |
| Waters | Xterra® Phenyl | ⊕ ⊕ | 0.683 |
| Waters | Spherisorb® ODS-1 | ⊕ ⊕ | 0.682 |
| ThermoHypersil | Hypersil® Prism C18 RPN | ⊕ ⊕ | 0.678 |
| ThermoHypersil | Fluophase® PFP F | ⊕ ⊕ | 0.675 |
| Agilent | Zorbax XDB-Phenyl | ⊕ ⊕ | 0.665 |
| Waters | Symmetry® 300 C4 | ⊕ ⊕ | 0.659 |
| Waters | Xterra® C8 RP | ⊕ ⊕ | 0.657 |
| Grace (Alltech) | Alltima™ HP C18 EPS | ⊕ ⊕ | 0.655 |
| Supelco | Discovery BIO Wide pore C5 | ⊕ ⊕ | 0.654 |
| Agilent | Zorbax Bonus RP EP | ⊕ ⊕ | 0.654 |
| Phenomenex | Synergi™ Polar-RP C18 | ⊕ ⊕ | 0.654 |
| MacMod/ACT | ACE® Phenyl | ⊕ ⊕ | 0.647 |
| Grace (Jones) | Genesis® 120 C4 EC | ⊕ ⊕ | 0.646 |
| ThermoHypersil | Hypersil® Prism C18 RP | ⊕ ⊕ | 0.645 |
| ThermoHypersil | BetaMax® Acid EP | ⊕ ⊕ | 0.635 |
| Supelco | Discovery HS F5 F | ⊕ ⊕ | 0.631 |
| Agilent | Zorbax SB-Phenyl | ⊕ ⊕ | 0.623 |
| Grace (Alltech) | Platinum™ EPS C18 | ⊕ ⊕ | 0.619 |
| Grace | VisionHT™ C18-P | ⊕ ⊕ | 0.619 |
| Grace (Alltech) | Prevail™ C8 | ⊕ ⊕ | 0.617 |
| Grace (Jones) | Genesis® 300 C4 C4 | ⊕ ⊕ | 0.615 |
| Grace (Jones) | Genesis® Phenyl | ⊕ ⊕ | 0.609 |
| Agilent | Zorbax StableBond 80A C3 | ⊕ ⊕ | 0.601 |
| MacMod/ACT | ACE® Phenyl-300 | ⊕ ⊕ | 0.599 |
| Agilent | Zorbax SB-AQ EP | ⊕ ⊕ | 0.593 |
| ZirChrom | ZirChrom®-PS EP | ⊕ ⊕ | 0.589 |
| Waters | MicroBondapak Phenyl | ⊕ ⊕ | 0.585 |
| Grace (Alltech) | Platinum™ C8 | ⊕ ⊕ | 0.584 |
| Grace (Alltech) | Platinum™ EPS C8 300 | ⊕ ⊕ | 0.584 |
| ThermoHypersil | BetaBasic® Phenyl | ⊕ ⊕ | 0.582 |
| Macherey Nagel | Nucleosil® C8 | ⊕ ⊕ | 0.575 |
| Macherey Nagel | EC Nucleosil® 100-5 Protect 1 EP | ⊕ ⊕ | 0.544 |
| Bischoff | Prontosil™ 120-5-C8 ace-EPS | ⊕ ⊕ | 0.532 |
| Phenomenex | Prodigy™ Phenyl-3 | ⊕ ⊕ | 0.529 |
| Agilent | Zorbax StableBond 300A C3 | ⊕ ⊕ | 0.526 |
| Grace (Alltech) | Alltima™ HP C18 Amide | ⊕ ⊕ | 0.497 |
| ThermoHypersil | BetaMax® Base EP | ⊕ ⊕ | 0.470 |
| Grace (Alltech) | Platinum™ EPS C8 | ⊕ ⊕ | 0.420 |
| ThermoHypersil | Hypurity® Advance | ⊕ ⊕ | 0.412 |

References:

1. "The "Hydrophobic-subtraction" Model of Reversed-phase Column Selectivity", L.R. Snyder, J.W. Dolan and P.W. Carr, *J. Chromatogr. A*, 1060 (2004) 77-116.
2. "A New Look at the Selectivity of Reversed-phase HPLC Columns", L.R. Snyder, J.W. Dolan and P.W. Carr, *Anal. Chem.*, 79 (2007) 3255-3262.

Grace® HPLC Packing Material Specifications

hplc columns | introduction

| Columns for Small Molecules | | | | | | | | | | | |
|--------------------------------------------|--------------------------------------|---------------|----------------|---------------|------------|----------------------|----------------------|------------|-------------|------------|----|
| Brand | Phase | Base Material | Particle Shape | Particle Size | Pore Size | Surface Area | Carbon Load | Phase Type | End-capped? | USP L-code | |
| Adsorbosil® <i>Alltech</i> | C18 | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | 15% | Polymeric | Yes | L1 | |
| | C8 | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | 10% | Polymeric | Yes | L7 | |
| | C2 | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | No | Polymeric | No | L16 | |
| | CN | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | — | Polymeric | Yes | L10 | |
| | NH ₂ | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | — | Polymeric | No | L8 | |
| Adsorbosphere™ <i>Alltech</i> | Silica | Silica | Irregular | 5, 10µm | 60Å | 450m ² /g | — | Polymeric | No | L3 | |
| | C18 | Silica | Spherical | 3, 5, 10µm | 80Å | 200m ² /g | 12% | Monomeric | Yes | L1 | |
| | C18 HS | Silica | Spherical | 3, 5µm | 60Å | 350m ² /g | 20% | Monomeric | Yes | L1 | |
| | C18 UHS | Silica | Spherical | 5, 10µm | 60Å | 500m ² /g | 30% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5, 10µm | 80Å | 200m ² /g | 8% | Monomeric | Yes | L7 | |
| | Phenyl | Silica | Spherical | 5µm | 80Å | 200m ² /g | 5% | Monomeric | Yes | L11 | |
| | Cyano | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | Monomeric | Yes | L10 | |
| | Cyano-AQ | Silica | Spherical | 5µm | 120Å | 170m ² /g | — | Polymeric | No | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 3, 5µm | 80Å | 200m ² /g | — | Polymeric | No | L8 | |
| | Silica | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | — | No | L3 | |
| | SAX | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | Monomeric | No | — | |
| | SCX | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | Monomeric | Yes | — | |
| Adsorbosphere™ XL <i>Alltech</i> | C18 | Silica | Spherical | 3, 5µm | 90Å | 200m ² /g | 11% | Monomeric | Yes | L1 | |
| | C18-B | Silica | Spherical | 5µm | 90Å | 200m ² /g | 12% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5µm | 90Å | 200m ² /g | 6% | Monomeric | Yes | L7 | |
| | C1 (TMS) | Silica | Spherical | 5µm | 90Å | 200m ² /g | — | Monomeric | Yes | L13 | |
| | Silica | Silica | Spherical | 5µm | 90Å | 200m ² /g | — | — | No | L3 | |
| | SAX | Silica | Spherical | 5, 10µm | 90Å | 200m ² /g | — | Monomeric | Yes | — | |
| | SCX | Silica | Spherical | 5, 10µm | 90Å | 200m ² /g | — | Monomeric | Yes | — | |
| Allsphere™ <i>Alltech</i> | ODS-1 | Silica | Spherical | 5µm | 80Å | 220m ² /g | 7% | Monomeric | Partial | L1 | |
| | ODS-2 | Silica | Spherical | 3, 5µm | 80Å | 220m ² /g | 12% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5µm | 80Å | 220m ² /g | 6% | Monomeric | Yes | L7 | |
| | C6 | Silica | Spherical | 5µm | 80Å | 220m ² /g | 4% | Monomeric | Yes | L15 | |
| | C1 (TMS) | Silica | Spherical | 5µm | 80Å | 220m ² /g | 3% | Monomeric | No | L13 | |
| | Phenyl | Silica | Spherical | 5µm | 80Å | 220m ² /g | 3% | Monomeric | Yes | L11 | |
| | Cyano | Silica | Spherical | 5µm | 80Å | 220m ² /g | 3.5% | Monomeric | No | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 5µm | 80Å | 220m ² /g | 3% | Monomeric | No | L8 | |
| | Silica | Silica | Spherical | 3, 5µm | 80Å | 220m ² /g | — | — | No | L3 | |
| | SAX | Silica | Spherical | 5µm | 100Å | 220m ² /g | 4% | Monomeric | No | — | |
| | SCX | Silica | Spherical | 5µm | 100Å | 220m ² /g | 4% | Monomeric | No | — | |
| | Alltima™ HP <i>Alltech</i> | C18 | Silica | Spherical | 3, 5µm | 190Å | 200m ² /g | 12% | Monomeric | Yes | L1 |
| C18 EPS | | Silica | Spherical | 3, 5µm | 190Å | 200m ² /g | 4% | Monomeric | No | L1 | |
| C18 HiLoad | | Silica | Spherical | 3, 5µm | 100Å | 450m ² /g | 24% | Monomeric | Yes | L1 | |
| C18 AQ | | Silica | Spherical | 3, 5µm | 100Å | 450m ² /g | 20% | Monomeric | Yes | L1 | |
| C18 Amide | | Silica | Spherical | 3, 5µm | 190Å | 200m ² /g | 12% | Monomeric | Yes | L1 | |
| C8 | | Silica | Spherical | 3, 5µm | 190Å | 200m ² /g | 8% | Monomeric | Yes | L7 | |
| Cyano | | Silica | Spherical | 3, 5µm | 190Å | 200m ² /g | 4% | Monomeric | Yes | L10 | |
| Silica | | Silica | Spherical | 3, 5µm | 100Å | 450m ² /g | — | — | No | L3 | |
| HILIC | | Silica | Spherical | 1.5, 3, 5µm | 120Å | 230m ² /g | — | — | No | L3 | |
| C18 | | Silica | Spherical | 3, 5, 10µm | 100Å | 340m ² /g | 16% | Polymeric | Yes | L1 | |
| Alltima™ <i>Alltech</i> | C18 LL | Silica | Spherical | 5µm | 100Å | 340m ² /g | 9% | Polymeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5, 10µm | 100Å | 340m ² /g | 9% | Polymeric | Yes | L7 | |
| | Phenyl | Silica | Spherical | 3, 5µm | 100Å | 340m ² /g | 7.5% | Polymeric | Yes | L11 | |
| | Cyano | Silica | Spherical | 3, 5µm | 100Å | 340m ² /g | — | Polymeric | Yes | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 3, 5µm | 100Å | 340m ² /g | — | Polymeric | No | L8 | |
| | Silica | Silica | Spherical | 3, 5, 10µm | 100Å | 340m ² /g | — | — | No | L3 | |
| | AlphaBond™ <i>Alltech</i> | C18 | Silica | Irregular | 5, 10µm | 125Å | 300m ² /g | 10% | Monomeric | Yes | L1 |
| C8 | | Silica | Irregular | 10µm | 125Å | 300m ² /g | — | Monomeric | Yes | L7 | |
| Phenyl | | Silica | Irregular | 10µm | 125Å | 300m ² /g | — | Monomeric | Yes | L11 | |
| Cyano | | Silica | Irregular | 10µm | 125Å | 300m ² /g | — | Monomeric | Yes | L10 | |
| Amino (NH ₂) | | Silica | Irregular | 10µm | 125Å | 300m ² /g | — | Polymeric | No | L8 | |
| Silica | | Silica | Irregular | 10µm | 125Å | 300m ² /g | — | — | No | L3 | |
| Apex™ I JONES | | ODS | Silica | Spherical | 3, 5, 10µm | 100Å | 170m ² /g | 10% | Polymeric | Yes | L1 |
| | C8 | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 7% | Monomeric | No | L7 | |
| | C8(EC) | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 7% | Monomeric | Yes | L7 | |
| | C1 | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 2.5% | Monomeric | Yes | L13 | |
| | Phe | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 3% | Monomeric | No | L11 | |
| | Basic ODS | Silica | Spherical | 5µm | 100Å | 200m ² /g | 12% | Monomeric | Yes | L1 | |
| | PAH | Silica | Spherical | 5µm | 100Å | 170m ² /g | — | Monomeric | Yes | — | |
| | CN | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 4% | Monomeric | No | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 2% | Monomeric | No | L8 | |
| | Carbohydrate | Silica | Spherical | 5µm | 100Å | 170m ² /g | — | Monomeric | † | — | |
| | Silica | Silica | Spherical | 3, 5, 10µm | 100Å | 170m ² /g | — | — | No | L3 | |
| | Apex™ II JONES | ODS | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 10.5% | Monomeric | Yes | L1 |
| | | Diol | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 3.5% | Monomeric | No | — |
| Amino (NH ₂) | | Silica | Spherical | 3, 5µm | 100Å | 170m ² /g | 2% | Monomeric | No | L8 | |
| Others—as Apex I | | Silica | Spherical | 5µm | 100Å | 170m ² /g | — | — | No | — | |
| Apex™ Prepsil JONES | ODS | Silica | Spherical | 8, 15µm | 130Å | 170m ² /g | 10% | Polymeric | Yes | L1 | |
| | C8 | Silica | Spherical | 8, 15µm | 130Å | 170m ² /g | 7% | Monomeric | No | L7 | |
| | C8(EC) | Silica | Spherical | 8, 15µm | 130Å | 170m ² /g | 7% | Monomeric | Yes | L7 | |
| | C2 | Silica | Spherical | 8µm | 130Å | 170m ² /g | 2.8% | Monomeric | No | L30 | |
| | CN | Silica | Spherical | 8µm | 130Å | 170m ² /g | 4% | Monomeric | Yes | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 8µm | 130Å | 170m ² /g | 2% | Monomeric | Yes | L8 | |
| | Silica | Silica | Spherical | 8, 15µm | 130Å | 170m ² /g | — | — | No | L3 | |
| | Diol | Silica | Spherical | 8µm | 130Å | 170m ² /g | 3.2% | Monomeric | No | L20 | |

†Proprietary.

| Columns for Small Molecules (continued) | | | | | | | | | | | |
|-----------------------------------------|-------------------------------------------------|---------------|----------------|--------------------|-------------------------------------|-------------------------------------|----------------------|------------|-------------|------------|----|
| Brand | Phase | Base Material | Particle Shape | Particle Size | Pore Size | Surface Area | Carbon Load | Phase Type | End-capped? | USP L-code | |
| Apollo™ <i>Alltech</i> | C18 | Silica | Spherical | 5µm | 100Å | 340m ² /g | 15% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 5µm | 100Å | 340m ² /g | 9% | Monomeric | Yes | L7 | |
| | Phenyl | Silica | Spherical | 5µm | 100Å | 340m ² /g | 8% | Monomeric | Yes | L11 | |
| | Silica | Silica | Spherical | 5µm | 100Å | 340m ² /g | — | — | No | L3 | |
| Brava™ <i>Alltech</i> | C18 BDS | Silica | Spherical | 3, 5µm | 145Å | 185m ² /g | 8.5% | Monomeric | Yes | L1 | |
| | C18 ODS | Silica | Spherical | 3, 5µm | 130Å | 195m ² /g | 8.5% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5µm | 130Å | 195m ² /g | 6% | Monomeric | Yes | L7 | |
| | C8 BDS | Silica | Spherical | 3, 5µm | 145Å | 185m ² /g | 5.5% | Monomeric | Yes | L7 | |
| | Phenyl | Silica | Spherical | 5µm | 130Å | 195m ² /g | — | Monomeric | No | L11 | |
| | Cyano | Silica | Spherical | 5µm | 130Å | 195m ² /g | — | Monomeric | No | L10 | |
| | Cyano BDS | Silica | Spherical | 5µm | 145Å | 185m ² /g | — | Monomeric | No | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 5µm | 130Å | 195m ² /g | — | Monomeric | No | L8 | |
| | Silica | Silica | Spherical | 5µm | 130Å | 195m ² /g | — | — | No | L3 | |
| | Carbohydrate <i>Alltech</i> | Amino | Silica | Irregular | 10µm | 80Å | 550m ² /g | — | Polymeric | No | — |
| Cation | | Polymer | Spherical | 10µm | — | — | — | — | No | — | |
| Denali® <i>VDAC</i> | 238DE C18 | Silica | Spherical | 3, 5, 10, 15, 20µm | 120Å | 280-340m ² /g | 20% | Monomeric | Yes | L1 | |
| Econosil™* <i>Alltech</i> | C18 | Silica | Irregular | 3, 5, 10µm | 80Å | 200m ² /g | 10% | Monomeric | Yes | L1 | |
| | C8 | Silica | Irregular | 3, 5, 10µm | 80Å | 200m ² /g | 5% | Monomeric | Yes | L7 | |
| | CN | Silica | Irregular | 5, 10µm | 80Å | 200m ² /g | — | Monomeric | Yes | L10 | |
| | NH ₂ | Silica | Irregular | 5, 10µm | 80Å | 200m ² /g | — | Polymeric | No | L8 | |
| | Silica | Silica | Irregular | 3, 5, 10µm | 80Å | 200m ² /g | — | — | No | L3 | |
| Econosphere™* <i>Alltech</i> | C18 | Silica | Spherical | 3, 5, 10µm | 80Å | 200m ² /g | 10% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 5, 10µm | 80Å | 200m ² /g | 5% | Monomeric | Yes | L7 | |
| | Cyano | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | Monomeric | Yes | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 5µm | 80Å | 200m ² /g | — | Polymeric | No | L8 | |
| | Silica | Silica | Spherical | 3, 5, 10µm | 80Å | 200m ² /g | — | — | No | L3 | |
| Genesis® 120 <i>JONES</i> | C18 | Silica | Spherical | 3, 4, 7, 15µm | 120Å | 300m ² /g | 18% | Monomeric | Yes | L1 | |
| | C18 AQ | Silica | Spherical | 4, 7µm | 120Å | 300m ² /g | 15% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 3, 4, 7, 15µm | 120Å | 300m ² /g | 11% | Monomeric | No | L7 | |
| | C8(EC) | Silica | Spherical | 3, 4, 7, 15µm | 120Å | 300m ² /g | 11% | Monomeric | Yes | L7 | |
| | C4 | Silica | Spherical | 4µm | 120Å | 300m ² /g | 6% | Monomeric | Yes | L26 | |
| | Phenyl | Silica | Spherical | 4µm | 120Å | 300m ² /g | 9% | Monomeric | Yes | L11 | |
| | CN | Silica | Spherical | 3, 4µm | 120Å | 300m ² /g | 7% | Monomeric | Yes | L10 | |
| | Amino (NH ₂) | Silica | Spherical | 3, 4µm | 120Å | 300m ² /g | 3.5% | Polymeric | No | L8 | |
| | Carbohydrate | Silica | Spherical | 4µm | 120Å | 300m ² /g | — | Monomeric | No | — | |
| | CN-TCA | Silica | Spherical | 4µm | 120Å | 300m ² /g | 7% | Monomeric | Yes | — | |
| | Petro-XP | Silica | Spherical | 4µm | 120Å | 300m ² /g | — | Monomeric | No | — | |
| | Silica | Silica | Spherical | 3, 4, 7, 15µm | 120Å | 300m ² /g | — | — | No | L3 | |
| | Grace Alpha® <i>GRACE</i> | C18 | Silica | Spherical | 5, 10, 15, 20µm | 120Å | 325m ² /g | 15% | Monomeric | Yes | L1 |
| | | C8 | Silica | Spherical | 5, 10, 15, 20µm | 120Å | 325m ² /g | 10% | Monomeric | No | L7 |
| Silica | | Silica | Spherical | 5, 10, 15, 20µm | 120Å | 325m ² /g | — | — | No | L3 | |
| GraceSmart™ <i>GRACE</i> | C18 | Silica | Spherical | 3, 5µm | 120Å | 220m ² /g | 10% | Monomeric | Yes | L1 | |
| Grom™ Sil <i>GROM</i> | ODS-0 AB (acid/base deactivated) | Silica | Spherical | 1.5, 3, 5, 10µm | 100Å | 200m ² /g | 11% | Monomeric | Yes | L1 | |
| | ODS-2 FE (fully endcapped) | Silica | Spherical | 1.5, 3, 5, 10µm | 80, 100, 300Å | 220, 200, 100m ² /g | 12, 11, 6% | Monomeric | Yes | L1 | |
| | ODS-3 CP (encapsulated) | Silica | Spherical | 3, 5, 7, 10µm | 120, 300Å | 320, 170m ² /g | 15, 6% | Polymeric | No | L1 | |
| | ODS-4 HE (hydrophilic endcapping) | Silica | Spherical | 3, 4, 5, 7, 10µm | 120, 200Å | 300, 200m ² /g | 16, 11% | Monomeric | Yes | L1 | |
| | ODS-5 ST (standard) | Silica | Spherical | 3, 4, 5, 7, 10µm | 60, 120, 200, 300Å | 580, 300, 200, 150m ² /g | 22, 17, 12, 7% | Monomeric | Yes | L1 | |
| | ODS-6 NE (non endcapped) | Silica | Irregular | 3, 5µm | 120Å | 300m ² /g | 17% | Monomeric | No | L1 | |
| | ODS-7 pH (pH-stable) | Silica | Irregular | 4µm | 80Å | 510m ² /g | 22% | Polymeric | No | L1 | |
| | Octyl-1 B (base deactivated) | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | 6.5% | Monomeric | Yes | L7 | |
| | Octyl-2 AB (acid/base deactivated) | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | 5% | Monomeric | Yes | L7 | |
| | Octyl-3 BA (for bases) | Silica | Spherical | 3, 5µm | 120Å | 300m ² /g | 9% | Monomeric | Yes | L7 | |
| | Octyl-4 FE (fully endcapped) | Silica | Spherical | 3, 5, 10µm | 80, 100, 300Å | 220, 200, 100m ² /g | 6.6, 6, 3% | Monomeric | Yes | L7 | |
| | Octyl-5 CP (encapsulated) | Silica | Spherical | 3, 5, 7, 10µm | 120, 300Å | 320, 170m ² /g | 10, 5.5% | Polymeric | No | L7 | |
| | Octyl-6 MB (monomer binding) | Silica | Spherical | 3, 5, 10µm | 120, 200, 300Å | 300, 200, 150m ² /g | 10, 7, 4% | Monomeric | Yes | L7 | |
| | Hexyl-1 MB (monomeric bonding) | Silica | Spherical | 5µm | 80, 100Å | 220, 200m ² /g | 4, 4% | Monomeric | Yes | — | |
| | Phenyl-1 FE (fully endcapped) | Silica | Spherical | 3, 5, 10µm | 120, 300Å | 300, 150m ² /g | 9, 5% | Monomeric | Yes | L11 | |
| | Phenyl-2 CP (encapsulated) | Silica | Irregular | 5µm | 120, 300Å | 320, 170m ² /g | 7, 4% | Polymeric | No | L11 | |
| | Phenyl-3 PE (partially endcapped) | Silica | Spherical | 3, 5, 10µm | 80, 100Å | 220, 200m ² /g | 6.6, 6% | Monomeric | Yes | L11 | |
| | Butyl-1 ST (standard) | Silica | Spherical | 3, 5µm | 120, 300Å | 300, 150m ² /g | 7, 2.5% | Monomeric | No | L26 | |
| | Butyl-2 FE (fully endcapped) | Silica | Spherical | 3, 5µm | 300Å | 100m ² /g | 1, 5% | Monomeric | No | L26 | |
| | TMS-1 ST (standard) | Silica | Spherical | 3, 5µm | 120, 300Å | 300, 150m ² /g | 4% | Monomeric | Yes | L13 | |
| | TMS-2 CP (encapsulated) | Silica | Spherical | 3, 5µm | 120, 300Å | 320, 170m ² /g | 3% | Polymeric | No | L13 | |
| | Cyan-1 ST (standard) | Silica | Spherical | 3, 5µm | 120, 300Å | 300, 150m ² /g | 4.8% | Monomeric | Yes | — | |
| | Cyan-2 PR (cyanopropyl) | Silica | Spherical | 3, 5µm | 80, 100Å | 220, 200m ² /g | 3.5% | Monomeric | Yes | — | |
| | Cyan-3 CP (encapsulated) | Silica | Spherical | 5µm | 120Å | 320m ² /g | 4% | Polymeric | No | — | |
| | Amino-1 PR (NH-propyl) | Silica | Spherical | 3, 5, 10µm | 80, 100Å | 220, 200m ² /g | 2% | Monomeric | Yes | L8 | |
| | Amino-2 PA (cross linked Poly-NH ₂) | Silica | Spherical | 5µm | 120Å | 300m ² /g | — | Polymeric | No | L8 | |
| | Amino-3 CP (encapsulated NH-residues) | Silica | Irregular | 5µm | 80Å | 420m ² /g | — | Monomeric | Yes | L8 | |
| | Amino-4 PR (propylamine bonded to silica) | Silica | Irregular | 3, 7µm | 300Å | 100m ² /g | — | Monomeric | No | L8 | |
| | Diol | Silica | Spherical | 5, 10µm | 60, 120, 200, 300Å | 580, 300, 200, 150m ² /g | — | Monomeric | No | L20 | |
| | Normal Phase-1 ST (standard silica) | Silica | Spherical | 3, 5, 10µm | 80, 100, 1000Å | 220, 200m ² /g | — | — | No | L3 | |
| | Normal Phase-2 SP (spherical silica) | Silica | Spherical | 3, 5, 10µm | 60, 120, 200, 1000Å | 580, 300, 200m ² /g | — | — | No | L3 | |
| | Normal Phase-3 PV (polyvinylalcohol) | Silica | Spherical | 5µm | 120Å | 300m ² /g | — | Polymeric | No | L3 | |
| SEC (size exclusion chromatography) | Silica | Spherical | 5, 10µm | 60, 120, 200, 300Å | 580, 300, 200, 150m ² /g | — | — | No | — | | |
| Strong Anion-1 | Silica | Spherical | 5, 10µm | 80, 100Å | 220, 200m ² /g | — | — | No | — | | |
| Weak Anion-2 (ion exchange) | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | — | No | — | | |
| Strong Cation-1 (ion exchange) | Silica | Spherical | 5, 10µm | 80, 100Å | 220, 200m ² /g | — | — | No | — | | |

*Available only online.

| Columns for Small Molecules (continued) | | | | | | | | | | |
|-----------------------------------------|--------------------------------------|---------------|----------------|----------------|-----------|---------------------------|--------------------------|------------|-------------|------------|
| Brand | Phase | Base Material | Particle Shape | Particle Size | Pore Size | Surface Area | Carbon Load | Phase Type | End-capped? | USP L-code |
| Grom™ Sil (cont.) CRoM | Weak Cation-2 (ion exchange) | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | — | No | — |
| | HIC (hydrophobic interaction chrom.) | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | — | No | — |
| Grom™ Sapphire CRoM | C18 | Silica | Spherical | 3, 5, 10µm | 65, 110Å | 500, 270m ² /g | 23, 16% | Monomeric | Yes | L1 |
| | C8 | Silica | Spherical | 3, 5, 10µm | 65, 110Å | 500, 270m ² /g | 15, 10% | Monomeric | Yes | L7 |
| | C4 | Silica | Spherical | 3, 5, 10µm | 65, 110Å | 500, 270m ² /g | 10.5, 7% | Monomeric | Yes | L26 |
| | Silica | Silica | Spherical | 3, 5, 10µm | 65, 110Å | 500, 270m ² /g | — | — | No | L3 |
| Mixed Mode Alltech | C18/Cation | Silica | Spherical | 5, 7µm | 100Å | 350m ² /g | — | Polymeric | No | — |
| | C8/Anion | Silica | Spherical | 7µm | 100Å | 350m ² /g | — | Polymeric | No | — |
| | C8/Cation | Silica | Spherical | 5µm | 100Å | 350m ² /g | — | Polymeric | No | — |
| Platinum™ Alltech | C18 | Silica | Spherical | 1.5, 3, 5µm | 100Å | 200m ² /g | 6% | Monomeric | Yes | L1 |
| | C18 EPS | Silica | Spherical | 1.5, 3, 5µm | 100Å | 200m ² /g | 5% | Monomeric | No | L1 |
| | C8 | Silica | Spherical | 1.5, 3, 5µm | 100Å | 200m ² /g | 4% | Monomeric | Yes | L7 |
| | C8 EPS | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | 2.5% | Monomeric | No | L7 |
| | Phenyl | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | — | Monomeric | Yes | L11 |
| | Cyano | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | — | Monomeric | No | L10 |
| | Amino (NH ₂) | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | — | Monomeric | No | L8 |
| | Silica | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | — | — | No | L3 |
| | SAX | Silica | Spherical | 3, 5µm | 100Å | 200m ² /g | — | Monomeric | No | — |
| | Prevail™ Alltech | C18 Select | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | 17% | Monomeric | Yes |
| C18 | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | 15% | Monomeric | Yes | L1 |
| C8 | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | 8% | Monomeric | Yes | L7 |
| Phenyl | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | 7% | Monomeric | Yes | L11 |
| Cyano | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | — | Monomeric | Yes | L10 |
| Amino (NH ₂) | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | — | Monomeric | No | L8 |
| Silica | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | — | — | No | L3 |
| Organic Acid | | Silica | Spherical | 3, 5µm | 110Å | 350m ² /g | — | Monomeric | Yes | — |
| Carbohydrate ES | | Polymer | Spherical | 5µm | — | — | — | — | No | — |
| Vydac® SP VYDAC | | 101SP Sil | Silica | Spheroidal | 5, 10µm | 100Å | 250–350m ² /g | — | unbonded | No |
| | 201SP C18 | Silica | Spheroidal | 3, 5, 10, 15µm | 90Å | 250–350m ² /g | 13% | Monomeric | Yes | L1 |
| | 208SP C8 | Silica | Spheroidal | 5, 10, 15µm | 90Å | 250–350m ² /g | 9% | Monomeric | Yes | L7 |
| VisionHT™ GRACE | C18 | Silica | Spherical | 1.5µm | 100Å | 200m ² /g | 6% | Monomeric | Yes | L1 |
| | C-18-B | Silica | Spherical | 1.5µm | 120Å | 220m ² /g | 5.5% | Monomeric | † | L1 |
| | C18-P | Silica | Spherical | 1.5µm | 100Å | 200m ² /g | 5% | Monomeric | No | L1 |
| | C18-HL | Silica | Spherical | 1.5µm | 120Å | 220m ² /g | 11% | Polymeric | † | L1 |
| | HILIC, Silica | Silica | Spherical | 1.5µm | 120Å | 220m ² /g | — | — | No | L3 |

| Columns for Large Molecules | | | | | | | | | | | |
|-------------------------------------------|------------------------------|------------------------|----------------|-----------------------------|-----------------------|-------------------------|-------------------------|------------|-------------|------------|----|
| Brand | Phase | Base Material | Particle Shape | Particle Size | Pore Size | Surface Area | Carbon Load | Phase Type | End-capped? | USP L-code | |
| Genesis® 300 JONES | C18 | Silica | Spherical | 4, 7µm | 300Å | 120m ² /g | 10% | Monomeric | Yes | L1 | |
| | C4 | Silica | Spherical | 4, 7µm | 300Å | 120m ² /g | 3% | Monomeric | Yes | L26 | |
| | CN | Silica | Spherical | 4µm | 300Å | 120m ² /g | 3.3% | Monomeric | Yes | L11 | |
| Macrosphere™ 300 Alltech | C18 | Silica | Spherical | 5, 7µm | 300Å | 100m ² /g | 10% | Monomeric | Yes | L1 | |
| | C8 | Silica | Spherical | 5, 7µm | 300Å | 100m ² /g | 2.2% | Monomeric | Yes | L7 | |
| | C4 | Silica | Spherical | 5, 7µm | 300Å | 100m ² /g | — | Monomeric | Yes | L26 | |
| | SAX | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | Monomeric | No | — | |
| | WAX | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | Monomeric | No | — | |
| | SCX | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | Polymeric | No | — | |
| | WCX | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | Polymeric | No | — | |
| Macrosphere™ GPC Alltech | GPC 60 | Silica | Spherical | 7µm | 60Å | 450m ² /g | — | Polymeric | No | L25 | |
| | GPC 100 | Silica | Spherical | 7µm | 100Å | 350m ² /g | — | Polymeric | No | — | |
| | GPC 150 | Silica | Spherical | 7µm | 150Å | 200m ² /g | — | Polymeric | No | — | |
| | GPC 300 | Silica | Spherical | 7µm | 300Å | 100m ² /g | — | Polymeric | No | — | |
| ProSphere™ Alltech | C18 | Silica | Spherical | 3, 5, 10µm | 300Å | 120m ² /g | 9% | Monomeric | Yes | L1 | |
| | C18-AQ | Silica | Spherical | 3, 5µm | 100Å | 450m ² /g | 18% | Monomeric | Yes | L1 | |
| | C4 | Silica | Spherical | 3, 5, 10µm | 300Å | 120m ² /g | 3% | Monomeric | Yes | L26 | |
| | Size-Exclusion 125 | Silica | Spherical | 4, 5µm | 125Å | — | — | — | No | — | |
| | Size-Exclusion 250 | Silica | Spherical | 4, 5µm | 250Å | — | — | — | No | — | |
| | Size-Exclusion 450 | Silica | Spherical | 8µm | 450Å | — | — | — | No | — | |
| | P-HR (reversed phase) | Polymer | Spherical | 4µm | 140Å | — | — | — | No | — | |
| | C18 ProZap™ | Silica | Spherical | 1.5µm | 500Å | 59m ² /g | 3% | Monomeric | Yes | L1 | |
| ProZap™ VYDAC | C18 ProZap™ | Silica | Spherical | 1.5µm | 500Å | 59m ² /g | 3% | Monomeric | Yes | L1 | |
| | Vydac® ATP VYDAC | Silica | Spheroidal | 5, 10–15µm | 300Å | 70–110m ² /g | 3% | Polymeric | Yes | L26 | |
| Vydac® TP VYDAC | 101TP Sil | Silica | Spheroidal | 5, 10, 10–15, 15–20µm | 300Å | 70–110m ² /g | — | unbonded | No | L3 | |
| | 201TP C18 | Silica | Spheroidal | 5, 7, 10, 10–15, 15–20µm | 300Å | 70–90m ² /g | 8% | Polymeric | No | L1 | |
| | 202TP C18 | Silica | Spheroidal | 3, 5, 10µm | 300Å | 60–90m ² /g | 9% | Polymeric | No | L1 | |
| | 208TP C8 | Silica | Spheroidal | 3, 5, 7, 10, 10–15, 15–20µm | 300Å | 60–110m ² /g | 5% | Polymeric | Yes | L7 | |
| | 214TP C4 | Silica | Spheroidal | 3, 5, 7, 10, 10–15, 15–20µm | 300Å | 60–110m ² /g | 3% | Polymeric | Yes | L26 | |
| | 218TP C18 | Silica | Spheroidal | 3, 5, 7, 10, 10–15, 15–20µm | 300Å | 60–110m ² /g | 8% | Polymeric | Yes | L1 | |
| | 219TP Di-Phe | Silica | Spheroidal | 3, 5, 7, 10, 10–15, 15–20µm | 300Å | 60–110m ² /g | 4% | Polymeric | Yes | — | |
| | 238TP C18 | Silica | Spheroidal | 3, 5, 7, 10, 10–15, 15–20µm | 300Å | 60–110m ² /g | 4% | Monomeric | Yes | L1 | |
| | Everest® VYDAC | 238EV C18 | Silica | Spherical | 5, 10, 10–15, 15–20µm | 300Å | 70–110m ² /g | 6% | Monomeric | Yes | L1 |
| | | Vydac® MS VYDAC | Silica | Spheroidal | 5µm | 300Å | 70m ² /g | 5% | Polymeric | Yes | L7 |
| Vydac® MS VYDAC | 214MS C4 | Silica | Spheroidal | 5µm | 300Å | 70–110m ² /g | 3% | Polymeric | Yes | L26 | |
| | 218MS C18 | Silica | Spheroidal | 3, 5, 10, 10–15µm | 300Å | 60–110m ² /g | 8% | Polymeric | Yes | L1 | |
| | 238MS C18 | Silica | Spheroidal | 5µm | 300Å | 70m ² /g | 4% | Monomeric | Yes | L1 | |
| | 219MS Di-Phe | Silica | Spheroidal | 5µm | 300Å | 70m ² /g | 4% | Polymeric | Yes | — | |

*Product information is available at www.discoverysciences.com. †Proprietary.

Grace® HPLC Hardware Formats

Grace Davison Discovery Sciences provides hardware formats to maximize the performance of the phases in our column families for all applications. Choose a hardware format based on your method requirements for speed, sensitivity, and resolution.

Ultra High-Pressure Hardware



7135

Ideally suited for ultra high-pressure systems, this ultra-low volume hardware is pressure rated to 18,000psig, and packed with 1.5µm media to maximize speed and efficiency.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|-------------------|-----------------|----------------|---------------------|
| 20, 30, 50, 100mm | 1.0, 2.0mm | 18,000psig | 316 Stainless Steel |

Expedite™ MS Hardware



6113

High-speed, low-volume columns for rapid resolution and high-throughput LC/MS applications.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|----------|-----------------|----------------|---------------------|
| 10, 20mm | 2.1, 4.6mm | 10,000psig | 316 Stainless Steel |

Rocket™ Hardware



6904

High-speed, high-resolution columns for high-throughput analysis. Large 7mm i.d. balances column volume with system volume to deliver excellent peak shapes on conventional HPLC instrumentation. Large i.d. also allows fast mobile-phase flow rates which minimizes peak broadening.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|----------|-----------------|----------------|---------------------|
| 33, 53mm | 7mm | 10,000psig | 316 Stainless Steel |

Solvent-Reducer and Microbore Hardware



7345

Smaller diameter columns reduce solvent consumption and increase sensitivity when compared to standard 4.6mm i.d. columns. Use with standard HPLC instrumentation or with MS and ELS detectors.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|-------------------------------|-----------------|----------------|---------------------|
| 50, 100, 150, 150, 250, 300mm | 1.0, 2.1, 3.0mm | 10,000psig | 316 Stainless Steel |

Capillary Hardware



7344

For LC/MS and other high-sensitivity and sample-limited applications.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|--------------------------|-------------------------------|----------------|--------------------------------------|
| 50, 100, 150, 250, 300mm | 0.075, 0.150, 0.3, 0.5, 0.8mm | 5000psig | 316 Stainless Steel and fused silica |

Analytical Hardware



7345

4.6mm i.d. columns for standard HPLC instrumentation, the most commonly used. Analytical columns have industry standard port configurations. Some column families are also available with port configurations for Waters® endfittings.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|------------------------------|-----------------|----------------|---------------------|
| 30, 50, 100, 150, 250, 300mm | 4.6mm | 10,000psig | 316 Stainless Steel |

Metal-Free Hardware



7341

Mechanically strong, metal-free columns offer biocompatibility, and chemical resistance that ion chromatography and biotechnology applications demand.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|--------------------------|-----------------|----------------|-----------------|
| 50, 100, 150, 250, 300mm | 2.1, 4.6mm | 5000psig | PEEK |

Capillary Guards



6931

Capillary guards offer zero dead-volume with finger-tight connections to maintain column performance. Use them to protective your capillary investment, as enrichment columns, or as short analytical columns.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|-------------|-----------------|----------------|--------------------------------------|
| 5, 10, 20mm | 0.3, 0.5mm | 5000psig | 316 Stainless Steel and fused silica |

All-Guard™ Guard Cartridges



6617

Guard system with reusable holder and disposable guard columns to protect analytical columns.

| Lengths | Inner Diameters | Pressure Limit | Wetted Surfaces |
|---------|-----------------|----------------|---------------------|
| 7.5mm | 2.1, 3.0, 4.6mm | 5000psig | 316 Stainless Steel |

more info

Looking for preparative HPLC hardware? See page 150.

High Throughput HPLC Introduction

The Benefits of Grace® Small Particles and Formats

Today's laboratories are under greater demand to analyze more samples in less time. To help meet this demand, major product advances have been made to decrease analysis time and increase HPLC throughput. HPLC system advances have greatly reduced the influence of extra column volumes and have extended the range of flow rates and pressure capabilities. In response to this, Grace now offers a wide range of column chemistries and formats appropriate for all types of high throughput systems.

The use of smaller particles offers two main improvements to the chromatographic separation—increased resolution and speed. Resolution is directly proportional to the square root of column efficiency, therefore the higher the efficiency, the narrower the peaks and the greater the resolution between them. Increased speed comes from higher mobile phase flow rates that can be used without loss in efficiency, and higher flow rates mean faster analysis times.

Understanding the type of system currently in use for high throughput separations and then matching the right column configuration is critical to achieving the best high throughput separation. Here we outline the three types of systems currently in use today for high throughput and the recommended HPLC column format for use with that system.

tech tip

Converting Methods From Traditional Column Formats to High Throughput Columns

| Convert Between Standard HPLC and VisionHT™ Columns | |
|-----------------------------------------------------|-----------|
| | Flow Rate |
| Standard HPLC (4.6mm) | 1.0X |
| VisionHT™ (2.0mm) | 2.3X |

When adjusting between standard LC conditions to VisionHT™ columns convert flow rates accordingly and then increase flow rate for faster analysis.

| Convert Between Standard HPLC Columns and Rocket™ Columns | |
|-----------------------------------------------------------|-----------|
| | Flow Rate |
| Standard Analytical (4.6mm) | 1.0X |
| Rocket™ Column (7.0mm) | 2.3X |

Use this conversion of flow rate to transfer methods between Rocket™ column or VisionHT™ columns. Backpressure on standard LC systems should be considered.

System Type 1 Ultra High-Pressure LC System

(>10,000psig pressure limitation)

Examples: Agilent 1200, Waters® Acquity®, Thermo Accela™, Jasco XLC.


Speed from Ultra High-Pressure Systems:

Representing the latest in LC instrumentation technology, ultra high-pressure systems theoretically have the potential to deliver the fastest separations. They have minimal system volume and offer a pressure limit upwards of 12,000psig. This allows the use of columns with sub 2µm particles and 2–7 times traditional flow rates. Sub 2µm particles extend the working range of acceptable mobile phase linear velocities without sacrificing efficiency. Therefore, you can push flow rates and still get equal or better performance.

Suggested Column Format:


VisionHT™ columns are designed for microbore and ultra high-pressure LC systems that have small system volumes to limit peak broadening from sample diffusion. VisionHT™ columns are packed with highly efficient 1.5µm phases that yield taller peaks and increase sensitivity. The column hardware incorporates a low dead volume design to minimize sample diffusion, and maintain peak integrity and efficiency. 12,000psig high-pressure stability allows fast flow rates, decreasing run times 10 fold.

System Type 1
Recommended Column Format: VisionHT™
Length: 20, 30, 50, 100mm
i.d.: 1.0 and 2.0mm




**VisionHT™ Columns for
Ultra High-Pressure LC**

CHROM
10724



- 1. Sulfadiazine
- 2. Sulfathiazole
- 3. Sulfamerazine
- 4. Sulfamethazine
- 5. Sulfamethoxazole
- 6. Sulfaisoxazole
- 7. Sulfaquinoxilone
- 8. Sulfamethoxine

System Agilent 1200
Mobile Phase: A: 0.1% Formic Acid
 B: Methanol
Detector: UV at 280nm
Temperature: Ambient



MZ-Analystechnik GmbH
 Wöhlerstraße 2-6 • D-55120 Mainz
 Tel +49 6131 68 66 19
 Fax +49 6131 68 66 20
 e-mail: info@mz-at.de
 www.mz-at.de

more info

See pages 34–37 for more information on VisionHT™ columns.

System Type 2 Low Volume, High Throughput (HTP) LC System

(*<10,000psig Pressure Limitation*)


Examples: Shimadzu Prominence® UFLC, Hitachi Ultra, LC Packings UltiMate®, Surveyor Plus

Speed from Low-Volume HTP LC Systems:

These systems concentrate on reducing the cycle time between injections thus allowing more injections per time frame. Typically they employ high-speed gradient pumps, fast autosamplers and a quick detector sampling rates. They are not pressure rated to the extremes of ultra high-pressure systems, but the low system volume allows for short columns of narrow i.d.

Suggested Column Format:

Expedite™ and VisionHT™ columns are designed for microbore high throughput LC systems that have small system volumes. These formats incorporate a low dead volume design so that sample bands do not diffuse within the column hardware, maintaining peak integrity and efficiency. Expedite™ columns are packed with highly efficient 1.5µm or 3µm in very short column lengths to minimize backpressure and reduce analysis times. Though not as fast as ultra high-pressure LC systems these systems and columns balance both speed and backpressure.



more info

Expedite™ columns are available in most popular Alltech® brand phases. See page 24–25 for overview of phases and formats offered.

System Type 3 Traditional LC System

(*<5000psi Pressure Limitation*)


Examples: Agilent 1100, Waters® Alliance®, Thermo Surveyor®, Dionex Ultimate®, Shimadzu Prominence®, Hitachi LaChrom®

Speed from Traditional LC Systems:

To get speed from a traditional LC system, you need to consider the 5000psig pressure limitation, and the typically “large” 2mL system volume. High throughput columns on this system need to deliver a highly efficient separation typically achieved with a small particle packing (≤3µm), but without the high backpressures. The “large” system volume also needs to be balanced with an equally large column volume or the separation will be plagued with extra-column effects.

Suggested Column Format:

Rocket™ columns provide low backpressure and fast analysis times while preserving column efficiency. They are available with both 1.5µm and 3µm packing materials for use on standard HPLC systems with backpressure limits less than 5000psig. The 7mm i.d. allows faster flow rate that “sweep” the extra system volume faster and reduce peak broadening. This larger diameter also means a larger column volume to system volume ratio to minimize the efficiency loss from extra system volume. This benefit is more pronounced over 2.1 and 1mm i.d. columns that have a smaller ratio than 4.6mm i.d. columns and require much lower flow rates for acceptable backpressures. Low flow rates allow more time for sample diffusion within the standard HPLC’s system volume to further degrade the column’s efficiency.



more info

Rocket™ columns are available in most popular Alltech® brand phases. See page 24–25 for overview of phases and formats offered.