

PFAS Analysis According to EPA 533

MACHEREY-NAGEL application department · Dr. H. R. Wollseifen, T. Kretschmer, L. Emmerich

Abstract

This application note describes the determination of per- and polyfluoroalkyl substances (PFAS) according to EPA directive 533. It demonstrates the enrichment of PFAS using the SPE column CHROMABOND® HR-XAW and subsequent chromatographic separation on a NUCLEODUR® PFAS analytical HPLC column.

Introduction

Per- and polyfluoroalkyl substances (PFAS) are a group of anthropogenic chemicals widely used as additives in consumer products like fire-fighting foam, fiber coating, cookware, paper finishing, food packaging (e. g. pizza cartons, paper cups), building material, (e. g. water resistant lacquer). These persistent and bioaccumulative, anthropogenic pollutants are characterized by a linear aliphatic backbone, a high degree of fluorination and often feature a carboxylic or sulfonic acid functionality. PFAS present entail numerous analytical challenges, including their widespread presence in a variety of environmental samples, occurrence of isomers for some compounds and precursor transformations that may occur during preservation and storage of the samples. There is also evidence that exposure to PFAS can lead to adverse human health effects.

This is the reason why authorities in the US published variety of laws and regulations to protect public health and the environment [1]:

- Safe Drinking Water Act
- Toxic Substances Control Act (TSCA)
- Comprehensive Environmental Response, Compensation and Liability Act
- Clean Air Act

To protect environment and human health, the environmental protection agency (EPA) has published an action plan for identifying and for understanding PFAS e.g. new approaches to address current PFAS contamination, to prevent future contamination and to effectively communicate with the public about PFAS [2].

There is a need for robust and fast analytical methods to ensure accurate quantitation of low levels. Therefore, EPA recommends two methods to detect PFAS from drinking water by solid phase extraction (SPE) and liquid chromatography/tandem mass spectrometry (LC/MS/MS) [3, 4].

This work presents the analysis of PFAS from drinking water according to directive EPA method 533. It shows high recoveries using a weak anion exchanger based on a polystyrene-divinylbenzene copolymer (PS /DVB) in SPE columns for the enrichment of PFAS from water samples. These mixed-mode SPE phases successfully combine several interactions like ionic, hydrophobic and π - π . The extracts are finally analyzed by HPLC-MS/MS.




Figure 1: Classes of per- and polyfluoroalkyl substances (PFAS).

PFAS Analysis According to EPA 533

Sample pretreatment for solid phase extraction (SPE)

- Samples are preserved, collected and stored as presented in polypropylene bottles.
- Add ammonium acetate (1.0 g/L) to the sample. Ammonium acetate will sequester free chlorine to form chloramine.
- Verify that the sample containing 1 g/L ammonium acetate has a pH between 6.0 and 8.0. Acetic acid may be added as needed to adjust the pH.
- Add 20 µL organic standard solution* ($\beta = 25.0$ ng/mL in methanol for each compound) to the 250 mL water sample.

* Contains native and isotopically labeled per- and polyfluoroalkyl substances.

Solid phase extraction

Column: CHROMABOND® HR- XAW, 85 µm, 6 mL, 500 mg, (REF 730745)

CHROMABOND® HR-XAW, 45 µm, 3 mL, 200 mg, (REF 730748P45)

Conditioning: Rinse each cartridge with 10 mL methanol, 10 mL of aqueous 0.1 M phosphate buffer. Close the valve and add 2–3 mL of phosphate buffer (pH 7.0)* to the cartridge reservoir and fill the remaining volume with reagent water.

* Mix 500 mL of 0.1 M dibasic sodium phosphate with approximately 275 mL of 0.1 M monobasic sodium phosphate. Verify that the solution pH is approximately 7.0.

Sample application: Attach the reservoir cartridges, turn on the vacuum, and begin adding the 250 mL water sample with a flow rate of 5 mL/min to the cartridge. Do not allow the cartridge to run dry before all the sample has passed through.

Sample bottle and cartridge rinse: After the entire sample has passed through the cartridge, rinse the sample bottles and the transfer cartridges with aliquots of 1 g/L ammonium acetate in water and draw each aliquot through the SPE columns. Add 1 mL of methanol to the sample bottle and draw through the transfer cartridge and SPE cartridge. Draw air or nitrogen through the cartridge for 5 min at high vacuum (15–20 inches Hg).

Sample bottle and cartridge elution: Rinse the sample bottles and the transfer cartridge with 5 mL of methanol with 2 % ammonium hydroxide (v/v) and elute the analytes from the cartridges by pulling the 5 mL of methanol with 2 % ammonium hydroxide (v/v) through the SPE column. Use a low vacuum such that the solvent exits the cartridge in a dropwise fashion. Repeat sample bottle rinse and cartridge elution with a second 5 mL aliquot of methanol with 2 % ammonium hydroxide (v/v).

Eluent exchange: Evaporate eluate to dryness at 40 °C under a stream of nitrogen and dissolve residue in 0.5 mL methanol.

Analysis by HPLC-MS/MS

Chromatographic conditions:

Delay column: EC 50/2 NUCLEODUR® PFAS Delay , 5 µm (REF 760673.20)

Analytical column: EC 100/2 NUCLEODUR® PFAS, 3 µm (REF 760666.20)

Eluent A: 5 mM ammonium acetate in water

Eluent B: 5 mM ammonium acetate in methanol

Gradient: hold 40 % B for 1.0 min, in 8 min from 40 % B to 95 % B, hold 95 % B for 3.0 min, in 0.1 min to 40 % B, hold 40 % B for 2.9 min

Flow rate: 0.3 mL/min

Temperature: 40 °C

Injection volume: 2 µL

Sample solution: Mixture of PFAS in methanol, concentration 1 ng/mL for each compound

MS conditions :

AB Sciex QTRAP 5500

| | |
|------------------------------|------------------------------------|
| Acquisition mode: SRM | Ion spray voltage: – 4500 V |
| Interface: ESI | Temperature: 400 °C |
| Polarity: negative | Ion source gas 1: 50 psig |
| Curtain gas: 30 psig | Ion source gas 2: 60 psig |
| Collision Gas: medium | Detection window: 60 s |



PFAS Analysis According to EPA 533

MRM transitions

| Abbreviation | Analyte | CASRN | Q ₁ Mass (DA) | Q ₃ Mass (DA) | RT (min) |
|---------------------------|---|-------------|--------------------------|--------------------------|----------|
| NFDHA | Nonafluoro-3,6-dioxaheptanoic acid | 151772-58-6 | 201.027 | 85.0 | 5.19 |
| PFBA | Perfluorobutanoic acid | 375-22-4 | 212.904 | 168.8 | 2.01 |
| PFMPA | Perfluoro-3-methoxypropanoic acid | 377-73-1 | 228.933 | 85.0 | 2.64 |
| PFPeA | Perfluoropentanoic acid | 2706-90-3 | 262.880 | 219.0 | 3.90 |
| PFMBA | Perfluoro-4-methoxybutanoic acid | 863090-89-5 | 279.184 | 84.9 | 7.97 |
| HFPO-DA | Hexafluoropropylene oxide dimer acid | 13252-13-6 | 284.991 | 168.7 | 5.77 |
| PFBS | Perfluorobutanesulfonic acid | 375-73-5 | 298.933 | 98.9 | 4.20 |
| PFHxA | Perfluorohexanoic acid | 307-24-4 | 312.911 | 268.8 | 5.40 |
| PFEESA | Perfluoro(2-ethoxyethane)sulfonic acid | 113507-82-7 | 315.118 | 135.1 | 4.85 |
| 4:2FTS | 1H,1H, 2H, 2H-Perfluorohexane sulfonic acid | 757124-72-4 | 326.940 | 306.9 | 5.27 |
| PFPeS | Perfluoropentanesulfonic acid | 2706-91-4 | 348.845 | 80.0 | 5.54 |
| PFHpA | Perfluoroheptanoic acid | 375-85-9 | 362.931 | 318.8 | 6.45 |
| ADONA | 4,8-Dioxa-3H-perfluorononanoic acid | 919005-14-4 | 376.901 | 250.7 | 6.58 |
| PFHxS | Perfluorohexanesulfonic acid | 355-46-4 | 398.942 | 79.8 | 6.49 |
| PFOA | Perfluorooctanoic acid | 335-67-1 | 412.910 | 369.0 | 7.26 |
| 6:2FTS | 1H,1H, 2H, 2H-Perfluorooctane sulfonic acid | 27619-97-2 | 426.927 | 406.9 | 7.24 |
| PFHpS | Perfluoroheptanesulfonic acid | 375-92-8 | 448.929 | 79.8 | 7.26 |
| PFNA | Perfluorononanoic acid | 375-95-1 | 462.893 | 418.9 | 7.92 |
| PFOS | Perfluorooctanesulfonic acid | 1763-23-1 | 498.836 | 79.9 | 7.89 |
| 8:2FTS | 1H,1H, 2H, 2H-Perfluorodecane sulfonic acid | 39108-34-4 | 526.821 | 506.8 | 8.50 |
| PFDA | Perfluorodecanoic acid | 335-76-2 | 512.841 | 468.9 | 8.49 |
| 9Cl-PF ₃ ONS | 9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid | 756426-58-1 | 530.752 | 350.7 | 8.25 |
| PFUnA | Perfluoroundecanoic acid | 2058-94-8 | 562.801 | 518.9 | 8.95 |
| PFDoA | Perfluorododecanoic acid | 307-55-1 | 612.787 | 568.9 | 9.33 |
| 11Cl-PF ₃ OUdS | 11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid | 763051-92-9 | 630.738 | 451.0 | 9.15 |

Table 1: MRM transitions and retention times of PFAS according to EPA 533.




Figure 2: Chromatogram of PFAS according to EPA 533 on NUCLEODUR® PFAS EC 100/2 mm column ($\beta = 1.0$ ng/mL for each compound).

PFAS Analysis According to EPA 533

Recovery rates

| Native Analyte | Abbreviation | CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg | | CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg | |
|---|---------------------------|---|-------|---|-------|
| | | Recovery % | RSD % | Recovery % | RSD % |
| 11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid | 11Cl-PF ₃ OUdS | 91.9 | 6.5 | 98.6 | 4.3 |
| 1H,1H, 2H, 2H-Perfluorohexane sulfonic acid | 4:2FTS | 96.7 | 2.5 | 99.2 | 5.8 |
| 1H,1H, 2H, 2H-Perfluorooctane sulfonic acid | 6:2FTS | 60.5 | 5.1 | 55.1 | 2.4 |
| 1H,1H, 2H, 2H-Perfluorodecane sulfonic acid | 8:2FTS | 79.8 | 7.4 | 83.4 | 8.7 |
| 9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid | 9Cl-PF ₃ ONS | 100.4 | 5.4 | 95.2 | 5.8 |
| 4,8-Dioxa-3H-perfluorononanoic acid | ADONA | 98.5 | 3.4 | 104.5 | 2.7 |
| Perfluoro-1-octanesulfonamide | FOSA | 56.9 | 6.8 | 57.0 | 4.0 |
| Hexafluoropropylene oxide dimer acid | HFPO-DA | 99.9 | 5.5 | 107.7 | 1.0 |
| N-ethyl perfluorooctanesulfonamidoacetic acid | NEtFOSAA | 84.5 | 17.2 | 101.2 | 1.7 |
| Nonafluoro-3,6-dioxaheptanoic acid | NFDHA | 100.1 | 8.4 | 115.6 | 13.0 |
| N-methyl perfluorooctanesulfonamidoacetic acid | NMeFOSAA | 83.0 | 14.7 | 100.3 | 8.1 |
| Perfluorobutanoic acid | PFBA | 67.7 | 6.8 | 81.9 | 3.2 |
| Perfluorobutanesulfonic acid | PFBS | 102.4 | 5.3 | 102.5 | 3.3 |
| Perfluorodecanoic acid | PFDA | 89.9 | 3.0 | 94.7 | 6.0 |
| Perfluorododecanoic acid | PFDoA | 79.9 | 5.3 | 77.0 | 12.3 |
| Perfluoro-1-decanesulfonate | PFDS | 95.1 | 9.2 | 88.5 | 2.0 |
| Perfluoro(2-ethoxyethane)sulfonic acid | PFEESA | 96.8 | 4.7 | 101.5 | 4.9 |
| Perfluoroheptanoic acid | PFHpA | 98.1 | 4.9 | 99.9 | 3.4 |
| Perfluoroheptanesulfonic acid | PFHpS | 99.0 | 2.4 | 101.0 | 1.6 |
| Perfluorohexanoic acid | PFHxA | 98.3 | 2.8 | 96.0 | 3.1 |
| Perfluorohexanesulfonic acid | PFHxS | 103.4 | 6.2 | 100.5 | 4.4 |
| Perfluoro-4-methoxybutanoic acid | PFMBA | 94.5 | 2.4 | 92.1 | 1.6 |
| Perfluoro-3-methoxypropanoic acid | PFMPA | 94.4 | 14.3 | 82.2 | 6.3 |
| Perfluorononanoic acid | PFNA | 94.9 | 8.2 | 97.1 | 3.7 |
| Perfluoro-1-nonanesulfonate | PFNS | 100.5 | 3.4 | 87.3 | 0.9 |
| Perfluorooctanoic acid | PFOA | 95.3 | 5.7 | 82.5 | 7.5 |
| Perfluorooctanesulfonic acid | PFOS | 100.0 | 8.7 | 100.8 | 3.2 |
| Perfluoropentanoic acid | PFPeA | 91.6 | 4.4 | 88.6 | 5.1 |
| Perfluoropentanesulfonic acid | PFPeS | 102.7 | 4.6 | 96.1 | 6.9 |
| Perfluorotetradecanoic acid | PFTA | 87.1 | 17.9 | 98.6 | 1.0 |
| Perfluorotridecanoic acid | PFTTrDA | 93.0 | 12.9 | 85.3 | 8.9 |
| Perfluoroundecanoic acid | PFUnA | 97.6 | 13.0 | 89.8 | 9.2 |
| Perfluoro-1-butanefulfonamide | FBSA | 46.5 | 27.1 | 80.1 | 9.3 |
| Perfluoro-1-hexanesulfonamide | FHxSA | 48.6 | 14.6 | 59.1 | 8.7 |



PFAS Analysis According to EPA 533

| Internal Standards | Abbreviation | CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg | | CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg | |
|--|--------------|---|-------|---|-------|
| | | Recovery % | RSD % | Recovery % | RSD % |
| <i>N</i> -methyl- <i>d</i> ₃ -perfluoro-1-octanesulfonamidoacetic acid | d3-NMeFOSAA | 95.7 | 15.3 | 95.7 | 2.3 |
| <i>N</i> -ethyl- <i>d</i> ₅ -perfluoro-1-octanesulfonamidoacetic acid | d5-NEtFOSAA | 85.8 | 14.0 | 104.7 | 5.8 |
| Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-hexane sulfonate (4:2) | M2-4:2FTS | 91.7 | 4.7 | 105.2 | 9.1 |
| Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-octane sulfonate (6:2) | M2-6:2FTS | 105.9 | 9.4 | 111.6 | 4.1 |
| Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-decane sulfonate (8:2) | M2-8:2FTS | 104.0 | 7.7 | 105.4 | 5.0 |
| Perfluoro- <i>n</i> -[1,2- ¹³ C ₂]tetradecanoic acid | M2PFTA | 86.3 | 20.4 | 100.8 | 4.4 |
| Sodium perfluoro-1-[2,3,4- ¹³ C ₃]-butanesulfonate | M3PFBS | 96.2 | 2.5 | 99.5 | 3.3 |
| Sodium perfluoro-1-[1,2,3- ¹³ C ₃]-hexanesulfonate | M3PFHxS | 98.7 | 3.2 | 96.1 | 7.5 |
| Perfluoro- <i>n</i> -[1,2,3,4- ¹³ C ₄]heptanoic acid | M4PFHpA | 95.6 | 3.0 | 94.3 | 4.3 |
| Perfluoro- <i>n</i> -[1,2,3,4,6- ¹³ C ₅]hexanoic acid | M5PFHxA | 96.0 | 3.2 | 93.4 | 3.6 |
| Perfluoro- <i>n</i> -[¹³ C ₅]pentanoic acid | M5PFPeA | 97.8 | 3.3 | 91.2 | 2.6 |
| Perfluoro- <i>n</i> -[1,2,3,4,5,6- ¹³ C ₆]decanoic acid | M6PFDA | 92.3 | 7.1 | 80.7 | 9.1 |
| Perfluoro- <i>n</i> -[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid | M7PFUnA | 97.9 | 7.1 | 91.0 | 3.1 |
| Perfluoro-1-[¹³ C ₈]octanesulfonamide | M8FOSA | 56.0 | 6.1 | 53.3 | 2.4 |
| Perfluoro- <i>n</i> -[¹³ C ₈]octanoic acid | M8PFOA | 90.7 | 4.6 | 85.3 | 5.7 |
| Sodium perfluoro-1-[¹³ C ₈]octanesulfonate | M8PFOS | 99.9 | 4.3 | 96.8 | 6.3 |
| Perfluoro- <i>n</i> -[¹³ C ₉]nonanoic acid | M9PFNA | 92.5 | 4.1 | 84.7 | 8.3 |
| Perfluoro- <i>n</i> -[1,2- ¹³ C ₂]dodecanoic acid | MPFDoA | 93.7 | 12.1 | 85.1 | 3.1 |
| Perfluoro- <i>n</i> -[¹³ C ₄]butanoic acid | MPFBA | 88.6 | 6.4 | 79.9 | 4.4 |

Table 2: Recovery rates for the presented SPE method for drinking water using SPE column, CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg and SPE column, CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg.



PFAS Analysis According to EPA 533




Figure 3: Recovery rates of PFAS according to EPA 533 method for drinking water using SPE column, CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg (spiked with 0.5 ng for each compound in 0.25 L water).




Figure 4: Recovery rates of PFAS according to EPA 533 method for drinking water using SPE column, CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg (spiked with 0.5 ng for each compound in 0.25 L water).




Figure 5: Recovery rates of further PFAS* according to EPA 533 method for drinking water using SPE column, CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg (spiked with 0.5 ng for each compound in 0.25 L water). * not specified in EPA method 533




Figure 6: Recovery rates of further PFAS* according to EPA 533 method for drinking water using SPE column, CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg (spiked with 0.5 ng for each compound in 0.25 L water). * not specified in EPA method 533



PFAS Analysis According to EPA 533

Conclusion

This application note presents the reliable and successful determination of per- and polyfluoroalkyl substances (PFAS) from drinking water according to EPA method 533. By using CHROMABOND® HR-XAW SPE columns, it was possible to achieve high recovery rates for PFAS from drinking water with good reproducibility.

The used SPE phase, a weak anion exchanger based on a polystyrene-divinylbenzene copolymer (PS/DVB), is excellently suited for the enrichment of PFAS from water samples. This mixed-mode SPE phase successfully combines the required retention modes, ion exchange and reversed phase.

The work shows that both tested SPE column dimensions, CHROMABOND® HR-XAW, 45 µm, 3 mL/200 mg and CHROMABOND® HR-XAW, 85 µm, 6 mL/500 mg, fulfill the analytical requirements of EPA method 533. Most of the PFAS show recovery rates between 80 % and 120 %. Further tests showed that good results could also be obtained for several other analytes with the before-mentioned SPE products. Due to the missing acid functionality, the recovery rates of the group of perfluorooctane sulfonamides, presented in figures 5 and 6, are smaller than the average. The use of SPE phases with stronger hydrophobic interaction could lead to better enrichment. We recommend to use CHROMABOND® PFAS for this issue.

The chromatographic separation of PFAS was performed by using NUCLEODUR® PFAS HPLC column. This phase is specially suitable for the analysis of PFAS compounds. It shows high retention for polar PFAS, high MS intensity and excellent batch-to-batch reproducibility. Contaminants from an LC system, especially PFOA, are separated from sample analytes by implementing an isolator column, NUCLEODUR® PFAS Delay.

In summary, the presented application shows that the utilized SPE and HPLC products allow a reliable and successful determination of per- and polyfluoroalkyl substances (PFAS) from drinking according to EPA method 533.

References


- [1] United States Environmental Protection Agency, PFAS Laws and Regulations, <https://www.epa.gov/pfas/pfas-laws-and-regulations>
- [2] EPA's Per- and Polyfluoroalkyl Substances (PFAS) Action Plan, EPA 823R18004, February 2019, www.epa.gov/pfas
- [3] Method 537.1: Determination of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS). U.S. Environmental Protection Agency, Office of Research and Development, National Center for Environmental Assessment, Washington, DC, 2018.
- [4] Method 533: Determination of Per- and Polyfluoroalkyl Substances in Drinking Water by Isotope Dilution Anion Exchange Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry. U.S. Environmental Protection Agency, Office of Research and Development, National Center for Environmental Assessment, Washington, DC, 2019.

Product information

The following MACHEREY-NAGEL products have been used in this application note:

| | |
|---------------|---|
| REF 760673.20 | EC 50/2 NUCLEODUR® PFAS Delay |
| REF 760666.20 | EC 100/2 NUCLEODUR® PFAS |
| REF 730748P45 | CHROMABOND® HR-XAW, 45µm, 3 mL, 200 mg |
| REF 730745 | CHROMABOND® HR- XAW 85 µm, 6 mL, 500 mg |
| REF 730382 | Reservoir columns, 70 mL with adaptor for 1, 3, 6 mL CHROMABOND® SPE PP columns |
| REF 702402 | Screw closure, N 9, PP, blue, c. hole, Sili. w./ Polyimide orange, 1.0mm, flourine-free |
| REF 702009 | Screw neck vial, N 9, 11.6 x 32.0 mm, 0.3 mL, inner cone, PP tr. |

MACHEREY-NAGEL shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
 Barcelona-Allee 17 • D-55129 Mainz
 Tel +49 6131 880 96-0
 Fax +49 6131 880 96-20
 e-mail: info@mz-at.de
www.mz-at.de