

Chromatography

Flash chromatography
outperforming the standard

CHROMABOND[®]
Flash cartridges
from the silica experts



MN products for Flash chromatography

Contents

| | | |
|--|--|-----------|
| CHROMABOND® Flash solutions | separations from 10 mg – 160 g | 3 |
| CHROMABOND® Flash RS cartridges | for Isco® instruments | 4 |
| CHROMABOND® Flash BT cartridges | for Biotage® instruments | 5 |
| CHROMABOND® Flash DL cartridges | direct load cartridges | 6 |
| CHROMABOND® Flash FM cartridges | for FlashMaster™ instruments | 7 |
| MN Flash adsorbents | the core of CHROMABOND® Flash cartridges | 8 |
| Technical support for Flash RS and Flash BT | useful technical data | 9 |
| Flash method development | TLC upscaling | 10 |
| CHROMABOND® Flash accessories | for connecting Flash cartridges | 11 |



CHROMABOND® Flash cartridges

ideal for Flash separations from 10 mg up to 160 g

- **Convenient operation and reliable upscaling**

Complete program of ready-to-use Flash cartridges for

- Isco Companion® and other Teledyne Isco CombiFlash® systems
- Biotage® Isolera™, Biotage® FlashMaster™
- or as stand-alone version for all pump/detector combinations, e.g., from Biotage®, Büchi

- **Enhanced flexibility**

- All common RP and NP phases available on request
- Adsorbent weights from 4 g to 1600 g


- **Outstanding price-performance ratio**

- **Increased analytical safety**

- Organic solvent resistant, low bleed polypropylene cartridges, thick column walls, one piece body, sophisticated length-to-diameter ratio for high plate numbers and excellent separation efficiencies
- Distribution of eluent stream via highly porous frits
- High pressure stability of 15 bar (12 bar for cartridges > 200 g), good reproducibility

- **High quality standard**

All flash cartridges and adsorbents undergo comprehensive during- and after-production quality assurance measures to ensure that the products conform to the specification.



CHROMABOND® Flash RS cartridges

CHROMABOND® Flash RS

solutions for Isco® Flash instruments

- Heavy-duty polypropylene cartridges designed for use in Teledyne Isco CombiFlash® systems (Companion®, Rf etc.) without additional connectors or capillaries
- *Column connection:*
cartridges up to RS 330: **female Luer lock inlet** and **male Luer outlet**
RS 800 and RS 1600: Maxi Luers
- Using the CHROMABOND® Flash Starter Kit, REF 730798 or the CHROMABOND® Flash Stand Alone Kit, REF 732903 (see page 11) CHROMABOND® RS cartridges can also be used as stand alone system with any pump/detector/fraction collector combination (except RS 800 and RS 1600 with Maxi Luers).

CHROMABOND® Flash RS columns with Luer exit · Ordering information

Filled with standard silica, unmodified (SiOH) or endcapped octadecyl modified (C₁₈ ec);
45 µm, specific surface 500 m²/g, pH stability 2–8

For RS columns filled with other modified silica phases and aluminium oxide see www.mn-net.com

| Designation | Column length [cm] | ID [mm] | Adsorbent weight [g] | Pack of [columns] | REF |
|--|--------------------|---------|----------------------|-------------------|--------|
| CHROMABOND® Flash RS 4 SiOH | 9.8 | 12.4 | 4 | 20 | 732800 |
| CHROMABOND® Flash RS 15 SiOH | 11.6 | 21.2 | 15 | 20 | 732801 |
| CHROMABOND® Flash RS 25 SiOH | 16.5 | 21.2 | 25 | 15 | 732802 |
| CHROMABOND® Flash RS 40 SiOH | 17.1 | 26.4 | 40 | 15 | 732803 |
| CHROMABOND® Flash RS 80 SiOH | 24.0 | 30.8 | 80 | 12 | 732804 |
| CHROMABOND® Flash RS 120 SiOH | 25.5 | 36.0 | 120 | 10 | 732805 |
| CHROMABOND® Flash RS 200 SiOH | 20.0 | 60.0 | 200 | 6 | 732806 |
| CHROMABOND® Flash RS 330 SiOH | 27.0 | 60.0 | 330 | 4 | 732807 |
| CHROMABOND® Flash RS 800 SiOH | 38.5 | 82.0 | 800 | 2 | 732808 |
| CHROMABOND® Flash RS 1600 SiOH | 43.0 | 104.0 | 1600 | 2 | 732809 |
| CHROMABOND® Flash RS 4 C ₁₈ ec | 9.8 | 12.4 | 4.3 | 2 | 732810 |
| CHROMABOND® Flash RS 15 C ₁₈ ec | 11.6 | 21.2 | 16.4 | 1 | 732811 |
| CHROMABOND® Flash RS 25 C ₁₈ ec | 16.5 | 21.2 | 26 | 1 | 732812 |
| CHROMABOND® Flash RS 40 C ₁₈ ec | 17.1 | 26.4 | 43 | 1 | 732813 |
| CHROMABOND® Flash RS 80 C ₁₈ ec | 24.0 | 30.8 | 86 | 1 | 732814 |
| CHROMABOND® Flash RS 120 C ₁₈ ec | 25.5 | 36.0 | 130 | 1 | 732815 |
| CHROMABOND® Flash RS 200 C ₁₈ ec | 20.0 | 60.0 | 220 | 1 | 732816 |
| CHROMABOND® Flash RS 330 C ₁₈ ec | 27.0 | 60.0 | 360 | 1 | 732817 |
| CHROMABOND® Flash RS 800 C ₁₈ ec | 38.5 | 82.0 | 880 | 1 | 732818 |
| CHROMABOND® Flash RS 1600 C ₁₈ ec | 43.0 | 104.0 | 1760 | 1 | 732819 |

On request, all column types listed above can be packed with any adsorbent from our program of CHROMABOND® or spherical NUCLEODUR® phases (please note that other packings often result in differing adsorbent weights).



CHROMABOND® Flash BT

solutions for Biotage® Flash instruments

- Heavy-duty polypropylene cartridges designed for use in the Biotage® Isolera™ systems without additional connectors or capillaries
- *Column connection:*
female Luer lock inlet and **male Luer lock outlet**
- Using the CHROMABOND® Flash Starter Kit, REF 730798 or the CHROMABOND® Flash Stand Alone Kit, REF 732903 (see page 11) CHROMABOND® Flash BT cartridges can also be used as stand alone system with any pump / detector / fraction collector combination.



CHROMABOND® Flash BT columns with Luer lock exit · Ordering information

Filled with unmodified standard silica, 45 µm, specific surface 500 m²/g, pH stability 2–8

| Designation | Column length [cm] | ID [mm] | Adsorbent weight [g] | Pack of [columns] | REF |
|-------------------------------|--------------------|---------|----------------------|-------------------|---------------|
| CHROMABOND® Flash BT 4 SiOH | 9.8 | 12.4 | 4 | 20 | 732960 |
| CHROMABOND® Flash BT 15 SiOH | 11.6 | 21.2 | 15 | 20 | 732961 |
| CHROMABOND® Flash BT 25 SiOH | 16.5 | 21.2 | 25 | 15 | 732962 |
| CHROMABOND® Flash BT 40 SiOH | 17.1 | 26.4 | 40 | 15 | 732963 |
| CHROMABOND® Flash BT 80 SiOH | 24.0 | 30.8 | 80 | 12 | 732964 |
| CHROMABOND® Flash BT 120 SiOH | 25.5 | 36.0 | 120 | 10 | 732965 |
| CHROMABOND® Flash BT 200 SiOH | 20.0 | 60.0 | 200 | 6 | 732966 |
| CHROMABOND® Flash BT 330 SiOH | 27.0 | 60.0 | 330 | 4 | 732967 |

On request all column types listed above can be packed with any adsorbent from our program of CHROMABOND® or spherical NUCLEODUR® phases (please note that other packings often result in differing adsorbent weights).


Partly filled CHROMABOND® Flash BT cartridges (e.g., filled up to 80%) are available on request. By removal of the top cap the sample can be applied directly on to the cartridges (see page 6).



CHROMABOND® Flash DL cartridges

CHROMABOND® Flash DL cartridges solutions for direct loading

- Suitable as **solid injection** system
- For individual **self-filling** and packing of flash cartridges
- *Column connection:*
female Luer lock inlet and **male Luer lock outlet**.
Each cartridge comes with 3 filter elements: one already inserted, two more filters aside.



CHROMABOND® Flash DL cartridge filled with sample on CHROMABOND® XTR on top of CHROMABOND® Flash RS or BT silica cartridge

CHROMABOND® Flash BT cartridge partly filled with silica topped with sample on CHROMABOND® XTR

Options for solid injection

The sample is dissolved in a suitable solvent and adsorbed onto CHROMABOND® XTR (diatomaceous earth, see page 11). After removal/evaporation of the residual solvent, the adsorbent is put on top of a partly filled CHROMABOND® Flash BT cartridge or into an empty CHROMABOND® Flash DL cartridge.

CHROMABOND® Flash DL empty cartridges · Ordering information

| Designation | Column length [cm] | ID [mm] | for adsorbent weight [g] | | Volume [mL] | Empty columns Packing unit | REF | PE filter elements | |
|--------------------------|--------------------|---------|--------------------------|------------|-------------|----------------------------|---------------|--------------------|-----------------|
| | | | SiOH | Kieselguhr | | | | Packing unit | REF |
| CHROMABOND® Flash DL 4 | 9.8 | 12.4 | 4 | 3 | 8 | 50 | 732980 | 250 | 732980FE |
| CHROMABOND® Flash DL 15 | 11.6 | 21.2 | 15 | 10 | 30 | 50 | 732981 | 250 | 732981FE |
| CHROMABOND® Flash DL 25 | 16.5 | 21.2 | 25 | 15 | 45 | 50 | 732982 | 250 | 732982FE |
| CHROMABOND® Flash DL 40 | 17.1 | 26.4 | 40 | 30 | 75 | 20 | 732983 | 250 | 732983FE |
| CHROMABOND® Flash DL 80 | 24.0 | 30.8 | 80 | 60 | 160 | 20 | 732984 | 250 | 732984FE |
| CHROMABOND® Flash DL 120 | 25.5 | 36.0 | 120 | 80 | 220 | 20 | 732985 | 250 | 732985FE |
| CHROMABOND® Flash DL 200 | 20.0 | 60.0 | 200 | 150 | 410 | 10 | 732986 | 100 | 732986FE |
| CHROMABOND® Flash DL 330 | 27.0 | 60.0 | 330 | 250 | 600 | 10 | 732987 | 100 | 732987FE |



CHROMABOND® Flash FM

solutions for FlashMaster™ instruments

- Polypropylene cartridges designed for use in the Biotage® FlashMaster™ systems without additional connectors or capillaries
- *Column connection:*
open-tubular inlet and **male Luer outlet**

CHROMABOND® Flash FM columns · Ordering information

Filled with standard silica, unmodified (SiOH) or endcapped octadecyl modified (C₁₈ ec); 45 µm, specific surface 500 m²/g, pH stability 2–8

| Designation | Column length [cm] | ID [mm] | Adsorbent weight [g] | Pack of [columns] | REF |
|--|--------------------|---------|----------------------|-------------------|---------------|
| CHROMABOND® Flash FM 15/2 SiOH | 9.0 | 15.8 | 2.0 | 50 | 730881 |
| CHROMABOND® Flash FM 25/5 SiOH | 10.0 | 20.5 | 5.0 | 50 | 730891 |
| CHROMABOND® Flash FM 25/10 SiOH | 10.0 | 20.5 | 10.0 | 50 | 730666 |
| CHROMABOND® Flash FM 70/10 SiOH | 15.4 | 26.8 | 10.0 | 30 | 730885 |
| CHROMABOND® Flash FM 70/20 SiOH | 15.4 | 26.8 | 20.0 | 30 | 730915 |
| CHROMABOND® Flash FM 70/25 SiOH | 15.4 | 26.8 | 25.0 | 30 | 730892 |
| CHROMABOND® Flash FM 150/25 SiOH | 17.0 | 38.2 | 25.0 | 20 | 730667 |
| CHROMABOND® Flash FM 150/50 SiOH | 17.0 | 38.2 | 50.0 | 20 | 730887 |
| CHROMABOND® Flash FM 150/70 SiOH | 17.0 | 38.2 | 70.0 | 10 | 730880 |
| CHROMABOND® Flash FM 15/2 C ₁₈ ec | 9.0 | 15.8 | 2.0 | 50 | 730890 |
| CHROMABOND® Flash FM 25/5 C ₁₈ ec | 10.0 | 20.5 | 5.0 | 20 | 730884 |
| CHROMABOND® Flash FM 70/10 C ₁₈ ec | 15.4 | 26.8 | 10.0 | 20 | 730886 |
| CHROMABOND® Flash FM 150/50 C ₁₈ ec | 17.0 | 38.2 | 50.0 | 10 | 730888 |

Custom filling sizes are available on request.

On request all column types listed above can be packed with any adsorbent from our program of CHROMABOND® or spherical NUCLEODUR® phases (please note that other packings often result in differing adsorbent weights).



MN Flash adsorbents

- Flash columns and cartridges from MACHEREY-NAGEL are available with all CHROMABOND® SPE /Flash packings (more than 40 phases, e.g., C₁₈, C₈, OH, Alox)
Additionally you can choose from our range of POLYGOPREP silica packings in particle sizes from 20 to 130 µm and pore sizes from 60 to 4000 Å.
- For high performance Flash separations you can order columns custom-packed with spherical NUCLEODUR® featuring very high separation efficiency and extremely increased column lifetime (particle size > 20 µm)

a unique variety of phases

Technical silica information

Specification of modified and plain silica:
acid-washed irregular silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8

Additionally available silicas and particle sizes:

- Irregular POLYGOPREP silica with particle sizes of 20 to 130 µm and pore sizes of 60 to 4000 Å
- Spherical high performance silica (NUCLEODUR®, 110 Å) with particle sizes of 20 or 30 µm for high separation efficiency



Comparison of separation efficiency and price of irregular versus spherical silica


Separation efficiency and reproducibility

Our optimized automatic packing process leads to an excellent packing quality, irrespective of the phase or particle size distribution (normal phase or reversed phase, spherical or irregular particles).

MACHEREY-NAGEL, as a manufacturer of silicas, has decades of experience in the production of first class separation phases and columns. This leads to highest separation efficiencies of the columns, a constant back pressure (via controlled narrow particle size distribution) and good reproducibilities from cartridge to cartridge.

The separation efficiency is in the first place not influenced by the dimension or the geometry of the Flash RS cartridges. The chromatograms below show an identical resolution and peak form for different column dimensions, when flow and sample amount is adjusted correctly. This is advantageous for optimization and upscaling experiments.

Resolution and peak shape for different column dimensions



Loadability

- Highest loadability at the best possible separation efficiency due to narrow particle size distribution, excellent packing quality and optimized stationary phases (acid washed silica, reduced particulate matter)
- Large range of different cartridge lengths and diameters allows optimum in loadability for a given sample amount

Rule of thumb for the loadability

| Separation | Loadability | g sample / g adsorbent |
|------------|-------------|------------------------|
| difficult | low | ≤ 1 % |
| easy | high | ≥ 10 % |

Loadability table for CHROMABOND® Flash RS and BT

| SiOH cartridge | Average loadability per cartridge [g] | |
|----------------|---------------------------------------|-----------------|
| | difficult separation | easy separation |
| RS/BT 4 | 0.04 | 0.4 |
| RS/BT 15 | 0.15 | 1.5 |
| RS/BT 25 | 0.25 | 2.5 |
| RS/BT 40 | 0.4 | 4 |
| RS/BT 80 | 0.8 | 8 |
| RS/BT 120 | 1.2 | 12 |
| RS/BT 200 | 2 | 20 |
| RS/BT 330 | 3.3 | 33 |
| RS 800 | 8 | 80 |
| RS 1600 | 16 | 160 |

Back pressure and pressure stability

The back pressure always depends on flow rate and viscosity of the eluent mixture, column length and diameter as well as on the particle size. The high performance CHROMABOND® Flash cartridges up to 200 g silica are stable up to 15 bar (> 200 g: 12 bar). We recommend using a pressure guard, because short time pressure peaks (viscosity of eluent or gradient changes) can exceed the pressure limit.

Back pressure of CHROMABOND® Flash SiOH cartridges at typical flow rates (eluent hexane – ethyl acetate 9:1 or 8:2)

| Cartridge | Flow rate | | | | | | |
|-----------|-----------------------------|-----------------|-----------------|-----------------|----------------|----------------|----------------|
| | 20 mL/min | 40 mL/min | 80 mL/min | 120 mL/min | 160 mL/min | 200 mL/min | 240 mL/min |
| RS/BT 4 | 0.75 bar | 1.5 bar | | | | | |
| RS/BT 15 | 0.25 bar | 0.75 bar | 1.5 bar | 2.0 bar | | | |
| RS/BT 25 | 0.5 bar | 1.0 bar | 1.75 bar | 3.0 bar | 4.0 bar | 5.0 bar | |
| RS/BT 40 | | 0.75 bar | 1.5 bar | 2.25 bar | 3.0 bar | 3.25 bar | 3.5 bar |
| RS/BT 80 | | | 1.5 bar | 2.5 bar | 3.0 bar | 3.5 bar | 4.0 bar |
| RS/BT 120 | | | 1.0 bar | 1.5 bar | 2.0 bar | 2.5 bar | 3.0 bar |
| RS/BT 200 | | | 1.0 bar | 1.5 bar | 2.0 bar | 2.5 bar | 3.0 bar |
| RS/BT 330 | (typical flow rates) | | 1.5 bar | 2.25 bar | 3.0 bar | 3.5 bar | 4.0 bar |

Conditioning volumes for CHROMABOND® Flash RS cartridges (normally 1.5 column volumes of the eluent)

| Cartridge | Volume of eluent for conditioning | Cartridge | Volume of eluent for conditioning |
|-----------|-----------------------------------|-----------|-----------------------------------|
| RS/BT 4 | 20 mL | RS/BT 120 | 440 mL |
| RS/BT 15 | 60 mL | RS/BT 200 | 750 mL |
| RS/BT 25 | 90 mL | RS/BT 330 | 1100 mL |
| RS/BT 40 | 140 mL | RS 800 | 2900 mL |
| RS/BT 80 | 280 mL | RS 1600 | 5000 mL |

Upscaling of the optimum flow rate


This depends on the eluent and the separation problem.

For flash cartridges the upscaling relation is simple:

The flow rate is proportional to the silica weight in g (for equal eluent polarity), e.g.,

4 g silica → optimum flow: ~ 6–12 mL/min

40 g silica → optimum flow: ~ 60–120 mL/min




MN TLC and Flash products

- Same selectivity and easy upscaling from TLC to Flash separations
- Saving time and money, because expensive optimizations are not required

TLC is often used for the development of a selective and reproducible method in Flash chromatography, because it is often necessary to test a large number of eluent and/or adsorbent combinations. MN TLC plates and sheets are coated with the same base silica, which is used in our CHROMABOND® Flash cartridges. This is an important prerequisite for the reproducible transfer of a TLC separation to the Flash column, because the parameters are identical in both systems.

TLC screening

For TLC separation you should start with an unmodified silica and a nonpolar eluent of low viscosity (e.g., mixtures of *n*-hexane – ethyl acetate or *n*-hexane – acetone). By changing the composition of the eluent the R_f value of the TLC separation is adjusted to approx. 0.3. Increasing polarity of the eluent decreases the R_f values. The difference in R_f values between the substances to be separated should be at least 0.1 to allow a reliable separation in the subsequent flash chromatography. Variation of the eluent components (e.g., acetone, dichloromethane) can be used to enhance the separation by eluent-specific selectivity.



Selection from our program of TLC plates

SIL G · unmodified standard silica layers on glass plates

| Plate size [cm] | 2.5 x 7.5 | 5 x 10 | 5 x 20 | 10 x 10 | 10 x 20 | 20 x 20 | Thickness of layer | Fluorescent indicator |
|----------------------------|------------|--------|--------|---------|---------|---------|--------------------|-----------------------|
| Pack of [plates] | 100 | 50 | 100 | 25 | 50 | 25 | | |
| SIL G-25 | | 809017 | 809011 | | 809012 | 809013 | 0.25 mm | – |
| SIL G-25 UV ₂₅₄ | 809028.100 | 809027 | 809021 | 809020 | 809022 | 809023 | 0.25 mm | UV ₂₅₄ |

ALUGRAM® Xtra SIL G · unmodified standard silica layers on aluminium sheets

| Plate size [cm] | 2.5 x 7.5 | 4 x 8 | 5 x 7.5 | 5 x 10 | 5 x 20 | 10 x 20 | 20 x 20 | Thickness of layer | Fluorescent indicator |
|-------------------------|-----------|--------|-----------|--------|--------|---------|---------|--------------------|-----------------------|
| Pack of [plates] | 200 | 50 | 20 | 50 | 50 | 20 | 25 | | |
| SIL G | | | 818230.20 | 818261 | 818232 | | 818233 | 0.20 mm | – |
| SIL G/UV ₂₅₄ | 818329 | 818331 | 818330.20 | 818360 | 818332 | 818362 | 818333 | 0.20 mm | UV ₂₅₄ |

CHROMABOND® XTR

adsorbent for solid injection

- Base material coarse-grained kieselguhr (also known as diatomaceous earth, hydromatrix, celite)
Large pore size, high pore volume, constantly high batch-to-batch quality, pH working range 1–13

Ordering information

| Description | Pack of | REF |
|---------------------------|---------|--------------------|
| CHROMABOND® XTR adsorbent | 500 g | 730595.500 |
| | 1 kg | 730595.1000 |
| | 5 kg | 730595.5000 |

Accessories for CHROMABOND® Flash columns



CHROMABOND® Flash connecting kits · Ordering information

CHROMABOND® Flash connecting kits allow to use CHROMABOND® Flash RS and BT cartridges as stand-alone system with any pump/detection/fraction collector combination.

| Description | Pack of | REF |
|---|---------|---------------|
| CHROMABOND® Flash Starter Kit consists of 1/8" PTFE tubing, 1.5 mm ID, 3 m long; 5 x 1/4"-28 PP nuts; 5 x 1/8" ETFE ferrules; 5 x 1/4"-28 nylon unions; 2 x 1/4"-28 PP Luer lock, female; 1 x 1/4"-28 PP Luer lock, male; 1 x 1/4"-28 PP Luer tip, male | 1 kit | 730798 |
| CHROMABOND® Flash Stand Alone Kit, Luer consists of 1 x 1/4"-28 PP Luer lock, female; 1 x 1/4"-28 PP Luer lock, male; 2 x 1/8" ETFE ferrules; 2 x 1/4"-28 nylon unions; 2 x 1/4"-28 PP nuts | 1 kit | 732903 |

Trademarks: CHROMABOND®, NUCLEODUR® (MACHEREY-NAGEL GmbH & Co. KG, Germany); Biotage®, Isolera™, FlashMaster™ (Biotage AB, Sweden); Isco®, Companion®, CombiFlash® (Teledyne Isco Inc., USA)



HPLC



GC



TLC




SPE and Flash



Syringe filters




Vials and caps



... we Meet your Needs

Local distributor



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

Germany and international:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

Switzerland:
MACHEREY-NAGEL AG
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

France:
MACHEREY-NAGEL EURL
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

USA:
MACHEREY-NAGEL Inc.
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com



Since 1911

KATEN2001.06 Flash brochure EN1/2002.2014 PD - Printed in Germany
Image credits © catolia, Malytschyls Viktor, Gerhard Seybert - fotocoll.com