


### Rapid Protein Separations on BIOshell™ A400

<b>column:</b>	BIOshell™ A400 Protein C4, 5 cm x 2.1 mm I.D., 3.4 µm
<b>mobile phase:</b>	[A] 75:25 (0.1% TFA in water):(0.1% TFA in acetonitrile); [B] 25:75 (0.1% TFA in water):(0.1% TFA in acetonitrile)
<b>gradient:</b>	12 to 100% B in 1 min; held at 100% B for 1 min
<b>flow rate:</b>	0.4 mL/min
<b>column temp.:</b>	90 °C
<b>detector:</b>	UV, 215 nm
<b>injection:</b>	1 µL
<b>sample:</b>	Protein mix, varied concentration, water (0.1% TFA)

Due to the shallow, porous shell, protein separations can be performed in less than one minute with the BIOshell™ A400 column.

Peak Name	kDa (native)
1. lysozyme; chicken	14.3
2. haptoglobin, phenotype 1-1; human	96
3. glutamate dehydrogenase; bovine	332
4. β-galactosidase; <i>E. coli</i>	465
5. lactate dehydrogenase contaminant	
6. lactate dehydrogenase; rabbit	140





### Protein Stability Can be Measured with BIOshell™ A400

<b>column:</b>	BIOshell™ A400 Protein C4, 10 cm x 2.1 mm I.D., 3.4 µm
<b>mobile phase:</b>	[A] 0.1% TFA in water; [B] 0.1% TFA in acetonitrile
<b>gradient:</b>	55 to 70% B in 5 min
<b>flow rate:</b>	0.2 mL/min
<b>column temp.:</b>	60 °C
<b>detector:</b>	UV, 215 nm
<b>injection:</b>	1 µL
<b>sample:</b>	Oxidized filgrastim, 98 µg/mL, 0.1% aqueous hydrogen peroxide

Proteins can be oxidized overtime if not stored properly. Oxidized proteins can lead to lower efficacy or toxicity. The BIOshell™ A400 Protein C4 resolved the oxidation products of Filgrastim in under five minutes..

- |                       |                       |
|-----------------------|-----------------------|
| 1. Hydrogen Peroxide  | 4. Oxidized Product 3 |
| 2. Oxidized Product 1 | 5. Oxidized Product 4 |
| 3. Oxidized Product 2 | 6. Native Filgrastim  |



#### AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0, Fax +49 6131 880 96-20  
e-mail: info@mz-at.de, www.mz-at.de