


HPLC Analysis of Acidic Herbicides in Water on a Polymeric C18 Column after SPE using Supelclean™ ENVI™-Carb

using Zymark AutoTrace Extraction WorkStation 1.20


sample preparation SPE (Solid Phase Extraction)
 sample/matrix Fresh 1 L water samples, dechlorinated with sodium thiosulfate when necessary, at ambient temperature and pH.
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/6 mL (57092)
 condition 10 mL DI water at 20 mL/min
 sample addition 0.9 L water sample at 20 mL/min
 drying 10 min using clean nitrogen
 washing 10 mL DI water at 20 mL/min
 elution 10 mL 0.1% phosphoric acid in methylene chloride:acetonitrile (80:20) at 5 mL/min
 column polymeric-coated silica-based PAH specialty column, 20 cm × 3 mm I.D., 5 µm (Supelco equivalent, SUPELCOSIL LC-PAH, available upon request)
 mobile phase gradient, (A): 0.05% phosphoric acid in DI water; (B): acetonitrile
 flow rate 0.5 mL/min
 column temp. 50 °C
 detector photodiode array- peak width: 0.053 min, sampling interval: 0.320 sec, monitor 210 nm & 225 nm
 injection 10 µL of extract (4-5 ppb each analyte in water)
 Application No. 796-0150



Air Monitoring and Industrial Hygiene

EPA Method TO-11/IP-6A and ASTM® D5197: HPLC Analysis of Aldehydes and Ketones on SUPELCOSIL™ LC-18 after Collection/Desorption using LpDNPH

sample/matrix 5 µg each carbonyl-DNPH on LpDNPH cartridge.
 Cartridge eluted with 5.0 mL acetonitrile
 adsorbent tube LpDNPH Cartridge (21024-U)
 column SUPELCOSIL LC-18, 25 cm × 4.6 mm I.D., 5 µm particles (58298)
 mobile phase (A): acetonitrile:tetrahydrofuran:water, 30:10:60;
 (B): acetonitrile:water, 60:40 gradient 0% B for 1 min, linear gradient to 100% B over 10 min
 flow rate 1.5 mL/min
 detector VIS, 360 nm
 injection 25 µL of extract
 Application No. 795-0298



**OSHA Method 42/47 and ASTM® D5836: HPLC
Analysis of Isocyanates on SUPELCOSIL™ LC-8 after
Collection/Desorption using ORBO™-80**

sample/matrix . . . calibration standard: isocyanate derivatives, 5 µg/mL in acetonitrile
 blank: 1-(2-pyridyl)piperazine (1-2PP) coated glass fiber filter
 desorbed in acetonitrile/DMSO, 9:1

adsorbent tube ORBO-80 (20811)

column SUPELCOSIL LC-8, 25 cm × 4.6 mm I.D., 5 µm particles (58297)


mobile phase 0.05 M ammonium acetate in water:acetonitrile,
 70:30 (pH 6-6.2 with acetic acid)

flow rate linear gradient, 0.7 mL/min to 2 mL/min in 15 min

detector UV, 254 nm

injection 10 µL

Application No. 796-0376



**HPLC Analysis of Fatty Acid Methyl Esters (FAMES)
on SUPELCOSIL™ LC-18**

column SUPELCOSIL LC-18, 25 cm × 4.6 mm I.D., 5 µm particles (58298)


mobile phase (A) acetonitrile; (B) acetone; (59.0:41.0, v/v, A:B)

flow rate 1 mL/min

detector RI

injection 10 µL of 9% C8 to C18 saturated and unsaturated FAMES in mobile phase

Application No. 797-0501



HPLC Analysis of Metabolites of 7,12-Dimethylbenz[a]anthracene on SUPELCOSIL™ LC-18

column SUPELCOSIL LC-18, 25 cm × 4.6 mm I.D., 5 µm particles (58298)
 mobile phase . . . (A) methanol:(B) water, (50:50, A:B), 10 min to 100:0 at 2.5% /min
 flow rate 0.8 mL/min
 detector fluorescence
 sample rat liver, 9000 × g supernatant fraction from Aroclor-treated animals
 Application No. 713-1043

- | | |
|---|---|
| 1. 2-OH-DMBA(t-diol) | 11. 7-OHM-12-MBA-4-ol |
| 2. 7-OHM-12-MBA-(t-10,11-diol) | 12. 7-M-12-OHMBA-(3-ol) |
| 3. 7-M-12-OHMBA-(t-8,9-diol) | 13. 7-M-12-OHMBA-2-ol |
| 4. 7-OHM-12-MBA-(t-8,9-diol) | 14. 7-OHM-12-MBA |
| 5. 7-OHM-(t-3,4-diol) | 15. 7-M-12-OHMBA |
| 6. Mixed diols, incl. DMBA-(t-8,9-diol & t-10, 11-diol) | 16. DMBA-2-ol |
| 7. Position of 7,12-bis-OHMBA | 17. DMBA-3-ol |
| 8. 7-OHM-12-MBA-2-ol | 18. DMBA-4-ol |
| 9. 7-OHM-12-MBA-(3-ol) | 19. DMBA |
| 10. 7-M-12-OHMNA-4-ol | (parentheses = tentative identifications) |




Figure provided by Drs. J. Milner and J. Grunau, University of Illinois, Urbana, Illinois, USA.


Cosmetics, Personal Care, and Cleaning Products

HPLC Analysis of Antioxidants on SUPELCOSIL™ LC-18

column SUPELCOSIL LC-18, 15 cm × 4.6 mm I.D., 5 µm particles (58230-U)
 mobile phase . . . (A) 5% acetic acid in deionized water; (B) acetonitrile:methanol (1:1)
 70% A/30% B to 100% B, linear gradient over 10 min, hold 10min
 flow rate 2 mL/min
 detector UV, 280 nm
 injection 10µL, 20 µg/mL each antioxidant
 Application No. 795-0438

Gradient Program Time (min)	%A	%B
0	70	30
10	0	100
20	70	30


1. Propyl gallate
2. THBP
3. TBHQ
4. NDGA
5. BHA
6. Ionox 100
7. Octyl gallate
8. BHT
9. Dodecyl gallate



HPLC Analysis of Acidic Herbicides in Water on a Polymeric C18 Column after SPE using Supelclean™ ENVI™-Carb


using Zymark AutoTrace Extraction WorkStation 1.20

sample preparation SPE (Solid Phase Extraction)
 sample/matrix Fresh 1 L water samples, dechlorinated with sodium thiosulfate when necessary, at ambient temperature and pH.
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/6 mL (57092)
 condition 10 mL DI water at 20 mL/min
 sample addition 0.9 L water sample at 20 mL/min
 drying 10 min using clean nitrogen
 washing 10 mL DI water at 20 mL/min
 elution 10 mL 0.1% phosphoric acid in methylene chloride:acetonitrile (80:20) at 5 mL/min
 column polymeric-coated silica-based PAH specialty column, 20 cm × 3 mm I.D., 5 µm (Supelco equivalent, SUPELCOSIL LC-PAH, available upon request)
 mobile phase gradient, (A): 0.05% phosphoric acid in DI water; (B): acetonitrile
 flow rate 0.5 mL/min
 column temp. 50 °C
 detector photodiode array- peak width: 0.053 min, sampling interval: 0.320 sec, monitor 210 nm & 225 nm
 injection 10 µL of extract (4-5 ppb each analyte in water)
 Application No. 796-0150




HPLC Analysis of Carbamate and Urea Pesticides on SUPELCOSIL™ LC-8

column SUPELCOSIL LC-8, 15 cm × 4.6 mm I.D., 5 µm particles (58220-U)
 mobile phase (A) acetonitrile; (B) water
 gradient 18:82 to 65:35 in 9 min, hold 3 min
 flow rate 2 mL/min
 column temp. 35 °C
 detector UV, 240 nm
 injection 10 µL
 Application No. 794-0803



HPLC Analysis of Carbamate and Urea Pesticides on SUPELCOSIL™ LC-8 after SPME using 60 µm PDMS/DVB Fiber

sample/matrix 3 mL water containing 8ng/mL of each analyte in 10% NaCl
 SPME fiber PDMS/DVB, 60 µm (57317)
 extraction immersion, 40 min, rapid stirring
 desorption process static, 5 min in acetonitrile:water (65:35);
 dynamic, valve open during run
 column SUPELCOSIL LC-8, 15 cm × 4.6 mm I.D., 5 µm particles (58220-U)
 mobile phase (A) acetonitrile; (B) water
 gradient 18:82 to 65:35 in 9 min, hold 3 min
 flow rate 2.0 mL/min
 column temp. 35 °C
 detector UV, 240 nm
 Application No. 797-0049



HPLC Analysis of Metabolites of 7,12-Dimethylbenz[a]anthracene on SUPELCOSIL™ LC-18

column SUPELCOSIL LC-18, 25 cm × 4.6 mm I.D., 5 μm particles (58298)
 mobile phase . . . (A) methanol:(B) water, (50:50, A:B), 10 min to 100:0 at 2.5% /min
 flow rate 0.8 mL/min
 detector fluorescence
 sample rat liver, 9000 × g supernatant fraction from Aroclor-treated animals
 Application No. 713-1043

1. 2-OH-DMBA(t-diol)
2. 7-OHM-12-MBA-(t-10,11-diol)
3. 7-M-12-OHMBA-(t-8,9-diol)
4. 7-OHM-12-MBA-(t-8,9-diol)
5. 7-OHM-(t-3,4-diol)
6. Mixed diols, incl. DMBA-(t-8,9-diol & t-10, 11-diol)
7. Position of 7,12-bis-OHMBA
8. 7-OHM-12-MBA-2-ol
9. 7-OHM-12-MBA-(3-ol)
10. 7-M-12-OHMNA-4-ol
11. 7-OHM-12-MBA-4-ol
12. 7-M-12-OHMBA-(3-ol)
13. 7-M-12-OHMBA-2-ol
14. 7-OHM-12-MBA
15. 7-M-12-OHMBA
16. DMBA-2-ol
17. DMBA-3-ol
18. DMBA-4-ol
19. DMBA
(parentheses = tentative identifications)





Figure provided by Drs. J. Milner and J. Grunau, University of Illinois, Urbana, Illinois, USA.

HPLC Analysis of Nonvolatile Pesticides in Water on SUPELCOSIL™ LC-18-DB after SPE using Supelclean™ ENVI-Carb™

sample preparation SPE (Solid Phase Extraction)
 sample/matrix water spiked with pesticides at 10-50 μg/L
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/3 mL (57088)
 condition 5 mL methylene chloride:methanol (80:20); 1 mL methanol;
 10 mL 2% acetic acid in water (keep bed moist until sample addition)
 sample addition 100 mL sample at 5 mL/min
 drying 1 minute with vacuum
 elution 0.8 mL methanol; 2 x 35 mL methylene chloride:methanol (80:20)
 eluate post-treatment dry eluate to 500 μL under gentle nitrogen stream,
 reconstitute to 1 mL with methanol
 column SUPELCOSIL LC-18-DB, 25 cm × 4.6 mm I.D., 5 μm particles (58355-U)
 mobile phase (A) water:acetonitrile (90:10), (B) acetonitrile;
 gradient 80% A for 5 min then to 30% A over 30 min
 flow rate 1.5 mL/min
 detector UV, 220 nm
 injection 20 μL
 Application No. 85-438




HPLC Analysis of PAHs in Water on SUPELCOSIL™ LC-PAH after SPME using 100 μm PDMS Fiber

sample/matrix 5 μL PAH mix (Cat. No. 4-8743) in 5 mL water
 SPME fiber polydimethylsiloxane, 100 μm (57301)
 extraction immersion, 30 min (rapid stirring)
 desorption process static, 200 μL acetonitrile:water, 40:60, 2 min
 column SUPELCOSIL LC-PAH, 15 cm × 4.6 mm I.D., 5 μm particles (58318)
 gradient 0 min - 50% A; 5 min - 50% A; 30 min - 100% A; 45 min - 100% A
 flow rate 0-2 min: 0.2 mL/min 2-45 min: 1.0 mL/min
 mobile phase (A) acetonitrile:(B) water
 detector UV, 254 nm
 Application No. 796-0086

1. Naphthalene 1,000 ng/mL
2. Acenaphthylene 2,000 ng/mL
3. Acenaphthene 1,000 ng/mL
4. Fluorene 200 ng/mL
5. Phenanthrene 100 ng/mL
6. Anthracene 100 ng/mL
7. Fluoranthene 200 ng/mL
8. Pyrene 100 ng/mL
9. Benzo(a)anthracene 100 ng/mL
10. Chrysene 100 ng/mL
11. Benzo(b)fluoranthene 200 ng/mL
12. Benzo(k)fluoranthene 100 ng/mL
13. Benzo(a)pyrene 100 ng/mL
14. Dibenzo(a,h)anthracene 200 ng/mL
15. Benzo(ghi)perylene 200 ng/mL
16. Indeno(1,2,3-cd)pyrene 100 ng/mL

Gradient Program Time (min)	% ACN
0	50
5	50
30	100
45	100


Flow increased at 2.0 min



HPLC Analysis of PAHs on SUPELCOTM LC-PAH

PAHs are traditionally evaluated using an acetonitrile/water gradient. Recently customers are looking for alternatives to this method using a methanol/water gradient due to the risk of another acetonitrile shortage, solvent costs, and toxicity to the environment. This analysis demonstrates the separation using a methanol gradient.

column SUPELCOTM LC-PAH, 15 cm x 4.6 mm I.D., 5 µm particles (58318)
 mobile phase (A) water, (B) methanol
 gradient 55 to 100% B in 5 min; held at 100%B for 10 min
 flow rate 1.5 mL/min
 column temp. 35 °C
 detector UV, 254 nm
 injection 5 µL
 sample EPA 610 mix
 Application No. **G005559**




HPLC Analysis of Paraquat and Diquat on SUPELCOTM LC-18 after SPE using ENVITM-8 DSK

Refer to US EPA Method 549.1 for full details

Sample Pre-treatment:

250 mL drinking water. Adjust water sample to pH to 10.5 ± 0.2 with 10% sodium hydroxide or 10% hydrochloric acid

sample preparation SPE (Solid Phase Extraction)
 sample/matrix drinking water adjusted to pH to 10.
 sample preparation ENVI-8 DSK, 47 mm (57172)
 SPE tube/cartridge ENVI-8 DSK SPE Disk, 47 mm (57172)
 condition 10 mL acetonitrile; 2 × 10 mL reagent water; 10 mL conditioning solvent A (5 g cetyl trimethyl ammonium bromide and 5 mL conc. ammonium hydroxide in 500 mL DI water, dilute to 1 L); (2 × 10 mL reagent water; 10 mL conditioning solvent B (10 g hexanesulfonic acid, sodium salt and 10 mL conc. ammonium hydroxide in 250 mL DI water, dilute to 500 mL))
 sample addition 250 mL water 100 mL/min
 elution 0.5-1.0 mL acetonitrile (to cover/solvate disk); 2 × 4 mL eluting solution (13.5 mL orthophosphoric acid and 10.3 mL diethylamine in 500 mL DI water, dilute to 1 L)
 column SUPELCOTM LC-18, 15 cm × 4.6 mm I.D., 5 µm particles (58230-U)
 mobile phase 3.5 mL triethylamine and 1.0 g 1-hexane-sulfonic acid, sodium salt to 800 mL DI H₂O add orthophosphoric acid to pH 2.5, dilute to 1 L
 flow rate 1.0 mL/min
 column temp. 35 °C
 detector photodiode array, diquat - 308 nm, paraquat - 257 nm
 injection 100 µL
 Application No. **794-0715**




HPLC Analysis of Triazine Herbicides in Grass on SUPELCOSIL™ LC-8-DB after SPE using Supelclean™ LC-SCX

Sample Pre-treatment:


5 g fresh grass clippings spiked with 2 ppm each herbicide. Add 4 g anhydrous sodium sulfate and 20 mL methylene chloride:acetone (80:20). Shake 20 min and allow mixture to stand 1 min.

sample preparation SPE (Solid Phase Extraction)
 sample/matrix fresh grass clippings spiked with herbicides at 2 ppm
 SPE tube/cartridge Supelclean LC-SCX, 500 mg/3 mL (57018)
 condition 1 mL methylene chloride
 sample addition 2 mL grass extract. Wash with 2 x 2 mL acetonitrile.
 drying Dry packing for 5 min. under nitrogen purge
 washing Dry packing for 5 min under nitrogen purge; wash with 2 x 2 mL DI water
 elution 1.5 mL methanol
 eluate post-treatment dilute to 2 mL with DI water
 column SUPELCOSIL LC-8-DB, 15 cm x 4.6 mm I.D.,
 5 µm with guard column (58347)
 mobile phase (A) acetonitrile; (B) water; (45:55, A:B)
 flow rate 1.5 mL/min
 detector UV, 254 nm
 injection 100 µL
 Application No. [84-109](#)




HPLC Analysis of Antioxidants on SUPELCOSIL™ LC-18

column SUPELCOSIL LC-18, 15 cm x 4.6 mm I.D., 5 µm particles (58230-U)
 mobile phase (A) 5% acetic acid in deionized water; (B) acetonitrile:methanol (1:1)
 70% A/30% B to 100% B, linear gradient over 10 min, hold 10min
 flow rate 2 mL/min
 detector UV, 280 nm
 injection 10µL, 20 µg/mL each antioxidant
 Application No. [795-0438](#)




HPLC Analysis of Fatty Acid Methyl Esters (FAMES) on SUPELCOTM LC-18

column SUPELCOTM LC-18, 25 cm × 4.6 mm I.D., 5 μm particles (58298)
 mobile phase (A) acetonitrile; (B) acetone; (59.0:41.0, v/v, A:B)
 flow rate 1 mL/min
 detector RI
 injection . . . 10 μL of 9% C8 to C18 saturated and unsaturated FAMES in mobile phase
 Application No. 797-0501




HPLC Analysis of Free Fatty Acids on SUPELCOTM LC-8

column SUPELCOTM LC-8, 15 cm × 4.6 mm I.D., 3 μm particles (58983)
 mobile phase acetonitrile:tetrahydrofuran: 0.1% phosphoric acid (50.4:21.6:28)
 flow rate 1 mL/min
 column temp. 35 °C
 detector UV, 215 nm
 injection 10 μL
 Application No. 713-0969



HPLC Analysis of Sugars on SUPELCOTM LC-NH2

column SUPELCOTM LC-NH2, 25 cm × 4.6 mm I.D., 5 μm particles (58338)
 mobile phase (A) water; (B) acetonitrile; (25:75, A:B)
 flow rate 1 mL/min
 temp. ambient
 detector refractive index (RI)
 injection 2 μL
 sample 10 mg/mL in water
 Application No. G006195




HPLC Analysis of Hippuric and Methylhippuric Acids in Urine on SUPELCOSIL™ LC-18

Urinary excretion of hippuric acid and m-or p-methylhippuric acid in the urine is a test of exposure to toluene and m-or p-xylene vapors.

column SUPELCOSIL LC-18, 7.5 cm × 4.6 mm I.D., 3 μm particles (58984)
 mobile phase (A) tetrahydrofuran:10 mM potassium phosphate, pH 3.0 (B) with phosphoric acid (3:97, A:B)

flow rate 1.5 mL/min
 detector UV, 260 nm
 injection 10 μL spiked urine (300 μg/mL HA, o-MHA; 150 μg/mL m-MHA, p-MHA)
 Application No. 713-0961A



HPLC Analysis of Metabolites of 7,12-Dimethylbenz[a]anthracene on SUPELCOSIL™ LC-18

column SUPELCOSIL LC-18, 25 cm × 4.6 mm I.D., 5 μm particles (58298)
 mobile phase (A) methanol:(B) water, (50:50, A:B), 10 min to 100:0 at 2.5% /min
 flow rate 0.8 mL/min
 detector fluorescence
 sample rat liver, 9000 × g supernatant fraction from Aroclor-treated animals
 Application No. 713-1043

1. 2-OH-DMBA(t-diol)
2. 7-OHM-12-MBA-(t-10,11-diol)
3. 7-M-12-OHMBA-(t-8,9-diol)
4. 7-OHM-12-MBA-(t-8,9-diol)
5. 7-OHM-(t-3,4-diol)
6. Mixed diols, incl. DMBA-(t-8,9-diol & t-10, 11-diol)
7. Position of 7,12-bis-OHMBA
8. 7-OHM-12-MBA-2-ol
9. 7-OHM-12-MBA-(3-ol)
10. 7-M-12-OHMNA-4-ol
11. 7-OHM-12-MBA-4-ol
12. 7-M-12-OHMBA-(3-ol)
13. 7-M-12-OHMBA-2-ol
14. 7-OHM-12-MBA
15. 7-M-12-OHMBA
16. DMBA-2-ol
17. DMBA-3-ol
18. DMBA-4-ol
19. DMBA
(parentheses = tentative identifications)




Figure provided by Drs. J. Milner and J. Grunau, University of Illinois, Urbana, Illinois, USA.

HPLC Analysis of Deoxyribonucleosides and Ribonucleosides on SUPELCOSIL™ LC-18-S

column SUPELCOSIL LC-18-S, 15 cm × 4.6 mm I.D., 5 μm particles (58931)
 mobile phase 0.05 M K₂HPO₄/KH₂PO₄, pH 4.0:methanol (A) 97.5:2.5; (B) 80:20
 flow rate 1 mL/min
 column temp. 30 °C
 detector UV, 254 nm
 sample nucleoside standards in water
 Application No. 713-0955





Figure provided by Dr. C. W. Gehrke and Mr. K. C. Kuo, University of Missouri-Columbia, Experimental Station Chemical Laboratories, Columbia, MO USA.

HPLC Analysis of Nucleotides on SUPELCOTM LC-18-T

column SUPELCOTM LC-18-T, 15 cm × 4.6 mm I.D., 3 μm particles (58970-U)
 mobile phase . . . (A) 0.1M potassium phosphate buffer/4mM tetrabutylammonium
 hydrogen sulfate, pH6.0; (B) A:methanol (70:30), pH 7.2
 flow rate 1.5 mL/min
 detector UV, 254 nm
 injection 20 μL deionized water (25-530 μM each analyte)
 Application No. 713-0951



HPLC Analysis of tRNA Hydrolysate on SUPELCOTM LC-18-S

column SUPELCOTM LC-18-S, 15 cm × 4.6 mm I.D., 5 μm particles (58931)
 mobile phase . . . (A) 0.1M potassium phosphate, dibasic (pH 5.3):methanol, 97.5:2.5;
 (B) 0.01M potassium phosphate, dibasic (pH 5.1):methanol, 80:20; (C) 0.01M
 potassium phosphate, dibasic (pH 4.9):acetonitrile, 65:35
 flow rate 1 mL/min
 column temp. 26 °C
 detector UV, 210 nm, 254 nm, 280 nm
 injection 50 μL hydrolysate buffer (pH 7.8) containing
 5 μg bovine tRNA^{Leu} hydrolysate
 Application No. 713-0501





Figure provided by Dr. C.W. Gehrke and Mr. K.C. Kuo of the Interdisciplinary Chromatography-Mass Spectrometry Laboratories, University of Missouri (Columbia, MO, USA).


HPLC Analysis of Anxiolytic Drugs on SUPEL COSIL™ LC-8

column SUPEL COSIL LC-8, 15 cm × 4.6 mm I.D., 5 µm particles (58220-U)
 mobile phase methanol:acetonitrile:0.005 M KH₂PO₄ and 0.1 M ammonium acetate
 buffer (pH 6.0 with glacial acetic acid),26.5:16.5:57 (v/v/v)
 flow rate 2 mL/min
 column temp. 30 °C
 detector UV, 245 nm
 injection 10 µL
 Application No. 713-1390




HPLC Analysis of Vitamins, Fat Soluble (A and E), on SUPEL COSIL™ LC-18

column SUPEL COSIL LC-18, 25 cm × 4.6 mm I.D., 5 µm particles (58298)
 mobile phase (A) methanol: (B) deionized water (98:2, A:B)
 flow rate 2 mL/min
 detector UV, 325 nm (retinol, retinyl acetate) or
 290 nm (tocopherols, tocopherol acetate)
 injection 20 µL
 sample 0.2-1 mg/mL each analyte in methanol
 Application No. 712-0679



HPLC Analysis of Acidic Herbicides in Water on a Polymeric C18 Column after SPE using Supelclean™ ENVI™-Carb

using Zymark AutoTrace Extraction WorkStation 1.20
 sample preparation SPE (Solid Phase Extraction)
 sample/matrix Fresh 1 L water samples, dechlorinated with sodium thiosulfate
 when necessary, at ambient temperature and pH.
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/6 mL (57092)
 condition 10 mL DI water at 20 mL/min
 sample addition 0.9 L water sample at 20 mL/min
 drying 10 min using clean nitrogen
 washing 10 mL DI water at 20 mL/min
 elution 10 mL 0.1% phosphoric acid in methylene
 chloride:acetonitrile (80:20) at 5 mL/min
 column polymeric-coated silica-based PAH specialty column, 20 cm × 3 mm I.D.,
 5 µm (Supelco equivalent, SUPEL COSIL LC-PAH, available upon request)
 mobile phase gradient, (A): 0.05% phosphoric acid in DI water; (B): acetonitrile
 flow rate 0.5 mL/min
 column temp. 50 °C
 detector photodiode array- peak width: 0.053 min,
 sampling interval: 0.320 sec, monitor 210 nm & 225 nm
 injection 10 µL of extract (4-5 ppb each analyte in water)
 Application No. 796-0150



Figures provided by A. Lichtman, Nassau County Dept. of Health, Hempstead, New York, USA. Request Application Note No. T196100.

HPLC Analysis of Paraquat and Diquat on SUPELCOSIL™ LC-18 after SPE using ENVI™-8 DSK

Refer to US EPA Method 549.1 for full details

Sample Pre-treatment:

250 mL drinking water. Adjust water sample to pH to 10.5 ± 0.2 with 10% sodium hydroxide or 10% hydrochloric acid

sample preparation SPE (Solid Phase Extraction)

sample/matrix drinking water adjusted to pH to 10.

SPE tube/cartridge ENVI-8 DSK SPE Disk, 47 mm (57172)

sample preparation ENVI-8 DSK, 47 mm (57172)

condition 10 mL acetonitrile; 2 × 10 mL reagent water; 10 mL conditioning solvent

A (5 g cetyl trimethyl ammonium bromide and 5 mL conc. ammonium

hydroxide in 500 mL DI water, dilute to 1 L);

(2 × 10 mL reagent water; 10 mL conditioning

solvent B (10 g hexanesulfonic acid, sodium salt and 10 mL conc. ammonium

hydroxide in 250 mL DI water, dilute to 500 mL))

sample addition 250 mL water 100 mL/min

elution 0.5-1.0 mL acetonitrile (to cover/solvate disk); 2 × 4 mL eluting solution

(13.5 mL orthophosphoric acid and 10.3 mL diethylamine

in 500 mL DI water, dilute to 1 L)

column SUPELCOSIL LC-18, 15 cm × 4.6 mm I.D., 5 μm particles (58230-U)

mobile phase 3.5 mL triethylamine and 1.0 g 1-hexane-sulfonic acid, sodium salt to

800 mL DI H₂O add orthophosphoric acid to pH 2.5, dilute to 1 L


flow rate 1.0 mL/min

column temp. 35 °C

detector photodiode array, diquat - 308 nm, paraquat - 257 nm

injection 100 μL

Application No. 794-0715



HPLC Analysis of Triazine Herbicides in Grass on SUPELCOSIL™ LC-8-DB after SPE using Supelclean™ LC-SCX

Sample Pre-treatment:

5 g fresh grass clippings spiked with 2 ppm each herbicide. Add 4 g anhydrous sodium sulfate and 20 mL methylene chloride:acetone (80:20). Shake 20 min and allow mixture to stand 1 min.

sample preparation SPE (Solid Phase Extraction)

sample/matrix fresh grass clippings spiked with herbicides at 2 ppm

SPE tube/cartridge Supelclean LC-SCX, 500 mg/3 mL (57018)

condition 1 mL methylene chloride

sample addition 2 mL grass extract. Wash with 2 × 2 mL acetonitrile.

drying Dry packing for 5 min. under nitrogen purge

washing Dry packing for 5 min under nitrogen purge; wash with 2 × 2 mL DI water

elution 1.5 mL methanol

eluate post-treatment dilute to 2 mL with DI water

column SUPELCOSIL LC-8-DB, 15 cm × 4.6 mm I.D.,

5 μm with guard column (58347)


mobile phase (A) acetonitrile; (B) water; (45:55, A:B)

flow rate 1.5 mL/min

detector UV, 254 nm

injection 100 μL


Application No. 84-109



HPLC Analysis of Acidic Herbicides in Water on a Polymeric C18 Column after SPE using Supelclean™ ENVI™-Carb


using Zymark AutoTrace Extraction WorkStation 1.20

sample preparation SPE (Solid Phase Extraction)
 sample/matrix Fresh 1 L water samples, dechlorinated with sodium thiosulfate when necessary, at ambient temperature and pH.
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/6 mL (57092)
 condition 10 mL DI water at 20 mL/min
 sample addition 0.9 L water sample at 20 mL/min
 drying 10 min using clean nitrogen
 washing 10 mL DI water at 20 mL/min
 elution 10 mL 0.1% phosphoric acid in methylene chloride:acetonitrile (80:20) at 5 mL/min
 column polymeric-coated silica-based PAH specialty column, 20 cm × 3 mm I.D., 5 μm (Supelco equivalent, SUPELCOSIL LC-PAH, available upon request)
 mobile phase gradient, (A): 0.05% phosphoric acid in DI water; (B): acetonitrile
 flow rate 0.5 mL/min
 column temp. 50 °C
 detector photodiode array- peak width: 0.053 min, sampling interval: 0.320 sec, monitor 210 nm & 225 nm
 injection 10 μL of extract (4-5 ppb each analyte in water)
 Application No. 796-0150



HPLC Analysis of Nonvolatile Pesticides in Water on SUPELCOSIL™ LC-18-DB after SPE using Supelclean™ ENVI™-Carb™

sample preparation SPE (Solid Phase Extraction)
 sample/matrix water spiked with pesticides at 10-50 μg/L
 SPE tube/cartridge Supelclean ENVI-Carb, 250 mg/3 mL (57088)
 condition 5 mL methylene chloride:methanol (80:20); 1 mL methanol; 10 mL 2% acetic acid in water (keep bed moist until sample addition)
 sample addition 100 mL sample at 5 mL/min
 drying 1 minute with vacuum
 elution 0.8 mL methanol; 2 x 35 mL methylene chloride:methanol (80:20)
 eluate post-treatment dry eluate to 500 μL under gentle nitrogen stream, reconstitute to 1 mL with methanol
 column SUPELCOSIL LC-18-DB, 25 cm × 4.6 mm I.D., 5 μm particles (58355-U)
 mobile phase (A) water:acetonitrile (90:10), (B) acetonitrile; gradient 80% A for 5 min then to 30% A over 30 min
 flow rate 1.5 mL/min
 detector UV, 220 nm
 injection 20 μL
 Application No. 85-438




HPLC Analysis of Paraquat and Diquat on SUPELCOSIL™ LC-18 after SPE using ENVI™-8 DSK

Refer to US EPA Method 549.1 for full details

Sample Pre-treatment:

250 mL drinking water. Adjust water sample to pH to 10.5 ± 0.2 with 10% sodium hydroxide or 10% hydrochloric acid

sample preparation SPE (Solid Phase Extraction)
 sample/matrix drinking water adjusted to pH to 10.
 SPE tube/cartridge ENVI-8 DSK SPE Disk, 47 mm (57172)
 sample preparation ENVI-8 DSK, 47 mm (57172)
 condition 10 mL acetonitrile; 2 × 10 mL reagent water; 10 mL conditioning solvent
 A (5 g cetyl trimethyl ammonium bromide and 5 mL conc. ammonium hydroxide in 500 mL DI water, dilute to 1 L);
 (2 × 10 mL reagent water; 10 mL conditioning solvent B (10 g hexanesulfonic acid, sodium salt and 10 mL conc. ammonium hydroxide in 250 mL DI water, dilute to 500 mL))
 sample addition 250 mL water 100 mL/min
 elution 0.5-1.0 mL acetonitrile (to cover/solvate disk);
 2 × 4 mL eluting solution (13.5 mL orthophosphoric acid and 10.3 mL diethylamine in 500 mL DI water, dilute to 1 L)
 column SUPELCOSIL LC-18, 15 cm × 4.6 mm I.D., 5 μm particles (58230-U)
 mobile phase 3.5 mL triethylamine and 1.0 g 1-hexane-sulfonic acid, sodium salt to 800 mL DI H₂O add orthophosphoric acid to pH 2.5, dilute to 1 L
 flow rate 1.0 mL/min
 column temp. 35 °C
 detector photodiode array, diquat - 308 nm, paraquat - 257 nm
 injection 100 μL
 Application No. 794-0715



HPLC Analysis of Triazine Herbicides in Grass on SUPELCOSIL™ LC-8-DB after SPE using Supelclean™ LC-SCX

Sample Pre-treatment:

5 g fresh grass clippings spiked with 2 ppm each herbicide. Add 4 g anhydrous sodium sulfate and 20 mL methylene chloride:acetone (80:20). Shake 20 min and allow mixture to stand 1 min.

sample preparation SPE (Solid Phase Extraction)
 sample/matrix fresh grass clippings spiked with herbicides at 2 ppm
 SPE tube/cartridge Supelclean LC-SCX, 500 mg/3 mL (57018)
 condition 1 mL methylene chloride
 sample addition 2 mL grass extract. Wash with 2 x 2 mL acetonitrile.
 drying Dry packing for 5 min. under nitrogen purge
 washing Dry packing for 5 min under nitrogen purge; wash with 2 x 2 mL DI water
 elution 1.5 mL methanol
 eluate post-treatment dilute to 2 mL with DI water
 column SUPELCOSIL LC-8-DB, 15 cm × 4.6 mm I.D.,
 5 μm with guard column (58347)
 mobile phase (A) acetonitrile; (B) water; (45:55, A:B)
 flow rate 1.5 mL/min
 detector UV, 254 nm
 injection 100 μL
 Application No. 84-109

