



# Glyphosate: Herbicide

## ANP Retention of an Extremely Polar Compound



**Notes:** Glyphosate is a nonselective herbicide and is used for the control of a wide range of weeds. It is strongly retained on soil components and due to its long half life and solubility in water can be detected long after application or far away from the application site.

### Method Conditions

**Column:** Cogent Diamond Hydride™, 4µm, 100Å

**Catalog No.:** 70000-15P-2

**Dimensions:** 2.1 x 150 mm

**Solvents:** A: DI H<sub>2</sub>O + 5 mM ammonium acetate  
B: 90% acetonitrile/ 10% DI H<sub>2</sub>O/ 10 mM ammonium acetate

Gradient:	time (min.)	%B
	0	80
	1	80
	1.1	5
	5	5
	6	80

**Post Time:** 5 min

**Flow rate:** 0.5 mL/min


**Detection:** ESI - neg - Agilent 6210 MSD TOF mass spectrometer

**Sample:** Sample stock solution was purchased from Sigma (1000 mcg/mL). Sample for analysis was made by diluting the stock 1:100 in 30:70 solution A and B

**Peaks:** 1. Glyphosate: 168.0067 m/z (M-H)<sup>-</sup>  
Figure A: injection #1, RT = 3.365 min  
Figure B: injection #5, RT = 3.366 min

### Discussion

A reliable method for the determination of glyphosate is presented. Analysis was performed using a Cogent Diamond Hydride HPLC column which provides very reproducible retention and fast equilibration even when a gradient analysis is used. The use of LC-MS detection allows avoiding time consuming derivatization of this compound which is lacking a chromophore for UV detection. The method shown, with ANP retention avoids derivatization, which is required when an ordinary C18 column is used.



**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20  
e-mail: info@mz-at.de  
www.mz-at.de