




COSMOSIL

Pentabromobenzyl Group Bonded HPLC Column COSMOSIL PBr

- **Separate hydrophilic compounds in reversed-phase conditions**
- **Simple mobile phase condition compared to HILIC**
- **More sample loading capacity than HILIC**
- **Alternative selectivity to C₁₈ column**





Separation of hydrophilic compounds in reversed-phase conditions

Hydrophilic interaction chromatography (HILIC) is an increasingly popular analytical technique. However, it is often difficult to develop a robust method due to users' less familiarity with HILIC conditions. Furthermore, high concentration of acetonitrile used in HILIC mobile phase makes it extremely sensitive to samples' water concentration; injecting samples in high water concentration often results in poor peak shapes. COSMOSIL PBr enables separation of hydrophilic compounds in reversed-phase conditions, maintaining sharp peak shapes even with aqueous samples.

COSMOSIL Application Data


Column: HILIC
Column size: 4.6mm I.D.-150mm
Mobile phase: Acetonitrile/ H₂O = 90/10
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV260nm

Sample: 1; Uracil (0.1mg/ml)
2; Uridine (0.2mg/ml)
Sam. Solution H₂O



Column: PBr
Column size: 4.6mm I.D.-150mm
Mobile phase: H₂O
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV260nm

Sample: 1; Uracil (0.1mg/ml)
2; Uridine (0.2mg/ml)
Sam. Solution H₂O



NACALAI TESQUE, INC



Comparison with C₁₈ Columns

COSMOSIL PBr retains hydrophilic compounds stronger than C₁₈ columns under the same reversed-phase conditions.

COSMOSIL Application Data

Column: 5C₁₈-PAQ
Column size: 4.6mm I.D.-250mm
Mobile phase: 20mmol/l Phosphate buffer(pH7.0)
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV225nm

Sample: 1; Ammelide (0.1mg/ml)
2; Cyanuric Acid (0.5mg/ml)
3; Ammeline (0.1mg/ml)
4; Melamine (0.1mg/ml)
Inj. Vol.: 1.0µl




T₀; Retention time of unretained peak

NACALAI TESQUE, INC

Separation Mechanism

In reversed-phase chromatography, compounds are separated by difference in hydrophobicity (the number of CH₂ bases, see chromatogram on the left (1)). Compounds with little hydrophobicity (the number of OCH₂CH₂, see chromatogram on the right (2)), are not retained by a C₁₈ column. However, these compounds can easily be separated by dispersion force interaction using the COSMOSIL PBr column.



► Dispersion force (instantaneous dipole–induced dipole force)

London Dispersion force is a weak intermolecular force that results from dipoles temporarily induced from random unsymmetrical electron positions in two adjacent atoms, also known as “instantaneous dipole–induced dipole force”. It is present in all molecules regardless if they are polar or non-polar. Compounds with high polarizability have stronger dispersion force.

► Compounds with stronger dispersion force


- Larger and heavier molecules
- Molecules with larger and heavier atoms (e.g. from weakest to strongest in halogens, F₂, Cl₂, Br₂, and I₂)
- Molecules with delocalized electrons and resonance (e.g. aromatic compounds)



COSMOSIL PBr column is packed with pentabromobenzyl bonded silica that enables separation by dispersion force interaction.

Suitable Compounds to Analyze by COSMOSIL PBr

COSMOSIL PBr retains cyclic compounds or amide compounds more strongly than C₁₈ columns.



(Mobile phase)

(1) Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 10/90

(2) 20mmol/l Phosphate buffer(pH2.5)

(3) 20mmol/l Phosphate buffer(pH7.0)

(4) Methanol/ H₂O = 10/90

Applications

Gentamicin


COSMOSIL Application Data

Column: PBr
 Column size: 2.0mm I.D.-150mm
 Mobile phase: 0.1% Pentafluoropropionic Acid
 -Acetonitrile/ H₂O = 10/90

Flow rate: 0.2 ml/min
 Temperature: 40°C
 Detection: ESI-MS, Positive, SIM

Sample: Gentamicin
 (C1a, C2, C2a, C1)

Inj. Vol.: 1.0µl



Nucleic Acid Base


COSMOSIL Application Data

Column: PBr
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/ 20mmol/l Phosphate
 buffer(pH7.0) = 10/90

Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260nm

Sample: 1; Cytosine (0.05mg/ml)
 2; Uracil (0.05mg/ml)
 3; Thymine (0.05mg/ml)
 4; Guanine (0.05mg/ml)
 5; Adenine (0.05mg/ml)

Inj. Vol.: 1.0µl



Coenzyme A


COSMOSIL Application Data

Column: PBr
 Column size: 4.6mm I.D.-150mm
 Mobile phase: Methanol/ 20mmol/l Phosphate
 buffer(pH7.0) = 20/80

Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260nm

Sample: 1; Malonyl Coenzyme A (0.25mg/ml)
 2; Coenzyme A (0.25mg/ml)
 3; Acetyl Coenzyme A (0.25mg/ml)

Inj. Vol.: 2.0µl



Methylimidazole Isomers


COSMOSIL Application Data

Column: PBr
 Column size: 4.6mm I.D.-250mm
 Mobile phase: 20mmol/l Phosphate buffer(pH2.5)

Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220nm

Sample: 1; *l*-Methylimidazole (0.1mg/ml)
 2; 2-Methylimidazole (0.1mg/ml)
 3; 4-Methylimidazole (0.1mg/ml)

Inj. Vol.: 1.0µl



Column Stability Test


COSMOSIL PBr column remains stable after 1000 injections under 100% aqueous condition.

COSMOSIL Application Data

Column: PBr
 Column size: 4.6mm I.D.-150mm
 Mobile phase: 20mmol/l Phosphate buffer(pH2.5)
 Flow rate: 1.0 ml/min
 Temperature: 40°C
 Detection: UV280nm

Sample: 1; *L*-Noradrenaline (0.5mg/ml)
 2; *L*-Adrenaline (0.5mg/ml)
 3; Dopamine (0.5mg/ml)
 4; *L*-DOPA (0.5mg/ml)

Inj. Vol.: 1.0µl




Specifications

Packing Material	5PBr
Silica gel	High-purity porous spherical silica
Average particle size	5 µm
Average pore size	approx. 120 Å
Specific surface area	approx. 300 m ² /g
Stationary phase	Pentabromobenzyl group
Bonding type	Monomeric
Endcapping treatment	Yes
Carbon load	approx. 8%
Usable pH range	2-7.5

Applications for Surfactants and Halogen Compounds


Surfactants

COSMOSIL PBr shows excellent separation for surfactants compared to C₁₈ or PFP columns.



Halogen Compounds

COSMOSIL PBr shows high selectivity for halogen compounds.



Ordering Informations

Column Size I.D. x Length (mm)	Product Number	Column Size I.D. x Length (mm)	Product Number	Column Size I.D. x Length (mm)	Product Number
2.0 x 50	12943-61	10 x 50	13253-71	4.6 mm I.D. x 10 mm cartridge*	12444-14
2.0 x 100	13245-81	10 x 100	13254-61	10 x 20	12396-41
2.0 x 150	12392-81	10 x 150	13255-51	20 x 20	13256-41
2.0 x 250	13247-61	10 x 250	12397-31	4.6 mm I.D. Cartridge Holder	38009-79
3.0 x 50	12592-61	20 x 50	13257-31	*Cartridge Holder is required.	
3.0 x 100	13249-41	20 x 100	13258-21		
3.0 x 150	13250-01	20 x 150	13259-11		
3.0 x 250	13251-91	20 x 250	12398-21		
4.6 x 50	13252-81	28 x 100	13260-71		
4.6 x 150	12394-61	28 x 150	13261-61		
4.6 x 250	12395-51	28 x 250	13262-51		

For research use only, not intended for diagnostic or drug use.