



## Contact Us

 [info@nanochrom.com](mailto:info@nanochrom.com)

 [www.nanochrom.com](http://www.nanochrom.com)

 Building NE-37, 999 Jinjihu Avenue  
Suzhou Industrial Park, Jiangsu Province, China

# NANOCHROM



BioCore Columns  
for antibody and protein separation



Website




WeChat Official Account

# Introduction

Antibody-based biologics have been a fast-growing area in the bio-therapeutic market for their specificity and efficacy, for the diagnosis and treatment of a wide range of diseases, including auto-immune, cardiovascular, infectious, cancer and inflammation. In the development and production of therapeutic antibodies, the content of impurities, structural variants and post-translational modification variants must be monitored, characterized and quantified, are critical quality attributes (CQAs) to prove the stability and effectiveness of final products. Antibody-based biologics are often of complex microheterogeneity by nature. Thus, their quality control and stability evaluation are highly challenging tasks.

HPLC is an important analytical technique for characterizing and quantifying impurities and variants in antibodies and related substances. Due to the complex nature of these molecules, a variety of separation modes are employed for a thorough characterization, including size exclusion chromatography (SEC) for aggregates and fragments, ion exchange chromatography (IEC) for charged variants, hydrophobic interaction chromatography (HIC) for DAR analysis in ADC, reversed phase chromatography (RPC) for exact mass determination for an antibody and its subunits, hydrophilic interaction chromatography (HILIC) for glycans, and protein A affinity chromatography for titer analysis.



# BioCore Columns Family

The BioCore family is designed to chromatographically analyze antibodies, proteins and related substances by liquid chromatography, and serves customers in bio-technology, biopharmaceutical and academic research.

The BioCore family consists of a suite of columns in different separation modes, including BioCore SEC size exclusion columns for aggregates and fragments, BioCore WCX/SCX/WAX/SAX ion exchange columns for charged variants, BioCore HIC hydrophobic interaction columns for DAR analysis in ADC, BioCore RP columns for exact mass determination of intact antibody and subunits, BioCore Glycan HILIC columns for glycan analysis, and BioCore Protein A affinity columns for titer analysis.

	SEC	IEC	HIC	RPC	HILIC	Affinity
Product	BioCore SEC-150 BioCore SEC-300 BioCore SEC-500	BioCore WCX BioCore SCX BioCore WAX BioCore SAX	BioCore HIC	BioCore RP	BioCore Glycan	BioCore Protein A
Application	Small-molecule drugs, peptides, proteins, oligos, glycans, etc.	Charged variants in mAbs, bi-specific antibodies, ADCs and proteins	MABs and ADCs	Intact proteins and protein fragments	N-glycans of proteins	MABs and Fc fusion proteins

BioCore columns employ innovative particle technology, advanced column chemistry and rugged packing method to ensure superior performance, quality and ruggedness. All BioCore columns are based on monodispersed spherical particles via innovative production processes, which provide excellent column efficiency and consistency.

BioCore media utilize uniquely designed surface modification and grafting methodology to ensure media bio-compatibility and desired selectivity for separating antibodies, proteins and related substances. Each batch of BioCore media is produced in accordance with a strict quality management standard and tested with relevant biological molecules to ensure performance and batch-to-batch consistency.

Each BioCore column is produced by a well-developed packing method and tested individually to ensure quality and consistency. A certificate of assurance of the media and a column quality assurance report are shipped with every shipped BioCore column.



# BioCore SEC Columns

BioCore SEC is a family of high performance, size exclusion chromatography columns, designed for separating antibodies, proteins, peptides, oligonucleotides and related substances, having a broad application range in bio-tech, biopharmaceutical and academic research.

## Main Features

- Innovative particle technology: monodispersed particles for high efficiency, high mechanical strength for better column lifetime, and high pore volume for high resolution
- Advanced column chemistry for minimal secondary interaction
- Multiple pore size options for broad application range
- Robust column packing for good column lifetime
- Good column-to-column consistency



## Specification


Product Name	SEC-150	SEC-300	SEC-500
Functional Group	Diol		
Substrate	Monodispersed, high pore volume, porous silica particles		
Particle Size	1.8, 3 & 5 $\mu$ m		
Pore Size	150 Å	300 Å	500 Å
Pressure Limit	5000 psi for 1.8 $\mu$ m 2500 psi for 3 $\mu$ m 1500 psi for 5 $\mu$ m		
Temperature Limit	40 °C		
pH Range	2-8		
Calibration Curve (PEG)	500-15,000	1,000-50,000	5,000-200,000
Calibration Curve (Glucan)	1,000-30,000	2,000-100,000	20,000-500,000
Calibration Curve (Globular Protein)	5,000-150,000	10,000-750,000	20,000-1,500,000
Application	Peptides, heparin, glycans, small oligos and small proteins	Aggregates and fragments in antibodies and proteins, DNA/RNA	High order aggregates in antibodies, larger proteins and DNA/RNA

# Application

Size exclusion chromatography (SEC) is an important analytical technique for determining aggregate and fragment content in antibody biologics, which is one of the most frequently run assays in mAb drug development and drug manufacturing.

Figure 1 illustrates the separation of a mixture of a mAb (~150 kDa), its aggregate (tetramer, ~900 kDa) and fragments (HC, ~50 kDa, and LC, ~25 kDa) on a 5  $\mu$ m, 7.8x300 mm BioCore SEC-300 column, demonstrating its good suitability for simultaneous separation of mAb, related aggregates and fragments.


Separation of MAb, Aggregate, Heavy Chain and Light Chain



Columns: **BioCore SEC-300**, 5  $\mu$ m  
 Dimension: 7.8 x 300 mm  
 Mobile Phase: 300 mM NaCl in 50 mM phosphate buffer, pH 6.8  
 Flow Rate: 0.7 mL/min  
 Temperature: 25 °C  
 Injection: 10  $\mu$ L  
 Detection: UV 214 nm  
 Sample: mAb  
 Peaks: 1. Aggregate (M.W.=900 kDa)  
 2. MAb (M.W.=150 kDa)  
 3. Heavy Chain (M.W.=50 kDa)  
 4. Light Chain (M.W.=25 kDa)

For routine analysis, SEC columns packed with 5  $\mu$ m media are often used. As shown in Figure 2, a 5  $\mu$ m, 4.6x300 mm SEC-300 column provides suffice separation between the main peak of trastuzumab biosimilar and its aggregates but failing to resolve the fragment adjacent to the main mAb peak. When a 4.6x300 mm SEC-300 column packed with 3  $\mu$ m media is used, both the resolutions between the mAb and its aggregates and among its aggregates are improved significantly. Moreover, a fragment peak is partially resolved from the main mAb peak. The separation between the mAb and the fragment can be further improved on a 4.6x300 mm SEC-300 column packed with 1.8  $\mu$ m media.


Separation of MAb, Aggregates and Fragments (from HPLC to UHPLC)



Column: **BioCore SEC-300**, 5  $\mu$ m  
 Black: **BioCore SEC-300**, 5  $\mu$ m  
 Blue: **BioCore SEC-300**, 3  $\mu$ m  
 Red: **BioCore SEC-300**, 1.8  $\mu$ m  
 Dimension: 4.6 x 300 mm  
 Mobile Phase: 90/10 v/v 50 mM phosphate buffer, pH6.8/MeCN  
 Flow Rate: 0.25 mL/min  
 Temperature: 30 °C  
 Injection: 5  $\mu$ L  
 Detection: UV 280 nm  
 Sample: Trastuzumab Biosimilar (5 mg/mL)  
 Peaks: 1-3. Aggregates  
 4. MAb  
 5. Fragment

Particle Size	N(4)	Rs(1,2)	Rs(2,3)	Rs(3,4)	Rs(4,5)
1.8 $\mu$ m	22668	2.21	1.95	3.83	2.33
3 $\mu$ m	17006	2.84	1.88	3.58	1.06
5 $\mu$ m	9616	2.29	0.94	1.68	/


Very often, fast analysis is desired for better productivity. Figure 3 compares one 1.8  $\mu\text{m}$ , 4.6 $\times$ 150 mm BioCore SEC-300 with one 5  $\mu\text{m}$ , 4.6 $\times$ 300 mm BioCore SEC-300 column. While the resolutions between aggregates and mAb are comparable, the total analysis time on the 1.8  $\mu\text{m}$  4.6 $\times$ 150 mm BioCore SEC-300 (1.8  $\mu\text{m}$  is two times faster than that on the 5  $\mu\text{m}$ , 4.6 $\times$ 300 mm column, demonstrating that 1.8  $\mu\text{m}$  UHPLC SEC columns can not only provide superior resolution power, but also support fast analysis.



**Figure 3**

Column: **BioCore SEC-300**, 5  $\mu\text{m}$   
 Black: **BioCore SEC-300**, 5  $\mu\text{m}$   
 Blue: **BioCore SEC-300**, 1.8  $\mu\text{m}$   
 Dimension: Black: 4.6 $\times$ 300 mm; Blue: 4.6 $\times$ 150 mm  
 Mobile Phase: 300 mM NaCl in 50 mM phosphate buffer, pH6.8  
 Flow Rate: 0.35 mL/min  
 Temperature: 30  $^{\circ}\text{C}$   
 Injection: 5  $\mu\text{L}$   
 Detection: UV 280 nm  
 Sample: Trastuzumab Biosimilar (5 mg/mL)  
 Peaks: 1-3. Aggregates; 4. MAb

BioCore SEC-150 columns are based on 150  $\text{\AA}$  pore silica particles, designed for separating peptides and small proteins. Terlipressin is a peptide-like pro-drug to improve kidney function in adults. Terlipressin and impurity can be well separated on a BioCore SEC-150 column, as shown in Figure 5.




**Figure 5**

Column: **BioCore SEC-150**, 5  $\mu\text{m}$   
 Dimension: 7.8 $\times$ 300 mm  
 Mobile Phase: 100 mM  $\text{Na}_2\text{SO}_4$  in 100mM phosphate buffer  
 Flow Rate: 0.5 mL/min  
 Temperature: 30  $^{\circ}\text{C}$   
 Injection: 10  $\mu\text{L}$   
 Detection: UV 274 nm  
 Peaks: 1. Impurity U; 2. Terlipressin

R.T. (min)	Theoretical Plate (USP)	Tailing Factor (USP)	Resolution (USP)
23.080	29589	1.19	3.63

Compared with the BioCore SEC-300, BioCore SEC-500 columns employ larger pore ( $\sim$ 500  $\text{\AA}$ ) silica particles, thus are suitable for separating larger proteins. IgM is the largest immunoglobulin by size and exists mainly as pentamer in human serum. IgM polymer assembly depend on noncovalent interaction between subunit and disulfide bridging. Figure 4 shows the separation of IgM and its aggregates on a BioCore SEC-500 column.



**Figure 4**

Column: **BioCore SEC-500**, 3  $\mu\text{m}$   
 Dimension: 4.6 $\times$ 150 mm  
 Mobile Phase: 10/90 v/v MeCN/300 mM NaCl in 50 mM phosphate buffer, pH6.8  
 Flow Rate: 0.21 mL/min  
 Temperature: 30  $^{\circ}\text{C}$   
 Injection: 1  $\mu\text{L}$   
 Detection: UV 280 nm  
 Sample: IgM in Human Serum  
 Peak: 1. IgM

### Ordering Information

Particle Size ( $\mu\text{m}$ )	Column Dimension L $\times$ ID (mm)	Product Name		
		BioCore SEC-150	BioCore SEC-300	BioCore SEC-500
5	300 $\times$ 4.6	B213-050015-04630S	B213-050030-04630S	B213-050050-04630S
	150 $\times$ 4.6	B213-050015-04615S	B213-050030-04615S	B213-050050-04615S
	50 $\times$ 4.6	B213-050015-04605S	B213-050030-04605S	B213-050050-04605S
	300 $\times$ 7.8	B213-050015-07830S	B213-050030-07830S	B213-050050-07830S
	150 $\times$ 7.8	B213-050015-07815S	B213-050030-07815S	B213-050050-07815S
3	300 $\times$ 4.6	B213-030015-04630S	B213-030030-04630S	B213-030050-04630S
	150 $\times$ 4.6	B213-030015-04615S	B213-030030-04615S	B213-030050-04615S
	50 $\times$ 4.6	B213-030015-04605S	B213-030030-04605S	B213-030050-04605S
	300 $\times$ 7.8	B213-030015-07830S	B213-030030-07830S	B213-030050-07830S
	150 $\times$ 7.8	B213-030015-07815S	B213-030030-07815S	B213-030050-07815S
1.8	300 $\times$ 4.6	B213-018015-04630S	B213-018030-04630S	B213-018050-04630S
	150 $\times$ 4.6	B213-018015-04615S	B213-018030-04615S	B213-018050-04615S



# BioCore Ion-Exchange Columns

BioCore ion-exchange columns, including BioCore WCX, BioCore SCX, BioCore WAX, BioCoreSAX, are designed for separating charged variants in antibodies, proteins, related substances, and provide a broad application range in the areas including bio-tech, biopharmaceutical and academic research.

## Main Features

- Optimal selectivity for separating charged variants in antibodies and proteins
- Good peak shape and low carryover
- High column efficiency for better resolution power
- Excellent chemical and mechanical stability
- Good column-to-column consistency

## Specification

Product Name	BioCore WCX	BioCore SCX	BioCore WAX	BioCore SAX
Functional Group	Carboxylate	Sulfonate	Tertiary Amine	Quaternary Ammonium
Substrate	Monodispersed, spherical, nonporous PS/DVB particles			
Particle Size	5 & 10 µm			
Pore Size	Nonporous			
Pressure Limit	4500 psi for 10 µm 5000 psi for 5 µm			
Temperature Limit	60 °C			
pH Range	2-12			
Colum Dimension	4.6x250 mm 4.6x150 mm 4.6x50 mm			
Application	Charged variants in antibodies and proteins			

# Application




Charge variants are often considered as a critical quality attribute (CQA) in biopharmaceutical development and need to be closely monitored to ensure the safety and efficacy of biotherapeutics. Depending on the specific conditions of manufacturing or the intrinsic properties of a mAb, charge variants can contain both acidic or basic variants. Because additional charge variation may occur during process development and formulation, charge variant analysis is used across the biopharmaceutical pipeline, from discovery to manufacturing and quality control (QC). LC offers a robust, high-resolution and high-throughput analytical method, thus widely used in biopharmaceutical analysis.

IEX is a commonly used separation mode, which separates molecules based on the difference in the ionic groups on the protein surface. In IEX chromatography, elution is performed using either a salt- or pH-based gradient, which elutes the bound protein by disrupting the electrostatic interactions between the protein and the stationary phase. When traditional salt-based gradients have been employed, significant method development is required to get optimal resolution of charge variants contained in different mAb samples. pH-based gradients have gained popularity, primarily for their more generic applicability across different mAbs.

WCX is the most commonly used IEX mode for separating charged variants in high pI antibodies. Figures 6-8 illustrate three examples of charged variant separation on three mAbs (IgG1, IgG2 and IgG4) using BioCore WCX columns. Figure 9 provides a comparison of charged variant separation using salt-gradient and pH-gradient methods on a BioCore WCX column. Figure 10 shows the charged variant separation of a bispecific antibody on a BioCore WCX column. Figure 11 gives an example of charged variant separation of a pegylated Fab and its native Fab.

IgG1 on BioCore WCX




Column: **BioCore WCX**, 10 µm  
 Dimension: 4.6 × 250 mm  
 Mobile Phase: A) 20 mM MES, pH6.5  
 B) 150 mM NaCl in 20 mM MES, pH6.5  
 Gradient: 

t(min)	%A	%B
-15	95	5
0	95	5
0.1	95	5
40	80	20
40.1	0	100
43	0	100

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30 °C  
 Injection: 10 µL  
 Detection: UV 280 nm  
 Sample: IgG1 (~2.5 mg/mL in mobile phase A)

IgG2 on BioCore WCX




Column: **BioCore WCX**, 5 µm  
 Dimension: 4.6 × 150 mm  
 Mobile phase: A) 20 mM phosphate buffer, pH6.5  
 B) 300 mM NaCl in 20 mM phosphate buffer, pH6.5  
 Gradient: 

t (min)	%A	%B
-15	85	15
0	85	15
0.1	85	15
20	70	30
20.1	0	100
23	0	100

  
 Flow rate: 1.0 mL/min  
 Temperature: 30 °C  
 Injection: 25 µL  
 Detection: UV 280 nm  
 Sample: IgG2 (~1 mg/mL in mobile phase A)

## IgG4 on BioCore WCX




Column: **BioCore WCX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150mm  
 Mobile phase: A) 20 mM phosphate buffer, pH6.5  
 B) 300 mM NaCl in 20 mM phosphate buffer, pH6.5  
 Gradient: 

t(min)	%A	%B
-15	95	5
0	95	5
0.1	95	5
20	80	20
20.1	0	100
23	0	100

  
 Flow rate: 1 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 25  $\mu$ L  
 Detection: UV 280 nm  
 Sample: IgG4 (1mg/mL in H<sub>2</sub>O)

The BioCore SCX offers selectivity different from the BioCore WCX, thus can complement BioCore WCX for separating mAb charged variants (Figure 12). In addition, The charged variants of a fusion protein (Figure 13) can be resolved well on a BioCore SCX column.

## IgG4 on BioCore SCX




Column: **BioCore SCX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM MES, pH6.5  
 B) 300 mM NaCl in 20 mM MES, pH6.5  
 Gradient: 

t(min)	%A	%B
-15	93	7
0	93	7
20	79	21
20.1	0	100
23	0	100

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 10  $\mu$ L  
 Detection: UV 280 nm  
 Sample: IgG4 (1.0 mg/mL in H<sub>2</sub>O)

## Fusion Protein on BioCore SCX




Column: **BioCore SCX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM MES, pH6.0  
 B) 300 mM NaCl in 20 mM MES, pH6.0  
 Gradient: 

t(min)	%A	%B
-15	72	28
0	72	28
25	45	55
25.1	0	100
28	0	100

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 10  $\mu$ L  
 Detection: UV 280 nm  
 Sample: Fusion Protein (2.5 mg/mL in mobile phase A)

## Salt-Gradient vs pH-Gradient on BioCore WCX



Column: **BioCore WCX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  250 mm  
 Mobile Phase: A) 10 mM Tris, pH8.5  
 B) 100 mM NaCl in 10 mM Tris, pH8.5  
 Salt gradient: 


t (min)	%C	%D
-15	100	0
0	100	0
2	100	0
32	60	40
32.1	0	100
35	0	100

  
 Mobile Phase: A) pH gradient buffer A, pH5.0  
 B) pH gradient buffer B, pH10.8  
 pH gradient: 

t (min)	%A	%B
0	40	60
2	40	60
60	2	98
64	2	98
65	40	60
72	40	60

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 10  $\mu$ L  
 Detection: UV 280 nm  
 Sample: Trastuzumab Biosimilar (5 mg/mL)

## Bispecific Antibody on BioCore WCX




Column: **BioCore WCX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM phosphate buffer, pH6.5  
 B) 300 mM NaCl in 20 mM phosphate buffer, pH6.5  
 Gradient: 

t (min)	%A	%B
0	100	0
40	65	35
40.1	0	100
45	0	100
45.1	100	0
60	100	0

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 20  $\mu$ L  
 Detection: UV 280 nm  
 Sample: Bispecific Antibody (~5.0 mg/mL in mobile phase A)


BioCore anion-exchange columns are often employed to separate charged variants in low pI proteins, such as ovalbumin using BioCore WAX, (Figure 14) as well as a glycoprotein (Figure 15), a recombinant fusion protein (Figure 16), and a hGH-L-vFc protein (Figure 17) using the BioCore SAX.



Column: **BioCore WAX**, 5  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM Tris, pH8.0  
 B) 500mM NaCl in 20 mM Tris, pH8.0  
 Gradient: 

t(min)	%A	%B
0	1	0
0.1	1	0
15	50	50
15.1	99	100
20	99	100
20.1	1	0
30	1	0


  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 10  $\mu$ L  
 Detection: UV 280 nm  
 Sample: Ovalbumin (5 mg/mL)



Column: **BioCore SAX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  250 mm  
 Mobile Phase: A) 20 mM phosphate buffer, pH3.0  
 B) 300 mM NaCl in 20 mM phosphate buffer, pH3.0  
 Gradient: 

t(min)	%A	%B
-15	100	0
0	100	0
20	0	100
23	0	100


  
 Flow Rate: 1.0 mL/min  
 Injection: 5  $\mu$ L  
 Temperature: 30  $^{\circ}$ C  
 Detection: UV 280 nm  
 Sample: API (40 mg/mL)  
 Injection (10 mg/mL)



Column: **BioCore SAX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM MES, pH6.5  
 B) 300 mM NaCl in 20 mM MES, pH6.5  
 Gradient: 

t(min)	%A	%B
-15	70	30
0	70	30
20	40	60
20.1	0	100
23	0	100

  
 Flow Rate: 0.8 mL/min  
 Injection: 10  $\mu$ L  
 Temperature: 20  $^{\circ}$ C  
 Detection: UV 280 nm  
 Sample: Recombinant Fusion Protein (1 mg/mL in H<sub>2</sub>O)



Column: **BioCore SAX**, 10  $\mu$ m  
 Dimension: 4.6  $\times$  150 mm  
 Mobile Phase: A) 20 mM phosphate buffer, pH6.5  
 B) 300 mM NaCl in 20 mM phosphate buffer, pH6.5  
 Gradient: 

t (min)	%A	%B
0	80	20
20	50	50
20.1	0	100
25	0	100
25.1	80	20
35	80	20

  
 Flow Rate: 1.0 mL/min  
 Temperature: 30  $^{\circ}$ C  
 Injection: 5  $\mu$ L  
 Detection: UV 220 nm  
 Sample: HGH-L-vFc Protein (pI=5.9-6.2, 8 mg/mL)

## Ordering Information

Particle Size ( $\mu$ m)	Column Dimension L $\times$ ID (mm)	Product Name			
		BioCore WCX	BioCore SCX	BioCore WAX	BioCore SAX
10	250 $\times$ 4.6	B311-100000-04625P	B411-100000-04625P	B511-100000-04625P	B611-100000-04625P
	150 $\times$ 4.6	B311-100000-04615P	B411-100000-04615P	B511-100000-04615P	B611-100000-04615P
	50 $\times$ 4.6	B311-100000-04605P	B411-100000-04605P	B511-100000-04605P	B611-100000-04605P
5	250 $\times$ 4.6	B311-050000-04625P	B411-050000-04625P	B511-050000-04625P	B611-050000-04625P
	150 $\times$ 4.6	B311-050000-04615P	B411-050000-04615P	B511-050000-04615P	B611-050000-04615P
	50 $\times$ 4.6	B311-050000-04605P	B411-050000-04605P	B511-050000-04605P	B611-050000-04605P

# BioCore HIC-Butyl

BioCore HIC-Butyl is a family of high-performance, hydrophobic interaction chromatography columns that separate monoclonal antibodies (mAbs) and antibody-drug conjugates (ADCs) based on their differences in surface hydrophobicity.

## Main Features

- Optimal selectivity for the DAR analysis in ADCs
- Minimal undesired interactions for low carryover
- Excellent mechanical strength for column robustness
- Good column-to-column consistency


## Specification

Product Name	BioCore HIC-Butyl
Functional Group	Butyl
Substrate	Monodispersed, wide-pore silica particles
Particle Size	5 $\mu\text{m}$
Pore Size	1000 $\text{\AA}$
Pressure Limit	6000 psi
Temperature Limit	60 $^{\circ}\text{C}$
pH Range	2-8
Colum Dimension	4.6x250 mm 4.6x100 mm 4.6x50 mm
Application	DAR analysis in ADCs, oxidation variants in mAbs

# Application



Hydrophobic interaction chromatography (HIC) is a traditional technique used for the separation, purification, and characterization of proteins. As the number of antibody-drug conjugates (ADCs) continues to increase in development of bio-therapeutics, HIC and RP methods, both utilizing changes in hydrophobicity for separation, are often used for ADC characterization and analysis. Unlike RP technique, HIC uniquely allows for protein analysis under mild nondenaturing conditions that preserve the native structure and activity of the molecules, thus analysis of the ADC in its native form is advantageous. Figure 18 presents a HIC method for the drug antibody ratio (DAR) analysis of a cysteine-conjugated antibody using an ammonium sulfate buffer system on a BioCore HIC-Butyl column. Figure 19 demonstrates the suitability of the BioCore HIC-Butyl for separating variants in a bispecific antibody.



## Ordering Information

Product Name	Particle Size ( $\mu\text{m}$ )	Column Dimension L x ID (mm)	Part Number
BioCore HIC-Butyl	5	250x4.6	B713-050100-04625S
		100x4.6	B713-050100-04610S
		50x4.6	B713-050100-04605S



# BioCore RP-Butyl

BioCore RP-Butyl is a family of high-performance, reversed-phase columns, designed for the separation and determination of exact mass of antibodies / proteins and related substances.

## Main Features

- High column efficiency and low carryover
- Excellent mechanical strength
- Good MS compatibility
- Good column-to-column consistency

## Specification

Product Name	BioCore RP-Butyl
Functional Group	Butyl
Substrate	Monodispersed, spherical PS/DVB particles
Particle Size	5 $\mu\text{m}$
Pore Size	Nonporous
Pressure Limit	4500 psi
Temperature Limit	100 $^{\circ}\text{C}$
pH Range	2-12
Column Dimension	3.0 $\times$ 100 mm 3.0 $\times$ 50 mm
Application	Intact protein and protein fragment analysis

## Ordering Information

Product Name	Particle Size ( $\mu\text{m}$ )	Column Dimension L $\times$ ID (mm)	Part Number
BioCore RP-Butyl	5	150 $\times$ 4.6	B821-050000-04615S
		100 $\times$ 4.6	B821-050000-04610S
		50 $\times$ 4.6	B821-050000-04605S
		150 $\times$ 3.0	B821-050000-03015S
		100 $\times$ 3.0	B821-050000-03010S
		50 $\times$ 3.0	B821-050000-03005S
		150 $\times$ 2.1	B821-050000-02115S
		100 $\times$ 2.1	B821-050000-02110S
		50 $\times$ 2.1	B821-050000-02105S

# Application



Reversed-phase HPLC is one of most important techniques for protein separations and the method of choice for peptide separation. In protein characterization, RP-HPLC is commonly used for separating variants in proteins or antibodies by their difference in hydrophobicity. When combined with high resolution mass spectrometry (MS), RP-HPLC can be utilized for determination of the exact mass for the whole proteins or their subunits (e.g., heavy chain, light chain, Fab, Fc). Figures 20 and 21 provide examples of exact mass determination by LC-MS using a BioCore RP-Butyl column.

Exact Mass of a Monoclonal Antibody






Figure 20

Column: **BioCore RP-Butyl**, 5  $\mu\text{m}$   
 Dimension: 3.0  $\times$  50 mm  
 Mobile Phase: A) 0.1% HCOOH in H<sub>2</sub>O  
 B) 0.1% HCOOH in MeCN  
 Gradient: 

t(min)	%A	%B
0	95	5
1	95	5
1.1	95	5
1.2	5	95
4	5	95
4.1	95	5
7	95	5

  
 Flow Rate: 0.6 mL/min  
 Temperature: 60  $^{\circ}\text{C}$   
 Injection:  $\sim$ 1  $\mu\text{g}$   
 Detection: QE  
 Sample: Monoclonal Antibody

Exact Mass of a Recombinant Protein






Figure 21

Column: **BioCore RP-Butyl**, 5  $\mu\text{m}$   
 Dimension: 3.0  $\times$  50 mm  
 Mobile Phase: A) 0.1% HCOOH in H<sub>2</sub>O  
 B) 0.1% HCOOH in MeCN  
 Gradient: 

t(min)	%A	%B
0	95	5
1	95	5
1.1	95	5
1.2	5	95
4	5	95
4.1	95	5
7	95	5

  
 Flow Rate: 0.6 mL/min  
 Temperature: 60  $^{\circ}\text{C}$   
 Injection:  $\sim$ 1  $\mu\text{g}$   
 Detection: QE  
 Sample: Recombinant Protein

# BioCore Glycan

BioCore Glycan is a family of high-performance HILIC columns designed for profiling N-glycans present in proteins/antibodies and related substances.

## Main Features

- Desired selectivity for separating fluorescently labeled N-glycans in proteins
- High resolution and stability
- Good MS compatibility
- Good column-to-column consistency

## Specification

Product Name	BioCore Glycan
Functional Group	Polyamide
Substrate	Monodispersed, spherical silica particles
Particle Size	3 µm
Pore Size	180 Å
Pressure Limit	6000 psi
Temperature Limit	80 °C
pH Range	2-9
Column Dimension	2.1×150 mm 2.1×100 mm 3.0×150 mm 3.0×100 mm
Application	N-glycans analysis in antibodies and proteins

## Ordering Information

Product Name	Particle Size (µm)	Column Dimension L x ID (mm)	Part Number
BioCore Glycan	3	150×2.1	B913-030018-02115S
		100×2.1	B913-030018-02110S
		150×3.0	B913-030018-03015S
		100×3.0	B913-030018-03010S
		150×4.6	B913-030018-04615S
		100×4.6	B913-030018-04610S


# Application



Biopharmaceuticals are highly complex molecules with remarkable heterogeneity. Protein glycosylation is an inherent source of this heterogeneity, which affects the safety, efficacy, and half-life of therapeutic glycoproteins. Therefore, analysis of the glycan pattern is an important issue for characterization and quality control in the biopharmaceutical industry.

For the analysis of protein N-glycans, the sample-preparation procedure consists of the release of the N-glycans by PNGase-F, followed by fluorescence labeling and removal of excess label. Subsequently, labeled glycans are usually analyzed by hydrophilic-interaction liquid chromatography (HILIC) with a fluorescence detector.


Figures 22 and 23 demonstrate the N-glycan profiling with 2-AB labeling of Avastin and a recombinant protein obtained from a BioCore Glycan column.



Column: **BioCore Glycan**, 3 µm  
 Dimension: 3.0 × 100 mm  
 Mobile Phase: A) 50 mM CH<sub>3</sub>COONH<sub>4</sub> in H<sub>2</sub>O, pH4.4  
 B) MeCN  
 Gradient:
 

t (min)	%A	%B
0	25	75
2	25	75
37	46	54
37.1	25	75
40	25	75

 Flow Rate: 0.4 mL/min  
 Temperature: 65 °C  
 Injection: 30 µL  
 Detection: Ex/Em= 250 nm/428 nm  
 Sample: N-Glycans of Avastin



Column: **BioCore Glycan**, 3 µm  
 Dimension: 3.0 × 100 mm  
 Mobile Phase: A) 50 mM CH<sub>3</sub>COONH<sub>4</sub> in H<sub>2</sub>O, pH4.4  
 B) MeCN  
 Gradient:
 

t (min)	%A	%B
0	25	75
2	25	75
37	46	54
37.1	25	75
40	25	75

 Flow Rate: 0.4 mL/min  
 Temperature: 65 °C  
 Injection: 30 µL  
 Detection: Ex/Em= 250 nm/428 nm  
 Sample: N-Glycans of a Recombinant Protein

# BioCore Protein A

BioCore Protein A is a family of high-performance affinity chromatography columns, designed for fast titer analysis of monoclonal antibodies (mAb) and Fc fusion proteins.

## Main Features

- High specificity for efficient capture of a broad selection of antibodies and antibody fragments
- High dynamic binding capacity for a wide linear range
- High mechanical strength for faster analysis and better column lifetime
- Low ligand leakage for higher purity

## Specification


Product Name	BioCore Protein A
Functional Group	pH stable rProtein A
Substrate	Monodispersed, porous, spherical PS/DVB particles
Particle Size	15 µm
Pore Size	1000 Å
Pressure Limit	1450 psi
Temperature Limit	2-40 °C
pH Range	2-12
Dynamic Binding Capacity	≥20 mg/mL (IgG)
Linear Range (≥0.99)	0-200 µg (2.1X30 mm)

# Application



In early development of recombinant monoclonal antibodies (mAbs) for biotherapeutics, a large number of harvest cell culture (HCC) samples need to be screened for immunoglobulin G (IgG) titer. Affinity chromatography employing protein A ligand is often used for the mAb concentration determination, as well as for downstream purification.


Titer analysis is required for accurate determination of amount of mAbs in the sample. Figure 24 shows a 20 µL injection of an antibody sample onto the BioCore Protein A column. The unbound material elutes first (large peak) at pH7.0. Then, the antibody is released using a low pH wash (pH 2.5). The BioCore Protein A column gives a sharp peak for the mAb with fast, efficient elution. The complete cycle time, including equilibration, is 2.5 minutes. The mAb titer is determined by back calculating the integrated IgG peak area against a previously generated calibration curve (see Figure 25). The BioCore Protein A provides accurate titer analysis with linearity over a wide concentration range (0.2 mg/mL to 6.25 mg/mL).



Column: **BioCore Protein A**, 15 µm  
 Dimension: 2.1×30 mm  
 Mobile Phase: 150 mM NaCl in 50 mM phosphate buffer, pH7.0  
 150 mM NaCl in 50 mM phosphate buffer, pH2.5  
 Gradient:

t(min)	%A	%B
0	100	0
0.5	100	0
0.51	0	100
1.5	0	100
1.51	100	0
2.5	100	0

Flow Rate: 2.0 mL/min  
 Temperature: 30 °C  
 Injection: 20 µL  
 Detection: UV 280 nm  
 Sample: Harvest Cell Culture (IgG~2.4 mg/mL)  
 Peaks: 1. Unbound Components  
 2. IgG




Column: **BioCore Protein A**, 15 µm  
 Dimension: 2.1×30 mm  
 Mobile Phase: A) 150 mM NaCl in 50 mM phosphate buffer, pH 7.0  
 B) 150 mM NaCl in 50 mM phosphate buffer, pH 2.5  
 Gradient:

t (min)	%A	%B
0	100	0
0.5	100	0
0.51	0	100
1.5	0	100
1.51	100	0
2.5	100	0

Flow Rate: 2.0 mL/min  
 Temperature: 30 °C  
 Injection: 20 µL  
 Detection: UV 280 nm  
 Sample: IgG (0.2-6.25 mg/mL in H<sub>2</sub>O)

The BioCore Protein A column illustrates satisfactory ruggedness when tested continuously for 2,000 cycles and every hundred cycles a set of calibration standards (from 0.2 mg/mL to 6.25 mg/mL) are analyzed. As shown in Figure 26, the retention time, peak area, and peak width of IgG remain virtually unchanged.



Column: **BioCore Protein A**, 15  $\mu\text{m}$   
 Dimension: 2.1  $\times$  30 mm  
 Mobile Phase : 150 mM NaCl in 50 mM phosphate buffer, pH7.0  
 150 mM NaCl in 50 mM phosphate buffer, pH2.5  
 Gradient:

t(min)	%A	%B
0	100	0
0.5	100	0
0.51	0	100
1.5	0	100
1.51	100	0
2.5	100	0

Flow Rate: 2.0 mL/min  
 Temperature: 30  $^{\circ}\text{C}$   
 Injection: 20  $\mu\text{L}$   
 Detection: UV 280 nm  
 Sample: IgG (1.0 mg/mL in  $\text{H}_2\text{O}$ )

## Ordering Information

Product Name	Particle Size ( $\mu\text{m}$ )	Column Dimension L x ID (mm)	Part Number
BioCore Protein A	15	100 $\times$ 2.1	B111-150100-02110S
		50 $\times$ 2.1	B111-150100-02105S
		30 $\times$ 2.1	B111-150100-02103S
		100 $\times$ 4.6	B111-150100-04610S
		50 $\times$ 4.6	B111-150100-04605S
		30 $\times$ 4.6	B111-150100-04603S