

The moment you adopt our Kromasil High Performance Concept, you join thousands of chromatographers who share a common goal: to achieve better separations when analyzing or isolating pharmaceuticals or other substances.

Not only will you benefit from our patented silica technology, but you gain a strong partner with a reliable track record in the field of silica products. For the past 60 years, Eka Chemicals has pioneered new types of silica. Our long experience in the field of silica chemistry is the secret behind the development of Kromasil, and the success of our Separation Products Group.


Kromasil is available in bulk, or in high-pressure slurry-packed columns. The development, production and marketing of Kromasil are ISO 9001 certified.

Eka Chemicals is a global company with 3,000 people in 30 countries. It is a business unit within AkzoNobel, one of the world's largest chemical groups, with more than 67,000 employees in 80 countries.


PUBBankv_1105

© Eka Chemicals AB 2011 This publication may not be reproduced in any way without the consent of Eka Chemicals AB.

Polysaccharide-based chiral HPLC phases for better performance



Kromasil[®] AmyCoat[™]
Kromasil[®] CelluCoat[™]




AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0, Fax +49 6131 880 96-20
e-mail: info@mz-at.de, www.mz-at.de

www.kromasil.com

Kromasil[®]



MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

AUTHORIZED DISTRIBUTOR

Designed to stretch the limits.

Kromasil® AmyCoat™ and Kromasil® CelluCoat™ stretch the limits for analytical and industrial chiral chromatography. The silica is based on an in-house developed matrix and coated with a functionalized amylose or cellulose selector. High resolution, excellent selectivity, no pressure limits* and stable performance when switching between compatible mobile phases are some important benefits.

ANALYTICAL CHROMATOGRAPHY

NO PRESSURE LIMITS ▶ HIGH FLOW RATES ▶ FASTER ANALYSES

SMALL PARTICLES (3 μm) ▶ HIGH RESOLUTION ▶ BETTER ANALYSES

UNIQUE COATING METHOD ▶ EXCELLENT SELECTIVITY ▶ BETTER ANALYSES

PREPARATIVE CHROMATOGRAPHY

UNIQUE SILICA MATRIX ▶ MECHANICALLY STRONG ▶ LONGER LIFETIME

WIDE RANGE OF PARTICLES ▶ CONSISTENT PERFORMANCE ▶ EASIER TO SCALE UP

UNIQUE COATING METHOD ▶ HIGHER LOADABILITY ▶ HIGHER PRODUCTIVITY

* Kromasil AmyCoat and Kromasil CelluCoat can withstand flow rates equivalent to pressures of up to 400 bar—i.e. the approximate limit for most HPLC systems.



Peak performance. Expect nothing less.

To speed up and simplify method development, we have removed some of the restrictions for coated polysaccharide phases. In analytical scale chromatography, 3 µm particles and the absence of pressure limits allow faster chromatography with better results than ever before.

BETTER RESULTS

Both Kromasil AmyCoat and Kromasil CelluCoat show great enantioselectivity for many different racemates. In a separate section of the brochure, you will find an application guide showing exactly the performance levels you can expect.

By having access to 3 µm particles, it enables a higher plate count and resolution for analytical chromatography. Combined with great selectivity, this facilitates the separation of the enantiomers.


FASTER RESULTS

Not only can you get better results with Kromasil AmyCoat and Kromasil CelluCoat you also get them faster. Thanks to a product characteristic like the absence of pressure limits, you can run analytical chromatography at very high flow rates and save time.


The high selectivity and 3 µm particles are also important factors in speeding up the method. You achieve an improved resolution, which, in turn, allows a higher flow rate.

Kromasil AmyCoat

Selectivity and resolution comparison—Kromasil AmyCoat 3 µm and Chiralpak AD-H (5 µm).




Difference in resolution—Kromasil AmyCoat 3 µm vs. 5 µm.



Baseline separation by using 3 µm particles


Fast analytical chromatography




Baseline separation in less than 1 minute

Kromasil CelluCoat

Selectivity and resolution comparison—Kromasil CelluCoat 3 µm and Chiralcel OD-H (5 µm).



Fast analytical chromatography



Baseline separation in less than 50 seconds



Kromasil AmyCoat and Kromasil CelluCoat are available in particle sizes down to 3 µm.



Kromasil AmyCoat and Kromasil CelluCoat —faster results, better results.



Makes everyday work so much easier.

Kromasil AmyCoat and Kromasil CelluCoat let you perform method development without interference from restrictive parameters. Forget pressure limits, long equilibration times and degradation in performance when switching from one compatible mobile phase to another—they simply do not exist anymore.

NO PRESSURE LIMITS

User-friendliness, as we see it, is the ability to perform method development without being restricted by various parameters. One well-known restriction for coated polysaccharide phases is the general pressure limit over the bed. Kromasil AmyCoat and Kromasil CelluCoat withstand flow rates equivalent to pressures up to 400 bar—which is about the limit for the HPLC system itself. This makes it possible to run analytical chromatography a lot faster.

STABLE PERFORMANCE

From a stability perspective, you can run Kromasil AmyCoat and Kromasil CelluCoat using compatible normal, polar organic and reversed mobile phases. Switching between compatible normal to polar organic mobile phases will not lead to any degradation in performance, no need for solvent dedicated columns.

To show the stability in column efficiency for Kromasil CelluCoat, we have performed mobile phase switches and visualized them after each switch with a Knox plot.


SHORT EQUILIBRATION TIMES

Column equilibration is a time-consuming activity when running chiral chromatography. Higher flow rates allow Kromasil AmyCoat and Kromasil CelluCoat to show substantially shorter equilibration times than other coated phases on the market.


In general, long equilibration times are most pronounced when switching mobile phases containing basic additives to acidic additives or the other way around. The test with a Kromasil CelluCoat 3 μm column switching between two compatible mobile phases shows how short equilibration times you can get. There are also two tests illustrating short equilibration times and additive switch for both Kromasil AmyCoat and Kromasil CelluCoat with absolutely no sign of memory effects.

Kromasil AmyCoat

Stable performance—No pressure limits—Freedom to switch solvents




Short equilibration times—Freedom to switch additives




Kromasil CelluCoat


Stable performance—No pressure limits—Freedom to switch solvents



Short equilibration times—Freedom to switch solvents



Short equilibration times—Freedom to switch additives



Works all the way.

Kromasil products are well known for their ability to work along the whole spectrum from analytical to industrial scale chromatography. Kromasil AmyCoat and Kromasil CelluCoat are no exception.


SIMPLIFIES METHOD DEVELOPMENT

With particle sizes from 3 μm to 10 μm giving identical selectivity, Kromasil AmyCoat and Kromasil CelluCoat make it easy to scale up while retaining excellent performance.

As for all Kromasil products, you can perform the required method development in analytical scale columns and then scale up to a larger column.

You can, for example, use 3 μm particles in an analytical scale column and translate all the generated data to a larger column packed with 10 μm particles.

If you plan to scale up the process, we suggest you use an analytical column packed with 10 μm particles right from the start.




N.B. The relative back-pressure increase is a measure of the degree of degradation of the material after repeated packings. Actual particle size for Chiralcel and Chiralpak is about 3 times larger than that for Kromasil, which explains the difference in initial back pressure (back-pressure is inversely proportional to the square of the particle size).

The test was designed to exert greater than normal mechanical stress on the chiral stationary phases, and is performed at a packing pressure above the maximum 50 bar recommended by the manufacturer of Chiralcel OD and Chiralpak AD.

FULLY BACK-INTEGRATED

We manufacture the super wide pore silica for Kromasil polysaccharide products and perform all subsequent steps leading to the final product.

This means you get a product with full traceability. Every manufacturing step is ensured through our detailed quality system, and the final product is never released until it has passed our rigorous quality control tests. Having all the manufacturing steps in-house also gives us a strong position in terms of production capacity.



• Kromasil 3-CelluCoat $\alpha = 2.4$ $R_s = 13.9$ • Kromasil 5-CelluCoat $\alpha = 2.4$ $R_s = 11.4$ • Kromasil 10-CelluCoat $\alpha = 2.4$ $R_s = 10.0$


Common conditions
Column size: 4.6 x 150 mm
Mobile phase: Heptane/2-Propanol (90/10)
Flow rate: 0.5 ml/min
Temperature: 25 °C
Solute: trans-Stilbene oxide
Detection: UV 229 nm

MECHANICALLY STRONG


Mechanical strength is an important product lifetime parameter. Kromasil AmyCoat and Kromasil CelluCoat have a mechanically strong spherical silica which withstands repeated cycles of packing.

Kromasil AmyCoat

Preparative applications




EASY TO SCALE UP




Kromasil CelluCoat

Loadability comparison



Preparative applications




Taking a closer look.


Kromasil AmyCoat

CHIRAL SELECTOR

The coated selector is tris-(3,5-dimethylphenyl) carbamoyl amylose.




PARTICLE SIZES



Kromasil AmyCoat is available in 3 and 5 µm as pre-packed columns, and in 10 µm as pre-packed column and bulk material.

COMPATIBLE MOBILE PHASES



Before testing conditions not mentioned, please contact the Kromasil application group, kromasil@akzonobel.com

PRODUCT ASSORTMENT*

	3 µm - AmyCoat	5 µm - AmyCoat	10 µm - AmyCoat
4.6 x 50 mm	Kromasil 3-AmyCoat 4.6 x 50 mm	Kromasil 5-AmyCoat 4.6 x 50 mm	—
4.6 x 150 mm	Kromasil 3-AmyCoat 4.6 x 150 mm	Kromasil 5-AmyCoat 4.6 x 150 mm	Kromasil 10-AmyCoat 4.6 x 150 mm
4.6 x 250 mm	Kromasil 3-AmyCoat 4.6 x 250 mm	Kromasil 5-AmyCoat 4.6 x 250 mm	Kromasil 10-AmyCoat 4.6 x 250 mm
10 x 250 mm	—	Kromasil-5-AmyCoat 10 x 250 mm	Kromasil-10-AmyCoat 10 x 250 mm
21.2 x 250 mm	—	Kromasil 5-AmyCoat 21.2 x 250 mm	Kromasil 10-AmyCoat 21.2 x 250 mm

* Other column dimensions available upon request.


PRODUCT CODES

Kromasil X-AmyCoat followed by column diameter and length. X indicates particle size: 3, 5 or 10 µm. For example: *Kromasil 3-AmyCoat 4.6 x 150 mm*


Kromasil CelluCoat

CHIRAL SELECTOR

The coated selector is tris-(3,5-dimethylphenyl) carbamoyl cellulose.




PARTICLE SIZES



Kromasil CelluCoat is available in 3 and 5 µm as pre-packed columns, and in 10 µm as pre-packed column and bulk material.

COMPATIBLE MOBILE PHASES



Before testing conditions not mentioned, please contact the Kromasil application group, kromasil@akzonobel.com

PRODUCT ASSORTMENT*

	3 µm - CelluCoat	5 µm - CelluCoat	10 µm - CelluCoat
4.6 x 50 mm	Kromasil 3-CelluCoat 4.6 x 50 mm	Kromasil 5-CelluCoat 4.6 x 50 mm	—
4.6 x 150 mm	Kromasil 3-CelluCoat 4.6 x 150 mm	Kromasil 5-CelluCoat 4.6 x 150 mm	Kromasil 10-CelluCoat 4.6 x 150 mm
4.6 x 250 mm	Kromasil 3-CelluCoat 4.6 x 250 mm	Kromasil 5-CelluCoat 4.6 x 250 mm	Kromasil 10-CelluCoat 4.6 x 250 mm
10 x 250 mm	—	Kromasil-5-CelluCoat 10 x 250 mm	Kromasil-10-CelluCoat 10 x 250 mm
21.2 x 250 mm	—	Kromasil 5-CelluCoat 21.2 x 250 mm	Kromasil 10-CelluCoat 21.2 x 250 mm

* Other column dimensions available upon request.

PRODUCT CODES

Kromasil X-CelluCoat followed by column diameter and length. X indicates particle size: 3, 5 or 10 µm. For example: *Kromasil 3-CelluCoat 4.6 x 150 mm*