


**Analysis of MAb with Charge Variance on Antibodix™ NP10 (10 μm, 4.6x250 mm)**

Mobile Phases:	A, 10 mM phosphate buffer, pH 7.5; B, A + 100 mM NaCl
Gradient:	15-55% B (30 min)
Flow Rate:	0.8 mL/min
Detection:	UV 214 nm
Temperature:	Ambient
Concentration:	5 mg/mL
Injection Volume:	5 μL
Sample:	Monoclonal antibody

Keywords: Weak cation exchange, Antibodix, monoclonal antibody, charge variance