




Analysis of MAb-X22 on Antibodix™ NP10 (10 µm, 4.6x250 mm) - Impact of Mobile Phase pH



Column: Antibodix™-NP10 (10 µm, 4.6x250 mm)
Mobile Phases: A, 10 mM phosphate buffer; B, A + 100 mM NaCl
PH: A, 7.5 B, 7.0
Gradients: 15-65%B in 60 min
Flow Rate: 0.8 mL/min
Detection: UV 214 nm
Temperature: 25 °C
Concentration: 1.5 mg/mL
Injection Volume: 10 µL
Sample: MAb-X22

Keywords: Weak cation exchange, Antibodix, monoclonal antibody, mobile phase pH