

# High Resolution and Fast Analysis of N-linked Sialylated Glycan by Sepax Proteomix SAX

Proteomix SAX-NP3, 3  $\mu\text{m}$ , 4.6 x 50 mm (PN: 403NP3-4605)

Proteomix SAX-NP5, 5  $\mu\text{m}$ , 2.1 x 250 mm PEEK (PN: 403NP5P-2125)



# 2-AB-labeled Fetuin N-linked Sialylated Glycan Separation

## Comparison - Proteomix SAX - 3 $\mu\text{m}$ , 4.6 x 50 mm vs. 5 $\mu\text{m}$ , 2.1 x 250 mm


**Column:** Proteomix SAX-NP3, 3  $\mu\text{m}$ , 4.6 x 50 mm (PN: 403NP3-4605), UPLC, 430 bar, Flow Rate: 1 mL/min

Proteomix SAX-NP5, 5  $\mu\text{m}$ , 2.1 x 250 mm PEEK (PN: 403NP5P-2125), HPLC, 225 bar, Flow Rate: 0.2 mL/min

**Mobile phase A:** H<sub>2</sub>O + 20% ACN, B: 250 mM Ammonium Formate pH 4.5 + 20% ACN

**Detector:** Excitation 260 nm, Emission 430 nm; PMT 10 **Column temperature:** Ambient

**Injection Volume:** 3.0  $\mu\text{L}$  and 1.5  $\mu\text{L}$  **Sample:** N-Linked Glycans Bovine Fetuin



Zoom-in view presented in following slides

[www.sepax-tech.com](http://www.sepax-tech.com)



# 2-AB-labeled Fetuin N-linked Sialylated Glycan Separation on Proteomix


## SAX 5 $\mu\text{m}$ , 2.1 x 250 mm - High Resolution

**Column:** Proteomix SAX-NP5, 2.1 x 250mm PEEK (403NP5P-2125)

**Mobile phase A:** H<sub>2</sub>O + 20% ACN, B: 250 mM Ammonium Formate pH 4.5 + 20% ACN

**Flow Rate** 0.2 mL/min; **Detector:** Excitation 260 nm, Emission 430 nm; PMT 10; **Column temperature:** Ambient

**Injection Volume:** 1.5  $\mu\text{L}$  **Sample:** N-Linked Glycans Bovine Fetuin



# 2-AB-labeled Fetuin N-linked Sialylated Glycan Separation on Proteomix SAX, 3 $\mu\text{m}$ , UPLC – FASTER Analysis

**Column:** Proteomix SAX-NP3, 4.6 x 50mm (PN: 403NP3-4605)

**Mobile phase A:** H<sub>2</sub>O + 20% ACN, B: 250 mM Ammonium Formate pH 4.5 + 20% ACN

**Flow Rate:** 1.0 mL/min; **Detector:** Excitation 260 nm, Emission 430 nm; PMT 10 **Column temperature:** Ambient

**Injection Volume:** 3.0  $\mu\text{l}$ ; **Sample:** N-Linked Glycans Bovine Fetuin (6.7 pmol/ $\mu\text{L}$ )

