

Polymer-based Hydrophilic Interaction Chromatography (HILIC) Columns (HILICpak)

Features

- VG-50**
- Suitable for saccharides analysis using HILIC mode
 - Recovers reducing saccharides with high ratio
 - Polymer-based packing material provides excellent chemical stability and minimum deterioration over extended time period
 - Easily regenerated by washing in an alkaline solution
 - Appropriate for evaporative light scattering detector, corona charged aerosol detector, and LC/MS

- VT-50**
- Suitable for anionic substances analysis using HILIC mode
 - Use of some eluents add ion exchange mode
 - Polymer-based packing material provides excellent chemical stability and minimum deterioration over extended time period
 - Suitable for LC/MS analysis


- New VC-50**
- Modified carboxyl group is suitable for cationic substance analysis including amines
 - The dominant separation mode is RP or IEX rather than HILIC mode

- New VN-50**
- The modified diol groups on the packing material create the HILIC mode
 - Suitable for oligosaccharide separation which is not possible by SEC column or conventional HILIC columns

Standard columns

VG-50

(Housing Material: SUS)



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17• D-55129 Mainz
Tel +49 6131 880 96-0, Fax +49 6131 880 96-20
e-mail: info@mz-at.de, www.mz-at.de

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630200	HILICpak VG-50 4D	≥ 5,500	Amino	5	100	4.6 x 150	H ₂ O/CH ₃ CN=20/80
F7630100	HILICpak VG-50 4E	≥ 7,500	Amino	5	100	4.6 x 250	H ₂ O/CH ₃ CN=20/80
F6711100	HILICpak VG-50G 4A	(guard column)	Amino	5	100	4.6 x 10	H ₂ O/CH ₃ CN=20/80

Base Material: Polyvinyl alcohol

(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630300	HILICpak VG-50 2D	≥ 3,500	Amino	5	100	2.0 x 150	H ₂ O/CH ₃ CN=15/85
F6711200	HILICpak VG-50G 2A	(guard column)	Amino	5	100	2.0 x 10	H ₂ O/CH ₃ CN=15/85

Base Material: Polyvinyl alcohol

VT-50

(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630400	HILICpak VT-50 2D	≥ 4,500	Quaternary ammonium	5	100	2.0 x 150	25mM HCOONH ₄ aq. /CH ₃ CN=15/85
F6711300	HILICpak VT-50G 2A	(guard column)	Quaternary ammonium	5	100	2.0 x 10	25mM HCOONH ₄ aq. /CH ₃ CN=15/85

Base Material: Polyvinyl alcohol

VC-50

(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630700	New HILICpak VC-50 2D	≥ 3,500	Carboxyl	5	100	2.0 x 150	H ₂ O
F6711600	New HILICpak VC-50G 2A	(guard column)	Carboxyl	5	100	2.0 x 10	H ₂ O






Base Material: Polyvinyl alcohol

VN-50

(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630600	New HILICpak VN-50 2D	≥ 3,500	Diol	5	100	2.0 x 150	H ₂ O/CH ₃ CN=25/75
F6711500	New HILICpak VN-50G 2A	(guard column)	Diol	5	100	2.0 x 10	H ₂ O/CH ₃ CN=25/75
F7630500	New HILICpak VN-50 4D	≥ 10,000	Diol	5	100	4.6 x 150	H ₂ O/CH ₃ CN=25/75
F6711400	New HILICpak VN-50G 4A	(guard column)	Diol	5	100	4.6 x 10	H ₂ O/CH ₃ CN=25/75

Base Material: Polyvinyl alcohol



● Polymer-based Hydrophilic Interaction Chromatography (HILIC) Columns (Asahipak)

Features

- NH2P**
- Suitable for saccharides analysis using HILIC mode
 - Polymer-based packing material provides excellent chemical stability and minimum deterioration over extended time period
 - Easily regenerated by washing in an alkaline solution
 - Appropriate for evaporative light scattering detector, corona charged aerosol detector, and LC/MS
 - Fulfills USP L82 requirements

- NH2P-40**
- Provides higher theoretical plate number than NH2P-50 series

■ Standard columns

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (μm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630005	Asahipak NH2P-50 4B	≥ 1,500	Amino	5	100	4.6 x 50	CH ₃ CN
F7630002	Asahipak NH2P-50 4D	≥ 5,500	Amino	5	100	4.6 x 150	CH ₃ CN
F7630001	Asahipak NH2P-50 4E	≥ 7,500	Amino	5	100	4.6 x 250	CH ₃ CN
F6710016	Asahipak NH2P-50G 4A	(guard column)	Amino	5	–	4.6 x 10	CH ₃ CN
F7630006	Asahipak NH2P-50 2D	≥ 3,500	Amino	5	100	2.0 x 150	CH ₃ CN
F6713000	Asahipak NH2P-50G 2A	(guard column)	Amino	5	–	2.0 x 10	CH ₃ CN
F7630007	Asahipak NH2P-40 3E	≥ 8,500	Amino	4	100	3.0 x 250	CH ₃ CN
F6710030	Asahipak NH2P-50G 3A	(guard column)	Amino	5	–	3.0 x 10	CH ₃ CN
F7630008	Asahipak NH2P-40 2B	≥ 2,000	Amino	4	100	2.0 x 50	CH ₃ CN
F7630009	Asahipak NH2P-40 2D	≥ 5,500	Amino	4	100	2.0 x 150	CH ₃ CN
F7630010	Asahipak NH2P-40 2E	≥ 7,000	Amino	4	100	2.0 x 250	CH ₃ CN
F6710100	Asahipak NH2P-LF	(line filter)	Amino	–	–	8.0 x 75	CH ₃ CN

Base Material: Polyvinyl alcohol









■ 3mm I.D columns [Customized columns]

Product Code	Product Name	Functional Group	Particle Size (μm)	Pore Size (Å)	Column Size (mm) I.D. x Length
F7630011	Asahipak NH2P-40 3B	Amino	4	100	3.0 x 50
F7630012	Asahipak NH2P-40 3D	Amino	4	100	3.0 x 150

Base Material: Polyvinyl alcohol

■ Preparative columns *Preparative columns are made to order.

Product Code	Product Name	Plate Number (TP/column)	Particle Size (μm)	Column Size (mm) I.D. x Length	Standard Column
F6830001	Asahipak NH2P-50 10E	≥ 10,000	5	10.0 x 250	NH2P-50
F6710016	Asahipak NH2P-50G 4A	(guard column)	5	4.6 x 10	(guard column)
F6830031	Asahipak NH2P-90 20F	≥ 10,000	9	20.0 x 300	NH2P-50
F6710017	Asahipak NH2P-130G 7B	(guard column)	13	7.5 x 50	(guard column)



Ligand Exchange Chromatography Columns

Lists summarizing elution volumes of various saccharides using Shodex columns is available. Please refer to our website (<http://www.shodex.com/en/>) or technical notebook (No.2 and 3).

Features

- SC1011**
 - Separates saccharides by combination of ligand exchange and size exclusion modes
 - SC1821**
 - Three types of counter ions are available: Ca²⁺, Pb²⁺, and Na⁺
 - SP0810**
 - Only water is required for the analysis of neutral sugars
 - KS-801 to 802**
 - SC1011 and SC1821 fulfill USP L19 and L22 requirements
 - SP0810 fulfills USP L22 and L34 requirements
 - KS-801 and KS-802 fulfill USP L22 and L58 requirements
-
- KS-803 to 806**
 - Suitable for separation of polysaccharides by size exclusion mode
 - Can be used in combination with other columns e.g., KS-802 and KS-801
 - Only water is required for the analysis of neutral sugars
 - Fulfills USP L22 and L58 requirements
-
- DC-613**
 - Separates elements by combination of ligand exchange and HILIC modes
 - SZ5532**
 - DC-613 can analyze sugars without removing sodium salts in the sample
 - SC1211**
 - SZ5532 is recommended for the separation of disaccharides or trisaccharides
 - SC1211 is suitable for separating sugar alcohols
 - DC-613 fulfills USP L22 and L58 requirements
 - SZ5532 fulfills USP L22 requirements
 - SC1211 fulfills USP L19 and L22 requirements
-
- SC1011-7F**
 - Fulfills mannitol analysis requirements of JP, USP, and EP methods
 - Ca²⁺ modified ligand exchange chromatography column
 - Only water is required for the analysis of neutral sugars
 - Fulfills USP L19 and L22 requirements

Standard columns

[Ligand exchange and size exclusion]

Product Code	Product Name	Plate Number (TP/column)	Functional Group (Counter Ion)	Exclusion Limit (Pullulan)	Particle Size (μm)	Column Size (mm) I.D. x Length	Shipping Solvent
F6378102	SUGAR SC1011	≥ 13,000	Sulfo (Ca ²⁺)	1,000	6	8.0 x 300	H ₂ O
F6378103	SUGAR SC1821	≥ 13,000	Sulfo (Ca ²⁺)	10,000	6	8.0 x 300	H ₂ O
F6700090	SUGAR SC-G 6B (SUGAR SC-LG)	(guard column)	Sulfo (Ca ²⁺)	–	10	6.0 x 50	H ₂ O
F6378105	SUGAR SP0810	≥ 11,000	Sulfo (Pb ²⁺)	1,000	7	8.0 x 300	H ₂ O
F6700081	SUGAR SP-G 6B (SUGAR SP-G)	(guard column)	Sulfo (Pb ²⁺)	–	10	6.0 x 50	H ₂ O
F6378010	SUGAR KS-801	≥ 17,000	Sulfo (Na ⁺)	1,000	6	8.0 x 300	H ₂ O
F6378020	SUGAR KS-802	≥ 17,000	Sulfo (Na ⁺)	10,000	6	8.0 x 300	H ₂ O
F6378025	SUGAR KS-803	≥ 17,000	Sulfo (Na ⁺)	50,000	6	8.0 x 300	H ₂ O
F6378035	SUGAR KS-804	≥ 17,000	Sulfo (Na ⁺)	400,000	7	8.0 x 300	H ₂ O
F6378050	SUGAR KS-805	≥ 9,000	Sulfo (Na ⁺)	5,000,000	17	8.0 x 300	H ₂ O
F6378060	SUGAR KS-806	≥ 9,000	Sulfo (Na ⁺)	*(50,000,000)	17	8.0 x 300	H ₂ O
F6700020	SUGAR KS-G 6B (SUGAR KS-G)	(guard column)	Sulfo (Na ⁺)	–	10	6.0 x 50	H ₂ O

*() Estimated value Base Material: Styrene divinylbenzene copolymer

[Ligand exchange and HILIC]

Product Code	Product Name	Plate Number (TP/column)	Functional Group (Counter Ion)	Particle Size (μm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7001003	RSpak DC-613	≥ 5,500	Sulfo (Na ⁺)	6	100	6.0 x 150	H ₂ O/CH ₃ CN=30/70
F6700170	RSpak DC-G 4A (RSpak DC-G)	(guard column)	Sulfo (Na ⁺)	10	–	4.6 x 10	H ₂ O/CH ₃ CN=30/70
F7001300	SUGAR SZ5532	≥ 5,500	Sulfo (Zn ²⁺)	6	–	6.0 x 150	H ₂ O/CH ₃ CN=30/70
F6700110	SUGAR SZ-G	(guard column)	Sulfo (Zn ²⁺)	6	–	4.6 x 10	H ₂ O/CH ₃ CN=30/70
F7001400	SUGAR SC1211	≥ 5,500	Sulfo (Ca ²⁺)	6	50	6.0 x 250	H ₂ O/CH ₃ CN=75/25
F6700120	SUGAR SC1211G 4A (SUGAR SC-G)	(guard column)	Sulfo (Ca ²⁺)	10	–	4.6 x 10	H ₂ O/CH ₃ CN=75/25

Base Material: Styrene divinylbenzene copolymer

■ For mannitol analysis following JP, USP, and EP methods

Product Code	Product Name	Functional Group (Counter Ion)	Particle Size (μm)	Column Size (mm) I.D. x Length	Shipping Solvent
F6379300	EP SC1011-7F	Sulfo (Ca ²⁺)	8	7.8 × 300	H ₂ O
F6700090	SUGAR SC-G 6B (SUGAR SC-LG) (guard column)	Sulfo (Ca ²⁺)	10	6.0 × 50	H ₂ O
F6379230	USPpak MN-431	Sulfo (Ca ²⁺)	8	4.0 × 250	H ₂ O

*See page 79 for USP40-NF35 Column List.

Base Material: Styrene divinylbenzene copolymer

■ Preparative columns *Preparative columns are made to order.

[Ligand exchange and size exclusion]

Product Code	Product Name	Plate Number (TP/column)	Particle Size (μm)	Column Size (mm) I.D. x Length	Standard Column
F6502007	SUGAR KS-2001	≥ 7,000	13	20.0 × 300	KS-801
F6502008	SUGAR KS-2002	≥ 7,000	13	20.0 × 300	KS-802
F6502009	SUGAR KS-2003	≥ 8,000	13	20.0 × 300	KS-803
F6502010	SUGAR KS-2004	≥ 6,000	18	20.0 × 300	KS-804
F6502011	SUGAR KS-2005	≥ 6,000	18	20.0 × 300	KS-805
F6502012	SUGAR KS-2006	≥ 6,000	18	20.0 × 300	KS-806
F6700002	SUGAR KS-G 8B (SUGAR KS-LG)	(guard column)	13	8.0 × 50	(guard column)

[Ligand exchange and HILIC]

Product Code	Product Name	Plate Number (TP/column)	Particle Size (μm)	Column Size (mm) I.D. x Length	Standard Column
F6514013	RSPak DC-2013	≥ 6,000	10	20.0 × 300	DC-613
F6700402	RSPak DC-G 8B (RSPak DC-LG)	(guard column)	10	8.0 × 50	(guard column)

Elution volumes of saccharides analyzed by Shodex columns

[Partial list only; refer to our website for complete list]

Substances	Elution Volume (mL)					
	SP0810	SC1011	KS-801	SZ5532	NH2P-50 4E	SC1211
Arabinose	10.42	8.91	8.21	5.11	6.18	5.56
D-Arabitol	15.86	11.33	7.63	7.27	6.29	8.16
Dulcitol	20.18	12.76	7.40	9.46	7.45	11.28
meso-Erythritol	12.70	10.09	7.86	5.73	5.43	6.27
D(-)-Fructose	11.05	8.85	7.71	5.37	6.75	5.90
D(+)-Fucose	10.48	8.84	8.09	4.50	5.43	4.96
D(+)-Galactose	9.74	7.98	7.58	6.46	8.10	4.98
Gentiobiose	7.22	6.08	5.75	10.50	16.36	*
Glucose	8.63	7.30	7.17	5.87	8.61	4.76
myo-Inositol	12.77	8.86	7.99	12.63	9.96	7.87
Isomaltose	7.68	6.26	5.95	10.57	15.18	*
Isomaltotriose	7.09	5.75	5.34	21.17	27.55	*
1-Kestose	6.79	5.75	5.26	13.09	20.11	*
Kojibiose	7.56	6.21	5.88	9.65	14.82	*
Lactitol	13.27	8.09	6.13	16.35	11.82	6.67
Lactose	8.05	6.51	5.99	10.12	13.27	4.07
Lactulose	9.13	6.99	6.19	9.16	10.72	4.65
Maltitol	12.23	8.26	6.03	13.04	11.82	6.77
Maltose	7.85	6.34	5.94	8.67	14.24	*
Maltotriose	7.48	5.89	5.38	13.79	24.96	*
Mannitol	15.80	11.10	7.23	8.75	7.39	9.03

Substances	Elution Volume (mL)					
	SP0810	SC1011	KS-801	SZ5532	NH2P-50 4E	SC1211
D-Mannose	10.72	8.17	7.64	5.83	7.84	5.01
Melibiose	8.16	6.45	5.98	11.69	14.70	4.23
Nystose	6.38	5.45	4.93	20.05	31.90	*
Palatinin	2peaks	2peaks	5.90	2peaks	12.73	2peaks
Palatinose	7.84	6.45	5.89	8.08	12.12	3.99
Panose	7.14	5.78	5.32	16.87	25.60	*
D(+)-Raffinose	7.14	5.78	5.29	16.36	20.25	*
Rhamnose	9.77	8.23	7.37	3.93	5.52	4.43
D(-)-Ribose	19.35	13.66	9.04	4.82	5.45	8.64
D(-)-Sorbitol	21.61	13.31	7.42	9.79	7.09	11.88
Sorbose	9.67	8.03	7.38	5.12	7.35	4.92
Stachyose	6.82	5.57	4.97	—	36.22	*
Sucrose	7.54	6.29	5.87	7.91	11.87	*
α-D-Talose	21.33	12.59	8.76	5.69	6.47	8.51
Trehalose	7.62	6.27	5.78	10.85	13.25	*
Trehalulose	8.92	6.95	6.10	9.54	11.68	4.78
Xylitol	19.87	13.14	7.94	7.77	6.10	10.16
Xylobiose	8.16	6.68	6.40	5.65	9.05	*
D(+)-Xylose	9.21	7.90	7.71	4.55	6.58	4.48
D-Xylose	10.64	9.02	8.04	4.06	5.41	5.07

(-)→Not detected (-)→Overlap with solvent peak

(-)→Not detected (-)→Overlap with solvent peak

Column : SUGAR SP0810, SC1011, KS-801
 Eluent : H₂O
 Flow rate : 1.0mL/min
 Detector : RI
 Column temp. : 80°C

Column : SUGAR SC1211
 Eluent : H₂O/CH₃CN=65/35
 Flow rate : 1.0mL/min
 Detector : RI
 Column temp. : 70°C


Column : SUGAR SZ5532
 Eluent : H₂O/CH₃CN=25/75
 Flow rate : 1.0mL/min
 Detector : RI
 Column temp. : 60°C

Column : Asahipak NH2P-50 4E
 Eluent : H₂O/CH₃CN=25/75
 Flow rate : 1.0mL/min
 Detector : RI
 Column temp. : 30°C

Saccharides anomer separation


Saccharides may present their anomers at lower temperatures. By setting the SUGAR series columns at higher temperatures will prevent the anomer separation and this results in providing better chromatograms of each saccharide.

Sample :
0.5% each, 10 μ L




Column : Shodex SUGAR SC1011
Eluent : H₂O
Flow rate : 0.7mL/min
Detector : RI
Column temp.: 29°C, 70°C

Calibration curves for KS-800 series using pullulan



Column : Shodex SUGAR KS-800 series
Eluent : H₂O
Detector : RI
Column temp. : 80°C


Hydrolyzed starch



Column : Shodex SUGAR KS-802 x 2
Eluent : H₂O
Flow rate : 1.0mL/min
Detector : RI
Column temp. : 80°C

Biomass sugars


Sample : 5 μ L
1. Cellobiose 1.0%
2. Glucose 1.5%
3. Xylose 0.5%
4. Galactose 0.5%
5. Arabinose 0.5%
6. Mannose 1.0%



Column : Shodex SUGAR SP0810
Eluent : H₂O
Flow rate : 0.6mL/min
Detector : RI
Column temp. : 85°C

Ketohexoses


Sample : 0.025% each, 10 μ L
1. Sorbose
2. Fructose
3. Tagatose
4. Psicose



Column : Shodex SUGAR SP0810
Eluent : H₂O
Flow rate : 1.0mL/min
Detector : RI
Column temp. : 80°C

Pinitol


Sample : 0.1% each, 20 μ L
1. Sucrose
2. Glucose
3. Pinitol
4. Fructose
5. *chiro*-Inositol
6. *myo*-Inositol



Column : Shodex SUGAR SP0810
Eluent : H₂O
Flow rate : 0.8mL/min
Detector : RI
Column temp. : 85°C

Oligosaccharides in soybean


Sample : 0.1% each, 20 μ L
1. Verbascose
2. Stachyose
3. Raffinose
4. Sucrose
5. Pinitol










Column : Shodex SUGAR KS-802 + KS-801
Eluent : H₂O
Flow rate : 0.6mL/min
Detector : RI
Column temp. : 85°C

Saccharides related to raffinose biosynthesis

Sample : 0.1% each, 20 μ L
1. Verbascose
2. Sucrose
3. Galactinol
4. Galactose
5. *myo*-Inositol



Column : Shodex SUGAR SC1011
Eluent : H₂O
Flow rate : 1.0mL/min
Detector : RI
Column temp. : 80°C



● Ion Exclusion Chromatography Columns

Features

- SH1011** • Columns for simultaneous analysis of saccharides and organic acids
- SH1821** • Separates neutral sugars in size exclusion mode and organic acids in ion exclusion mode
 • Suitable for the analysis of uronic and aldonic acids
 • Fulfills USP L17 and L22 requirements
-
- KC-811** • Columns suitable for the analysis of organic acids
 • Separates compounds by ion exclusion mode and reversed phase mode
 • Highly selective when used with post column method
 • KC-811 6E is suitable for the analysis of cyanide ions and cyanogen chloride in accordance with the Japanese Water Supply Act
 • Fulfills USP L17 and L22 requirements

■ Standard columns

[For simultaneous analysis of saccharides and organic acids]

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Exclusion Limit (Pullulan)	Particle Size (µm)	Column Size (mm) I.D. x Length	Shipping Solvent
F6378100	SUGAR SH1011	≥ 17,000	Sulfo	1,000	6	8.0 x 300	H ₂ O
F6378101	SUGAR SH1821	≥ 17,000	Sulfo	10,000	6	8.0 x 300	H ₂ O
F6700080	SUGAR SH-G	(guard column)	Sulfo	–	10	6.0 x 50	H ₂ O

Base Material: Styrene divinylbenzene copolymer

[For organic acids, cyanide ions and cyanogen chloride]

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Column Size (mm) I.D. x Length	Shipping Solvent
F6378030	RSpak KC-811	≥ 17,000	Sulfo	6	8.0 x 300	0.1% H ₃ PO ₄ aq.
F6378033	RSpak KC-811 6E	≥ 13,000	Sulfo	6	6.0 x 250	0.1% H ₃ PO ₄ aq.
F6700030	RSpak KC-G 6B (RSpak KC-G)	(guard column)	Sulfo	10	6.0 x 50	0.1% H ₃ PO ₄ aq.
F6700010	RSpak KC-G 8B (RSpak KC-LG)	(guard column)	Sulfo	13	8.0 x 50	0.1% H ₃ PO ₄ aq.

*Use KC-G 8B for samples with relatively high impurity and KC-G 6B for samples with relatively low impurity.

Base Material: Styrene divinylbenzene copolymer

■ Preparative columns *Preparative columns are made to order.

Product Code	Product Name	Plate Number (TP/column)	Particle Size (µm)	Column Size (mm) I.D. x Length	Standard Column
F6505012	RSpak KC-2011	≥ 8,000	13	20.0 x 300	KC-811
F6700010	RSpak KC-G 8B (RSpak KC-LG)	(guard column)	13	8.0 x 50	(guard column)

