


# Shodex Technical Training

## NEW Columns

### 2017



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20  
e-mail: [info@mz-at.de](mailto:info@mz-at.de)  
[www.mz-at.de](http://www.mz-at.de)

**Shodex™**



**HPLC Columns**

New

## New Products

### HILIC Columns

- HILICpak VN Series **Polymer-based**
- HILICpak VC Series **Polymer-based**

### SEC Columns

- PROTEIN LW-803 **Silica-based**
- OHpak LB-800 Series **Polymer-based**
- GPC HK-404L **Polymer-based**



# Background - Materials




## Housing Materials

- Steel (SUS)
- PEEK (Polyether Ether Keton)



## Base Materials (Stationary Phase)

unmodified



modified




**FG (Functional Group):**

e.g. -OH, -NH<sub>2</sub>, -CN, -SO<sub>3</sub><sup>-</sup>, -NR<sub>3</sub><sup>+</sup>, -C18, -C8, etc.

# Background – Comparison Silica vs. Polymer

## Base Materials (Stationary Phase)




- Monomeric alkylsilane
- Polymeric alkylsilane
- Endcapping



- Unstable under basic conditions (pH > 7)
- Basic compounds (e.g. amines) are absorbed
- Highly polar compounds are absorbed

# Background – Comparison Silica vs. Polymer



Polymer-based

## Advantages

- Extended pH range (2-13)
- 2-3 times longer lifetime (chemical resistance)
- Almost **no bleeding** (MS, LS & CAD detection)
- **Low salt** concentrations (MS detection)
- Higher **resolution**
- Better **retention** of highly polar compounds
- Simple analysis of **basic compounds**

**FG** = Functional Group

**Polymer** = Polyvinyl alcohol, Polymethacrylate, Polyhydroxymethacrylate or Polystyrene divinylbenzene copolymer

## 2. HILIC

**Product:** Asahipak NH2P HILICpak VG-50 & VT-50

**Application:** small polar molecules, saccharides

**Fields:** pharma, food, beverages, life science

Technical  
Notebook  
No. 2  
HILIC





# HILIC Chromatography - Background

## HILIC: Hydrophilic Interaction Liquid Chromatography

- Suitable for polar analytes (like in **IEC** chromatography)
- Polar stationary phase (like in **NP** chromatography)
- Polar mobile phase, e.g. water/ACN mixtures  
(like in **RP** chromatography, but with higher ACN content (> 95%))

### Functional groups:

- **Amino (-NH<sub>2</sub>)**
- Diol
- Zwitterionic
- Ionic FG
- Pure silica




# HILIC Chromatography - Background

When **HILIC** should be used?

- Analysis of highly polar compounds, saccharides, etc.
- $\text{LogP} < 0$

highly polar  
stationary phase



**Interactions:**


- Hydrophilic partitioning
- Electrostatic interaction
- Hydrogen bonding



# HILIC Chromatography – Stationary Phases


## Asahipak NH2P Serie:

Polyvinyl alcohol




## HILICpak VG Serie:

Polyvinyl alcohol




## HILICpak VT Serie:

Polyvinyl alcohol



## HILICpak VC Serie:


Polyvinyl alcohol



New


## HILICpak VN Serie:

Polyvinyl alcohol



New

# HILICpak VC and VN Series



FG: -COO<sup>-</sup>, Diol

## VC-50

(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630700	<b>HILICpak VC-50 2D</b>	≥ 3,500	Carboxyl	5	100	<b>2.0 x 150</b>	H <sub>2</sub> O
F6711600	<b>HILICpak VC-50G 2A</b>	Guard column	Carboxyl	5	-	<b>2.0 x 10</b>	H <sub>2</sub> O

Base Material : Polyvinyl alcohol

## VN-50


(Housing Material: PEEK)

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7630500	<b>HILICpak VN-50 4D</b>	≥ 10,000	Diol	5	100	<b>4.6 x 150</b>	H <sub>2</sub> O/CH <sub>3</sub> CN=25/75
F6711400	<b>HILICpak VN-50G 4A</b>	-	Diol	5	-	<b>4.6 x 10</b>	H <sub>2</sub> O/CH <sub>3</sub> CN=25/75
F7630600	<b>HILICpak VN-50 2D</b>	≥ 3,500	Diol	5	100	<b>2.0 x 150</b>	H <sub>2</sub> O/CH <sub>3</sub> CN=25/75
F6711500	<b>HILICpak VN-50G 2A</b>	-	Diol	5	-	<b>2.0 x 10</b>	H <sub>2</sub> O/CH <sub>3</sub> CN=25/75

Base Material : Polyvinyl alcohol

New

# HILIC Chromatraphy - Properties of polymeric HILIC Columns




**FG:** -NH<sub>2</sub>, -NR<sub>2</sub>,  
NR<sub>3</sub><sup>+</sup>, -COO<sup>-</sup>, Diol

## Polymeric HILIC Columns

- pH 2 - 13 (depends on column series)
- Regeneration with NaOH Solution
- suitable for ELSD, CAD and LC/MS
- USP L82 (Asahipak NH2P)
- combination of HILIC and IEC possible
- max. pressure up to 15 MPa (150 bar)
- Gradient is possible

New

# HILIC Chromatography - Properties of polymeric HILIC Columns




FG: -NH<sub>2</sub>, -NR<sub>2</sub>,  
NR<sub>3</sub><sup>+</sup>, -COO<sup>-</sup>, Diol

Polymeric HILIC Columns	Application Overview
<b>Asahipak NH2P Series</b> (FG: polyamine, -NH <sub>2</sub> )	<b>Saccharide analysis</b> by HILIC
<b>HILICpak VG Series</b> (FG: -NR <sub>2</sub> ; R: hydrophilic rest)	<b>Saccharide analysis</b> (improved recovery of reducing sugars)
<b>HILICpak VT Series</b> (FG: -NR <sub>3</sub> <sup>+</sup> ; R: hydrophilic rest)	high-sensitivity analysis of <b>anionic compounds</b>
<b>HILICpak VC Series</b> <span>New</span> (FG: -COO <sup>-</sup> )	high-sensitivity analysis of <b>cationic compounds</b> (e.g. choline, acetylcholine, amino acids and amino-based neurotransmitters)
<b>HILICpak VN Series</b> <span>New</span> (FG: Diol)	high-sensitivity analysis of <b>oligosaccharides</b> (e.g. hydrolysis fragments of cellulose)


New

## Asahipak NH2P Series - Applications



**Sample** : 50 ng/mL each, 5  $\mu\text{L}$

1. *meso*-Erythritol;
2. Xylitol;
3. Fructose;
4. Glucose;
5. Sucrose;
6. Lactose;
7. Maltose;
8. Raffinose



Column: Shodex **Asahipak NH2P-40 2D** (2.0 mm I.D. x 150 mm)

Eluent: (A) 0.05%  $\text{NH}_3$  aq. / (B)  $\text{CH}_3\text{CN}$  Isocratic ; (B%) 75%

Flow rate: 0.2 mL/min

Detector: ESI-MS/MS (MRM Negative)


( $m/z$ ) 121 $\rightarrow$ 82(-) for *meso*-Erythritol, 151 $\rightarrow$ 59(-) for Xylitol,

179 $\rightarrow$ 89(-) for Fructose and Glucose, 341 $\rightarrow$ 59(-) for Sucrose,


341 $\rightarrow$ 161(-) for Lactose and Maltose, 503 $\rightarrow$ 89(-) for Raffinose

Column temp. : 40  $^\circ\text{C}$

# HILICpak VG Series - Applications



**Sample:** 50 ng/mL each  
(in H<sub>2</sub>O/CH<sub>3</sub>CN=1/9), 5 µL



Column : Shodex **HILICpak VG-50 2D** (2.0 mm I.D. x 150 mm)

Eluent : (A); 0.1% NH<sub>3</sub> aq. (B); CH<sub>3</sub>CN


Linear gradient ; (B%) 95% (0 to 5min), 95% to 80% (5 to 15min), 80% (15 to 20min)

Flow rate : 0.3 mL/min


Detector : ESI-MS/MS (MRM Negative)

Column temp. : 60 °C

# HILICpak VT Series - Applications



Sample: 5  $\mu\text{L}$   
1  $\mu\text{g/mL}$  each



Column : Shodex **HILICpak VT-50 2D** (2.0 mm I.D. x 150 mm)

Eluent :  $\text{H}_2\text{O}/1\% \text{HCOOH aq.}/\text{CH}_3\text{CN} = 70/20/10$


Flow rate : 0.3 mL/min

Detector : ESI-MS (SIM)

Column temp. : 40 °C


New

# HILICpak VC Series - Applications



## Sample:

0.1  $\mu\text{M}$  each (in  $\text{H}_2\text{O}$ ), 20  $\mu\text{L}$



Column : **Shodex HILICpak VC-50 2D** (2.0 mm I.D. x 150 mm)

Eluent : (A) 200 mM  $\text{HCOOH}$  aq, (B)  $\text{CH}_3\text{CN}$

Linear gradient (High pressure) ; B% = 60% (0 to 5min) 60% to 10% (5 to 6min) 10% (6 to 20min)

Flow rate : 0.3 mL/min



Detector : ESI-MS/MS (MRM Positive)

Column temp. : 40°C



New

# HILICpak VN Series - Applications



Sample: 100 µg/mL , 1 µL

Column: Shodex Asahipak **VN-50 2D**

Eluent: 100 mM HCOONH<sub>4</sub> aq. 30% / Acetonitrile 70% Isocratic

Flow rate: 0.2 mL/min

Column temp.: 40 °C

Detector ESI-MS (SIM)


# Summary

Application	Shodex Column Name	Separation Technique
Saccharides, organic acids, polar compounds	• <b>Asahipak NH2P Series</b>	HILIC
	• <b>HILICpak VG and VT Series</b>	
	• <b>HILICpak VC and VN Series</b> <span>New</span>	
Proteins & peptides	• <b>PROTEIN KW and LW Series</b> <span>New</span>	size exclusion (SEC)
	• <b>OHpak SB and LB Series</b> <span>New</span>	size exclusion (SEC)
Synthetic polymers, plastics	• <b>GPC KF and HK Series</b> <span>New</span>	size exclusion (SEC)


# Shodex Information



Column Catalogue



Detector Catalogue



Technical Notebooks



[www.shodex.de](http://www.shodex.de)  
Applications



## Technical Support

[info@shodex.de](mailto:info@shodex.de)


+49 (0)89/ 93 99 62 37



## Purchase, Quotes

[order@shodex.de](mailto:order@shodex.de)

+49 (0)89/ 93 99 62 34



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20  
e-mail: [info@mz-at.de](mailto:info@mz-at.de)  
[www.mz-at.de](http://www.mz-at.de)