



SiliCycle[®]

Analytical Chemistry Catalog



Analytical Products



SiliaPrepTM

Silica-Based SPE Cartridges and Well Plates

p. 9



SiliaPrepXTM

Polymeric SPE Cartridges and Well Plates

p. 27



SiliaPrepTM

QuEChERS

p. 55



SiliaPlateTM

TLC Plates

p. 63



SiliaSphereTM

Spherical Silica Gels

p. 77



SiliaChrom[®]

HPLC Columns

p. 87



SiliaChrom[®]

HPLC Technical Section

p. 124



SiliCycle[®]

Consumables - Syringe Filters, Vials & Caps

p. 133

Contact Us

p. 161



About SiliCycle

We provide solutions to the global chemical industry.

Founded in 1995, SiliCycle® Inc. is a worldwide leader in the development, the manufacturing and the commercialization of high value silica-based and specialty products for chromatography, analytical and organic chemistry.

Our business extends to more than fifty countries and our customer portfolio includes companies in the pharmaceutical and biopharmaceutical industries, contract research and manufacturing organizations as well as university laboratories, hospital research centers, agriculture and food, environmental, petrochemicals, and industrial process companies.

Since we moved to our new state of the art plant in March 2009, SiliCycle has successfully passed 55 audits, and never failed any. As a certified ISO 9001:2008 company, all procedures and all employees are in line to assure you ultimate quality and an unbeatable customer service.

At SiliCycle, we are at the forefront of the chromatography industry, owing to the extraordinary purity of our silica gels and polymeric sorbents, combined with our capacity to rapidly adapt our products to meet the specific requirements of scientists worldwide.

We lead the way in offering innovative products such as SiliaChrom® HPLC columns, SiliaSphere™ spherical silica gels, SiliaPrep™ silica-based and SiliaPrepX™ polymeric SPE cartridges and Well Plates, SiliaPlate™ TLC plates, SiliaMetS® metal scavengers, SiliaCat® heterogeneous catalysts, SiliaBond® functionalized silica gels, SiliaFlash® irregular silica gels, IMPAQ® angular silica gels, SiliaSep™ flash cartridges, QuEChERS, and lab consumables such as syringe & membrane filters, vials & caps.

We offer a wide variety of first-rate *UltraPure* products. Our automated manufacturing process, which includes acid washing and multiple analyses, is continuously optimized to ensure high purity and a low percentage of fine particles, thereby guaranteeing optimal performance. With our multi-ton manufacturing capability, we are your partner of choice for all your analysis, metal removal, catalysis, synthesis, and purification requirements.

We are committed to providing with the highest quality products and services in the industry.



Table of Contents

About SiliCycle.....	1
Analytical Chemistry Fields.....	4
Word from the President.....	5
Quality Commitment for SiliCycle Products.....	6
Analytical Products.....	7
SiliaPrep™ Silica-Based SPE Cartridges and Well Plates.....	9
Product Selection Guide by Sample Properties.....	11
Product Selection Guide by Manufacturer.....	12
Standard Method Development Procedure.....	14
SiliaPrep Reversed-Phases.....	15
SiliaPrep Normal Phases.....	19
SiliaPrep Ion Exchange Phases.....	21
SiliaPrep Mixed-Mode and Special Phases.....	23
SiliaPrep ClanDRUG.....	24
SiliaPrep CleanENVI.....	26
SiliaPrepX™ Polymeric SPE Cartridges and Well Plates.....	27
Determination of the Optimal SiliaPrepX Polymeric Phase.....	29
Product Selection Guide by Manufacturer.....	29
SiliaPrepX HLB Phase.....	30
SiliaPrepX DVB Phase.....	30
SiliaPrepX SCX Phase.....	31
SiliaPrepX SAX Phase.....	31
SiliaPrepX WCX Phase.....	32
SiliaPrepX WAX Phase.....	32
SiliaPrepX Applications.....	33
SiliaPrepMB™ SPE for MiniBlock® Systems.....	44
Webinar on Applications Developed by SiliCycle.....	45
SiliaPrep™ Accessories.....	46
SiliaPrep Adapters and Vacuum Adapters.....	46
SiliaPrep SPE Vacuum Manifolds.....	47
SiliaPrep 96 Well Collection Plates.....	49
SiliaPrep Phase Separator Columns.....	51
SiliaPrep™ Tips SPE Micro-Cartridges.....	52
SiliaPrep Tips General Experimental Procedure.....	53
SiliaPrep Tips Sorbent Descriptions & Ordering Information.....	54
SiliaQuick™ QuEChERS.....	55
SiliaQuick QuEChERS for Pesticides Residue Analysis.....	56
Schematic Flow Chart of the Most Used QuEChERS Technique.....	58
How to Choose the Proper SiliaQuick QuEChERS Kit.....	59
SiliaPlate™ TLC Plates.....	63
SiliaPlate Feature and Benefits.....	64
SiliaPlate Ordering Information.....	68

Silia <i>Plate</i> TLC Accessories	70
Thin Layer Chromatography Practical Guide	71
Silia<i>Sphere</i>™ Spherical Silica Gels	77
Features and Benefits of Silia <i>Sphere</i> Spherical Silica Gels	78
Silia <i>Sphere</i> Characteristics & Ordering Information	81
SiliCycle; Experts in Functionalization	82
Bonded Silia <i>Sphere</i> Spherical Silica Gels Ordering Information	85
Silia<i>Chrom</i>® HPLC Columns	87
SiliCycle; Experts in HPLC Column Packing	89
Silia <i>Chrom</i> Selection Guide for Small Molecules (MW < 2,000 Da)	91
Silia <i>Chrom</i> HPLC columns Portfolio	92
Silia <i>Chrom</i> HPLC Selection Guide by Manufacturer	94
Silia <i>Chrom</i> HPLC Selection Guide by USP Code	96
How to Choose the Right Silia <i>Chrom</i> C18 Phase	97
Silia <i>Chrom</i> Reversed-Phase HPLC Column Character Evaluation	98
Silia <i>Chrom</i> dt C18 & Silia <i>Chrom</i> AQ C18	101
Silia <i>Chrom</i> SB C18 and C8	107
Silia <i>Chrom</i> XT C18 & XT C18 Fidelity	109
Silia <i>Chrom</i> XDB C18 & XDB C8	111
Silia <i>Chrom</i> XDB1 Family	112
Silia <i>Chrom</i> XDB2 C18	114
Silia <i>Chrom</i> HILIC	115
Silia <i>Chrom</i> SCX-SAX & Others Silia <i>Chrom</i> Products	116
Silia<i>Chrom</i> HPLC Columns for Biochromatography	117
Silia<i>Chrom</i> GF Phases for Size Exclusion Chromatography	118
Silia<i>Chrom</i> IEC Phases for Ion Exchange Chromatography	119
Silia<i>Chrom</i> Chiral Phases for Chiral Chromatography	120
Silia<i>Chrom</i> Phases for Supercritical Fluid Chromatography (SFC)	121
Silia<i>Chrom</i> Guard Cartridges & Holders	122
Silia<i>Chrom</i> HPLC Technical Section	124
Silia <i>Chrom</i> Cleaning and Regeneration Procedures	124
Silia <i>Chrom</i> Suggested Storage Conditions	125
Important HPLC Definitions and Equations	126
HPLC Method Scaling Up or Scaling Down Theory	127
How to Select the Right Silia <i>Chrom</i> HPLC Column	128
Acceptable Modifications to an HPLC Validated Method	130
Silia<i>Chrom</i> HPLC Column Storage Cabinet	131
SiliCycle Consumables	133
SiliCycle Syringe Filters	134
SiliCycle Membrane Filters	138
SiliCycle Vials & Caps	140
Contact Us - Order Now	161

Analytical Chemistry Fields

SiliCycle has developed products that are used in many field of the analytical industry to help customers for their analytical needs. The SiliCycle Analytical Chemistry Catalog is designed in same way that scientists are developing their applications, starting by the sample treatment (*extraction, purification, enrichment, filtration*), to the final analysis (*determination, recovery, yield and selectivity*) including the use of consumables and accessories. In order to facilitate selection of the best product for your requirements, SiliCycle has introduced icons representing each field of the analytical industry.



Bio-Pharma



Pharmaceutical

Biotechnological & Pharmaceutical

Products and applications for each step of the drug discovery & development, purification, characterization, manufacturing and quality control of small pharmaceutical molecules to large peptides and proteins.



Food

Food & Beverage

Products and applications available for food & beverage industry including the food safety testing, fragrance & flavor, quality control testing of intermediate and final products, nutraceutical and natural products analysis.



Forensic

Forensic

Products and applications used for forensic analysis, clinical study and toxicology testing from the preparation of the sample through the analysis.



Environment

Environment

Products and applications covering environmental testing of broad range of matrices such as water, waste water, soil, sludge and air.



Energy

Energy

Products and applications covering the petrochemical, biodiesel and alternative fuels development, testing and analysis.

Word from the President



Dear Colleague,

We are pleased to present you our New SiliCycle Analytical Chemistry Catalog.

The importance of analytical chemistry has never been greater than it is today. Therefore, we have created this new catalog as an essential tool in providing solutions to today's demand for safe food, pure water, safe consumer goods, and safe APIs. Whether you come from the pharmaceutical or biopharmaceutical industry, from agriculture and food, from petrochemicals, environmental industry, quality assurance, quality control or any other analytical lab, this catalog is meant for you.

As a key feature of our new analytical catalog, you will find for the first time in SiliCycle's history, the new SiliaPrepX™ family of polymeric SPE cartridges and well plates. This addition of our new polymeric sorbents marks an important milestone for SiliCycle. It is part of our quest to offer you the most appropriate selection of high quality products providing solutions to the most challenging analytical applications. You will also find other novelties such as QuEChERS, syringe filters, vials & caps, and lots of other consumables... Of course, included within this new catalog, you will find our silica-based best-selling products such as our SiliaChrom® HPLC columns, SiliaPrep™ SPE cartridges and well plates, SiliaSphere™ spherical silica gels, and our SiliaPlate™ TLC plates.

In May 2012, we acquired Chromatography Sciences Company (CSC) Inc. Founded in 1980, CSC was a Canadian pioneer in the manufacturing of HPLC columns and the marketing of other analytical products for the market of research laboratories in North America. With this acquisition, all manufacturing operations, equipments, and know-how of CSC were transferred to SiliCycle state of the art facility, in Quebec City. The key personnel of CSC, including its President, Mr. Denis Boudriau, also joined the SiliCycle team.

We are confident that you will find herein the perfect fit for your day-to-day work. We invite you to visit our Multi-Currency eCommerce website at www.SiliCycle.com for a secure, fast and easy ordering experience, and to get complementary information with regards to our full product lines and services.

You may also contact our highly skilled representatives and knowledgeable technical support people who are available to assist you in application development, and in finding the right solution to any questions you may encounter in your work.

Finally, with over 18 years of market leadership as a worldwide provider of the highest quality products and services, we remain committed to offering you the best and most diversified product lines for analytical and organic chemistry, as well as chromatography, purification and sample clean-up.

To remain at the forefront of the industry, we have increased our presence worldwide in the past few years. Follow us and meet us through our numerous participations in the major trade shows and conferences around the world. It's always a pleasure to meet our fellow colleagues.

Thank you for your confidence and support.

Hugo St-Laurent

President & CEO

Quality Commitment for SiliCycle Products

Quality assurance

Having rigorous quality controls (QC) with high standards does not guarantee absolute satisfaction of the customer. This is why SiliCycle created a Quality Assurance department with a clear goal: always ensure that products are consistently produced and controlled to the quality standards appropriate to their intended use. SiliCycle bases its quality management system (QMS) on the ISO standard.

SiliCycle is ISO 9001:2008 certified. This registration shows that we constantly improve the effectiveness of our quality management system; we follow our policies and fulfill our objectives which lead to customer satisfaction.

You can be sure of the outstanding quality of SiliCycle's silica gels because of the tightly controlled manufacturing conditions at our new state of the art facility. Our tight control of every manufacturing process step, affords identical and reproducible properties (*chemical, physical and structural*) and ensures constant chromatographic selectivities.

Furthermore, our stringent Quality Control and Quality Assurance ensures high performance with no scale-up limitations. Every product meets our quality specifications and is shipped with a Certificate of Analysis (CoA). Individual data sheets are also available directly from our website.

Audits

For many years, SiliCycle products have been used by major players in the industry (*pharmaceutical, biotechnology, etc.*) who are regulated by strict rules (*GMP for example*). SiliCycle has been audited by several customers and successfully passed each one.

Bare Silica Gel

The backbone of most of SiliCycle's products is SiliaFlash F60 (40-63 μm , 60 \AA) silica gel. It provides superior performance for chromatographic applications due to its narrow particle size distribution and high purity.

Before functionalization, every silica is rigorously characterized and analyzed by the procedures below to ensure lot-to-lot reproducibility.

Quality Control	
Type of Analysis	Performed by:
Bare Silica Gel	
Carbon, nitrogen & sulfur content	Elemental analyzer
Total trace metal	ICP-OES
Surface area & porosity	Nitrogen adsorption analyzer
Particle size distribution	Laser light diffraction
Tapped density analysis	Density measurement
Water content	Moisture balance
pH	pH-meter
Functionalized Silica Gel	
Residual solvent content	Moisture balance
Specific reactivity analysis	GC-FID, GC-MS, LC-MS/MS, ICP-OES
Organic function signature	Infrared spectroscopy
Purity analysis	GC-MS


Functionalized Silica Gel

The process for functionalizing the silica is highly dependent on the group being attached. However, it is still possible to functionalize 90% of the surface, verified by ^{29}Si MAS NMR. The remaining 10% of the surface may be endcapped to provide a completely inert support. After being functionalized, the product is submitted to further analysis and quality control as outlined below.





Analytical Products





SiliaPrepTM

Silica-Based SPE Cartridges and Well Plates



SiliaPrep™ SPE Cartridges and Well Plates

Using SiliaPrep SPE Cartridges and Well Plates guarantees the following benefits:

- Choice of a wide variety of SiliaBond high-quality functionalized silica gels.
- Excellent separation (*tight particle size distribution and no fines*).
- High recovery and yield.
- Less time and solvent required for conditioning the sorbent.
- Reproducible flow rates from lot-to-lot.
- Excellent packing and storage qualities.



SiliaPrep Solid-Phase Extraction SPE Cartridges and Well Plates

Solid-phase extraction (SPE) is designed for rapid sample preparation and purification prior to chromatographic analysis. You can optimize your SPE protocols by using SiliCycle SiliaPrep SPE Cartridges and Well Plates.

SiliCycle offers products to meet your specific purification needs. SiliaPrep products are available in different formats including SPE cartridges and 48- & 96-well plates, with different sorbents (*SiliaFlash and SiliaBond*), and in bed weights up to 10 grams (*>10 g are also available in SiliaSep OT formats*)

The well plates are used in high throughput drug discovery and screening, metabolic pharmacokinetic applications, and for automated methods such as a multiprobe approach.

By using SiliaPrep products you will generate higher purity samples and reduce the number of false positives in your screenings, giving you higher quality data. SiliaPrep cartridges are packed with fines-free SiliaFlash silica gel sorbents.

Sorbent Specifications


SiliaPrep products are packed with SiliCycle's SiliaFlash UltraPure silica gels to provide superior performance for all types of applications. This is due to the narrow particle size distribution and high purity. Although the standard products included in this catalog are made of SiliaFlash F60 (40-63 μm , 60 \AA), custom products are available with any type of silica (*irregular, spherical and IMPAQ, etc, in various pore and particle sizes*) offered in our catalog or website and in any format on a custom order basis. Contact us for more information.

Plastic Device Specifications

Standard SiliaPrep cartridges are made with flanged polypropylene (PP) tubes and 20 μm polyethylene (PE) frits. Other plastic materials (*Teflon®, HDPE, etc.*), frit porosity (10 μm), and cartridge rim's (*flangeless*) are available on a custom order basis.



Product Selection Guide by Sample Properties



« I had a difficult time purifying a compound having a basic center by the conventional chromatography on silica gel. Then, I could purify the compound quickly and cleanly with the SiliaPrep SCX cartridge. »

Sangdon Han, Ph.D. from Arena Pharmaceuticals, San Diego, CA, USA

Product Selection Guide by Manufacturer

Product Selection Guide by Manufacturer				
SiliaCycle SiliaPrep	SiliaCycle Part Number	Agilent Bond Elut®	Biotage Isolute®	Macherey-Nagel Chromabond®
Non Polar Phases				
SiliaPrep C18 <i>nec</i> (23 %)	SPE-R30130B-xxx		C18	
SiliaPrep C18 (17 %)	SPE-R31930B-xxx	C18	C18 (EC)	C18 ec
SiliaPrep C18 <i>nec</i> (17 %)	SPE-R35530B-xxx	C18 OH		C18
SiliaPrep C18 WPD	SPE-R33229G-xxx		MFC18	C18 ec f
SiliaPrep C8	SPE-R31030B-xxx		C8 (EC)	
SiliaPrep C8 <i>nec</i>	SPE-R31130B-xxx		C8	C8
SiliaPrep Cyclohexyl	SPE-R61530B-xxx	CH	CH (EC)	C ₆ H ₁₁ ec
SiliaPrep Phenyl	SPE-R34030B-xxx	PH	PH (EC)	C ₆ H ₅
Polar Phases				
SiliaPrep Silica	SPE-R10030B-xxx	SI	SI	SiOH
SiliaPrep Silica WPD	SPE-R10029G-xxx			
SiliaPrep Cyano	SPE-R38030B-xxx	Cyano	CN (EC)	CN
SiliaPrep Diol <i>nec</i>	SPE-R35030B-xxx	Diol (2OH) ^b	DIOL	OH
SiliaPrep Florisil	SPE-AUT-0014-xxx	Florisil	FL	Florisil
SiliaPrep Florisil PR	SPE-AUT-0015-xxx			
SiliaPrep Alumina Acidic	SPE-AUT-0053-xxx	Alumina A (AL-A)	AL-A	Alox A
SiliaPrep Alumina Neutral	SPE-AUT-0054-xxx	Alumina N (AL-N)	AL-N	Alox N
SiliaPrep Alumina Basic	SPE-AUT-0055-xxx	Alumina B (AL-B)	AL-B	Alox B
Ion Exchange Phases				
SiliaPrep SAX <i>nec</i>	SPE-R66530B-xxx	SAX ^b	SAX	SB
SiliaPrep SAX-2 <i>nec</i>	SPE-R66430B-xxx	PRS ^b	PE-AX	
SiliaPrep SCX	SPE-R60530B-xxx	SCX ^b	SCX-3 ^b	SA
SiliaPrep SCX-2	SPE-R51230B-xxx		SCX-2 ^b	PSA
SiliaPrep WAX	SPE-R52030B-xxx	NH ₂ ^b	NH ₂	NH ₂
SiliaPrep Diamine (WAX-2)	SPE-R49030B-xxx	PSA ^b	Diamino	Diamino
SiliaPrep WCX	SPE-R70030B-xxx	CBA	CBA ^b	PCA
Mixed-Mode and Specialty Phases				
SiliaPrep C8/SAX-2 <i>nec</i>	SPM-R026630B-xxx	Certify II	HAX	
SiliaPrep SCX-2/SAX <i>nec</i>	SPM-R802830B-xxx	AccuCAT		
SiliaPrep CleanDRUG	SPEC-R651230B-xxx	Certify ^b	HCX ^d	Drug 1
SiliaPrep CleanENVI	SPEC-R31930B-xxx			C18 PAH
SiliaPrep Activated Carbon	SPE-AUT-0110-xxx	Carbon		
SiliaPrep DL AC/WAX	SP2-R11098-xxx			
SiliaPrep DL AC/Diamine	SP2-R11007-xxx			
SiliaPrep PCB <i>nec</i>	SP2-R00650030B-xxx			SA/SiOH

^a Mallinkrodt Baker, ^b Non-encapped, ^c Encapped, ^d Ion exchange phase is non-encapped xxx = Formats



Avantor Performance Material® Bakerbond®	Phenomenex Strata®	Supelco Discovery® and SupelClean®	Thermo Scientific HyperSep	Waters Sep-Pak®
Octadecyl (C18)	C18-E	DSC-18 and ENVI-18	C18	tC18
Light Load Octadecyl	C18-U			
	C18-T			C18
Octyl (C8)	C8	DSC-8 and ENVI-8	C8	C8
Cyclohexyl (C ₆ H ₁₁)				
Phenyl (C ₆ H ₅)	Phenyl	DSC-Ph and LC-Ph	Phenyl	
Silica gel (SiOH)	Silica (Si-1)	Silica	Silica	Silica
Cyano (CN)	Cyano (CN) ^b	DSC-CN and LC-CN	Cyano	Cyanopropyl
Diol (COHCOH)		DSC-Diol and LC-Diol	Diol	Diol ^b
Florisil (Mg ₂ SiO ₃)		ENVI-Florisil	Florisil	Florisil
	Florisil (FL-PR)			
		LC-Alumina-A		Alumina A
Alumina Neutral	Alumina-N (AL-N)	LC-Alumina-N		Alumina N
		LC-Alumina-B		Alumina B
Quaternary Amine	SAX ^b	DSC-SAX and LC-SAX	SAX	Accell Plus QMA
Aromatic Sulfonic Acid	SCX ^b	DSC-SCX and LC-SCX	SCX	
Amino (NH ₂)	NH ₂ /WAX ^b	DSC-NH ₂ and LC-NH ₂ ^b	Aminopropyl	Aminopropyl
Diamino (NH ₂ NH)		PSA		PSA
Carboxylic Acid (COOH)	WCX ^b	DSC-WCX & LC-WCX		Accell Plus CM
	Screen-A	DSC-MCAX	Verify AX	
			Verify CX	
		ENVI-Carb		AC2
		ENVI-Carb/NH ₂		Carbon Black/Amino
		ENVI-CarbII/PSA		Carbon Black/PSA

All SiliCycle products are endcapped unless noted by « nec » (*non-endcapped*)

Standard Method Development Procedure

Solid-phase extraction methodology will vary depending on the sorbent (*normal, reversed, ion exchange*). Here, we propose generic methods for each mode based on sample and sorbent properties. However, procedures can be slightly different from one sample to another.

Standard Method Development Procedure			
Procedure Step	Reversed-Phase	Ion Exchange Phase	Normal Phase
Analyte properties	Non-polar, uncharged or neutralized, hydrophobic	Ionized or charged	Slightly to moderately polar, uncharged
Matrix sample properties	Organic solvents and aqueous (<i>buffer</i>)	Aqueous (<i>buffer</i>) and pH-adjusted solutions	Organic solvents
Conditioning step	Water-miscible organic solvents	Water-miscible organic solvents or aqueous buffered solution	Sample solvent or methanol
Sample loading	Dissolve analyte in highly polar solvents	Dissolve analyte in highly polar solvents	Dissolve analyte in low polar solvents
Washing	Aqueous or buffered solution and polar solvents	Aqueous solutions containing salts	Non-polar solvents
Elution	Polar or non-polar organic solvents	Polar solvents, may contain acids or bases	Mixture of non-polar (5 - 50%) and polar solvents

Suggested Elution Solvents				
Reversed-Phase	Polarity	Ion Exchange Phase	Polarity	Normal Phase
THF Acetone Ethyl Acetate Acetonitrile Methanol	Low ↓ High	For complete ionization, sample should be adjusted 2 pH units above or below the analyte pKa. pH can be used to neutralize analyte or sorbent. Use 2% strong acid or base in acetonitrile or methanol.	Low ↓ High	Hexane CH ₂ Cl ₂ THF Acetone Acetonitrile




SiliaPrep Reversed-Phases

Description

SiliaPrep C18

SiliCycle recently developed an innovative C18 phase characterized by a homogeneous coverage of the silane on the surface. This strongly hydrophobic and non-polar sorbent is used to extract acidic, neutral and basic compounds from aqueous solutions, various organic compounds from water, and drugs and metabolites from physiological fluids.




- SiliCycle Sorbent Number: R31930B
- Loading: 17% C
- Endcapping: Yes
- Silica type: 60 Å, 500 m²/g, 40 - 63 μm

Description

SiliaPrep C18 WPD

This strongly hydrophobic, non-polar and high-loading capacity sorbent is similar to SiliaPrep C18 but can accommodate larger molecules and untreated matrices.




- SiliCycle Sorbent Number: R33229G
- Loading: 13% C
- Endcapping: Yes
- Silica type: 125 Å, 300 m²/g, 37 - 55 μm

Description

SiliaPrep C18 nec

This strongly hydrophobic and non-polar sorbent is similar to SiliaPrep C18, but presents higher retention and polar selectivity for basic compounds. Unreacted surface OH's can be used for soft condition catch and release purification of glucuronides.



- SiliCycle Sorbent Number: R35530B
- Loading: 17 %C
- Endcapping: No
- Silica type: 60 Å, 500 m²/g, 40 - 63 μm



SiliaPrep Reversed-Phases C18

SiliaPrep Reversed-Phases C18 SPE Formats				
Formats	Qty/Box	SiliaPrep C18	SiliaPrep C18 WPD	SiliaPrep C18 nec
SiliaPrep SPE Cartridges				
1 mL/50 mg	100	SPE-R31930B-01B	SPE-R33229G-01B	SPE-R35530B-01B
1 mL/100 mg	100	SPE-R31930B-01C	SPE-R33229G-01C	SPE-R35530B-01C
3 mL/200 mg	50	SPE-R31930B-03G	SPE-R33229G-03G	SPE-R35530B-03G
3 mL/500 mg	50	SPE-R31930B-03P	SPE-R33229G-03P	SPE-R35530B-03P
6 mL/500 mg	50	SPE-R31930B-06P	SPE-R33229G-06P	SPE-R35530B-06P
6 mL/1 g	50	SPE-R31930B-06S	SPE-R33229G-06S	SPE-R35530B-06S
6 mL/2 g	50	SPE-R31930B-06U	SPE-R33229G-06U	SPE-R35530B-06U
12 mL/2 g	20	SPE-R31930B-12U	SPE-R33229G-12U	SPE-R35530B-12U
25 mL/5 g*	20	SPE-R31930B-20X	SPE-R33229G-20X	SPE-R35530B-20X
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL/200 mg	50	SPC-R31930B-10G	SPC-R33229G-10G	SPC-R35530B-10G
10 mL/500 mg	50	SPC-R31930B-10P	SPC-R33229G-10P	SPC-R35530B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R31930B-P	SPS-R33229G-P	SPS-R35530B-P
1,000 mg	50	SPS-R31930B-S	SPS-R33229G-S	SPS-R35530B-S
SiliaPrep 96-Well Plates				
2 mL/50 mg	1	96W-R31930B-B	96W-R33229G-B	96W-R35530B-B
2 mL/100 mg	1	96W-R31930B-C	96W-R33229G-C	96W-R35530B-C


*Commercialized under SiliaSep OT branding



Determination of Testosterone in Human Urine

General Procedure

1. Mini-SiliaPrep C18 WPD (PN: SPS-R33229G-P) was conditioned with 5 mL of methanol and 5 mL of H₂O.
2. The urine sample (2 mL) was then slowly aspirated through the cartridge.
3. Cartridge was washed with 5 mL of H₂O and 5 mL of hexane.
4. Analyte was eluted with 5 mL of methanol.
5. The sample was evaporated under a nitrogen stream for 30 min at 40°C.
6. The analyte was derivatized using 800 µL of Girard-P (100 mM ammonium acetate buffer, pH 4.2) and 200 µL of methanol maintained at room temperature for 12 h.
7. Quantification was done using LC-MS/MS apparatus.

Testosterone Recovery		
Testosterone	Recovery (%) ^a	
	lot #1	lot #2
	94 ± 2	96 ± 1

^aMean Recovery n = 3, 250 ng/mL



Δ^9 -Tetrahydrocannabinol in Human Plasma

SiliaPrep C18 3 mL/500 mg

SiliCycle PN: SPE-R31930B-03P

Sample Preparation

- Mix 250 μ L of plasma with 1 mL of phosphate buffer (0.1M pH 6.0)

Conditioning Step

- 3 mL of MeOH, 3 mL of HCl 1M and 3 mL of H₂O

Loading Step

- Pass the treated sample through the cartridge

Washing Step

- 2 mL of H₂O
- 1 mL of acetic acid 1M
- 2 mL of (20/80) MeOH/H₂O (v/v)

Elution Step

- 3 mL of (50/50) CH₂Cl₂/Acetone (v/v)

Evaporation Step

- Evaporate under a stream of nitrogen (10 min at 40°C)

Derivatization Step

- Mix under vortex 100 μ L of carbonate buffer 0.1M with 200 μ L of dansyl chloride solution for 1 min (1 mg/mL in acetone)
- Incubate 40 min at 40°C

Liquid-liquid Extraction

- Add 2 mL of 1-chlorobutane
- Centrifugate at 3000 rpm for 5 min

Flash/Freeze Recuperation Step


- Flash/freeze the excess of water from the organic phase in a bath of dry ice/acetone for 3 min

Reconstitution Step


- Evaporate under a stream of nitrogen (10 min at 40°C)
- Reconstitute with 200 μ L of (80/20) ACN/H₂O, 0.1% formic acid (v/v)

Chromatographic Conditions:


Column: SiliaChrom dt C18, 2.5 μ m
 Column Size: 3.0 x 30 mm
 SiliCycle PN: H141802E-H030
 Mobile Phase: MPA: 1 mM ammonium formate in (10/90) H₂O/ACN, 0.1% formic acid (v/v)
 MPB: 1 mM ammonium formate in (90/10) H₂O/ACN, 0.1% formic acid (v/v)
 Temperature: 23°C
 Flow Rate: 1.000 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 325°C, ESI⁺, MRM SCAN
 Injection Volume: 5 μ L



Δ^9 -Tetrahydrocannabinol (THC)




11-nor-9-Carboxy- Δ^9 -Tetrahydrocannabinol (THC-COOH)



11-nor-9-Hydroxy- Δ^9 -Tetrahydrocannabinol (THC-OH)

Recovery Results (n = 6)



Gradient		
Time (min)	MPA (%)	MPB (%)
0	10	90
1.00	10	90
1.01	0	100
3.50	0	100
3.51	10	90
5.00	10	90

SiliaPrep Reversed-Phases

Description

SiliaPrep C8 and SiliaPrep C8 nec

A moderately hydrophobic and non-polar sorbent used to extract extremely non-polar compounds. This phase is more selective than SiliaPrep C18 for large compounds such as PAH, vitamin D, and oils as well as greasy compounds.

- SiliCycle Sorbent Number: R31030B and R31130B (*nec*)
- Loading: 12% C
- Endcapping: Yes (*R31030B*) and No (*R31130B*)
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep Phenyl

A moderately hydrophobic and non-polar sorbent used to extract non-polar compounds with different selectivities through π - π interactions including aromatic compounds and other non-polar phases.

- SiliCycle Sorbent Number: R34030B
- Loading: 9% C
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep Cyano

A moderately polar sorbent used as a normal phase (*less polar compared to silica*) to extract acidic, basic and neutral compounds from aqueous solutions. It is also used as a reversed-phase (*less hydrophobic than C8 and C18*).

- SiliCycle Sorbent Number: R38030B
- Loading: 7% C
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

SiliaPrep Reversed-Phases SPE Formats

Formats	Qty/Box	SiliaPrep C8	SiliaPrep C8 nec	SiliaPrep Phenyl	SiliaPrep Cyano
SiliaPrep SPE Cartridges					
1 mL/50 mg	100	SPE-R31030B-01B	SPE-R31130B-01B	SPE-R34030B-01B	SPE-R38030B-01B
1 mL/100 mg	100	SPE-R31030B-01C	SPE-R31130B-01C	SPE-R34030B-01C	SPE-R38030B-01C
3 mL/200 mg	50	SPE-R31030B-03G	SPE-R31130B-03G	SPE-R34030B-03G	SPE-R38030B-03G
3 mL/500 mg	50	SPE-R31030B-03P	SPE-R31130B-03P	SPE-R34030B-03P	SPE-R38030B-03P
6 mL/500 mg	50	SPE-R31030B-06P	SPE-R31130B-06P	SPE-R34030B-06P	SPE-R38030B-06P
6 mL/1 g	50	SPE-R31030B-06S	SPE-R31130B-06S	SPE-R34030B-06S	SPE-R38030B-06S
6 mL/2 g	50	SPE-R31030B-06U	SPE-R31130B-06U	SPE-R34030B-06U	SPE-R38030B-06U
12 mL/2 g	20	SPE-R31030B-12U	SPE-R31130B-12U	SPE-R34030B-12U	SPE-R38030B-12U
25 mL/5 g*	20	SPE-R31030B-20X	SPE-R31130B-20X	SPE-R34030B-20X	SPE-R38030B-20X
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL/200 mg	50	SPC-R31030B-10G	SPC-R31130B-10G	SPC-R34030B-10G	SPC-R38030B-10G
10 mL/500 mg	50	SPC-R31030B-10P	SPC-R31130B-10P	SPC-R34030B-10P	SPC-R38030B-10P
SiliaPrep 96-Well Plates					
2 mL/50 mg	1	96W-R31030B-B	96W-R31130B-B	96W-R34030B-B	96W-R38030B-B
2 mL/100 mg	1	96W-R31030B-C	96W-R31130B-C	96W-R34030B-C	96W-R38030B-C

*Commercialized under SiliaSep OT branding



SiliaPrep Normal Phases

Description

SiliaPrep Silica

The most polar sorbent, which presents a slightly acidic character and is used to extract various compounds from non-polar solvents through hydrogen bonding.

- SiliCycle Sorbent Number: R10030B
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep Silica WPD

The Silica WPD sorbent is used for the same application as the Silica sorbent but can accommodate larger molecules and untreated matrices.

- SiliCycle Sorbent Number: R10029G
- Silica Type: 125 Å, 300 m²/g, 37 - 55 µm

Description

SiliaPrep Diol *nec*

Moderate polar sorbent presenting neutral character used to extract polar compounds from non-polar solvents and structural isomers. Alternative to silica when the acidic character is problematic.

- SiliCycle Sorbent Number: R35030B
- Loading: 8% C
- Endcapping: No
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

SiliaPrep Normal Phases SPE Formats

Formats	Qty/Box	SiliaPrep Silica	SiliaPrep Silica WPD	SiliaPrep Diol <i>nec</i>
SiliaPrep SPE Cartridges				
1 mL/50 mg	100	SPE-R10030B-01B	SPE-R10029G-01B	SPE-R35030B-01B
1 mL/100 mg	100	SPE-R10030B-01C	SPE-R10029G-01C	SPE-R35030B-01C
3 mL/200 mg	50	SPE-R10030B-03G	SPE-R10029G-03G	SPE-R35030B-03G
3 mL/500 mg	50	SPE-R10030B-03P	SPE-R10029G-03P	SPE-R35030B-03P
6 mL/500 mg	50	SPE-R10030B-06P	SPE-R10029G-06P	SPE-R35030B-06P
6 mL/1 g	50	SPE-R10030B-06S	SPE-R10029G-06S	SPE-R35030B-06S
6 mL/2 g	50	SPE-R10030B-06U	SPE-R10029G-06U	SPE-R35030B-06U
12 mL/2 g	20	FLH-R10030B-15U	FLH-R10029G-15U	SPE-R35030B-12U
25 mL/5 g*	20	FLH-R10030B-25X	FLH-R10029G-25X	SPE-R35030B-20X
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL/200 mg	50	SPC-R10030B-10G	SPC-R10029G-10G	SPC-R35030B-10G
10 mL/500 mg	50	SPC-R10030B-10P	SPC-R10029G-10P	SPC-R35030B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R10030B-P	SPS-R10029G-P	SPS-R35030B-P
1,000 mg	50	SPS-R10030B-S	SPS-R10029G-S	SPS-R35030B-S
SiliaPrep 96-Well Plates				
2 mL/50 mg	1	96W-R10030B-B	96W-R10029G-B	96W-R35030B-B
2 mL/100 mg	1	96W-R10030B-C	96W-R10029G-C	96W-R35030B-C

*Commercialized under SiliaSep OT branding

SiliaPrep Normal Phases

Description

SiliaPrep Florisil and SiliaPrep Florisil PR

A polar sorbent (MgO_3Si) presenting a basic character used to extract non-polar to moderately polar compounds from non-polar solvents. The magnesium ion allows retention of chlorinated pesticides, polychlorinated biphenyl (PCB's) and polysaccharides.

- SiliCycle Sorbent Number: AUT-0014
AUT-0015 (PR)
- Florisil Type: 75 - 150 μm
- Florisil PR Type: 150 - 200 μm

Description

SiliaPrep Alumina-Acidic, Neutral and Basic

Alumina can present either cationic, neutral and acidic character. It is used in a similar fashion as for the SiliaPrep Silica. The difference is that Alumina is more stable at high pH than silica. These sorbents present favorable retention of aromatic compounds, aliphatic amines and compounds containing electronegative functions.

- SiliCycle Sorbent Number: Acidic: AUT-0053
Neutral: AUT-0054, Basic: AUT-0055
- Alumina Type: 60 \AA , 0.9 g/mL, 50 - 200 μm

SiliaPrep Normal Phases SPE Formats

Formats	Qty/Box	SiliaPrep Florisil	SiliaPrep Florisil PR	SiliaPrep Acidic Alumina	SiliaPrep Neutral Alumina	SiliaPrep Basic Alumina
SiliaPrep SPE Cartridges						
1 mL/50 mg	100	SPE-AUT-0014-01B	SPE-AUT-0015-01B	SPE-AUT-0053-01B	SPE-AUT-0054-01B	SPE-AUT-0055-01B
1 mL/100 mg	100	SPE-AUT-0014-01C	SPE-AUT-0015-01C	SPE-AUT-0053-01C	SPE-AUT-0054-01C	SPE-AUT-0055-01C
3 mL/200 mg	50	SPE-AUT-0014-03G	SPE-AUT-0015-03G	SPE-AUT-0053-03G	SPE-AUT-0054-03G	SPE-AUT-0055-03G
3 mL/500 mg	50	SPE-AUT-0014-03P	SPE-AUT-0015-03P	SPE-AUT-0053-03P	SPE-AUT-0054-03P	SPE-AUT-0055-03P
6 mL/500 mg	50	SPE-AUT-0014-06P	SPE-AUT-0015-06P	SPE-AUT-0053-06P	SPE-AUT-0054-06P	SPE-AUT-0055-06P
6 mL/1 g	50	SPE-AUT-0014-06S	SPE-AUT-0015-06S	SPE-AUT-0053-06S	SPE-AUT-0054-06S	SPE-AUT-0055-06S
6 mL/2 g	50	SPE-AUT-0014-06U	SPE-AUT-0015-06U	SPE-AUT-0053-06U	SPE-AUT-0054-06U	SPE-AUT-0055-06U
12 mL/2 g	20	SPE-AUT-0014-12U	SPE-AUT-0015-12U	SPE-AUT-0053-12U	SPE-AUT-0054-12U	SPE-AUT-0055-12U
25 mL/5 g*	20	SPE-AUT-0014-20X	SPE-AUT-0015-20X	SPE-AUT-0053-20X	SPE-AUT-0054-20X	SPE-AUT-0055-20X
SiliaPrep Large Reservoir Volume SPE Cartridges						
10 mL/200 mg	50	SPC-AUT-0014-10G	SPC-AUT-0015-10G	SPC-AUT-0053-10G	SPC-AUT-0054-10G	SPC-AUT-0055-10G
10 mL/500 mg	50	SPC-AUT-0014-10P	SPC-AUT-0015-10P	SPC-AUT-0053-10P	SPC-AUT-0054-10P	SPC-AUT-0055-10P
Mini-SiliaPrep SPE Cartridges						
500 mg	50	SPS-AUT-0014-P	SPS-AUT-0015-P	SPS-AUT-0053-P	SPS-AUT-0054-P	SPS-AUT-0055-P
1,000 mg	50	SPS-AUT-0014-S	SPS-AUT-0015-S	SPS-AUT-0053-S	SPS-AUT-0054-S	SPS-AUT-0055-S
SiliaPrep 96-Well Plates						
2 mL/50 mg	1	96W-AUT-0014-B	96W-AUT-0015-B	n/a	n/a	n/a
2 mL/100 mg	1	96W-AUT-0014-C	96W-AUT-0015-C	n/a	n/a	n/a

*Commercialized under SiliaSep OT branding



SiliaPrep Ion Exchange Phases

Description

SiliaPrep TMA Chloride *nec* (Si-SAX)

Strong anion exchanger sorbent positively charged under all conditions. Used to extract acidic molecules (pK_a 3 - 5).

- SiliCycle Sorbent Number: R66530B
- Loading: ≥ 1.01 mmol/g (or meq/g)
- Endcapping: No
- Silica Type: 60 Å, 500 m²/g, 40 - 63 μm

Description

SiliaPrep TMA Acetate *nec* (Si-SAX-2)

Strong anion exchanger (low-selectivity acetate counter ion) sorbent positively charged under all conditions. Used to extract acidic molecules (pK_a 3 - 5).

- SiliCycle Sorbent Number: R66430B
- Loading: ≥ 0.71 mmol/g (or meq/g)
- Endcapping: No
- Silica Type: 60 Å, 500 m²/g, 40 - 63 μm

Description

SiliaPrep Amine (Si-WAX)

A weak anion exchanger used instead of a strong anion exchanger for strong anions, thus avoiding irreversible retention (acidic molecules $pK_a < 3$). This sorbent is used in different applications such as the separation of peptides, drugs and metabolites from physiological fluids, poly- and monosaccharides and structural isomers.

- SiliCycle Sorbent Number: R52030B
- Loading: ≥ 1.20 mmol/g (or meq/g)
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 μm

SiliaPrep Ion Exchange Phases SPE Formats

Formats	Qty/Box	SiliaPrep TMA Chloride <i>nec</i>	SiliaPrep TMA Acetate <i>nec</i>	SiliaPrep Amine
SiliaPrep SPE Cartridges				
1 mL/50 mg	100	SPE-R66530B-01B	SPE-R66430B-01B	SPE-R52030B-01B
1 mL/100 mg	100	SPE-R66530B-01C	SPE-R66430B-01C	SPE-R52030B-01C
3 mL/200 mg	50	SPE-R66530B-03G	SPE-R66430B-03G	SPE-R52030B-03G
3 mL/500 mg	50	SPE-R66530B-03P	SPE-R66430B-03P	SPE-R52030B-03P
6 mL/500 mg	50	SPE-R66530B-06P	SPE-R66430B-06P	SPE-R52030B-06P
6 mL/1 g	50	SPE-R66530B-06S	SPE-R66430B-06S	SPE-R52030B-06S
6 mL/2 g	50	SPE-R66530B-06U	SPE-R66430B-06U	SPE-R52030B-06U
12 mL/2 g	20	SPE-R66530B-12U	SPE-R66430B-12U	SPE-R52030B-12U
*25 mL/5 g	20	SPE-R66530B-20X	SPE-R66430B-20X	SPE-R52030B-20X
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL/200 mg	50	SPC-R66530B-10G	SPC-R66430B-10G	SPC-R52030B-10G
10 mL/500 mg	50	SPC-R66530B-10P	SPC-R66430B-10P	SPC-R52030B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R66530B-P	SPS-R66430B-P	SPS-R52030B-P
1,000 mg	50	SPS-R66530B-S	SPS-R66430B-S	SPS-R52030B-S
SiliaPrep 96-Well Plates				
2 mL/50 mg	1	96W-R66530B-B	96W-R66430B-B	96W-R52030B-B
2 mL/100 mg	1	96W-R66530B-C	96W-R66430B-C	96W-R52030B-C

*Commercialized under SiliaSep OT branding

SiliaPrep Ion Exchange Phases

Description

SiliaPrep Tonic Acid (Si-SCX)

Strong cation exchanger sorbent positively charged under all conditions. Used to extract basic molecules (pK_a 7 - 10).

- SiliCycle Sorbent Number: R60530B
- Loading: ≥ 0.54 mmol/g (or meq/g)
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep Propylsulfonic Acid (Si-SCX-2)

Strong cation exchanger sorbent positively charged under all conditions. Used to extract basic molecules (pK_a 7 - 10).

- SiliCycle Sorbent Number: R51230B
- Loading: ≥ 0.63 mmol/g (or meq/g)
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep Carboxylic Acid (Si-WCX)

A weak cation exchanger sorbent used to extract strong basic compounds ($pK_a > 9$).

- SiliCycle Sorbent Number: R70030B
- Loading: ≥ 0.92 mmol/g (or meq/g)
- Endcapping: Yes
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

SiliaPrep Ion Exchange Phases SPE Formats

Formats	Qty/Box	SiliaPrep Tonic Acid	SiliaPrep Propylsulfonic Acid	SiliaPrep Carboxylic Acid
SiliaPrep SPE Cartridges				
1 mL/50 mg	100	SPE-R60530B-01B	SPE-R51230B-01B	SPE-R70030B-01B
1 mL/100 mg	100	SPE-R60530B-01C	SPE-R51230B-01C	SPE-R70030B-01C
3 mL/200 mg	50	SPE-R60530B-03G	SPE-R51230B-03G	SPE-R70030B-03G
3 mL/500 mg	50	SPE-R60530B-03P	SPE-R51230B-03P	SPE-R70030B-03P
6 mL/500 mg	50	SPE-R60530B-06P	SPE-R51230B-06P	SPE-R70030B-06P
6 mL/1 g	50	SPE-R60530B-06S	SPE-R51230B-06S	SPE-R70030B-06S
6 mL/2 g	50	SPE-R60530B-06U	SPE-R51230B-06U	SPE-R70030B-06U
12 mL/2 g	20	SPE-R60530B-12U	SPE-R51230B-12U	SPE-R70030B-12U
25 mL/5 g*	20	SPE-R60530B-20X	SPE-R51230B-20X	SPE-R70030B-20X
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL/200 mg	50	SPC-R60530B-10G	SPC-R51230B-10G	SPC-R70030B-10G
10 mL/500 mg	50	SPC-R60530B-10P	SPC-R51230B-10P	SPC-R70030B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R60530B-P	SPS-R51230B-P	SPS-R70030B-P
1,000 mg	50	SPS-R60530B-S	SPS-R51230B-S	SPS-R70030B-S
SiliaPrep 96-Well Plates				
2 mL/50 mg	1	96W-R60530B-B	96W-R51230B-B	96W-R70030B-B
2 mL/100 mg	1	96W-R60530B-C	96W-R51230B-C	96W-R70030B-C

*Commercialized under SiliaSep OT branding



SiliaPrep Mixed-Mode and Specialty Phases

Description

SiliaPrep C8/SAX-2 *nec*

Mixed-mode sorbent designed to extract or isolate acidic and neutral drugs and metabolites from physiological fluids.

- SiliCycle Sorbent Number: R661230B
- Loading: 11% C
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep SCX-2/SAX *nec*

This mixed-mode sorbent is typically used for the separation of acidic and basic molecules from non-ionizable molecules.

- SiliCycle Sorbent Number: R802830B
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm

Description

SiliaPrep PCB *nec*

This special phase is specially designed for extraction of PCB's from waste oil (*hexane extract*).

- SiliCycle Sorbent Number: R00650030B
- Endcapping: No
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm



SiliaPrep Mixed-Mode and Specialty Phases SPE Formats

Formats	Qty/Box	SiliaPrep C8/SAX-2	SiliaPrep SCX-2/SAX	SiliaPrep PCB <i>nec</i>
SiliaPrep SPE Cartridges				
1 mL/50 mg	100	SPM-R661230B-01B	SPM-R802830B-01B	n/a
1 mL/100 mg	100	SPM-R661230B-01C	SPM-R802830B-01C	n/a
3 mL/200 mg	50	SPM-R661230B-03G	SPM-R802830B-03G	n/a
3 mL/500 mg	50	SPM-R661230B-03P	SPM-R802830B-03P	n/a
6 mL/500 mg	50	SPM-R661230B-06P	SPM-R802830B-06P	n/a
6 mL/1 g	50	SPM-R661230B-06S	SPM-R802830B-06S	SP2-R00650030B-06S
6 mL/2 g	50	SPM-R661230B-06U	SPM-R802830B-06U	n/a
12 mL/2 g	20	SPM-R661230B-12U	SPM-R802830B-12U	n/a
25 mL/5 g*	20	SPM-R661230B-20X	SPM-R802830B-20X	n/a
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL/200 mg	50	SPC-R661230B-10G	SPC-R802830B-10G	n/a
10 mL/500 mg	50	SPC-R661230B-10P	SPC-R802830B-10P	n/a

*Commercialized under SiliaSep OT branding

SiliaPrep CleanDRUG

Description

SiliaPrep CleanDRUG

SiliaPrep CleanDRUG is designed to extract specific analytes with more reproducibility and efficacy when using sensitive detectors. This product was developed, tested, and quality controlled for drugs of abuse applications.

- SiliCycle Sorbent Number: R651230B
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm



Easy SPE Method for Drugs of Abuse Determination in Human Urine

General Procedure

1. Sample (0.5 mL) is mixed with 2.5 mL of aqueous H₂SO₄ (0.1 M).
2. SiliaPrep CleanDRUG (3 mL/200 mg cartridges) is conditioned with 2 column volumes of methanol, then 2 column volumes of aqueous H₂SO₄ (0.1 M).
3. Slowly force or aspirate the sample of urine through the cartridge.
4. Wash the cartridge with 3 mL of phosphate buffer (KH₂PO₄/K₂HPO₄ pH = 7.0), then with 3 mL of aqueous H₂SO₄ (0.1 M), and finally with 3 mL of methanol.
5. Analyte is eluted with 2 x 3 mL of aqueous NH₄OH (5% v/v).
6. Sample is evaporated under a nitrogen stream and, reconstituted with distilled water and methanol (9:1 v/v). Finally, the quantification is done using LC-MS apparatus.

SiliaPrep CleanDRUG SPE Formats

Formats	Qty/Box	SiliaPrep Product Number
SiliaPrep SPE Cartridges		
1 mL/50 mg	100	SPEC-R651230B-01B
1 mL/100 mg	100	SPEC-R651230B-01C
3 mL/200 mg	50	SPEC-R651230B-03G
3 mL/500 mg	50	SPEC-R651230B-03P
6 mL/500 mg	50	SPEC-R651230B-06P
6 mL/1 g	50	SPEC-R651230B-06S
6 mL/2 g	50	SPEC-R651230B-06U
12 mL/2 g	20	SPEC-R651230B-12U
25 mL/5 g*	20	SPEC-R651230B-20X

*Commercialized under SiliaSep OT branding

Drugs of Abuse Recovery

Drugs	MDMA	MDEA	Amphetamine
Recovery (%) ^a	96	98	99

^aMean Recovery n = 2, 10 ng/mL to 100 ng/mL



Fentanyl and Norfentanyl in Urine

SiliaPrep CleanDRUG 1 mL/100 mg

SiliCycle PN: SPEC-R651230B-01C

Sample Preparation

- Spike 200 μL of urine and 600 μL of sodium acetate in H_2O (100 mM, pH 6.0) with 40 μL of internal standard (fentanyl- d_5 and norfentanyl- d_5 , 200 ng/mL in MeOH)

Conditioning Step

- 1 mL of MeOH, 1 mL of H_2O and 1 mL of sodium acetate in H_2O (100 mM, pH 6.0)

Loading Step

- Pass the treated sample through the cartridge

Washing Step

- 1 mL of H_2O
- 1 mL of MeOH


Elution Step

- 1 mL of (78/20/2) EtOAc/IPA/ NH_4OH (v/v)


Spotting Step*

- Spot 2 μL in a LazWell™ plate

*Organic phase can be evaporated and reconstituted



Fentanyl



Norfentanyl

LDTD-MS/MS Conditions:

Detector: Phytronix LDTD System on Thermo Vantage Mass Spectrometer
 Gas Flow: 3 L/min
 Mode: MRM, ESI⁺

Laser Pattern	
Time (s)	Power (%)
0	0
2.0	0
5.0	45
7.0	45
7.1	0
8.0	0

MRM Transition			
Drug	MRM Transition	CE	S-Lens
Fentanyl	337 \rightarrow 188	22	120
Fentanyl- d_5	342 \rightarrow 188	22	120
Norfentanyl	233 \rightarrow 150	15	85
Norfentanyl- d_5	238 \rightarrow 155	15	85

Accuracy and Precision Results						
Parameters	Fentanyl			Norfentanyl		
	QC Low	QC Med	QC High	QC Low	QC Med	QC High
Concentration (ng/mL)	25	100	500	25	100	500
N	12	12	12	12	12	12
Mean (ng/mL)	26.38	95.25	481.44	37.58	93.17	489.69
% RSD	2.0	3.9	1.4	15.6	8.1	6.2
% Nominal	105.5	95.2	96.3	110.3	93.2	97.9

SiliaPrep CleanENVI

Description

SiliaPrep CleanENVI

SiliaPrep CleanENVI is designed for typical environmental samples such as PAH's, PCB's, herbicides and herbicides from water or waste water.

- SiliCycle Sorbent Number: R31930B
- Silica Type: 60 Å, 500 m²/g, 40 - 63 µm



Easy SPE Method of Pesticides Determination from Drinking Water

General Procedure




1. SiliaPrep CleanENVI (6 mL/500 mg cartridge) is conditioned with 2 column volumes of methanol, then 2 column volumes of distilled water.
2. Slowly force or aspirate 10 mL of drinking water through the cartridge.
3. Wash the cartridge with 2 column volumes of distilled water (2 x 5 mL).
4. Analyte is eluted with 2 x 3 mL acetone.
5. Sample is evaporated under a nitrogen stream and, reconstituted with distilled water and methanol (1:1 v/v). Finally, the quantification is done using LC-MS apparatus.

SiliaPrep CleanENVI SPE Formats

Formats	Qty/Box	SiliaPrep Product Number
SiliaPrep SPE Cartridges		
1 mL/50 mg	100	SPEC-R31930B-01B
1 mL/100 mg	100	SPEC-R31930B-01C
3 mL/200 mg	50	SPEC-R31930B-03G
3 mL/500 mg	50	SPEC-R31930B-03P
6 mL/500 mg	50	SPEC-R31930B-06P
6 mL/1 g	50	SPEC-R31930B-06S
6 mL/2 g	50	SPEC-R31930B-06U
12 mL/2 g	20	SPEC-R31930B-12U
25 mL/5 g*	20	SPEC-R31930B-20X

*Commercialized under SiliaSep OT branding

Pesticides Recovery

Pesticides	 Atrazine	 Simazine	 Alachlor
Recovery (%) ^a	95	96	86

^aMean Recovery n = 2, 10 ng/mL to 100 ng/mL



SiliaPrepX™

Polymeric SPE Cartridges
and Well Plates



SiliaPrepX™ Polymeric SPE Cartridges and Well Plates

Using SiliaPrepX Polymeric SPE Cartridges and Well Plates guarantees the following benefits:

- High quality products covering a broad spectrum of applications for a wide range of complex matrices.
- Exceptional lot-to-lot reproducibility.
- High recovery and yield.
- Very good separation and flow characteristics (*high loadability and great elution properties allowing low solvent consumption and economical analysis*).
- No contamination from the matrix (*clean extract reducing the ion suppression and increasing the selectivity for LC-MS/MS applications*).



SiliaPrepX Polymeric; An Addition to Our Silica-Based SPE Cartridges

At SiliCycle, as leaders in the industry, we are committed to offer the best and most diversified portfolio for analytical, chromatography and organic chemistry. Although we are primarily renowned for *UltraPure* silica gels, we have created the SiliaPrepX family of polymeric SPE cartridges and well plates to cover the whole spectrum of your solid-phase extraction needs. This new family of polymeric sorbents includes all the phases currently used in the field of sample purification, namely Hydrophilic-Lipophilic Balance (*HLB*), Divinylbenzene (*DVB*), Strong Cation Exchanger (*SCX*), Strong Anion Exchanger (*SAX*), Weak Cation Exchanger (*WCX*) and Weak Anion Exchanger (*WAX*). This complete range of sorbents allows the treatment of the most common matrices, such as human and animal biological fluids, waste waters, petrochemical residues, toxicological residues, food and beverage.

Our new SiliaPrepX polymeric products are made using state of the art technology that provides the highest quality and lot-to-lot reproducibility. In addition, we conduct strict quality controls and analysis during the manufacturing process to remove any impurity or defect that could alter our products.

Plastic Device Specifications


Standard SiliaPrepX cartridges are made with flanged polypropylene (*PP*) tubes and 20 μm polyethylene (*PE*) frits. Other plastic materials (*Teflon*®, *HDPE*, etc.), frit porosity (10 μm), and/or cartridge rim's (*flangeless*) are also available on a custom order basis.



Determination of the Optimal SiliaPrepX Polymeric Phase by a Simple & Logical Method

Follow the simple and logical steps outlined in the flow chart beside to determine the optimal SiliaPrepX polymeric phase to use, as well as to achieve high recovery and clean extract.

- Determine the classification of the analyte (*neutral, acidic or basic compound*)
- Determine the pK_a of the analyte
- Select the SiliaPrepX phase
- Apply the indicated treatment
- Determine recovery by LC analysis



*Polar compounds as organic acids and bases can also be eluted after the Wash/Elution 1.

Product Selection Guide by Manufacturer

Product Selection Guide by Manufacturers					
SiliCycle®	Waters®	Phenomenex®	Agilent®	Supelco®	Macherey-Nagel®
SiliaPrepX HLB	Oasis® HLB	Strata™-X	Plexa® or Nexus®	Supel™ Select HLB	Chromabond® HR-X
SiliaPrepX DVB	Oasis® HLB	Strata™-XL	Plexa® or SimpliQ® PS-DVB	Supel™ Select HLB	Chromabond® HR-P
SiliaPrepX SCX	Oasis® MCX	Strata™-X-C	Plexa® PCX or SimpliQ® SCX	Supel™ Select SCX	Chromabond® HR-XC
SiliaPrepX SAX	Oasis® MAX	Strata™-X-A	Plexa® PAX or SimpliQ® SAX	Supel™ Select SAX	Chromabond® HR-XA
SiliaPrepX WCX	Oasis® WCX	Strata™-X-CW	Nexus® WCX or SimpliQ® WCX	n/a	Chromabond® HR-XCW
SiliaPrepX WAX	Oasis® WAX	Strata™-X-AW	SimpliQ® WAX	n/a	Chromabond® HR-XAW

SiliaPrepX HLB Phase

Description

SiliaPrepX HLB is a wettable copolymer presenting a Hydrophilic-Lipophilic Balance (HLB) allowing a strong retention for neutral, acidic and basic compounds and a higher stability in organic solvents.

Typical Applications

- Drugs and metabolites in biological fluids
- API from tablets, creams, in waste water & drinking water
- Environmental analysis: trace of PAHs, pesticides, herbicides, phenols and PCB in water
- Antibiotics and pesticides in food & beverage

- Particle Size: 40 μm
- Pore Size: 110 \AA
- Surface Area: 850 m^2/g
- pH Stability: 0 to 14

SiliaPrepX DVB Phase

Description

SiliaPrepX DVB is a polystyrene-divinylbenzene copolymer presenting a high hydrophobicity used as a reversed-phase for the extraction of neutral, acidic and basic compounds in viscous matrices or for post synthesis clean-up.

Typical Applications

- Drugs & metabolites in biological fluids
- API from tablets, creams, in waste water & drinking water
- Environmental analysis: trace of PAHs, pesticides, herbicides, phenols & PCBs in water

SiliaPrepX HLB General Extraction Procedure

Conditioning step	1 x Column volume of CH_3OH
Equilibration step	1 x Column volume of H_2O
Loading step	Dilute sample (<i>with H_2O</i>)
Washing step	1 x Column volume of 5% CH_3OH in H_2O
Elution step	1 x Column of CH_3OH

Note: This procedure is a convenient starting point for method development (*format 1 mL/30 mg*). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.

- Particle Size: 85 μm
- Pore Size: 60 \AA
- Surface Area: 950 m^2/g
- pH Stability: 0 to 14

SiliaPrepX DVB General Extraction Procedure

Conditioning step	1 x Column volume of CH_3OH
Equilibration step	1 x Column volume of H_2O
Loading step	Dilute sample (<i>with H_2O</i>)
Washing step	1 x Column volume of 5% CH_3OH in H_2O
Elution step	1 x Column of CH_3OH

Note: This procedure is a convenient starting point for method development (*format 1 mL/30 mg*). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.

SiliaPrepX HLB and DVB Polymeric Formats

Formats	Qty/Box	SiliaPrepX HLB	SiliaPrepX DVB
SiliaPrepX Polymeric SPE Cartridges			
1 mL/30 mg	100	SPE-P0002-01AA	SPE-P0001-01AA
3 mL/60 mg	50	SPE-P0002-03BB	SPE-P0001-03BB
6 mL/100 mg	30	SPE-P0002-06C	SPE-P0001-06C
6 mL/200 mg	30	SPE-P0002-06G	SPE-P0001-06G
6 mL/500 mg	30	SPE-P0002-06P	SPE-P0001-06P
Custom formats available on request			
SiliaPrepX Polymeric 96-Well Plates			
2 mL/10 mg	1	96W-P0002-1A	96W-P0001-1A
2 mL/30 mg	1	96W-P0002-AA	96W-P0001-AA



SiliaPrepX SCX Phase

Description

SiliaPrepX SCX is a polystyrene-divinylbenzene copolymer functionalized by a strong cation exchanger presenting a high selectivity for bases (pK_a 2 - 10). It is highly stable in organic solvents.

Typical Applications

- Basic drugs from biological fluids & tissues
- Pesticides, herbicides, fungicides & melamine from food & beverage

- Particle Size: 85 μ m
- Pore Size: 60 Å
- Surface Area: 800 m²/g
- pH Stability: 0 to 14
- Ionic Capacity: \geq 0.85 mmol/g (or meq/g)

SiliaPrepX SCX General Extraction Procedure

Conditioning step	1 x Column volume of CH ₃ OH
Equilibration step	1 x Column volume of H ₂ O
Loading step	Dilute sample with 1% AcOH in H ₂ O (pH 4-5)
Washing step 1	1 x Column volume of H ₂ O
Washing step 2	1 x Column volume of CH ₃ OH
Elution step	1 x Column of 5% NH ₄ OH in CH ₃ OH

Note: This procedure is a convenient starting point for method development (format 1 mL/30 mg). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.

SiliaPrepX SAX Phase

Description

SiliaPrepX SAX is a polystyrene-divinylbenzene copolymer functionalized by a strong anion exchanger presenting a high selectivity for acids (pK_a 2 - 8). It is highly stable in organic solvents.

Typical Applications

- Acidic compounds & metabolites from biological fluids & tissues
- Food additives & contaminants
- Phenolic acids
- Acidic herbicides

- Particle Size: 85 μ m
- Pore Size: 60 Å
- Surface Area: 900 m²/g
- pH Stability: 0 to 14
- Ionic Capacity: \geq 0.25 mmol/g (or meq/g)

SiliaPrepX SAX General Extraction Procedure

Conditioning step	1 x Column volume of CH ₃ OH
Equilibration step	1 x Column volume of H ₂ O
Loading step	Dilute sample with 5% NH ₄ OH in H ₂ O (pH 7-8)
Washing step 1	1 x Column volume of H ₂ O
Washing step 2	1 x Column volume of CH ₃ OH
Elution step	1 x Column of 2% HCO ₂ H in CH ₃ OH

Note: This procedure is a convenient starting point for method development (format 1 mL/30 mg). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.

SiliaPrepX Strong Exchangers Polymeric Formats

Formats	Qty/Box	SiliaPrepX SCX	SiliaPrepX SAX
---------	---------	----------------	----------------

SiliaPrepX Polymeric SPE Cartridges			
1 mL/30 mg	100	SPE-P0005-01AA	SPE-P0010-01AA
3 mL/60 mg	50	SPE-P0005-03BB	SPE-P0010-03BB
6 mL/100 mg	30	SPE-P0005-06C	SPE-P0010-06C
6 mL/200 mg	30	SPE-P0005-06G	SPE-P0010-06G
6 mL/500 mg	30	SPE-P0005-06P	SPE-P0010-06P

Custom formats available on request

SiliaPrepX Polymeric 96-Well Plates			
2 mL/10 mg	1	96W-P0005-1A	96W-P0010-1A
2 mL/30 mg	1	96W-P0005-AA	96W-P0010-AA

SiliaPrepX WCX Phase

Description

SiliaPrepX WCX is a polystyrene-divinylbenzene copolymer functionalized by a weak cation exchanger used to catch and release strong basic compounds ($pK_a > 10$). It is highly stable in organic solvents.

Typical Applications

- Strong basic compounds from biological fluids & tissues
- Streptomycin from food

- Particle Size: 85 μm
- Pore Size: 60 \AA
- Surface Area: 800 m^2/g
- pH Stability: 0 to 14
- Ionic Capacity: $\geq 0.65 \text{ mmol/g}$ (or meq/g)

SiliaPrepX WAX Phase

Description

SiliaPrepX WAX is a polystyrene-divinylbenzene copolymer functionalized by a weak anion exchanger used to catch and release strong acidic compounds ($pK_a < 2$). It is highly stable in organic solvents.

Typical Applications

- Strong acidic compounds & metabolites from biological fluids & tissues
- Sulfonates & perfluorinated surfactants

- Particle Size: 85 μm
- Pore Size: 60 \AA
- Surface Area: 800 m^2/g
- pH Stability: 0 to 14
- Ionic Capacity: $\geq 0.65 \text{ mmol/g}$ (or meq/g)

SiliaPrepX Weak Exchangers Polymeric Formats

Formats	Qty/Box	SiliaPrepX WCX	SiliaPrepX WAX
SiliaPrepX Polymeric SPE Cartridges			
1 mL/30 mg	100	SPE-P0015-01AA	SPE-P0020-01AA
3 mL/60 mg	50	SPE-P0015-03BB	SPE-P0020-03BB
6 mL/100 mg	30	SPE-P0015-06C	SPE-P0020-06C
6 mL/200 mg	30	SPE-P0015-06G	SPE-P0020-06G
6 mL/500 mg	30	SPE-P0015-06P	SPE-P0020-06P

Custom formats available on request

SiliaPrepX Polymeric 96-Well Plates

2 mL/10 mg	1	96W-P0015-1A	96W-P0020-1A
2 mL/30 mg	1	96W-P0015-AA	96W-P0020-AA

SiliaPrepX WCX General Extraction Procedure

Conditioning step	1 x Column volume of CH_3OH
Equilibration step	1 x Column volume of H_2O
Loading step	Dilute sample with 5% NH_4OH in H_2O (pH 7-8)
Washing step 1	1 x Column volume of H_2O
Washing step 2	1 x Column volume of CH_3OH
Elution step	1 x Column of 2% HCO_2H in CH_3OH

Note: This procedure is a convenient starting point for method development (format 1 mL/30 mg). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.

SiliaPrepX WAX General Extraction Procedure

Conditioning step	1 x Column volume of CH_3OH
Equilibration step	1 x Column volume of H_2O
Loading step	Dilute sample with 1% AcOH in H_2O (pH 4-5)
Washing step 1	1 x Column volume of H_2O
Washing step 2	1 x Column volume of CH_3OH
Elution step	1 x Column of 5% NH_4OH in CH_3OH

Note: This procedure is a convenient starting point for method development (format 1 mL/30 mg). In general, SPE protocols tend to be very specific to each molecule. Further optimization may be required to tailor the method to your application needs.



Marbofloxacin & Sarafloxacin in Salmon

Marbofloxacin and Sarafloxacin are fluoroquinolone antibiotics used in veterinary medicine for infection treatment of skin and respiratory system. The excessive use of fluoroquinolones in food from animals has led to fluoroquinolone-resistant *Salmonella* causing several human health diseases. In 1998, the World Health Organization (**WHO**) has encouraged the promotion of prudent use of quinolones in veterinary medicine in order to minimize the emergence of antimicrobial resistance.

SiliaPrepX SCX 3 mL/60 mg
SiliCycle PN: SPE-P0005-03BB

Sample Preparation

- Add 2 g of salmon and 15 mL of 3% H₃PO₄ aqueous solution in a 50 mL tube
- Shake the tube in a horizontal position for 15 min
- Add 5 mL of hexane and vortex for 2 min
- Centrifuge at 3,000 rpm for 5 min
- Recuperate the aqueous phase from the gelled organic phase by filtration

Conditioning Step

- 3 mL of MeOH, 3 mL of HCl 1M and 3 mL of H₂O

Loading Step

- Pass 3 mL of the filtered sample through the cartridge

Washing Step


- 2 mL of HCl 2M
- 1 mL of MeOH

Elution Step


- 3 mL of 10% NH₄OH in MeOH

Reconstitution Step

- Evaporate under a stream of nitrogen (20 min at 40°C)
- Reconstitute with 1 mL of (20/80) MeOH/H₂O (v/v)



Recovery Results (n = 3)



Chromatographic Conditions:

Column: SiliaChrom dt C18, 3 μm
 Column Size: 3.0 x 30 mm
 SiliCycle PN: H141803E-H030
 Mobile Phase: MPA: 1 mM ammonium formate in (10/90) MeOH/H₂O, 0.1% formic acid (v/v)
 MPB: 1 mM ammonium formate in (90/10) MeOH/H₂O, 0.1% formic acid (v/v)
 Temperature: 23°C
 Flow Rate: 0.800 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 400°C, ESI⁺
 MRM Transition: Marbofloxacin 363.1 → 320.2
 Sarafloxacin 386.3 → 368.3
 Injection Volume: 5 μL

Gradient		
Time (min)	MPA (%)	MPB (%)
0	90	10
0.50	90	10
0.51	65	35
1.50	65	35
1.51	0	100
2.25	0	100
2.26	90	10
4.00	90	10



Acrylamide Determination in Fried Potato Chips

Acrylamide is a chemical contaminant produced during the cooking step at high temperature of foods containing high levels of carbohydrate and low levels of protein. High levels of acrylamide have been detected in french fries and potato chips as well as other fried foods. Acrylamide is suspected to be a cancer agent. Health Canada is currently collecting information on the properties and prevalence of this substance in order to make their assessment. Other governmental health organizations have added acrylamide to the list of potentially cancerogenic substances.

SiliaPrepX HLB 6 mL/200 mg
SiliCycle PN: SPE-P0002-06G

SiliaPrepX SCX 3 mL/60 mg
SiliCycle PN: SPE-P0005-03BB

Sample Preparation

- First extraction:
 - In a PTFE tube, add 1 g of potato chips, 8 mL of 4M NaCl aqueous solution and vortex 1 min
 - Incubate 30 min at 60°C (*vortex 10 sec each 10 min*)
 - Centrifugate for 10 min at 4,500 rpm and collect the supernatant
- Second extraction: repeat previous 3 steps with same potato chips
- Add 1 mL of Cirraz 1^a solution and 1 mL of Cirraz 2^b solution
- Vortex 1 min and centrifugate at 4,500 rpm for 5 min
 - ^aCirraz 1 solution: 15 g K₃[Fe(CN)₆] in 100 mL H₂O
 - ^bCirraz 2 solution: 30 g Zn(O₂CCH₃)₂ in 100 mL H₂O

Chromatographic Conditions:

Column: SiliaChrom dt C18, 2.5 μm
 Column Size: 2.0 x 50 mm
 SiliCycle PN: H141802E-E050
 Mobile Phase: 1 mM ammonium formate in (2/98) MeOH/H₂O, 0.1% formic acid (v/v)
 Temperature: 23°C
 Flow Rate: 0.600 mL/min
 MS Splitting Flow: 0.300 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 400°C, ESI⁺
 MRM Transition: 71.9 → 55.1
 Injection Volume: 10 μL

Step 1: SiliaPrepX HLB 6 mL/200 mg

Conditioning Step

- 3 mL of MeOH and 3 mL of H₂O

Loading Step

- Pass 1.5 mL of the treated sample through the cartridge

Washing Step

- 1.5 mL H₂O

Elution Step

- 3 mL 1% HCO₂H in MeOH

Step 2: SiliaPrepX SCX 3 mL/60 mg

Conditioning Step

- 3 mL of MeOH

Loading Step


- Pass the treated sample from step 1 through the cartridge by gravity and collect this fraction

Washing Step

- 1 mL of MeOH and mix the collected fraction

Reconstitution Step

- Evaporate to dryness with a rotary evaporator at 45°C
- Reconstitute with 2 mL of (5/95) MeOH/H₂O (v/v)






Food



Acrylamide Determination in Fried Potato Chips (con't)

Calibration Curve (100 - 2,000 $\mu\text{g}/\text{Kg}$, $n = 3$)



Method Precision Results ($n = 3$)

Acrylamide added ($\mu\text{g}/\text{Kg}$)	Theoretical Concentration ($\mu\text{g}/\text{Kg}$)	Measured Concentration ($\mu\text{g}/\text{Kg}$)
0	200	195 \pm 15
100	300	308 \pm 19
200	400	394 \pm 11
500	700	699 \pm 11
1,000	1,200	1,260 \pm 28
2,000	2,200	2,116 \pm 47

**HIGH
PRECISION**

Method Accuracy Results ($n = 3$)

Level of Concentration ($\mu\text{g}/\text{Kg}$)	Accuracy (%)	CV (%)
100	102 \pm 6	4.9
500	99 \pm 2	1.0
2,000	96 \pm 2	3.2

**GREAT
ACCURACY**



Determination of Carbendazim in Orange Juice

Carbendazim is a fungicide used in some countries to preserve agricultural crops. The use of this fungicide on oranges and orange juice concentrates is not approved by the Environmental Protection Agency (**EPA**) and is illegal in the United States. Recently, the Food and Drug Administration (**FDA**) has received reports that low levels of the fungicide have been found in some orange juice that contain imported orange juice concentrates.

SiliaPrepX SCX 3 mL/60 mg
SiliCycle PN: SPE-P0005-03BB

Sample Preparation

- Centrifugate 5 mL of orange juice 5 min at 3,000 rpm
- Sample 1 mL of the supernatant
- Add 2 mL of acetic acid 10% and vortex 1 min

Conditioning Step

- 3 mL of MeOH and 3 mL of acetic acid 10%

Loading Step

- Pass 3 mL of the treated sample through the cartridge

Washing Step


- 2 mL of acetic acid 10%
- 2 mL of MeOH

Elution Step

- 3 mL of 5% NH₄OH in MeOH

Reconstitution Step

- Evaporate under a stream of nitrogen (10 min at 40°C)
- Reconstitute with 3 mL of (10/90) MeOH/H₂O (v/v)




Carbendazim

Chromatographic Conditions:

Column: SiliaChrom dt C18, 2.5 μm
 Column Size: 3.0 x 30 mm
 SiliCycle PN: H141802E-H030
 Mobile Phase: 1 mM ammonium formate in (20/80) MeOH/H₂O, 0.1% formic acid (v/v)
 Temperature: 23°C
 Flow Rate: 0.800 mL/min
 MS Splitting Flow: 0.300 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 400°C, ESI⁺
 MRM Transition: 192.1 → 160.2
 Injection Volume: 5 μL

Recovery Results at 0.5 & 100 ng/mL (n = 3)



Method Accuracy Results (n = 3)

Level of Concentration	Concentration (ng/mL)	Accuracy (%)
LLQC	0.5	89 ± 4
QC1 (3x LLQC)	1.5	106 ± 2
QC2 (30% ULQC)	30	89 ± 1
QC3 (70% ULQC)	70	94 ± 1
ULQC	100	107 ± 1

Ion Suppression Measured at 100 ng/ml (n = 3)

Brand	Ion Suppression (%)
SiliaPrepX SCX	-11
Plexa PCX	-11
Oasis MCX	-10
Strata-X-C	-11



Amphetamine Quantification in Human Urine

SiliaPrepX HLB 3 mL/60 mg
SiliCycle PN: SPE-P0002-03BB

Sample Preparation

- 10 mL of urine is treated with 100 μ L of TFA

Conditioning Step

- 3 mL of MeOH and 3 mL of H₂O

Loading Step

- Pass 1 mL of the treated sample through the cartridge

Washing Step







- 3 mL (5/95) MeOH/H₂O, 2% NH₄OH (v/v)
- 3 mL (20/80) MeOH/H₂O, 2% NH₄OH (v/v)
- 1 mL (80/20) MeOH/H₂O (v/v)

Elution Step


- 3 mL MeOH
- 3 mL 2% formic acid in MeOH

Reconstitution Step


- Evaporate under a stream of nitrogen (20 min at 40°C)
- Reconstitute with 3 mL of the mobile phase solution



Recovery Results 5 ng/mL (n = 3)



Recovery Results 250 ng/mL (n = 3)



Chromatographic Conditions:

Column: SiliaChrom XT C18, 5 μ m
 Column Size: 4.6 x 50 mm
 SiliCycle PN: H171805H-N050
 Mobile Phase: 1 mM ammonium formate in (70/30) MeOH/H₂O, 0.1% NaOH (v/v)
 Temperature: 23°C
 Flow Rate: 1.000 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 350°C, ESI+, MRM SCAN
 Injection Volume: 5 μ L

Ion Suppression Measured at 250 ng/mL (n = 3)

Compounds	SiliaPrepX HLB (%)	Plexa* (%)	Oasis* HLB (%)	Strata*-X (%)
Amphetamine	-9	-4	-4	7
MDA	-7	-5	-5	3
MDEA	-6	-5	-4	-2
MDMA	-12	-11	-12	-10
Methamphetamine	-7	-3	-4	7
Phentermine	11	15	15	21



Sulfonamides, Tetracyclines & Pyrimethamines in Milk

SiliaPrepX HLB 3 mL/60 mg
SiliCycle PN: SPE-P0002-03BB

SiliaPrepX DVB 3 mL/60 mg
SiliCycle PN: SPE-P0001-03BB

Sample Preparation

- Vortex (2 min) 250 μ L of 20% trichloroacetic acid (TCA) in H₂O (w/v) with 600 μ L of bovine milk
- Add 2.5 mL of McIlvain buffer (vortex 3 min)
- pH adjustment of the solution at 5.5 with NaOH 1.0 M
- Centrifuge at 3,000 rpm for 5 min

Conditioning Step

- 3 mL of MeOH and 3 mL of H₂O

Loading Step

- Pass 1 mL of the treated sample through the cartridge

Washing Step









- 2 x 3 mL (10/90) MeOH/buffer ammonium acetate pH 5.5 (v/v) and dry the cartridge

Elution Step

- 3 mL of MeOH

Reconstitution Step

- Evaporate under a stream of nitrogen (40 min at 40°C)
- Reconstitute with 300 μ L of (90/10) MeOH/H₂O (v/v)




Chromatographic Conditions:

Column: SiliaChrom dt C18, 2.5 μ m
 Column Size: 3.0 x 30 mm
 SiliCycle PN: H141802E-H030
 Mobile Phase: MPA 1mM ammonium formate in (90/10) MeOH/H₂O, 0.1% formic acid (v/v)
 MPB 1mM ammonium formate in (10/90) MeOH/H₂O, 0.1% formic acid (v/v)
 Temperature: 23°C
 Flow Rate: 1.000 mL/min
 Detector: Sciex API 3000
 Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
 Turbo Ion Spray Heater Temperature: 375°C, ESI⁺, MRM SCAN
 Injection Volume: 10 μ L


Gradient		
Time (min)	MPA (%)	MPB (%)
0	90	10
1.25	90	10
2.51	60	40
3.00	60	40
3.01	90	10
5.00	90	10

Sulfonamides, Tetracyclines & Pyrimethamines in Milk (con't)

Recovery Results 10 ng/mL ($n = 3$)



Recovery Results 100 ng/mL ($n = 3$)



Ion Suppression Measured at 100 ng/mL ($n = 3$)

Compounds	SiliaPrepX HLB (%)	SiliaPrepX DVB (%)	Plexa* (%)	Oasis* HLB (%)	Strata*-X (%)
Sulfathiazole	-1	-1	1	-2	-2
Sulfadiazine	-1	-2	-2	-3	2
Sulfamethoxy-pyridazine	-6	-3	-5	-13	-6
Sulfamethazole	7	8	1	7	7
Sulfamethazine	5	6	1	5	7
Pyrimethamine	-12	-9	-4	-7	-7
Tetracycline	-5	-1	-1	-5	-4
Oxytetracycline	1	-5	-1	-3	2



Determination of Trace Pesticides in Water

SiliaPrepX HLB 6 mL/200 mg
SiliCycle PN: SPE-P0002-06G

Conditioning Step

- 6 mL of MeOH and 3 mL of H₂O (HPLC grade)

Loading Step

- Pass 100 mL of drinking water through the cartridge

Washing Step











- 6 mL of H₂O (HPLC grade)

Elution Step

- 2 x 6 mL of MeOH

Reconstitution Step

- Evaporate under a stream of nitrogen (10 min at 50°C)
- Reconstitute with 2 x 3 mL of MeOH



Chromatographic Conditions:


Column: SiliaChrom dt C18, 3 µm
Column Size: 3.0 x 30 mm
SiliCycle PN: H141803E-H030
Mobile Phase: MPA: 1 mM ammonium formate in (10/90) MeOH/H₂O, 0.1% formic acid (v/v)
MPB: 1 mM ammonium formate in (90/10) MeOH/H₂O, 0.1% formic acid (v/v)
Temperature: 23°C
Flow Rate: 0.800 mL/min
MS Splitting Flow: 0.300 mL/min
Detector: Sciex API 3000
Turbo Ion Spray Heater Gas Flow: 8,000 cc/min
Turbo Ion Spray Heater Temperature: 375°C, ESI⁺, MRM
Injection Volume: 5 µL

Gradient & MRM Transition

Time (min)	MPA (%)	MPB (%)
0.00	100	0
10.00	0	100
10.01	100	0
12.00	100	0

Analyte	MRM Transition
Atrazine	216.1 → 174.1
Benalaxyl	326.2 → 148.2
Carbendazim	192.1 → 160.2
Chloroxuron	291.1 → 72.0
Imazalil	297.0 → 159.1
Methalaxyl	280.2 → 220.3
Myclobutanil	289.1 → 70.0
Propoxur	210.1 → 111.2
Simazine	202.1 → 174.2
Thiambazole	202.2 → 175.1

Recovery Results at 10 & 1,000 pg/mL (n = 3)



■ SiliaPrepX HLB
■ Waters Oasis® HLB



Caffeine, Cotinine & Nicotine in Human Urine

SiliaPrepX HLB 3 mL/60 mg

SiliCycle PN: SPE-P0002-03BB

Sample Preparation

- Mix 500 μ L of urine with 1.5 mL of NaOH 0.1 M

Conditioning Step

- 3 mL of MeOH and 3 mL of H₂O

Loading Step

- Pass 1 mL of the treated sample through the cartridge

Washing Step




- 3 mL H₂O and dry the cartridge

Elution Step

- 3 mL of MeOH

Reconstitution Step

- Evaporate under a stream of nitrogen (20 min at 40 °C)
- Reconstitute with 1.5 mL (80/20) MeOH/H₂O (v/v)



Chromatographic Conditions:

Column: SiliaChrom XT C18, 5 μ m

Column Size: 4.6 x 50 mm

SiliCycle PN: H171805H-N050

Mobile Phase: MPA 1 mM ammonium formate in (90/10) MeOH/H₂O, 0.1% formic acid (v/v)
MPB 1 mM ammonium formate in (10/90) MeOH/H₂O, 0.1% formic acid (v/v)

Temperature: 23°C

MS Splitting Flow: 0.25 mL/min

Flow Rate: 1.000 mL/min

Detector: Sciex API 3000

Turbo Ion Spray Heater Gas Flow: 8,000 cc/min



Turbo Ion Spray Heater Temperature: 375°C, ESI⁺, MRM SCAN

Injection Volume: 10 μ L

Gradient		
Time (min)	MPA (%)	MPB (%)
0	75	25
1.25	75	25
1.26	10	90
2.50	10	90
2.51	100	0
3.00	100	0
3.01	0	100
5.00	0	100

Caffeine, Cotinine & Nicotine in Human Urine (*con't*)

GREAT RECOVERY



LOW ION SUPPRESSION

Ion Suppression Measured at 100 ng/mL ($n = 3$)



Compounds	SiliaPrepX HLB (%)	Plexa® (%)	Oasis® HLB (%)	Strata®-X (%)
Caffeine	-10	-10	-6	-6
Cotinine	-13	-13	-14	-14
Nicotine	8	9	8	8



Caffeine, Cotinine & Nicotine in Human Urine (*con't*)

SiliaPrepX HLB


Lot-to-lot Reproducibility 100 ng/mL ($n = 3$)



SiliaPrepMB™ SPE for MiniBlock® Systems

SiliCycle Products and Mettler-Toledo MiniBlock - An Ideal Partnership

- The productivity enhancement of MiniBlock combined with the cutting-edge technology available from SiliCycle enable chemists to design reactions that eliminate tedious work-up and purification issues.
- The MiniBlock is compatible with the full range of SiliCycle products from the synthesis through the purification.
- SiliCycle has the exclusive distribution of the MiniBlock product line in North America.



Mettler-Toledo MiniBlock

The MiniBlock is an easy to use reaction block designed for parallel synthesis and screening. The unique valve body design of the MiniBlock enables processes where filtration is critical, including solid-phase organic synthesis, use of scavenger resins with solution phase synthesis and parallel purification via Solid-Phase Extraction (SPE).

MiniBlock Reactors

Patented reactor with built-in valve design. Available in 48, 24, 12, and 6-position arrays for reaction vessel volumes respectively of 4mL, 10mL, 20mL and 40mL.

Shaking and Washing Station

High performance orbital shaker with integrated basins for wash and rinse capability. Customized and configured to provide vigorous vortex mixing for up to 2 (*compact*) and 6 (*high capacity*) MiniBlocks.



Parallel Synthesis & Purification

MiniBlock is ideal for parallel synthesis and post-reaction cleanup using SiliaPrepMB SPE prepacked with either our chromatographic and ion exchange phases, our silica-supported reagents and catalysts, or our metal or organic scavengers. This is achieved by transferring reaction products from one block to a second MiniBlock or by passing through an SPE filter plate.






Webinar on Applications Developed by SiliCycle

Reductive Amination



A library of 24 amines (secondary and tertiary) was synthesized in parallel using a SiliaBond Cyanoborohydride reagent. The reductive amination reaction produced the desired products through screening reaction conditions in parallel (time of reaction, solvents) and improved yields from minimal effort. In addition workup was also minimized and performed in parallel via filtration of the reaction mixtures.



Carbon-Carbon Coupling Reactions




Carbon-carbon coupling reactions have generated a lot of interest over the past few years to make novel structures and expand the possibilities that chemists have at their disposal. This webinar will present results from Suzuki, Sonogashira, and Heck coupling reactions run in parallel using novel SiliaCat heterogeneous catalysts for screening and library generation.



Metal Scavenging



SiliaMetS Metal Scavengers are known as an efficient solution for removing residual trace metals within valuable materials, such as Active Pharmaceutical Ingredients (APIs). Metal scavenging results depend on quite a few parameters such as: temperature, nature of the catalyst, solvent, etc. MiniBlock can be used in conjunction with SiliaMetS Metal Scavengers to quickly determine the most efficient conditions with limited time and resources.



Note: Scan QR Codes to view the full webinar hosted by Mettler-Toledo.

All SiliCycle SiliaBond Functionalized Silica Gels, SiliaMetS Metal Scavengers & SiliaCat Heterogeneous Catalysts are available in prepacked SiliaPrepMB SPE Cartridges in two different formats (4.5 mL and 10.0 mL, rimless).

*Contact us for ordering your SiliaPrepMB products or to get more details.

SiliaPrep™ Accessories

SiliaPrep Accessories; the #1 solution to simplify your solid-phase extractions.

- Great complement to our SiliaPrep & SiliaPrepX SPE Cartridges and Well Plates products.
- Wide variety of accessories available to increase your productivity thereby save time and money.



Maximize your Productivity with SiliaPrep Accessories

SiliCycle offers various accessories for SPE Cartridges and Well Plates to simplify method development and expedite highthroughput analysis.

SiliaPrep Adapters

Enable cartridge stacking or easy SPE cartridge connection with syringe or gas lines (*for positive pressure*).

AUT-0172 SiliaPrep Adapter for 1, 3, 6 & 12 mL SPE (10/box)

AUT-0173 SiliaPrep Adapter for 25 & 60 mL SPE (10/box)



AUT-0172



AUT-0173

SiliaPrep Vacuum Adapters

Fast, user friendly, and economical adapters for SPE cartridges. Only a vacuum source is needed.

SiliaPrep Vacuum Adapter - Flasks

Joint	PN	Description
20/40	AUT-0043	20/40 - SiliaPrep Vacuum Adapter
19/22	AUT-0044	19/22 - SiliaPrep Vacuum Adapter
14/22	AUT-0045	14/22 - SiliaPrep Vacuum Adapter

Note: One unit per box.

SiliaPrep Vacuum Adapter - Screw Thread Vials

Thread	PN	Description
22/400	AUT-0046	22/400 Vial - SiliaPrep Vacuum Adapter Without Vial Connector
22/400	AUT-0047	22/400 Vial - SiliaPrep Vacuum Adapter With Vial Connector

Note: One unit per box.



AUT-0043

AUT-0044

AUT-0045



AUT-0046

AUT-0047



SiliaPrep Empty Tubes

Looking to pack your own SPE cartridges using our bulk sorbent, use our SiliaPrep Empty Tubes with frits for this purpose.

SililaPrep Empty Tubes	
Formats	Description
SIM-0007-001	Empty 1 mL SPE tube with 2 frits (100/box)
SIM-0008-003	Empty 3 mL SPE tube with 2 frits (100/box)
SIM-0002-006	Empty 6 mL SPE tube with 2 frits (100/box)
SIM-0003-012	Empty 12 mL SPE tube with 2 frits (100/box)
SIM-0004-020	Empty 25 mL SPE tube with 2 frits (100/box)
SIM-0006-060	Empty 60 mL SPE tube with 2 frits (100/box)
SIM-0009-150	Empty 150 mL SPE tube with 2 frits (20/box)

SiliaPrep SPE Vacuum Manifolds

Run multiple samples simultaneously with a controlled flow rate for higher reproducibility with SiliaPrep SPE Vacuum Manifolds. These manifolds are available in 12, 16 and 24 position configurations and allow consistent extraction and no possibility of cross-contamination from one sample to another.

The design consists in a clear glass chamber equipped with replaceable individual stopcocks (*also known as control valves*) and solvent guide needles. The adjustable racks allow the use of a wide variety of collection vessels including 13 and 16 mm test tubes, autosampler & scintillations vials, volumetric and Erlenmeyer flasks.

Simply apply a vacuum source to elute sample through a cartridge or a disk directly to the collection vessel of choice.

Complete Set Includes

- Glass chamber, vacuum gauge & bleed valve
- Cover, gasket, male and female luer fittings
- Individual stopcocks and needles
- Collection racks (6) with supporting legs, retaining clips, shelves and posts



SiliaPrep SPE Vacuum Manifolds (complete set)

Product Number	Description
AUT-0128-12	SiliaPrep SPE Vacuum Manifold - 12 positions
AUT-0128-16	SiliaPrep SPE Vacuum Manifold - 16 positions
AUT-0129-24	SiliaPrep SPE Vacuum Manifold - 24 positions

SiliaPrep Flash Cartridge Vacuum Manifold

The latest addition to our SPE manifold portfolio is the SiliaPrep Flash Cartridge Vacuum Manifold for large volume samples. This 10 ports manifold can handle 150 mL flash column and collection vessel up to 15 cm long. The complete set comes with the same equipment as the SiliaPrep SPE Vacuum Manifold excepts that it has only four collection racks.



SiliaPrep Flash Cartridge Vacuum Manifold (complete set)

Product Number	Description
AUT-0130-10	SiliaPrep Flash Cartridge Vacuum Manifold - 10 positions

SiliaPrep Vacuum Manifold Accessories



Various replacement parts are available for each SiliaPrep Vacuum Manifold offers by SiliCycle.

SiliaPrep Vacuum Manifold Accessories				
Description	10 positions Vacuum Manifold	12 positions Vacuum Manifold	16 positions Vacuum Manifold	24 positions Vacuum Manifold
SiliaPrep Vacuum Manifold Complete Set	AUT-0130-10 (1/box)	AUT-0128-12 (1/box)	AUT-0128-16 (1/box)	AUT-0129-24 (1/box)
Glass Chamber [Dimensions: Length x Width x Heigh]	AUT 0162 (1/box) [12" x 5.25" x 12"]	AUT-0163 (1/box) [7" x 5.25" x 7"]	AUT-0184 (1/box) [12" x 5.25" x 7"]	AUT-0185 (1/box) [12" x 5.25" x 7"]
Vacuum Gauge, Valve & Glass Chamber Kit	AUT-0186 (1/box)	AUT-0187 (1/box)	AUT-0188 (1/box)	AUT-0189 (1/box)
Top Cover Gasket	AUT-0190 (2/box)	AUT-0174 (2/box)	AUT-0175 (2/box)	AUT-0193 (2/box)
Polypropylene Stopcocks	AUT-0194 (10/box)	AUT-0146 (12/box)	AUT-0195 (16/box)	AUT-0147 (24/box)
Teflon® Stopcocks	AUT-0149-25 (25/box) or AUT-0149-50 (50/box)			
Polypropylene Needles	AUT- 0196 (10/box)	AUT-0154 (12/box)	AUT-0197 (16/box)	AUT-0155 (24/box)
Stainless Steel Needles	AUT-0198 (10/box)	AUT-0178 (12/box)	AUT-0199 (16/box)	AUT-0179 (24/box)
Teflon® Needles	AUT-0200 (100/box)			
Collection Racks Kit (<i>supporting legs, retaining clips, shelves and posts included</i>)	AUT-0201 (1/box)	AUT-0202 (1/box)	AUT-0203 (1/box)	AUT-0204 (1/box)
Plate for 13 mm Test Tubes	-	AUT-0205 (1/box)	AUT-0206 (1/box)	AUT-0207 (1/box)
Plate for 16 mm Test Tubes	-	AUT-0208 (1/box)	AUT-0209 (1/box)	AUT-0210 (1/box)
Plate for 19 mm Test Tubes	AUT-0211 (1/box)	-	-	-
Plate for 25 mm Test Tubes	AUT-0212 (1/box)	-	-	-
Plate for Autosampler Vials	-	AUT-0213 (1/box)	-	-
Plate for Volumetric Flasks	AUT-0214 (1/box)	-	-	-

SiliaPrep Waste Containers

Disposable solvent resistant polypropylene containers are available for the 12 port manifolds. These waste containers greatly simplify sample preparation, solvent disposal and clean-up. Depending on the nature of the solvent used, the waste container can be reused many times prior to discarding.



AUT-0176 (10/box)

SiliaPrep Drying Manifold Covers

SiliaPrep Drying Manifold Covers can be used to concentrate samples with a flow of air or gaz (*nitrogen*). These covers are available for the 12, 16 and 24 ports vacuum manifolds.

SiliaPrep Drying Manifold Covers (1/box)	
Product Number	Description
AUT-0215-12	SiliaPrep Drying Manifold Cover - 12 positions
AUT-0215-16	SiliaPrep Drying Manifold Cover - 16 positions
AUT-0215-24	SiliaPrep Drying Manifold Cover - 24 positions



AUT-0215-12






SiliaPrep 96-Well Collection Plates

NEW



SiliCycle offers SiliaPrep 96-Well Collection Plates made from polypropylene with extremely low extractable levels. These collection plates are available with square deep shape in both 1.0 mL and 2.0 mL well volume and with round bottom in 1 mL only. Cap mats are available for all these collection plates (See page 50 for more details).

SiliaPrep 96-Well Collection Plates	
Product Number	Description
96W-0009	SiliaPrep 96 Well Collection Plate Square Bottom, 2 mL (50/box)
96W-0010	SiliaPrep 96 Well Collection Plate Square Bottom, 1 mL (50/box)
96W-0011	SiliaPrep 96 Well Collection Plate Round Bottom, 1 mL (50/box)

96-Well Collection Plates Square Shape



96-Well Collection Plates Round Shape



SiliaPrep Disposable Reservoir Trays for 96-Well Plates

SiliCycle offers SiliaPrep Disposable Reservoir Trays to collect waste solvents used during activation, loading and washing steps. These disposable trays are made in PVC and are compatible with all manifolds used with well plates.



SiliaPrep Disposable Reservoir Trays	
Product Number	Description
96M-0012	SiliaPrep Disposable Reservoir Trays (25/box)





SiliaPrep 96-Well Plate Cap Mats



SiliCycle offers SiliaPrep 96-Well Plate Cap Mats compatible with most 96-Well Plate available on the market. These cap mats are made from a premium silicone quality with a PTFE coating for ultra low bleed. Slit and 384 well plate cap mats are available under request. Contact us for more details.



SiliaPrep 96 Well Plate Square
Silicone/PTFE Cap Mats



SiliaPrep 96 Well Plate Round
Silicone/PTFE Cap Mats

SiliaPrep 96-Well Plate Cap Mats Ordering Information

SiliaPrep 96-Well Plate Cap Mats			
Well Shape	Quantity per box	Product Number	Description
Square	5 / box	96M-0001S	SiliaPrep 96-Well Plate Square Silicone/PTFE Cap Mats (use with 96W-0009 & 96W-0010 collection plate)
	25 / box	96M-0001S-25	
	50 / box	96M-0001S-50	
	100 / box	96M-0001S-100	
Round	5 / box	96M-0001R	SiliaPrep 96-Well Plate Round Silicone/PTFE Cap Mats (use with 96W-0011 collection plate)
	25 / box	96M-0001R-25	
	50 / box	96M-0001R-50	
	100 / box	96M-0001R-100	

* Contact us if you are looking for a cap mat not listed inside this table.

SiliaPrep Phase Separator Cartridges

NEW

SiliCycle offers SiliaPrep Phase Separator Cartridges to separate aqueous phase from chlorinated solvents under gravity. These ready-to-use cartridges are fitted with a proprietary hydrophobic frit and are a great alternative to liquid-liquid extraction, the most popular technique to do this separation. However, this last method is time consuming, requires the use of a glass funnel which needs to be washed between each separate extraction and is not suitable for multiple extractions. SiliaPrep Phase Separator Cartridges solve these drawbacks and offer many advantages such as:

Why choose SiliaPrep Phase Separator Cartridges

- Ease of use
- Efficient and cost saving
- Comply with “Green Chemistry” philosophy
- Compatible with automated systems


SiliaPrep Phase Separator Cartridges

Product Number	Description
PS-012	SiliaPrep Phase Separator Column 12 mL (100/box)
PS-060	SiliaPrep Phase Separator Column 60 mL (50/box)
PS-150	SiliaPrep Phase Separator Column 150 mL (25/box)

Typical Experimental Procedure

1. Select the appropriate size of SiliaPrep Phase Separator Cartridge to hold the entire sample volume (*both aqueous and chlorinated phases*).
2. Connect the SiliaPrep Phase Separator Cartridge on a vacuum manifold. Ensure the collection vessel volume is sufficient enough to recuperate entirely the organic layer.
(Note: Do not connect the manifold to a vacuum source.)
3. Transfer the sample mixture to be separated on the top of the SiliaPrep Phase Separator Cartridge.
4. After a few seconds (*under gravity*), the water immiscible chlorinated solvent will start to pass through the frit and is collected in the suitable vial already placed inside the manifold.
5. The proprietary frit used in the SiliaPrep Phase Separator Cartridge allows the aqueous layer to be left on the column for at least 48 hours

without passing through the frit.



SiliaPrep Phase Separator Typical Experimental Procedure

Important Advices

- **Process under gravity only - Do not apply vacuum or positive pressure**
The SiliaPrep Phase Separator Cartridges are designed to be used under gravity only. The use of vacuum or positive pressure source can yield to a lost in the separation efficiency.
- **Biphasic or two phase system required**
The sample to be separated needs to contain water and a water immiscible solvent with **greater** density than water to form the lower layer. Most common solvents are dichloromethane, chloroform and chlorinated solvents. Furthermore, try to minimize the presence of water miscible solvent (*i.e. methanol, ethanol or acetone*) which can cause problem to obtain a real biphasic system and consequently, the phase separator may not work effectively.

* To obtain a most efficient compound partition between the aqueous and the organic layer, a liquid-liquid extraction can be done prior to use the phase separator column.

SiliaPrep™ Tips Micro-SPE Cartridges

Using SiliaPrep Tips Micro-SPE Cartridges guarantees the following benefits:


- Simple, fast analyte retention & elution with a minimal loss.
- Sorbents are directly embedded on inner cartridge wall which reduces the risk of contamination.
- Exceptional binding capacity.
- No back-pressure.



SiliaPrep Tips for Micro Sample Preparation

SiliaPrep Tips Micro-SPE Cartridges are designed for micro-purification and micro-extraction of femtomole (*fmol*) to picomole (*pmol*) quantities of analytes prior to the analysis by mass spectrometric and/or chromatographic techniques. The constant improvement of these techniques of analysis has allowed scientists to decrease the limit of quantification of several applications. This lower limit has pushed SPE manufacturers to design new SPE cartridges accepting small volumes of analytes.

These Tips are specially designed to achieve extraction and purification of small molecules, peptides, phosphopeptides and proteins. SiliaPrep Tips Micro-SPE Cartridges are available in 3 different cartridge formats based on the binding capacity of each embedded sorbent. They are packed with our SiliaBond functionalized silica gels and specialty phases to cover the broadest spectrum of applications requiring small volume of analytes. The SiliaBond phases are embedded directly in the inner surface of the tip to provide consistent flow rates. Finally, no glue has been used during the packing procedure in order to prevent any contamination of the analyte.





SiliaPrep Tips General Experimental Procedure

The following lines present the general experimental procedure for the purification and enrichment of small molecules, peptides and proteins using SiliaPrep Tips Micro-SPE Cartridges.

Conditioning Step:

Attach the SiliaPrep Tips to a micropipette and aspirate/expel 5 times the elution solution and 3 times the binding solution.

A) Loading Step:

Aspirate/expel 20 to 50 times the sample to allow the compounds to adsorb to the sorbent.

B) Washing Step:

Aspirate/expel 10 times the binding solution and discard the expelled solution each time.

C) Elution Step:

Aspirate/expel 10 times the elution solution and collect the expelled solution in a suitable clean tube. Repeat with a fresh portion of elution solution if you want to be sure to collect all of the adsorbed compounds. (Note : repeat 3-5 times for the carbon black sorbent.)


Binding Solution:

0.1% formic acid or 0.05% trifluoroacetic acid (TFA).


Elution Solution:

0.1% formic acid or 0.05% trifluoroacetic acid (TFA) plus $\pm 60\%$ of acetonitrile, propanol or methanol depending on the compound polarities.


A) Loading Step



B) Washing Step



C) Elution Step



SiliaPrep Tips Micro-SPE Cartridges Specifications

Tip Volume (μL)	Sample Volume (μL)	Binding Capacity (μg)	Sorbent Weight (μg)
1 - 10	0.5 - 10	1	30
10 - 200	2 - 25	2.5	75
10 - 200	5 - 50	15	400

SiliaPrep Tips Sorbent Descriptions

The table below presents the sorbent descriptions and characteristics available for SiliaPrep Tips Micro-SPE Cartridges.

SiliaPrep Tips Sorbent Descriptions and Ordering Information					
SiliaPrep Tips Sorbent	Description	Tip Formats (μL/μg)	Binding Capacity (μg)	Tips/box	Product Number
C18	This phase presents the highest hydrophobic character of the SiliaPrep Tips sorbent mainly used for small peptides or protein purification, enrichment or desalting analysis.	10 / 30	1	96	SPET-C18-T1
		200 / 75	2.5	96	SPET-C18-T2
		200 / 400	15	96	SPET-C18-T3
C8	The C8 phase is the mid-level hydrophobic sorbent of the reversed-phase family. This phase is mainly used for the sample treatment of proteins and peptides requiring a lower hydrophobic capacity than the C18 sorbent.	10 / 30	1	96	SPET-C8-T1
		200 / 75	2.5	96	SPET-C8-T2
		200 / 400	15	96	SPET-C8-T3
C4	The C4 phase presents the lowest hydrophobic character of the SiliaPrep Tips sorbent mainly used for protein purification, enrichment or desalting analysis.	10 / 30	1	96	SPET-C4-T1
		200 / 75	2.5	96	SPET-C4-T2
		200 / 400	15	96	SPET-C4-T3
Carbon Black	The carbon black sorbent presents both a hydrophilic and hydrophobic character. This phase is mainly used for purification of oligosaccharides and other macromolecules containing sugar moieties.	10 / 30	1	96	SPET-CB-T1
		200 / 75	2.5	96	SPET-CB-T2
		200 / 400	15	96	SPET-CB-T3
TiO ₂	The TiO ₂ sorbent is mainly used for phosphopeptide enrichment and presents high selectivity for multiple phosphylated peptides.	10 / 30	1	96	SPET-TI-T1
		200 / 75	2.5	96	SPET-TI-T2
		200 / 400	15	96	SPET-TI-T3
ZrO ₂	The ZrO ₂ sorbent is mainly used for phosphopeptide enrichment and presents high selectivity for mono-phosphylated peptides.	10 / 30	1	96	SPET-ZR-T1
		200 / 75	2.5	96	SPET-ZR-T2
		200 / 400	15	96	SPET-ZR-T3
TiO ₂ /ZrO ₂	The literature suggests only 30% overlap in phosphopeptides isolated by TiO ₂ versus ZrO ₂ . The ZrO ₂ -TiO ₂ sorbent is an excellent alternative for the enrichment of a broad spectrum of phosphopeptides.	10 / 30	1	96	SPET-TIZR-T1
		200 / 75	2.5	96	SPET-TIZR-T2
		200 / 400	15	96	SPET-TIZR-T3





Silia*Quick*TM
QuEChERS



SiliaQuick™ QuEChERS



Food



Environment

Using SiliaPrep QuEChERS ensures the following benefits:

- Clean extracts from pure products.
- High recovery and lot-to-lot reproducibility.
- Great variety of QuEChERS to cover the full spectrum of food applications.
- Reduction of analysis cost.




SiliaQuick QuEChERS for Pesticide Residue Analysis

The QuEChERS technique was developed in 2003 by USDA scientists to simplify and accelerate the analysis of pesticides in various fruit and vegetable samples. The name **QuEChERS** is formed by an acronym of the properties that are observed with this technique: **Q**uick, **E**asy, **C**heap **E**ffective, **R**ugged and **S**afe. The QuEChERS method has gained in popularity to become the most valuable alternative for the determination of traces of analytes in a high throughput environment. Presently, scientists have expanded the use of this method to the analysis of a vast array of pesticides, herbicides, fungicides and other compounds present in all food and beverage matrices.

The QuEChERS technique can be summarized as a three-step methodology, starting with a liquid extraction, followed by a dispersive solid-phase extraction clean-up, and completed by a LC or GC analysis. The first step is to carry out the extraction of compounds of interest from food or beverage matrices through a solvent (*acetonitrile*). The dispersive solid-phase extraction clean-up is designed to remove specific undesired compounds such as sugars, lipids, organic acids, proteins, pigments and excess water from the final solution. The analysis step consists of a simple injection into a LC-MS/MS or GC-MS instrument to determine the analyte recovery.

Step 1

Liquid Extraction




Step 2

Dispersive SPE Clean-up



Step 3

LC or GC Analysis





SiliaQuick QuEChERS for Food Sample Treatments

SiliaQuick QuEChERS are designed to ensure the ultimate performance in pesticide analysis.


- Quick:** Pre-packed liquid extraction kits and dispersive solid-phase extraction clean-up kits contain the right amount of salts and/or sorbents to suit the specific food matrices, hence eliminating the sample preparation measurement step.
- Cheap:** No specialized equipment or glassware is required to achieve the pesticide residue analysis.
- Effective:** General procedure for all food and beverage matrices allowing a significant reduction of the analysis cost.
- Rugged:** Useful for the treatment of complex food matrices such as fish, meat or nuts without the requirement of additional treatments.
- Safe:** Limited time of contact with dangerous compounds and solvents.

Extraction and Dispersive Reagents

The following table presents each extraction and dispersive reagent and their specific functions in the QuEChERS technique.

Extraction and Dispersive Reagents	
Extraction Reagents	Specific Function
Magnesium Sulfate Anhydrous (MgSO ₄)	Facilitates solvent partitioning.
Acetic Acid	Used for pH adjustment.
Acetonitrile	Solvent providing the best characteristics for extracting a wide variety of pesticides. Amenable for both LC and GC analysis.
Buffers	Maintain optimal pH and prevent pH degradation of sensitive analytes.
Sodium Chloride (NaCl)	Limits the amount of polar interferences.
Dispersive Reagents	Specific Function
SiliaMetS Diamine	Removes sugars, fatty acids, organic acids, lipids, and some pigments. Sterols and additional lipids can also be removed in combination with SiliaBond C18.
SiliaBond Amine	Removes sugars and fatty acids as well as the SiliaMetS Diamine but is less likely to catalyze degradation of base sensitive analytes.
SiliaBond C18	Removes long chain, non-polar compounds, and sterols.
Carbon Black	Removes pigments, polyphenols, and other polar compounds.
Magnesium Sulfate Anhydrous (MgSO ₄)	Removes residual water from the organic phase.

Schematic Flow Chart of the Most Used QuEChERS Technique





How to Choose the Proper SiliaQuick QuEChERS Kit

Step 1: For Liquid Extraction

The table below presents the SiliaQuick QuEChERS liquid extraction kits specially pre-packed with anhydrous salts and/or sorbents to suit the QuEChERS technique of your choice.

SiliaPrep QuEChERS Liquid Extraction Kits			
QuEChERS Method	Content	Units/box	Product Number without 50 mL Tube
Original QuEChERS	4 g magnesium sulfate anhydrous 1 g sodium chloride	100	QE-0001-50T
Buffered AOAC 2007.01	6 g magnesium sulfate anhydrous 1.5 g sodium acetate	100	QE-0002-50T
Buffered EN 15662	4 g magnesium sulfate anhydrous 1 g sodium chloride 1.5 g sodium citrate dibasic sesquihydrate 0.5 g sodium citrate tribasic dihydrate	100	QE-0003-50T

SiliaQuick QuEChERS Troubleshooting

Poor recovery of pesticide compounds:

- Each sample has to be at the minimum 80% hydrated to perform optimal liquid extraction.
- For base sensitive compounds use buffered method.
- Always mix the sample with the solvent first to reduce the exothermic reaction between the magnesium sulfate and water.
- Add an analyte protector like toluene or sorbitol to prevent loss of thermally unstable pesticides in the GC inlet.
- Add formic acid after the dispersive SPE clean-up step to limit the degradation of base sensitive compounds prior the LC analysis.

Step 2: For Dispersive Solid-Phase Extraction Clean-Up

The following table presents the SiliaQuick QuEChERS dispersive solid-phase extraction clean-up kits to match your food matrices. It is recommended to use the 2 mL dispersive tube for an extract volume of 1 mL and the 15 mL dispersive tube for extract volumes higher than 3 mL.

SiliaQuick QuEChERS Dispersive Solid-Phase Extraction Kits								
Matrix	Method	Tube (mL)	Units/box	Content (mg)				Product Number
				MgSO ₄	PSA	CB	C18	
General Fruits & Vegetables	AOAC 2007.01	2	100	150	50	-	-	QD-1000-2T
	EN 15662	2	100	150	25	-	-	QD-1001-2T
	AOAC 2007.01	15	50	1,200	400	-	-	QD-2000-15T
	EN 15662	15	50	900	150	-	-	QD-2001-15T
Pigmented Fruits & Vegetables	AOAC 2007.01	2	100	150	50	50	-	QD-1002-2T
	EN 15662	2	100	150	25	2.5	-	QD-1003-2T
	AOAC 2007.01	15	50	1,200	400	400	-	QD-2002-15T
	EN 15662	15	50	900	150	15	-	QD-2003-15T
Highly Pigmented & Fatty Fruits and Vegetables	AOAC 2007.01	2	100	150	50	50	50	QD-1004-2T
	EN 15662	2	100	150	25	7.5		QD-1005-2T
	AOAC 2007.01	15	50	1,200	400	400	400	QD-2004-15T
	EN 15662	15	50	900	150	45		QD-2005-15T
Fatty and Waxed Fruits & Vegetables	AOAC 2007.01	2	100	150	50	-	50	QD-1006-2T
	EN 15662	2	100	150	25	-	25	QD-1007-2T
	AOAC 2007.01	15	50	1,200	400	-	400	QD-2006-15T
	EN 15662	15	50	900	150	-	150	QD-2007-15T

MgSO₄ = Magnesium sulfate anhydrous, PSA = SiliaMetS Diamine, CB = Carbon Black, and C18 = SiliaBond C18






























Choose your SiliaQuick QuEChERS Dispersive SPE Clean-Up Kit by Food Type

The SiliaQuick QuEChERS dispersive solid-phase extraction clean-up kits are assembled to match food matrices to the right method.

SiliaQuick QuEChERS Dispersive Solid-Phase Extraction Kits				
Food Matrices	General Fruits & Vegetables	Pigmented Fruits & Vegetables	Highly Pigmented and Fatty Fruits & Vegetables	Fatty and Waxed Fruits & Vegetables
Root and Tuber Vegetables				
Beets				
Carrot				
Radish				
Potato				
Fruiting Vegetables				
Eggplant				
Cucumber				
Pepper (<i>green or red</i>)				
Pumpkin				
Tomato				
Cabbage				
Broccoli				
Brussels sprouts				
Cauliflower				
Stem Vegetables				
Aparagus				
Celery				
Leek				
Rhubarb				
Leafy Vegetables				
Lettuce				
Basil				
Parsley				
Spinach				
Leek Plants				
Garlic				
Onion				
Shallot				


Choose your SiliaQuick QuEChERS Dispersive SPE Clean-Up Kits by Food Types (*con't*)

SiliaQuick QuEChERS Solid-Phase Extraction Kits				
Food Matrices	General Fruits & Vegetables	Pigmented Fruits & Vegetables	Highly Pigmented and Fatty Fruits & Vegetables	Fatty and Waxed Fruits & Vegetables
Small Fruits				
Blackberry				
Blueberry				
Grapes (<i>red</i>)				
Cranberry				
Strawberry				
Pome Fruits				
Apple				
Pear				
Quince				
Citrus Fruits				
Grapefruit				
Lemon & Lime				
Orange				
Tangerine				
Stone Fruits				
Apricot				
Cherry				
Peach				
Plum				
Other Fruits				
Avocado				
Banana				
Mango				
Pineapple				
Other				
Cereals (<i>wheat, corn, rice</i>)				
Coffee beans				
Tea Leaves				



SiliaPlateTM

TLC Plates



Thin Layer Chromatography (TLC)

SiliCycle is your partner of choice for your purification and chromatography needs.

- Optimize your separation conditions by using the same silica gel as in your flash columns and cartridges.
- Made with an extra hard layer that ensures the plates don't lose silica on rubbing and heating.
- The consistent thickness of our SiliaPlate™ ensures lot-to-lot reproducibility.



Introduction to Thin Layer Chromatography (TLC)

Thin layer chromatography (TLC) is a quick, simple and inexpensive analytical technique frequently used in various laboratories. It is used for reaction monitoring, compound purity evaluation as well as a rapid and cost-efficient selection and optimization of chromatographic conditions prior to purification by flash chromatography or HPLC.

Besides speed and low cost, TLC analysis presents other non-negligible advantages like the small quantity of compound required and high sample throughput capability (*up to 20 samples simultaneously*).

SiliaPlate Features and Benefits

For over 18 years, SiliCycle has been offering a wide selection of TLC plates in various sizes (*plate size, thickness*) and chemistries (*10% Silver Nitrate, CN, C18, NH₂*). SiliaPlate represents an efficient and economical alternative to other TLC plate manufacturers while demonstrating high separation power, which is due to the narrow particle size distribution silica gel used for manufacturing.

The extraordinary silica layer hardness combined to a homogeneous coating and layer thickness allow excellent separation. Each TLC batch is chemically and physically controlled by our Quality Control department to ensure lot-to-lot and layer-to-layer reproducibility.

« Many products have been successfully purified with the silica gel. We have had problems with other companies' TLC plates not running the same as their silica gel, but everything was fixed when we switched over to all SiliCycle products.. »

William Nguyen from Stanford University, Stanford, CA, USA



Types of plates available (TLC/HPTLC/PLC)

SiliaCycle offers different types of plates for thin layer chromatography applications: classical TLC, high performance TLC (*also called HPTLC*) and preparative TLC (*PLC*). The plate types are selected based on the type of analysis required and the available budget.

Differences between classical TLC, HPTLC and PLC			
Properties	Classical TLC	HPTLC	Preparative PLC
Application	Quick, inexpensive, flexible and portable separations	Highly sophisticated separation problems, complex samples	Purification on a TLC plate
Analysis	Qualitative	Qualitative & Quantitative	Quantitative
Detection	UV - Stains	Instrumented analysis (<i>use of scanners for detection</i>)	UV
Price	Lower prices than HPTLC	Higher prices than TLC	-
Distribution [Mean Particle Size]	5 - 20 μm [10 - 12 μm]	4 - 8 μm [5 - 6 μm]	5 - 40 μm [25 μm]
Layer Thickness	250 μm	200 μm	0.5 mm, 1 mm, 2 mm
Typical Sample Volume	1 - 5 μL	0.1 - 0.5 μL	5 - 20 μL

TLC Backings

TLC plates are available with different backings (*also called supports*): rigid (*glass-backed*) or flexible sheets (*aluminum & plastic-backed*). Glass backed plates are the most frequently used due to the ease of handling, transparency (*spot can be seen on both sides*) as well as the chemical resistance and inertness of the support. However, glass plates also present certain disadvantages like superior fragility and higher weight over flexible backings. On the other hand, aluminum and plastic backings also offer both pros and cons as presented in the table below.

TLC Backings Comparison			
Properties	Glass	Aluminum	Plastic
Approximate Thickness	1.5 mm	1.5 mm	2 mm
Total Weight	High	Low	Medium
Heating Stability	High	High	Below 175°C
Fragility	High	Low	Low
Scissors Cut	Impossible	Easily	Possible
Chemical Resistance Against			
Mineral Acids	High	Low	High
Bases (<i>ammoniac</i>)	High	Low	High

Available Sorbents

Various adsorbents can be used for TLC coating; silica, alumina, florisil, etc. However, silica gel is probably the most versatile since it covers almost all types of separation (*if the right solvent system is selected*). More than 80% of all purifications are performed using silica gel as the adsorbent.

The particle size distribution used for the silica is related to the nature of the plate. For standard TLC, silica gel with a mean particle size of 10 - 14 μm is used compared to HPTLC where a smaller particle size is required. In both cases, pore diameter is always 60 Å. Some functionalized silica gels like reversed-phase (*C18, C8, Amine, Cyano, Diol, ...*) and specialty (*Silver Nitrate*) plates can also be used as TLC adsorbent for particular needs.

The two most popular modes of separation employed in TLC are normal and reversed phases. In normal phase separation, the mobile phase is less polar than the stationary phase. Inversely, in reversed mode, the mobile phase (*usually a mixture of water and organic solvent*) is more polar than the stationary phase (*C18*).

Layer Thickness

The layer thickness is related to the nature of the analysis (*analytical or preparative*) as well as the performance of the plate (*TLC or HPLTC*). The most common layer thicknesses are 150 μm (*HPTLC plates*), 200 - 250 μm (*analytical TLC plates*) and 500 - 2,000 μm (*preparative TLC plates*).

Binder & UV Indicator


All standard SiliaPlate products are made with a Gypsum binder and have an UV indicator (*F254*). Contact us for custom products.

Plate Size

SiliaPlate TLC plates are available in the following standard sizes depending on the coating used: 20 x 20 cm, 10 x 20 cm, 5 x 20 cm, 5 x 10 cm & 10 x 10 cm. Also for your convenience, SiliCycle provides ready to use micro TLC plates in 2.5 x 10 cm, 2.5 x 7.5 cm & 2.5 x 5 cm formats.

An interesting compromise between standard and micro plate sizes is our Scored SiliaPlate (*glass backing*). Two different formats are available and possible cut combinations are shown in the image below.

- 20 x 20 cm plates scored to get four 5 x 20 cm plates (*or multiple of 5 cm width*).
- 10 x 20 cm plates scored to get seven 2.5 x 10 cm plates (*or multiple of 2.5 cm width*).






SiliaPlate TLC Plates Portfolio

SiliCycle offers the possibility to analyze reactions on thin layer chromatography support and rapidly develop optimized purification conditions for efficient transfer to flash columns packed with the same SiliaFlash silica support. Maximize the benefits by using our *UltraPure* SiliaPlate TLC plates with an extra hard layer of silica. For your convenience, SiliCycle offers different sizes, choice of backings, reversed-phase & specialty plates. Contact us for more information.

Various attribute combinations are possible with SiliaPlate TLC plates and are summarized in the table below.

SiliaPlate TLC Plates Variety of Attributes			
Properties	TLC		HPTLC 
	Analytical	Preparative	
Backing Available			
Glass	Yes	Yes	Yes
Aluminum	Yes	No	No
Plastic	Yes	No	No
Adsorbent Available			
Bare silica	Yes	Yes	Yes
Silica - functionalized	No	Yes	Yes
Silica Specifications			
Mean Particle Size	10 - 14 μm	20 - 25 μm	$\leq 10 \mu\text{m}$
Mean Pore Diameter	60 \AA	60 \AA	60 \AA
Type of Plate Available			
Scored Plate Available	Yes	Yes	No
Channeled Plate Available	Yes	No	No
Layer Thickness	Glass: 250 μm Flexible: 200 μm	Glass: 500 μm , 1,000 μm Flexible: 1,500 μm , 2,000 μm	Glass: 150 μm
Typical Possible Plate Size*	2.5 x 5 cm; 2.5 x 7.5 cm; 2.5 x 10 cm; 5 x 10 cm; 5 x 20 cm; 10 x 20; 20 x 20 cm	20 x 20 cm	2.5 x 5 cm; 2.5 x 7.5 cm; 2.5 x 10 cm; 5 x 10 cm; 5 x 20 cm; 10 x 20; 20 x 20 cm

*For the glass-backing TLC plates.

« We had tried working with TLC plates of another brand and realized that the SiliCycle brand was the most durable and long-lasting as well as clear when visualizing with UV light so we switched back. »

Jessica Kisunzu from UC Berkeley, Berkeley, CA, USA

SiliaPlate Ordering Information.

SiliaPlate TLC with Glass Backing				
SiliCycle PN	Product Name	Plate Size (cm)	Thickness (µm)	#/box
Analytical SiliaPlate Glass				
TLG-R10011B-423	Micro SiliaPlate Glass	2.5 x 5	250	25
TLG-R10011B-124	Micro SiliaPlate Glass	2.5 x 7.5	250	100
TLG-R10011B-2575B	Micro SiliaPlate Glass (<i>bulk</i>)	2.5 x 7.5	250	384
TLG-R10011B-624	Micro SiliaPlate Glass	2.5 x 10	250	100
TLG-R10011B-527	SiliaPlate Glass	5 x 10	250	200
TLG-R10011B-424	SiliaPlate Glass	5 x 20	250	100
TLG-R10011B-723	SiliaPlate Glass	10 x 20	250	25
TLG-R10011B-2020	SiliaPlate Glass	20 x 20	250	20
TLG-R10011B-323	SiliaPlate Glass	20 x 20	250	25
Scored Analytical SiliaPlate Glass				
TLGSR10011B-723	SiliaPlate Glass (<i>scored</i>)	10 x 20	250	25
TLGSR10011B-423	SiliaPlate Glass (<i>scored</i>)	20 x 20	250	25
Channeled Analytical SiliaPlate Glass (with Preadsorbent Zone)				
TLGCZ-R10011B-323	Channeled SiliaPlate Glass	20 x 20	250	25
Preparative SiliaPlate Prep (Glass Preparative)				
TLG-R10011B-333	SiliaPlate Prep	20 x 20	500	25
TLG-R10011B-341	SiliaPlate Prep	20 x 20	1,000	25
TLG-R10011B-353	SiliaPlate Prep	20 x 20	2,000	25
Scored SiliaPlate Prep (Glass Preparative)				
TLGSR10011B-333	SiliaPlate Prep Glass (<i>scored</i>)	20 x 20	500	25
TLGSR10011B-341	SiliaPlate Prep Glass (<i>scored</i>)	20 x 20	1,000	25
TLGSR10011B-350	SiliaPlate Prep Glass (<i>scored</i>)	20 x 20	2,000	25
SiliaPlate Prep C18 (Glass Preparative)				
TLG-R30411B-341	SiliaPlate C18 Prep Glass	20 x 20	1,000	25

SiliaPlate TLC with Flexible Backings				
SiliCycle PN	Product Name	Plate Size (cm)	Thickness (µm)	#/box
SiliaPlate Al (Aluminum)				
TLA-R10011B-2575	Micro SiliaPlate Aluminum	2.5 x 7.5	200	200
TLA-R10011B-323	SiliaPlate Aluminum	20 x 20	200	25
SiliaPlate Al C18 (Aluminum)				
TLA-R30411B-303	SiliaPlate Aluminum C18	20 x 20	200	25
SiliaPlate PI (Plastic)				
TLP-R31001B-2575	Micro SiliaPlate Plastic	2.5 x 7.5	200	200
TLP-R31001B-323	SiliaPlate Plastic	20 x 20	200	25



NEW

SiliaPlate HPTLC Silica with Glass Backing (Thickness: 150 microns, 25 plates/box)


SiliCycle PN	Plate Size (cm)	SiliCycle PN	Plate Size (cm)
SiliaPlate Silica HPTLC			
HPTLG-R10011B-1010	10 x 10	HPTLG-R10011B-2020	20 x 20
Scored SiliaPlate Silica HPTLC			
HPTLGSR10011B-1010	10 x 10 (to 5 x 5 cm)	HPTLGSR10011B-1020	10 x 20

Functionalized SiliaPlate HPTLC (Thickness: 150 microns, 25 plates/box)

SiliCycle PN	Plate Size (cm)	SiliCycle PN	Plate Size (cm)
SiliaPlate C18 HPTLC			
TLG-R30411B-213	10 x 10	TLG-R30411B-303	20 x 20
SiliaPlate C8 HPTLC			
TLG-R31011B-203	10 x 10	TLG-R31011B-303	20 x 20
SiliaPlate C2 HPTLC			
TLG-R32611B-203	10 x 10	TLG-R32611B-303	20 x 20
SiliaPlate NH₂ (Amine) HPTLC			
TLG-R52011B-203	10 x 10	TLG-R52011B-303	20 x 20
SiliaPlate CN (Cyano) HPTLC			
TLG-R38011B-203	10 x 10	TLG-R38011B-303	20 x 20
SiliaPlate Diol HPTLC			
TLG-R35011B-203	10 x 10	TLG-R35011B-303	20 x 20
SiliaPlate Ag (Silver Nitrate 10% impregnated) HPTLC			
TLG-R23511B-423	5 x 20	TLG-R23511B-303	20 x 20

Trial Package of Functionalized SiliaPlate HPTLC: TLG-R1234511B-723

[5 plates (10 x 20 cm) of each SiliaPlate C18, C8, C2, NH₂ & CN scored to 2.5 x 10 cm]



SiliaPlate™ TLC Accessories

SiliaPlate™ TLC Developing Chamber

The most commonly used accessories to develop a TLC plate.

AUT-0160 SiliaPlate™ Cylinder Micro TLC Developing Chamber (1/box)

AUT-0161 SiliaPlate™ Rectangular TLC Developing Chamber (1/box)

Other SiliaPlate™ TLC Accessories

AUT-0162 SiliaPlate™ TLC Cutter

AUT-0163 SiliaPlate™ TLC Spotting Capillary Tubes

AUT-0164 SiliaPlate™ TLC Spotting Guide



An Ideal Partnership Between SiliCycle and Society AR2i

SiliCycle has entered into an exclusive world strategic specialty distribution partnership with the Society AR2i specialized in the conception and the manufacturing of innovative devices in the field of Thin-Layer Chromatography since 1994.

Chromimage® Documentation

- Perform qualitative analysis on TLC plates in a few minutes.
- Detect and numerize your TLC plates under UV254 nm and visible mode.
- Classify and archive your TLC analysis under several storage formats (*jpg, eps, pdf, etc.*).
- Suitable for reading 10 x 10 cm, 10 x 20 cm and 20 x 20 cm plates.



PN: AUT-0165

Derivapress® System

It's simple as opening and closing a book: the Derivapress immersion derivatization system provides a cost-effective, efficient and safe alternative to perfect this essential stage of TLC and to move towards densitometric measurements like quantitative and semi-quantitative TLC.

Furthermore Derivapress complies with the GLP requirements and can be used in 21 CFR Part 11 work environments.



PN: AUT-0166



Thin Layer Chromatography Practical Guide

Select a Stationary Phase

As almost 80% of all separations can be performed using silica gel plates, it is suggested to try using this coating. However, for acid sensitive compounds, alumina is probably a better choice (*useful for amine purification*). If you are working with highly polar compounds, reversed-phase mode is more suitable.

Select a Mobile Phase (*Solvent Systems*)

The selection of the mobile phase (*also called solvent system or eluent*) is perhaps the most important parameter to achieve efficient thin layer chromatography separation. It is based on the compound's solubility with the solvent and the difference in the affinity for the mobile phase versus the stationary adsorbent (*silica*).

In normal phase chromatography, where non-polar solvents such as hexane or pentane are used, non-polar compounds will move up the plate while most polar compounds will stay on the baseline. Inversely, polar solvents will allow polar compounds to move off the origin. The most suitable solvent system is the one that moves all components off the baseline with R_f values between 0.15 and 0.85 (*ideally, close to 0.2 - 0.4*).

Remember that it is not always possible in TLC but should be possible in flash chromatography where solvent gradients can be used.

For most applications, a common solvent system to start with is 1:1 Ethylacetate (*EtOAc*) / Hexane. Varying the ratio can have a pronounced effect on the R_f . If it is not working, then try: Methanol (*MeOH*) / Dichloromethane (*DCM*) (1:99 - 10:90); or toluene with acetone, EtOAc, or DCM.

Remember: To increase the compound's R_f , increase the polarity of the mobile phase; increase the ratio of the polar solvent or choose another solvent. Inversely, to decrease R_f , decrease the polarity of the eluent.

Rules of Thumb

- **Standard compounds (*most popular solvent system*):** 10 - 50% EtOAc/Hexane
- **Polar compounds:** 100% EtOAc or 5 - 10% MeOH/DCM
- **Non-polar compounds:** 5% EtOAc (*or ether*) / Hexane or 100% Hexane
- **For basic compounds:** (*amine or nitrogen containing*), it could be useful or required to add a small quantity of triethylamine (*Et₃N*) to the solvent mixture (0.1 - 2.0% but typical quantity is 0.1%) or 1 - 10% ammonia (*NH₃*) in MeOH/DCM.
- **For acidic compounds:** it could be useful to add acetic (*AcOH*) or formic acid (*FA*) to the solvent mixture (0.1 - 2.0%).

Reversed-phase mode

In reversed-phase chromatography, the typical solvent systems are:

- Mixtures of water or aqueous buffers and water miscible organic solvents such as acetonitrile (*ACN*), methanol, and tetrahydrofuran (*THF*). Other solvents can be used such as ethanol (*EtOH*) & isopropanol (*IPA*).
- If needed, to improve peak shape in flash chromatography, 0.1% of acetic, formic or trifluoroacetic acid (*TFA*) can be added to the solvent system.

« Have given your products to other folks within organisation and used it myself with great success (both the Prep SPE, HPLC columns, TLC plates, and silica gel).. »

Kerry M. Keertikar, Merck Research Labs, Kenilworth, NJ, USA

TLC Plate Preparation

Using a pencil, lightly draw a straight-line parallel to the width of the plate at about 1 cm from the base end of the plate. Sample application will be done on this line called baseline or origin.

Note: Never use a pen because ink can move with some solvents used as eluent.

Sample preparation


Thorough sample preparation is a prerequisite for an optimal and efficient TLC separation. Typical sample preparation processes could consist in a sample crushing, filtration, extraction or concentration of the product of interest.

Sample Application

Sample preparation will differ depending on the nature of the plate (*analytical or preparative*). For analytical plates, because thin layer chromatography is extremely sensitive, it is really important to apply a small quantity using a glass capillary (*or a micro pipette*) to get optimal resolution. For preparative plates, apply a series of small adjacent spots to form a band or a streak using a glass capillary (*or a microliter syringe*). In both cases, a spotting guide can be used to facilitate sample application.

Co-spotting

For analytical chromatography, co-spotting is frequently used for similar polarity products. This consists to apply on the same spot, the starting material and reaction mixture as shown by the image below.



TLC Plate Development

The most commonly used method to perform thin layer chromatography separation is to place vertically the TLC plate inside a sealed developing chamber to ensure solvent saturation. Place approximately 0.5 cm of the suitable solvent system inside the chamber. Slowly place the TLC inside the chamber and allow the eluent to travel up the plate until it gets to 1 cm from the top of the plate. Immediately remove the plate and draw a line along the solvent front.

Note: for optimal solvent saturation, a filter paper can be added inside the TLC chamber. This also prevents eluent evaporation. The solvent level needs to be below the baseline; otherwise the spots will be dissolved.



TLC Plate Visualization

If components of the reaction are colored, no visualization method is required (*spots can be seen directly on the silica layer*). However, most of the time it is not the case, therefore one of the methods described below should be used to reveal the spots.

Non-destructive methods

As a general visualization procedure, before treating the TLC plate with any destructive methods, UV-active compounds can be viewed under an ultraviolet lamp (*usually for polyconjugated compounds like benzophenones and anthracenes*). Furthermore, an iodine chamber can be useful for thiols, phosphines, and alkenes but it works in about 50% of cases for alkanes. It is recommended to circle the spots with a pencil on the TLC plate prior to visualization by destructive methods.

Destructive methods

For compounds that are not UV-active, there are several varieties of stains that can be used depending on the nature of the compound of interest. To use a stain, simply dip the TLC plate into the staining solution as quickly as possible, and then immediately absorb the excess stain with paper and heat carefully with a heat gun or on a hot plate at 110°C until spots are revealed. See next two pages.

Chromatogram Interpretation

Retention factor (*R_f*) definition


Retention factor analysis is used to evaluate if the solvent system is adequate. *R_f* is defined as the distance traveled by the compound divided by the distance traveled by the solvent front. This means: the larger the *R_f* value of a compound, the larger is the distance traveled by the compound. In other words, when comparing *R_f* values of various compounds under identical chromatography conditions, the compound with the larger *R_f* is less polar because it interacts less strongly with the polar adsorbent on the plate.

Remember, a good solvent system is one that moves all components off the baseline with *R_f* values between 0.15 and 0.85 (*ideal R_f is 0.2 - 0.4*). Otherwise, when possible, it is preferable to choose another solvent system.

$$\text{Retention factor (Rf)} = \frac{\text{distance traveled by the compound}}{\text{distance traveled by the solvent front}}$$

R_f calculation based on the example shown here:

$$R_f = 4.0 \text{ cm} / 5.5 \text{ cm} = 0.73$$



Prediction of Column Volumes (CV)

TLC data can be used to predict column elution based on the relationship between the retention factor and the column volume. CV is the number of column volumes required to elute the component from the column regardless of column dimensions [(*bed volume*) - (*volume of packing*)].

$$CV = 1 / R_f \quad \& \quad \Delta CV = 1 / R_{f_1} - 1 / R_{f_2}$$

The greater the ΔCV , the greater will be the separation and resolution between the spots (easier separation). A bigger ΔCV will therefore allow more sample to be loaded onto the column.

Described below, are the most frequently used TLC visualization methods (also called stains) in alphabetical order.

Stains for Thin Layer Chromatography			
Name	Visualization of...	Stain Recipe	Comments
<i>p</i>-Anisaldehyde #1	Universal Stain <i>Good for nucleophiles and oxygenated compounds</i>	Prepare stain as follows <ul style="list-style-type: none"> • 2 mL of glacial acetic acid • 5 mL of <i>p</i>-anisaldehyde • 7 mL of conc. sulfuric acid • 185 mL 95% ethanol <i>Tip: Add dropwise the acid at the end and stir vigorously.</i>	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: Orange to pink Appropriate Storage <ul style="list-style-type: none"> • Aluminum wrapped at 0°C
N.B.: Tends to be insensitive to alkenes, alkynes and aromatic compounds unless other functional groups are present.			
<i>p</i>-Anisaldehyde #2	Acronycine Cineoles Terpenes	Prepare stain as follows [1:10:20:80] <ul style="list-style-type: none"> • <i>p</i>-anisaldehyde • perchloric acid • acetone • water 	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: Orange to pink Appropriate Storage <ul style="list-style-type: none"> • Aluminum wrapped at 0°C
Bromocresol Green	Acidic groups (pK_a < 5) Carboxylic Acids	Prepare stain as follows <ul style="list-style-type: none"> • 0.04 g of bromocresol green • 100 mL of 95% ethanol • 0.1 M solution of sodium hydroxide <i>Tip: Add the base slowly at the end until the solution turns pale blue.</i>	Visualization Colors <ul style="list-style-type: none"> • Spots: Yellow to green • BG: Blue Appropriate Storage <ul style="list-style-type: none"> • Aluminum wrapped at 0°C Heating NOT required
Cerium Molybdate (CAM or Hanessian's Stain)	Universal Stain <i>Good for peptides</i>	Prepare stain as follows <ul style="list-style-type: none"> • 12 g of ammonium molybdate • 0.5 g of ceric ammonium molybdate • 15 mL of conc. sulfuric acid • 235 mL of water 	Visualization Colors <ul style="list-style-type: none"> • Spots: Blue • BG: White Appropriate Storage <ul style="list-style-type: none"> • Aluminum wrapped
N.B.: Highly sensitive stain; very low concentration of product may appear as a significant impurity.			
Cerium Sulfate (Ce(SO₄)₂)	Difficultly stainable compounds	Prepare stain as follows <ul style="list-style-type: none"> • 15% aqueous sulphuric acid saturated with ceric sulfate 	Visualization Colors <ul style="list-style-type: none"> • Spots: Black • BG: Yellow to white
Chromic Acid	Difficultly stainable compounds	Prepare stain as follows <ul style="list-style-type: none"> • 2.5 g of potassium chromate • 100 mL of 20% sulfuric acid in water 	Visualization Colors <ul style="list-style-type: none"> • Spots: Orange to green • BG: Yellow to red
Cobalt Chloride (CoCl₂)	Universal stain <i>Used in conjunction with PMA when this one is not effective enough</i>	Prepare stain as follows <ul style="list-style-type: none"> • 2 g of cobalt chloride • 100 mL of water • 10 mL of conc. sulfuric acid <i>Tip: simply dip PMA treated plate in CoCl₂ solution.</i>	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: Pink Heating NOT required
<i>p</i>-Dimethylamino-benzaldehyde (PDAB or Ehrlich's Reagent)	Amines Indoles	Prepare stain as follows <ul style="list-style-type: none"> • 0.5 g of <i>p</i>-dimethylamin-benzaldehyde • 10 mL of conc. hydrochloric acid • 40 mL of acetone (or 95% ethanol) 	Visualization Colors <ul style="list-style-type: none"> • Blue



Stains for Thin Layer Chromatography (Con't)

Name	Visualization of...	Stain Recipe	Comments
2,4-Dinitrophenyl-hydrazine (DNP)	Aldehydes Ketones	Prepare stain as follows <ul style="list-style-type: none"> • 12 g of 2,4-dinitrophenylhydrazine • 60 mL of conc. sulfuric acid • 80 mL of water • 200 mL of 95% ethanol 	Visualization Colors <ul style="list-style-type: none"> • Spots: Yellow to red • BG: Light orange DO NOT HEAT dipped plate
Dragendorff Reagent	Nitrogenous Compounds <i>Alkaloids, amines, organics bases, etc.</i> Phenols	Prepare stain as follows Solution A <ul style="list-style-type: none"> • 1.7 g of bismuth nitrate • 80 mL of water • 20 mL of acetic acid Solution B <ul style="list-style-type: none"> • 40 g of potassium iodide • 100 mL of water <i>Tip: mix 5 mL of each solution A and B to a solution of 20 mL of acetic acid in 70 mL of water.</i>	Visualization Colors <ul style="list-style-type: none"> • Spots: Orange to red • BG: Yellow Appropriate Storage <ul style="list-style-type: none"> • Aluminum wrapped Stain Shelf-Life <ul style="list-style-type: none"> • One or two weeks • Solutions A and B are long term storable DO NOT HEAT dipped plate
Ferric Chloride (FeCl ₃)	Phenols	Prepare stain as follows <ul style="list-style-type: none"> • 2 g of ferric chloride • 102 mL of 0.5N hydrochloric acid 	Visualization Colors <ul style="list-style-type: none"> • Spots: Red • BG: Yellow
Iodine	Unsaturated & Aromatic Compounds	Prepare stain as follows <ul style="list-style-type: none"> • Iodine crystals in an amber bottle 	Visualization Colors <ul style="list-style-type: none"> • Spots: Dark brown • BG: Light brown
N.B.: iodine stain can be removed by heating.			
Morin Hydrate (Hydroxy Flavone)	Universal stain <i>Fluorescently active</i>	Prepare stain as follows <ul style="list-style-type: none"> • 0.1% of morin hydrate in methanol <i>Tip: by weight.</i>	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: White
Ninhydrin (Indanetrione Hydrate)	Amino Acids Amino Sugars Amines	Prepare stain as follows <ul style="list-style-type: none"> • 1.5 g of ninhydrin • 3 mL acetic acid • 100 mL of <i>n</i>-butanol 	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: White
Phosphomolybdic Acid (PMA)	Universal stain <i>Very effective against dilute sample</i>	Prepare stain as follows <ul style="list-style-type: none"> • 10% of PMA solution in ethanol • or 10 g of PMA in 100 mL ethanol 	Visualization Colors <ul style="list-style-type: none"> • Spots: Dark green to black • BG: Light green
Potassium Permanganate (KMnO ₄)	Olefins Readily oxidized groups <i>Alcohols, aldehydes, alkenes, alkynes, etc.</i>	Prepare stain as follows <ul style="list-style-type: none"> • 1.5 g of potassium permanganate • 10 g of potassium carbonate • 1.25 mL of 10% sodium hydroxyde • 200 mL of water 	Visualization Colors <ul style="list-style-type: none"> • Spots: Yellow to light brown • BG: Purple to pink Stain Shelf-Life <ul style="list-style-type: none"> • Three months
N.B.: Can be used for detection of alcohols, amines, sulfides and mercaptans groups when gently heated.			
Vanillin	Universal stain <i>Very effective for some polarity products (Rf)</i>	Prepare stain as follows <ul style="list-style-type: none"> • 15 g of vanillin • 250 mL of 95% ethanol • 2.5 mL of conc. sulfuric acid 	Visualization Colors <ul style="list-style-type: none"> • Spots: Various colors • BG: Light tan

(N.B.: Shaded lines refer to "universal stains". Occasionally, spots can be seen more clearly from glass side with glass backed TLC plate. Otherwise mentioned, stains are long-term stable when stored in a tightly-closed container to prevent solvent evaporation. "BG" stands for "background".)

SiliaPlate TLC Troubleshooting

Problem: Streaking or elongated spot rather than a defined spot?

Possible Solutions:

- Sample was overloaded: run the TLC again using a more diluted solution of your sample.
- In presence of a base sensitive compound: try to add acetic or formic acid to the eluent (0.1 - 2.0%).
- In presence of an acid sensitive compound: try to add triethylamine to the eluent (0.1 - 2.0%) or 1 - 10% ammonia in MeOH/DCM. If it is not working use Alumina as TLC coating.
- In presence of too highly polar compounds: try using a specialized silica TLC plate like reversed-phase (C18 for example).

Problem: Unable to see any spots on the TLC?

Possible Solutions:

- If you have not been able to visualize any spots on your TLC using UV light, try another method; maybe your compound is not UV-active.
- Maybe your sample is too diluted. Try to apply several times your sample on the same spot (*do not forget to dry solvent between each application for optimal results*) or to concentrate your solution.
- Make sure the solvent level inside the tank is lower than the spotting line to avoid sample dissolution by the eluent.

Problem: How to monitor a reaction in presence of similar Rf's for both starting materials and product of interest?

Possible Solutions:

- Try the co-spotting method.
- Try to visualize the plate using anisaldehyde or molybdene. Spot color or brightness differ for two compounds when using these stains.
- If none of the two previous solutions work, change solvent systems (*use another class of solvent*).

Tips: *in chromatography, there are three classes of solvent systems providing significantly different results:*

- 1: Mixture of polar/hydrocarbon solvents (*i.e.: EtOAc/Hexane; Ether/Petroleum ether*).
- 2: Mixture of polar/dichloromethane solvents (*examples of polar solvent: Ether, EtOAc, MeOH*).
- 3: Mixture of polar/benzene (*or toluene*) solvents (*examples of polar solvent: Ether, EtOAc, MeOH*).

Problem: Compounds stay too close to the baseline or solvent front.


Possible Solutions:

- Too close to the baseline: your eluent is not polar enough; increase the proportion of polar solvent in the same solvent system or chose a more polar solvent.
- Too close to the solvent front: inversely, your eluent is too polar; decrease the proportion of polar solvent in the same solvent system or chose a less polar solvent.



SiliaSphere™

Spherical Silica Gels



SiliaSphere Spherical Silica Gel

SiliCycle is your partner of choice for your purification and chromatography needs. Here is why:

- Recognized worldwide as a leader with outstanding quality spherical silica gels for analytical and preparative chromatography.
- One of the best silica on the market with unbeatable performance.
- Broadest portfolio of chromatographic phases to suit all purification needs.
- Wide range of formats for all purification scales bulk & HPLC columns.



SiliCycle is a World Class Silica Expert

At SiliCycle, we offer two grades of spherical silica gel; SiliaSphere Monodispersed Spherical Silica Gel for analytical needs and SiliaSphere PC for preparative ones. Both grades are manufactured under the same well controlled processes to ensure constant reproducibility and easy scale-up.

With pore diameters ranging from 60 to 1,000 Å and particle sizes from 1.8 to 200 µm, we offer products to meet all your separation requirements. This is one of the most reliable portfolios of spherical silica gels for high pressure chromatography. SiliaSphere are ideal for both analytical and preparative chromatography, from laboratory to pilot-plant processes and production scales. Furthermore, the excellent properties of SiliaSphere make them the packing of choice for high performance liquid chromatography (HPLC), Supercritical Fluid Chromatography (SFC) and simulated moving bed (SMB) applications.

Features and Benefits of SiliaSphere Spherical Silica Gels

Features and Benefits of SiliaSphere Spherical Silica Gels	
Features	Benefits
High purity silica gels	Consistency, reliability & reproducibility
Perfect spherical shape free of any cavities or cracks	Ease of column packing and high resolution
Exceptional narrow particle size distribution	Optimal separation and resolution
Strong mechanical stability	Low back-pressure without surface abrasion
Same well controlled processes for all SiliaSphere	Easy scalability
Availability in bulk quantities at affordable price	On-time delivery



SiliaSphere as a Silica Matrix

For over 18 years, SiliCycle has been dedicated to silica gel manufacturing and we have achieved a strong know-how and expertise in this field. Furthermore, to support the increasing demand on our spherical silicas, we have developed an optimized and highly controlled large-scale production process for all of our SiliaSphere products without decreasing the quality of the silica.

Particle shape, pore and particle size distributions, silica gel purity, and surface properties, all have their influence on chromatographic performance. Therefore, in order to develop the most efficient process, all previously mentioned parameters need to be evaluated and optimized to ensure batch-to-batch reproducibility.


SiliaSphere is manufactured from an organic form of silicon (*alkoxydes*). This ensures very low metal content as the starting material is purified by distillation. Deionized water is used to hydrolyze the silicon alkoxydes. Careful monitoring and control of the parameters that induce precipitation afford spherical silica gels with the desired characteristics. SiliaSphere products are characterized by a very low metal content and exceptional stability even at extreme pHs. Furthermore, our manufacturing process ensures quality and reproducibility in pore size, surface area and particle sizes and morphology for all SiliaSphere products.

SiliaSphere, the right choice for superior...

- Chromatographic performance
- Loading capacity
- Reproducibility
- Chemical & physical stability

Perfectly Spherical Particle Shape


The perfectly spherical shape of SiliaSphere silicas, combined to their smooth surfaces free of any cracks, cavities and fines make them the packing of choice for chromatography. The SiliaSphere sphericity exceeds or compares favorably to well-known brands in spherical silica gel as demonstrated by the scanning electron microscope (SEM) picture.



SiliaSphere

Narrow Pore Size Distribution


The right pore diameter to use is related to the type of molecules present in the sample to be purified. Typically, for small molecules, 100 - 120 Å pore size is recommended. However, if higher resolution is required, then 60 - 80 Å will be more suitable. For large molecules, such as peptides and proteins, 300 Å or higher is recommended.



Tight Particle Size Distribution

The importance of the particle size distribution varies depending on the type of chromatography being done. For instance, it is very important when using HPLC that the particle size distribution of the spherical particles being used be very narrow. Tight particle size distribution yields greater column performance (*separation*), better peak shape, lower back-pressure as well as higher packing stability.

SiliaSphere offers one of the narrowest particle size distributions available on the market. To evaluate our particle size distribution, we use the D90/D10 ratio.




High Surface Area and Pore Volume

Our optimized manufacturing process ensures high specific surface area for greater loading capacity with an uniform and reproducible surface coverage.

Low Trace Metal Content & High Purity

SiliCycle's proprietary technology generates a silica gel with one of the lowest trace metal content on the market today. Our low trace metal content ensures you will get optimal performance chromatography. Tight control of trace metals in every batch also improves the reproducibility and reduces risks of interaction between metals and desired compounds. Low metal content limiting any unwanted metal ion solute interactions – providing symmetrical peaks with little or no tailing.

To probe low metal content in SiliaSphere silicas, we ran the following chromatographic test and we compared our products with a well-known competitor. With SiliaSphere, you can be assured that peak tailing or missing peaks are not coming from our silica.



- | | |
|----------------------------|---|
| 1. Toluene: | Column Packing Efficiency => Neutral solute |
| 2. Acetoacetanilide: | Metal Content => Can interact with heavy metals resulting in tailing or missing peaks |
| 3. Catechol: | Metal Content => Can interact with heavy metals resulting in tailing or missing peaks |
| 4. 3-Pyridyl-acetonitrile: | Silica's Acidic Character => Strongly sorbs to acidic sites |




Scalability

SiliaSphere are available in bulk for a wide range of HPLC purification; from laboratory to plant scale. All SiliaSphere products are manufactured under tightly controlled manufacturing processes, and stringent quality control ensures the highest quality and reproducibility possible. Scaling-up is extremely straight-forward with SiliaSphere silicas and performance will remain the throughout the range of particle sizes.

Lot-to-lot reproducibility

SiliCycle is committed to high quality standards and always strives to provide defect-free products and lot-to-lot reproducibility. In doing so, all products are manufactured in an ISO 9001:2008 compliant facility and subjected to a stringent quality control. Every lot needs to meet the established specifications to be released. This ensures reliable & reproducible chromatographic performances.

The graph at right presents results obtained for the molecular loading based on sulfur content as well as surface coverage of the last 5 lots of SiliaSphere 10 μm, 120 Å produced which shows high product manufacturing reproducibility.



SiliaSphere Characteristics & Ordering Information

SiliaSphere Monodispersed Characteristics

Properties / Pore Diameter	60 & 80 Å	100 Å	120 Å	300 Å	1,000 Å
Specific Surface Area (m ² /g)	≥ 450	≥ 400	≥ 300	≥ 80	≥ 20
Pore Volume (mL/g)	0.85 - 1.15	0.85 - 1.15	0.85 - 1.15	0.75 - 1.05	0.75 - 1.05
pH (5% w/w)	4 - 7	4 - 7	4 - 7	4 - 7	4 - 7
Available Particle Sizes (μm)	1.8', 2.2', 3, 5, 10, 15, 20	1.8, 2.2, 3, 5, 7, 10, 15, 20	1.8, 2.2, 3, 5, 7, 10, 15, 20	3, 5, 10, 15	10, 15

*Not available in 60 Å

SiliaSphere Monodispersed Spherical Silica Gels Product Number

Particle Size (μm)	Pore Diameter (Å)					
	60 Å	80 Å	100 Å	120 Å	300 Å	1,000 Å
BARE SiliaSphere Monodispersed Spherical Silica						
1.8	n/a	S10001F-A	S10001E-A	S10001G	n/a	n/a
2.2	n/a	S10002F-A	S10002E-A	S10002G	n/a	n/a
3	S10003B	S10003F-A	S10003E-A	S10003G-A	S10003M	n/a
5	S10005B	S10005F-A	S10005E-A	S10005G-A	S10005M	n/a
7	n/a	n/a	S10006E-A	S10006G-A	n/a	n/a
10	S10007B	S10007F-A	S10007E-A	S10007G-A	S10007M	S10007T
15	S10008B	S10008F-A	S10008E-A	S10008G-A	S10008M	S10008T
20	n/a	S10009F-A	S10009E-A	S10009G-A	n/a	n/a

Formats : 100 g, 500 g, 1 kg, 5 kg, 10 kg, 25 kg, ... up to multi-ton scale. Refer to page 75 for bonded phases information.

SiliCycle; Experts in Functionalization

Every year, various companies release new chromatographic phases to the market and claim that their phases are the most efficient. Choosing the right phase along with the right supplier is becoming a challenging task. For the past 18 years, SiliCycle has been dedicated to the development of silica-based products and has acquired extensive experience in grafting technology. This outstanding expertise, combined to tightly controlled proprietary functionalization processes, allows SiliCycle to be the supplier of choice for all chromatographic applications.

Grafting Process

Silica surface possesses active silanols (*Si-OH, free OH groups of the silica*), which permits the modification of the surface chemistry by grafting silane moieties. This property allows the control of the surface polarity useful in separation techniques. Various types of silanes can be grafted on the surface to afford monomeric or polymeric bonded phases which both have pros and cons as described below.

Main Differences Between Monomeric and Polymeric Bonded Phases

Monomeric Functionalization

By grafting a monofunctional alkylsilane reagent, only one bond can be formed with the silica surface. This type of grafting is called monomeric. The dimethyl groups help to protect the surface by creating steric hindrance, which, at the same time, prevent from reaching the highest possible silane density. The residual silanol groups are inhibited by the grafting of a capping agent. Usually, even after the endcapping, a small portion of the initial silanols is still present due to steric hindrance.


Monomeric phases present a very high stability, batch-to-batch reproducibility and good hydrophobic properties. The fact that the silane possesses only one bond with the surface makes this phase less stable at low pH, which may lead to silane hydrolysis and consequently leaching. For low pH, the polymeric phase is preferred.

Polymeric Functionalization

Historically, trichloroalkylsilanes were the first to be employed for economic reasons. By grafting a di- or trifunctional alkylsilane reagent, it is possible to form multidimensional bondings with the surface and also between silane molecules. This grafting method is called polymeric functionalization. The silica surface is more hydrophobic, has greater stability in strong acidic condition (*pH 2-3*) and has a longer lifetime. However, polymeric phases present major disadvantages over monomeric ones: lower homogeneous surface coverage due to cross-polymerization reactions, poorer batch-to-batch reproducibility leading to variation in retention even for the same molecule. (*Note: see picture next page*).



Monomeric vs Polymeric Bonded Phases



Homogeneous Surface Coverage

SiliCycle developed a new and innovative grafting technique which is characterized by a homogeneous coverage of the alkyl chains on the surface. This proprietary process can be used with all silane types and ensures greater chemical stability as well as better performance due to the greater homogeneity of the surface coverage.

Endcapping

When functionalizing silica, it is impossible to react with every silanol group, so endcapping technology is often used to prevent peak tailing caused by non-specific interactions and thus improve separations. Furthermore, more sophisticated methods lead to strong layer protection offering very high sorbent durability in harsh conditions.

The endcapping step can be done using various methodologies. The easiest one is to treat the surface with a small silylating agent, such as trimethylchlorosilane (*TMSCl*). However, at SiliCycle, we always try to improve and control this critical step to afford highly deactivated silanol phases. For some phases, we use the conventional single endcapping step technique, for others we use our proprietary endcapping processes which can include multistep methods, use of specific silylating agents or other special treatments.

Most Popular Bonded Phases

In liquid chromatography, there are various modes of operation possible based on the interaction mechanism of the solute with the stationary phase (*sorbent*). Most known separation modes are summarized in the table below of which reversed-phase is the most popular one.

SiliaSphere Most Popular Bonded Phases			
Mode	Normal (NP)	Reversed (RP)	Ion Exchange (IEX)
Mode Mechanism	Polar or hydrophilic	Non-polar or lipophilic	Ionic
Typical Stationary Phase	Bare silica or functionalized silica (<i>Amine, Cyano or Diol</i>)	Functionalized silica (<i>mostly C18, C8, C4, Cyano, Phenyl and PFP</i>)	Ionic functionalized silica (<i>SAX, SCX</i>) or polymer
Stationary Phase Polarity	Polar	Non-polar	Anionic or cationic exchanger
Typical Mobile Phase	Non-polar organic solvents such as hexane, dichloromethane, THF	Mixtures of water or aqueous buffers and organic solvents (<i>mostly ACN, MeOH, THF</i>). Ion pairing agents can also be added.	Water, buffers; acid; base
Stationary Phase Polarity	Non-polar	Polar	Buffer or ionic

Typical Applications of Most Common Bonded Phases

Typical Applications of Most Common Bonded Phases					
Sorbent Phase	Functional Group	Mode			Typical Applications
		NP	RP	IEX	
Silica	- OH	X			Most polar sorbent with a slight acidic character used for purification of non-polar and non-ionic compounds.
C18	- (CH ₂) ₁₇ CH ₃		X		Great start for method development. Presents the maximum retention of non-polar compounds. Usually used for peptides, pesticides, PCBs, PAHs, drugs, etc.
C8	- (CH ₂) ₇ CH ₃		X		Presents less retention compared to C18. Usually used for highly hydrophobic pesticides, small peptides and heavy drugs.
C4	- (CH ₂) ₃ CH ₃		X		Presents less retention compared to C18 and C8. Widely used for molecules with large hydrophilic regions such as peptides, proteins and zwitterions (<i>in 300 Å</i>).
C1	- (CH ₃) ₃		X		Lower retention compared to other reversed-phases used for the purification of polar and non-polar pharmaceutical products, highly hydrophobic molecules.
Cyano (CN)	- (CH ₂) ₂ CN	X	X		Normal: less polar sorbent compared to silica used for the purification of polar organic compounds. Reversed: Moderate non-polar sorbent with less hydrophobicity than C18 or C8. Purification of cyclosporine and carbohydrates.
Phenyl (PHE)	- C ₆ H ₅		X		Moderate non-polar sorbent with different selectivity for aromatic compounds compared to other non-polar sorbents.
Pentafluorophenyl (PFP)	- (CH ₂) ₃ C ₆ F ₅		X		For a new selectivity approach or the purification of conjugated compounds (isomers).
Amine (NH ₂)	- (CH ₂) ₃ NH ₂	X		X	Normal: polar sorbent with a basic character with less retention and a different selectivity for acidic/basic compounds compared to silica. Ion Exchange: A weak anion exchanger with pK _a of 9.8. At pH 7.8 or below, the functional groups are positively charged. It facilitates the rapid release of very strong anions such as sulfonic acids that may be retained irreversibly on SAX.
Diol	- (CH ₂) ₃ OCH ₂ CH(OH)CH ₂ OH	X			Moderate polar sorbent with a neutral character used to extract polar compounds. Alternative to silica when acidic character is problematic.
Toxic Acid (SCX)	- (CH ₂) ₂ C ₆ H ₄ SO ₃ H			X	Due to the very low pK _a (< 7), this silica is a strong cation exchanger. The most common use is likely for catch and release purification.
TMA Chloride (SAX)	- (CH ₂) ₃ N ⁺ (CH ₃) ₃ Cl ⁻			X	The quaternary amine is permanently charged and commonly used for the extraction of weak cations that may not bind strongly enough to weaker anion exchangers.
TMA Acetate (SAX-2)	- (CH ₂) ₃ N ⁺ (CH ₃) ₃ (CO ₂ CH ₃) ⁻			X	The acetate counter ion is easier to exchange compared to the chloride ion. It is used for compounds with pK _a < 5, such as carboxylic acids or to selectively purify acidic compounds or remove acidic impurities from reaction mixtures.



Bonded SiliaSphere Spherical Silica Gels Ordering Information

The table below presents the most popular SiliaSphere bonded phases available from SiliCycle. Please note that we can provide phases below on any of our other silica gels and all other sorbents, contact us.

Bonded SiliaSphere Spherical Silica Gels Ordering Information						
Particle Size (µm)	Pore Diameter (Å)					
	60 Å	80 Å	100 Å	120 Å	300 Å	1,000 Å
C18 SiliaSphere Monodispersed Spherical Silica (endcapped)						
1.8	n/a	S03201F-A	S03201E-A	S03201G	n/a	n/a
2.2	n/a	S03202F-A	S03202E-A	S03202G	n/a	n/a
3	S03203B	S03203F-A	S03203E-A	S03203G-A	S03203M	n/a
5	S03205B	S03205F-A	S03205E-A	S03205G-A	S03205M	n/a
7	n/a	n/a	S03206E-A	S03206G-A	n/a	n/a
10	S03207B	S03207F-A	S03207E-A	S03207G-A	S03207M	S03207T
15	S03208B	S03208F-A	S03208E-A	S03208G-A	S03208M	S03208T
20	n/a	S03209F-A	S03209E-A	S03209G-A	n/a	n/a

Formats : 100 g, 500 g, 1 kg, 5 kg, 10 kg, 25 kg, ... up to multi-ton scale.

Bonded SiliaSphere Spherical Silica Gels Ordering Information						
Particle Size (µm)	Pore Diameter (Å)					
	60 Å	80 Å	100 Å	120 Å	300 Å	1,000 Å
C8 SiliaSphere Monodispersed Spherical Silica (endcapped)						
1.8	n/a	S30801F-A	S30801E-A	S30801G	n/a	n/a
2.2	n/a	S30802F-A	S30802E-A	S30802G	n/a	n/a
3	S30803B	S30803F-A	S30803E-A	S30803G-A	S30803M	n/a
5	S30805B	S30805F-A	S30805E-A	S30805G-A	S30805M	n/a
7	n/a	n/a	S30806E-A	S30806G-A	n/a	n/a
10	S30807B	S30807F-A	S30807E-A	S30807G-A	S30807M	S30807T
15	S30808B	S30808F-A	S30808E-A	S30808G-A	S30808M	S30808T
20	n/a	S30809F-A	S30809E-A	S30809G-A	n/a	n/a

Formats : 100 g, 500 g, 1 kg, 5 kg, 10 kg, 25 k g, ... up to multi-ton scale.

« I was able to separate two very close diastereomers that were inseparable using reverse phase HPLC. »

Jeff Warrington from Cytokinetics, South San Francisco, CA, USA

Bonded SiliaSphere Spherical Silica Gels Ordering Information

The table below presents the most popular SiliaSphere bonded phases available from SiliCycle. Please note that we can provide phases below on any of our other silica gels and all other sorbents, please contact us.

Bonded SiliaSphere Spherical Silica Gels Ordering Information						
Particle Size (µm)	Pore Diameter (Å)					
	60 Å	80 Å	100 Å	120 Å	300 Å	1,000 Å
Amine SiliaSphere Monodispersed Spherical Silica (endcapped)						
1.8	n/a	S52001F-A	S52001E-A	S52001G	n/a	n/a
2.2	n/a	S52002F-A	S52002E-A	S52002G	n/a	n/a
3	S52003B	S52003F-A	S52003E-A	S52003G-A	S52003M	n/a
5	S52005B	S52005F-A	S52005E-A	S52005G-A	S52005M	n/a
7	n/a	n/a	S52006E-A	S52006G-A	n/a	n/a
10	S52007B	S52007F-A	S52007E-A	S52007G-A	S52007M	S52007T
15	S52008B	S52008F-A	S52008E-A	S52008G-A	S52008M	S52008T
20	n/a	S52009F-A	S52009E-A	S52009G-A	n/a	n/a

Formats : 100 g, 500 g, 1 kg, 5 kg, 10 kg, 25 kg, ... up to multi-ton scale.

Bonded SiliaSphere Spherical Silica Gels Ordering Information						
Particle Size (µm)	Pore Diameter (Å)					
	60 Å	80 Å	100 Å	120 Å	300 Å	1,000 Å
Cyano SiliaSphere Monodispersed Spherical Silica (endcapped)						
1.8	n/a	S38001F-A	S38001E-A	S38001G	n/a	n/a
2.2	n/a	S38002F-A	S38002E-A	S38002G	n/a	n/a
3	S38003B	S38003F-A	S38003E-A	S38003G-A	S38003M	n/a
5	S38005B	S38005F-A	S38005E-A	S38005G-A	S38005M	n/a
7	n/a	n/a	S38006E-A	S38006G-A	n/a	n/a
10	S38007B	S38007F-A	S38007E-A	S38007G-A	S38007M	S38007T
15	S38008B	S38008F-A	S38008E-A	S38008G-A	S38008M	S38008T
20	n/a	S38009F-A	S38009E-A	S38009G-A	n/a	n/a

Formats : 100 g, 500 g, 1 kg, 5 kg, 10 kg, 25 kg, ... up to multi-ton scale.

« SiliCycle provided improved peak shape for an acidic compounds improving our yield off the column. »

Sean Brown from Amgen, South San Francisco, CA, USA



SiliaChrom[®]

HPLC Columns



SiliaChrom HPLC Columns

Using SiliaChrom HPLC Columns in chromatographic applications ensures the following:

- Excellent column efficiency.
- Long lifetime and column-to-column reproducibility.
- Broad pH range from 0.8 to 12.
- Compatibility with 100% aqueous and organic mobile phases.
- High surface coverage presenting no bleeding for LC-MS applications.



Presentation of the SiliaChrom HPLC Column Series

SiliCycle manufactures a variety of HPLC columns for reversed and normal phase applications. The SiliaChrom series contain more than 40 different phases, and we continue to develop additional, unique and powerful HPLC sorbents. Most of the SiliaChrom are silica-based products. You can be assured of the quality, from raw material synthesis through to the packing process.

We pack bonded phases in a wide range of column dimensions, including standard narrow bore and analytical columns in lengths of 20 to 250 mm, internal diameters (*ID*) of 2.0 - 4.6 mm, with particle sizes of 2.5, 3.0, 5.0, 10.0 or 20.0 μm . Also, preparative and semi-preparative HPLC columns are available, in 10, 20, 30 and 50 mm ID with lengths of 50, 100, 150 and 250 mm with particle sizes of up to 20 μm . These columns exhibit superior performances for any type of compound. The SiliaChrom

series, with its unique sol-gel process technology, offers the total solution for HPLC end-users: broad pH range (0.5 - 12), compatibility with 100% aqueous and organic mobile phases, low bleeding for LC-MS, high surface coverage, and excellent column efficiency. All columns are packed using a consistent proprietary packing process to achieve uniform and stable bed for long lifetime and column-to-column reproducibility.

The SiliaChrom HPLC portfolio offers a broad variety of separations for various types of chromatography such as biochromatography of large molecules, size exclusion chromatography for large proteins and peptides, chiral chromatography for enantioselective separations and supercritical fluid chromatography for API separations. The following pages will highlight SiliaChrom phases that can be used for these applications.



SiliCycle; Experts in HPLC Column Packing

Superior HPLC columns can be produced only with excellent packing materials and excellent packing techniques. SiliaChrom columns are made from extremely pure silicas and are well known for their high efficiency and high resolution capacity. Based on spherical, totally porous silica, SiliaChrom columns provide enhanced chemical and mechanical stability as well as very high loading capacity and full end-capping.

All SiliaChrom columns are packed using a proprietary slurry packing process to achieve uniform and column-to-column reproducibility. SiliaChrom columns have good selectivity, good asymmetry and long lifetime for HPLC separation of acidic, neutral and basic organic compounds, polar or non-polar.

Standard HPLC Columns

SiliaChrom HPLC columns are available in Narrow Bore, Analytical, Semi-Preparative, and Preparative formats.

SiliaChrom HPLC Standard Column Dimensions								
		SiliaChrom HPLC Column Length (mm)						
		20	30	50	100	150	200	250
Particle Size (µm)		2.5, 3 & 5	2.5, 3 & 5	2.5, 3, 5, 7 & 10	2.5, 3, 5, 7, 10 & 20	2.5, 3, 5, 7, 10 & 20	5 & 10	2.5, 3, 5, 7, 10 & 20
SiliaChrom HPLC Column Internal Diameter (mm)	2.0	S	S	S	S	S	C	C
	2.1	S	S	S	S	S	C	C
	3.0	C	S	S	S	S	C	C
	4.6	C	S	S	S	S	S	S
	10	C	C	C	S	S	S	S
	20	C	C	C	C	S	S	S
	30	C	C	C	S	S	S	S
	50	C	C	C	S	S	S	S
	100	C	C	C	C	C	C	C

S = Standard C = Custom

« SiliCycle has been able to repeatedly come through and produce high quality semi-prep HPLC columns (50+ mm ID) for several different projects that we have done. For many of these projects prices is not the driving force, the timing is.»

Jason Blanchard from Ricerca Biosciences, Concord, OH, USA

Column Packing Reproducibility

SiliCycle is recognized for its strong expertise in column packing technology. All SiliChrom columns are packed using a consistent packing methodology to achieve an extremely stable and uniform column packing bed leading to high column lifetime and column-to-column reproducibility. To prove this, we packed and tested several analytical columns 4.6 x 250 mm using the same SiliSphere C18 3 μm , 100 Å for reproducibility and high efficiency evaluation.


Chromatographic conditions

- **Sample mixture in mobile phase:** Uracil / Phenol / Nitrobenzene / Naphtalene
- **Injection volume:** 2 μL
- **Temperature:** 30°C
- **Flow rate:** 0.8 mL/min
- **Mobile phase:** 15% Water, 85% Methanol

Observed Column Parameters for Naphtalene			
Column Number	Retention Time (min)	Theoretical Plates Number / meter	Tailing Factor
1	9.148	28,481	1.01
2	9.382	28,391	1.00
3	9.398	28,712	1.00
4	8.998	28,150	1.01
5	9.307	28,393	1.00
6	9.307	28,267	1.03
7	9.015	28,153	1.04
8	9.373	28,801	1.06
9	9.298	28,357	1.00
10	9.298	28,206	1.04
Average	9.252	28,391	1.02
Standard Deviation	0.147	222	0.02
Relative Standard Deviation	1.589	0.783	2.14

SiliaChrom Selection Guide for Small Molecules

(Molecular Weight < 2,000 Dalton)



SiliaChrom HPLC columns Portfolio

How to build your Part Number

SiliaChrom HPLC columns are available in Narrow Bore, Analytical, Semi-Preparative, and Preparative sizes.

Below is an example of a SiliaChrom product number that shows you the way they are structured;

The product numbers start with the **phase** code, followed by the **particle size**, the **pore size**, the **internal diameter**, and finally the **length** codes.

Note: For Guard Columns, add the letter "G" between the "H" and the phase code.

Example;

SiliaChrom dt C18, 3 μm , 100 \AA , 4.6 mm x 150 mm = H141803E-N150

Particle Size		Pore Size		Internal Diameter			Column Length	
μm	Code	\AA	Code	Type of Columns	mm	Code	mm	Code
2.5	02	100	E	Narrow Bore	2.0	E	10	010
3.0	03	150	H	Narrow Bore	2.1	G	20	020
5.0	05	300	M	Narrow Bore	3.0	H	30	030
7.0	06			Analytical	4.6	N	50	050
10	07			Semi-Preparative	10	Q	100	100
20	09			Preparative	20	Y	150	150
				Preparative	30	V	200	200
				Preparative	50	W	250	250
				Preparative	100	X		

Labels for the diagram:

- Particle Size (points to 03)
- Pore Size (points to E)
- Internal Diameter (points to N)
- Column Length (points to 150)

*You may also find and buy your SiliaChrom online at www.silicycle.com/products/siliachrom-hplc-columns

phase code

SiliaChrom HPLC column Characteristics

SiliaChrom	Particle Size (µm)	Pore Size (Å)	Specific Surface Area (m ² /g)	Carbon Load (%)	pH Range	USP Code	T Limit* (°C)	Pressure Limit (psi)	Phase Code
SiliaChrom AQ C18	3, 5, 10	100	380	18	1.5 - 9.0	L1	60	5,000	H1518
SiliaChrom AQ C8	3, 5, 10	100	380	14	1.5 - 9.0	L7	60	5,000	H1508
SiliaChrom dt C18	2.5, 3, 5, 10	100	410 - 440	18	1.5 - 9.0	L1	60	5,000	H1418
SiliaChrom dt Si	3, 5, 10	100	410 - 440	n/a	1.0 - 8.0	L3	45	4,500	H1430
SiliaChrom XT C18	3, 5, 10	150	200	15	1.5 - 12.0	L1	60	5,000	H1718
SiliaChrom XT Fidelity C18	3, 5, 10	100	380	21	1.5 - 12.0	L1	60	5,000	HF1718
SiliaChrom SB C18	3, 5, 10	150	200	12	0.5 - 7.5	L1	60	4,500	H1018
SiliaChrom SB C18-300	5	300	80	5	0.5 - 7.5	L1	60	4,500	H1018
SiliaChrom SB C8	5	150	200	7	1.0 - 7.5	L7	60	4,500	H1008
SiliaChrom SB C8-300	5	300	80	3	1.0 - 7.5	L7	60	4,500	H1008
SiliaChrom XDB C18	5	150	200	15	1.5 - 9.0	L1	60	5,500	H1118
SiliaChrom XDB C8	5	150	200	8	1.5 - 9.0	L7	60	5,500	H1108
SiliaChrom XDB Si	5	150	200	n/a	1.0 - 8.0	L3	45	4,000	H1100
SiliaChrom XDB1 C18	3, 5	100	380 - 400	22	1.5 - 10.0	L1	60	5,500	H1218
SiliaChrom XDB1 C18-300	5, 10	300	80	8	1.5 - 9.0	L1	60	5,500	H1218
SiliaChrom XDB1 C8	5, 10	100	380 - 400	14	1.5 - 8.5	L7	60	5,500	H1208
SiliaChrom XDB1 C8-300	5	300	80	4	1.5 - 8.5	L7	60	5,500	H1208
SiliaChrom XDB1 C4	5	100	380 - 400	7	1.5 - 8.5	L26	60	5,500	H1204
SiliaChrom XDB1 C4-300	3, 5, 10	300	80	3	2.0 - 8.0	L26	60	5,500	H1204
SiliaChrom XDB1 C1	5	100	380 - 400	3	1.5 - 8.5	L13	60	5,500	H1201
SiliaChrom XDB1 C1-300	5	300	80	1	2.0 - 8.0	L13	60	5,500	H1201
SiliaChrom XDB1 CN	5, 10	100	380 - 400	5	2.0 - 8.5	L10	60	5,500	H1220
SiliaChrom XDB1 CN-300	5	300	80	3.5	2.0 - 8.0	L10	60	5,500	H1220
SiliaChrom XDB1 Amino	5, 10	100	380 - 400	7	2.0 - 8.5	L8	45	5,500	H1260
SiliaChrom XDB1 Amino-300	5	300	80	3.5	2.0 - 8.0	L8	45	5,500	H1260
SiliaChrom XDB1 Phenyl	5	100	380 - 400	12	1.5 - 9.0	L11	60	4,000	H1240
SiliaChrom XDB1 Phenyl-300	5	300	80	4.5	2.0 - 8.0	L11	60	4,000	H1240
SiliaChrom XDB1 Diol	5	100	380 - 400	5	2.0 - 8.0	L20	45	4,000	H1250
SiliaChrom XDB1 Diol-300	5	300	80	1	2.0 - 8.0	L20	45	4,000	H1250
SiliaChrom XDB1 Si	3, 5, 10	100	380 - 400	n/a	1.0 - 8.0	L3	45	4,000	H1230
SiliaChrom XDB1 Si-300	3, 5, 10	300	80	n/a	2.0 - 8.0	L3	45	4,000	H1230
SiliaChrom XDB2 C18	3, 5, 10	100	380	18	1.5 - 9.0	L1	60	5,000	H1318
SiliaChrom SCX	3, 5, 10	150	200	10	2.0 - 8.5	L9	45	5,000	H1800
SiliaChrom SCX-300	3, 5	300	80	3.5	2.0 - 8.0	L9	45	5,000	H1800
SiliaChrom SAX	3, 5, 10	100	380	6	2.0 - 8.5	L14	45	5,000	H1900
SiliaChrom SAX-300	3, 5	300	80	1	2.0 - 8.0	L14	45	5,000	H1900
SiliaChrom HILIC	3, 5, 10	100	380	8	2.0 - 8.0	-	60	5,000	H1600
SiliaChrom HILIC-300	5	300	80	2.5	2.0 - 8.0	-	60	5,000	H1600
SiliaChrom RPC	5, 7, 10, 20	n/a	750	polymer	1.0 - 14.0	L21	-	-	H920
SiliaChrom GF	5, 10	100	340	5	2.0 - 8.0	-	45	4,000	H900
SiliaChrom GF-300	5, 10	300	80	1	2.0 - 8.0	-	45	4,000	H900
SiliaChrom GF AMIDE	5, 10	100	340	5	2.0 - 8.0	-	60	4,000	H901
SiliaChrom GF AMIDE-300	5, 10	300	80	1	2.0 - 8.0	-	60	4,000	H901
SiliaChrom Cellulose T-DPC	5, 10	-	-	n/a	2.0 - 7.0	L40	40	< 700	H800
SiliaChrom Cellulose T-MB	5, 10	-	-	n/a	2.0 - 7.0	-	40	< 700	H820
SiliaChrom Amylose T-DPC	5, 10	-	-	n/a	2.0 - 7.0	L51	40	< 700	H810

*At pH range 5.0 - 7.5

SiliaChrom HPLC Selection Guide by Manufacturer

SiliaChrom HPLC Selection Guide by Manufacturer				
SiliaChrom HPLC Column	Agilent®	Eka Chemicals® (AksoNobel)	EMD Merck®	Phenomenex®
Reversed-Phases				
SiliaChrom dt C18 SiliaChrom AQ C18	Zorbax® SB Aq			Synergy™ Hydro RP Synergy™ Fusion RP
SiliaChrom XT C18	Zorbax® Extend C18	Kromasil® Eternity		Gemini® C18
SiliaChrom XT Fidelity C18	Zorbax® Extend C18	Kromasil® Eternity		Gemini®-NX C18
SiliaChrom SB C18	Zorbax® SB C18			
SiliaChrom SB C8	Zorbax® SB C8			
SiliaChrom XDB C18	Pursuit™ C18 Zorbax® XDB C18			
SiliaChrom XDB C8	Pursuit™ C8 Zorbax® XDB C8			
SiliaChrom XDB1 C18	Pursuit™ XRS C18	Kromasil® C18	LiChrospher® RP18e	Luna® C18
SiliaChrom XDB1 C8	Pursuit™ XRS C8	Kromasil® C8	LiChrospher® RP8	Luna® C8
SiliaChrom XDB1 Phenyl	Zorbax® SB Phenyl	Kromasil® Phenyl		
SiliaChrom XDB1 CN	Zorbax® SB CN			Luna® CN
SiliaChrom XDB2 C18	Zorbax® Eclipse Plus C18 Zorbax® Rx C18			Luna® C18(2)
Normal Phases				
SiliaChrom XDB1 Si	Zorbax® SIL Pursuit™ XRS Si	Kromasil® Si	LiChrospher® Si 100	Luna® Silica
SiliaChrom XDB1 Diol		Kromasil® Diol	LiChrospher® Diol	Luna® Diol
SiliaChrom XDB1 Amino	Zorbax® SB NH ₂	Kromasil® NH ₂		Luna® NH ₂
Ion Exchange Phases				
SiliaChrom SCX				Luna® SCX
SiliaChrom SAX	Agilent® SB-AX			Luna® SAX

SiliaChrom HPLC Chiral Selection Guide by Manufacturer				
SiliaChrom Chiral Column	Daicel	Eka Chemicals (AksoNobel)	Phenomenex	Supelco
SiliaChrom Cellulose T-DPC	ChiralCell® OD	Kromasil® CelluCoat	Lux® Cellulose-1	Astec™ Cellulose DMP
SiliaChrom Cellulose T-MB	ChiralCell® OJ		Lux® Cellulose-3	
SiliaChrom Amylose T-DPC	Chiralpak™ AD			



Supelco®	Thermo Fisher Scientific®	YMC®	Waters®	Others
	Acclaim® Polar Advantage Hypersil™ GOLD aQ C18	YMC™-PACK ODS -AQ	Atlantis® T3 Symmetry™ Shiels C18	Inertsil® ODS-3 ACE AQ C18
		YMC™ Triart C18	XTerra® C18	Nucleodur® C18 HTec
		YMC™ Triart C18	XBridge™ C18	
			X-Select™ CSH C18	
Discovery® C18 SUPELCOSIL™ LC-18-DB		YMC™-PACK ODS-A		Pinnacle™ DB C18
		YMC™-PACK C8		Pinnacle™ DB C8
Ascentis® C18		YMC™-PACK Pro RS	Sunfire™ C18 Symmetry™ C18	Ace® C18 HL Alltima™ HP C18 HiLoad
Ascentis® C8	Acclaim® C8	YMC™-PACK Pro C8	Sunfire™ C8 Symmetry™ C8	Ultra™ C8 ProntoSIL™ C8 SH
Ascentis® Phenyl	Hypersil® Phenyl	YMC™-PACK Ph		Ace® Phenyl ProntoSIL™ Phenyl
Ascentis® Cyano	Hypersil® Cyano	YMC™-PACK CN		ACE® CN
	Acclaim® 120 C18	YMC™-PACK Pro C18	SunFire™ C18	Ace® C18 Pinnacle™ II C18
Ascentis® Si		YMC™-PACK SIL	SunFire™ Si	Partasil™ Silica Nucleodur® SiOH
		YMC™-PACK Diol NP		ProntoSIL® Diol
	Acclaim® WAX	YMC™-PACK NH ₂	Spherisorb® Amino	Ultra® Amino ProntoSIL® Amino E
	Hypersil® SCX		Spherisorb® SCX	Nucleosil® SCX Partisil® SCX
	Hypersil® SAX		Spherisorb® Amino	Nucleosil® SAX Partisil® SAX

« Your analytical HPLC columns are the best!
I improved my separations with several of my methods and
got an award for my efforts. Thanks for the great products!!! »

Cliff Klimas from Bristol-Myers-Squibb, Pennington, NJ, USA

SiliaChrom HPLC Selection Guide by USP Code

SiliaChrom HPLC Selection Guide by USP Code			
USP Code	Packing Type	Description	SiliaChrom HPLC Columns
L1	Bonding: Octadecyl (C18) Support Type: Silica Particle size: 1.5 - 10 µm	Octadecyl silane chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.5 to 10 µm in diameter, or a monolithic rod	SiliaChrom dt C18 SiliaChrom AQ C18 SiliaChrom XT C18 SiliaChrom XT Fidelity C18 SiliaChrom SB C18 SiliaChrom XDB C18 SiliaChrom XDB1 C18 SiliaChrom XDB2 C18
L3	Bonding: Silica Support Type: Silica Particle size: 1.5 - 10 µm	Porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	SiliaChrom dt Si SiliaChrom XDB Si SiliaChrom XDB1 Si
L7	Bonding: Octyl (C8) Support Type: Silica Particle size: 1.5 - 10 µm	Octylsilane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	SiliaChrom AQ C8 SiliaChrom SB C8 SiliaChrom XDB C8 SiliaChrom XDB1 C8
L8	Bonding: Amine (NH ₂) Support Type: Silica Particle size: 1.5 - 10 µm	An essentially monomolecular layer of aminopropylsilane chemically bonded to totally porous silica gel support, 1.5 to 10 µm in diameter.	SiliaChrom XDB1 Amino
L9	Bonding: Strong cation exchange Support Type: Silica Particle size: 3 - 10 µm	Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter.	SiliaChrom SCX
L10	Bonding: Nitrile (CN) Support Type: Silica Particle size: 1.5 - 10 µm	Nitrile groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter.	SiliaChrom XDB1 CN
L11	Bonding: Phenyl Support Type: Silica Particle size: 1.5 - 10 µm	Phenyl groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter.	SiliaChrom XDB1 Phenyl
L13	Bonding: TMS (C1) Support Type: Silica Particle size: 3 - 10 µm	Trimethylsilane chemically bonded to porous silica particles, 3 to 10 µm in diameter.	SiliaChrom XDB1 C1
L14	Bonding: Strong anion exchange Support Type: Silica Particle size: 5 - 10 µm	Silica gel having a chemically bonded, strongly basic quaternary ammonium anion-exchange coating, 5 to 10 µm in diameter.	SiliaChrom SAX
L17	Bonding: Strong cation exchange Support Type: Polymer Particle size: 6 - 12 µm	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 6 to 12 µm in diameter.	SiliaChrom IEC SC-H
L21	Bonding: N/A Support Type: Polymer Particle size: 3 - 30 µm	A rigid, spherical styrene-divinylbenzene copolymer, 3 to 30 µm in diameter.	SiliaChrom RPC
L22	Bonding: Strong cation exchange Support Type: Polymer Particle size: ~10 µm	A cation-exchange resin made of porous polystyrene gel with sulfonic acid groups, about 10 µm in size.	SiliaChrom IEC SC-M
L26	Bonding: Butyl (C4) Support Type: Silica Particle size: 1.5 - 10 µm	Butyl silane chemically bonded to totally porous silica particles, 1.5 to 10 µm in diameter.	SiliaChrom XDB1 C4
L40	Bonding: Chiral Support Type: Silica Particle size: 5 - 20 µm	Cellulose tris-3,5-dimethylphenylcarbamate coated porous silica particles, 5 to 20 µm in diameter.	SiliaChrom Chiral Cellulose T-DPC
L42	Bonding: Mixed-mode C18/C8 Support Type: Silica Particle size: ~5 µm	Octylsilane and octadecylsilane groups chemically bonded to porous silica particles, 5 µm in diameter.	SiliaChrom C18/C8
L51	Bonding: Chiral Support Type: Silica Particle size: 5 - 10 µm	Amylose tris-3,5-dimethylphenylcarbamate-coated, porous, spherical, silica particles, 5 to 10 µm in diameter.	SiliaChrom Chiral Amylose T-DPC



How to Choose the Right SiliaChrom C18 Phase

C18 reversed-phase is the most used sorbent for HPLC applications. SiliCycle has developed over the years several C18 phases for specific analytes and/or matrices. The table below presents all SiliaChrom C18 phases available in the SiliCycle portfolio including a short description and characteristics. This table will help you choose the right SiliaChrom C18 phase based on your separation needs.

SiliaChrom C18 Reversed-Phase Characteristics						
SiliaChrom Phases	Description	%C	Pore Size (Å)	Surface Area (m ² /g)	pH Stability Range	Phase Description
SiliaChrom dt C18 SiliaChrom AQ C18	Universal 100% aqueous compatible C18 column. Most versatile column of the SiliCycle portfolio. Great retention for hydrophilic compounds. High sensitivity for LC-MS analysis. Same C18 functionalization but the SiliaChrom dt C18 is free of metal content.	18	100	410 - 440 380	1.5 - 9.0	Page 101
SiliaChrom SB C18	Column designed for extremely low pH conditions Compatibility with 100% aqueous mobile phase. Great sensitivity for LC-MS.	12	150	200	0.5 - 7.5	Page 107
SiliaChrom XT C18	High stability under high pH conditions Ideal for basic compounds.	15	150	380	1.5 - 12.0	Page 109
SiliaChrom XT Fidelity C18	Excellent stability under extreme pH and temperature conditions Ideal HPLC column for either metabolic or metabolite analysis.	21	150	380	1.5 - 12.0	Page 109
SiliaChrom XDB1 C18	Highest level of hydrophobicity of the SiliCycle C18 phases Designed for dirty samples. Oldest C18 phase technology.	22	100	380 - 400	1.5 - 10.0	Page 112
SiliaChrom XDB2 C18	Mid-level hydrophobicity and most popular phase for QC analysis Typical average value of carbon loading.	18	100	380 - 400	1.5 - 9.0	Page 114
SiliaChrom XDB C18	Lowest level of hydrophobicity of the SiliCycle C18 phase Ideal for separation of highly hydrophobic molecules such as fatty acids, barbiturates, fat-soluble vitamins & steroids.	15	150	200	1.5 - 9.0	Page 111
SiliaChrom XDB1 C18-300	Highest level of hydrophobicity for a C18 with wide pore size Designed for biochromatography applications (<i>peptides, proteines or nucleic acids</i>).	8	300	80	1.5 - 9.0	Page 112

« Needed a set of columns that work with a wide pH range.
SiliaChrom XT columns did the trick. »


Victor Nicolaev from Sanofi, Oro Valley, AZ, USA

SiliaChrom Reversed-Phase HPLC Column Character Evaluation

Our SiliaChrom HPLC columns are evaluated by USP and NIST tests for classification purpose and based on the selectivity chart. These tests allow the characterization and the comparison of various HPLC columns in order to determine the following parameters: void volume, retention capacity of hydrophobic compound, selectivity, efficiency and silanol activity. To run this test, we use a mixture of the five organic compounds listed below. Furthermore, we used the same test for side-by-side comparison on various SiliaChrom C18 columns against three well-known suppliers¹.

Reaction Mixture


- **Uracil (1)**: void volume marker (T_0)
- **Toluene (2)**: retention capacity of hydrophobic compounds (k'_{Tol})
- **Ethylbenzene (3)**: marker for the calculation of column efficiency for hydrophobic compounds ($k'_{Ethylbenzene}$)
- **Ratio Toluene/Ethylbenzene**: determination of selectivity ($\alpha_{Ethylbenzene/Toluene}$)
- **Amitriptyline (4)**: activity towards bases (*silanol activity evaluation*)
- **Quinizarin (5)**: activity towards chelating reagents (*metal contamination evaluation*)




Uracil (1)




Toluene (2)



Ethylbenzene (3)



Amitriptyline (4)



Quinizarin (5)

Description of the Column Aspects Evaluated

Hydrophobicity is measured by the retention factor of the hydrophobic analyte (*ethylbenzene*) using the following equation:

$$k' = \frac{(T_R - T_0)}{T_0} = \frac{\text{Ethylbenzene retention time} - \text{Uracil retention time (Void volume)}}{\text{Uracil retention time (Void volume)}}$$

Selectivity (α) is measured by the retention factor ratio between two similar compounds, ethylbenzene (k_2) and toluene (k_1):

$$\alpha = k_2/k_1$$

Column Efficiency is usually measured by the plate count (N) obtained for the ethylbenzene peak.

Chelating Tailing Factor – Metal Content is measured by the quinizarin peak symmetry. A symmetric peak shape indicates low activity toward chelating agent (*absence of metals*) and an asymmetric peak shape indicates the presence of metals by peak tailing (*high activity toward chelating reagents*).


Tailing Factor of Amitriptyline (Amitr.) – Silanol Activity is measured by the peak symmetry of amitriptyline (*basic compound*). Important silanol activity is often associated with peak tailing or an asymmetric peak. In other words, a highly deactivated column will have a lower peak asymmetry.

Chromatographic conditions:

- **HPLC System:** Thermo Surveyor with PDA
- **HPLC Software:** Xcalibur handling version 2.0
- **Column Size:** All HPLC columns: 4.6 x 150 mm, 5 μ m
- **Mobile Phase:** Methanol/buffer (80/20, v/v)
- **Buffer:** 20 mM of phosphate buffer adjusted at pH=7.0
- **Temperature:** 30°C
- **Flow rate:** 1.000 mL/min
- **Temperature:** UV scan (PDA), Total scan 200-600 nm

(¹Pharmacopeial Forum, Vol. 31(2) March-Apr. 2005, p.637)


SiliaChrom C18 HPLC Columns Versus the Competition




Uracil (1)




Toluene (2)




Ethylbenzene (3)



Amitriptyline (4)



Quinizarin (5)



C18 Column Character Evaluation Comparison

HPLC Columns	Hydrophobicity		Selectivity
	k' Toluene	k' Ethylbenz.	α Ethylbenz./Tol.
SiliaChrom XDB1 C18	2.14	3.09	1.44
SiliaChrom XDB2 C18	0.61	2.22	1.41
Phenomenex Luna C18	1.50	2.13	1.42
Agilent Zorbax SB C18	1.38	2.01	1.45
Waters Sunfire C18	1.72	2.45	1.43

C18 Column Character Evaluation Comparison

HPLC Columns	Efficiency		Metal Content	Silanol Activity
	N (/ meter) Ethylbenz.	R	A_s Quinizarin	TF Amitr.
SiliaChrom XDB1 C18	45,000	4.73	1.09	1.65
SiliaChrom XDB2 C18	28,000	3,30	1.10	1.18
Phenomenex Luna C18	22,000	2.90	1.23	1.20
Agilent Zorbax SB C18	25,000	3.20	1.10	1.55
Waters Sunfire C18	35,500	3.90	1.80	1.12

Results interpretation:


SiliaChrom columns compared advantageously over the competition; they present high column performances and with our wide portfolio, you can select the most suitable phase depending on the compound's nature. For example:

Basic analytes: SiliaChrom XDB2
Less polar analytes: SiliaChrom XDB1

The SiliaChrom columns performed very well compared to the competition. The Phenomenex Luna C18 has a lower efficiency as shown by its plate count (N). The Agilent Zorbax SB C18 column, shows peak tailing for amitriptyline which means that there are still some free OH presents on the surface (*activity towards bases*). Finally, the Waters Sunfire C18 seems to have high metal impurities as shown by the peak asymmetry obtained for quinizarin.

Full Range of Selectivity with SiliaChrom C18 HPLC Columns


Our most popular SiliaChrom reversed-phase C18 HPLC columns were evaluated by USP and NIST tests for classification purpose based on the selectivity chart. Select the most suitable SiliaChrom C18 based on your sample's properties.




Uracil (1)




Toluene (2)



Ethylbenzene (3)



Amitriptyline (4)



Quinizarin (5)

100% water compatible

Universal SiliaChrom C18

Lowest level of hydrophobicity


Ideal for low pH conditions

Highest level of hydrophobicity

Ideal for strong hydrophobic analytes

Ideal for dirty samples

Ideal for low interactions with basic analytes



SiliaChrom dt C18


SiliaChrom AQ C8 and C18

Description

Universal 100% aqueous compatible HPLC columns **SiliaChrom dt** adsorbent presents an optimum ratio of C18 short TMS chains and some free silanol groups. This new technology shows good peak shapes for any type of molecule (*acid, neutral and base*). The silica framework is exempt of any metal permitting a high sensitivity for LC-MS applications.

SiliaChrom AQ presents the same modified surface chemistry as dt but the silica framework contain low level of metal. C8 and C18 functions are available.

Structure



For C18 R = (CH₂)₁₇CH₃

SiliaChrom dt Purity: 99.9999% SiO₂ (no metal content)

SiliaChrom AQ Purity: 99.999% SiO₂

Sorbent Characteristics

- **Pore Size:** 100 Å
- **Specific Surface Area:** SiliaChrom dt C18 410 - 440 m²/g
SiliaChrom AQ C8 & C18 380 m²/g
- **Particle Sizes Available:** 2.5, 3, 5 and 10 μm
- **USP Code:** SiliaChrom dt C18 and AQ C18: L1
SiliaChrom AQ C8: L7
- **Typical Carbon Loading:** SiliaChrom dt & AQ C18: 18%
SiliaChrom AQ C8 14%

SiliaChrom dt and AQ Main Characteristics


- Enhanced retention of hydrophilic molecules
- Inertness for acidic and basic analytes
- Compatible from 100% aqueous mobile phase to 100% organic
- Exceptional stability from pH 1.5 to 9.0
- Good tolerance to direct injection of biological matrix (*dirty samples*)
- Reduces the need for mobile phase modifiers
- Low bleeding and high sensitivity for LC-MS
- Partially endcapped



Forensic

Separation of Catecholamines in Acidic Mobile Phase

Catecholamines are hydrophilic compounds with acidic functions. The mobile phase needs to be acidic to have the catecholamines under the molecular configuration and use the sorbent hydrophilic character to drive the separation.



Chromatographic conditions

- **Column:** SiliaChrom dt C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H141805E-N150
- **Mobile phase:** 1% Acetic Acid in water
- **Temperature:** 23°C
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 265 nm
- **Injection volume:** 5 μL


« Polar metabolites separation is very challenging.
Using SiliaChrom, dt C18 in normal phase solved the problem. »

Huns Nejad from BASF, Research Triangle Park, NC, USA




Assay for QC Testing of Blood Pressure and Cholesterol Medication


The SiliaChrom dt C18 presents a high lot-to-lot reproducibility, which makes it an excellent choice for quality control analysis in pharmaceutical laboratories.



A: Aspirine




B: Valsartan



C: Atorvastatin

Chromatographic conditions

- **Column:** SiliaChrom dt C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H141805E-N150
- **Mobile phase:** Methanol/H₂O (70/30), 0.1% (v/v) formic acid
- **Temperature:** 30°C
- **Flow rate:** 0.800 mL/min
- **Detector:** UV at 280 nm
- **Injection volume:** 10 μL




Peak Shape Evaluation for Zwitterion Fluoroquinolones

The SiliaChrom dt C18 presents a high separation capacity for zwitterion analysis.

Chromatographic conditions

- **Column:** SiliaChrom dt C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H141805E-N150
- **Mobile phase:** 2.5 mM potassium phosphate monobasic (adjust to pH 2.5 with H₃PO₄)/ethanol (68/32)
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 275 nm
- **Injection volume:** 10 μL



Peak Shape Results

Product	Asymmetry (USP) SiliaChrom dt C18	Asymmetry (USP) Atlantis dC18
Marbofloxacin	1.11	1.29
Sarafloxacin	1.08	1.14



Ropinirole and Amitriptyline Detection in Human Plasma

SiliaChrom dt C18 presents low bleeding and is excellent for dirty samples. Partial endcapping allows for some interactions with free silanol groups. The use of SiliaPrep CleanDRUG prior to injection onto the column will insure a very clean sample which results in very low ionic suppression when using in LC-MS/MS analysis. Another big advantage is the high selectivity of SiliaChrom dt C18 for all concentration levels.

Chromatographic conditions


- **Column:** SiliaChrom dt C18, 2.5 μm
- **Column size:** 3.0 x 30 mm
- **SiliCycle PN:** H141802E-H030
Sample preparation by SPE
SiliaPrep CleanDRUG 3 mL/200 mg
PN: SPEC-R651230B-03G
- **Mobile phase:**
MPA: 1 mM ammonium formate in (ACN/water, 10/90), 0.1% formic acid (v/v)
MPB: 1 mM ammonium formate in (ACN/water, 90/10), 0.1% formic acid (v/v)

Gradient		
Time (min)	MPA (%)	MPB (%)
0.00 - 0.20	85	15
0.21 - 1.20	50	50
1.21 - 1.60	0	100
1.61 - 3.50	85	15


- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **MS splitting flow:** 0.30 mL/min
- **Injection volume:** 5 μL

Tandem mass spectroscopy conditions




- **Detector:** Sciex API 3000, Applied Biosystem
- **Ion Source:** Positive Electrospray (ESI+)
- **Turbolon Ion Spray heater gas flow:** 8000 cc/min
- **Turbolon Ion Spray heater temperature:** 375°C
- **MRM Transition:** Ropinirole: m/z (261.2 \rightarrow 114.2)
Amitriptyline: m/z (278.4 \rightarrow 233.1)



Ropinirole




Amitriptyline




SiliaChrom AQ C18 is Highly Efficient for Basic Compounds

Amitriptyline, a strong basic compound, can be adsorbed on residual silanols on the surface of the packing material. With the traditional endcapping technique, this results in poor peak shapes. SiliCycle has developed a new method of silanol deactivation to eliminate the peak tailing from adsorption of compounds on residual silanol groups. This enables highly qualitative and quantitative analysis of strong basic compounds.




Chromatographic conditions

- **Column:** SiliaChrom AQ C18, 5 μ m
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H151805E-N150
- **Mobile phase:** 80/20 methanol/
20 mM potassium phosphate pH 7.00
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 254 nm
- **Injection volume:** 1 μ L




Evaluation of Resolution and Peak Shape

The SiliaChrom AQ C18 column is universal, efficient even for mixtures of basic and acidic compounds.



Chromatographic conditions

- **Column:**
SiliaChrom AQ C18, 5 μ m
Phenomenex Luna, C18 (2) 5 μ m
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H151805E-N150
- **Mobile phase:**
MPA: 5 mM potassium phosphate monobasic
(adjust to pH 2.5 with H_3PO_4)/ACN (90/10)
MPB: 5 mM potassium phosphate monobasic
(adjust to pH 2.5 with H_3PO_4)/ACN (10/90)
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 254 nm
- **Injection volume:** 5 μ L






Retention Capacity of DMSO on SiliaChrom AQ C18

DMSO (*Dimethylsulfoxide*) is an excellent solvent to solubilize most compounds. Unfortunately, this solvent is not volatile and with some C18 columns, the DMSO can interact with the stationary phase and decrease the selectivity. In this case, the only way to inhibit this effect is to use preparative chromatography. In this study, we show that DMSO does not interact with our SiliaChrom AQ C18. A linear gradient has been used from a highly aqueous mobile phase to a highly organic phase.

Chromatographic conditions

- **Column:** SiliaChrom AQ C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H151805E-N150
- **Mobile phase:** MPA 0.1% formic acid in water
MPB 0.1% formic acid in ACN
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 254 nm
- **Reconstitution solution (REC):** DMSO

Gradient		
Time (min)	% MPA	% MPB
0	90	10
9	10	90
10	10	90
11	90	10



Statistic Analysis Results

Conditions	A_{DMSO}	Tr_{DMSO} (min)	K'_{DMSO}	W_{DMSO}	$Tr_{\text{diphenhydramine}}$ (min)	$Tr_{\text{amitriptyline}}$ (min)
0% DMSO 5 μL	-	-	-	-	6.14	6.63
100% DMSO 5 μL	2.29	1.80	0.09	0.3	6.15	6.64
100% DMSO 20 μL	4.10	1.78	0.08	0.5	6.16	6.64

Conclusion: The study shows that DMSO does not interact with the SiliaChrom AQ C18. No specific retention is observed. The SiliaChrom AQ C18 is an excellent choice to purify components contaminated with DMSO.



Pharmaceutical



Forensic

Dewetting Phenomena








The dewetting phenomena is the formation of drops on the solid surface caused by hydrophobic repulsions of highly hydrophobic sorbents. This phenomena is illustrated by the following scheme.

General procedure

- The mixture of catecholamines is eluted on the column
- The flow is then stopped
- The column is stored in this condition during 18 h
- The mixture is then re-injected after a reconditioning step

Chromatographic conditions

- **Column:** SiliaChrom AQ C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H151805E-N150
- **Mobile phase:** 1% acetic acid in water
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 265 nm
- **Injection volume:** 5 μL



Conclusion: A small decrease in retention time is observed, but is not significant. The displacement has been resolved after the reconditioning step. The SiliaChrom AQ C18 does not present the dewetting phenomena.




SililaChrom SB C18 and C8

Description

SililaChrom SB C18 and **C8** surfaces are treated with an organic form of silicon to increase the number of silanol groups on the surface. After this step, the surface is bonded with a silane containing two functions. One function is a protecting group that shields the area and protects the surface from an acid attack from the mobile phase. The H_3O^+ ion does not have access to the surface to break the O-Si bond (*steric effect*). The other function is the linear hydrophobic chain with 18 or 8 carbons.

Structure



For C18 R = $(CH_2)_{17}CH_3$
For C8 R = $(CH_2)_7CH_3$

SililaChrom SB C18


SililaChrom SB C8

Sorbent Characteristics

- Pore Size: 150 Å
- Specific Surface Area: 200 - 220 m²/g
- Particle Sizes Available: 3, 5 and 10 µm
- USP Code: SililaChrom SB C18: L1
SililaChrom SB C8: L7
- Typical Carbon Loading: SililaChrom SB C18: 12%
SililaChrom SB C8: 7%

SililaChrom SB Main Characteristics

- Extremely low pH limits (0.5 - 7.5)
- Extremely low bleeding for LC-MS applications under acidic conditions
- Compatible with mobile phases ranging 100% water to 100% organic
- Non endcapped




Stability of SiliaChrom SB C18 at Low pH Conditions


Acidic mobile phases have widespread applications in the reversed phase HPLC separation of many important pharmaceutical and environmental compounds. Analytes such as pharmaceuticals and biomolecules often show peak shape, retention and selectivity changes when the mobile phase pH is changed from neutral to acidic pH (*pH 1.0*). In fact, lowering the pH helps to suppress silanol interactions between basic compounds and the residual surface silanols, thus resulting in less tailing and better retention of acidic compounds (*pKa lower than 2*).

Chromatographic conditions


- **Column:** SiliaChrom SB C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H101805H-N150
- **Mobile phase:** 2% TFA in ACN/water (60/40)
Solution pH: 1.00
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 270 nm
- **Injection volume:** 10 μL



SiliaChrom SB C18 (Ethylbenzene)			
Time (hour)	RT (min)	TF (USP)	N (USP)
0	5.91	1.01	14,014
24	5.89	1.02	14,085
48	5.77	1.02	14,023
72	5.83	1.02	14,076
96	5.85	1.01	14,087
120	5.84	1.02	14,050
Mean	5.85	1.02	14,056
RSD (%)	0.84	0.51	0.23



SiliaChrom SB C18 before



SiliaChrom SB C18 after

No column degradation under extreme pH conditions

The HPLC column was used under extreme pH conditions and, even after 5 days of continuous injections, the number of theoretical plates (*N*), the tailing factor (*TF*) and the retention time (*RT*) are comparable. The sorbent kept its chemical and structural integrity, which we have proven with similar chromatograms and scanning electron microscope pictures (*SEM*) before and after 120 hours of use.

In conclusion, our SiliaChrom SB C18 and SB C8 columns are stable at low pH conditions.



SiliaChrom XT C18 and XT C18 Fidelity


Description

SiliaChrom XT C18 and **XT C18 Fidelity** are suitable with low or high pH conditions. The key is to have a hybrid surface to reduce the solubility of silica at high pH. In fact, the SiliaChrom XT C18 and the XT C18 Fidelity silica are coated with a monomeric methyltriethoxysilane/tetraethoxysilane prepolymer, followed by a special thermic treatment to get a rigid surface that is less soluble than untreated silica itself at high pH.

The SiliaChrom XT C18 column is designed for applications at very high pH (*up to 12*) at room temperature but is also suitable for low pH (*down to 1.5*).

The SiliaChrom XT C18 Fidelity is used at high pH conditions and offers a higher thermal stability. The only difference between SiliaChrom XT C18 and the XT C18 Fidelity is the carbon loading. The SiliaChrom XT C18 Fidelity (21% C) presents a higher hydrophobic capacity than the SiliaChrom XT C18 (15% C).

Structure



SiliaChrom XT C18 and XT C18 Fidelity

Sorbent Characteristics

- Pore Size: SiliaChrom XT C18: 150 Å
SiliaChrom XT C18 Fidelity: 100 Å
- Specific Surface Area: 380 m²/g
- Particle Sizes Available: 3, 5 and 10 μm
- USP Code: L1
- Typical Carbon Loading: SiliaChrom XT C18: 15%
SiliaChrom XT C18 Fidelity: 21%

SiliaChrom XT Main Characteristics

- Excellent durability at high pH (*up to 12*)
- Ideal for basic compounds
- High thermal stability
- Ideal for auto-purification (*Prep. LC-MS*)
- Double endcapped
- Best HPLC columns for either metabolic or metabolite studies

« The high quality nature of the HPLC columns and plates from Silicycle has allowed us to achieve a level of reproducibility with our compound libraries that would be unheard of with any other production line. »

Steven Marois from Boston University CMLD, Boston, MA, USA




Stability of SiliaChrom XT C18 Fidelity at High pH Conditions


For some applications, it is necessary to work at high pH to increase the selectivity or to optimize peak shape. This is the case with basic organic compounds ($pK_a > 9.0$). It is the reason why it is important to have chromatographic phases stable at alkaline pH. This study demonstrates the stability of the SiliaChrom XT C18 Fidelity at high pH.

Chromatographic conditions


- **Column:** SiliaChrom XT C18 Fidelity, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** HF171805H-N150
- **Mobile phase:** 0.2% TEA in ACN/water (55/45) (v/v)
Solution pH: 11.5
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 270 nm



SiliaChrom XT C18 Fidelity (Ethylbenzene)			
Time (hour)	RT (min)	TF (USP)	N (USP)
0	13.35	1.01	13,623
24	13.29	1.01	13,648
48	13.27	1.01	13,689
72	13.25	1.00	13,604
96	13.24	1.00	13,649
120	13.28	1.00	13,582
Mean	13.28	1.01	13,633
RSD (%)	0.29	0.54	0.28



SiliaChrom XT C18 Fidelity before



SiliaChrom XT C18 Fidelity after

The HPLC column was used under extreme pH conditions, and even after 5 days of continuous injections, the number of theoretical plates (N), the tailing factor (TF) and the retention times (RT) remain constant. The sorbent kept its chemical and structural integrity, which we have proven with similar chromatograms and scanning electron microscope (SEM) pictures before and after 120 hours of use.

In conclusion, our SiliaChrom XT C18 Fidelity column is stable at high pH conditions.


SiliaChrom XDB C18 and C8

Description

SiliaChrom XDB C18 and **C8** are made of a special silica with a larger pore size and lower surface area for the separation of large hydrophobic molecules. The relatively low surface area allows a shorter retention time for such compounds.

SiliaChrom XDB phases are ideal for separation of barbiturates, fat-soluble vitamins, fatty acids and steroids.

Structure



For C18 R = (CH₂)₁₇CH₃
For C8 R = (CH₂)₇CH₃

SiliaChrom XDB C18

SiliaChrom XDB C8

Sorbent Characteristics

- Pore Size: 150 Å
- Specific Surface Area: 200 m²/g
- Particle Sizes Available: 3, 5 and 10 μm
- USP Code: SiliaChrom SB C18 L1
SiliaChrom SB C8 L7
- Typical Carbon Loading: SiliaChrom XDB C18: 15%
SiliaChrom XDB C8: 8%

SiliaChrom XDB C18 Main Characteristics

- Better choice for molecules > 500 Dalton
- High Loading capacity
- Wide pH range: 1.5 to 9.0
- Double endcapped




Environment

Resolution and Peak Shape of a Highly Hydrophobic Domestic Insecticide

This application illustrates the high separation efficiency of the SiliaChrom XDB C18 for very hydrophobic compounds.

Chromatographic conditions

- Column: SiliaChrom XDB C18, 5 μm
- Column size: 4.6 x 150 mm
- SiliCycle PN: H111805H-N150
- Mobile phase: ACN/water (90/10)
- Temperature: 23°C
- Flow rate: 1.000 mL/min
- Detector: UV at 235 nm
- Injection Volume: 1 μL



Column Performance Results

Compounds	Retention Time (min)	Peak Asymmetry Factor (USP)	Theoretical Plates (USP)
Uracil	1.72	1.26	5,936
Bioresmethrin	4.34	1.03	14,090

SiliaChrom XDB1 Family


Description

SiliaChrom XDB1 phases have a wider range of polarity than other SiliCycle HPLC phases (*C18 to normal phase*). These phases have the maximum bonding density regardless of the compound's polarity. This allows for the least amount of interaction between the analytes and the surface OH's. These phases are not recommended for samples containing highly hydrophobic compounds.

All **SiliaChrom XDB1** are available in 3, 5 and 10 μm except the Diol-300 which is not available in 3 μm .

The **SiliaChrom XDB1 C18** is designed for maximum hydrophobicity and efficiency for dirty samples.

Structure



SiliaChrom XDB1 C18

Sorbent Characteristics

See table next page.

SiliaChrom XDB1 Family Main Characteristics

- Better choice for molecules > 500 Dalton
- High loading capacity
- Double endcapped

Highly Base Deactivated C18




1. Uracil




2. Toluene



3. Ethylbenzene




4. Quinizarin



5. Amitriptyline

Chromatographic conditions

- **Column:** SiliaChrom XDB1 C18, 5 μm
- **Column size:** 4.6 x 150 mm
- **SiliCycle PN:** H121805E-N150
- **Mobile phase:** MeOH/20 mM potassium phosphate monobasic pH = 7.00 (80/20)
- **Temperature:** 23°C
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 254 nm
- **Injection Volume:** 1 μL



Column Performance Results

Compounds	Retention Time (min)	Peak Asymmetry Factor (USP)	Theoretical Plates (USP)
Uracil	1.49	1.27	3,778
Toluene	4.86	1.09	12,144
Ethylbenzene	6.40	1.02	13,026
Quinizarin	12.24	1.07	11,525
Amitriptyline	13.66	1.76	8,190

SiliaChrom XDB1 Sorbent Characteristics

SiliaChrom XDB1 Sorbent Characteristics						
SiliaChrom Phases	Description	USP Code	%C	Pore Size (Å)	Surface Area (m ² /g)	pH Stability Range
Reversed-Phases						
SiliaChrom XDB1 C18	Designed for maximum hydrophobicity and efficiency for dirty samples.	L1	22	100	380 - 400	1.5 - 10.0
SiliaChrom XDB1 C18-300		L1	8	300	80	1.5 - 9.0
SiliaChrom XDB1 C8	Exceptionally stable with high bonding coverage and low silanol activity.	L7	14	100	380 - 400	1.5 - 8.5
SiliaChrom XDB1 C8-300		L7	4	300	80	1.5 - 8.5
SiliaChrom XDB1 C4		L26	7	100	380 - 400	1.5 - 8.5
SiliaChrom XDB1 C4-300		L26	3	300	80	2.0 - 8.0
SiliaChrom XDB1 C1		L13	3	100	380 - 400	1.5 - 8.5
SiliaChrom XDB1 C1-300		L13	1	300	80	2.0 - 8.0
SiliaChrom XDB1 CN	Maximum hydrophobicity and works in normal and reversed-phase conditions.	L10	5	100	380 - 400	2.0 - 8.5
SiliaChrom XDB1 CN-300		L10	3.5	300	80	2.0 - 8.0
SiliaChrom XDB1 Phenyl	Highly retentive phase for aromatic and unsaturated compounds.	L11	12	100	380 - 400	1.5 - 9.0
SiliaChrom XDB1 Phenyl-300		L11	4.5	300	80	2.0 - 8.0
Normal Phases						
SiliaChrom XDB1 Si	Designed for normal phase conditions, presents a high surface area and a low metal content.	L3	n/a	100	380 - 400	1.0 - 8.0
SiliaChrom XDB1 Si-300		L3	n/a	300	80	2.0 - 8.0
SiliaChrom XDB1 Diol	Excellent for normal phase applications with the highest hydrophobic activity.	n/a	5	100	380 - 400	2.0 - 8.0
SiliaChrom XDB1 Diol-300		n/a	1	300	80	2.0 - 8.0
SiliaChrom XDB1 Amino	Superior general purpose amino phase. Recommended for normal phase analysis and excellent for sugar analysis.	L8	6	100	380 - 400	2.0 - 8.5
SiliaChrom XDB1 Amino-300		L8	2.5	300	80	2.0 - 8.0

« I have successfully used regular HPLC Analytical Columns for some analytical purpose, it works perfectly and accomodate good separation. »


Xiaohai Li from Scripps Reserch Institute, Jupiter, FL, USA

SiliaChrom XDB2 C18

Description

SiliaChrom XDB2 C18 is designed to be a mid-hydrophobic C18 phase with 18% of carbon loading, like most of the popular reversed-phase HPLC columns on the market. This phase demonstrates a balanced hydrophobic adsorption in order to avoid excessive retention of hydrophobic compounds.

Structure



SiliaChrom XDB2 C18


Sorbent Characteristics

- Pore Size: 100 Å
- Specific Surface Area: 380 - 400 m²/g
- Particle Sizes Available: 3, 5 and 10 μm
- USP Code: L1
- Typical Carbon Loading: 18%
- pH Stability: 1.5 - 9.0

SiliaChrom XDB2 C18 Main Characteristics

- Great column-to-column and batch-to-batch reproducibility (*popular for QC/QA laboratory*)
- Typical average value for carbon loading (**18%**)
- Good peak shape for basic, acidic and neutral analytes
- Stronger separation power for isomers
- Double endcapped

Highly Base Deactivated C18




1. Uracil




2. Toluene



3. Ethylbenzene




4. Quinizarin



5. Amitriptyline

Chromatographic conditions

- Column: SiliaChrom XDB2 C18, 5 μm
- Column size: 4.6 x 150 mm
- SiliCycle PN: H131805E-N150
- Mobile phase: MeOH/20 mM potassium phosphate monobasic pH = 7.00 (80/20)
- Temperature: 23°C
- Flow rate: 1.000 mL/min
- Detector: UV at 254 nm
- Injection Volume: 1 μL



Column Performance Results


Compounds	Retention Time (min)	Peak Asymmetry Factor (USP)	Theoretical Plates (USP)
Uracil	1.61	1.24	4 618
Toluene	4.73	1.04	12 858
Ethylbenzene	6.19	1.00	13 633
Quinizarin	11.18	1.03	12 277
Amitriptyline	13.53	1.29	9 451

SiliaChrom HILIC

Description

SiliaChrom HILIC (*hydrophilic interaction chromatography*) HPLC columns are designed to retain highly polar analytes. SiliaChrom HILIC has a selectivity that is complementary to reversed-phase columns. In fact, it has a higher retention for hydrophilic compounds in HILIC mode. HILIC sorbent is more stable and offers higher reproducibility than normal phase silica or amino columns. This phase is ideal for MedChem laboratories and is approved for SFC applications.

Structure




SiliaChrom HILIC




Food



SiliaChrom HILIC: Separation of Vitamin B Complex and Vitamin C



A. Thiamine (B1)




B. Pyridoxine (B6)



D. Riboflavin (B2)

Chromatographic conditions

- **Column:** SiliaChrom HILIC, 5 μ m
- **Column size:** 4.6 x 200 mm
- **SiliCycle PN:** H131805E-N150
- **Mobile phase:** 0.1% TFA in water/0.1% in ACN (90/10)
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 280 nm




SiliaChrom SCX-SAX

Description

SiliaChrom SCX provides excellent resolution and peak shape for cationic analytes. The benzene sulfonic acid function of the SiliaChrom SCX is providing the cationic phase and also the $\pi - \pi$ (*aromatic*) interaction. The SiliaChrom SCX is used for specific analysis of amino acids, anilines, drug salts, inorganic cations, and nucleosides.

Structure



Sorbent Characteristics

- **Pore Size:** SiliaChrom SCX: 150 Å
SiliaChrom SAX: 100 Å
- **Specific Surface Area:** SiliaChrom SCX: 200 m²/g
SiliaChrom SAX: 380 m²/g
- **Particle Sizes Available:** 3, 5 and 10 μm
- **USP Code:** SiliaChrom SCX L9
SiliaChrom SAX L14
- **Typical Carbon Loading:** SiliaChrom SCX: 10%
SiliaChrom SAX: 6%

Other SiliaChrom Products

Apart from the classic stationary phases, SiliCycle has also developed specific HPLC columns based on a silica matrix like our mixed-mode HPLC columns.

Mixed-Mode SiliaChrom HPLC Columns


Conjugate two surface function chemistries to optimize your separation in a single experiment. SiliCycle offers the following SiliaChrom Mixed-Mode HPLC columns:

- SiliaChrom C18/C8
- SiliaChrom C18/Amide
- SiliaChrom C18/Phenyl
- SiliaChrom C18/CN
- SiliaChrom C18/SCX
- SiliaChrom C18/SAX
- SiliaChrom C18/Nitrophenyl

Description

SiliaChrom SAX provides excellent resolution and peak shape for anionic analytes. SiliaChrom SAX presents propyltrimethyl ammonium chloride functions allowing ion exchange interactions to achieve effective ion chromatography. SiliaChrom SAX is used for specific analysis of pesticides, herbicides, inorganic anions and biological species such as nucleotides and glucosinolates

Structure



SiliaChrom SCX and SAX Main Characteristics

- **Narrow peak shape**
- **Rapid equilibration**
- **Compatible with organic modifiers**
- **Provides high efficiency and rapid separations**
- **Endcapped**



SiliaChrom HPLC Columns for Biochromatography

The rapid progress in the areas of genomics, proteomics, metabolomics and other biotechnology sectors has pushed scientists to develop innovative and efficient chromatographic methods. These methods have opened the way to better understanding of biomolecules and now offer impactful solutions effective at each level of the development of new commercial biopharmaceutical ingredients. Sorbent materials used in biochromatography and small molecule chromatography are similar but they require specific characteristics such as wide pore sizes and/or precise chemical resistance.

Separation and determination of peptides, proteins and nucleic acids can be done via different chromatography techniques. This section will highlight the SiliaChrom HPLC columns used in each following technique:

- Reversed-phase biochromatography for molecular weights (MW) < 5,000 Da
- Reversed-phase biochromatography for MW between 5,000 and 100,000 Da
- Ion exchange chromatography (*IEC*)
- Size exclusion chromatography (*SEC*)

SiliaChrom Reversed-Phases for Biochromatography ($MW < 5,000$ Da)


SiliaChrom Reversed-Phases for Biochromatography ($MW < 5,000$ Da)					
SiliaChrom Phases	Pore Size (Å)	%C	pH Stability Range	Characteristics	Phase Description
XT C18	150	21	1.5 - 12.0	Superior separation of basic & hydrophobic compounds Excellent peak shape in every condition. Excellent durability.	Page 109
XT C18 Fidelity	150	15	1.5 - 12.0		Page 109
dt C18	100	18	1.5 - 9.0	Superior separation of hydrophilic compounds. Mobile phase compatibility 100% aqueous to 100% organic. Inert and stable for acidic & basic analytes.	Page 101
RPC	n/a	Polymer	1.0 - 14.0	Guarantees chemical stability between pH 1.0 to 14.0 Basic compounds are well separated without peak tailing.	Page 117

Polymeric-Based SiliaChrom RPC

Description

SiliaChrom RPC phase is a hydrophobic copolymer based on polystyrene and divinylbenzene. The macroporous RPC reversed-phase resins are available in different particle sizes within a very narrow size distribution. The chemically inert polymer matrix of the SiliaChrom RPC guaranteed chemical stability and allows for use with applications in the pH range from 1 to 14. The capacity factor (K') values measured for aromatic and conjugated molecules on RPC columns are high due to the very pure uniform hydrophobic surface. The high efficiency and high selectivity of SiliaChrom RPC columns allow the separation of analytes in minutes. Even basic substances are separated efficiently without any peak tailing.

Structure



Phase Code: H920

SiliaChrom Reversed-Phases for Biochromatography (MW 5,000 - 100,000 Da)

SiliaChrom Reversed Phases for Biochromatography (MW 5,000 - 100,000)					
SiliaChrom Phases	Pore Size (Å)	%C	pH Stability Range	Characteristics	Phase Description
XDB1 C18-300	300	8	1.5 - 9.0	SiliaChrom C18 phase with wide pore size specially designed for peptide & protein separation	Page 112
XDB1 C8-300	300	4.5	1.5 - 8.0	SiliaChrom C8 phase with wide pore diameter presenting lower hydrophobicity than C18	Page 112
XDB1 C4-300	300	3	2.0 - 8.0	SiliaChrom C4 phase with wide pore diameter presenting lower hydrophobicity than C8 ideal for protein separation	Page 113
XDB1 CN-300	300	3.5	2.0 - 8.0	This Cyano phase provides the maximum hydrophobicity for normal phase analysis conditions	Page 113
XDB1 Phenyl-300	300	4.5	2.0 - 8.0	Reversed-phase permitting π - π interactions Excellent for aromatic and unsaturated compounds	Page 113

SiliaChrom GF Phases for Size Exclusion Chromatography



Size exclusion chromatography (*SEC*) also known as gel permeation chromatography (*GPC*) or gel filtration chromatography, separates molecules according to their size (*or, more accurately, according to their hydrodynamic diameter or hydrodynamic volume*). Smaller molecules are able to enter the pores of the media and, are therefore trapped and removed from the flow of the mobile phase. The average residence time in the pores depends upon the effective size of the analyte and the pore size itself. Larger molecules are excluded with essentially no retention. SiliaChrom GF column series are an appropriate set of phases to be used for size exclusion chromatography with silica-based material in normal phase conditions.

SiliaChrom GF Phases for Size Exclusion Chromatography					
SiliaChrom Phases	Functional Group	Pore Size (Å)	pH Stability Range	Separation of molecules with molecular weights between:	Phase Code
GF	Diol	100	2.0 - 8.0	5,000 and 100,000 Da	H900
GF-300	Diol	300	2.0 - 8.0	50,000 and 1,000,000 Da	H900
GF AMIDE	Amide	100	2.0 - 8.0	5,000 and 100,000 Da	H901
GF AMIDE-300	Amide	300	2.0 - 8.0	50,000 and 1,000,000 Da	H901

SiliaChrom IEC Phases for Ion Exchange Chromatography

SiliaChrom IEC series are composed of polystyrene polymer-based packing bearing different functionalities such as weak or strong cationic and anionic functions. SiliaChrom IEC phases are compatible with most mobile phases and samples with a pH range from 1 to 14. Polymer-based columns tend to have lower efficiencies for small molecules compared to silica-based columns due to their smaller surface area.

Nevertheless, SiliaChrom IEC packings are a good alternative for samples that require a mobile phase pH outside the normal operating range of standard silica-based columns. SiliaChrom IEC columns are generally used for ion exchange separations, and are also useful for non-aqueous gel permeation chromatography size exclusion analysis and ion exclusion analysis of organic acids and carbohydrates.

SiliaChrom IEC Phases for Ion Exchange Chromatography

SiliaChrom Phases	Functional Group	%C	pH Stability Range	Characteristics	Phase Code
SiliaChrom IEC SA	Dimethylammonium Chloride	8	2.0 - 8.0	Strong anion exchanger	H950
SiliaChrom IEC SC	Sulfonic Acid	4.5	2.0 - 8.0	Strong cation exchanger	H930
SiliaChrom IEC WA	Amino	3	2.0 - 8.0	Weak anion exchanger	H960
SiliaChrom IEC WC	Carboxylic Acid	4.5	2.0 - 8.0	Weak cation exchanger	H940

Each SiliaChrom IEC phases is available in particule size 5, 7, 10 and 20 μm





SiliaChrom Chiral Phases for Chiral Chromatography

Pharmaceutical Bio-Pharma

SiliaChrom Chiral Phases


SiliaChrom Chiral coated polysaccharide stationary phases are made with a spherical high quality silica support physically coated with a polymeric chiral selector such as amylose or cellulose derivatives. Due to the coated nature of these supports, solvents should be carefully selected for normal phase conditions.

Description

SiliaChrom Chiral Amylose T-DPC

Amylose tris-(3,5-dimethylphenylcarbamate) coated on a spherical silica support (*USP L51*). SiliaChrom Chiral Amylose T-DPC is used for chiral separation of alkaloids, tropines, amines, and beta blockers.

Structure




SiliaChrom Chiral Amylose T-DPC Phase Code: H810

Description

SiliaChrom Chiral Cellulose T-DPC

Cellulose tris-(3,5-dimethylphenylcarbamate) coated on a spherical silica support (*USP L40*). SiliaChrom Chiral Cellulose T-DPC is the most popular phase for chiral separation of alkaloids, tropines, amines, and beta blockers.

Structure




SiliaChrom Chiral Cellulose T-DPC Phase Code: H800

Description

SiliaChrom Chiral Cellulose T-MB

Cellulose tris-(4-methylbenzoate) coated on a spherical silica support. SiliaChrom Chiral Cellulose T-MB is used for chiral separation of aryl methyl esters and aryl methoxy esters.

Structure




SiliaChrom Chiral Cellulose T-MB Phase Code: H820

SiliaChrom Chiral Amylose T-DPC Enantiomeric Separation of L and D-val PMB

Chromatographic conditions

- **Column:** SiliaChrom Chiral Amylose T-DPC, 5 μm
- **Column size:** 4.6 x 250 mm
- **SiliCycle PN:** H81005T-N250
- **Mobile phase:** Hexane/Isopropanol (80/20)
- **Flow rate:** 1.000 mL/min
- **Detector:** UV at 254 nm



SiliaChrom Phases for Supercritical Fluid Chromatography (SFC)

Supercritical Fluid Chromatography (SFC) is a globally accepted powerful «green» chromatographic technique for separation of enantiomeric compounds and complex mixtures. For decades, it has been the preferred technique for preparative chromatography. The recent advances in preparative and analytical equipment for SFC coupled with the industry demand for reliable rapid analysis chromatography has created the need for a dependable source for SFC columns and necessary technical support. SFC is a chromatographic technique where the main component of the mobile phase is Carbon Dioxide (CO₂). A CO₂ based mobile phase composition is a «green» alternative to conventional HPLC mobile phases. The use of CO₂ based mobile phases enables the use of high performance preparative columns (10 to 50 mm ID) with a variety of particle sizes from 3 to 10 µm. Many SFC separations have successfully utilized stationary phases from normal phase HPLC such as unmodified silica, Diol, Amino and Cyano without the need for special packing techniques or hardware. The low viscosity of supercritical CO₂ allows separations to occur 3 to 5 times faster with 70 to 90% less in solvent usage than those for normal phase HPLC. Speed of the SFC separations, conservation of organic solvents and more concentrated product fractions make SFC a desirable preparative chromatographic technique for purifying chemical mixtures.



Pharmaceutical

SiliaChrom Phases for Supercritical Fluid Chromatography

SiliaChrom Phases	Pore Size (Å)	Carbon Loading %	Particle Size (µm)	Phase Description
SiliaChrom XDB1 Si	100	-	3, 5, 10	Page 113
SiliaChrom XDB1 Diol	100	5	3, 5, 10	Page 113
SiliaChrom XDB1 Amino	100	6	3, 5, 10	Page 113
SiliaChrom XDB1 CN	100	5	3, 5, 10	Page 113
SiliaChrom Hilic	100	8	3, 5, 10	Page 115

Hydrophobicity
-
↓
↓
+

SiliaChrom Guard Columns and Holders

SiliaChrom HPLC Guard Columns are designed to effectively protect both analytical and preparative HPLC columns. The usage of this shorter column is highly recommended to prolong column lifetime and does not alter the chromatography. SiliaChrom Guard Columns are cost effective and easy to use as a pre-filter to remove contaminants prior to injection. In liquid chromatography, contaminants introduced into the column can cause:

- Higher backpressure
- Resolution loss
- Baseline noise or drift
- Peak shape changes
- Irreversible damages (*column + system*)

SiliaChrom Guard Columns Packing and Dimensions

For optimal results and maximal protection, it is recommended to always use a guard column packed with the same packing material than the HPLC column. However, only the same chemistry is really needed. Particle size can be different but it is highly recommended to match the characteristics of the HPLC column used.

SiliaChrom Guard Columns are available in two different lengths (*10 and 20 mm*) and four internal diameters (*ID: 2.1, 4.6, 10 and 20 mm*). In most cases, a 10 mm length would be enough but if the sample contains important quantity of impurities, the 20 mm would then be more suitable.

The Guard Column internal diameter should be the same as the HPLC column or one size smaller. Never use a guard column with a larger ID than to the HPLC column (*risk of efficiency loss*).

SiliaChrom Guard Columns and HPLC Column Combinations					
		SiliaChrom Guard Cartridges Internal Diameter (mm)			
		2.1	4.6	10	20
SiliaChrom HPLC Column Internal Diameter (mm)	2.0	r			
	2.1	r			
	3.0	r			
	4.6		r		
	10		O	r	
	20			O	r
	30				r
	50				r

X = Preferred O = Possible

SiliaChrom Guard Holders




To use a SiliaChrom Guard Column you need to purchase the appropriate holder:

SiliaChrom Guard Holders					
Product Number		SiliaChrom Guard Cartridges Internal Diameter (mm)			
		2.1	4.6	10	20
Holders	HDW-000	r			
	HDW-001		r		
	HDW-002			r	
	HDW-003				r

Installation Procedure

1. If a new capillary tubing has been installed or if the LC system has not been operated for some time, flush the lines free of particulate before attaching the SiliaChrom Guard Column.
2. Insert the stainless column fitting (B) into the metal housing male connector (A) of the SiliaChrom Guard Holder.
3. Insert the SiliaChrom Guard Column (C) into the metal housing male connector (A) of the SiliaChrom Guard Holder. Make sure that the flat side of the stainless column fitting (B) is placed in front the SiliaChrom Guard Column frit (C).
4. Insert the stainless column fitting (D) into the metal housing female connector (E) of the SiliaChrom Guard Holder.
5. Finger tight both parts of the assembled SiliaChrom Guard Holder until leak free.
6. Connect the assembled SiliaChrom Guard Holder into the male fitting of the LC system tubing.
7. Once you have connected the SiliaChrom Guard Holder to the system and the LC column, connect the LC column to the detector and start pumping the working mobile phase at a low flow rate to equilibrate both the Guard Column and the LC column.
8. Gradually increase the flow rate to working condition and check for leaks. If leaking still occurs after persistent hand tightening, replace the leaking fitting.

- A Metal housing male connector
- B Stainless column fitting
- C SiliaChrom Guard Column
- D Stainless column fitting
- E Metal housing female connector



SiliaChrom Cleaning and Regeneration Procedures

If adequate care is taken, it is possible to maintain column efficiency and reliability over an extended period of time. This section is intended to give information on the different procedures to help extend HPLC column lifetime.

Difference between cleaning and regeneration

We usually make the assumption that, after a separation, all the material initially present in the column or cartridge has been eluted. After a run, the column is simply washed with 2-3 column volumes of the initial solvent mixture before starting a new separation. However, some impurities that are strongly retained on the column will accumulate at the inlet, if the mobile-phase composition is not strong enough to elute them during a regular run. Some non-negligible problems can arise when this happens: loss of performance, back-pressure build up, peak tailing, retention time shift or baseline drift.

To avoid this, it is highly recommended to perform regular cleaning of the column before any of these symptoms occurs. This process is simple and does not require modification of the usual chromatographic set up. When cleaning is not sufficient, a more thorough treatment, i.e. regeneration, may be necessary to avoid discarding the column.

Suggested Cleaning Procedure

The more you use a cleaning procedure, the less rigorous conditions be necessary. Cleaning should be performed after running a known “dirty” sample and prior to column storage using lower flow rate than usual (*typically from 20% to 50%*).

Column volume estimation is done using the following equation:

$$\text{Column Volume (packing's volume included) in mL} = \pi * [\text{Column Radius in cm}]^2 * [\text{Column Length in cm}]$$

SiliaChrom Suggested Cleaning Procedure	
SiliaChrom HPLC Column	Suggested Procedure
Reversed-Phase Columns (C18, C8, C4, Amine, Cyano, Phenyl, etc.)	<ul style="list-style-type: none"> - Water/ACN (95/5) to remove buffer - Water/ACN (5/95) - Mobile phase used during the separation
Normal Phase Columns (Amine, Cyano, Diol, etc.)	<ul style="list-style-type: none"> - MeOH/CHCl₃ (50/50) - Ethyl Acetate - Mobile phase used during the separation <p>Note: Never use water.</p>
Unbonded Silica Columns (Silica)	<ul style="list-style-type: none"> - Hexane - Isopropanol - Methylene Chloride - Mobile phase used during the separation
Ion Exchange Columns (SCX, SAX, etc.)	<ul style="list-style-type: none"> - 5 mM Phosphate Buffer pH 7.00 - Acetic Acid/Water (10/90) - Water - Methanol - Water



SiliaChrom Suggested Storage Conditions


When SiliaChrom HPLC Columns are not used for an extended period of time, do not allow high aqueous or high salt mobile phases to remain in the column. Remove aqueous buffers remaining in the column by washing with 20-30 column volumes of a 50% methanol or acetonitrile aqueous solution, followed by 20 column volumes of organic solvent such as methanol or acetonitrile.

Each column is shipped with two removable column end plugs to prevent the drying of the column bed. Always put these plugs back on tightly before column storage or when column is not being used.

SiliaChrom Suggested Storage Conditions		
SiliaChrom HPLC Columns		Recommended Storage Solvent
SiliaChrom AQ C18 SiliaChrom AQ C8 SiliaChrom dt C18 SiliaChrom SB C18 SiliaChrom SB C18-300 SiliaChrom SB C8 SiliaChrom SB C8-300 SiliaChrom XDB C18 SiliaChrom XDB C8 SiliaChrom XDB1 C18 SiliaChrom XDB1 C18-300 SiliaChrom XDB1 C8 SiliaChrom XDB1 C8-300 SiliaChrom XDB1 C4 SiliaChrom XDB1 C4-300 SiliaChrom XDB1 C1 SiliaChrom XDB1 C1-300 SiliaChrom XDB1 CN SiliaChrom XDB1 CN-300 SiliaChrom XDB1 Phenyl SiliaChrom XDB1 Phenyl-300 SiliaChrom XDB2 C18	SiliaChrom XT C18 SiliaChrom XT Fidelity C18 SiliaChrom C18/C8 SiliaChrom C18/Amide SiliaChrom C8/Amide SiliaChrom C18/Phenyl SiliaChrom C18/CN SiliaChrom C18/SCX SiliaChrom C18/SAX SiliaChrom C18/Nitrophenyl SiliaChrom Hilic SiliaChrom Hilic-300 SiliaChrom SCX SiliaChrom SCX-300 SiliaChrom SAX SiliaChrom SAX-300 SiliaChrom GF-300 SiliaChrom GF Amide SiliaChrom GF Amide-300	Methanol or Acetonitrile
SiliaChrom XDB1 Amino SiliaChrom XDB1 Amino-300		Butyl Chloride/Methanol
SiliaChrom Chiral Cellulose T-DPC SiliaChrom Chiral Cellulose T-MB	SiliaChrom Chiral Amylose T-DPC	Hexane/Isopropyl Alcohol (90/10)
SiliaChrom dt Si SiliaChrom XDB Si SiliaChrom XDB1 Si SiliaChrom XDB1 Si-300	SiliaChrom XDB1 Diol SiliaChrom XDB1 Diol-300	Isooctane/Ethanol

Important HPLC Definitions and Equations

Typical Chromatogram in liquid chromatography



Capacity Factor or Retention Factor (k') is measured by the retention factor of the analyte compared to an unretained peak (*void volume marker*) using the following equation:

$$k' = \frac{(T_R - T_0)}{T_0}$$

Where:

T_R : is the retention time of the analyte

T_0 : is the retention time of the unretained product

Efficiency (N) is usually measured by the plate count (N or also called *theoretical plate number*) using various equations. The most popular ones are:

By USP (*United States Pharmacopeia*)

$$N = 16 \times \left[\frac{t}{W} \right]^2$$

Where:

N : is the number of theoretical plates

t : is the retention time of the analyte

W : is the width at the base of the analyte

By DAB (*German Pharmacopeia*)

$$N = 5.54 \times \left[\frac{t}{W_{0.5}} \right]^2$$

Where:

N : is the number of theoretical plates

t : is the retention time of the analyte

$W_{0.5}$: is the width-at-half-height of the analyte

Selectivity (α) is measured by the retention factor ratio between two similar compounds.

$$\alpha = \frac{k'_2}{k'_1}$$

Where:

K'_1 : is the retention factor of product #1

K'_2 : is the retention factor of product #2

Separation's difficulty based on the selectivity value.
If the selectivity is:

- ≥ 2 : Easy separation
- 1.5 - 2: Possible separation*
- 1.2 - 1.5: Difficult separation
- ≤ 1.2 : Very difficult separation**

* Method adjustment could be required

** Selectivity's optimization may be required



Important HPLC Definitions and Equations (con't)

Resolution (R) can be expressed using the two following equations

$$R = \frac{\sqrt{N}}{4} \times \left(\frac{\alpha - 1}{\alpha}\right) \times \left(\frac{1 + k'_2}{k'_2}\right)$$

Where:

N: is the number of theoretical plates

α: is the selectivity

K₂': is the retention factor of product #2

$$R = \frac{2(t_2 - t_1)}{W_2 + W_1}$$

Where:

T₁: is the retention time of the product #1

T₂: is the retention time of the product #2

W₁: is the width at the base of the product #1

W₂: is the width at the base of the product #2

Summary of Influencing Factors in HPLC

To choose the most suitable HPLC column, various parameters need to be taken into account: the desired selectivity and the sample load as well as the efficiency and the resolution. All these parameters are influenced by different factors in HPLC summarized in the table below.

Liquid Chromatography Influencing Factors			
Properties	Typical Parameters	Affected Influencing Factors	Limitations
Chromatographic Conditions	Solvent	Retention, Efficiency	Back-pressure & phase stability
	pH	Selectivity, Resolution & Retention	Phase stability
	Flow Rate	Analysis Time, Efficiency & Resolution	Back-pressure & phase stability
Packing Characteristics	Chemistry (SiO ₂ , C18, etc.)	Selectivity, Resolution & Retention	Solvent used
	Pore Size (Å)	Sample Load & Selectivity	Size of the molecule
	Particle Size (µm)	Back-pressure, Efficiency & Resolution	Back-pressure & flow rate
HPLC Column Dimensions	Internal Diameter	Sample Load & Sensitivity	Back-pressure & flow rate
	Length	Analysis Time & Resolution	Back-pressure & analysis time too long

HPLC Method Scaling Up or Scaling Down Theory

When your experimental conditions are well optimized to get the most suitable purification, it is possible to scale up/down your method by keeping the same particle size and sorbent using these two equations:

Adjustment of the Sample Load

$$x_2 = \frac{x_1 \times r_2^2 \times C_L}{r_1^2} \quad \text{where} \quad \left[C_L = \frac{L_2}{L_1} \right]$$

Where:

x₁: is the maximum sample load in initial column

x₂: is the maximum sample load in final column

r₁: is the radius of the initial column

r₂: is the radius of the final column

L₁: is the length of the initial column

L₂: is the length of the final column

Adjustment of the Flow Rate

$$V_2 = \frac{V_1 \times r_2^2}{r_1^2}$$

Where:

V₁: is the flow rate use with the initial column

V₂: is the flow rate use with the final column

r₁: is the radius of the initial column

r₂: is the radius of the final column

How to Select the Right SiliaChrom HPLC Column

To select the right HPLC Column to use in your method development, read the section below to select the most appropriate SiliaChrom HPLC column to try first. However, before going forward in the selection, you need to have an idea of the sample quantity you need to purify as well as the liquid chromatography equipment available.

Remember: Resolution $R = \frac{\sqrt{N}}{4} \times \left(\frac{\alpha - 1}{\alpha}\right) \times \left(\frac{1 + k'_2}{k'_2}\right)$

Step 1. Find the Desired Selectivity by Selecting the Chemistry

When selecting an HPLC column, the most important factor is the selectivity in order to achieve an optimal resolution. A good knowledge of the composition of the sample mixture is crucial to select the most suitable chromatography mode to use in order to have good interactions between the sorbent and the compounds.

In liquid chromatography, there are various modes of operation possible based on the interaction mechanism of the solute with the stationary phase. Please refer you to previous sections to choose the most suitable phases to get optimal separation results.

Step 2. Select the Pore Diameter

To select the right pore diameter to use, find out the molecular weight of the solute. Typically, for small molecules, 100 - 150 Å pore size is recommended (*molecular weights below 5,000 Da*). For large molecules, such as peptides and proteins, 300 Å or higher is recommended.

Step 3. Find the Desired Efficiency & Resolution

Once you found the right selectivity, the second step is;

Be able to separate your sample with the shortest possible analysis time WITH optimal efficiency.

Two factors can influence the efficiency of a chromatography:

1. The particle size: influence on the resolution and back-pressure
2. The column dimensions (*internal diameter & length*): influence on the resolution and the sample load

Step 3.1. Select the Particle Size

For analytical applications, different particle sizes are available. The most common one being the 5 µm due to a good price/performance ratio. However, if you require a better separation and want to decrease analysis time, then 3 µm would be a better choice. Keep in mind that with a smaller particle size the backpressure will be higher.

For preparative applications, a larger particle size is usually used (*most frequently used is 10 µm*) with a larger column diameter (≥ 20 mm).



How to Select the Right SiliaChrom HPLC Column (con't)

Step 3.2 Select the Column Dimensions (*Influence on the Resolution*)

For analytical applications, the most often recommended format for initial trial is the 4.6 x 150 mm. Then, if you need more resolution, look at: decreasing the internal diameter or increasing the column length.

3.2.1 Select the Internal Diameter (*Influence on the Sample Load*)

With smaller internal diameters, you reduce solvent consumption due to lower flow rate required but increase analysis time. Furthermore, loading capacity is decreased as the diameter decreases. The table below identifies typical applications associated with typical internal diameters used in HPLC.

Select the Internal Diameter (ID)				
Type of columns	ID (mm)	Typical Sample Load	Typical Flow Rate	Typical Applications
Narrow Bore	2.1	0.04 - 1.5 mg	0.1 - 0.3 mL/min	Used with low sample volumes or when more sensitivity and selectivity are needed over 3 mm ID.
	3.0	0.08 - 3.0 mg	0.2 - 0.6 mL/min	Used to reduce flow rate and solvent consumption over 4.6 mm ID. It is gaining popularity.
Analytical	4.6	0.2 - 7.0 mg	0.5 - 1.5 mL/min	This is the most common ID used for traditional quantitative analysis.
Semi-Preparative	10	0.95 - 33.0 mg	2.5 - 7.0 mL/min	Used for small-scale (<i>mg</i>) preparative purifications.
Preparative	20	4.0 - 132.0 mg	9.0 - 28.0 mL/min	Used for large-scale (<i>hundreds of mg to gram</i>) purifications. The higher the diameter, the greater the loading capacity.
	30	8.5 - 297.0 mg	20.0 - 60.0 mL/min	
	50	24.0 - 800.0 mg	60.0 - 175.0 mL/min	
	100	96.0 - 3,200.0 mg	240.0 - 700.0 mL/min	

3.2.2 Select the Column Length (*Influence on the Resolution*)

The rule of thumb is that in presence of the same packing, longer columns provide better resolution and efficiency over shorter ones but with longer retention times and higher pressure. In general, it is preferable to try using the shortest column length possible. If the resolution is not good enough, increase the column length or use a smaller particle size with the same length. The table below presents the most suitable length/particle combinations.

Select the Column Length		
Length (mm)	Most Suitable Particle Size (µm)	Typical Applications
30 & 50	3 µm or smaller	Used to reduce flow rate and solvent consumption over 100 & 150 mm lengths.
100 & 150	3 or 5 µm	These are the most common lengths used for traditional quantitative analysis.
200 & 250	5 µm or larger	For difficult separations or for higher resolution.

Acceptable Modifications to an HPLC Validated Method

Even if you are using an FDA validated or a USP recommended method, some operating conditions can be adjusted if the modifications respect the acceptable specifications proposed by Pharmacopeias¹⁻³ and the FDA⁴. A side-by-side comparison of both the original and the adjusted method needs to be performed to demonstrate that the method's accuracy and precision is not affected by these modifications.

Acceptable Modifications to an HPLC Validated Method		
Parameters	Allowable modification	Examples of possible modifications
Mobile phase pH	± 0.2 units	Validated pH: 7.0 Allowed pH range: 6.8 - 7.2
Concentration of salts in buffer	± 10%	Validated concentration: 20 mM Allowed concentration range: 18 - 22 mM
Ratio of components in mobile phase	Only the minor components can be adjusted by ± 30% or ± 2% absolute (<i>i.e.</i> : in regards to the total mobile phase), whichever is the larger but should never exceed ± 10% absolute or removed totally.	Binary mixtures: Validated ratio: 50/50 Allowed ratio: 40/60 to 60/40 Validated ratio: 95/5 Allowed ratio: 93.5/6.5 to 96.5/3.5 Ternary mixtures: Validated ratio: 60/35/5 Allowed % of the 1 st component: 60% Allowed % of the 2 nd component: 25 - 45% Allowed % of the 3 rd component: 3.5 - 6.5% The total of the three components together need to be 100%.
Wavelength of UV detector	No modification allowed.	n/a
Column length	± 70%	Validated length: 150 mm Allowed length range: 45 - 255 mm
Column inner diameter	± 50%	Validated inner diameter: 4.6 mm Allowed inner diameter range: 2.3 - 10.6 mm
Flow rate	± 50%	Validated flow rate: 1.00 mL/min Allowed flow rate range: 0.5 - 1.5 mL/min
Injection volume	May be increased to as much as 2 times if no adverse effects on LOD and repeatability.	n/a
Particle size	No increase permitted. May be decreased by as much as 50%.	Validated particle size: 5 µm Allowed particle size range: 2.5 - 5 µm
Column temperature	± 20%	Validated temperature: 23°C Allowed length range: 18.4 - 27.6°C

¹ USP. USP 32-NF 27, Chromatography <621>. Rockville, MD: USP; 2009:227.

² USP. Second Supplement to USP 32-NF 27. Rockville, MD: USP; 2009:4147.

³ USP. USP 32-NF 27, Verification of Compendial Procedures <1226>. Rockville, MD: USP; 2009:736.

⁴ ORA Laboratory Procedure, Food and Drug Administration, modification criteria.

NEW

SiliaChrom HPLC Column Storage Cabinet

Protect your HPLC Column Inventory with the SiliaChrom Column Storage Cabinet

The SiliaChrom Column Storage Cabinet has been designed to safely store your HPLC column investment. Poor column storage can lead to reduced column performance and decreased product life.

The SiliaChrom Column Storage Cabinet is a bench top storage unit of solid steel construction with chrome-plated D-ring handles for added resistance. Up to 30 columns of 300 mm long (*or shorter*) can be stored in 5 separate drawers. Each drawer has a 6 position secure molded foam insert providing several storage possibilities. The foam insert can easily be customized to accommodate shorter columns, guard cartridges or HPLC tools and fittings. Each storage cabinet is stackable and supplied with rubber mounts. With the addition of a mounting bracket (*sold separately*), your SiliaChrom Column Storage Cabinet can be expanded to a multi-unit storage device.

The SiliaChrom Column Storage Cabinet is a cost-effective, expandable solution to conveniently index and store your HPLC column inventory.

Using the SiliaChrom HPLC Column Storage Cabinet ensures the following benefits

- Easy column identification.
- No more misplaced or lost columns.
- Increases column lifetime.
- Saves time and storage space.

SiliaChrom Column Storage Cabinet Characteristics

Column Storage Cabinet PN	AUT-0167
Mounting Bracket PN	AUT-0168
Dimensions (W x H x D)	279 x 325 x 408 mm
Drawer Height	51 mm
Column Formats	From 20 to 300 mm length
Column Storage Cabinet Weight	12 Kg





SiliCycle[®]

Consumables



SiliCycle® Syringe Filters



Using SiliCycle Syringe Filters guarantees the following benefits:


- Suitable for all laboratory filtrations.
- Increases column and apparatus life.
- Consistent and reproducible analysis.
- Limits system down time.



SiliCycle Syringe Filters Are Suitable for your Laboratory Filtration Needs

Syringe filters are offered by SiliCycle to help scientists with their laboratory filtration needs. They are compatible either with aqueous, organic or inorganic solutions. SiliCycle Syringe Filters are available in a wide variety of membranes (*Nylon, PTFE, PVDF, PES, and RC*) with polypropylene housing. They are the perfect choice for many applications in the fields of environmental, pharmaceutical, biotechnological, forensic and food & beverage laboratories.

SiliCycle Syringe Filters are specially designed to offer efficient and superior flow rate for any type of solution. These high quality products provide you with optimal filtration and particulate-free samples prior to injection. This extends apparatus lifetime which decreases overall analysis cost.



Nylon Syringe Filter 0.45 µm 25 mm
(PN: SF-NYL25-45)

A: Identification

Membrane type and pore size are clearly indentified on each syringe filter.

B: Luer-Lok™ female connector

Secures the connections to prevent «blow off».

C: Distribution ring

Generates even distribution of the flow rate during the elution.

D: Luer male connector

Can be easily connected to an automated system.

Field of Applications

Environmental

Environmental analysis of water, waste water, soil and sludge can be easily filtered using the SiliCycle Syringe Filter portfolio.

Pharmaceutical & Biotechnological

At each step of the drug discovery process, API or target compounds have to be isolated, purified and filtered prior to analysis. SiliCycle Syringe Filters achieve easy purification of complex physiological matrices.

Forensic

Elimination of particulates before injection to a sub-micron device is necessary prior to each clinical or toxicological analysis (*HPLC, GC or Mass Spectrometer*).

Food & beverage

Improving the detection limits for pesticides, herbicides, flavours and fragrances analysis requires high quality syringe filter products because of the presence of different types of particulates in the samples.



Find the Perfect Syringe Filter in 3 Simple Steps

Step 1: Selection of the Membrane by Application

Step 2: Selection of the Membrane Porosity Based on the Sample Nature

Step 3: Selection of the Membrane Diameter Based on Sample Volume

Step 1: Selection of the Membrane by Application

In the table below you can find the most popular applications where scientists are using SiliCycle Syringe Filters in their laboratories. The table below presents recommended membranes and alternative membranes for each application type. Finally, your choice of membrane can be validated by reviewing the description of each membrane in the SiliCycle portfolio.

Membrane Selection by Application		
Applications	Recommended Membranes	Alternative Membranes
Sample preparation prior to HPLC & GC analysis	Nylon	PTFE, PES or RC
Protein analysis or biological solutions	PVDF or PES	Nylon or RC
High particulate loads	Nylon	RC or PTFE
Pure organic and aggressive solvents	PTFE	RC or Nylon
Environmental analysis	PTFE or Nylon	RC
Clinical or toxicological analysis	PES	RC
Food & beverage analysis	PTFE or Nylon	RC
Sample preparation prior to ionic chromatography analysis	PES	PTFE or RC
Capillary electrophoresis analysis	PES	RC
Sample preparation prior to ICP-MS or AAS analysis (<i>trace metals</i>)	PES or Nylon	RC

Membrane Descriptions

Nylon:

Hydrophilic membrane is working well for general filtrations, aqueous & mixed organic solutions, medical assays and HPLC sample preparations.

Polytetrafluoroethylene (PTFE):

Hydrophobic polytetrafluoroethylene membrane is an excellent media for filtration of strong acids and aggressive organic solvents. Wetting the membrane with alcohol and water gives a hydrophilic characteristic.

Polyvinylidene Fluoride (PVDF):

Hydrophilic polyvinylidene fluoride membrane presents low extractables and provides good filtration of aqueous solutions, organic solvents and biological solutions.

Polyethersulfone (PES):

Hydrophilic polyethersulfone membrane is mainly used for life-science applications (*biological & pharmaceutical*) because of the ultra-low protein binding characteristic of the membrane.

Regenerated Cellulose (RC):

Regenerated cellulose is a universal hydrophilic membrane used in chromatography for filtration of aqueous samples and solvents. This membrane is also used for filtration of biomolecules because of the ultra-low binding capability of the media.

Step 2: Selection of the Membrane Porosity Based on the Sample Nature

Pore Size of 0.45 μm

This porosity is recommended for filtration of viscous solutions or solutions containing high levels of particulate matter. Generally used for aqueous or mixed organic solutions prior to injection on an HPLC column packed with $> 3 \mu\text{m}$ particles.

Pore Size of 0.20 μm

Generally used for aqueous or mixed organic solution prior to injection on an HPLC or UHPLC column packed with $\leq 3 \mu\text{m}$ particles. This porosity is recommended for particulate-sensitive methods.

Step 3: Selection of the Membrane Diameter Based on Sample Volume

4 mm Membrane

Use this diameter when the volume to filtrate is less than 1 mL.

13 mm Membrane

Use this diameter when the volume to filtrate is between 1 to 10 mL.

25 mm Membrane

Use this diameter when the volume to filtrate is between 10 to 100 mL.



SiliCycle Syringe Filters Chemical Compatibility Chart

SiliCycle Syringe Filters Chemical Compatibility Chart											
Chemical	Nylon	PTFE	PVDF	PES	RC	Chemical	Nylon	PTFE	PVDF	PES	RC
Acids						Halogenated Hydrocarbons					
Acetic, Glacial	Yellow	Green	Green	Red	Green	Carbon Tetrachloride	Green	Green	Green	Green	Green
Acetic, 25%	Green	Green	Green	Green	Green	Chloroform	Green	Green	Green	Red	Green
Formic, 25%	Red	Green	Green	-	Green	Methylene Chloride	Yellow	Green	Green	Red	Green
Hydrochloric, 25%	Red	Green	Green	Green	Red	Ketones					
Nitric, 25%	Red	Green	Green	Green	Red	Acetone	Green	Green	Green	Red	Green
Phosphoric, 25%	Red	Green	Green	-	Yellow	Cyclohexanone	Green	Green	Green	Red	Green
Sulfuric, 25%	Red	Green	Green	Green	Red	Methyl Ethyl Ketone	Green	Green	Yellow	Red	Green
Trichloroacetic, 10%	Red	Green	-	-	Green	Methyl Isobutyl Ketone	-	Green	Red	Red	Green
Bases						Esters					
Ammonium Hydroxide, 25%	Green	Green	Green	Green	Yellow	Amyl Acetate	Green	Green	Yellow	Red	Green
Sodium Hydroxide, 3 N	Green	Green	Green	Green	Yellow	Butyl Acetate	Green	Green	Yellow	Red	Green
Alcohols						Ethyl Acetate, Methyl Acetate					
Benzyl Alcohol	Green	Green	Green	Red	Green	Propyl Acetate	Green	Green	Yellow	Red	Green
Ethanol, 70%	Yellow	Green	Green	Green	Green	Organic Oxides					
Ethanol, 98%	Green	Green	Green	Green	Green	Dimethylsulfoxide	-	Green	Red	Red	Green
Ethylene Glycol	Green	Green	Green	Green	Green	Dioxane	Green	Green	Red	Red	Green
Glycerol	Green	Green	Green	Green	Green	Tetrahydrofuran	Green	Green	Red	Red	Green
Isopropanol	Green	Green	Green	Green	Green	Ethyl Ether	Green	Green	Green	Green	Green
Methanol, 98%	Green	Green	Green	Green	Green	Amides & Amines					
Hydrocarbons						Acetonitrile					
Hexane	Green	Green	Green	Yellow	Green	Aniline	-	Green	Red	-	Green
Xylene	Green	Green	Green	Green	Green	Diethylacetamide	Green	Green	Red	-	Green
Toluene, Benzene	Green	Green	Green	Red	Green	Dimethyl Formamide	Yellow	Red	Red	Red	Yellow
Legend						Pyridine					
Compatible						Green					
Not Compatible						Red					
Limited Compatibility (small volumes or short contact time)						Yellow					
						Miscellaneous					
						Hydrogen Peroxide, 30%					
						Phenol, Aqueous, 10%					



Typical Experimental Procedure

Loading

Fill the syringe with the liquid sample and allow a small volume of air to enter the Syringe. This small volume of air is used as a purge to minimize fluid retention when expelling the liquid sample from the syringe.

Assembly

Twist the SiliCycle Syringe Filter Luer-Lok™ onto the syringe. Make sure that the Luer-Lok™ is matching with the syringe to prevent the filter from coming off during the filtration.

Filtration

Direct the SiliCycle Syringe Filter tip into the collection vessel and apply gentle pressure onto the syringe plunger. Push the liquid sample and the air volume through the SiliCycle Syringe Filter to maximize sample recovery.

SiliCycle Syringe Filters Ordering Information

SiliCycle Syringe Filter Ordering Information					
Membrane Type	Quantity per box	13 mm Diameter Sample volume: 1 - 10 mL		25 mm Diameter Sample volume: 10 - 100 mL	
		0.20 µm	0.45 µm	0.20 µm	0.45 µm
Nylon	100 / box	SF-NYL13-20	SF-NYL13-45	SF-NYL25-20	SF-NYL25-45
	500 / box	SF-NYL13-20-L	SF-NYL13-45-L	SF-NYL25-20-L	SF-NYL25-45-L
	1,000 / box	SF-NYL13-20-M	SF-NYL13-45-M	SF-NYL25-20-M	SF-NYL25-45-M
Polytetrafluoroethylene (PTFE)	100 / box	SF-PTF13-20	SF-PTF13-45	SF-PTF25-20	SF-PTF25-45
	500 / box	SF-PTF13-20-L	SF-PTF13-45-L	SF-PTF25-20-L	SF-PTF25-45-L
	1,000 / box	SF-PTF13-20-M	SF-PTF13-45-M	SF-PTF25-20-M	SF-PTF25-45-M
Polyvinylidene Fluoride (PVDF)	100 / box	SF-PVD13-20	SF-PVD13-45	SF-PVD25-20	SF-PVD25-45
	500 / box	SF-PVD13-20-L	SF-PVD13-45-L	SF-PVD25-20-L	SF-PVD25-45-L
	1,000 / box	SF-PVD13-20-M	SF-PVD13-45-M	SF-PVD25-20-M	SF-PVD25-45-M
Polyethersulfone (PES)	100 / box	SF-PES13-20	SF-PES13-45	SF-PES25-20	SF-PES25-45
	500 / box	SF-PES13-20-L	SF-PES13-45-L	SF-PES25-20-L	SF-PES25-45-L
	1,000 / box	SF-PES13-20-M	SF-PES13-45-M	SF-PES25-20-M	SF-PES25-45-M
Regenerated Cellulose (RC)	100 / box	SF-RC13-20	SF-RC13-45	SF-RC25-20	SF-RC25-45
	500 / box	SF-RC13-20-L	SF-RC13-45-L	SF-RC25-20-L	SF-RC25-45-L
	1,000 / box	SF-RC13-20-M	SF-RC13-45-M	SF-RC25-20-M	SF-RC25-45-M

* Contact us for 4 mm membrane diameter, 1 mm pre-filter, and sterile syringe filters.

SiliCycle® Membrane Filters



Advantages of using SiliCycle Membrane Filters for your filtrations:

- Broad portfolio to cover all purification needs: 25 mm & 47 mm diameter available in 0.20 µm & 0.45 µm pore size.
- Many choices of membrane types to suit all fields of applications.



Fast and Cost Effective Separations with SiliCycle Membrane Filters

Membrane filters consist in a microporous films having a specific porosity which can be used to retain components (*particles and microorganisms*) with larger pore size compared to the membrane. Table below presents an overview of SiliCycle Membrane Filters specifications and typical applications.

SiliCycle Membrane Filters Overview

Characteristics	Nylon	Polytetra-fluoroethylene (PTFE)	Polyvinylidene Fluoride (PVDF)	Polyethersulfone (PES)	Regenerated Cellulose (RC)
Membrane Specifications					
• Hydrophilic / Hydrophobic	Hydrophilic	Hydrophobic*	Hydrophobic	Hydrophilic	Hydrophilic
• Extractable level	Extremely low	Low	Low	Low	Low
Chemical & Thermal Properties					
• Chemical resistance	Very good	Higher resistance	Good	Lower resistance	Good
• pH stability range	3 - 12	1 - 14	1 - 14	1 - 14	3 - 12
• Autoclavable	Yes	Yes	Yes	Yes	Yes
• Protein binding capacity	Medium	Low	Very low	Medium	Very low
Typical Applications: Sample Preparation prior to...					
• GC analysis	Preferred	Alternative	-	Alternative	Alternative
• HPLC analysis	Preferred	-	-	-	-
• ICP-MS or AAS analysis	Preferred	-	-	Preferred	Alternative
• Ion chromatography	-	-	-	Preferred	-
Typical Applications: Solvent Filtration & Sterilization of...					
• Agressive solutions	Alternative	Preferred	-	-	Alternative
• Aqueous solutions	Preferred	-	-	-	Alternative
• Organic solutions	Preferred	-	-	-	Alternative
Typical Applications: Filtration & Sterilization for...					
• Biomolecules analysis	Alternative	-	Preferred	Preferred	Alternative
• Clinical & forensic analysis	-	-	-	Preferred	Alternative
• Environmental analysis	Preferred	Preferred	-	-	Alternative
• Food & beverage analysis	Preferred	Preferred	-	-	Alternative
Venting Applications	-	Preferred	-	-	-

* Hydrophilic character in presence of alcohol and water.



SiliCycle Membrane Filters Ordering Information

SiliCycle Membrane Filters Ordering Information (100/box)				
Membrane Type	25 mm Membrane Diameter		47 mm Membrane Diameter	
	0.20 μm	0.45 μm	0.20 μm	0.45 μm
Nylon	MF-NYL25-20	MF-NYL25-45	MF-NYL47-20	MF-NYL47-45
Polytetrafluoroethylene (PTFE)	MF-PTF25-20	MF-PTF25-45	MF-PTF47-20	MF-PTF47-45
Polyvinylidene Fluoride (PVDF)	MF-PVD25-20	MF-PVD25-45	MF-PVD47-20	MF-PVD47-45
Polyethersulfone (PES)	MF-PES25-20	MF-PES25-45	MF-PES47-20	MF-PES47-45
Regenerated Cellulose (RC)	MF-RC25-20	MF-RC25-45	MF-RC47-20	MF-RC47-45

* Contact us if you are looking for a membrane not listed inside this table.



25 mm
Membrane Diameter



47 mm
Membrane Diameter

SiliCycle® Vials & Caps



Using SiliCycle Vials & Caps ensures the following benefits:

- Compatible with most autosamplers.
- Tight seal each and every time.
- Affordable price and excellent quality.
- Wide range of products available.



SiliCycle Vials & Caps Are Adapted to your Sample Requirements

Vials & Caps are new consumables that SiliCycle now offers to help customers with their day-to-day sample handling needs. SiliCycle Vials & Caps are compatible with most autosampler systems. They respect general characteristics offering the right product with optimal results. These characteristics are; sample volume, volatility, stability, chemical compatibility, and the type of autosampler (*robotic specifications*).

The SiliCycle Vials & Caps portfolio also offers products that are adapted to all types of samples and for specific storage conditions. Product characteristics such as the type of glass, type of septum and type of vial closure can be easily be selected to fit your needs.

Type of Glass

Clear Glass

SiliCycle uses USP Type 1 glass which is a borosilicate presenting the lowest leaching characteristics. This glass has an expansion coefficient of 33. It is a clear glass that does not offer protection from light.

Amber Glass

SiliCycle also uses USP Type 1 glass with an expansion coefficient of 51. The amber glass protects the sample from light.

Type of Closure

Snap-Top Vial

This type of vial closure is an extension of the crimp-top vial system but the plastic cap is stretched over the rim of the vial to seal the septum. This is the least reproducible sealed system. This type of vial is not recommended for volatile samples.

Crimp-Top Vial

This type of vial prevents evaporation of the sample because the seal is squeezed between the glass vial rim and the aluminum cap. A crimper is required to seal the cap on the vial.

Screw-Top Vial

This type of vial is the universal and most used vial. Users need simply to screw the cap to seal the septum between the cap and the glass rim. No tool is required.



SiliCycle Snap-Top Vials and Snap Caps

SiliCycle snap-top vials and caps are designed to rapidly seal sample prior to analysis. Because the snap sealing technology is less secure than the screw and the crimp technology, it is highly recommended to use these vials for short time storage and/or for low volatile samples. Snap-top vials are compatible with almost all autosamplers including autosamplers with robotic arms (*refer to the Autosampler Compatibility Table page 155*).

Snap-Top Vials Using 11 mm Snap Caps



2SP-C11-C



2SP-A11-C



2SP-CW11-C



2SP-AW11-C



1-5HSP-C11-C



0-3SP-C11-D



0-3SP-A11-D

SiliCycle Snap-Top Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Snap-Top Large Opening Clear Vial	12 x 32	2.0	2SP-C11-C	2SP-C11-M
Snap-Top Large Opening Amber Vial	12 x 32	2.0	2SP-A11-C	2SP-A11-M
Snap-Top Large Opening Clear Vial with White Patch	12 x 32	2.0	2SP-CW11-C	2SP-CW11-M
Snap-Top Large Opening Amber Vial with White Patch	12 x 32	2.0	2SP-AW11-C	2SP-AW11-M
Snap-Top High Recovery Clear Vial	12 x 32	1.5	1-5HSP-C11-C	1-5HSP-C11-M
Snap-Top Fused Insert Clear Vial (0.3 mL)	12 x 32	0.3	500/box, 0-3SP-C11-D	
Snap-Top Fused Insert Amber Vial (0.3 mL)	12 x 32	0.3	500/box, 0-3SP-A11-D	



CSP11-MI-C



BSP11-MI-C



NSP11-MI-C



RSP11-MI-C



YSP11-MI-C



CSP11-JI-C



CSP11-KSI-C



BSP11-KSI-M

SiliCycle Snap Cap Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
11 mm Clear Snap Cap	PP	Red PTFE/White Silicone	0.040/1.0	CSP11-MI-C	CSP11-MI-M
11 mm Blue Snap Cap	PP	Red PTFE/White Silicone	0.040/1.0	BSP11-MI-C	BSP11-MI-M
11 mm Black Snap Cap	PP	Red PTFE/White Silicone	0.040/1.0	NSP11-MI-C	NSP11-MI-M
11 mm Red Snap Cap	PP	Red PTFE/White Silicone	0.040/1.0	RSP11-MI-C	RSP11-MI-M
11 mm Yellow Snap Cap	PP	Red PTFE/White Silicone	0.040/1.0	YSP11-MI-C	YSP11-MI-M
11 mm Clear Snap Cap	PP	Red PTFE/White Silicone/Red PTFE	0.040/1.0	CSP11-JI-C	CSP11-JI-M
11 mm Clear Snap Cap, Slit	PP	Blue PTFE/White Silicone	0.040/1.0	CSP11-KSI-C	CSP11-KSI-M
11 mm Blue Snap Cap, Slit	PP	Blue PTFE/White Silicone	0.040/1.0	BSP11-KSI-C	BSP11-KSI-M

PP = Polypropylene

SiliCycle Crimp-Top Vials and Crimp Caps

SiliCycle crimp-top vials and caps are designed to be the safest alternative to seal sample prior to analysis or storage. When properly assembled, these vials are the best alternative for long-term storage and/or for high volatile solvents. Crimp-top vials are compatible with almost all autosamplers available on the market, including autosamplers with robotic arms (refer to the *Autosampler Compatibility Table* page 155).

Crimp-Top Vials Using 8 mm Crimp Caps



1-2CP-C8-D 1CPT-C8-D 0-8CP-C8-D 0-8CPZ-C8-D 0-7CP-A8-D 0-7CP-C8-D

SiliCycle Crimp-Top Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Quantity per box	Product Number
Crimp-Top Clear Vial	8 x 40	1.2	500	1-2CP-C8-D
Crimp-Top Tapered Bottom Clear Vial	8 x 40	1.0	500	1CPT-C8-D
Crimp-Top Clear Vial	8 x 30	0.8	500	0-8CP-C8-D
Crimp-Top Clear Vial	7 x 40	0.8	500	0-8CPZ-C8-D
Crimp-Top Tapered bottom Amber Vial	7 x 40	0.7	500	0-7CP-A8-D
Crimp-Top Tapered bottom Clear Vial	7 x 40	0.7	500	0-7CP-C8-D



ALCP8-MI-C



ALCP8-JI-C



ALCP8-CI-C

SiliCycle Crimp Cap Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
8 mm Crimp Cap	Al	Red PTFE/White Silicone	0.040/1.0	ALCP8-MI-C	ALCP8-MI-M
8 mm Crimp Cap	Al	Red PTFE/White Silicone/Red PTFE	0.040/1.0	ALCP8-JI-C	ALCP8-JI-M
8 mm Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	ALCP8-CI-C	ALCP8-CI-M

Al = Aluminum



Crimp-Top Vials Using 11 mm Crimp Caps



SiliCycle Crimp-Top Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Crimp-Top Large Opening Clear Vial	12 x 32	2.0	2CP-C11-C	2CP-C11-M
Crimp-Top Large Opening Amber Vial	12 x 32	2.0	2CP-A11-C	2CP-A11-M
Crimp-Top Large Opening Clear Vial with White Patch	12 x 32	2.0	2CP-CW11-C	2CP-CW11-M
Crimp-Top Large Opening Amber Vial with White Patch	12 x 32	2.0	2CP-AW11-C	2CP-AW11-M
Crimp-Top Clear Vial	15 x 46	4.0	500/box, 4CP-C11-D	
Crimp-Top Clear Vial	12 x 40	2.5	500/box, 2-5CP-C11-D	
Crimp-Top Round Bottom Clear Vial with White Patch	12 x 32	2.0	500/box, 2CPR-CW11-D	
Crimp-Top High Recovery Clear Vial	12 x 32	1.5	100/box, 1-5HCP-C11-C	
Crimp-Top Tapered Bottom Clear Vial	12 x 32	1.1	500/box, 1-1CPT-C11-D	
Crimp-Top Tapered Bottom Amber Vial	12 x 32	1.1	500/box, 1-1CPT-A11-D	
Crimp-Top Tapered Bottom Clear Vial	10 x 32	0.9	500/box, 0-9CPT-C11-D	
Crimp-Top Fused Insert Clear Vial	12 x 32	0.9	500/box, 0-9CP-C11-D	
Crimp-Top Fused Insert Clear Vial	12 x 32	0.3	500/box, 0-3CP-C11-D	
Crimp-Top Fused Insert Amber Vial	12 x 32	0.3	500/box, 0-3CP-A11-D	



SiliCycle Crimp Cap Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
11 mm Aluminum Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	ALCP11-CI-C	ALCP11-CI-M
11 mm Blue Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	BCP11-CI-C	BCP11-CI-M
11 mm Red Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	RCP11-CI-C	RCP11-CI-M
11 mm Green Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	GCP11-CI-C	GCP11-CI-M
11 mm Yellow Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	YCP11-CI-C	YCP11-CI-M
11 mm Aluminum Crimp Cap	Al	Red PTFE/White Silicone	0.040/1.0	ALCP11-MI-C	ALCP11-MI-M
11 mm Aluminum Magnetic Crimp Cap	Al	Red PTFE/White Silicone	0.040/1.0	ALCP11-MIM-C	ALCP11-MIM-M
11 mm Aluminum Magnetic Crimp Cap	Al	Clear FEP/Orange Silicone	0.040/1.0	ALCP11-CIM-C	ALCP11-CIM-M
11 mm Aluminum Magnetic Crimp Cap	Al	Red PTFE/White Silicone/Red PTFE	0.040/1.0	ALCP11-JIM-C	ALCP11-JIM-M

Al = Aluminum

TEL.: 1 418 874.0054 FAX: 1 418 874.0355 TOLL-FREE: 1 877.SILICYCLE (NORTH AMERICA ONLY) WWW.SILICYCLE.COM INFO@SILICYCLE.COM

Crimp-Top Vials (*Headspace*) Using 20 mm Crimp Caps



6CPTR-C20-C 10CPTR-C20-C 20CPTR-C20-C 6CPTF-C20-C 10CPTF-C20-C 20CPTF-C20-C 6CP-C20-C 10CP-C20-C 20CP-C20-C

SiliCycle Crimp-Top (*Headspace*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Crimp-Top Tapered Top/Radius Bottom Clear Vial	22 x 38	6.0	6CPTR-C20-C	6CPTR-C20-M
Crimp-Top Tapered Top/Radius Bottom Clear Vial	22 x 46	10.0	10CPTR-C20-C	10CPTR-C20-M
Crimp-Top Tapered Top/Radius Bottom Clear Vial	22 x 75	20.0	20CPTR-C20-C	20CPTR-C20-M
Crimp-Top Tapered Top/Flat Bottom Clear Vial	22 x 38	6.0	6CPTF-C20-C	6CPTF-C20-M
Crimp-Top Tapered Top/Flat Bottom Clear Vial	22 x 46	10.0	10CPTF-C20-C	10CPTF-C20-M
Crimp-Top Tapered Top/Flat Bottom Clear Vial	22 x 75	20.0	20CPTF-C20-C	20CPTF-C20-M
Crimp-Top Standard Top/Flat Bottom Clear Vial	22 x 38	6.0	6CP-C20-C	6CP-C20-M
Crimp-Top Standard Top/Flat Bottom Clear Vial	22 x 46	10.0	10CP-C20-C	10CP-C20-M
Crimp-Top Standard Top/Flat Bottom Clear Vial	22 x 75	20.0	20CP-C20-C	20CP-C20-M



GOCP20-5DIM-C



GOCP20-8DIM-C



BMCP20-DIM-C



ALCP20-BI-C



O800-B0125-C



O800-N0125-C



ALCP20-M-M

SiliCycle Crimp Cap and Septum Products

Description	Cap	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
20 mm Gold Magnetic Crimp Cap (5 mm)	Al	Natural PTFE/White Ultra Low Bleed Silicone	0.125/3.0	GOCP20-5DIM-C	GOCP20-5DIM-M
20 mm Gold Magnetic Crimp Cap (8 mm)	Al	Natural PTFE/White Ultra Low Bleed Silicone	0.125/3.0	GOCP20-8DIM-C	GOCP20-8DIM-M
20 mm Bi-Metallic Magnetic Crimp Cap	Al/Ti	Natural PTFE/White Ultra Low Bleed Silicone	0.125/3.0	BMCP20-DIM-C	BMCP20-DIM-M
20 mm Aluminum Crimp Cap	Al	Natural PTFE/White Ultra Low Bleed Silicone	0.125/3.0	ALCP20-BI-C	ALCP20-BI-M
20 mm Gold Magnetic Crimp Cap (5 mm)	Al	White PTFE/White Silicone	0.125/3.0	GOCP20-5MIM-C	GOCP20-5MIM-M
20 mm Gold Magnetic Crimp Cap (8 mm)	Al	White PTFE/White Silicone	0.125/3.0	GOCP20-8MIM-C	GOCP20-8MIM-M
20 mm Aluminum Crimp Cap	Al	White PTFE/White Ultra Low Bleed Silicone	0.125/3.0	ALCP20-NI-C	ALCP20-NI-M
20 mm Bi-Metallic Magnetic Crimp Cap	Al/Ti	White PTFE/White Ultra Low Bleed Silicone	0.125/3.0	BMCP20-NIM-C	BMCP20-NIM-M
20 mm Gold Magnetic Crimp Cap (5 mm)	Al	White PTFE/White Ultra Low Bleed Silicone	0.040/1.0	GOCP20-5NIM-C	GOCP20-5NIM-M
20 mm Gold Magnetic Crimp Cap (8 mm)	Al	White PTFE/White Ultra Low Bleed Silicone	0.040/1.0	GOCP20-8NIM-C	GOCP20-8NIM-M
0.800" Septum	-	Natural PTFE/White Ultra Low Bleed Silicone	0.040/1.0	O800-B0125-C	O800-B0125-M
0.800" Septum	-	White PTFE/White Ultra Low Bleed Silicone	0.040/1.0	O800-N0125-C	O800-N0125-M
20 mm Gold Magnetic Crimp Cap (5 mm)	Al	-	-	-	GOCP20-5M-M
20 mm Gold Magnetic Crimp Cap (8 mm)	Al	-	-	-	GOCP20-8M-M
20 mm Aluminum Crimp Cap	Al	-	-	-	ALCP20-M-M

Al = Aluminum, Al/Ti = Aluminum/Tin Plate



SiliCycle Screw-Top Vials and Screw Caps

SiliCycle screw-top vials and caps are developed to be easily reusable. No tools are required to securely seal the sample comparatively to crimp-top vials and caps. Screw-top technology provides low evaporation and rapid alternative for volatile sample storage and/or analysis. Screw-top vials are compatible with almost all autosamplers available on the market (refer to the *Autosampler Compatibility Table* page 155).

Screw-Top Vials (*Standard Opening*) Using 8 mm Screw Caps



2SW-C8-C



2SW-A8-C



2SW-CW8-C



2SW-AW8-C

SiliCycle Screw-Top (*Standard Opening*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Screw-Top Standard Opening Clear Vial	12 x 32	2.0	2SW-C8-C	2SW-C8-M
Screw-Top Standard Opening Amber Vial	12 x 32	2.0	2SW-A8-C	2SW-A8-M
Screw-Top Standard Opening Clear Vial with White Patch	12 x 32	2.0	2SW-CW8-C	2SW-CW8-M
Screw-Top Standard Opening Amber Vial with White Patch	12 x 32	2.0	2SW-AW8-C	2SW-AW8-M



NSW8-HI-C



NSW8-HSI-C



NSW8-KI-C



NSW8-KSI-C



NSW8-JI-C



NSW8-C

SiliCycle Screw Cap and Septum Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
8 mm Black Screw Cap	PP	Red PTFE/White Silicone	0.060/1.5	NSW8-HI-C	NSW8-HI-M
8 mm Black Screw Cap, Slit	PP	Red PTFE/White Silicone	0.060/1.5	NSW8-HSI-C	NSW8-HSI-M
8 mm Black Screw Cap	PP	Blue PTFE/White Silicone	0.060/1.5	NSW8-KI-C	NSW8-KI-M
8 mm Black Screw Cap, Slit	PP	Blue PTFE/White Silicone	0.060/1.5	NSW8-KSI-C	NSW8-KSI-M
8 mm Black Screw Cap	PP	Red PTFE/White Silicone/Red PTFE	0.060/1.5	NSW8-JI-C	NSW8-JI-M
0.321" Septum	-	Red PTFE/White Silicone	0.060/1.5	0321-H0060-C	0321-H0060-M
0.321" Slit Septum	-	Red PTFE/White Silicone	0.060/1.5	0321-HS0060-C	0321-HS0060-M
0.321" Septum	-	Blue PTFE/White Silicone	0.060/1.5	0321-K0060-C	0321-K0060-M
0.321" Slit Septum	-	Blue PTFE/White Silicone	0.060/1.5	0321-KS0060-C	0321-KS0060-M
0.321" Septum	-	Red PTFE/White Silicone/Red PTFE	0.040/1.0	0321-J0040-C	0321-J0040-M
0.321" Septum	-	Red PTFE/White Silicone/Red PTFE	0.060/1.5	0321-J0060-C	0321-J0060-M
8 mm Black Screw Cap	PP	-	-	NSW8-C	NSW8-M

PP= Polypropylene

Screw-Top Vials (*Large Opening*) Using 9 mm Screw Caps



SiliCycle Screw-Top (*Large Opening*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Screw-Top Large Opening Clear Vial	12 x 32	2.0	2SW-C9-C	2SW-C9-M
Screw-Top Large Opening Amber Vial	12 x 32	2.0	2SW-A9-C	2SW-A9-M
Screw-Top Large Opening Clear Vial with White Patch	12 x 32	2.0	2SW-CW9-C	2SW-CW9-M
Screw-Top Large Opening Amber Vial with White Patch	12 x 32	2.0	2SW-AW9-C	2SW-AW9-M
Screw-Top Fused Insert Polypropylene Vial	12 x 32	0.25	0-25SW-PP9-C	0-25SW-PP9-M
Screw-Top Fused Insert Polypropylene Vial	12 x 32	0.7	0-7SW-PP9-C	0-7SW-PP9-M
Screw-Top Large Opening Clear Vial	16 x 46	4.0	500/box, 4SW-C9-D	
Screw-Top High Recovery Clear Vial	12 x 32	1.5	500/box, 1-5HSW-C9-D	
Screw-Top High Recovery Amber Vial	12 x 32	1.5	500/box, 1-5HSW-A9-D	
Screw-Top Fused Insert Clear Vial	12 x 32	0.9	500/box, 0-9SW-C9-D	
Screw-Top Fused Insert Clear Vial	12 x 32	0.3	500/box, 0-3SW-C9-D	
Screw-Top Fused Insert Amber Vial	12 x 32	0.3	500/box, 0-3SW-A9-D	



SiliCycle Screw Cap and Septum Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
9 mm Blue Open-Top Screw Cap	PP	-	-	BSW9-C	BSW9-M
9 mm Black Open-Top Screw Cap	PP	-	-	NSW9-C	NSW9-M
9 mm Red Open-Top Screw Cap	PP	-	-	RSW9-C	RSW9-M
9 mm Green Open-Top Screw Cap	PP	-	-	GSW9-C	GSW9-M
9 mm Yellow Open-Top Screw Cap	PP	-	-	YSW9-C	YSW9-M
0.346" Septum	-	Red PTFE/White Silicone	0.040/1.0	0346-H0040-C	0346-H0040-M
0.346" Septum	-	White PTFE/Red Silicone	0.040/1.0	0346-P0040-C	0346-P0040-M
0.346" Septum, Slit	-	White PTFE/Red Silicone	0.040/1.0	0346-PS0040-C	0346-PS0040-M
0.346" Septum	-	Blue PTFE/White Silicone	0.040/1.0	0346-K0040-C	0346-K0040-M
0.346" Septum	-	Red PTFE/White Silicone/Red PTFE	0.040/1.0	0346-J0040-C	0346-J0040-M

PP = Polypropylene



Screw-Top Vials (Large Opening) Using 9 mm Screw Caps (con't)



BSW9-HI-C



BSW9-HSI-C



BSW9-HB-C



NSW9-HI-C



RSW9-HI-C



GSW9-HI-C



YSW9-HI-C



BSW9-OB-C



BSW9-BI-C



BSW9-PI-C



NSW9-PI-C



NSW9-PSI-C



BSW9-KI-C



BSW9-KSI-C



BSW9-KSB-C



NSW9-KI-C



NSW9-KSI-C



BSW9-JI-C



BSW9-EI-C



NSW9-EI-C



NSW9-ESI-C

SiliCycle Screw-Top (Large Opening) Cap Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
9 mm Blue Screw Cap	PP	Red PTFE/White Silicone	0.040/1.0	BSW9-HI-C	BSW9-HI-M
9 mm Blue Screw Cap, Slit	PP	Red PTFE/White Silicone	0.040/1.0	BSW9-HSI-C	BSW9-HSI-M
9 mm Blue Screw Cap, Bonded	PP	Red PTFE/White Silicone	0.040/1.0	BSW9-HB-C	BSW9-HB-M
9 mm Black Screw Cap	PP	Red PTFE/White Silicone	0.040/1.0	NSW9-HI-C	NSW9-HI-M
9 mm Red Screw Cap	PP	Red PTFE/White Silicone	0.040/1.0	RSW9-HI-C	RSW9-HI-M
9 mm Green Screw Cap	PP	Red PTFE/White Silicone	0.040/1.0	GSW9-HI-C	GSW9-HI-M
9 mm Yellow Screw Cap	PP	Red PTFE/White Silicone	0.040/1.0	YSW9-HI-C	YSW9-HI-M
9 mm Blue Screw Cap, Bonded	PP	White PTFE/Orange Silicone	0.040/1.0	BSW9-OB-C	BSW9-OB-M
9 mm Blue Screw Cap	PP	Clear PTFE/Orange Silicone	0.040/1.0	BSW9-BI-C	BSW9-BI-M
9 mm Blue Screw Cap	PP	White PTFE/Red Silicone	0.040/1.0	BSW9-PI-C	BSW9-PI-M
9 mm Blue Screw Cap, Slit	PP	White PTFE/Red Silicone	0.040/1.0	BSW9-PSI-C	BSW9-PSI-M
9 mm Back Screw Cap	PP	White PTFE/Red Silicone	0.040/1.0	NSW9-PI-C	NSW9-PI-M
9 mm Black Screw Cap, Slit	PP	White PTFE/Red Silicone	0.040/1.0	NSW9-PSI-C	NSW9-PSI-M
9 mm Blue Screw Cap	PP	Blue PTFE/White Silicone	0.040/1.0	BSW9-KI-C	BSW9-KI-M
9 mm Blue Screw Cap, Slit	PP	Blue PTFE/White Silicone	0.040/1.0	BSW9-KSI-C	BSW9-KSI-M
9 mm Blue Screw Cap, Slit, Bonded	PP	Blue PTFE/White Silicone	0.040/1.0	BSW9-KSB-C	BSW9-KSB-M
9 mm Black Screw Cap	PP	Blue PTFE/White Silicone	0.040/1.0	NSW9-KI-C	NSW9-KI-M
9 mm Black Screw Cap, Slit	PP	Blue PTFE/White Silicone	0.040/1.0	NSW9-KSI-C	NSW9-KSI-M
9 mm Blue Screw Cap	PP	Red PTFE/White Silicone/Red PTFE	0.040/1.0	BSW9-JI-C	BSW9-JI-M
9 mm Blue Screw Cap	PP	Natural PTFE/White Silicone	0.040/1.0	BSW9-EI-C	BSW9-EI-M
9 mm Black Screw Cap	PP	Natural PTFE/White Silicone	0.040/1.0	NSW9-EI-C	NSW9-EI-M
9 mm Black Screw Cap, Slit	PP	Natural PTFE/White Silicone	0.040/1.0	NSW9-ESI-C	NSW9-ESI-M

PP = Polypropylene

Screw-Top Vials (*Large Opening*) Using 10 mm Screw Caps



2SW-C10-C



2SW-A10-C



2SW-CW10-C



2SW-AW10-C

SiliCycle Screw-Top (*Large Opening*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Screw-Top Large Opening Clear Vial	12 x 32	2.0	2SW-C10-C	2SW-C10-M
Screw-Top Large Opening Amber Vial	12 x 32	2.0	2SW-A10-C	2SW-A10-M
Screw-Top Large Opening Clear Vial with White Patch	12 x 32	2.0	2SW-CW10-C	2SW-CW10-M
Screw-Top Large Opening Amber Vial with White Patch	12 x 32	2.0	2SW-AW10-C	2SW-AW10-M



NSW10-HI-C



NSW10-HSI-C



WSW10-HI-C



0374-H0060-C



0374-HS0060-C



NSW10-C



WSW10-C

SiliCycle Screw Cap and Septum Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
10 mm Black Screw Cap	PP	Red PTFE/White Silicone	0.060/1.5	NSW10-HI-C	NSW10-HI-M
10 mm Black Screw Cap, Slit	PP	Red PTFE/White Silicone	0.060/1.5	NSW10-HSI-C	NSW10-HSI-M
10 mm White Screw Cap	PP	Red PTFE/White Silicone	0.060/1.5	WSW10-HI-C	WSW10-HI-M
0.374" Septum	-	Red PTFE/White Silicone	0.060/1.5	0374-H0060-C	0374-H0060-M
0.374" Septum, Slit	-	Red PTFE/White Silicone	0.040/1.0	0374-HS0060-C	0374-HS0060-M
10 mm Black Screw Cap	PP	-	-	NSW10-C	NSW10-M
10 mm White Screw Cap	PP	-	-	WSW10-C	WSW10-M

PP= Polypropylene



Screw-Top Vials (*Large Opening*) Using 13 mm Screw Caps



4SW-C13-C



4SW-A13-C



4SW-CW13-C



4SW-AW13-C



3-5HSW-C13-E

SiliCycle Screw-Top (*Large Opening*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Screw-Top Large Opening Clear Vial	15 x 45	4.0	4SW-C13-C	4SW-C13-M
Screw-Top Large Opening Amber Vial	15 x 45	4.0	4SW-A13-C	4SW-A13-M
Screw-Top Large Opening Clear Vial with White Patch	15 x 45	4.0	4SW-CW13-C	4SW-CW13-M
Screw-Top Large Opening Amber Vial with White Patch	15 x 45	4.0	4SW-AW13-C	4SW-CW13-M
Screw-Top High Recovery Clear Vial	15 x 45	3.5	250/box, 3-5HSW-C13-E	



NSW13-HI-C



NSW13-HSI-C



NSW13-HB-C



NSW13-HSB-C



NSW13-ZHI-C



NSW13-ZHSI-C



NSW13-C



0479-H0060-C



0479-HS0060-C



0479-H0075-C



0479-HS0075-C

SiliCycle Screw Cap and Septum Products

Description	Cap Material	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
13 mm Black Screw Cap	PP	Red PTFE/White Silicone	0.060/1.5	NSW13-HI-C	NSW13-HI-M
13 mm Black Screw Cap, Slit	PP	Red PTFE/White Silicone	0.060/1.5	NSW13-HSI-C	NSW13-HSI-M
13 mm Black Screw Cap, Bonded	PP	Red PTFE/White Silicone	0.060/1.5	NSW13-HB-C	NSW13-HB-M
13 mm Black Screw Cap, Slit, Bonded	PP	Red PTFE/White Silicone	0.060/1.5	NSW13-HSB-C	NSW13-HSB-M
13 mm Black Screw Cap	PP	Red PTFE/White Silicone	0.075/1.9	NSW13-ZHI-C	NSW13-ZHI-M
13 mm Black Screw Cap, Slit	PP	Red PTFE/White Silicone	0.075/1.9	NSW13-ZHSI-C	NSW13-ZHSI-M
13 mm Black Screw Cap	PP	-	-	NSW13-C	WSW13-M
0.479" Septum	-	Red PTFE/White Silicone	0.060/1.5	0479-H0060-C	0479-H0060-M
0.479" Septum, Slit	-	Red PTFE/White Silicone	0.060/1.5	0479-HS0060-C	0479-HS0060-M
0.479" Septum	-	Red PTFE/White Silicone	0.075/1.9	0479-H0075-C	0479-H0075-M
0.479" Septum, Slit	-	Red PTFE/White Silicone	0.075/1.9	0479-HS0075-C	0479-HS0075-M

PP= Polypropylene

Screw-Top Vials (*Headspace*) Using 18 mm Screw Caps



10SW-C18-C



10SW-A18-C



20SW-C18-C



20SW-A18-C

SiliCycle Screw-Top (*Headspace*) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
Screw-Top Standard Opening Clear Vial	22 x 46	10.0	10SW-C18-C	10SW-C18-M
Screw-Top Standard Opening Amber Vial	22 x 46	10.0	10SW-A18-C	10SW-A18-M
Screw-Top Standard Opening Clear Vial	22 x 76	20.0	20SW-C18-C	20SW-C18-M
Screw-Top Standard Opening Amber Vial	22 x 76	20.0	20SW-A18-C	20SW-A18-M



ALSW18-LIM-C



ALSW18-IIM-C



ALSW18-ZQIM-C



ALSW18-QIM-C



0680-I0050-C



0680-L0060-C



0680-Q0125-C



0680-Q0050-C

SiliCycle Screw Cap (*Headspace*) and Septum Products

Description	Cap	Septum	Thickness (in/mm)	Product Number	
				100/box	1,000/box
18 mm Magnetic Aluminum Screw Cap	Al	Blue PTFE/White Ultra Low Bleed Silicone	0.060/1.5	ALSW18-LIM-C	ALSW18-LIM-M
18 mm Magnetic Aluminum Screw Cap	Al	Red PTFE/White Ultra Low Bleed Silicone	0.050/1.3	ALSW18-IIM-C	ALSW18-IIM-M
18 mm Magnetic Aluminum Screw Cap	Al	White PTFE/Translucent Blue Ultra Low Bleed Silicone	0.125/3.0	ALSW18-ZQIM-C	ALSW18-ZQIM-M
18 mm Magnetic Aluminum Screw Cap	Al	White PTFE/Translucent Blue Ultra Low Bleed Silicone	0.050/1.3	ALSW18-QIM-C	ALSW18-QIM-M
0.680" Septum	-	Blue PTFE/White Ultra Low Bleed Silicone	0.060/1.5	0680-L0060-C	0680-L0060-M
0.680" Septum	-	Red PTFE/White Ultra Low Bleed Silicone	0.050/1.3	0680-I0050-C	0680-I0050-M
0.680" Septum	-	White PTFE/Translucent Blue Ultra Low Bleed Silicone	0.125/3.0	0680-Q0125-C	0680-Q0125-M
0.680" Septum	-	White PTFE/Translucent Blue Ultra Low Bleed Silicone	0.050/1.3	0680-Q0050-C	0680-Q0050-M

Al = Aluminum



Insert Vials for Microsampling



0-3IF-C9-M



0-2IT-C8-D



0-2IT-C9-M



0-2IF-C8-M

SiliCycle Insert Vial Products

Description	Insert Vial	Vial Size OD x Height (mm)	Capacity (mL)	Quantity per box	Product Number
Flat Bottom Clear Vial Insert	9 mm Screw-Top 11 mm Snap-Top 11 mm Crimp-Top	6 x 31	0.3	1,000	0-3IF-C9-M
Tapered Bottom Clear Vial Insert	8 mm Screw-Top	5 x 30	0.2	500	0-2IT-C8-D
Flat Bottom Clear Vial Insert	9 mm Screw-Top 11 mm Snap-Top 11 mm Crimp-Top	6 x 30	0.2	1,000	0-2IT-C9-M
Tapered Bottom Clear Vial Insert	8 mm Screw-Top	5 x 31	0.2	1,000	0-2IF-C8-M

Shell (Neckless) Vials and Caps



1SLF-C8-D



2SLF-C12-D

SiliCycle Shell (Neckless) Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Quantity per box	Product Number
Shell Flat Bottom Clear Vial	8 x 40	1.0	500	1SLF-C8-D
Shell Flat Bottom Clear Vial	12 x 32	2.0	500	2SLF-C12-D



CSL8-M



CSL12-M

SiliCycle Cap Products for Shell Vials

Description	Cap Material	Product Number 1,000/box
8 mm Clear Cap for Shell Vial	PE	CSL8-M
12 mm Clear Cap for Shell Vial	PE	CSL12-M

PE= Polyethylene

EPA Type Screw-Top Vials and Screw Caps

SiliCycle EPA screw-top vials and caps are developed to meet USP and EPA standard for water samples and/or environmental testing. Neutral borosilicate type 1 glass clear and amber are used coupled to an ultra low or an extreme low bleed septa liner ensure tight sealed samples prior analysis.

EPA Screw-Top Vials Using 24 mm Screw Caps



SiliCycle EPA Type Screw-Top Vial Products

Description	Vial Size OD x Height (mm)	Capacity (mL)	Product Number	
			100/box	1,000/box
EPA Screw-Top Clear Vial	28 x 57	20.0	20SW-C24-C	20SW-C24-M
EPA Screw-Top Amber Vial	28 x 57	20.0	20SW-A24-C	20SW-A24-M
EPA Screw-Top Clear Vial	28 x 70	30.0	30SW-C24-C	30SW-C24-M
EPA Screw-Top Amber Vial	28 x 70	30.0	30SW-A24-C	30SW-A24-M
EPA Screw-Top Clear Vial	28 x 95	40.0	40SW-C24-C	40SW-C24-M
EPA Screw-Top Amber Vial	28 x 95	40.0	40SW-A24-C	40SW-A24-M

















SiliCycle EPA Type Screw Cap and Septum Products

Description	Cap	Septum	Thickness (in/mm)	Product Number	
				100/box	500/box
24 mm EPA Type Black Screw Cap	PP	Natural PTFE/Clear Extreme Low Bleed Silicone	0.125/3.0	NSW24-RI-C	NSW24-RI-D
24 mm EPA Type White Screw Cap	PP	Natural PTFE/Clear Extreme Low Bleed Silicone	0.125/3.0	WSW24-RI-C	WSW24-RI-D
24 mm EPA Type White Screw Cap	PP	Natural PTFE/Clear Ultra Low Bleed Silicone	0.125/3.0	WSW24-GI-C	WSW24-GI-D
24 mm EPA Type White Screw Cap, Bonded	PP	Natural PTFE/Clear Extreme Low Bleed Silicone	0.060/1.5	WSW24-RB-C	WSW24-RB-D
24 mm EPA Type Black Screw Cap, Bonded	PP	Natural PTFE/Clear Extreme Low Bleed Silicone	0.125/3.0	NSW24-RB-C	NSW24-RB-D
24 mm EPA Type White Screw Cap, Bonded	PP	Natural PTFE/Clear Extreme Low Bleed Silicone	0.125/3.0	WSW24-GB-C	WSW24-GB-D
0.880" Septum	-	Natural PTFE/Clear Extreme Low Bleed Silicone	0.060/1.5	0880-R0060-C	0880-R0060-D
0.880" Septum	-	Natural PTFE/Clear Extreme Low Bleed Silicone	0.125/3.0	0880-R0125-C	0880-R0125-D
0.880" Septum	-	Natural PTFE/Clear Ultra Low Bleed Silicone	0.125/3.0	0880-G0125-C	0880-G0125-D





PP= Polypropylene














SiliCycle Vials Comparison Chart (Actual Size)

					
Product No.	8-CPZ-C8-D	0-8CP-C8-D	0-7CP-A8-D	1-2CP-C8-D	1CPT-C8-D
Dimensions	7 x 40 mm	8 x 30 mm	7 x 40 mm	8 x 40 mm	8 x 40 mm
					
Part No.	2SW-C8-C	2SW-C9-C	4SW-C9-D	1-5HSW-C9-D	0-3CP-A11-D
Dimensions	12 x 32 mm	12 x 32 mm	15 x 45 mm	12 x 32 mm	12 x 32 mm
					
Part No.	0-3SW-C9-D	0-9SW-C9-D	2-5CP-C11-D	4CP-C11-D	
Dimensions	12 x 32 mm	12 x 32 mm	12 x 40 mm	15 x 46 mm	

SiliCycle Vials Comparison Chart (con't)

				
Part No.	0-2IT-C8-D	0-2IT-C9-M	0-2IF-C8-M	0-3IF-C9-M
Dimensions	5 x 30 mm	6 x 30 mm	5 x 31 mm	6 x 31 mm

					
Part No.	2CP-C11-C	2CPR-CW11-D	1-1CPT-C11-D	1-5HSP-C11-C	2SP-C11-C
Dimensions	12 x 32 mm	12 x 32 mm	12 x 32 mm	12 x 32 mm	12 x 32 mm

						
Part No.	0-9CPT-C11-D	0-9CP-C11-D	0-3SP-C11-D	4SW-C13-C	20CP-C20-C	10CP-C20-C
Dimensions	10 x 32 mm	12 x 32 mm	12 x 32 mm	15 x 45 mm	22 x 75 mm	22 x 46 mm



Instrument and Vial Compatibility Chart

The table below indicates the categories of vials that are compatible with various models of autosamplers.

*Certain autosamplers require the purchase of optional vial trays and, in few cases, programming upgrades to use all of the vials listed.

Autosampler Compatibility Table										
Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neckless)	Plate (page 49)
Agilent	1050, 1090		•		•					
	1050 (34 Pos. Tray), 1090 (34 Pos. Tray)	•								
	1100/1200		•		•					
	G1888A							•		
	7673A/7683A	•	•		•					
	7693A		•	•	•		•		•	
	HS7694						•	•		
	7697A							•		
	79855(A)		•		•					
	5880/5890		•		•					
	6850 (27 Pos. Tray)		•		•					
	6850 (22 Pos. Tray)							•		
	6890			•		•				
	CTC HTS+HTC PAL+CTC GC PAL	•	•		•					•
	CTC Combi PAL								•	
1100 Well-Plate/1100 Nanoflow		•	•	•					•	
1200 Well-plate/1200 SL plus		•	•	•					•	
AI	42 vial tray		•	•	•					
	60 vial tray	•	•	•	•					
	CTC A200S	•	•	•	•					
	Headspace							•		
AIM	CPS-100+CPS-200		•	•	•					
Alcott	708 AL, 728								•	
	738, 719 D/ D-PCS		•	•	•	•				
	719 AL		•	•	•	•	•	•		
Alpha M.O.S.	Prometheus/Fox/Kronos							•		
Antec Leyden	AS 100, 736 Unisampler, 738		•	•	•	•				
	Alexys		•	•	•					
ATAS GL	Focus		•	•	•			•		
Beckman	501, 502/502e, 507/507e	•	•	•	•	•				
	504	•								
	508 (System Gold)				•			•		
	Marathon, Promis		•	•	•					
	Triathlon, Standard Tray		•	•	•			•		
	Triathlon, LSV Tray	•						•		
	Triathlon, Super-LSV Tray							•		
Triathlon, Micro-Tray	•									
Bruker	LC51						•			
	Mapi1									•
Cambridge Scientific Instruments	205 Series, 300 Series		•	•	•	•	•			

- indicates that a cap having an outer flange is required for the vial to operate properly with the autosampler.
- indicates that the vials from this category are compatible with the autosampler in most configurations.
- indicates that a magnetic seal is required for use with the autosampler.

Autosampler Compatibility Table (con't)

Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neckless)	Plate (page 49)
Carlo Erba	AS100, A200LC, AS300	•	•	•	•				•	
	AS200, AS200S	•	•	•	•					
	AS800, 42 vial tray		•	•	•					
	AS800, 60 vial tray	•	•	•	•					
	HS250, 500, 800, 850							•		
Cecil Instruments	CE4800		•	•	•	•				
	AutoQuest		•	•	•	•				
CTC	A200S	•	•	•	•	•				
	A200 LC	•	•	•	•	•		•		
	HS 500							•		
CTC (LEAP)	LC PAL (216 Pos.)		•	•	•	•		•		
	HTX PAL, HTC PAL, HTS PAL (200 Pos. Tray), Combi PAL (200 Pos. Tray), GC PAL (200 Pos. Tray)	•								
	HTX PAL, HTC PAL, HTS PAL (54/98 Pos. Tray)	•	•	•	•	•		•		
	HTX PAL, HTC PAL, HTS PAL (32 Pos. Tray), Combi PAL (32 Pos. Tray), GC PAL (32 Pos. Tray), Combi PAL SPME Mode (32 Pos. Tray)							•		
	Combi PAL (98 Pos. Tray), GC PAL (98 Pos. Tray)	•	•		•					
	Combi PAL SPME Mode (98 Pos. Tray)		•		•					
DANI	ALS 39.80, ALS 86.80, ALS 1000		•		•					
	HS39.50, HS86.50							•		
	Master AS		•		•			•		
	Master DHS							•		
Dionex	Gina 50	•	•		•		•			
	AS 50	•	•	•	•	•				
	Summit ASI 100, Micro-Tray (192 Pos.)	•								
	Summit ASI 100, Analytical-Tray (117 Pos.)		•	•	•					
	Summit ASI 100, Semiprep.-Tray (63 Pos.)						•			
	Famos (LC Packings/Dionex)		•	•	•	•		•		
	UltiMate Analytical, cylindrical, WPS-3000 SL, 120 Pos. Rack (2ml)		•	•	•	•		•		
	UltiMate Analytical, conical, WPS-3000 SL, 120 (3x40) Pos. Rack (1.1ml=2ml w. Inserts)		•					•		
	UltiMate Micro, conical, WPS-3000 SL, 120 (3x40) Pos. Rack (250µl), UltiMate Nano/Cap/Micro, WPS-3000 SL, 216 (3x72) Pos. Rack (1.2ml)	•							•	
	UltiMate Semipreparative, WPS-3000 SL, 66 (3x22) Pos. Rack (4ml)						•	•		
	AS 40						•		•	
AS-HV			•							
D-Star	DAS 10		•	•						
Dynatech	42 vial tray		•	•	•					
	60 vial tray	•	•	•	•					
	LC2000	•								
	GC111, GC311	•	•	•						
	LC-241	•	•	•						
Eksigent	NanoLC-AS1		•	•						
ESA	540-MT/540		•	•	•					•
EST	LC-241plus		•	•						
EST Analytical	Cobra L/S GC Autosampler; 120 vial tray		•	•	•	•				
	Cobra L/S GC Autosampler; 60 vial tray, Markelov HS9000							•		
Finnigan	A200S	•	•	•	•					

Autosampler Compatibility Table (con't)

Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neck/less)	Plate (page 49)
Fisons	AS100, A200LC, AS300	•	•	•	•				•	
	AS200	•	•	•	•					
	AS200S	•	•	•	•					
	AS800, 42 vial tray		•	•	•					
	AS800, 60 vial tray	•	•	•	•					
	HS250, HS500, HS800, HS 850							•		
GBC	Avanta Ultra Z		•		•				•	
	LC 1650		•	•						
GE Healthcare	Ettan A-905		•		•	•				
GE Instruments	Sievers 900							•		
Gerstel	MPS	•	•	•	•		•	•		•
Gilson	201/202, 221/222, 231/401/232/402, Aspec, Aspec Xli, Aspec XL4			•	•				•	
	221XL/222XL, 223, 231XL/232XL/233XL	•								
	Nano Injektor			•	•					
	235/235P/SP 235/SP 235P	•		•	•					
Gynkotek	Gina 50	•	•		•		•			
HTA	HT200H							•		
	HT250D, HT280T, HT300L		•	•	•	•		•		
	HT300A, HT310A		•	•	•	•				
ICI	LC1600	•	•							
IMT GmbH	PTA3000							•		
Jasco	AS 2055/AS 2055 (i), AS 2057/AS 2057 (i), AS 2059	•	•	•	•	•				
	851/AS-950/AS-1550/AS-1555			•						
	AS-2059/AS-2059Plus			•						•
	AS-2059-SF/X-LC	•		•						•
Knauer	K-3800 (Basic Marathon), Smartline K-3950, PLATINblue AS-1		•	•	•			•		
Konik -Tech	Robokrom Static HS							•		
	Robokrom HRGC	•	•							
	Robokrom HPLC		•	•	•	•				
Kontron	MSI 660			•			•			
	360, 460	•	•	•	•					
	360/460/560/565	•	•	•	•					
LDC	713-60	•								
	Marathon, Promis		•	•	•					
Metrohm	Triathlon		•	•						
PerkinElmer	Series 200, 25 vial tray, ISS-225, 25 vial tray							•		
	Series 200, 85 vial tray, ISS-100, 85 vial tray, ISS-200, 85 vial tray, ISS-225, 85 vial tray		•			•		•		
	Series 200, 81/100 vial tray, Integral 4000, ISS-100, 100 vial tray, ISS-200, 100 vial tray		•			•				
	Series 200, 205 vial tray	•	•			•				
	Series 200, 225 vial tray	•								
	AI-1	•	•							
	AS-100/AS-100B	•	•							
	AS2000/AS2000B	•	•			•				
	AS-300, AS8300, Autosystem	•	•							
	HS 6, HS40/HS100/101							•		
	TurboMatrix HS16/HS40/HS40 XL/ HS40 Trap/HS110/ HS110 Trap							•		
	ISS-200, 145 vial tray	•								
	ISS-225, 205 vial tray	•	•			•				
	ISS-225, 100 vial tray + 80 vial tray		•			•				
	LC 600, 42 vial tray	•								
	LC 600, 60 vial tray		•			•				
Clarus 400, 500, 600		•								

Autosampler Compatibility Table (con't)

Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neckless)	Plate (page 49)
Pharmacia	LKB 2157-010		•	•	•					
	LKB 2157-020	•	•							
	Akta A-900		•	•						
Polymer Laboratories	PL-AS RT		•	•	•	•	•			
	GPC 110/210		•	•						
Quma Elektronik	QHSS-40							•		
Sedere	-		•		•					
Selerity	3100		•	•						
Sepiatech	Sepmatix									•
SGE	LS-3200	•								
Shimadzu	AOC-5000	•	•		•			•		
	AOC-14/1400, AOC-17, AOC-20/20i/20s 150 Pos. Tray		•	•	•	•	•			
	AOC-20/20i/20s 96 Pos. Tray						•			
	LC-20A		•	•	•	•	•			
	SIL-2AS, SIL-6A, SIL-10A/SIL-10AF/SIL-10AP/SIL-10Ai/SIL-10AxL/Rack S 100 Pos.	•	•	•	•	•	•		•	
	SIL-6B/SIL-7A/SIL-8A/SIL-9A		•	•	•	•	•		•	
	SIL-10A/SIL-10AF/SIL-10AP/SIL-10Ai/SIL-10AxL/Rack L 80 Pos.						•		•	
	SIL-10A/SIL-10AF/SIL-10AP/SIL-10Ai/SIL-10AxL/Rack MTP2 192 Pos., SIL-10HTA/SIL-10HTC 350 pos. Tray								•	
	SIL-10HTA/SIL-10HTC 140 Pos. Tray		•	•	•	•			•	
	SIL-10HTA/SIL-10HTC 100 Pos. Tray						•		•	
	SIL-10ADvp		•	•	•	•	•		•	
	HTA 200 H							•		
	SIL-20A (Prominence) 105 vial tray/SIL-20AC (Prominence) 70 vial tray	•	•	•	•	•				
	SIL-20A/Sil-20AC (Prominence) 175 vial tray								•	
	SIL-20A/Sil-20AC (Prominence) 50 vial tray, LC2010C + LC2010A 100 Pos. Tray						•		•	
	LC2010C + LC2010A 350 Pos. Tray								•	
	LC2010C + LC2010A 140 Pos. Tray		•	•	•	•			•	
HSS-2B								•		
Spark	Marathon Basic, Standard 96 Pos. Tray, Midas, Large Capacity 96 Pos. Tray, Promis, SPH 125		•	•	•					
	Marathon Basic Prep King Size 48 Pos. Tray, Midas, Large Volume 24 Pos. Tray							•		
	Midas, Standard 84 Pos. Tray, Alias		•	•	•			•		
	Triathlon, Standard 96 Tray		•	•	•				•	
	Triathlon, LSV 72 Pos. Tray						•			
	Triathlon, Super-LSV 32 Pos. Tray							•		
	Triathlon, Micro 160 Pos. Tray	•								
	Endurance 48 Pos. Tray, Reliance 48 Pos. Tray		•	•	•					
	Integrity		•	•	•					•
	Prospekt 2		•	•						
Reliance/Symbiosis Pharma		•	•						•	
Symbiosis Pico									•	
Spectra-Physics	8875, 8880		•	•	•					
	SpectraSYSTEM AS1000, AS3000, AS3500	•	•	•	•				•	
Sykam	S 5200		•		•					
Talbot	ASI		•		•					
Teledyne Tekmar	7000/7000HT/7050							•		
	HT3A							•		

- indicates that a cap having an outer flange is required for the vial to operate properly with the autosampler.
- indicates that the vials from this category are compatible with the autosampler in most configurations.
- indicates that a magnetic seal is required for use with the autosampler.

Autosampler Compatibility Table (con't)

Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neckless)	Plate (page 49)
Thermo Scientific	AS1000 (Trace GC), AS200, AS2000 90 vial tray (Trace GC)	•	•	•	•					
	AS300	•	•	•	•				•	
	AS2000 30 vial tray							•		
	AI3000 (II)/AS3000 (II) AS3500 (Trace GC + Focus GC)	•	•		•			•		
	A200LC, AS 100	•	•	•	•				•	
	SpectraSYSTEM AS 1000, AS 3000, AS 3500	•	•	•	•				•	
	A200S	•	•	•	•					
	AS800, 42 vial tray		•	•	•					
	AS800, 60 vial tray	•	•	•	•					
	HS250, HS500, HS800, HS 850, HS2000								•	
	TriPlus (=GC PAL) (AS+ Duo)	•	•	•	•				•	
	TriPlus HS, TriPlus SPME								•	
	Surveyor (Surveyor Plus)	•	•	•	•				•	•
	Accela High Speed LC Autosampler (200 Pos.)	•	•	•	•					
Accela Open Autosampler (342 Pos)	•	•	•	•					•	
Tosoh	AS 8010		•		•					
	TSK-6080				•		•			
Tracor	770/771/772		•	•	•					
Unicam	4247, 4710		•	•	•					
	4700 (GC)	•								
	4700 (LC)	•		•	•					
	LC-XP		•	•	•		•			
	S4/S8	•								
Varian	ProStar 400, Standard 96 Pos. Tray, ProStar 410, Large Capacity 96 Pos. Tray		•	•	•	•				
	ProStar 400, King Size 48 Pos. Tray, ProStar 410, Large Volume 24 Pos. Tray							•		
	ProStar 410, Standard 84 Pos. Tray		•	•	•	•		•		
	ProStar 420, Standard 96 Pos. Tray		•	•	•	•			•	
	ProStar 420, LSV 72 Pos. Tray	•					•			
	ProStar 420, Super-LSV 32 Pos. Tray							•		
	ProStar 420, Micro 160 Pos. Tray	•								
	ProStar 430, 48 Pos. Tray		•	•	•					
	8035			•	•					
	8000, 8100		•	•	•					
	8200		•	•	•	•				
	8400 (100 Pos.), 8410-Autoinjector (10 x 2ml; 6 x 5ml; 5 x 10ml)		•	•	•			•		
	CP-910, 911, 912		•	•	•					
	CP-940, 941		•							
	LC 9100/LC 9095/LC 9090		•		•					
	COMBI PAL (200 Pos. Tray) GC PAL (200 pos. Tray)	•							•	
	COMBI PAL (98 Pos. Tray) GC PAL (98 Pos. Tray)	•	•		•				•	
	COMBI PAL SPME mode (98 Pos. Tray)		•		•				•	
	COMBI PAL (32 Pos. Tray) GC PAL (32 Pos. Tray), COMBI PAL SPME mode (32 Pos. Tray)								•	
	Genesis								•	
	Marathon Basic, Standard 96 Pos. Tray		•	•	•					
	Marathon Basic, Prep, King Size 48 Pos. Tray								•	
	Vista			•	•					
	CP-9020/CP-9025, CP-9060								•	
	CP-9010		•	•	•					
	CP-8410/8034/8035/8100/8200/9095/9100		•	•	•					
	920-LC/940-LC		•	•						

Autosampler Compatibility Table (con't)

Manufacturer	Model	8 mm Crimp	11 mm Crimp and Snap	8 mm Screw	9 mm Screw	10 mm Screw	13 mm Screw	18 mm Screw (Headspace)	Shell Vials (Neckless)	Plate (page 49)
Viscotek	GPC Autosampler			•	•	•				
VWR(Merck)/Hitachi	L2200 (LaChrom Elite)/L2200-U (LaChrom Ultra) (200 Pos. Tray), L7200 (LaChrom) (80 Pos. Tray)/L7250(LaChrom) (Pos. Tray)		•	•	•					
	L2200 (LaChrom Elite) (128 Pos. Tray)						•			
	L7250 (LaChrom) (Rack Holder for combination Racks)	•	•	•	•		•			
	655-A40 (108 Pos. Tray), L-9100, AS 2000 (50 Pos. Tray), AS 4000 (150 Pos. Tray)		•	•	•					
	AS 4000 (198 Pos. Tray)	•								
	AS 6000	•	•	•	•					
Waters	Acquity Sample Organizer		•		•					•
	Acquity/CapLC/Waters/Nano Acquity		•		•					•
	Alliance HTS									•
	Model 2767		•	•						•
	Model 2707		•	•						•
	Model 2777		•	•						•
	ACQUITY™ UPLC Systems				•		•			
	Wisp 48 position						•		•	
	Wisp 96 position, 717, 96 Position Carousel								•	
	717, 48 Position Carousel						•		•	
	Alliance®, Alliance HT Syst.		•		•	•				
	Alliance® GPC 2000						•	•		
	Alliance® 2790/2795, Alliance 2690/2695		•		•	•				

- indicates that a cap having an outer flange is required for the vial to operate properly with the autosampler.
- indicates that the vials from this category are compatible with the autosampler in most configurations.
- indicates that a magnetic seal is required for use with the autosampler.



Contact Us
Order Now

Ordering Information

Quote Form

General Information

Company Name:

SiliCycle Customer Number (*if known*):

Contact Name:

E-mail:

Tel: Fax:

Address Information

.....

Department:

Address:

ZIP/Postal Code:

Country:

Quote Request		
Product Number	Product Description	Requested Quantity

Please confirm this request by:

- FAX
- E-mail
- Phone

Comments

.....

.....

.....

.....

Please Copy and Fax this form to SiliCycle Inc. at **418-874-0355**



SiliCycle®



Information Request Form

- | | | |
|--|---|---|
| <input type="checkbox"/> SiliaPrep™ Silica-Based SPE | <input type="checkbox"/> SiliaPlate™ TLC Plates | <input type="checkbox"/> SiliCycle Syringe Filters |
| <input type="checkbox"/> SiliaPrepX™ Polymeric SPE | <input type="checkbox"/> SiliaPlate™ TLC Accessories | <input type="checkbox"/> SiliCycle Membrane Filters |
| <input type="checkbox"/> SiliaPrepMB™ SPE for MiniBlock® | <input type="checkbox"/> SiliaSphere™ Spherical Silica Gels | <input type="checkbox"/> SiliCycle Vials & Caps |
| <input type="checkbox"/> SiliaPrep™ Accessories | <input type="checkbox"/> SiliaChrom® HPLC Column | <input type="checkbox"/> Other: _____ |
| <input type="checkbox"/> SiliaPrep™ Tips Micro-SPE | <input type="checkbox"/> SiliaChrom® Guard Cartridge | |
| <input type="checkbox"/> SiliaPrep™ QuEChERS | | |

How Can We Help You?

.....

.....

.....

.....

.....

.....

.....

General Information

Company Name:

SiliCycle Customer Number (if known):

Contact Name:

E-mail:

Tel: Fax:

Address Information

.....

Department:

Address:

ZIP/Postal Code:

Country:

Please Copy and Fax this form to SiliCycle Inc. at **418-874-0355**

Order Form

General Information

Company Name:
 SiliCycle Client Number: Purchase Order Number:
 Contact Name:
 E-mail address:
 Tel: Fax:

Payment Details

Verification Code*:

Credit card number: *The 3 last numbers in the back side of your card

Expiry date:

VISA Master Card Amex

Name as it appears on card:



Address Information for the Shipment

Department:
 Address:

 ZIP/Postal Code:
 Country:

SiliCycle Products				
Part Number	Product Description	Packaging (g, Kg, box)	Quantity	Price (USD)

Confirm this order by

- FAX
- E-mail
- Phone

Comments

.....

..... Please Copy and Fax this form to SiliCycle Inc. at 418-874-0355



Ordering Information

How to order

You can order any SiliCycle product on-line through the new SiliCycle e-commerce website at www.SiliCycle.com.

Orders can also be placed by phone, fax, mail or e-mail. You will find an order form on page 166 in the present catalog for fax (1 418.874.0355) and mail (*SiliCycle headquarters address at the bottom of this page*) orders. If you prefer, you can reach us by e-mail (info@silicycle.com) or by phone (1 418.874.0054 or Toll free for North America only 1 877.745.4292). Please have the following information on hand:

- Your name
- Company name, billing and shipping address
- Purchase order number
- Credit card information
- Product number and description
- Size, quantity and unit of measure
- E.I.N. or F.I.N. (*for United States clients only*)

Technical Support


At SiliCycle, we are committed to providing the best technical support possible. Our worldwide Technical Support Group is comprised of a team of highly qualified M.Sc., Engineers and PhD Chemists, Technical Support Professionals and Service Coordinators who are prepared to troubleshoot, answer questions, and provide solutions for your service and applications needs.

In order to better respond to your technical inquiries, feel free to contact us in three different ways:

- E-mail: support@silicycle.com
- Phone: International 1 418.874.0054
North America 1 877.745.4292 (*Toll-Free*)
- Online forum at www.silicycle.com

SiliCycle headquarters address:

2500, Parc-Technologique Blvd
Quebec City, Quebec
G1P 4S6, CANADA



Terms and Conditions

General

Unless otherwise stated, all transactions are expressly subject to these Terms and Conditions. Modifications or additions will be recognized only if accepted in writing by an officer of SiliCycle Inc. (*hereinafter named SiliCycle*), or an officially designated representative. Provisions of Buyer's Purchase Order or other documents that add to or differ from these Terms and Conditions are expressly rejected. No waiver of these Terms and Conditions or acceptance of others shall be construed as failure of the Company to raise objections.

Privacy Policy

Because your clientele is our most vital asset, we take privacy very seriously and won't share your personal information with anyone. Your information is used only to personalize your profile and to facilitate the transaction. You can change or update your information at any time.

Quotation and Published Prices

Quotations automatically expire 30 calendar days from the date issued unless otherwise stated. Quotes are subject to withdrawal with notice within that period. Prices shown on the published price lists and other published literature issued by SiliCycle are not unconditional offers to sell, and are subject to change without notice.

Warranty

SiliCycle guarantees to the original Buyer that the products sold conform to the composition and purity described therein at the time of their shipment. The Buyer's sole remedy in the event that SiliCycle fails to meet said warranty shall be the replacement of the unused portion of the product(s), or if approved by SiliCycle, a refund (*at the purchase price*) provided that the Buyer returns the alleged non-conforming product(s) within 30 days after reception of product(s). SiliCycle makes no other guarantee of suitability for a particular purpose or of the merchantability in the use or handling of the product, and does not accept any liability for consequential, special, indirect or incidental damages resulting therefrom.

Changes

The Buyer may, with the express written consent of SiliCycle, make changes in the specifications for products or work covered by the contract. In such an event, the contract price and delivery dates shall be equitably adjusted. SiliCycle shall be entitled to payment for reasonable profit plus costs and expenses incurred by work and materials rendered unnecessary as a result of such changes and for work and materials required to effect said changes.

If the Buyer has made a mistake on his/her purchase order, and the material has already been shipped and received, SiliCycle may approve the exchange of said material (*if price is identical*); however the Buyer will be responsible for all shipping costs. See return authorization policy section on the next page to obtain a return merchandise authorization form prior to returning goods.

Cancellation

Undelivered parts of any order may be cancelled by the Buyer only with the written approval of SiliCycle. If the Buyer makes an assignment for the benefit of creditors, or in the event that SiliCycle, for any reason feels insecure about Buyer's willingness or ability to perform, SiliCycle shall have the unconditional right to cancel the sales transaction or demand full or partial payment.

In the event of any cancellation of this order by either party, the Buyer shall pay SiliCycle for reasonable costs and expenses incurred by the SiliCycle prior to receipt of the cancellation notice, plus SiliCycle's usual rate of profit for similar work.

Taxes

The Company's prices do not include any applicable sales, goods and services, use, excise or similar taxes and the amount of any such tax SiliCycle may be required to pay or collect will be added to each invoice and paid by the Buyer.

Terms of Payment

All merchandise purchased remains the property of SiliCycle until such time as all invoices for the merchandise have been paid in full. Except for purchases paid online, or unless explicitly stated elsewhere in writing, terms are cash net 30 days from date of invoice. Additional fees of 2% per month (*26.8% per year*) will accrue on all accounts past due. If any payment is in default, and it becomes necessary to hire a recovery agency or lawyer, the client accepts to pay, in addition to the outstanding balance, recovery fees equal to 20% of the balance in capital and interests. By reason of the financial condition of Buyer or otherwise, SiliCycle may require full or partial payment in advance.

Certain orders may require a deposit or progressive payments as referenced in the quote. Such deposits may be increased upon receipt of purchase order based upon the Buyer's most current credit rating. Subject to the warranties stated in this policy, all sales are final without right of return.



Return Policy

Our Customer Service Department is available to assist you at any time should a problem arise with your order. Please make sure to inspect your packages immediately upon receipt and notify us within the next two (2) business days of any damage and/or discrepancies. Should a product be sent to you incorrectly, as the result of an error on our part, we will take quick and appropriate action to correct the problem at no charge to you.

In order to maintain the quality of our products and continue to provide competitive prices, some products may not be returned for credit. SiliCycle will not grant credit for:

- (i) Shelf-worn, used or defaced products;
- (ii) Scavengers, reagents, catalysts, or any other bounded silica whose containers have been opened;
- (iii) Products that are personalized or customized;
- (iv) Refrigerated or temperature-controlled products;
- (v) Products that have been discontinued;
- (vi) Products not directly purchased from SiliCycle

Products sold in distribution by SiliCycle will be subject to the Terms and Conditions Policy of the respective manufacturer.

Prior to any return, an authorization and a return material authorization (*RMA*) number must be obtained from our Customer Service Department. Shipping instructions will also be provided at this point. The RMA will ensure the safe and proper handling of material; it should therefore be referenced on all shipping labels.

The Buyer has 30 days from the issuance of the RMA to return the goods. Returns made without an authorization number will not be accepted and will be returned to the Buyer.

Returns are subject to a 40% restocking and/or disposal fee.

Shipping Policy

SiliCycle uses a two-day or five-day delivery (*or equivalent*) depending on weight and availability of product. Standard overnight delivery can also be arranged. Freight charges are prepaid and added to the invoice unless special instructions are requested by the customer. These conditions apply to all North American shipments. International delivery delays will vary according to orders and destination countries.

Delivery

Delivery dates indicated in the contract documents are approximate and based on prompt receipt of all necessary information regarding the product covered by the contract. SiliCycle will use reasonable efforts to meet the indicated delivery dates, but cannot be held responsible for its failure to do so.

In the event of any delivery delay caused by the Buyer, SiliCycle will store and handle all items ordered at Buyer's risk and will invoice Buyer for the unpaid portion of the contract price, plus storage, insurance, and handling charges on or after the date on which the product is ready for delivery. The invoice will be payable in full within 30 days from the invoice date, unless otherwise expressly agreed to in writing by SiliCycle.

SiliCycle will not hold orders unless specifically approved. SiliCycle has the right to make partial shipments and bill for those shipments; the buyer will make payment in accordance with the terms mentioned in this policy.

Shipping and Handling Charges

Shipping charges plus the applicable company handling charges will be prepaid and billed as a separate item on the product invoice. Title to the product and risk of loss shall pass to Buyer upon delivery to a carrier.

Application

All products are sold for laboratory or manufacturing uses. Only professional laboratory staff should handle the chemicals.

Product Numbers Index

0-25SW-PP9-C	146	0880-R0125-C	152	6CPTF-C20-C	144	ALCP11-CI-M	143	AUT-0188	48
0-2IF-C8-M	151	0880-R0125-D	152	6CPTR-C20-C	144	ALCP11-CIM-C	143	AUT-0189	48
0-2IT-C8-D	151	0-8CP-C8-D	142	96M-0001R	48	ALCP11-CIM-M	143	AUT-0190	48
0-2IT-C9-M	151	0-8CPZ-C8-D	142	96M-0001R-100	48	ALCP11-JIM-C	143	AUT-0193	48
0321-H0060-C	145	0-9CP-C11-D	142	96M-0001R-25	48	ALCP11-JIM-M	143	AUT-0194	48
0321-H0060-M	145	0-9CPT-C11-D	142	96M-0001R-50	48	ALCP11-MI-C	143	AUT-0195	48
0321-HS0060-C	145	0-9SW-C9-D	146	96M-0001S	48	ALCP11-MI-M	143	AUT-0196	48
0321-HS0060-M	145	10CP-C20-C	144	96M-0001S-100	48	ALCP11-MIM-C	143	AUT-0197	48
0321-J0040-C	145	10CPTF-C20-C	144	96M-0001S-25	48	ALCP11-MIM-M	143	AUT-0198	48
0321-J0040-M	145	10CPTR-C20-C	144	96M-0001S-50	48	ALCP20-BI-C	144	AUT-0199	48
0321-J0060-C	145	10SW-A18-C	150	96M-0012	47	ALCP20-BI-M	144	AUT-0200	48
0321-J0060-M	145	10SW-C18-C	150	96W-0009	47	ALCP20-M-M	144	AUT-0201	48
0321-K0060-C	145	1-1CPT-A11-D	143	96W-0010	47	ALCP20-NI-C	144	AUT-0202	48
0321-K0060-M	145	1-1CPT-C11-D	143	96W-0011	47	ALCP20-NI-M	144	AUT-0203	48
0321-KS0060-C	145	1-2CP-C8-D	142	96W-AUT-0014-B	18	ALCP8-CI-C	142	AUT-0204	48
0321-KS0060-M	145	1-5HCP-C11-D	143	96W-AUT-0014-C	18	ALCP8-CI-M	142	AUT-0205	48
0346-H0040-C	146	1-5HSP-C11-D	141	96W-AUT-0015-B	18	ALCP8-JI-C	142	AUT-0206	48
0346-H0040-M	146	1-5HSW-A9-C	146	96W-AUT-0015-C	18	ALCP8-JI-M	142	AUT-0207	48
0346-J0040-C	146	1-5HSW-C9-C	146	96W-P0001-1A	28	ALCP8-MI-C	142	AUT-0208	48
0346-J0040-M	146	1CPT-C8-D	142	96W-P0001-AA	28	ALCP8-MI-M	142	AUT-0209	48
0346-K0040-C	146	1SLF-C8-D	151	96W-P0002-1A	28	ALSW18-IIM-C	150	AUT-0210	48
0346-K0040-M	146	20CP-C20-C	144	96W-P0002-AA	28	ALSW18-IIM-M	150	AUT-0211	48
0346-P0040-C	146	20CPTF-C20-C	144	96W-P0005-1A	29	ALSW18-LIM-C	150	AUT-0212	48
0346-P0040-M	146	20CPTR-C20-C	144	96W-P0005-AA	29	ALSW18-LIM-M	150	AUT-0213	48
0346-PS0040-C	146	20SW-A18-C	152	96W-P0010-1A	29	ALSW18-QIM-C	150	AUT-0214	48
0346-PS0040-M	146	20SW-A24-C	152	96W-P0010-AA	29	ALSW18-QIM-M	150	AUT-0215-12	48
0374-H0060-C	148	20SW-C18-C	152	96W-P0015-1A	30	ALSW18-ZQIM-C	150	AUT-0215-16	48
0374-H0060-M	148	20SW-C24-C	152	96W-P0015-AA	30	ALSW18-ZQIM-M	150	AUT-0215-24	48
0374-HS0060-C	148	2-5CP-C11-C	143	96W-P0020-1A	30	AUT-0043	44	BCP11-CI-C	143
0374-HS0060-M	148	2CP-A11-C	143	96W-P0020-AA	30	AUT-0044	44	BCP11-CI-M	143
0-3CP-A11-D	143	2CP-AW11-C	143	96W-R10029G-B	17	AUT-0045	44	BMCP20-DIM-C	144
0-3CP-C11-D	143	2CP-C11-C	143	96W-R10029G-C	17	AUT-0046	44	BMCP20-DIM-M	144
0-3IF-C9-M	151	2CP-CW11-C	143	96W-R10030B-B	17	AUT-0047	44	BMCP20-NIM-C	144
0-3SP-A11-D	141	2CPR-CW11-D	143	96W-R10030B-C	17	AUT-0128-12	45-46	BMCP20-NIM-M	144
0-3SP-C11-D	141	2SLF-C12-D	151	96W-R31030B-B	16	AUT-0128-16	45-46	BSP11-AI-C	141
0-3SW-A9-D	146	2SP-A11-C	141	96W-R31030B-C	16	AUT-0129-24	45-46	BSP11-AI-M	141
0-3SW-C9-D	146	2SP-AW11-C	141	96W-R31130B-B	16	AUT-0130-10	45-46	BSP11-KSI-C	141
0479-H0060-C	149	2SP-C11-C	141	96W-R31130B-C	16	AUT-0146	48	BSP11-KSI-M	141
0479-H0060-M	149	2SP-CW11-C	141	96W-R31930B-B	14	AUT-0147	48	BSP11-MI-C	141
0479-H0075-C	149	2SW-A10-C	148	96W-R31930B-C	14	AUT-0149-25	48	BSP11-MI-M	141
0479-H0075-M	149	2SW-A8-C	145	96W-R33229G-B	14	AUT-0149-50	48	BSW9-BI-C	147
0479-HS0060-C	149	2SW-A9-C	146	96W-R33229G-C	14	AUT-0154	48	BSW9-BI-M	147
0479-HS0060-M	149	2SW-AW10-C	148	96W-R34030B-B	16	AUT-0155	48	BSW9-C	146
0479-HS0075-C	149	2SW-AW8-C	145	96W-R34030B-C	16	AUT-0160	48	BSW9-EI-C	147
0479-HS0075-M	149	2SW-AW9-C	146	96W-R35030B-B	17	AUT-0161	48	BSW9-EI-M	147
0680-I0050-C	150	2SW-C10-C	148	96W-R35030B-C	17	AUT-0162	48	BSW9-HB-C	147
0680-I0050-M	150	2SW-C8-C	145	96W-R35530B-B	14	AUT-0163	48	BSW9-HB-M	147
0680-L0060-C	150	2SW-C9-C	146	96W-R35530B-C	14	AUT-0164	48	BSW9-HI-C	147
0680-L0060-M	150	2SW-CW10-C	148	96W-R38030B-B	16	AUT-0167	48	BSW9-HI-M	147
0680-Q0050-C	150	2SW-CW8-C	145	96W-R38030B-C	16	AUT-0168	48	BSW9-HSI-C	147
0680-Q0050-M	150	2SW-CW9-C	146	96W-R51230B-B	20	AUT-0172	48	BSW9-HSI-M	147
0680-Q0125-C	150	30SW-A24-C	152	96W-R51230B-C	20	AUT-0173	48	BSW9-JI-C	147
0680-Q0125-M	150	30SW-C24-C	152	96W-R52030B-B	19	AUT-0174	48	BSW9-JI-M	147
0-7CP-A8-D	142	3-5HSW-C13-D	149	96W-R52030B-C	19	AUT-0175	48	BSW9-KI-C	147
0-7CP-C8-D	142	40SW-A24-C	152	96W-R60530B-B	20	AUT-0176	48	BSW9-KI-M	147
0-7SW-PP9-C	146	40SW-C24-C	152	96W-R60530B-C	20	AUT-0178	48	BSW9-KSB-C	147
0800-B0125-C	144	4CP-C11-C	143	96W-R66430B-B	19	AUT-0179	48	BSW9-KSB-M	147
0800-B0125-M	144	4SW-A13-C	149	96W-R66430B-C	19	AUT-0182	48	BSW9-KSI-C	147
0800-N0125-C	144	4SW-AW13-C	149	96W-R66530B-B	19	AUT-0183	48	BSW9-KSI-M	147
0800-N0125-M	144	4SW-C13-C	149	96W-R66530B-C	19	AUT-0184	48	BSW9-M	146
0880-G0125-C	152	4SW-C9-D	146	96W-R70030B-B	20	AUT-0185	48	BSW9-OB-C	147
0880-R0060-C	152	4SW-CW13-C	149	96W-R70030B-C	20	AUT-0186	48	BSW9-OB-M	147
0880-R0060-D	152	6CP-C20-C	144	ALCP11-CI-C	143	AUT-0187	48	BSW9-PI-C	147

BSW9-PI-M	147	NSP11MI-C	141	QD-1004-2T	58	S10003E-A	79	S38002G-A	83
BSW9-PSI-C	147	NSP11MI-M	141	QD-1005-2T	58	S10003F-A	79	S38003B	83
BSW9-PSI-M	147	NSW10-C	148	QD-1006-2T	58	S10003G-A	79	S38003E-A	83
CSL8-M	151	NSW10-HI-C	148	QD-1007-2T	58	S10003M	79	S38003F-A	83
CSL12-M	151	NSW10-HI-M	148	QD-2000-15T	58	S10005B	79	S38003G-A	83
CSP11-JI-C	141	NSW10-HSI-C	148	QD-2001-15T	58	S10005E-A	79	S38003M	83
CSP11-JI-M	141	NSW10-HSI-M	148	QD-2002-15T	58	S10005F-A	79	S38005B	83
CSP11-KSI-C	141	NSW10-M	148	QD-2003-15T	58	S10005G-A	79	S38005E-A	83
CSP11-KSI-M	141	NSW13-C	149	QD-2003-15T	58	S10005M	79	S38005F-A	83
CSP11-MI-C	141	NSW13-HB-C	149	QD-2004-15T	58	S10006E-A	79	S38005G-A	83
CSP11-MI-M	141	NSW13-HB-M	149	QD-2006-15T	58	S10006G-A	79	S38005M	83
FLH-R10029G-15U	17	NSW13-HI-C	149	QD-2007-15T	58	S10007B	79	S38006E-A	83
FLH-R10029G-25X	17	NSW13-HI-M	149	QE-0001-50T	57	S10007E-A	79	S38006G-A	83
FLH-R10030B-15U	17	NSW13-HSB-C	149	QE-0002-50T	57	S10007F-A	79	S38007B	83
FLH-R10030B-25X	17	NSW13-HSB-M	149	QE-0003-50T	57	S10007G-A	79	S38007E-A	83
GCP11-CI-C	143	NSW13-HSI-C	149	RCP11-CI-C	143	S10007M	79	S38007F-A	83
GCP11-CI-M	143	NSW13-HSI-M	149	RCP11-CI-M	143	S10007T	79	S38007G-A	83
GOCP20-5DIM-C	144	NSW13-M	149	RSP11-MI-C	141	S10008B	79	S38007M	83
GOCP20-5DIM-M	144	NSW13-ZHI-C	149	RSP11-MI-M	141	S10008E-A	79	S38007T	83
GOCP20-5MIM-C	144	NSW13-ZHI-M	149	RSW9-C	146	S10008F-A	79	S38008B	83
GOCP20-5MIM-M	144	NSW13-ZHSI-C	149	RSW9-HI-C	147	S10008G-A	79	S38008E-A	83
GOCP20-5M-M	144	NSW13-ZHSI-M	149	RSW9-HI-M	147	S10008M	79	S38008F-A	83
GOCP20-5NIM-C	144	NSW24-RB-C	152	RSW9-M	146	S10008T	79	S38008G-A	83
GOCP20-5NIM-M	144	NSW24-RB-D	152	S03201E-A	83	S10009E-A	79	S38008M	83
GOCP20-8DIM-C	144	NSW24-RI-C	152	S03201F-A	83	S10009F-A	79	S38008T	83
GOCP20-8DIM-M	144	NSW24-RI-D	152	S03201G-A	83	S10009G-A	79	S38009E-A	83
GOCP20-8MIM-C	144	NSW8-C	145	S03202E-A	83	S3801E-A	83	S38009F-A	83
GOCP20-8MIM-M	144	NSW8-HI-C	145	S03202F-A	83	S3801F-A	83	S38009G-A	83
GOCP20-8M-M	144	NSW8-HI-C	145	S03202G-A	83	S3801G-A	83	S52001E-A	83
GOCP20-8NIM-C	144	NSW8-HSI-C	145	S03203B	83	S3802E-A	83	S52001F-A	83
GOCP20-8NIM-M	144	NSW8-HSI-M	145	S03203E-A	83	S3802F-A	83	S52001G-A	83
GSW9-C	146	NSW8-JI-C	145	S03203F-A	83	S3802G-A	83	S52002E-A	83
GSW9-HI-C	147	NSW8-JI-M	145	S03203G-A	83	S3803B	83	S52002F-A	83
GSW9-HI-M	147	NSW8-KI-C	145	S03203M	83	S3803E-A	83	S52002G-A	83
GSW9-M	146	NSW8-KI-M	145	S03205B	83	S3803F-A	83	S52003B	83
HDW-000	123	NSW8-KSI-C	145	S03205E-A	83	S3803G-A	83	S52003E-A	83
HDW-001	123	NSW8-KSI-M	145	S03205F-A	83	S3803M	83	S52003F-A	83
HDW-002	123	NSW8-M	145	S03205G-A	83	S3805B	83	S52003G-A	83
HDW-003	123	NSW8-ZJI-C	145	S03205M	83	S3805E-A	83	S52003M	83
HPTLG-R10011B-1010	67	NSW8-ZJI-M	145	S03206E-A	83	S3805F-A	83	S52005B	83
HPTLG-R10011B-2020	67	NSW9-C	146	S03206G-A	83	S3805G-A	83	S52005E-A	83
HPTLGSR10011B-1010	67	NSW9-EI-C	147	S03207B	83	S3805M	83	S52005F-A	83
HPTLGSR10011B-1020	67	NSW9-EI-M	147	S03207E-A	83	S3806E-A	83	S52005G-A	83
MF-NYL25-20	139	NSW9-ESI-C	147	S03207F-A	83	S3806G-A	83	S52005M	83
MF-NYL25-45	139	NSW9-ESI-M	147	S03207G-A	83	S3807B	83	S52006E-A	83
MF-NYL47-20	139	NSW9-HI-C	147	S03207M	83	S3807E-A	83	S52006G-A	83
MF-NYL47-45	139	NSW9-HI-M	147	S03207T	83	S3807F-A	83	S52007B	83
MF-PES25-20	139	NSW9-KI-C	147	S03208B	83	S3807G-A	83	S52007E-A	83
MF-PES25-45	139	NSW9-KI-M	147	S03208E-A	83	S3807M	83	S52007F-A	83
MF-PES47-20	139	NSW9-KSI-C	147	S03208F-A	83	S3807T	83	S52007G-A	83
MF-PES47-45	139	NSW9-KSI-M	147	S03208G-A	83	S3808B	83	S52007M	83
MF-PTF25-20	139	NSW9-M	146	S03208M	83	S3808E-A	83	S52007T	83
MF-PTF25-45	139	NSW9-PI-C	147	S03208T	83	S3808F-A	83	S52008B	83
MF-PTF47-20	139	NSW9-PI-M	147	S03209E-A	83	S3808G-A	83	S52008E-A	83
MF-PTF47-45	139	NSW9-PSI-C	147	S03209F-A	83	S3808M	83	S52008F-A	83
MF-PVD25-20	139	NSW9-PSI-M	147	S03209G-A	83	S3808T	83	S52008G-A	83
MF-PVD25-45	139	PS-012	49	S10001E-A	79	S3809F-A	83	S52008M	83
MF-PVD47-20	139	PS-060	49	S10001F-A	79	S3809G-A	83	S52008T	83
MF-PVD47-45	139	PS-150	49	S10001G-A	79	S3801E-A	83	S52009E-A	83
MF-RC25-20	139	QD-1000-2T	58	S10002E-A	79	S38001F-A	83	S52009F-A	83
MF-RC25-45	139	QD-1001-2T	58	S10002F-A	79	S38001G-A	83	S52009G-A	83
MF-RC47-20	139	QD-1002-2T	58	S10002G-A	79	S38002E-A	83	SF-NYL13-20	137
MF-RC47-45	139	QD-1003-2T	58	S10003B	79	S38002F-A	83	SF-NYL13-20-L	137

Product Numbers Index (con't)

SF-NYL13-20-M	137	SIM-0008-003	45	SPE-AUT-0015-06U	18	SPE-PO010-01AA	29	SPE-R33229G-12U	14
SF-NYL13-45	137	SIM-0009-150	45	SPE-AUT-0015-12U	18	SPE-PO010-03BB	29	SPE-R33229G-20X	14
SF-NYL13-45-L	137	SPC-AUT-0014-10G	18	SPE-AUT-0015-20X	18	SPE-PO010-06C	29	SPE-R34030B-01B	16
SF-NYL13-45-M	137	SPC-AUT-0014-10P	18	SPE-AUT-0053-01B	18	SPE-PO010-06G	29	SPE-R34030B-01C	16
SF-NYL25-20	137	SPC-AUT-0015-10G	18	SPE-AUT-0053-01C	18	SPE-PO010-06P	29	SPE-R34030B-03G	16
SF-NYL25-20-L	137	SPC-AUT-0015-10P	18	SPE-AUT-0053-03G	18	SPE-PO015-01AA	30	SPE-R34030B-03P	16
SF-NYL25-20-M	137	SPC-AUT-0053-10G	18	SPE-AUT-0053-03P	18	SPE-PO015-03BB	30	SPE-R34030B-06P	16
SF-NYL25-45	137	SPC-AUT-0053-10P	18	SPE-AUT-0053-06P	18	SPE-PO015-06C	30	SPE-R34030B-06S	16
SF-NYL25-45-L	137	SPC-AUT-0054-10G	18	SPE-AUT-0053-06S	18	SPE-PO015-06G	30	SPE-R34030B-06U	16
SF-NYL25-45-M	137	SPC-AUT-0054-10P	18	SPE-AUT-0053-06U	18	SPE-PO015-06P	30	SPE-R34030B-12U	16
SF-PES13-20	137	SPC-AUT-0055-10G	18	SPE-AUT-0053-12U	18	SPE-PO020-01AA	30	SPE-R34030B-20X	16
SF-PES13-20-L	137	SPC-AUT-0055-10P	18	SPE-AUT-0053-20X	18	SPE-PO020-03BB	30	SPE-R35030B-01B	17
SF-PES13-20-M	137	SPC-R10029G-10G	17	SPE-AUT-0054-01B	18	SPE-PO020-06C	30	SPE-R35030B-01C	17
SF-PES13-45	137	SPC-R10029G-10P	17	SPE-AUT-0054-01C	18	SPE-PO020-06G	30	SPE-R35030B-03G	17
SF-PES13-45-L	137	SPC-R10030B-10G	17	SPE-AUT-0054-03G	18	SPE-PO020-06P	30	SPE-R35030B-03P	17
SF-PES13-45-M	137	SPC-R10030B-10P	17	SPE-AUT-0054-03P	18	SPE-R10029G-01B	17	SPE-R35030B-06P	17
SF-PES25-20	137	SPC-R31030B-10G	16	SPE-AUT-0054-06P	18	SPE-R10029G-01C	17	SPE-R35030B-06S	17
SF-PES25-20-L	137	SPC-R31030B-10P	16	SPE-AUT-0054-06S	18	SPE-R10029G-03G	17	SPE-R35030B-06U	17
SF-PES25-20-M	137	SPC-R31130B-10G	16	SPE-AUT-0054-06U	18	SPE-R10029G-03P	17	SPE-R35030B-12U	17
SF-PES25-45	137	SPC-R31130B-10P	16	SPE-AUT-0054-12U	18	SPE-R10029G-06P	17	SPE-R35030B-20X	17
SF-PES25-45-L	137	SPC-R31930B-10G	14	SPE-AUT-0054-20X	18	SPE-R10029G-06S	17	SPE-R35530B-01B	14
SF-PES25-45-M	137	SPC-R31930B-10P	14	SPE-AUT-0055-01B	18	SPE-R10029G-06U	17	SPE-R35530B-01C	14
SF-PTF13-20	137	SPC-R33229G-10G	14	SPE-AUT-0055-01C	18	SPE-R10030B-01B	17	SPE-R35530B-03G	14
SF-PTF13-20-L	137	SPC-R33229G-10P	14	SPE-AUT-0055-03G	18	SPE-R10030B-01C	17	SPE-R35530B-03P	14
SF-PTF13-20-M	137	SPC-R34030B-10G	16	SPE-AUT-0055-03P	18	SPE-R10030B-03G	17	SPE-R35530B-06P	14
SF-PTF13-45	137	SPC-R34030B-10P	16	SPE-AUT-0055-06P	18	SPE-R10030B-03P	17	SPE-R35530B-06S	14
SF-PTF13-45-L	137	SPC-R35030B-10G	17	SPE-AUT-0055-06S	18	SPE-R10030B-06P	17	SPE-R35530B-06U	14
SF-PTF13-45-M	137	SPC-R35030B-10P	17	SPE-AUT-0055-06U	18	SPE-R10030B-06S	17	SPE-R35530B-12U	14
SF-PTF25-20	137	SPC-R35530B-10G	14	SPE-AUT-0055-12U	18	SPE-R10030B-06U	17	SPE-R35530B-20X	14
SF-PTF25-20-L	137	SPC-R35530B-10P	14	SPE-AUT-0055-20X	18	SPE-R31030B-01B	16	SPE-R38030B-01B	16
SF-PTF25-20-M	137	SPC-R38030B-10G	16	SPEC-R31930B-01B	24	SPE-R31030B-01C	16	SPE-R38030B-01C	16
SF-PTF25-45	137	SPC-R38030B-10P	16	SPEC-R31930B-01C	24	SPE-R31030B-03G	16	SPE-R38030B-03P	16
SF-PTF25-45-L	137	SPC-R51230B-10G	18	SPEC-R31930B-03G	24	SPE-R31030B-03P	16	SPE-R38030B-06P	16
SF-PTF25-45-M	137	SPC-R51230B-10P	18	SPEC-R31930B-03P	24	SPE-R31030B-06P	16	SPE-R38030B-06S	16
SF-PVD13-20	137	SPC-R52030B-10G	19	SPEC-R31930B-06P	24	SPE-R31030B-06S	16	SPE-R38030B-06U	16
SF-PVD13-20-L	137	SPC-R52030B-10P	19	SPEC-R31930B-06S	24	SPE-R31030B-06U	16	SPE-R38030B-12U	16
SF-PVD13-20-M	137	SPC-R60530B-10G	18	SPEC-R31930B-06U	24	SPE-R31030B-12U	16	SPE-R38030B-20X	16
SF-PVD13-45	137	SPC-R60530B-10P	18	SPEC-R31930B-12U	24	SPE-R31030B-20X	16	SPE-R51230B-01B	18
SF-PVD13-45-L	137	SPC-R661230B-10G	21	SPEC-R31930B-20X	24	SPE-R31130B-01B	16	SPE-R51230B-01C	18
SF-PVD13-45-M	137	SPC-R661230B-10P	21	SPEC-R651230B-01B	22	SPE-R31130B-01C	16	SPE-R51230B-03G	18
SF-PVD25-20	137	SPC-R66430B-10G	19	SPEC-R651230B-01C	22	SPE-R31130B-03G	16	SPE-R51230B-03P	18
SF-PVD25-20-L	137	SPC-R66430B-10P	19	SPEC-R651230B-03G	22	SPE-R31130B-03P	16	SPE-R51230B-06P	18
SF-PVD25-20-M	137	SPC-R66530B-10G	19	SPEC-R651230B-03P	22	SPE-R31130B-06P	16	SPE-R51230B-06S	18
SF-PVD25-45	137	SPC-R66530B-10P	19	SPEC-R651230B-06P	22	SPE-R31130B-06S	16	SPE-R51230B-06U	18
SF-PVD25-45-L	137	SPC-R70030B-10G	20	SPEC-R651230B-06S	22	SPE-R31130B-06U	16	SPE-R51230B-12U	18
SF-PVD25-45-M	137	SPC-R70030B-10P	20	SPEC-R651230B-06U	22	SPE-R31130B-12U	16	SPE-R51230B-20X	18
SF-RC13-20	137	SPC-R802830B-10G	21	SPEC-R651230B-12U	22	SPE-R31130B-20X	16	SPE-R52030B-01B	19
SF-RC13-20-L	137	SPC-R802830B-10P	21	SPEC-R651230B-20X	22	SPE-R31930B-01B	14	SPE-R52030B-01C	19
SF-RC13-20-M	137	SPE-AUT-0014-01B	18	SPE-PO001-01AA	28	SPE-R31930B-01C	14	SPE-R52030B-03G	19
SF-RC13-45	137	SPE-AUT-0014-01C	18	SPE-PO001-03BB	28	SPE-R31930B-03G	14	SPE-R52030B-03P	19
SF-RC13-45-L	137	SPE-AUT-0014-03G	18	SPE-PO001-06C	28	SPE-R31930B-03P	14	SPE-R52030B-06P	19
SF-RC13-45-M	137	SPE-AUT-0014-03P	18	SPE-PO001-06G	28	SPE-R31930B-06P	14	SPE-R52030B-06S	19
SF-RC25-20	137	SPE-AUT-0014-06P	18	SPE-PO001-06P	28	SPE-R31930B-06S	14	SPE-R52030B-06U	19
SF-RC25-20-L	137	SPE-AUT-0014-06S	18	SPE-PO002-01AA	28	SPE-R31930B-06U	14	SPE-R52030B-12U	19
SF-RC25-20-M	137	SPE-AUT-0014-06U	18	SPE-PO002-03BB	28	SPE-R31930B-12U	14	SPE-R52030B-20X	19
SF-RC25-45	137	SPE-AUT-0014-12U	18	SPE-PO002-06C	28	SPE-R31930B-20X	14	SPE-R60530B-01B	20
SF-RC25-45-L	137	SPE-AUT-0014-20X	18	SPE-PO002-06G	28	SPE-R33229G-01B	14	SPE-R60530B-01C	20
SF-RC25-45-M	137	SPE-AUT-0015-01B	18	SPE-PO002-06P	28	SPE-R33229G-01C	14	SPE-R60530B-03G	20
SIM-0002-006	45	SPE-AUT-0015-01C	18	SPE-PO005-01AA	29	SPE-R33229G-03G	14	SPE-R60530B-03P	0
SIM-0003-012	45	SPE-AUT-0015-03G	18	SPE-PO005-03BB	29	SPE-R33229G-03P	14	SPE-R60530B-06P	20
SIM-0004-020	45	SPE-AUT-0015-03P	18	SPE-PO005-06C	29	SPE-R33229G-06P	14	SPE-R60530B-06S	20
SIM-0006-060	45	SPE-AUT-0015-06P	18	SPE-PO005-06G	29	SPE-R33229G-06S	14	SPE-R60530B-06U	20
SIM-0007-001	45	SPE-AUT-0015-06S	18	SPE-PO005-06P	29	SPE-R33229G-06U	14	SPE-R60530B-12U	20

SPE-R60530B-20X	20	SPM-R802830B-06P	21	TLG-R35011B-203	67
SPE-R66430B-01B	19	SPM-R802830B-06S	21	TLG-R35011B-303	67
SPE-R66430B-01C	19	SPM-R802830B-06U	21	TLG-R38011B-203	67
SPE-R66430B-03G	19	SPM-R802830B-20X	21	TLG-R38011B-303	67
SPE-R66430B-03P	19	SPS-AUT-0014-P	18	TLG-R52011B-203	67
SPE-R66430B-06P	19	SPS-AUT-0014-S	18	TLG-R52011B-303	67
SPE-R66430B-06S	19	SPS-AUT-0015-P	18	TLGSR10011B-333	66
SPE-R66430B-06U	19	SPS-AUT-0015-S	18	TLGSR10011B-341	66
SPE-R66430B-12U	19	SPS-AUT-0053-P	18	TLGSR10011B-350	66
SPE-R66430B-20X	19	SPS-AUT-0053-S	18	TLGSR10011B-350	66
SPE-R66530B-01B	19	SPS-AUT-0054-P	18	TLGSR10011B-353	66
SPE-R66530B-01C	19	SPS-AUT-0054-S	18	TLGSR10011B-423	66
SPE-R66530B-03G	19	SPS-AUT-0055-P	18	TLGSR10011B-723	66
SPE-R66530B-03P	19	SPS-AUT-0055-S	18	TLP-R31001B-2575	66
SPE-R66530B-06P	19	SPS-R10029G-P	17	TLP-R31001B-323	66
SPE-R66530B-06S	19	SPS-R10029G-S	17	WSW10-C	148
SPE-R66530B-06U	19	SPS-R10030B-P	17	WSW10-HI-C	148
SPE-R66530B-12U	19	SPS-R10030B-S	17	WSW10-HI-M	148
SPE-R66530B-20X	19	SPS-R31930B-P	14	WSW10-M	148
SPE-R70030B-01B	20	SPS-R31930B-S	14	WSW24-GB-C	152
SPE-R70030B-01C	20	SPS-R33229G-P	14	WSW24-GB-D	152
SPE-R70030B-03G	20	SPS-R33229G-S	14	WSW24-GI-C	152
SPE-R70030B-03P	20	SPS-R35030B-P	17	WSW24-GI-D	152
SPE-R70030B-06P	20	SPS-R35030B-S	17	WSW24-RB-C	152
SPE-R70030B-06S	20	SPS-R35530B-P	14	WSW24-RB-D	152
SPE-R70030B-06U	20	SPS-R35530B-S	14	WSW24-RI-C	152
SPE-R70030B-12	20	SPS-R51230B-P	18	WSW24-RI-D	152
SPE-R70030B-20X	20	SPS-R51230B-S	18	WSW24-ZRI-C	152
SPET-C18-T1	52	SPS-R52030B-P	19	WSW24-ZRI-D	152
SPET-C18-T2	52	SPS-R52030B-S	19	YCP11-CI-C	143
SPET-C18-T3	52	SPS-R60530B-P	18	YCP11-CI-M	143
SPET-C4-T1	52	SPS-R60530B-S	18	YSP11-MI-C	143
SPET-C4-T2	52	SPS-R66430B-P	19	YSP11-MI-M	143
SPET-C4-T3	52	SPS-R66430B-S	19	YSW9-C	146
SPET-C8-T1	52	SPS-R66530B-P	19	YSW9-HI-C	147
SPET-C8-T2	52	SPS-R66530B-S	19	YSW9-HI-M	147
SPET-C8-T3	52	SPS-R70030B-P	20	YSW9-M	146
SPET-CB-T1	52	SPS-R70030B-S	20		
SPET-CB-T2	52	TLA-R10011B-2575	66		
SPET-CB-T3	52	TLA-R10011B-323	66		
SPET-TI-T1	52	TLA-R30411B-303	66		
SPET-TI-T2	52	TLGCZ-R10011B-323	66		
SPET-TI-T3	52	TLG-R10011B-124	66		
SPET-TIZR-T1	52	TLG-R10011B-2020	66		
SPET-TIZR-T2	52	TLG-R10011B-2575B	66		
SPET-TIZR-T3	52	TLG-R10011B-323	66		
SPET-ZR-T1	52	TLG-R10011B-333	66		
SPET-ZR-T2	52	TLG-R10011B-341	66		
SPET-ZR-T3	52	TLG-R10011B-353	66		
SPM-R661230B-01B	21	TLG-R10011B-423	66		
SPM-R661230B-01C	21	TLG-R10011B-424	66		
SPM-R661230B-03G	21	TLG-R10011B-527	66		
SPM-R661230B-03P	21	TLG-R10011B-624	66		
SPM-R661230B-06P	21	TLG-R10011B-723	66		
SPM-R661230B-06S	21	TLG-R23511B-303	67		
SPM-R661230B-06U	21	TLG-R23511B-423	67		
SPM-R661230B-12U	21	TLG-R30411B-213	67		
SPM-R661230B-20X	21	TLG-R30411B-303	67		
SPM-R802830B-12U	21	TLG-R30411B-341	67		
SPM-R802830B-01B	21	TLG-R31011B-203	67		
SPM-R802830B-01C	21	TLG-R31011B-303	67		
SPM-R802830B-03G	21	TLG-R32611B-203	67		
SPM-R802830B-03P	21	TLG-R32611B-303	67		

Applications Index

SiliaPrep SPE Cartridges and Well Plates

Determination of Testosterone in Human Urine	16
Δ^9 -Tetrahydrocannabinol in Human Plasma	17
Easy SPE Method for Drugs of Abuse Determination in Human Urine	24
Fentanyl and Norfentanyl in Urine	25
Easy SPE Method of Pesticides Determination from Drinking Water	26

SiliaPrepX Polymeric SPE Cartridges and Well Plates

Marbofloxacin & Sarafloxacin in Salmon	33
Acrylamide Determination in Fried Potato Chips	34
Determination of Carbendazim in Orange Juice	36
Amphetamines Quantification in Human Urine	37
Sulfonamides, Tetracyclines & Pyrimethamines in Milk	38
Determination of Trace Pesticides in Water	40
Caffeine, Cotinine & Nicotine in Human Urine	41

SiliaPrepMB SPE for MiniBlock Systems

Webinar on Reductive Amination	45
Webinar on Carbon-Carbon Coupling Reactions	45
Webinar on Metal Scavenging	45

SiliaPlate TLC Plates

Thin Layer Chromatography Practical Guide	71
---	----

SiliaChrom HPLC Columns

Separation of Catecholamines in Acidic Mobile Phase	101
Assay for QC Testing of Blood Pressure and Cholesterol Medication	102
Peak Shape Evaluation for Zwitterion Fluoroquinolones	102
Ropinirole and Amitriptyline Detection in Human Plasma	103
SiliaChrom AQ C18 is Highly Efficient for Basic Compounds	104
Evaluation of Resolution and Peak Shape	104
Retention Capacity of DMSO on SiliaChrom AQ C18	105
Dewetting Phenomena	106
Stability of SiliaChrom SB C18 at Low pH Conditions	108
Stability of SiliaChrom XT C18 Fidelity at High pH Conditions	110
Resolution and Peak Shape of Highly Hydrophobic Domestic Insecticide	111
Highly Base Deactivated C18 (SiliaChrom XDB1 C18)	112
Highly Base Deactivated C18 (SiliaChrom XDB2 C18)	114
SiliaChrom HILIC: Separation of Vitamin B Complex and Vitamin C	115
SiliaChrom Amylose T-DPC Enantiomeric Separation of L and D-val PMB	120
How to Select the Right SiliaChrom HPLC Column	128

Solvent Properties and Miscibility Chart




Your Partner of Choice for HPLC Analysis



Solvent Properties and Miscibility Chart

Solvent Strength	Polarity Index	UV Cutoff (nm)	Refractive Index	Viscosity (cP, 20°C)	Boiling Point (°C)	Water Solubility (w/w%)	Solvent
0.01	1.0	200	1.391	0.50	99	0.0002	Isooctane
0.04	0.0	200	1.410	0.92	174	0.01	n-Decane
0.05	0.1	200	1.407	0.44	49	0.01	Cycloheptane
0.1	1.0	220	1.402	0.45	78	0.11	1-Cyclobutane
0.21	2.1	220	1.397	0.64	142	0.19	n-Butyl Ether
0.28	2.4	220	1.368	0.37	68	0.62	Isopropyl Ether
0.42	3.1	233	1.424	0.44	40	1.6	Methylene Chloride
0.43	4.2	334	1.396	0.51	117	-	Methyl Butyl Ketone
0.47	4.7	320	1.451	2.00	156	-	Cyclohexanone
0.55	5.5	210	1.402	1.72	125	Miscible	Methoxyethanol
0.6	4.5	260	1.362	0.37	57	-	Methyl Acetate
0.64	6.0	380	1.344	0.67	101	2.1	Nitromethane
0.65	6.5	268	1.438	0.84	166	Miscible	N,N'-Dimethylacetamide
0.69	6.0	265	1.447	1.65	182	-	N-Methylformamide
1.11	6.9	210	1.432	19.9	198	Miscible	Ethylene Glycol
2	6.0	260	1.372	1.26	118	Miscible	Acetic Acid
0.56	5.1	330	1.359	0.36	56	Miscible	Acetone
0.65	5.8	190	1.344	0.38	82	Miscible	Acetonitrile
-	2.7	238	1.501	0.65	80	0.18	Benzene
0.39	3.9	215	1.399	2.98	117	7.8	n-Butanol
-	4.0	254	1.394	0.73	126	0.43	Butyl Acetate
-	1.6	265	1.460	0.97	77	0.08	Carbon Tetrachloride
0.4	4.1	245	1.446	0.57	61	0.815	Chloroform
0.04	0.2	200	1.427	1.00	81	0.01	Cyclohexane
-	3.5	228	1.445	0.79	83	0.81	1,2-Dichloroethane
-	3.1	235	1.424	0.44	40	1.3	Dichloromethane
0.64	6.4	268	1.431	0.92	153	Miscible	N,N'-Dimethylformamide
0.62	7.2	270	1.478	2.24	189	Miscible	Dimethyl Sulphoxide
0.56	4.8	220	1.422	1.37	101	Miscible	Dioxane
0.88	4.3	210	1.361	1.20	79	Miscible	Ethanol
0.58	4.4	260	1.372	0.45	77	8.7	Ethyl Acetate
-	2.8	218	1.352	0.23	35	6.89	Diethyl Acetate
0.01	0.1	200	1.388	0.40	98	0.0004	n-Heptane
0.01	0.1	200	1.375	0.31	69	0.0012	n-Hexane
0.95	5.1	205	1.329	0.55	65	Miscible	Methanol
0.35	2.5	220	1.369	0.27	55	4.8	Methyl-t-Butyl Ether
0.51	4.7	329	1.379	0.43	80	24	Methyl Ethyl Ketone
-	0.0	190	1.358	0.23	36	0.004	Pentane
0.82	4.0	210	1.385	2.30	97	Miscible	n-Pentanol
0.82	3.9	205	1.378	2.40	82	Miscible	Isopropanol
-	2.2	220	1.368	0.37	68	-	Diisopropyl Ether
0.45	4.0	212	1.407	0.55	66	Miscible	Tetrahydrofuran
0.29	2.4	285	1.496	0.59	111	0.05	Toluene
-	1.0	273	1.477	0.57	87	0.11	Trichloroethylene
2	10.2	190	1.000	1.00	100	-	Water
0.26	2.5	288	1.506	0.81	144	0.018	O-Xylene

■ Immiscible (2 phases are produced when both solvents are mixed)



- O-Xylene
- Water
- Trichloroethylene
- Toluene
- Tetrahydrofuran
- Diisopropyl Ether
- Isopropanol
- n-Propanol
- Pentane
- Methyl Ethyl Ketone
- Methyl-t-Butyl Ether
- Methanol
- Hexane
- Heptane
- Diethyl Ether
- Ethyl Acetate
- Ethanol
- Dioxane
- Dimethyl Sulphoxide
- Dimethylformamide
- Dichloromethane
- 1,2-Dichloroethane
- Cyclohexane
- Chloroform
- Carbon Tetrachloride
- Butyl Acetate
- n-Butanol
- Benzene
- Acetonitrile
- Acetone
- Acetic Acid

References: CRC Handbook of Chemistry and Physics, 73rd Edition
 The HPLC Solvent Guide, 2nd Edition, Paul C Sadek
 HPLC Columns, Theory, Technology & Practice, Uwe D Neue
 High-Performance Liquid Chromatography, 5th Edition, Veronica R Meyer
 The Merck Index, 12th Edition



As a recognized industry leader in the development, manufacturing and commercialization of innovative silica gel products, and with multi-ton manufacturing capability, SiliCycle® is your partner of choice for all your **metal removal**, **catalysis**, **synthesis**, and **purification** requirements.

SiliaMetS®
Metal Scavengers



- SiliCycle is **THE** world leader and **THE** pioneer in metal scavenging
- Over 18 years of **know-how** in silica grafting and scavenging technology
- Strong, extensive, and confidential support and scientists to help you
- **Broadest** portfolio of scavengers (*wide variety of ligands*) and applications developed
 - Wide range of formats for all purification scales; from laboratory to plant scale purifications
- Cited in many external publications and patents by satisfied customers
 - Innovative **Metal Scavenging Screening Service** to develop and provide an on-hand solution
 - Cost efficient: low cost per gram of metal scavenged
 - Available in bulk quantities (*multi-ton scale*)

SiliaBond®
Functionalized Silica Gels



- SiliCycle is the leader in the development of functionalized silica gels; **SiliaBond** products
- **SiliaBond** are suitable for use in batch and in flow through reactions
- Broad application domains: **chromatographic** and **ion exchange** phases, selective **scavengers** for nucleophile, electrophile and genotoxic impurities and **supported reagents** for organic synthesis
 - **SiliaBond** offer several advantages compared to polymer-based equivalents: no swelling, no static charge, greater solvent compatibility, faster kinetics and easier scalability
 - Available in tonnage quantities

SiliaCat®
Heterogeneous Catalysts



- Compatible with a wide range of reactions like: **Suzuki, Stille, Heck, Sonogashira** and **Kumada** couplings, debenzoylation, hydrogenation, hydrosilylation and oxydation reactions
 - Silica-based catalysts for cleaner products and cleaner waste streams
 - Ideal for batch or flowthrough applications
- Reusable and leach proof
 - Cost efficient and highly stable
- Available from grams to multi-kilograms quantities

SiliaFlash®
Irregular Silica Gels



- SiliCycle; a Silica Expert recognized worldwide as a leader with an outstanding quality silica gel
- SiliCycle offers one of the largest selections of silicas; various pore and particle sizes available
- With **SiliaFlash**, no scale-up limitations; performance will remain the same with every particle size
 - **SiliaFlash's** Exceptional Characteristics:
 - Tight particle and pore size distributions
 - Very low level of fines
 - Neutral pH, low metal content, and controlled water level content
 - Lot-to-lot reproducibility

SiliaSep™
Flash Cartridges



- **SiliaSep** pre-packed flash cartridge improves purification efficiency by offering tight and homogeneous packing
- **SiliaSep** offers superior performances over competitive cartridges
- Two silica grades available; 40-63 microns (**SiliaSep**) & 15-40 microns (**SiliaSep** HP *High performance cartridges*)
 - Wide choice of cartridge sizes available; from 4 grams up to 1.6 kg
 - Universal luer fittings for compatibility with any flash systems
 - Wide variety of sorbents to address any separation need (*reversed & normal phases, ion exchangers and metal scavengers*)

SiliaPlate™
TLC Plates




- Made of **SiliaFlash** to optimize the conditions by using the same silica as in your purification
- Made with an extra hard layer to ensure no loss of silica
- Various types of plates available: TLC, HPTLC and preparative
- Different sizes and choice of backings
 - Different possible phases (*normal, reversed and more*)

SiliaPrep™ & SiliaPrepX™
Solid-Phase Extraction



- Complete portfolio of SPE phases (*normal, reversed, fluorinated, ion exchange & polymeric phases*)
- Wide choice of SPE cartridges & 96-well plates formats
 - **SiliaPrep** QuEChERS for pesticide multiresidues analysis
 - **SiliaPrep** Tips for micro solid-phase extraction

SiliaSphere™
Spherical Silica Gels



- Outstanding quality spherical silica gels for analytical and preparative chromatography
- One of the best spherical silica on the market with unbeatable performance
 - Broadest portfolio of chromatographic phases to suit all purification needs

SiliaChrom®
HPLC Columns



- Broad pH range (0.5 - 12.0)
 - High surface coverage
- Wide variety of chemistries available
 - Excellent column efficiency
 - Long column lifetime
 - Low bleed

SiliCycle
Consumables



- Broad syringe filter portfolio
- Wide variety of membranes available (*Nylon, PTFE, PVDF, PES & RC*)
 - Various Vials & Caps available for autosamplers and laboratory needs
 - Highest glass quality (USP Type 1)

Printed in Quebec, Canada, © SiliCycle Inc.

