

## Unparalleled performance of Advanced Electron Ionization GC-MS\_MS technology for the determination of nitrosamines in drinking water on TG-1701MS

#	Compound Name	CAS Number	Compound Class	RT (min)
1	NDMA-d6 (/Search?CompoundName=NDMA-d6)	17829-05-9 (/Search?CasNumber=17829-05-9)	Nitrosamine	4.6
2	NDMA (/Search?CompoundName=NDMA)	62-75-9 (/Search?CasNumber=62-75-9)	Nitrosamine	4.7
3	NMEA (/Search?CompoundName=NMEA)	10595-95-6 (/Search?CasNumber=10595-95-6)	Nitrosamine	5.3
4	NDEA-d10 (/Search?CompoundName=NDEA-d10)	1219794-54-3 (/Search?CasNumber=1219794-54-3)	Nitrosamine	5.7
5	NDEA (/Search?CompoundName=NDEA)	55-18-5 (/Search?CasNumber=55-18-5)	Nitrosamine	5.8
6	NDPA (/Search?CompoundName=NDPA)	621-64-7 (/Search?CasNumber=621-64-7)	Nitrosamine	7.0
7	NPYR (/Search?CompoundName=NPYR)	930-55-2 (/Search?CasNumber=930-55-2)	Nitrosamine	7.4
8	NPIP (/Search?CompoundName=NPIP)	100-75-4 (/Search?CasNumber=100-75-4)	Nitrosamine	7.6
9	NDBA (/Search?CompoundName=NDBA)	924-16-3 (/Search?CasNumber=924-16-3)	Nitrosamine	8.4

### Autosampler

<b>Manufacturer</b>	Thermo Fisher Scientific™
<b>Brand</b>	Thermo Scientific
<b>System</b>	Trace™
<b>Module Type</b>	Autosampler
<b>Instrument Module</b>	TriPlus RSH™
<b>Description</b>	The Thermo Scientific™ TriPlus RSH™ Autosampler expands automated capabilities beyond liquid, headspace, and SPME injections. With its advanced, built-in robotics, the system delivers exceptional precision, unprecedented sample handling flexibility, and reliable, weekend-long unattended operation.
<b>Catalog #</b>	1R77010-0600

### Mass Spectrometer

<b>Manufacturer</b>	Thermo Fisher Scientific™
<b>Brand</b>	Thermo Scientific™
<b>System</b>	TSQ™ 9000
<b>Module Type</b>	Mass Spectrometer
<b>Instrument Module</b>	TSQ™ 9000 AEI

## Gas Chromatograph

<b>Manufacturer</b>	Thermo Fisher Scientific™
<b>Brand</b>	Thermo Scientific™
<b>System</b>	TRACE 1310 GC™
<b>Module Type</b>	Gas Chromatograph
<b>Instrument Module</b>	TRACE 1310 GC™
<b>Description</b>	Gas Chromatograph Manufacturer Thermo Fisher Scientific™ Brand Thermo Scientific™ System TRACE 1310 GC™ Module Type Gas Chromatograph Instrument Module TRACE 1310 GC™ Description
<b>Catalog #</b>	14800302

## Columns

<b>Manufacturer</b>	Thermo Fisher Scientific™
<b>Brand</b>	Thermo Scientific™
<b>Model</b>	TG-1701MS
<b>Length</b>	30 m
<b>Inner Diameter</b>	0.25 mm
<b>Film Thickness</b>	0.5 µm
<b>Stationary Phase</b>	(14%-Cyanopropyl-phenyl)- methylpolysiloxane
<b>Catalog #</b>	26090-2230

### Auto-sampler parameters

Syringe volume (µL)	10
Sample volume (µL)	2.0
Plunger strokes	9
Air volume (µL)	1.0
Rinses	1
Post injection solvent washes	5

### GC Inlet Parameters

Carrier gas, Flow mode/ flow rate (mL/ min)	Helium, Constant flow / 1.3
Injection mode	Splitless with surge
Temperature (°C)	240
Split flow (mL/ min)	81
Splitless Time	1.0
Purge flow (mL/ min)	5
Surge pressure (psi)	25
Surge duration (min)	1.01
Vacuum compensation	On

### MS Conditions

Transfer line temperature (°C)	250
Ion source temperature (°C)	300
Source used	Thermo Scientific™ Advanced Electron Ionization (AEI)
Ionization Type	Electron Ionization (EI)
Acquisition Mode	Timed SRM
Tune Type	AEI SmartTune
Collision gas and pressure	Argon at 70 psi
Peak width	0.7 Da (both Q1 and Q3)

### Gradient GC Oven

Step#	Rate (C/min)	Temperature (C)	Hold Time (min)
1	0.0	35.0	1.00
2	25.0	130.0	0.00
3	20.0	250.0	2.00