

HIC

HYDROPHOBIC INTERACTION

CHROMATOGRAPHY

HIC PRODUCTS

- TSKgel Ether-5PW
- TSKgel Phenyl-5PW
- TSKgel Butyl-NPR

≡ TOSOH FACT

Tosoh Bioscience provides solutions for today's biological purification needs. In fact, some of the first commercial HIC products were manufactured by Tosoh. We take pride in our ability to design new products based on existing chemistries to solve specific customer applications.

We encourage you to have a confidential discussion with us about your specific needs. Whether it is a surface modification of an existing product or the creation of a new one, we encourage you to call on us to meet your needs for a customized solution.



INTRODUCTION TO TSKgel HIC COLUMNS

Hydrophobic Interaction Chromatography (HIC) is based on the interaction between hydrophobic groups on a protein and a hydrophobic ligand on the solid support. HIC offers a distinct advantage for easily denatured proteins; it can be run using moderate concentrations of ammonium sulfate, which favors the stability of many proteins.

The binding of proteins to a hydrophobic matrix is affected by a number of factors including (1) the type of ligand, (2) the ligand density on the solid support, (3) the backbone material of the matrix, (4) the hydrophobic nature of the protein, and (5) the type of salt used. All of these factors help to make HIC a powerful technique for the separation of biomolecules.


Tosoh Bioscience offers three different HIC column types in analytical format: TSKgel Phenyl-5PW, Ether-5PW and Butyl NPR. TSKgel Phenyl-5PW is also available in preparative column formats. See **FIGURE 1** for the structure of the HIC resins.

COLUMN SELECTION

The HIC packing materials are based on the polymeric TSKgel G5000PW size exclusion resin (a hydrophilic gel with an estimated protein exclusion limit of 5,000,000 Da) which is then derivatized with oligoethylene-glycol (Ether-5PW) or phenyl (Phenyl-5PW) groups. Columns, depending on diameter, are packed with 10, 13 or 20 μm particles.

TSKgel ETHER-5PW is less hydrophobic than TSKgel Phenyl-5PW. It displays weaker interaction and thus shorter retention times compared to Phenyl-5PW. TSKgel Ether-5PW is the best choice for the separation of very hydrophobic proteins such as membrane proteins or monoclonal antibodies.

FIGURE 1
Structure of TSKgel HIC resins



FEATURES

- Choice of three hydrophobic ligands (ether, phenyl or butyl)
- Rigid polymeric base resin
- Similar chemistry to TOYOPEARL resins
- TSKgel Phenyl-5PW offered in PEEK hardware
- Ether and Phenyl available in 2 mm ID format

The **TSKgel PHENYL-5PW** columns were the first commercially available, polymer-based columns for high performance HIC. These columns have been instrumental to the increase in popularity of this technique for analytical, preparative, and process scale separations of biopolymers. **FIGURE 2** compares the separation of standard proteins on the Ether, Phenyl, and Butyl supports under similar operating conditions.

The base material of TSKgel Butyl-NPR is of the same chemical composition as the G5000PW base material used to prepare Phenyl-5PW and Ether-5PW. The difference between the two packings is that the G5000PW packing is porous, whereas the base material of the TSKgel Butyl-NPR column consists of spherical 2.5 μm nonporous particles. Nonporous resins (NPR) are typically used for high-speed analytical applications.

TSKgel BUTYL-NPR is the least hydrophobic among the three TSKgel HIC columns and requires a higher salt concentration for binding. TSKgel Butyl-NPR columns provide fast and quantitative HIC, because smaller particles provide higher efficiency. By packing the 2.5 μm nonporous resin particles into shorter columns, typical analysis times are reduced to less than 10 minutes. Pore diffusion is often the rate-limiting step in the overall mass transport of large biomolecules through a porous column. Eliminating the pores provides higher resolution at higher flow rates. Another benefit of NPR resins is excellent mass recovery, allowing quantitation down to nanogram levels. These properties make TSKgel Butyl-NPR the preferred choice for process monitoring and quality control.

TSKgel HIC columns are compatible with water-soluble organic solvents at concentration below 50 % (20 % for Butyl-NPR).

TABLE I
Column selection for the TSKgel HIC columns

Sample	MW range (Da)	TSKgel Column
Peptides	< 10,000	Butyl-NPR
Medium to large proteins	> 10,000	Phenyl-5PW Ether-5PW Butyl-NPR
DNA, RNA, and PCR products	> 500,000	Phenyl-5PW Butyl-NPR
Oligonucleotides	> 10,000	Phenyl-5PW Butyl-NPR

BENEFITS

- Added flexibility during method development
- Wide pH range (2-12) enabling robust cleaning options
- Seamless scalability from analytical to preparative scale
- Eliminates undesirable interactions with column hardware
- LC-MS applications

HIC

SAMPLE CAPACITY

One definition of sample capacity is the amount of pure compound injected onto the column at which the peak width is 10% larger than the peak width under low loading conditions. Using this definition, the capacity of a 7.5 mm ID x 7.5 cm L TSKgel Phenyl-5PW column varies from 0.1 to 1 mg of protein. Resolution and peak width are dependent on sample loading, as shown in **FIGURE 3**. Therefore, sample loading should be kept within 0.1 - 0.5 mg in order to obtain the highest resolution.

Separations on TSKgel Ether-5PW columns usually take 30 - 60 minutes. 0.5 mg of pure protein can be purified from a 5 - 10 mg crude protein mixture using a 7.5 mm ID x 7.5 cm L column.

Since almost all of the surface area of a porous particle is inside the pores, the capacity of the 4.6 mm ID x 3.5 cm L TSKgel Butyl-NPR column is significantly less than that for the 7.5 mm ID x 7.5 cm L Phenyl-5PW column. Capacities for the Butyl-NPR column are 100 µg for crude sample and 2 µg for pure sample.

CHEMICAL STABILITY


TSKgel 5PW-type HIC columns are physically and chemically stable in water-soluble organic solvents (at < 50% methanol, ethanol, ACN, DMF, DMSO or < 30 % chloroform). Change the solvent gradually by reducing the flow rate (preferably with a gradient) because rapid change may cause degradation of column efficiency. Note: When changing to an organic solvent, reduce the salt concentration to prevent precipitation of the salt on the column. Also, chaotropic agents (urea, SDS, etc.) will reduce the adsorption of biomolecules; therefore, use low levels of these agents (<2 mol/L).

Polymer-based columns are stable when cleaning at alkaline pH. All TSKgel HIC columns can be routinely operated from pH 2-12. **Table II** shows that the phenyl groups on the TSKgel Phenyl-5PW are stable for more than 10 days upon exposure to 0.5 mol/L NaOH or 0.5 mol/L acetic acid.

TABLE II
Long-term exposure of TSKgel Phenyl-5PW to acid and base


Acid/base	Phenyl content (mmol/mL - resin)	
	Before exposure	After 10 days exposure
0.5 mol/L CH ₃ COOH	0.105	0.106
0.5 mol/L NaOH	0.105	0.104

FIGURE 2
Comparing conventional and nonporous HIC columns



Column: TSKgel Ether-5PW & TSKgel Phenyl-5PW, 7.5 mm ID x 7.5 cm L TSKgel Butyl-NPR, 4.6 mm ID x 3.5 cm L; Sample: 1. myoglobin, 2. ribonuclease A, 3. lysozyme, 4. α-chymotrypsin, 5. α-chymotrypsinogen; Injection: 5PW-type columns: 100 µL (50-100 µg), NPR-type column: 20 µL (1.5-40 µg); Elution: 60 min linear gradient from 1.8 mol/L to 0 mol/L (NH₄)₂SO₄ in 0.1 mol/L phosphate buffer, pH 7.0, for 5PW-type columns; 12 min linear gradient from 2.3mol/L to 0 mol/L (NH₄)₂SO₄ in 0.1 mol/L phosphate buffer, pH 7.0 for TSKgel Butyl-NPR; Flow Rate: 1.0mL/min; Detection: UV@280nm

FIGURE 3
Dependence of peak width on sample loading in the separation of proteins



Column: TSKgel Phenyl-5PW, 7.5 mm ID x 7.5 cm L; Sample: 1. myoglobin; 2. ribonuclease A; 3. ovalbumin; 4. α-chymotrypsin; concentration: 0.025 % to 1.6 %; Elution: 60 min linear gradient of (NH₄)₂SO₄ from 1.5 mol/L to 0 mol/L in 0.1 mol/L phosphate buffer (pH 7.0); Flow Rate: 0.5 mL/min; Temperature: 25 °C; Detection: UV@280 nm

APPLICATIONS - TSKgel HIC COLUMNS


MODULATION OF SELECTIVITY

The addition of organic solvents or chaotropic agents in the final buffer can improve separations. However, relative elution positions may change. Therefore, add chaotropic agent and organic solvent in small quantities. See **FIGURE 4** for the effect of chaotropic agents and organic solvents on the HIC separation of two different samples.

ANTIBIOTICS


The TSKgel Ether-5PW column was used to determine the relative purity of the antibiotic components C-1027 and C-1027-AG as shown in **FIGURE 5**. Antibiotic C-1027 is composed of a protein consisting of many hydrophobic and hydroxyamino acids with a non-protein chromophore. Antibiotic C-1027-AG is composed of the hydrophobic and hydroxyamino acids without the chromophore.

FIGURE 4 Effect of urea and isopropanol on the separation of commercial lipoxidase and a standard protein mixture



Column: TSKgel Phenyl-5PW, 7.5 mm ID x 7.5 cm L;
 Sample: A & B: commercial lipoxidase, C & D: protein mixture: 1. cytochrome C, 2. myoglobin, 3. ribonuclease A, 4. lysozyme, 5. α -chymotrypsinogen, 6. α -chymotrypsin; Elution: A: 60 min linear gradient from 0.1 mol/L phosphate buffer containing 1.5 mol/L (NH₄)₂SO₄ (pH 7.0) to 0.1 mol/L phosphate buffer (pH 7.0), B: 60 min linear gradient from 0.1 mol/L phosphate buffer containing 1.5 mol/L (NH₄)₂SO₄ (pH 7.0) to 0.1 mol/L phosphate buffer containing 2 mol/L urea (pH 7.0), C: 60 min linear gradient from 0.1 mol/L phosphate buffer containing 1.8 mol/L (NH₄)₂SO₄ (pH 7.0) to 0.1 mol/L phosphate buffer (pH 7.0), D: 60 min linear gradient from 0.1 mol/L phosphate buffer containing 1.8 mol/L (NH₄)₂SO₄ (pH 7.0) to 0.1 mol/L phosphate buffer (pH 7.0) containing 7% isopropanol; Flow Rate : A & B: 0.5 mL/min; C & D: 1.0 mL/min;
 Temp.: 25°C; Detection: UV@280nm

FIGURE 5 Purification of anti-tumor antibiotic



Column: TSKgel Ether-5PW, 7.5 mm ID x 7.5 cm L; Sample: C-1027, C-1027-AG concentration: 1 mg/mL; Injection: 20 μ L; Elution: linear gradient from 1.5 mol/L to 0 mol/L (NH₄)₂SO₄ in 0.1 mol/L phosphate buffer, pH 7.0; Flow Rate: 0.8 mL/min; Detection: UV@220nm


HIC

APPLICATIONS - TSKgel HIC COLUMNS

RNAs


FIGURE 6 illustrates the separation of 16S and 23S ribosomal RNA on a TSKgel Phenyl-5PW column. The approximate molecular weights of these RNAs are 560,000 and 1,100,000 Da, respectively.

FIGURE 6
Retain large RNAs on TSKgel Phenyl-5PW



Column: TSKgel Phenyl-5PW, 7.5 mm ID x 7.5 cm L;
Sample: 16S and 23S rRNA from E. coli, 0.05 mg in 0.1 mL; Elution: 0 min linear gradient from 2 mol/L to 0 mol/L $(\text{NH}_4)_2\text{SO}_4$ in 0.1 mol/L phosphate buffer, pH 7.0; Flow Rate: 60.5 mL/min; Detection: UV@280nm

FIGURE 7
Scale up to preparative separations



Column: TSKgel Phenyl-5PW, A.) 7.5 mm ID x 7.5 cm L and B.) 21.5 mm ID x 15 cm L; Sample: 1. myoglobin, 2. ribonuclease A, 3. lysozyme, 4. α -chymotrypsinogen, 5. α -chymotrypsin; Elution: 60 min linear gradient from 1.8 mol/L to 0 mol/L $(\text{NH}_4)_2\text{SO}_4$ in 0.1 mol/L phosphate buffer, pH 7.0; Flow Rate: 0.5 mL/min (7.5 mm ID) or 4 mL/min (21.5 mm ID); Detection: UV@280nm

PROTEINS

FIGURE 7 compares the resolution of standard proteins on analytical and preparative TSKgel Phenyl-5PW columns. Different flow rates compensated for the change in particle size and column dimensions. High resolution was obtained on both columns.

ANTIBODY FRAGMENTS

FIGURE 8 shows the separation of Fab and Fc fragments of an antibody on TSKgel Butyl-NPR. The appearance of additional Fc fragments is due to the oxidation of methionine residues by 0.10% t-butylhydroperoxide (tBHP). The numbers above the Fc peaks correspond to the number of oxidized residues in each fragment.


Visit our website:

www.tosohbioscience.com for additional applications, product specifications and literature.

Contact our Technical Service specialists to discuss your specific application: +49 (0)711 13257-57 or techsupport.tb@tosoh.com.

Please see next page for ordering information.

FIGURE 8
Separation of Fab and Fc fragments on TSKgel Butyl-NPR



Column: TSKgel Butyl-NPR, 4.6 mm ID x 3.5 cm L; Elution: Buffer A: 2 mol/L $(\text{NH}_4)_2\text{SO}_4$, 20 mmol/L Tris, pH 7, Buffer B: 20 mmol/L Tris, pH 7; Gradient: linear from 10 % B to 100 % B in 34 minutes; Flow rate: 1 mL/min; Temperature: 30°C



► ORDERING INFORMATION

Part #	Description	ID (mm)	Length (cm)	Particle size (μm)	Number theoretical	Flow rate (mL/min)		Maximum pressure drop (MPa)
						Range	Max plates	
TSKgel Glass columns								
14013	Ether-5PW Glass, 1000 Å	5.0	5.0	10.0	≥ 600	0.5 - 0.8	1.0	2.0
14014	Ether-5PW Glass, 1000 Å	8.0	7.5	10.0	$\geq 1,000$	0.5 - 1.0	1.2	2.0
13063	Phenyl-5PW Glass, 1000 Å	5.0	5.0	10.0	≥ 600	0.5 - 0.8	1.0	2.0
08804	Phenyl-5PW Glass, 1000 Å	8.0	7.5	10.0	$\geq 1,000$	0.5 - 1.0	1.2	2.0
18147	Butyl-TOYOPEARL PAK 650S,	8.0	7.5	30.0				
18148	Phenyl-TOYOPEARL PAK 650S,	8.0	7.5	30.0				
18149	Ether-TOYOPEARL PAK 650S,	8.0	7.5	30.0				

TSKgel Stainless Steel Columns

18760	Ether-5PW, 1000 Å	2.0	7.5	10.0	$\geq 1,000$	0.05 - 0.1	0.12	0.6
08641	Ether-5PW, 1000 Å	7.5	7.5	10.0	$\geq 1,000$	0.5 - 1.0	1.2	2.0
18759	Phenyl-5PW, 1000 Å	2.0	7.5	10.0	$\geq 1,000$	0.05 - 0.1	0.12	0.8
07573	Phenyl-5PW, 1000 Å	7.5	7.5	10.0	$\geq 1,000$	0.5 - 1.0	1.2	2.0
07656	Phenyl-5PW, 1000 Å	21.5	15.0	13.0	$\geq 3,000$	4.0 - 6.0	8.0	2.0
07938	Phenyl-5PW, 1000 Å	55.0	20.0	20.0	$\geq 1,500$	20.0 - 40.0	50.0	0.4
14947	Butyl-NPR, nonporous	4.6	3.5	2.5		0.5 - 1.0	1.2	20.0

TSKgel PEEK columns

20023	BioAssist Phenyl, 1000 Å	7.8	5	10.0	$\geq 1,000$	0.5 - 1.0	1.2	2.0
-------	--------------------------	-----	---	------	--------------	-----------	-----	-----

Guard column products

	ID (mm)	Length (cm)	Particle Size (μm)	
14025	Ether-5PW Guardgel Kit, Glass		20.0	For P/Ns 14013 and 14014
08643	Ether-5PW Guardgel Kit		20.0	For P/N 08641
07652	Phenyl-5PW Guardgel Kit		20.0	For P/N 07573
16095	Phenyl-5PW Prep Guardgel Kit		20.0	For P/N 07656
07936	Phenyl-5PW Guard column	45.0	5.0	20.0 For P/N 07938

HIC

**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de