



SEC COLUMNS

SEC
SIZE EX
CLUSION
CHROMATO
GRAPHY


TOSOH BIOSCIENCE

PRINCIPLES OF CHROMATOGRAPHY

The analysis, isolation, and purification of biomolecules can be accomplished by a number of chromatographic modes. Each mode is based on specific physical, chemical, or biological interactions between the sample biomolecule and packing material.

Tosoh Bioscience offers a comprehensive line of TSKgel and TOYOPEARL media and pre-packed TSKgel columns for all common modes of liquid chromatography.

The various modes of chromatography involve separations that are based on specific features of the target or sample, like size, charge, hydrophobicity, function or specific content of the molecule. To find out more about general principles of liquid chromatography and on how each of them works, visit us on our YouTube channel <https://youtu.be/E3z1wllmvHI>



SIZE EXCLUSION CHROMATOGRAPHY

ABOUT US

WITH A GLOBAL PERSPECTIVE.

Tosoh Bioscience is a leading manufacturer in the field of liquid chromatography. The portfolio of over 500 specialty products encompasses instruments for size exclusion/gel permeation chromatography and a comprehensive line of media and prepacked (U)HPLC columns for all common modes of liquid chromatography. Over the last 40 years, TSKgel SW columns have become the worldwide industry standard for size exclusion chromatography of biomolecules.

Tosoh manufacturing sites in Japan provide products to the sales and support subsidiaries in the U.S. and Europe, ensuring full global coverage. Our technical specialists in the European Headquarters provide assistance in developing HPLC applications or purification methods, in up-scaling, or packing process columns. We offer chromatographic workshops, on-site training, and are the sole sponsor of the HIC/IPC Bioseparation Conference series.



1
TOSOH BIOSCIENCE GMBH
IM LEUSCHNERPARK 4
64347 GRIESHEIM
GERMANY

T + 49 (0) 6155 70437 00
F + 49 (0) 6155 83579 00
INFO.TBG@TOSOH.COM
WWW.TOSOHBIOSCIENCE.DE

2
TOSOH BIOSCIENCE LLC
3604 HORIZON DRIVE,
SUITE 100
KING OF PRUSSIA, PA 19406, USA

T +1 484 805 1219
F +1 610 272 3028
INFO.TBL@TOSOH.COM
WWW.SEPARATIONS.US.TOSOHBIOSCIENCE.COM

3
TOSOH CORPORATION
3-8-2 SHIBA, MINATO-KU
TOKYO 105-8623
JAPAN

T +81 3 5427 5118
F +81 3 5427 5198
INFO@TOSOH.CO.JP
WWW.TOSOHBIOSCIENCE.COM




4
TOSOH BIOSCIENCE SHANGHAI CO. LTD.

ROOM 301, PLAZA B,
NO. 1289 YI SHAN ROAD
XU HUI DISTRICT
SHANGHAI, 200233, CHINA
T +86 21 3461 0856
F +86 21 3461 0858
INFO@TOSOH.COM.CN
WWW.SEPARATIONS.ASIA.TOSOHBIOSCIENCE.COM

5
TOSOH ASIA PTE. LTD.
63 MARKET STREET #10-03
BANK OF SINGAPORE CENTRE
SINGAPORE 048942, SINGAPORE

T +65 6226 5106
F +65 6226 5215
INFO.TSAS@TOSOH.COM
WWW.TOSOHASIA.COM



AUTHORIZED DISTRIBUTOR

MZ-Analysetechnik GmbH
Wöhlerstraße 2-6 • D-55120 Mainz
Tel +49 6131 68 66 19
Fax +49 6131 68 66 20
e-mail: info@mz-at.de
www.mz-at.de

TOSOH HISTORY

- 1935 TOYO SODA MANUFACTURING CO., LTD. IS FOUNDED
- 1936 OPERATION OF NANYO MANUFACTURING COMPLEX BEGINS
- 1971 SCIENTIFIC INSTRUMENTS DIVISION DEVELOPS FIRST GPC COLUMN USING TSKgel
- 1974 HPLC COLUMN PLANT IS COMPLETED
- 1979 TOSOH DEVELOPS TOYOPEARL MEDIA
- 1983 TOSOH DEVELOPS HYDROPHOBIC INTERACTION MEDIA
- 1987 TOSOH US STARTS OPERATING FROM MONTGOMERYVILLE
- 1989 TOSOH GMBH STARTS OPERATING FROM STUTTGART
- 1995 TOSOH NANYO GEL FACILITY RECEIVES ISO 9001
- 2000/2001 FORMER TOSOH US AND EUROPE OPERATIONS BECOME TOSOH BIOSEP, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2002/ 2003 ALL SCIENTIFIC & DIAGNOSTIC SYSTEM RELATED COMPANIES IN EUROPE AND THE US ARE UNIFIED UNDER THE NEW NAME TOSOH BIOSCIENCE
- 2008 EcoSEC, THE 7TH GENERATION GPC SYSTEM IS INTRODUCED
- 2009 TOSOH BIOSCIENCE GMBH CELEBRATES ITS 20TH ANNIVERSARY
- 2010 TOSOH CELEBRATES ITS 75TH YEAR IN BUSINESS AND CONTINUED RAPID EXPANSION IN CHINA
- 2011 TOSOH BIOSCIENCE CELEBRATES 40 YEARS OF OPERATION
- 2012 TOSOH RELEASES FIRST TOYOPEARL MIXED-MODE RESIN
- 2013 TOSOH RELEASES A HIGH CAPACITY PROTEIN A RESIN
- 2014 TOSOH BIOSCIENCE GMBH CELEBRATES ITS 25TH ANNIVERSARY
- 2015 TOSOH BIOSCIENCE GMBH MOVES TO GRIESHEIM, GERMANY



SEC SIZE EXCLUSION CHROMATOGRAPHY

Size exclusion chromatography (SEC) separates molecules based on their size, or more precisely, their hydrodynamic volume. It is usually applied to large molecules such as proteins or synthetic polymers. When an aqueous mobile phase is used, SEC is also referred to as gel filtration chromatography (GFC). When an organic eluent is applied, SEC is referred to as gel permeation chromatography (GPC). GPC is typically used to determine the molecular weight (MW) and the MW distribution of synthetic polymers while GFC is used to separate biopolymers based on their size.

Aqueous SEC is a popular technique for the separation and purification of proteins because of its effectiveness and non-denaturing mobile phase conditions. It is popular for the isolation of proteins, removal of aggregates, desalting or characterization of water-soluble polymers used in food products, paints, pharmaceutical formulations and the like. Stationary phases for aqueous SEC range from soft packing materials, such as dextran or agarose, over hydrophilic polymers to silica. Soft particles were employed as stationary phases for early GFC whereas today porous silica particles with high mechanical strength are applied for aqueous SEC in high performance liquid chromatography (HPLC).

Tosoh Bioscience offers a broad portfolio SEC columns packed with silica or polymer based porous beads. They are well suited for a wide range of applications in R&D, method development and quality control. TSKgel SW, SW_{XL}, SuperSW, UltraSW, and UP-SW are silica SEC phases with pore size distributions suited to protein separations. TSKgel SW-type packings feature low adsorption and well-defined pore size distribution. It is the leading SEC column series for (U)HPLC due to its excellent resolution.

Polymeric TSKgel PW and PW_{XL} columns are designed for GFC of water soluble organic polymers, polysaccharides, oligosaccharides, DNA and RNA. The TSKgel Alpha and SuperAW series, based on a unique hydrophilic, polyvinyl resin, is suited for SEC of water-soluble and polar organic-soluble polymers. TSKgel columns for gel permeation chromatography of organic soluble polymers are described in a separate brochure on TSKgel GPC columns.

Tosoh Corporation employs state-of-the-art manufacturing techniques that result in uniformly bonded packing materials with narrow pore size distributions and well-defined particle sizes to ensure high performance and efficiency.







SEC HOW IT WORKS

Size exclusion chromatography (SEC) is a method in which components of a mixture are separated according to their molecular size, based on the flow of the sample through a porous packing. In contrast to all other modes of liquid chromatography the prerequisite for SEC is that the analyte does not interact with the surface of the stationary phases. Differences in elution time are ideally based solely on the volume the analyte passes.

Large biomolecules that cannot penetrate the pores of the packing material elute first from the column. They are said to be excluded from the packing; they flow with the mobile phase in the interparticle space of the packed column. The exclusion limit characterizes the upper limit of molecular weight (or size), beyond which molecules will elute at the same retention volume called the exclusion or void volume of the column. Many SEC columns are referred to by their exclusion limit.

Smaller molecules can partially or completely enter the porous particles. Because these smaller molecules have to flow through the interparticle space, as well as through the pore volume, they will elute from the column after the excluded sample components. Molecules small enough to penetrate the whole pore system of the stationary phase will pass the entire pore and interparticle volume, and will elute late. Their retention volume is referred to as 'total permeation' in SEC, whereas it is interpreted as 'unretained peak' in conventional LC modes.

SEC is a very simple method for separating biomolecules, because it is not necessary to change the composition of the mobile phase during elution. However, the separation capacity of this method is limited. For a baseline separation it is necessary that the molecular weights of the molecules differ by at least 10 to 20 %.



SEC

TSKgel SEC COLUMNS



Tosoh Corporation has a proud history of innovation in size exclusion chromatography. TSKgel SEC columns are known worldwide for their reliability and suitability for the analysis of proteins, peptides and other biological macro-molecules. The complete TSKgel SW, PW, Alpha and SuperAW column lines consist of either silica based or polymer based packings, ranging in particle size from 2 μm to 20 μm . Columns are available in analytical through semi- preparative size, in stainless steel, PEEK or glass.

COLUMN SELECTION

The main criterion in choosing between the TSKgel SW, PW, Alpha and SuperAW SEC columns is the molecular weight of the sample and its solubility. The fact that the TSKgel SW columns are based on silica and the TSKgel PW, Alpha and SuperAW columns are derived from a hydrophilic polymer network has less impact on the separation than the particle and pore size differences.

TSKgel SW SERIES

Tosoh Bioscience TSKgel SW, SW_{XL}, SuperSW, UltraSW, and UP-SW series are silica SEC phases with pore size distributions suited to protein separations. A hydrophilic diol-type bonded phase shields the silica surface from interacting with protein samples. Due to their high resolving power, the TSKgel SW columns are suitable for the separation of mono-disperse biopolymers such as proteins and nucleic acids.

TSKgel SW-type packings feature low adsorption and well-defined pore size distribution. They are the leading SEC columns in bioanalysis due to their excellent resolution. The new TSKgel UP-SW columns offer the opportunity to transfer existing TSKgel SW series HPLC methods to UHPLC.

TSKgel PW SERIES

TSKgel PW and PW_{XL} columns are packed with hydrophilic, rigid polymethacrylate beads. They are commonly used for the separation of synthetic water soluble polymers because they exhibit a much larger separation range, better linearity of calibration curves, and less adsorption than the TSKgel SW columns. While a TSKgel SW column is typically the first column to try for biopolymers, TSKgel PW columns have demonstrated good results for smaller peptides (<1,000 Da), protein aggregates, DNA fragments, and viruses. TSKgel PW_{XL}-CP columns are especially suited for the separation of cationic polymers at low salt.

TSKgel AW/ALPHA SERIES

The TSKgel Alpha series columns offer a new alternative for performing SEC. Their compatibility with a wide range of solvents makes them useful for both GFC and GPC. TSKgel SuperAW columns are based on the same chemistry as Alpha columns but have smaller particle sizes and shorter, narrower column dimensions for high throughput applications.

TABLE 1


CHARACTERISTICS OF TSKgel SIZE EXCLUSION COLUMN LINES

Column line	TSKgel SW / SW _{XL} / SuperSW/UltraSW/UP-SW	TSKgel PW / PW _{XL}	TSKgel Alpha / SuperAW
Resin type	Silica	Polymethacrylate	highly crosslinked Polymethacrylate
No. of available pore sizes	3/3/2/1/1	7	5
PH stability	2.5 - 7.5	2.0 - 12.0	2.0 - 12.0
Solvent compatability	100% polar	50% polar	100% polar, and nonpolar
Max. temp.	30°C	80°C*	80°C
Pressure**(MPa)	1.0 - 34.0	1.0 - 4.0	2.0 - 4.0
Application focus	Proteins	Water-soluble polymers	Intermediate polar polymers

* Except for the TSKgel G-DNA-PW, which can be operated up to 50°C. When operating below 10°C, reduce the flow rate to ensure that the maximum pressure is not exceeded.

** Depends on column dimensions and particle size

Note: The operating conditions and specifications for each column are listed on the Operating Conditions and Specifications sheet (OCS) shipped with the column.



SEC TSKgel SEC COLUMN SELECTION

SAMPLE		COLUMN SELECTION		SELECTION CRITERIA	
		FIRST CHOICE	ALTERNATIVE		
Carbohydrates	polysaccharides	TSKgel GMPW _{XL} TSKgel SuperMultiporePW	TSKgel G5000PW _{XL} & TSKgel G3000PW _{XL}	large pore size, small particles, linear calibration curve, high resolving power	
	oligosaccharides	TSKgel G-Oligo-PW TSKgel SuperOligoPW	TSKgel G2500PW _{XL}	small particles, high resolving power	
Nucleic acids	DNA fragments	large	TSKgel G-DNA-PW or TSKgel G5000PW _{XL}	large pore size, small particles, high resolving power	
		medium and small	TSKgel G4000SW _{XL} , TSKgel BioAssist G4SW _{XL} TSKgel SuperSW3000 or TSKgel G3000SW _{XL}	suitable pore sizes	
	RNA		TSKgel G4000SW _{XL} TSKgel SuperSW3000 or TSKgel G3000SW _{XL}	TSKgel BioAssist G3SW _{XL}	suitable pore sizes
		oligonucleotides	TSKgel G2500PW _{XL}		small pore size, ionic interaction
Proteins	small to medium sized proteins	TSKgel UP-SW3000 TSKgel SuperSW3000 TSKgel G3000SW _{XL} TSKgel G4000SW _{XL} TSKgel SuperSW2000 or TSKgel G2000SW _{XL}	TSKgel G3000/G4000PW _{XL} TSKgel BioAssist G3SW _{XL} TSKgel BioAssist G4SW _{XL}	small particles small to medium range pore sizes	
		antibodies	TSKgel SuperSW mAB HR/HTP TSKgel UP-SW3000 TSKgel UltraSW Aggregate		fragments/monomer & dimer higher aggregates
	large proteins	low density lipoprotein	TSKgel G6000PW _{XL} or TSKgel G5000PW _{XL}		large pore sizes
		gelatin	TSKgel GMPW _{XL} TSKgel SuperMultiporePW-M TSKgel G3000SW _{XL}	TSKgel G5000PW _{XL} & G3000PW _{XL}	large pore size, linear calibration curve
Peptides	large	TSKgel SuperSW3000 TSKgel G3000SW _{XL} TSKgel BioAssist G3SW _{XL} or TSKgel G2000SW _{XL}	TSKgel SuperSW2000 / TSKgel G3000PW _{XL} TSKgel BioAssist G2SW _{XL}	small to medium range pore size, versatile	
	small	TSKgel G2500PW _{XL}	TSKgel SuperSW2000 / TSKgel G2000SW _{XL}	linear calibration curve, high resolving power	
Viruses		TSKgel G6000PW _{XL} or TSKgel G5000PW _{XL} TSKgel SuperMultiporePW-H		large pore size, high resolving power	
Synthetic polymers		TSKgel GMPW _{XL} or TSKgel Alpha-M TSKgel SuperMultiporePW	TSKgel G5000PW _{XL} & G3000PW _{XL} / TSKgel Alpha- 5000 & Alpha-3000	large pore size, low adsorption, linear calibration curve	
	cationic	TSKgel G3000PW _{XL} -CP TSKgel G5000PW _{XL} -CP TSKgel G6000PW _{XL} -CP		medium to large pore size, low adsorption, linear calibration curve	
Synthetic oligomers	nonionic	TSKgel G-Oligo-PW TSKgel G2500PW _{XL} or TSKgel Alpha-2500 TSKgel SuperOligoPW and TSKgel SuperMultiporePW-N	TSKgel G2500PW / TSKgel SuperAW2500	small pore size, high resolving power	
	anionic	TSKgel G2500PW _{XL} or TSKgel Alpha-2500	TSKgel G2500PW / TSKgel SuperAW2500	small pore size, ionic interaction	

SEC

TSKgel SW SERIES



TSKgel SW-type columns (SW, SW_{XL}, SuperSW, UltraSW, and UP-SW) are all based on spherical silica particles with very high internal pore volume. They are stable from pH 2.0 to 7.5 and have excellent solvent stability up to 100% polar organic solvents.


Three different pore sizes of the SW and SW_{XL} packings result in different exclusion limits for several sample types, as shown by the calibration curve in Figure 1. From this data, recommended separation ranges for globular proteins can be made for each column (see Table 2).

Different particle sizes, column dimensions, and column hardware materials are available. The resulting differences in column characteristics allow the scientist to select the appropriate column to his individual separation requirements. Latest innovations in SW technology comprise dedicated small particle columns for antibody analysis (Super SW mAb and UltraSW Aggregate) and UHPLC columns with 2 μm particle and 25 nm pore size (UP-SW3000).

HIGHLIGHTS

- Rigid spherical silica gel chemistry bonded with hydrophilic groups
- Well defined pore size distribution
- Low non specific adsorption
- SuperSW and UltraSW designed for IgG analysis
- UHPLC columns with 2 μm particles available

➤ **FIGURE 1**



PROTEIN CALIBRATION CURVES FOR TSKgel SW_{XL} COLUMNS

Column: TSKgel SW_{XL} columns, 5 or 8 μm, 7.8 mm ID x 30 cm L
 Sample: 1. thyroglobulin (660,000 Da); 2. IgG (160,000 Da); 3. BSA (67,000 Da); 4. ovalbumin (43,000 Da); 5. peroxidase (40,200 Da); 6. β-lactoglobulin (18,400 Da); 7. myoglobin (16,900 Da); 8. ribonuclease A (12,600 Da); 9. cytochrome C (12,400 Da); 10. glycine tetramer (246 Da)
 Mobile phase: 0.3 mol/L NaCl in 0.1 mol/L sodium phosphate buffer, pH 7.0; Detection: UV @ 220 nm



SEC TSKgel SW SERIES

TABLE 2

PROPERTIES AND SEPARATION RANGES FOR TSKgel SW TYPE PACKINGS

TSKgel COLUMN	ID (mm) X LENGTH (cm L)	PARTICLE SIZE (μm)	PORE SIZE (nm)	MIN. NO. THEORET. PLATES	MOLECULAR WEIGHT OF PROTEINS (Da)
SuperSW2000	4.6 x 30	4	12.5	30,000	5 x 10 ³ –1.5 x 10 ⁵
G2000SW _{XL}	7.8 x 30	5	12.5	20,000	5 x 10 ³ –1.5 x 10 ⁵
BioAssist G2SW _{XL}	7.8 x 30	5	12.5	20,000	5 x 10 ³ –1.5 x 10 ⁵
QC-PAK GFC 200	7.8 x 15	5	12.5	10,000	5 x 10 ³ –1.5 x 10 ⁵
G2000SW	7.5 x 30/60 21.5 x 30/60	10 13	12.5 12.5	10,000/20,000 10,000/20,000	5 x 10 ³ –1.5 x 10 ⁵ 5 x 10 ³ –1.5 x 10 ⁵
SuperSW3000	4.6 x 30	4	25	30,000	1 x 10 ⁴ –5 x 10 ⁵
UP-SW3000	4.6 x 15/30	2	25	25,000/45,000	1 x 10 ⁴ –5 x 10 ⁵
SuperSW mAb HTP	4.6 x 15	4	25	15,000	1 x 10 ⁴ –5 x 10 ⁵
SuperSW mAb HR	7.8 x 30	4	25	30,000	1 x 10 ⁴ –5 x 10 ⁵
G3000SW _{XL}	7.8 x 30	5	25	20,000	1 x 10 ⁴ –5 x 10 ⁵
BioAssist G3SW _{XL}	7.8 x 30	5	25	20,000	1 x 10 ⁴ –5 x 10 ⁵
QC-PAK GFC 300	7.8 x 15	5	25	10,000	1 x 10 ⁴ –5 x 10 ⁵
G3000SW	7.5 x 30/60	10	25	10,000/20,000	1 x 10 ⁴ –5 x 10 ⁵
UltraSW mAb Aggregate	7.8 x 30	3	30	35,000	1 x 10 ⁴ –2 x 10 ⁶
G4000SW _{XL}	7.8 x 30	8	45	16,000	2 x 10 ⁴ –7 x 10 ⁶
BioAssist G4SW _{XL}	7.8 x 30	8	45	16,000	2 x 10 ⁴ –7 x 10 ⁶
G4000SW	7.5 x 30/60 21.5 x 30/60	13 17	45 45	8,000/16,000 8,000/16,000	2 x 10 ⁴ –7 x 10 ⁶ 2 x 10 ⁴ –7 x 10 ⁶



SEC TSKgel SW SERIES APPLICATIONS

AGGREGATE ANALYSIS


Protein aggregation is a common issue encountered during expression, purification and formulation of protein biotherapeutics, which needs to be characterized and controlled during the development and production of protein pharmaceuticals such as monoclonal antibodies (mAbs). Even small amounts of aggregates can alter the therapeutic function. TSKgel G3000_{XL} columns are the industry standard for quality control of MABs by SEC. Besides the traditional detection of proteins using their UV absorption at 280 nm, multi angle light scattering (MALS) detection gains more and more interest in protein analysis. Being a universal detection method, MALS can deliver valuable additional information.

As it will also detect several other impurities, pure solvents and samples are of utmost importance. This also applies to the stationary phase, which should not generate interfering baseline noise under the conditions used for analysis. Figure 2 shows the analysis of MAB aggregates of a commercial monoclonal antibody with UV, refractive index (RI) and MALS detection. Separation was performed on a TSKgel G3000_{SWXL} column under standard conditions. With modern (U)HPLC instrumentation resolution can be further enhanced by using the latest generation of SW columns, TSKgel UltraSW Aggregate or TSKgel UP-SW3000 (see pages 10 ff)

When the analysis of proteins needs to be performed in a metal free environment, the BioAssist_{SW} series offers TSKgel _{SWXL} packings in PEEK housings, featuring the same performance as stainless steel columns. Figure 3 shows a typical separation performed with a BioAssist _{SW} PEEK column.

USE OF DETERGENTS

FIGURE 2




SEC-MALS-UV-RI ANALYSIS OF MAB AGGREGATES

Column: TSKgel G3000_{SWXL} column, 5 μm, 7.8 mm ID x 30 cm L, Sample: monoclonal antibody, Inj. volume: 20 μL, Mobile phase: phosphate buffered saline (PBS); Flow rate: 1 mL/min, Detection: MALS (red), refractive index (blue) & UV @ 280 nm (green) HPLC System: LC-20A prominence, Shimadzu, MALS detector: miniDAWN™ TREOS, Wyatt Techn. Corp.

Some SEC separations require denaturing conditions like sodiumdodecylsulfate (SDS) containing eluents. In other cases the formulations of biopharmaceuticals contain some detergents (e.g. Tween 20 or Triton). TSKgel _{SW} type columns can be operated under these conditions although certain amounts of the detergent will stick to the column, affecting column lifetime and the future use of the column. If analysis under denaturing conditions was performed once, the affected column should be used with detergent containing eluents only. Regular maintenance of the column, the use of guard columns and monitoring of the column status by analyzing control samples are recommended as well.

FIGURE 3



QC ANALYSIS OF AN ANTI-TSH ANTIBODY PURIFIED FROM CELL CULTURE SUPERNATANT

Column: TSKgel BioAssist G3_{SWXL}, 5 μm, 7.8 mm ID x 30 cm L, Mobile phase: 0.3 mol/L phosphate buffer, pH 7.0, Flow rate: 1.0 mL/min; Inj. volume: 50 μL



SEC TSKgel SuperSW & SuperSW mAb SERIES

Speed and resolution is an increasing demand in liquid chromatography. The need for high sensitivity applicable to trace analysis is increasing as sample size or sample concentrations become limited. To meet the needs of high sensitivity and high resolution protein analysis Tosoh Bioscience developed TSKgel SuperSW columns packed with 4 µm spherical silica particles. TSKgel SuperSW columns are available in two pore sizes, 125 Å and 250 Å, both featuring a minimum of 30,000 theoretical plates / column. Compared to the well established TSKgel SW_{XL} (5 µm) series, SuperSW columns show higher resolution due to a 50 percent increase in theoretical plate numbers (Table 3).

To further improve performance, TSKgel SuperSW media are packed into columns with smaller inner diameter (1.0, 2.0, 4.6 mm ID). The smaller diameters are one reason for increased peak heights. In addition, the high resolution of the 4 µm particles and accordingly smaller peak widths further increase peak height provided the HPLC system is optimized with regard to extra column dead volume.

Figure 4 demonstrates the superior sensitivity reached with TSKgel SuperSW2000 compared to a TSKgel G2000SW_{XL} column of the same length but larger inner diameter. TSKgel SuperSW can yield peak heights approximately 4 times that of TSKgel SW_{XL} due to downsizing in column diameter and increased theoretical plates.

HIGHLIGHTS

- SuperSW mAb columns tailored to antibody analysis
- High throughput and high resolution columns available
- 4 µm particle size featuring superior resolution and highest sensitivity
- Low non-specific adsorption
- High reproducibility due to well-defined pore size distribution
- 30,000 theoretical plates / column (30 cm L)


A new series of TSKgel SuperSW mAb columns was developed to improve IgG analysis in QC or research applications. They feature the SuperSW3000 stationary phase packed in dimensions that are either ideal for either reaching maximum resolution or speeding up analysis without compromising resolution much.

➤ **TABLE 3**

SPECIFICATIONS OF TSKgel SuperSW SERIES COMPARED TO TSKgel SW_{XL} SERIES

TSKgel COLUMN	PARTICLE SIZE (µm)	COLUMN SIZE (mm ID X cm L)	GUARANTEED THEOR. PLATES
TSKgel SuperSW2000	4	4.6 x 30	30,000
TSKgel SuperSW3000	4	4.6 x 30	30,000
TSKgel G2000SW _{XL}	5	7.8 x 30	20,000
TSKgel G3000SW _{XL}	5	7.8 x 30	20,000

➤ **FIGURE 4**



COMPARISON OF TSKgel SuperSW2000 AND TSKgel G2000SW_{XL} FOR THE SEPARATION OF PROTEINS

Column: A. TSKgel G2000SW_{XL}, 7.8 mm ID x 30 cm L; B. TSKgel SuperSW2000, 4.6 mm ID x 30 cm L
 Sample: 1. thyroglobulin (0.2 mg/mL); 2. albumin (1.0 mg/mL); 3. ribonuclease A (1.0 mg/mL); 4. p-aminobenzoic acid (0.01 mg/mL)
 Inj. volume: 5 µL, Mobile phase: 0.1 mol/L phosphate buffer + 0.1 mol/L Na₂SO₄ + 0.05% NaN₃ (pH 6.7), Flow rate: 0.35 mL/min (SuperSW2000), 1.0 mL/min (G2000SW_{XL}), Temp: 25°C; Detection: UV @ 280 nm

SEC TSKgel SuperSW SERIES



SEPARATION RANGE OF TSKgel SuperSW

The TSKgel SuperSW series has the same pore sizes as the conventional TSKgel SW_{XL} series with equivalent grade. Therefore it has similar calibration curves and separation ranges as well. Method transfer from conventional SEC to high resolution SEC is very straight forward. TSKgel SuperSW columns are available in two pore sizes, 125 Å (TSKgel SuperSW2000) and 250 Å (TSKgel SuperSW3000). Figure 5 shows the SEC calibration curves for standard proteins. In general, TSKgel SuperSW2000 is suited to separate proteins with molecular weights of 150 KDa or smaller. TSKgel SuperSW3000 can be used for the separation of proteins with molecular weights up to 500 KDa.

INCREASED DETECTION LIMIT

Table 4 shows the detection limits for some proteins. The high sensitivity allows for analysis of nanogram sample amounts. If sample amount is limited a reduction of column inner diameter can further enhance sensitivity. TSKgel SuperSW3000 columns are available with 4.6; 2 and 1 mm ID. Figure 6 shows the levels of sensitivity which can be reached with semi-micro or micro columns. When limited sample amount is an issue (e.g. in proteomics research) enhancing detection limits by using a micro column can increase the number of hits.


TABLE 4

DETECTION LIMIT FOR PROTEINS (S/N=3)

	TSKgel SuperSW	TSKgel SW _{XL}
FLOW CELL	STANDARD CELL (LOW DEAD VOLUME TYPE)	STANDARD CELL (LOW DEAD VOLUME TYPE)
Light path length	10 mm	10 mm
Thyro-globulin	70 ng	200 ng
γ-globulin	50 ng	100 ng
Bovine serum albumin	70 ng	200 ng
Ovalbumin	50 ng	100 ng
Myoglobin	15 ng	30 ng

Column: TSKgel SuperSW3000, 4.6 mm ID x 30 cm L; TSKgel G3000SW_{XL}, 7.8 mm ID x 30 cm L, Mobile phase: 0.2 mol/L phosphate buffer, pH 6.7, Detection: UV @ 220 nm


FIGURE 5



PROTEIN CALIBRATION CURVES FOR TSKgel SuperSW

Column: TSKgel SuperSW Series, 4.6 mm ID X 30 cm L, Sample: Standard proteins (5 μL, 0.1 g/L each); 1.thyroglobulin 2. γ-globulin 3. bovine serum albumin, 4. β-lactoglobulin 5. lysozyme 6. cytochrome C 7. glycine tetramer, Mobile phase: 0.2 mol/L phosphate buffer (pH 6.7), Flow rate: 0.35 mL/min; Detection: UV @ 280 nm

FIGURE 6



ESTIMATION OF SENSITIVITY

Column: TSKgel SuperSW3000, 1.0, 2.0, 4.6 mm ID x 30 cm L, Sample: 1. thyroglobulin (1.0 g/L), 2. γ-globulin (2.0 g/L), 3. ovalbumin (2.0 g/L), 4. ribonuclease A (3.0 g/L), 5. p-aminobenzoic acid (0.02 g/L) Mobile phase: 0.1 mol/L phosphate buffer + 0.1 mol/L Na₂SO₄ + 0.05% NaN₃, Flow rate: 16 μL/min (1 mm), 65 μL/min (2 mm); 350 μL/min (4.6 mm), Inj. volume: 0.2 μL; Temperature: 25 °C, Detection: UV @ 280 nm, cell vol. 2 μL (4.6 mm ID), 35 nL (1.0, 2.0 mm ID)



SEC - TSKgel SuperSW mAb HTP/HR FOR mAb APPLICATIONS


QC ANALYSIS OF ANTIBODIES

Thermally induced denaturation or aggregation of therapeutic antibodies can be a significant problem during different stages of its production and formulation, since aggregates affect the efficiency of the biopharmaceutical. Thus the quantification of aggregates is an important parameter in the quality control analysis of biopharmaceuticals. Using TSKgel SuperSW mAb columns the amounts of tri-, di- and monomers of monoclonal antibodies can be monitored.

TSKgel SuperSW mAb columns are specifically designed to facilitate analysis of antibodies. Each of the columns is tailored to a specific separation problem. TSKgel SuperSW mAb HR - "HR" indicating high resolution - delivers superior resolution over the whole range of typical mAb SEC analysis, from fragments to aggregates. It has the same column dimensions as the established TSKgel G3000SW_{XL} column.

Figure 7 shows a comparison of the analysis of a mouse-human chimeric IgG using these two columns. Resolution is improved with TSKgel SuperSW mAb HR.

➤ **FIGURE 7**




COMPARISON OF AGGREGATE ANALYSIS

Columns: A. TSKgel G3000SW_{XL}, B. TSKgel SuperSW mAb HR,
Dimension: 7.8 mm ID × 30 cm;
Eluent: 0.2 mol/L phosphate buffer (pH 6.7) + 0.05% Na₃
Flow rate: 0.8 mL/min; Detection: UV @ 280 nm; Temp.: 25°C;
Sample: monoclonal antibody, (mouse-human chimeric IgG, Erbitux), 10 µL

FAST ANALYSIS OF AGGREGATION

TSKgel SuperSW mAb HTP - "HTP" indicating high throughput - was developed to enable an easy transfer of established methods to fast (U)HPLC analysis. The SuperSW3000 stationary phase is packed into column hardware with smaller (4.6 mm) inner diameter. This column enables doubling the throughput without compromising resolution too much. Figure 8 shows the optimization of the separation of mAb Aggregates with regard to analysis time. Compared to the standard separation with maximum resolution using TSKgel SuperSW mAb HR (see Figure 7) the analysis time can be reduced by factor 3 when using TSKgel SuperSW mAb HTP and adjusting the flow rate accordingly. Resolution between IgG monomer and dimer is still sufficient for quantification.

➤ **FIGURE 8**



FAST ANALYSIS OF mAb AGGREGATION

Column: TSKgel SuperSW mAb HTP (4.6 mm ID × 15 cm)
Elution: 0.2 mol/L phosphate buffer (pH 6.7) + 0.05% Na₃
Flow rate: 0.50 mL/min, 0.35 mL/min; Detection: UV @ 280 nm
Temp.: 25°C; Sample: monoclonal antibody (mouse-human chimeric IgG, Erbitux), 5 µL

SEC TSKgel UltraSW COLUMN



TSKgel UltraSW AGGREGATE FEATURES

TSKgel UltraSW Aggregate provides a smaller particle size of 3 micron and a higher exclusion limit through slightly larger pores (30 nm) than TSKgel SuperSW mAb columns. This column was developed to offer a wider separation window in the molecular mass range of antibody aggregates.

Figure 9 shows the calibration curve of TSKgel UltraSW Aggregate compared to the calibration curves of TSKgel SuperSW mAb and G3000 SW_{XL}. The calibration curve of TSKgel SuperSW mAb HR is most similar to the one of TSKgel G3000SW_{XL}, the current industrial standard for antibody analysis. The curve for TSKgel UltraSW Aggregate shows a shallower slope in the high molecular weight region (MW > 300 000).


HIGHLIGHTS

- High resolution through 3 μm particle size
- Large pore size expands separation range for higher aggregates
- Dimensions optimized for mAb aggregate analysis

ANALYSIS OF mAb AGGREGATES

The wide separation window of TSKgel UltraSW Aggregate at higher molecular weights combined with the small particle size results in a superior separation of higher antibody aggregates. Resolution of dimer, trimer, and higher aggregates is higher than for TSKgel SuperSW mAb HR and TSKgel G3000 SW_{XL} (Figure10).


➤ **FIGURE 9**



CALIBRATION CURVES


Columns: A. TSKgel G3000SW_{XL}, B. TSKgel SuperSW mAb HR, C. TSKgel UltraSW Aggregate; Dimension: 7.8 mm ID x 30 cm; Eluent: 0.2 mol/L phosphate buffer (pH 6.7) + 0.05% NaN₃; Flow rate: 0.8 mL/min; Detection: UV @ 280 nm; Temp.: 25°C; Sample: monoclonal antibody, (mouse-human chimeric IgG, Erbitux), 10 μL

➤ **FIGURE 10**



COMPARISON OF AGGREGATE ANALYSIS

Columns: TSKgel UltraSW Aggregate (7.8 mm ID x 30 cm), TSKgel SuperSW mAb HR(7.8 mm ID x 30 cm), TSKgel G3000SW_{XL} (7.8 mm ID x 30 cm); Eluent: 0.2 mol/L phosphate buffer (pH 6.7) + 0.05% NaN₃; Flow rate: 1.0 mL/min; Temperature: 25°C; Detection: UV @ 280 nm; Injection vol.: 20 μL; Sample: IgG (human monoclonal, 2.5 g/L), 1 monomer, 2 dimer, 3 trimer, 4 tetramer, 5 multimers



SEC TSKgel UP-SW3000 SERIES


TSKgel UP-SW3000 FEATURES

TSKgel UP-SW3000 columns packed with 2 μm silica based particles are the latest addition to the popular TSKgel SW series, the gold standard for QC analysis of antibody therapeutics. The new silica-based UHPLC columns are based on the proven proprietary surface technology of the renowned TSKgel SW series and were engineered to facilitate the transfer of existing HPLC methods to UHPLC systems.

The new SEC columns can be used with modern, low extra column dead volume HPLC and UHPLC systems and are available in the lengths of 15 or 30 cm. The short one enables short analysis times; the long one provides higher resolution for mAb analysis. The lifetime of the columns can be improved when using the corresponding guard columns. A "direct connect" (DC) guard column is available to minimize extra column dead volume.

Figure 11 shows the calibration curve and the molecular weight range of the new 2 μm TSKgel UP-SW3000 compared to the 5 micron TSKgel G3000SW_{XL} and 4 micron TSKgel SuperSW3000. The calibration curves and mass ranges are almost identical which facilitates method transfer from existing methods.

➤ **FIGURE 11**



COMPARISON OF CALIBRATION CURVES

Columns : A: TSKgel UP-SW3000 (4.6 mm ID \times 30 cm, red)
 B: TSKgel SuperSW3000 (4.6 mm ID \times 30 cm, blue)
 C: TSKgel G3000SW_{XL} (7.8 mm ID \times 30 cm, grey)
 Eluent: 100 mmol/L phosphate buffer (pH 6.7) + 100 mmol/L sodium sulfate + 0.05% NaN₃; Flow rate: A & B: 0.35 mL/min; C: 1.0 mL/min
 Temperature: 25°C; Detection: UV @ 280 nm; Injection vol.: 10 μL
 Samples: 1. thyroglobulin (640,000 Da); (1' thyroglobulin aggregate);
 2. γ -globulin (155,000 Da); (2' γ -globulin dimer); 3. ovalbumin (47,000 Da);
 4. ribonuclease A (13,700 Da); 5. p-amino benzoic acid (137 Da)

HIGHLIGHTS


- Highest resolution through 2 μm particle size
- Easy transfer of existing methods
- Optimized for mAb quality control
- Short runtimes with 15 cm column

UHPLC ANALYSIS OF ANTIBODIES

The separation of an antibody sample on the new 2 mikron UP-SW3000 column compared to the competitor UHPLC column is depicted in figure 12. The two columns differ slightly in pore size and particle size (25 nm for 2 μm TSKgel versus 20 nm for the 1.7 μm material). This difference in pore sizes results in a better separation in the molecular weight range of antibodies, fragments, and aggregates on TSKgel UP-SW.

Based on the wider separation window the resolution between monomer and dimer as well as dimer and trimer is slightly higher with TSKgel UP-SW3000 although particle size is slightly larger than in the competitor column. Moreover, also the fragment peak is more clearly separated from the monomer peak.

➤ **FIGURE 12**



COMPARISON OF mAb ANALYSIS

Column	RS (peak 1/2)	RS (peak 2/3)
TSKgel UP-SW3000 2 μm	1.52	3.56
Competitor UHPLC-SEC 1.7 μm	1.25	3.47

SEC - TSKgel SuperSW, UltraSW AND UP-SW SYSTEM REQUIREMENTS



OPTIMIZATION OF HPLC EQUIPMENT

To benefit from the improved features of TSKgel SuperSW, UltraSW and UP-SW columns the HPLC system should be optimized and extra column peak broadening reduced. This means reduction of dead volume and adjustment of sample concentration and injection volume.

SYSTEM DEAD VOLUME

Key components of the HPLC system with regard to dead volume reduction are the void volume of tubings, the cell volume of the detector cell and the void volume of the injection unit. Modern UHPLC systems designed for use with sub 2 μm particles exhibit extremely small dead volumes and can be used for SEC analysis without modification.

VOID VOLUME OF THE TUBING

The volume of tubing from injector to column, column to detector influences the diffusion within the tubing and the column efficiency. Column efficiency starts deteriorating remarkably when the volume of the tubing exceeds 10 μL (e.g. 0.1 mm ID x 150 cm L). Shortening of tubings of 0.1 or 0.125 mm inner diameter is often better than using longer capillaries with smaller inner diameters. The backpressure increases with smaller inner diameters and the system becomes more susceptible towards clogging.

DETECTOR CELL VOLUME

The detector cell volume also contributes to the dead volume of the system and might impair peak resolution. For most separations with 4.6 mm ID TSKgel SuperSW columns a 8-10 μL standard detector cell might be sufficient but for semi-micro (2 mm ID) or micro columns (1 mm ID), and for TSKgel UP-SW we strongly recommend using semi-micro/micro detector cells.

INJECTOR

The maximum number of theoretical plates in isocratic separations can be reached when using a low diffusion type manual injector like the Rheodyne 8125. All kinds of automated HPLC injectors will deteriorate column efficiency to a certain extend but due to practical reasons, auto-samplers are nowadays standard. All the more it is important to select an auto-sampler capable of trace injection mode. Dead volume of the outlet capillary should be minimized to the utmost (as short as possible, 0.1 mm ID). Figure 13 shows the effect of injector tubings on column efficiency for a 1 mm ID column.

TSKgel SuperSW, UltraSW, AND UP-SW OPERATING CONDITIONS

For best results, it is recommended to use the following experimental conditions for TSKgel SuperSW, UltraSW, and UP-SW columns:

CONNECTIONS	The conventional 0.1 mm tubing may be used, but length should be kept as short as possible. Void volume between the column and detector cell should be less than 20 μL.
INJECTOR	Best results are obtained with a low diffusion type manual injector (Rheodyne 8152). Autosampler outlet void volume should be as low as possible.
SAMPLE VOLUME	Sample volume should be 10 μL or less. Sample load should be less than 100 μg (4.6 mm ID column).
GUARD COLUMN	A guard column or an inline filter is highly recommended to reduce clogging and contamination.
DETECTOR	
Flow Cell	For best results, use a flow cell with a maximum of 2 μL. The 2 μL flow cell will give the highest efficiencies. A 2-10 μL flow cell can be used for 4.6 mm ID columns. However, theoretical plates will be reduced.
Time Constant	A small time constant (less than 0.5 sec) is needed to achieve best column performance.
PUMP	A pump capable of accurately delivering a flow rate between 0.01 mL/min and 0.35 mL/min is recommended.



SEC - TSKgel SuperSW, UltraSW AND UP-SW SYSTEM REQUIREMENTS

SAMPLE LOAD AND INJECTION VOLUME

Although the efficiency of TSKgel Super/Ultra/UP-SW columns is high, it is obvious that it decreases at high sample loads. Figure 14 shows that sample load should not exceed 100 µg for a TSKgel SuperSW3000 column of 4.6 mm ID x 30 cm L. On the other hand the injection volume itself is a critical parameter. As for all HPLC applications injection volume should be as small as possible. If injection volume exceeds 20 µL on a 4.6 mm ID column, a considerable deterioration of column efficiency is observed for TSKgel SuperSW2000 (80 µL for TSKgel SuperSW3000). In general the sample load should be less than 100 µg in less than 10 µL injection volume for a 4.6 mm ID TSKgel SuperSW column.

FLOW RATE DEPENDENCE

The effect of flow rate on column efficiency depends on particle size of packing materials, sample molecular size, eluent viscosity, etc. The appropriate flow rate for TSKgel SuperSW columns is up to 0.4 mL/min for a 4.6 mm ID column, up to 75 µL/min for a 2 mm ID column, and up to 20 µL/min for a 1 mm ID column, respectively. If higher resolution is required the flow rate can be lowered.

MOBILE PHASE

The eluent plays an important role in SEC separations. When denaturing agents are used, the exclusion limits for proteins become smaller since they lose their compact globular structure. Proper selection of eluting conditions is necessary to maximize the molecular sieving mechanism and to minimize secondary effects, such as ionic and hydrophobic interactions between the sample and the column packing material. In general, the use of relatively high ionic strength buffers is recommended for most protein applications. A neutral salt is often added to increase ionic strength.

RECOVERY OF PROTEIN

TSKgel SuperSW series is capable of obtaining high protein recovery even in trace analysis with sample load of 1 µg or lower. Most proteins are recovered quantitatively with TSKgel SuperSW series, but it is important to make sure that samples in small concentrations are not adsorbed to the sample vial or to the HPLC system itself. Similar samples should be injected several times before measurement so that adsorption points within the system are inactivated in advance when trace analysis is performed.

FIGURE 13





FIGURE 14



INFLUENCE OF TUBING (INJECTOR TO COLUMN)


Column: TSKgel SuperSW3000 1.0 mm ID x 30 cm L, Mobile phase: 0.1 mol/L phosphate buffer + 0.1 mol/L Na₂SO₄ + 0.05 % NaN₃
 Flow rate: 16 µL/min; Inj. volume: 0.2 µL; Temp.: 25 °C, Detection: UV @ 280 nm, Sample: p-Aminobenzoic acid (20mg/L),
 Tubing: ID (mm) x L (cm), Vol.
 0.050 x 20, 393 nL; 0.050 x 40, 785 nL; 0.050 x 60, 1178 nL;
 0.075 x 20, 883 nL; 0.075 x 40, 1766 nL; 0.075 x 60, 2469 nL;
 0.130 x 20, 2653 nL; 0.130 x 40, 5307 nL; 0.130 x 60, 7960 nL

EFFECT OF SAMPLE LOAD

Column: TSKgel SuperSW series, 4.6 mm ID x 30 cm L; TSKgel SW_{XL} series, 7.8 mm ID x 30 cm L, Sample: Bovine serum albumin, Mobile phase: 0.2 mol/L phosphate buffer, pH 6.7, Flow rate: 0.35 mL/min (SuperSW series), 1.00 mL/min (SW_{XL} series), Temp.: 25 °C; Detection: UV @ 280 nm, micro flow cell

SEC

TSKgel SW SERIES ORDERING INFORMATION



ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (μm)	NUMBER THEORETICAL PLATES	FLOW RATE (mL/min)		MAXIMUM PRESSURE DROP (MPa)
						RANGE	MAX.	
STAINLESS STEEL COLUMNS								
0023448	UP-SW3000 -NEW-	4.6	30	2	≥ 45000	0.10-0.35		34
0023449	UP-SW3000 -NEW-	4.6	15	2	≥ 25000	0.10-0.35		25
0022856	UltraSW Aggregate -NEW-	7.8	30	4	≥ 35,000	0.5 - 1.0		12.0
0018674	SuperSW2000	4.6	30	4	≥ 30,000	0.1 - 0.35	0.4	12.0
0021845	SuperSW3000	1.0	30	4	≥ 18,000	0.016	0.02	12.0
0021485	SuperSW3000	2.0	30	4	≥ 25,000	0.065	0.075	12.0
0018675	SuperSW3000	4.6	30	4	≥ 30,000	0.1 - 0.35	0.4	12.0
0022855	SuperSW mAb HTP -NEW-	4.6	15	4	≥ 15,000	0.1 - 0.35		8.0
0022854	SuperSW mAb HR -NEW-	7.8	30	4	≥ 30,000	0.5 - 1.0		12.0
0008540	G2000SW _{XL}	7.8	30	5	≥ 20,000	0.5 - 1.0	1.2	7.0
0008541	G3000SW _{XL}	7.8	30	5	≥ 20,000	0.5 - 1.0	1.2	7.0
0008542	G4000SW _{XL}	7.8	30	8	≥ 16,000	0.5 - 1.0	1.2	3.5
0016215	QC-PAK GFC 200	7.8	15	5	≥ 10,000	0.5 - 1.0	1.2	4.0
0016049	QC-PAK GFC 300	7.8	15	5	≥ 10,000	0.5 - 1.0	1.2	4.0
0005788	G2000SW	7.5	30	10	≥ 10,000	0.5 - 1.0	1.2	2.0
0005789	G3000SW	7.5	30	10	≥ 10,000	0.5 - 1.0	1.2	2.5
0005790	G4000SW	7.5	30	13	≥ 8,000	0.5 - 1.0	1.2	1.5
0005102	G2000SW	7.5	60	10	≥ 20,000	0.5 - 1.0	1.2	4.0
0005103	G3000SW	7.5	60	10	≥ 20,000	0.5 - 1.0	1.2	5.0
0005104	G4000SW	7.5	60	13	≥ 16,000	0.5 - 1.0	1.2	3.0
0006727	G2000SW	21.5	30	13	≥ 10,000	3.0 - 6.0	8.0	1.0
0006728	G3000SW	21.5	30	13	≥ 10,000	3.0 - 6.0	8.0	1.5
0006729	G4000SW	21.5	30	17	≥ 8,000	3.0 - 6.0	8.0	1.0
0005146	G2000SW	21.5	60	13	≥ 20,000	3.0 - 6.0	8.0	2.0
0005147	G3000SW	21.5	60	13	≥ 20,000	3.0 - 6.0	8.0	3.0
0005148	G4000SW	21.5	60	17	≥ 16,000	3.0 - 6.0	8.0	2.0
PEEK COLUMNS								
0020027	BioAssist G2SW _{XL}	7.8	30	5	≥ 20,000	0.5 - 1.0	1.2	7.0
0020026	BioAssist G3SW _{XL}	7.8	30	5	≥ 20,000	0.5 - 1.0	1.2	7.0
0020025	BioAssist G4SW _{XL}	7.8	30	8	≥ 16,000	0.5 - 1.0	1.2	3.5
GLASS COLUMNS								
0008800	G3000SW, Glass	8.0	30	10	≥ 10,000	0.4 - 0.8	0.8	2.0
0008801	G4000SW, Glass	8.0	30	13	≥ 8,000	0.4 - 0.8	0.8	2.0



SEC TSKgel SW SERIES ORDERING INFORMATION

► ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (µm)	
0023450	TSKgel Guardcolumn UP-SW	4.6	2.0	2	
0023451	TSKgel Guardcolumn UP-SW DC	4.6	2.0	2	
0022857	SuperSW mAb Guardcolumn	6.0	4.0	4	
0022858	SuperSW mAb Guardcolumn	3.0	2.0	4	
0022859	UltraSW Guardcolumn	6.0	4.0	3	
0018762	SuperSW Guardcolumn	4.6	3.5	4	For P/N 0018674 and 0018675
0008543	SW _{XL} Guardcolumn	6.0	4.0	7	For all SW _{XL} and P/N 0016215 and 0016049
0005371	SW Guardcolumn	7.5	7.5	10	For all 7.5 mm ID SW columns
0005758	SW Guardcolumn	21.5	7.5	13	For all 21.5 mm ID SW columns
0008805	SW Guardcolumn, Glass	8.0	4.0	10	For all 8.0 mm ID SW and QC-PAK Glass column
0018008	BioAssist SW _{XL} Guardcolumn	6.0	4.0	7	For all BioAssist SW _{XL} columns



SEC

TSKgel PW SERIES



Polymeric TSKgel PW and high resolution TSKgel PW_{XL} columns are designed for SEC of water soluble organic polymers, polysaccharides, DNA and RNA. They are based on a hydrophilic polymethacrylate matrix. Stable from pH 2 to 12, TSKgel PW series columns can be used in mobile phases of water or buffer (up to 50% polar organic solvent). A large pore G6000PW phase is available in PEEK column hardware (TSKgel BioAssist G6PW) for ultra-low sample adsorption during virus analysis. The properties of all TSKgel PW columns are summarized in Table 5.

When the molecular weight range of the sample is broad or unknown, Tosoh Bioscience offers two mixed-bed columns: The TSKgel GMPW column and its high resolution counterpart, TSKgel GMPW_{XL}, are packed with the G2500, G3000 and G6000 PW or corresponding PW_{XL} resins.

The new generation of TSKgel SuperMultiporePW columns for semi-micro SEC provide near linear calibration curves. They are packed with spherical, mono-disperse particles incorporating a proprietary multi-pore particle technology. They are ideally suited to analyze water soluble polymers, such as polyvinylpyrrolidones or dextrans.

The TSKgel PW_{XL} product line also offers specialty columns for analyzing carbohydrate oligomers (TSKgel G-Oligo-PW) and DNA and RNA fragments of 500-5000 base pairs (TSKgel G-DNA-PW). The new SuperOligoPW semi-micro SEC column featuring a small particle size has been designed to enable fast analysis of oligosaccharides and other water soluble oligomers.

TSKgel PW_{XL}-CP columns have the same base matrix as the PW_{XL} columns and were specifically developed for the analysis of water-soluble cationic polymers.

HIGHLIGHTS

- Hydrophilic spherical polymethacrylate particles
- pH range of 2-12 with up to 50% polar organic solvent
- Six different TSKgel PW pore sizes
- Linear SEC column line incorporating proprietary multipore technology
- Specialty columns for challenging SEC separations

TABLE 5

PROPERTIES AND SEPARATION RANGES OF TSKgel PW, PW_{XL} AND PW_{XL}-CP COLUMNS

TSKgel COLUMN	PARTICLE SIZE (μm)	PORE SIZE (nm)	MW RANGE (PEG/PEO)
G2000PW	12	12.5	<2 × 10 ³
G2500PW	12, 17	<20	<3 × 10 ³
G3000PW	12, 17	20	<5 × 10 ⁴
G4000PW	17	50	<3 × 10 ⁵
G5000PW	17	100	<1 × 10 ⁶
G6000PW/ BioAssist G6PW	17	>100	<8 × 10 ⁶
GMPW	17	10-100	5 × 10 ² - 8 × 10 ⁶
G2500PW _{XL}	7	<20	<3 × 10 ³
G3000PW _{XL}	7	20	<5 × 10 ⁴
G4000PW _{XL}	10	<50	<3 × 10 ⁵
G5000PW _{XL}	10	100	<1 × 10 ⁶
G6000PW _{XL}	13	>100	<8 × 10 ⁶
G-DNA-PW	10	>100	<8 × 10 ⁶
GMPW _{XL}	13	10-100	5 × 10 ² - 8 × 10 ⁶
G-Oligo-PW	7	12.5	<5 × 10 ³
SuperMultiporePW-N	4	n/a	3 × 10 ² - 5 × 10 ⁴
SuperMultiporePW-M	5	n/a	5 × 10 ² - 1 × 10 ⁶
SuperMultiporePW-H	8 (6-10)	n/a	1 × 10 ³ - 1 × 10 ⁷
SuperOligoPW	3	n/a	1 × 10 ² - 3 × 10 ³
G3000PW _{XL} -CP	7	20	2 × 10 ² - 5 × 10 ⁴
G5000PW _{XL} -CP	10	100	4 × 10 ² - 5 × 10 ⁵
G6000PW _{XL} -CP	13	>100	1 × 10 ³ - 1 × 10 ⁷




SEC TSKgel PW SERIES

CALIBRATION CURVES

Figure 15 shows the calibration curves for polyethylene glycol (PEG) and oxides (PEO) for TSKgel PW and TSKgel PW_{XL} columns, respectively. In general silica based SW type columns are recommended for the analysis of proteins, but for special applications, e.g. at basic pH or for large molecular weight proteins, PW type columns can be applied (Figure 16). Figure 17 shows the near linear calibration curves for PEG/PEO on TSKgel SuperMultiporePW columns.


▶ **FIGURE 16**



PROTEIN CALIBRATION CURVES ON TSKgel PW_{XL} COLUMNS

Column: 1. TSKgel G3000PW_{XL}, 2. TSKgel G4000PW_{XL}, 3. TSKgel G5000PW_{XL}, 4. TSKgel G6000PW_{XL}, 5. TSKgel GMPW_{XL}
 Sample: a. thyroglobulin (660,000 Da), b. γ -globulin (150,000 Da), c. albumin (67,000 Da), d. ovalbumin (43,000 Da), e. β -lactoglobulin (36,000 Da), f. myoglobin (16,900 Da), g. cytochrome C (12,400 Da)
 Mobile phase: 0.2 M phosphate buffer (pH 6.8); Flow rate: 1.0 mL/min;
 Detection: UV @ 280 nm



▶ **FIGURE 17**



CALIBRATION CURVES FOR TSKgel SuperMultiporePW

Sample: PEO & PEG standards; Mobile phase: H₂O; Flow rate: 0.6 mL/min; Detection: RI; Temperature: 25 °C

▶ **FIGURE 15**



POLYETHYLENE GLYCOL AND OXIDE CALIBRATION CURVES ON TSKgel PW AND TSKgel PW_{XL} COLUMNS

Column: TSKgel PW columns: A. G2000PW, B. G2500PW, C. G3000PW, D. G4000PW, E. G5000PW, F. G6000PW, G. GMPW, all 7.5 mm ID x 60 cm L
 Mobile phase: distilled water; Flow rate: 1.0 mL/min; Detection: RI

TSKgel PW_{XL} columns: H. G2500PW_{XL}, J. G3000PW_{XL}, K. G4000PW_{XL}, L. G5000PW_{XL}, M. G6000PW_{XL}, N. GMPW_{XL}, all 7.8 mm ID x 30 cm L

SEC TSKgel PW SERIES APPLICATIONS



LARGE DNA FRAGMENTS


For the separation of large DNA fragments greater than 1,000 base pairs, a four-column system is typically required. Baseline resolution of DNA fragments up to 7,000 base pairs can be achieved, provided there is a two-fold difference in the chain length of the fragments. Figure 14A shows the elution of double stranded DNA fragments, obtained from pBR322 DNA cleaved by both Eco RI and Bst NI, on four TSKgel G-DNA-PW columns in series. The eluted peaks were collected and subjected to polyacrylamide gel electrophoresis, which showed almost complete separation of the 1060, 1857, and 4362 base pair fragments. Although lower flow rates typically yield better separations of most fragments, the resolution of the 1857 and 4362 base pair fragments was slightly greater at the higher flow rate, as shown in Figure 18B.

OLIGOMERS

The influence of particle size on resolution and analysis time can be seen in Figure 19. It compares the separation of PEG 200 on two TSKgel G-Oligo-PW columns in series with 7 μm beads and two newly developed TSKgel SuperOligoPW semi-micro columns with a 3 μm material. The TSKgel SuperOligoPW column is designed for high resolution separations water soluble oligomers. Figure 19 demonstrates excellent resolution of the PEG 200 obtained by using the smaller, 3 μm particle size packing in the TSKgel SuperOligoPW column.

➤ **FIGURE 18**


OLIGOMERS



SEPARATION OF LARGE DNA FRAGMENTS ON TSKgel G-DNA-PW

Column: TSKgel G-DNA-PW, 10 μm, 4 x 7.8 mm ID x 30 cm L
 Sample: 60 μL of Eco RI and Bst NI cleaved pBR322 DNA,
 Base pairs: a. 4362, b. 1857, c. 1060 & 928, d. 383, e. 121, f. 13
 Mobile phase: 0.3 M NaCl, 1 mM EDTA, in 0.1 M Tris-HCl, pH 7.5,
 Flow rate: A. 0.15 mL/min, B. 0.5 mL/min; Detection: UV @ 260 nm

➤ **FIGURE 19**



ANALYSIS OF PEG 200. COMPARISON BETWEEN TSKgel SuperOligoPW AND TSKgel G-OLIGO-PW

Column: A. TSKgel SuperOligoPW, 6.0 mm ID x 15 cm L x 2
 B. TSKgel G-Oligo-PW, 7.8 mm ID x 30 cm L x 2
 Mobile phase: H₂O; Flow rate: A: 0.6 mL/min, B: 1.0 mL/min
 Detection: RI; Temp.: 25°C; Inj. vol.: A: 20 μL, B: 100 μL



SEC TSKgel PW-CP SERIES


TSKgel PW_{XL}-CP size exclusion columns were specifically developed for the analysis of water soluble cationic polymers. Three columns are available within the TSKgel PW_{XL}-CP series, each with a different particle size, separation range and exclusion limit, allowing polymers within a wide molecular mass range to be separated and characterized.

When using conventional SEC columns the analysis of cationic polymers requires a high salt concentration in the mobile phase to prevent adsorption of the polymers onto the particles in SEC columns. The TSKgel PW_{XL}-CP columns eliminate ionic adsorption onto the particle by incorporating a cationic functionality on the particle surface. This modification results in high recovery for cationic polymers and enables elution under low salt conditions.

These columns show high theoretical plate numbers, linear calibration curves and high durability. The base resin is the same as that used in the TSKgel PW_{XL} columns. Figure 20 shows the calibration curves for PEG/PEO calibration curves obtained with TSKgel PW_{XL}-CP columns.

Figure 21 demonstrates that these SEC columns can be utilized for the analysis of a wide variety of cationic polymers. Various cationic polymers with different functional groups and molecular weights were injected on the three TSKgel PW_{XL}-CP columns (TSKgel G6000PW_{XL}-CP, G5000PW_{XL}-CP and G3000PW_{XL}-CP, connected in series).


FIGURE 20



CALIBRATION CURVE FOR TSKgel PW_{XL}-CP COLUMNS

Columns: TSKgel G3000PW_{XL}-CP, 7 μ m; TSKgel G5000PW_{XL}-CP, 10 μ m; TSKgel G6000PW_{XL}-CP, 13 μ m; Samples: polyethylene oxides (PEO) standards; polyethylene glycols (PEG) standards
Mobile phase: 0.1 mol/L NaNO₃; Flow rate: 1 mL/min; Detection: RI; Temp: 25 °C

FIGURE 21



ANALYSIS OF CATIONIC POLYMERS

Columns: TSKgel G3000PW_{XL}-CP, 7 μ m (7.8 mm ID x 30 cm L), TSKgel G5000PW_{XL}-CP, 10 μ m (7.8 mm ID x 30 cm L), TSKgel G6000PW_{XL}-CP, 13 μ m (7.8 mm ID x 30 cm L); Mobile phase: 0.1 mol/L NaNO₃; Flow rate: 1 mL/min; Detection: RI; Temperature: 25 °C; Sample Load: 3 g/L, 100 μ L

SEC TSKgel SuperMultiporePW SERIES



The new TSKgel SuperMultiporePW column line is incorporating Tosoh's proprietary multi-pore particle technology. These semi-micro SEC columns provide near linear calibration curves. They are ideally suited to analyze the molecular weight and the MW distribution of water soluble polymers, such as polyvinylpyrrolidones or dextrans.


TSKgel SuperMultiporePW columns are packed with spherical mono-disperse polymethacrylate particles, each containing a wide range of pore sizes. They belong to the semi-micro type of SEC columns (6 mm ID, 15 cm length) providing high theoretical plate numbers at half of the length of a conventional SEC column. The TSKgel SuperMultiporePW series comprises of three column types covering different molecular weight ranges (PW-N; PW-M, PW-H).

Multi-pore particle technology is the most elegant way to achieve near linear SEC calibration curves. It solves the known problem of peak disturbances/inflection points, which typically occur due to a mismatch of pore sizes when columns with different molecular weight ranges are coupled. Particles produced by multi-pore technology contain a broad range of pore sizes in a single polymeric bead. This innovative approach essentially creates a linear calibration curve within each particle (Figure 22).

Multi-pore, semi-micro SEC columns provide high resolution and smooth peak shapes without shoulders or inflection points. This leads to better accuracy and reproducibility when determining the molecular mass distribution of water soluble polymers.

Figure 23 shows the SEC analysis of a real sample - Polyvinylpyrrolidone (PVP) K-30 - on a series of conventional TSKgel G3000PW_{XL} and G5000PW_{XL} columns compared to the one obtained with a single TSKgel SuperMultiporePW-M semi-micro linear SEC column (MW range 600,000 – 1,500,000). On a series of conventional SEC columns the Polyvinylpyrrolidone peak shows an inflection point, which does not appear on the SuperMultiporePW-M column. Analysis is much faster and more sensitive when applying the new multi-pore packing.

FIGURE 23





ANALYSIS OF POLYVINYLPIRROLIDONE

Columns: TSKgel SuperMultiporePW-M, 6 mm ID x 15 cm L x 1 (red) TSKgel G3000PW_{XL} & G5000PW_{XL}, each 7.8 mm ID x 30 cm L in line (blue); Sample: Polyvinylpyrrolidone (K-30); Mobile phase: 0.1 mol/L NaNO₃; Flow rate: 0.6 mL/min; Detection: RI

FIGURE 22

STRATEGIES FOR WIDE RANGE SEPARATION USING SIZE EXCLUSION CHROMATOGRAPHY





SEC TSKgel PW SERIES ORDERING INFORMATION

► ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (µm)	NUMBER THEORETICAL PLATES	FLOW RATE (mL/min)		MAXIMUM PRESSURE DROP (MPa)
						RANGE	MAX.	
STAINLESS STEEL COLUMNS								
0022789	SuperMultiporePW-N	6.0	15	4	>16,000	0.3 - 0.6	0.6	4.5
0022790	SuperMultiporePW-M	6.0	15	5	>12,000	0.3 - 0.6	0.6	2.7
0022791	SuperMultiporePW-H	6.0	15	8 (6-10)	>7,000	0.3 - 0.6	0.6	0.9
0022792	SuperOligoPW	6.0	15	3	>16,000	0.3 - 0.6	0.6	5.0
0008031	G-Oligo-PW	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0008032	G-DNA-PW	7.8	30	10	≥ 10,000	0.2 - 0.5	0.6	2.0
0008020	G2500PW _{XL}	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0008021	G3000PW _{XL}	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0008022	G4000PW _{XL}	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	2.0
0008023	G5000PW _{XL}	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	2.0
0008024	G6000PW _{XL}	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
0008025	GMPW _{XL}	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
0021873	G3000PW _{XL} -CP	7.8	30	7	≥ 16,000	1.0 5.5		
0021874	G5000PW _{XL} -CP	7.8	30	10	≥ 10,000	1.0 2.5		
0021875	G6000PW _{XL} -CP	7.8	30	13	≥ 7,000	1.0 2.0		
0005761	G2000PW	7.5	30	12	≥ 5,000	0.5 - 1.0	1.2	2.0
0008028	G2500PW	7.5	30	12	≥ 5,000	0.5 - 1.0	1.2	2.0
0005762	G3000PW	7.5	30	12	≥ 5,000	0.5 - 1.0	1.2	2.0
0005763	G4000PW	7.5	30	17	≥ 3,000	0.5 - 1.0	1.2	1.0
0005764	G5000PW	7.5	30	17	≥ 3,000	0.5 - 1.0	1.2	1.0
0005765	G6000PW	7.5	30	17	≥ 3,000	0.5 - 1.0	1.2	1.0
0008026	GMPW	7.5	30	17	≥ 3,000	0.5 - 1.0	1.2	1.0
0005105	G2000PW	7.5	60	12	≥ 10,000	0.5 - 1.0	1.2	4.0
0008029	G2500PW	7.5	60	12	≥ 10,000	0.5 - 1.0	1.2	4.0
0005106	G3000PW	7.5	60	12	≥ 10,000	0.5 - 1.0	1.2	4.0
0005107	G4000PW	7.5	60	17	≥ 6,000	0.5 - 1.0	1.2	2.0
0005108	G5000PW	7.5	60	17	≥ 6,000	0.5 - 1.0	1.2	2.0
0005109	G6000PW	7.5	60	17	≥ 6,000	0.5 - 1.0	1.2	2.0
0008027	GMPW	7.5	60	17	≥ 6,000	0.5 - 1.0	1.2	2.0
0008030	G2500PW	21.5	60	17	≥ 10,000	1.6 - 6.0	8.0	2.0
PEEK								
0020024	BioAssist G6PW	7.8	30	17	≥ 3,000	0.5 - 1.0	1.2	10
GUARD COLUMNS								
0022793	SuperMP (PW)-N Guard column	4.6	3.5	4				
0022794	SuperMP (PW)-M Guard column	4.6	3.5	5				
0022795	SuperMP (PW)-H Guard column	4.6	3.5	8				
0022796	SuperOligoPW Guard column	4.6	3.5	3				
0008034	Oligo Guard column	6.0	4.0	13	For 7.8 mm ID G-Oligo-PW columns			
0008033	PW _{XL} Guard column	6.0	4.0	12	For 7.8 mm ID PW _{XL} & G-DNA-PW (TSKgel G3000PW packing)			
0021876	PW _{XL} -CP Guard column	6.0	4.0	13	For 7.8 mm ID PW _{XL} -CP columns			
0006763	PW-L Guard column	7.5	7.5	13	For 7.5 mm ID G1000PW & G2000PW (TSKgel G2000PW packing)			
0006762	PW-H Guard column	7.5	7.5	13	For 7.5 mm ID G2500PW through GMPW columns			
0006758	PW-H Guard column	21.5	7.5	17	For 21.5 mm ID G2500PW through G5000PW columns			



SEC TSKgel Alpha & SuperAW SERIES

ORDERING INFORMATION

► ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (µm)	NUMBER THEORETICAL PLATES	FLOW RATE (mL/min)		MAXIMUM PRESSURE DROP (MPa)
						RANGE	MAX.	
STAINLESS STEEL COLUMNS								
0018339	Alpha-2500	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018340	Alpha-3000	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018341	Alpha-4000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018342	Alpha-5000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018343	Alpha-6000	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
0018344	Alpha-M (mixed bed)	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
GUARD COLUMNS								
0018345	Alpha Guard column	6.0	4	13	For all Alpha columns			
VMPAK COLUMNS*								
0020011	VMpak-25	2.0	5	7	≥ 1,000	0.1 - 0.2	0.25	2.0
0020012	VMpak-25	2.0	15	7	≥ 3,000	0.1 - 0.2	0.25	6.0
STAINLESS STEEL COLUMNS								
0019315	SuperAW2500	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019316	SuperAW3000	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019317	SuperAW4000	6.0	15	6	≥ 10,000	0.3 - 0.6		4.0
0019318	SuperAW5000	6.0	15	7	>10,000	0.3 - 0.6		3.0
0019319	SuperAW6000	6.0	15	9	>7,000	0.3 - 0.6		2.0
0019320	SuperAWM-H	6.0	15	9	>7,000	0.3 - 0.6		2.0
GUARD COLUMNS								
0019321	SuperAW-L Guard Column		4.6	3.5	7	For SuperAW2500-4000 columns.		
0019322	SuperAW-H Guard Column		4.6	3.5	13	For SuperAW5000-AWM-H columns		

*TSKgel VMpak-25 series contains a similar packing as TSKgel Alpha-2500. It can be used for multimodal LC and LC-MS separations.



SEC TSKgel Alpha & SuperAW SERIES

ORDERING INFORMATION

► ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (µm)	NUMBER THEORETICAL PLATES	FLOW RATE (mL/min)		MAXIMUM PRESSURE DROP (MPa)
						RANGE	MAX.	
STAINLESS STEEL COLUMNS								
0018339	Alpha-2500	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018340	Alpha-3000	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018341	Alpha-4000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018342	Alpha-5000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018343	Alpha-6000	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
0018344	Alpha-M (mixed bed)	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
GUARD COLUMNS								
0018345	Alpha Guard column	6.0	4	13	For all Alpha columns			
VMPAK COLUMNS*								
0020011	VMpak-25	2.0	5	7	≥ 1,000	0.1 - 0.2	0.25	2.0
0020012	VMpak-25	2.0	15	7	≥ 3,000	0.1 - 0.2	0.25	6.0
STAINLESS STEEL COLUMNS								
0019315	SuperAW2500	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019316	SuperAW3000	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019317	SuperAW4000	6.0	15	6	≥ 10,000	0.3 - 0.6		4.0
0019318	SuperAW5000	6.0	15	7	>10,000	0.3 - 0.6		3.0
0019319	SuperAW6000	6.0	15	9	>7,000	0.3 - 0.6		2.0
0019320	SuperAWM-H	6.0	15	9	>7,000	0.3 - 0.6		2.0
GUARD COLUMNS								
0019321	SuperAW-L Guard Column		4.6	3.5	7	For SuperAW2500-4000 columns.		
0019322	SuperAW-H Guard Column		4.6	3.5	13	For SuperAW5000-AWM-H columns		

*TSKgel VMpak-25 series contains a similar packing as TSKgel Alpha-2500. It can be used for multimodal LC and LC-MS separations.

SEC TSKgel Alpha & SuperAW SERIES




The TSKgel Alpha and SuperAW column series offer a new alternative for performing SEC. The columns are packed with a hydrophilic, highly crosslinked polymer which is compatible to a wide range of solvents ranging from pure aqueous up to 100 % organic mobile phases (see Figure 24). Both series consist of six columns with different pore sizes, spanning a wide MW separation range from 100 to over 1,000,000 Da when using polyethylene glycol (PEG) as a standard. Exclusion limits for polyethylene oxides in water and other physical properties for the Alpha and SuperAW columns are listed in Table 6.

The TSKgel Alpha and SuperAW column series can be used for separations of synthetic polymers, oligomers, additives and detergents as well as for saccharides, nucleic acids and peptides. TSKgel SuperAW columns with reduced particle size and semi-micro column dimensions of 6 mm ID and 15 cm length provide short analysis times and higher resolution power. For samples with big differences in molecular weights, the mixed bed columns TSKgel Alpha-M and TSKgel SuperAWM-H show linear calibration curves over the whole range.

HIGHLIGHTS

- Unique hydrophilic polymethacrylate resin with rigid spherical beads
- Minimal swelling characteristics from 100% water to 100% non-polar solvents
- Excellent mechanical and chemical stability
- TSKgel SuperAW columns with reduced particle size and shorter columns length provide short analysis times and high resolution power

➤ **FIGURE 24**



SOLVENT COMPATABILITY OF TSKgel Alpha-3000 WITH ORGANIC SOLVENT

Conditions for solvent change: Flow rate: 1.0 mL/min; Temp.: 25 °C; Time for purge: 8 h; Conditions for TP measurement: Sample: ethylene glycol; Flow rate: 1.0 mL/min; Temp.: 25 °C; Detection: RI

➤ **TABLE 6**

PROPERTIES AND SEPARATION RANGES OF TSKgel Alpha AND SuperAW-SERIES

TSKgel COLUMN	ID (mm) X LENGTH (cm L)	PARTICLE SIZE (µm)	MIN NO. THEORET. PLATES	EXCLUSION LIMIT (PEO/H ₂ O)
Alpha-2500	7.8 x 30	7	16,000	5 x 10 ³
Alpha-3000	7.8 x 30	7	16,000	9 x 10 ⁴
Alpha-4000	7.8 x 30	10	10,000	4 x 10 ⁵
Alpha-5000	7.8 x 30	10	10,000	1 x 10 ⁶
Alpha-6000	7.8 x 30	13	7,000	>1 x 10 ⁷
Alpha-M	7.8 x 30	13	7,000	>1 x 10 ⁷
SuperAW2500	6.0 x 15	4	>16,000	5 x 10 ³
SuperAW3000	6.0 x 15	4	>16,000	9 x 10 ⁴
SuperAW4000	6.0 x 15	6	>10,000	1 x 10 ⁶
SuperAW5000	6.0 x 15	7	>10,000	1 x 10 ⁶
SuperAW6000	6.0 x 15	9	>6,000	1 x 10 ⁷
SuperAWM-H	6.0 x 15	9	>6,000	1 x 10 ⁷



SEC TSKgel Alpha & SuperAW SERIES

ORDERING INFORMATION

► ORDERING INFORMATION

PART #	DESCRIPTION	ID (mm)	LENGTH (cm)	PARTICLE SIZE (µm)	NUMBER THEORETICAL PLATES	FLOW RATE (mL/min)		MAXIMUM PRESSURE DROP (MPa)
						RANGE	MAX.	
STAINLESS STEEL COLUMNS								
0018339	Alpha-2500	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018340	Alpha-3000	7.8	30	7	≥ 16,000	0.5 - 0.8	1.0	4.0
0018341	Alpha-4000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018342	Alpha-5000	7.8	30	10	≥ 10,000	0.3 - 0.6	1.0	3.0
0018343	Alpha-6000	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
0018344	Alpha-M (mixed bed)	7.8	30	13	≥ 7,000	0.3 - 0.6	1.0	2.0
GUARD COLUMNS								
0018345	Alpha Guard column	6.0	4	13	For all Alpha columns			
VMPAK COLUMNS*								
0020011	VMpak-25	2.0	5	7	≥ 1,000	0.1 - 0.2	0.25	2.0
0020012	VMpak-25	2.0	15	7	≥ 3,000	0.1 - 0.2	0.25	6.0
STAINLESS STEEL COLUMNS								
0019315	SuperAW2500	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019316	SuperAW3000	6.0	15	4	≥ 16,000	0.3 - 0.6		6.0
0019317	SuperAW4000	6.0	15	6	≥ 10,000	0.3 - 0.6		4.0
0019318	SuperAW5000	6.0	15	7	>10,000	0.3 - 0.6		3.0
0019319	SuperAW6000	6.0	15	9	>7,000	0.3 - 0.6		2.0
0019320	SuperAWM-H	6.0	15	9	>7,000	0.3 - 0.6		2.0
GUARD COLUMNS								
0019321	SuperAW-L Guard Column		4.6	3.5	7	For SuperAW2500-4000 columns.		
0019322	SuperAW-H Guard Column		4.6	3.5	13	For SuperAW5000-AWM-H columns		

*TSKgel VMpak-25 series contains a similar packing as TSKgel Alpha-2500. It can be used for multimodal LC and LC-MS separations.

SEC

OPTIMIZING SEC



SAMPLE LOAD


As SEC is a partition chromatography, sample load on the column is limited. High sample loads distort peak shapes and cause an overall decrease in efficiency due to column overload. Optimal sample load highly depends on the sample properties (sample matrix) and the separation task. For analytical columns, sample concentrations of 1-20 mg/ml are recommended. Proteins can be loaded at higher concentrations and higher total loads than synthetic macromolecules. For preparative purposes for example, 100 mg of BSA can be loaded on two 21.5 mm ID x 60 cm L TSKgel G3000SW columns, but only 20 mg of PEG 7500.

Sample volume depends very much on the type of column. On TSKgel SuperSW columns for example, a 5 μ L injection volume ensures optimal results. Standard injection volumes for 7.5 and 7.8 mm ID columns are 20-100 μ L, whereas for preparative purposes on 21.5 mm ID columns, injection volumes may be raised up to 2 ml.

MOBILE PHASE

Proper selection of the mobile phase is necessary to maximize molecular sieving mechanism and to minimize secondary effects such as ionic and hydrophobic interaction between the sample and the column packing material. For each sample there will be an optimum buffer type and concentration that results in the highest resolution and recovery.

➤ **FIGURE 24**



INFLUENCE OF MOBILE PHASE

A: No ethanol in mobile phase; B: 10% ethanol in mobile phase
 Column: TSKgel G3000SWXL columns, 5 μ m, 7.8 mm ID x 30 cm L
 Sample: 10 mL PEG r-HuMGDF;
 — initial injection; after 150 injections
 Mobile phase: 0.1 M sodium phosphate, pH 6.9, 0.5 M NaCl,
 Flow rate: 0.7 mL/min; Detection: UV @ 220 nm

For TSKgel SW columns mobile phases a buffer concentration between 0.1 M and 0.5 M is recommended. Under low ionic strength (< 0.1 M), ionic interactions between the sample molecules and the silica surface may occur. Under conditions of high ionic strength (>1.0 M), hydrophobic interactions are more likely to occur. A neutral salt, such as sodium sulphate may be added to the buffer to increase buffer ionic strength. Also the ionic species of the buffer has an effect on the separation. As a good starting point, a 0.1 M sodium phosphate buffer together with 0.1 M sodium sulphate has proved to be of value.

As the polymeric TSKgel PW and Alpha-type resins carry less residual charged groups on the surface than silica gels, salt concentration of the mobile phase can be lower. Non-ionic, non-polar compounds such as polyethylene glycols can simply be analyzed with distilled water. For ionic polymeric compounds, a neutral salt such as sodium nitrate is added to the aqueous eluent. Generally, a concentration of 0.1 M to 0.2 M is sufficient to overcome undesirable ionic interactions.


If hydrophobic interaction occurs between the sample and the column matrix, a water soluble organic solvent can be added to the mobile phase. The addition of acetonitrile, acetone, ethanol or methanol up to a concentration of 20% may also prevent columns from fouling by suppressing interaction of hydrophobic impurities of the sample. An example is shown in Figure 24 with the analysis of a pegylated protein on a TSKgel G3000SWXL column. As pegylated products are more hydrophobic, they tend to interact with the column matrix. Over time the pegylated product can foul the column, which is indicated by shifts of retention time and decreasing separation performance. By adding 10% of ethanol to the elution buffer, this problem is overcome and no differences in performance at the first and the 150th injection are observed (courtesy of J.J. Ratto et al. Amgen Inc., 1996).

COLUMN PROTECTION

To protect the column and increase its lifetime, the use of a guard column is strongly recommended. Sample purity, sample load and the composition of the mobile phase have an influence on column lifetime. For information on TSKgel SEC columns for GPC analysis of organic polymers please refer to the TSKgel GPC column brochure.



TOSOH



A N A L Y S E N T E C H N I K

AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de

TOSOH BIOSCIENCE

Im Leuschnerpark 4 64347 Griesheim, Germany
Tel: +49 6155-7043700 Fax: +49 6155-8357900
info.tbg@tosoh.com www.tosohbioscience.de

B15L34A