

FORENSICS



CLEAN SCREEN XCEL®

EXCELLENCE JUST GOT BETTER



INNOVATION THROUGH CHEMISTRY

CLEAN SCREEN XCEL® COLUMNS

Clean Screen Xcel® solid phase extraction columns are designed to reduce the number of steps in the extraction. The result is a column that reduces sample prep times versus traditional SPE and minimizes the amount of solvent necessary. Additional advantages include reduced sample size and improved cleanliness and recovery versus dilute and shoot.

BENEFITS

- Sorbent conditioning step is eliminated
- Reduced sample size
- Increased sensitivity vs. dilute and shoot
- Decreased extraction steps
- Increased recovery values vs. dilute and shoot

CLEAN SCREEN XCEL® I:

Extracts a wide array of basic drugs including benzodiazepines and opiates.

Organic Loading = 12.4%

Average Pore Size = 60 Å

Surface Area = 500 m²/g

Pore Volume = 0.77 cm³/g

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	Part Number	
1	130	100	CSXCE111	
3	130	50	CSXCE103	
6	130	50	CSXCE106	
6	200	50	CSXCE206	
10	130	50	ZSXCE010	
WELL PLATES				
Number of Wells	Sorbent Amount (mg)	Units per Pack	Extended Drip Tip	Part Number
48	130	1	NO	WSH48XCE11
96	80	1	YES	WSH96XCE108-LD
96	130	1	NO	WSH96XCE11
96	130	1	YES	WSH96XCE11-LD

CLEAN SCREEN XCEL® II:

Designed solely for the extraction of the THC metabolite.

Organic Loading = 16.7%

Average Pore Size = 60 Å



Surface Area = 500 m²/g

Pore Volume = 0.77 cm³/g

COLUMNS				
Tube Volume (mL)	Sorbent Amount (mg)	Units per Pack	Part Number	
1	130	100	CSXCE211	
3	130	50	CSXCE2103	
6	130	50	CSXCE2106	
6	200	50	ZSXCE2010	
10	130	50	ZSXCE010	
WELL PLATES				
Number of Wells	Sorbent Amount (mg)	Units per Pack	Extended Drip Tip	Part Number
48	130	1	NO	WSH48XCE211
96	80	1	YES	WSH96XCE208-LD
96	130	1	NO	WSH96XCE211

CLEAN SCREEN XCEL®

uses less solvent and has shorter prep times than either traditional SPE or SLE extraction techniques.



SCREENING METHOD FOR 121 ACIDIC, NEUTRAL AND BASIC DRUG ANALYTES IN PLASMA, SERUM, URINE, OR TISSUE BY LC-MS/MS

Part

CSXCE111 – CLEAN SCREEN XCEL® I 130 mg, 1 mL Tube
 BETA-GLUC-10 – SELECTRAZYME® Beta-glucuronidase
 SLDA50ID21-5UM – SELECTRA® DA HPLC Column, 50 x 2.1 mm, 5 µm
 SLDAGDC21-5UM - SELECTRA® DA 5 µm Guard Cartridge
 SLGRDHLDR - Guard Cartridge Holder



CLEAN SCREEN XCEL® I Columns

1. PREPARE SAMPLE:

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards
 Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate
 Mix/vortex and let stand for 5 minutes
 Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex
 Sample pH should be 6.0 ± 0.5.
 Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate. Centrifuge for 10 minutes at 2000 rpm and discard pellet
NOTE: See Hydrolysis step if required

Hydrolysis (for urine samples only): To 1-2 mL of urine sample, add 1 mL of acetate buffer (pH 5.0) containing 5,000 units/mL Selectrazyme® β-glucuronidase.
 Optionally, add 1 mL of acetate buffer and 25-50 µL of concentrated β-glucuronidase.
 Vortex and heat for 1-2 hours at 65°C.
 (Hydroxylamine can be added to sample here if oxime derivative is preferred.) Allow sample to cool

2. APPLY SAMPLE:

Load sample directly to column without any preconditioning
 Pull sample through at a rate of 1-2 mL/ minute
 Dry column thoroughly under full vacuum or positive pressure for 1 minute

3. WASH 1 – ACIDIC & NEUTRAL COMPOUNDS (FRACTION 1):

Add 1 x 1 mL of DI H₂O
 Apply pressure to column for ~1 minute (either vacuum (10mm Hg) or positive pressure(~80-100psi). This ensures that the entire sample and any residual is pulled through to waste
 Add 1 x 1 mL of 0.1M Acetic Acid
 Apply pressure to column for ~1minute (either vacuum (10mm Hg) or positive pressure(~80-100psi).
 Add 1 x 2 mL Hexane to remove residual aqueous phase
 Dry column (5 minutes at full vacuum or pressure)

4. ELUTION 1 – ACIDIC & NEUTRAL COMPOUNDS (FRACTION 1):

Add 1 x 1 mL Ethyl Acetate: Hexane (50:50) Collect eluate at 1 to 2 mL/minute

5. DRY ELUTE:

Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C Reconstitute with 100 µL of Mobile Phase

6. WASH 2 - BASIC COMPOUNDS (FRACTION 2):

Add 1 x 1 mL of 2% Acetic Acid/98% Methanol
 Dry column 5 minutes at full vacuum (10mm Hg) or positive pressure (~80-100 psi)

7. ELUTION 2 - BASIC COMPOUNDS (FRACTION 2):

1 x 1 mL of CH₂Cl₂/ IPA/ Ammonium Hydroxide (78/20/2).

8. DRY ELUTE:

Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C. Take care not to overheat or over evaporate. Certain compounds are heat labile, such as the amphetamines and phencyclidine. Reconstitute with 100 µL of Mobile Phase

NOTES:


(1) Fraction 1 (Acid Neutrals) and Fraction 2 (Bases) can be combined together if need be. This is not generally recommended as the Acid/ Neutral fraction tends to be dirtier than the Basic one, so for more effective results, keep fractions separate.

(2) A keeper solvent such as DMF can be used to prevent the volatilization of amphetamines and phencyclidine. Use 30-50 µL of high purity DMF in the sample (Fraction 2) before evaporation.

(3) A 1% HCl in CH₃OH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs. Add 1 drop of the solution prior to evaporating then continue to dryness.

(4) The hexane wash step can be removed if user is looking to analyze for Parent THC.

(5) To extract the benzodiazepine group at higher recovery, following the elution of the Acidic/Neutral drugs, a second elution can be done prior to moving on to the second wash phase. The second elution solvent would consist of 98% Ethyl Acetate/ 2% Ammonium Hydroxide.



HPLC Conditions

Instrument: Shimadzu HPLC 20-AD

Detector: AB Sciex API 3200 Qtrap MS/MS

LC Column: UCT Selectra® DA HPLC Column 50 x 2.1mm, 5 µm

Ionization mode: ESI+

Injection Volume: 20 µL

Flow Rate: 0.6 mL/minute

Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in MeOH

Gradient:

Time	%A	%B
0.00	90	10
0.50	90	100
4.00	60	40
7.50	15	85
8.50	10	90
8.51	90	10
10.00	STOP	

Limits of detection for each compound class†

Representative Compounds	LOD Range (ng/mL)	Representative Compounds	LOD Range (ng/mL)
Amphetamine/Methamphetamine	10	Benzodiazepines (24)	5-10
MDMA/MDA/MDEA	5-10	Fentanyl/Norfentanyl	1-5
Opiates (9)*	2	Carisoprodol/Meprobamate	25-100
Methadone/EDDP	5-20	Cathinones/Phenethylamines (10)	10-20
Sympathomimetic Amines	20	Anti-Psychotics (7)	10-50
Meperidine/Normeperidine	20	Z Drugs/Insomnia Treatment (3)	20-25
Cocaine/BE/EME	5	Pain Management Compounds (14)	2-20
Tricyclic Antidepressants (16)	5-50	Caffeine/Theobromine/Theophylline	20-50
Gabapentin/Pregabalin	50	PCP/Ketamine/Norketamine	5-10
Acetaminophen	50	Anti-histamines	10-20

()* - Number of compounds analyzed in this class.

†For the full application note with the complete compound list and LOD's contact your regional representative or download from our website, www.unitedchem.com:

SCREENING METHOD FOR 121 ACIDIC, NEUTRAL AND BASIC DRUG ANALYTES IN PLASMA, SERUM, URINE, OR TISSUE BY LC-MS/MS.

CARBOXY-THC IN URINE BY LC-MS/MS OR GC-MS USING CLEAN SCREEN XCEL® II EXTRACTION COLUMN

Part

CSXCE2106 – CLEAN SCREEN XCEL® II 130 mg 6 mL Tube
 SPHPHO6001-5 – Select pH Buffer Pouches 100 mM Phosphate pH 7.0
 SLDA50ID21-5UM - Selectra® DA HPLC Column, 50 x 2.1 mm, 5 µm
 SMSTFA-1-1 – SELECTRA-SIL® MSTFA w/ 1% TMCS
 SBSTFA-1-1 – SELECTRA-SIL® BSTFA w/ 1% TMCS
 SLDAGDC21-5UM - SELECTRA® DA 5 µm Guard Cartridge
 SLGRDHLDR - Guard Cartridge Holder



CLEAN SCREEN XCEL® II Columns

1. PREPARE SAMPLE-ENZYME HYDROLYSIS OF GLUCURONIDES:

To 2 mL of urine add internal standard and 50 µL of 10 M NaOH
 Mix/vortex
 Hydrolyze for 15 minutes at 60-70 °C. Cool before proceeding
 Adjust sample pH to 7.0 with 50 µL of 1:1 H₂O: Glacial Acetic Acid.
 Add 200 µL pH 7.0 100mM Phosphate Buffer (pH should be ~7.0)

2. APPLY SAMPLE:

Load at 1 to 2 mL/minute
 Dry column (2 minutes at full vacuum or pressure)

3. WASH COLUMN:

1 x 2 mL Hexane
 Dry Column at full vacuum or pressure for 10 minutes

4. ELUTE ANALYTE:

1 x 3 mL Hexane/ Ethyl Acetate/ Glacial Acetic Acid (49:49:2)
 Collect eluate at 1 to 2 mL/minute

5. DRY ELUATE:

Evaporate to dryness at < 40 °C

6. RECONSTITUTE / DERIVATIZE:

- LC-MS/MS: Reconstitute sample in 100 µL of mobile phase Inject 20 µL.
- GC-MS: Dissolve residue in 50 µL of Ethyl Acetate and 50 µL MSTFA (with 1%TMCS)
 Overlay with N₂ and cap. Mix/vortex
 React 30 minutes at 70 °C; Cool and inject 1 µL

Alternate Derivatization


1. Form TMS Derivatives by adding 50 µL BSTFA w/ 1% TMCS and 50 µL of Ethyl Acetate; React 45 minutes at 70 °C

THC-delta-9-COOH

Comparison of Extraction Methods - XCEL II vs. Liquid/Liquid

Sample ID#	Liquid / Liquid Result (ng/mL)	2 Step XCEL I SPE Result (ng/mL)	Liquid / Liquid Area Response	2 Step XCEL I SPE Area Response
60 control	60.8	63.2	193689	521993
300 control	226	299	663859	1793767
Patient #1	157	171	366285	1210294
Patient #2	26	30.1	74689	283042
Patient #3	57	73	125662	497268
Patient #4	13.8	14	50238	168674
Patient #5	57.9	56.2	148767	469350
30 control	26.5	27	95809	215516

CARBOXY-THC



	Analyte	MRM Transitions		Relative Retention Time (minutes)
		Q1	Q3	
1.	THC-DELTA 9-COOH D9	352	308	2.44
2.	THC-DELTA 9-COOH	343	299	2.49

HPLC Conditions

Instrument: Shimadzu Prominence UFLC

Detector: API 3200 Qtrap MS/MS

HPLC Column: SELECTRA® DA HPLC Column 50 x 2.1 mm 5 µm

Polarity: Negative

Reconstitute: 100 µL

Injection Volume: 20 µL

Flow Rate: 0.5 mL/minute

Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in MeOH

Gradient:

Time	%A	%B
0.00	60	40
2.00	30	70
2.50	10	90
2.51	60	40
4.00	STOP	

SYMPATHOMIMETIC AMINES IN BLOOD, PLASMA/ SERUM, URINE, OR TISSUE BY LC-MS/MS CLEAN SCREEN XCEL® I EXTRACTION COLUMN

Part

CSXCE111 – CLEAN SCREEN XCEL® 130 mg, 1 mL Tube
 SPHPHO6001-5 – Select pH Buffer Pouches 100 mM Phosphate pH 6.0
 PFAA-0-1 – SELECTRA-SIL® PFAA
 SPFPOH-1 – SELECTRA-SIL® PFPOH
 SHFAA-0-1– SELECTRA-SIL® HFAA
 SBSTFA-1-1 – SELECTRA-SIL® BSTFA w/ 1% TMCS
 SLDA100ID21-3UM – SELECTRA® DA HPLC Column, 100 x 2.1 mm, 3 µm
 SLDAGDC20-3UM - SELECTRA® DA 3 µm Guard Cartridge
 SLGRDHLDR - Guard Cartridge Holder



Select pH Buffer Pouches

1. PREPARE SAMPLE:

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards
 Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate
 Mix/vortex and let stand for 5 minutes
 Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex
 Sample pH should be 6.0 ± 0.5.
 Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
 Centrifuge for 10 minutes at 2000 rpm and discard pellet

2. APPLY SAMPLE:

Load sample directly to column without any preconditioning.
 Pull sample through at a rate of 1-2 mL/ minute.
 Dry column thoroughly under full vacuum or positive pressure for 1 minute.

3 WASH:

1 x 3 mL 98% Methanol: 2% Acetic Acid
 Dry column thoroughly under full vacuum or positive pressure for a minimum of 5 minutes.

4. ELUTION:

1 x 3 mL CH₂Cl₂/ IPA/ NH₄OH (78:20:2)
 Collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily.
 Add IPA/ NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5. DRY ELUTE:

Add 50 µL of 1% HCl in CH₃OH to each tube
 Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.

NOTE: A 1% HCl in CH₃OH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs.


6. RECONSTITUTE / DERIVATIZE

- LC-MS/MS: Reconstitute sample in 100 µL of mobile phase Inject 5-20 µL.
- GC-MS: Fluoroacylate with PFAA (PFAA)
 Add 50 µL PFAA. Over lay with N₂ and cap
 *Improve derivatization by addition of PFPOH
 React 20 minutes at 70 °C. Evaporate to dryness <40 °C
 Reconstitute with 100 µL Ethyl Acetate.

NOTES:

1. It is important to dry the column thoroughly to achieve the highest recovery of all compounds. Any residual moisture will slow down the drying of the elution solvents prior to derivatization for GC/MS analysis, if being used. Also, any residual moisture could reduce the reactivity of the derivatization agent resulting in low GC/MS sensitivity and inconsistent results.
2. **Sodium periodate can be added to sample during preparation if oxidation is preferred.**
 To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standard(s). Add 2 mL of urine and 1 mL 0.35 M sodium periodate.
 Mix/vortex
 Incubate at room temp. for 20 min.
 Sample pH should be 6.0 ± 0.5.
 Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate
 Sample is now ready to be added to the extraction column
3. **Alternate Derivatization**
 1. Fluoroacylate with HFAA (HFAA)
 Add 50 µL HFAA. Over lay with N₂ and cap
 *Improved derivatization by addition of PFPOH
 React 20 minutes at 70 °C. Evaporate to dryness <40 °C
 Reconstitute with 100 µL Ethyl Acetate
 2. Form TMS Derivatives by adding 50 µL BSTFA w/ 1% TMCS and 50 µL of Ethyl Acetate; React 45 minutes at 70 °C

SYMPATHOMIMETIC AMINES



Analyte	MRM Transitions		Relative Retention Time (minutes)
	Q1	Q3	
1. Epedrine	166.0	115.0	4.45
2. Flephedrone	182.0	149.0	4.78
3. Methamphetamine	150.1	91.1	5.18
4. Methylone	208.0	132.0	5.26
5. MDA	180.1	77.0	5.27
6. MDMA	194.1	77.0	5.69
7. Butylone	222.0	131.0	5.84
8. Ethylone	222.0	131.0	6.36
9. Mephedrone	178.1	145.0	6.10
10. MDEA	208.1	105.0	6.48
11. Methcathinone	192.1	144.0	6.93
12. Ritalinic Acid	220.1	56.0	7.03
13. Methylphenidate	234.1	56.0	8.26
14. MDPV	276.1	126.0	8.67
15. Pyrovalerone	246.1	105.0	8.83

HPLC Conditions

Instrument: Thermo Scientific™ Dionex™ Ultimate™ 3000 LC system

Detector: TSQ Vantage™ tandem mass spectrometer

Ionization mode: ESI+

HPLC column: UCT Selectra® DA, 100 × 2.1 mm, 3 μm
(p/n: SLDA100ID21-3UM)

Guard column: UCT Selectra® DA, 10 × 2.0 mm, 3 μm,
(p/n: SLDAGDC21-3UM)

Guard column holder: p/n: SLDGRDHLDR

Column temp.: 40 °C

Injection volume: 10 μL

Flow rate: 300 μL/min

Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in MeOH

Gradient:

Time	%A	%B
0	98	2
1	65	35
5	65	35
7	0	100
10	0	100
10.2	98	2
15	98	2

EXTRACTION OF BENZODIAZEPINES FROM URINE, BLOOD, PLASMA/SERUM, TISSUE BY LC-MS/MS

Part

CSXCE106 CLEAN SCREEN XCEL® I 130 mg, 6 mL Tube
 SPHPHO6001-5 – Select pH Buffer Pouches 100 mM Phosphate pH 6.0
 SMTBSTFA-1-1 – SELECTRA-SIL® MTBSTFA w/ 1% TBDMCS
 SLDA100ID21-3UM – Selectra® DA HPLC Column 100 x 2.1mm 3 µm
 SLDAGDC21-3UM - SELECTRA® DA 3 µm Guard Cartridge
 SLGRDHLDR - Guard Cartridge Holder



SELECTRAZYME®
Beta-glucuronidase

1a. PREPARE SAMPLE (Blood, Plasma and Serum)

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards.
 Add 1-2 mL of blood, plasma/ serum, or 1 g (1:4) tissue homogenate.
 Mix/vortex and let stand for 5 minutes.
 Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex
 Sample pH should be 6.0 ± 0.5.
 Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
 Centrifuge for 10 minutes at 2000 rpm and discard pellet.

1b. PREPARE SAMPLE (Urine) ENZYME HYDROLYSIS OF GLUCURONIDES

To 1-2 mL of urine sample, add 1 mL of acetate buffer (pH 5.0) containing 5,000 units/mL of Selectrazyme® β-glucuronidase.
 Optionally, add 1 mL of acetate buffer and 25-50 µL of concentrated β-glucuronidase.
 Vortex and heat for 1-2 hours at 65°C.
 Allow sample to cool.
 Do not adjust pH~ sample is ready to be added to the extraction column.

2. APPLY SAMPLE

Load sample directly to column without any preconditioning.
 Pull sample through at a rate of 1-2 mL/ minute.
 Dry column thoroughly under full vacuum or positive pressure for 1 minute.

3. WASH

1 x 3mL 100 mM phosphate buffer (pH6).
 1 x 3 mL CH₂Cl₂
 Dry column thoroughly under full vacuum or positive pressure for a minimum of 5-10 minutes.

4. ELUTION

1 x 3 mL Ethyl Acetate:NH₄OH (98:2)
 Collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily.


5. DRY ELUTE

Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.

6. RECONSTITUTE / DERIVATIZE

- LC-MS/MS: Reconstitute sample in 100 µL of mobile phase Inject 10-20 µL.
- GC-MS: Dissolve residue in 50 µL of ACN and 50 µL MTBSTFA w/ 1%TBDMCS
 Overlay with N₂ and cap. Mix/vortex
 React 30 minutes at 70 °C; Cool and inject 1-2 µL

BENZODIAZEPINES



Analyte	MRM Transitions		Relative Retention Time (minutes)	% Recovery (Urine)	% Recovery (Blood)
	Q1	Q3			
1. 7-Aminoclonazepam	286.0	222.3	2.10	93%	95%
2. Midazolam	326.0	291.0	3.26	89%	89%
3. Lorazepam	321.0	303.3	4.73	93%	76%
4. Oxazepam	287.0	241.3	4.98	78%	96%
5. Clonazepam	316.1	270.2	5.05	78%	78%
6. Alpha-Hydroxy-Alprazolam	325.1	297.1	6.08	103%	82%
7. Nordiazepam	271.0	140.1	6.65	101%	80%
8. Temazepam	301.1	255.2	7.59	79%	68%
9. Alprazolam	309.1	205.3	8.91	90%	89%
10. Diazepam	285.1	193.2	10.98	100%	90%

HPLC Conditions

Instrument: Agilent 1200 Binary Pump SL

Detector: API 4000 Qtrap MS/MS

HPLC Column: SELECTRA® DA HPLC Column 100 x 2.1mm 3 µm

Polarity: Positive

Reconstitute: 100 µl

Injection Volume: 10 µL

Flow Rate: 0.3 mL/minute

Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in MeOH

Gradient:

Time	%A	%B
0.0	40	60
8.00	40	60
8.01	5	95
9.01	5	95
9.50	40	60
13.00	STOP	

BATH SALTS IN BLOOD, PLASMA/SERUM, URINE, OR TISSUE BY LC-MS/MS OR GC-MS CLEAN SCREEN XCEL® I EXTRACTION COLUMN

Part

CSXCE111 – CLEAN SCREEN XCEL® 130 mg, 1 mL Tube
 PFAA-0-1 – SELECTRA-SIL® PFAA
 SPFPOH-1 – SELECTRA-SIL® PFPOH
 SLDA100ID21-3UM – SELECTRA® DA HPLC Column 100 × 2.1 mm, 3 µm
 SLDAGDC21-3UM - SELECTRA® DA 3 µm Guard Cartridge
 SLGRDHLDR - Guard Cartridge Holder



SELECTRA-SIL®
Derivatizing Reagents

1. PREPARE SAMPLE

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards
 Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate
 Mix/vortex and let stand for 5 minutes
 Add 2 mL of 100 mM phosphate buffer (pH 6.0).
 Mix/vortex
 Sample pH should be 6.0 ± 0.5.
 Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
 Centrifuge for 10 minutes at 2000 rpm and discard pellet

2. APPLY SAMPLE

Load sample directly to column without any preconditioning.
 Pull sample through at a rate of 1-2 mL/ minute.
 Dry column thoroughly under full vacuum or positive pressure for 1 minute.

3. WASH

1 x 3 mL 98% Methanol: 2% Acetic Acid
 Dry column thoroughly under full vacuum or positive pressure for a minimum of 5 minutes.

4. ELUTION

1 x 3 mL CH₂Cl₂/ IPA/ NH₄OH (78:20:2)
 Collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily.
 Add IPA/ NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5. DRY ELUTE

Add 50 µL of 1% HCl in CH₃OH to each tube
 Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.


NOTE: A 1% HCl in CH₃OH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs.

6. RECONSTITUTE / DERIVATIZE

- LC-MS/MS: Reconstitute sample in 100 µL of mobile phase Inject 5 µL.
- GC-MS: Fluoroacylate with PFAA (PFAA)
 Add 50 µL PFAA. Over lay with N₂ and cap
 *Improved derivatization by addition of PFPOH
 React 20 minutes at 70 °C. Evaporate to dryness <40 °C
 Reconstitute with 100 µL Ethyl Acetate

NOTE: It is important to dry the column thoroughly to achieve the highest recovery of all compounds. Any residual moisture will slow down the drying of the elution solvents prior to derivatization for GC/MS analysis, if being used. Also, any residual moisture could reduce the reactivity of the derivatization agent resulting in low GC/MS sensitivity.

BATH SALTS



Analyte	MRM Transitions		Relative Retention Time (minutes)
	Q1	Q3	
1. Flephedrone	182.1	164.2	4.78
2. Methylone	208.1	160.1	5.26
3. Ethylone	222.0	131.0	5.84
4. Mephedrone	178.1	145.0	6.10
5. Butylone	222.0	131.0	6.36
6. Methcathinone	192.2	144.0	6.93
7. MDPV	276.2	126.1	8.67
8. Pyravalerone	246.2	105.2	8.83

HPLC Conditions

Instrument: Thermo Scientific™ Dionex™ Ultimate™ 3000 LC system

Detector: TSQ Vantage™ tandem mass spectrometer

Ionization mode: ESI+

HPLC column: Selectra® DA, 100 × 2.1 mm, 3 μm (p/n: SLDA100ID21-3UM)

Guard column: Selectra® DA, 10 × 2.0 mm, 3 μm, (p/n: SLDAGDC21-3UM)

Guard column holder: p/n: SLDGRDHLDLDR

Column temp.: 40 °C

Injection volume: 10 μL

Flow rate: 300 μL/min

Mobile Phase A: 0.1% formic acid in water

Mobile Phase B: 0.1% formic acid in MeOH

Gradient:

Time	%A	%B
0	98	2
1	65	35
5	65	35
7	0	100
10	0	100
10.2	98	2
15	98	2

EXTRACTION OF BASIC DRUGS AND METABOLITES FROM URINE/ BLOOD USING CLEAN SCREEN XCEL® I EXTRACTION COLUMN

Part #:

CSEXCE106 - CLEAN SCREEN XCEL® I 130 mg / 6 mL Tube
 SPACE5001-5 - SELECT pH Buffer Pouch, 100mM acetate pH 5.0
 SPHPHO6001-5 - SELECT pH Buffer Pouch, 100mM phosphate pH 6.0
 BETA-GLUC-10 - Selectrzyme® beta-glucuronidase, 10 mL

1a. Sample Preparation (blood)

To 1-2 mL blood add 2 ml of 100mM phosphate buffer pH 6.0. Add appropriate volume and concentration internal standards. Sample is ready to be added to extraction column.

1b. Sample Preparation (urine hydrolysis)

To 1-2 mL urine sample add 500 µL of 100mM acetate buffer pH 5.0 containing 5,000 units/mL β-glucuronidase. **Optionally**, add 500 µL of acetate buffer and 25 µL of concentrated β-glucuronidase. Add appropriate volume and concentration internal standards. Vortex and heat for 1-2 hours at 65 °C. Allow sample to cool. **Do not adjust pH- sample is ready to be added to extraction column.**

2. Applying Sample to Column

Load sample directly to column without any preconditioning. Pull sample through at a rate of 1-2 mL/ minute. Dry column thoroughly under vacuum (10 mm Hg) or positive pressure (~ 80-100 psi) for 1 minute.

3a. Wash (Blood Only)

Wash sample with 2 mL of 100mM phosphate buffer pH 6.0.

3b. Wash (Urine and Blood)

Wash sample with 1 mL of 2% glacial acetic acid/ 98 % methanol. Dry column thoroughly under vacuum (10 mm Hg) or positive pressure (~ 80-100 psi) for a minimum of 5 minutes.

NOTE 1: It is important to dry the column thoroughly to achieve the highest recovery of all compounds. Any residual moisture will slow down the drying of the elution solvents prior to derivatization for GC/MS analysis. Also, any residual moisture could reduce the reactivity of the derivatization agent resulting in low GC/MS sensitivity.

4. Elution

Elute samples with ~1-2 mL CH₂Cl₂ / IPA / ammonium hydroxide (78/20/2)
 Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.

(**Note:** For amphetamine group analysis, add 200 µL of 1% HCl in MeOH to eluate to minimize amphetamine group loss to volatilization.)

GC/MS Analysis

Derivatize compounds with appropriate derivatizing procedure or reconstitute in 100 µL ethyl acetate and inject 1-2 µL into the GC/ MS system for analysis.

LC/MS Analysis

Reconstitute in methanol or appropriate mobile phase.

Analyte	% Recovery	
	Urine	Blood
6-MAM	97	100
Amphetamine	112	104
Benzoylcegonine	54	58
Buprenorphine	92	92
Chlorpheniramine	85	94
Citalopram	99	95
Clonidine	85	98
Clozapine	98	94
Cocaine	103	96
Codeine	96	99
Cyclobenzaprine	102	92
Diphenhydramine	87	88
Doxepin	97	93
EDDP	85	88
Ephedrine	75	57
Fentanyl	88	89
Fluoxetine	97	92
Hydrocodone	98	94
Hydromorphone	91	91
Imipramine	97	97
MDA	95	84
MDEA	100	93
MDMA	99	90
Meperidine	82	95
Methadone	92	93
Methamphetamine	99	104
Morphine	92	85
Nalophrine	99	97
Naltrexone	99	90
Naxolone	99	94
Norfentanyl	84	85
Nortriptyline	79	89
Oxycodone	98	94
Oxymorphone	94	88
Paroxetine	90	92
Phencyclidine	90	93
Propoxyphene	90	95
Pseudoephedrine	79	56
Sertraline	75	70
Tramadol	89	81
Venlafaxine	102	93
Zolpidem	92	95

PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Terms of payment are net 30 days.

MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

SHIPMENTS

Normal processing is within 24 hours after receipt of an order. Unless special shipping requests have been made, our trained staff will send all orders by UPS Ground service. The appropriate shipping charges (freight & insurance costs) will be added to the invoice, unless otherwise instructed by the customer.

SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. These discounts apply to bonded phase extraction column purchases only. Please call a sales representative for more information on special pricing qualifications.

RETURN POLICY

Our Quality Manager will handle all returns. Before returning merchandise, please call to obtain a return authorization number from the quality manager. We will need to know the reason for the return, date of purchase, purchase order number and invoice number in order to issue a return authorization number. Return merchandise must be received before a credit can be issued. Returns will not be accepted after 90 days. A restocking fee of 25% of the price paid, or a minimum of \$25.00 (whichever is greater) will be charged on all returns.

WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period.

The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing. Secondly, the defective product must be returned within 30 days after proper approval from our Quality Manager. All returns must indicate the purchase order number, the lot number and the shipping date. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

Contact Us

Phone: 215.781.9255
800.385.3153
Fax: 215.785.1226

UCT, Inc.
2731 Bartram Rd.
Bristol, PA 19007

Email: info@unitedchem.com
Web: www.unitedchem.com

