


# ATS Series Auto-Sampler

## —Provide Efficient Purification

ATS Auto-Sampler is an efficient product for sample injection, which provides fast and reliable purification by connected with a prep LC system

- Zero Contamination  
Ensure complete cleaning with no residue
- Continuous Sampling  
Sample injection could be conducted up to a thousand times
- Excessive Large Sampling Volume  
1-25 ml for a single run
- Good Compatibility  
Compatible with the purification system from most of the manufactures in the market



AUTHORIZED DISTRIBUTOR



MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20  
e-mail: info@mz-at.de  
www.mz-at.de

# Still purify your sample by hand?

Free yourself by choosing a fully automatic purification kit!

Just build a connection between

**ATS Auto-Sampler----- Column Switcher----- CHEETAH™ System**



## Specifications

Sample channel	5
Cleaning channel	1
Sample loop	Standard: 2 mL
Single Sample Size	1-25 mL
Cleaning	Automatic/Manual
Language	Chinese/English
Temperature	20-35 °C
Size	361 x 208.5 x 344.4 mm (length x width x height )
Weight	11.4 kg
Power	1-120 VAC, 200-240 VAC, 50 Hz/60Hz

## Accessories

Cat. No	Discription	Quantity
DLH-1-1/16 (optional)	Sample loop, 1mL; 1/16 inch; stainless steel; Vaclo	1
DLH-2-1/16 (standard)	Sample loop, 2mL; 1/16 inch; stainless steel; Vaclo	1
DLH-5-1/16 (optional)	Sample loop, 5mL; 1/16 inch; stainless steel; Vaclo	1
DLH-10-1/16 (optional)	Sample loop, 10mL; 1/16 inch; stainless steel; Vaclo	1