

Conditions

Column: ACE Excel 2 SuperC18
Dimensions: 150 x 2.1 mm
Part Number: EXL-1011-1502U
Mobile Phase:
 A: 10 mM sodium dihydrogen phosphate pH 2.8 in H₂O
 B: 10 mM sodium dihydrogen phosphate pH 2.8 in H₂O/MeCN (20:80 v/v)
Gradient:

Time (mins)	%B
0	39.5
4	48.0

Flow Rate: 0.6 mL/min
Injection: 1 µL
Temperature: 50 °C
Detection: UV, 200 nm

Analytes

1. Rebaudioside D
2. Rebaudioside A
3. Stevioside
4. Rebaudioside F
5. Rebaudioside C
6. Dulcoside A
7. Rubusoside
8. Rebaudioside B
9. Steviolbioside

