





CHIRAL
TECHNOLOGIES INC

SUBSIDIARY OF DAICEL CHEMICAL INDUSTRIES, LTD.

CHIRALPAK® IC™

*Immobilized Polysaccharide
chiral stationary phases
and columns — the solution
you've been looking for.*





- *Authentic Daicel stationary phases and columns*
- *Unique chiral selector*
- *Outstanding solvent tolerance*
- *Solvent flexibility for the resolution of compounds with limited solubility*
- *Extended durability*
- *Excellent separation of enantiomers*



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

CHIRALPAK® IC™



For more information about all our proprietary columns and CSPs, call 1-800-6CHIRAL, or visit us online at www.chiraltech.com.

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
83020	CHIRALPAK IC	20	1 kg		CSP
83225	CHIRALPAK IC	20	4.6	250	Method Develop
83233	CHIRALPAK IC	20	10	100	SMB (set of 8)
83242	CHIRALPAK IC	20	21	50	Semi-Prep Guard
83245	CHIRALPAK IC	20	21	250	Semi-Prep
83256	CHIRALPAK IC	20	50	500	Preparative
83266	CHIRALPAK IC	20	100	500	Preparative
833C4	CHIRALPAK IC	5	0.3	150	Microflow
83311	CHIRALPAK IC	5	4.0	10	Guard Cartridges
83324	CHIRALPAK IC	5	4.6	150	Analytical
83325	CHIRALPAK IC	5	4.6	250	Analytical
83335	CHIRALPAK IC	5	10	250	Semi-Prep

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
83345	CHIRALPAK IC	5	20	250	Semi-Prep
83375	CHIRALPAK IC	5	30	250	Semi-Prep
83394	CHIRALPAK IC	5	2.1	150	Analytical
83422	CHIRALPAK IC	5	4.6	50	SFC Guard
83423	CHIRALPAK IC	5	4.6	100	SFC Screening
83435	CHIRALPAK IC	5	10	250	SFC Semi-Prep
83445	CHIRALPAK IC	5	21	250	SFC Semi-Prep
83455	CHIRALPAK IC	5	50	250	SFC Preparative
83475	CHIRALPAK IC	5	30	250	SFC Semi-Prep
83522	CHIRALPAK IC-3	3	4.6	50	Analytical
83523	CHIRALPAK IC-3	3	4.6	100	Analytical
83524	CHIRALPAK IC-3	3	4.6	150	Analytical