

ChiralTek Chiral Column Possible Compatible Table

ChiralTek	Daicel	Supelco	Phenomenex	Merck
AM-2/AMCE-2/AMCD-2/AM2AMx	IA		Lux i-Amylose-1	
CE-2/AMCE-2/CECD-2/CE2CEx	IB			
CE-4/AMCE-4/CE4CEx	IC		Lux i-Cellulose-5	
AM-4/AMCE-4/AM4AMx	IE			
AM-3/AMCE-3/AM3AMx	IF			
AM-5/AMCE-5/AM5AMx	IG			
AM-6/AMCE-6/AM6AMx	IH			
CE-7/ AMCE-7/CE7CEx	IJ			
AM-1				
AM-7				
AM-8				
AM-9				
AM-10				
AM-11				
AM-12				
CE-1				
CE-3				
CE-5				
CE-6				
CE-8				
CE-9				
CE-10				
CE-11				
CE-12				
CTCE-1	OC/OC-H			
CTCE-2	OD/OD-H		Lux Cellulose-1	
CTCE-7	OJ/OJ-H		Lux Cellulose-3	
CTCE-3	OZ/OZ-H		Lux Cellulose-2	
CTCE-9	OX/OX-H		Lux Cellulose-4	
CTCE-8	OG/OG-H			
CTAM-6	AS/AS-H			
CTAM-2	AD/AD-H		Lux Amylose-1	
CTAM-3	AZ/AZ-H			
CTAM-10	AY/AY-H		Lux Amylose-2	
CTCE-2R	OD-R/OD-RH			
CTCE-7R	OJ-R/OJ-RH			
CTCE-3R	OZ-R/OZ-RH			
CTCE-9R	OX-R/OX-RH			
CTAM-6R	AS-R/AS-RH			
CTAM-2R	AD-R/AD-RH			
CTAM-3R	AZ-R/AZ-RH			
CTAM-10R	AY-R/AY-RH			
CTCE-1R	OC-R/OC-RH			
CTCE-8R	OG-R/OG-RH			
CF6	CR(+)			
CF6	CR(-)			
β CD-0		Astec Cyclobond I 2000		
γ CD-0		Astec Cyclobond II		
β CD-2		Astec [®] CYCLOBOND I 2000 DMP		
β CD-1/ β CD-2/ β CD-3/ β CD-4/ β CD-5/ β CD-6/ β CD-7/ β CD-8/ β CD-9/ β CD-10/ β CD-11/ β CD-12/				ChiralDex

PS: x=1-12, this table only shows possible compatible chiral columns from ChiralTek and other suppliers.