


RiboSep RNA SemiPrep Column

P/N: RPC-99-3010

Processing Method: RNASep_30x100 (new)
 Application: Oligo_HC
 Sample Name: pUC18

Volume: 100.0 ul
 Column size: 30mm x 100mm
 Column S/N: RNA-07302020-1


UV @ 260nm



Flowrate: 3.0mL/min

Pump A Pressure Profile: Max: 1180.8 Min: 1081.2

Temperature Profile : Max: 50.0 Min: 49.9



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
 Barcelona-Allee 17 • D-55129 Mainz
 Tel +49 6131 880 96-0
 Fax +49 6131 880 96-20
 e-mail: info@mz-at.de
 www.mz-at.de

RT (min)	Name	k'	Asym	N (USP)	Res (USP)	Alpha	S/N	Noise (uV)
14.59	80 bp	11.58	---	67374	---	---	---	---
16.26	102 bp	13.01	1.41	56625	6.70	1.12	---	---
21.66	174 bp	17.68	2.13*	115609	20.48	1.36	---	---
25.12	257 bp	20.66	---	177587	14.03	1.17	---	---
25.46	267 bp	20.95	1.32	177033	1.42	1.01	---	---
26.30	298 bp	21.68	1.92*	158300	3.32	1.03	---	---
28.90	434 bp	23.91	---	190783	9.80	1.10	---	---
29.31	458 bp	24.27	1.48	188065	1.56	1.02	---	---
30.71	587 bp	25.47	2.23*	154942	4.78	1.05	---	---

