



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH  
Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0  
Fax +49 6131 880 96-20  
e-mail: info@mz-at.de  
www.mz-at.de

# Transgenomic RPsep Fluoro Silica Columns (pentafluorophenylpropyl phase)

Finally! An improved Silica-based line of Reverse Phase Columns for those who need increased selectivity and reproducibility in their assays. Transgenomic is introducing special reverse phase column types to enhance your everyday separation needs, as well as to provide more power towards R&D projects.

Let us help you choose the perfect column for your separation needs

## Benefits

- LC-MS Compatible
- Can be used with Reversed or Normal Phase and HILIC Mode
- Special Selectivity compared to Normal C18 phase

## Uses

- Pharmaceuticals, natural products, research & method development
- Strong interaction with Basic, Halogenated analytes and Aromatics such as Fluorines and Amines



**Beta Blockers, Fluorinated Corticosteroids, Nucleosides, Nucleotides, Purines, Pyrimidines, Taxanes and Tricyclic Antidepressants**

Cat No.	Sphere	Length	ID		Cat No.	Sphere	Length	ID
SIL-99-1101	3µm	50mm	2.1mm		SIL-99-1106	5µm	50mm	2.1mm
SIL-99-1102	3µm	100mm	2.1mm		SIL-99-1107	5µm	100mm	2.1mm
SIL-99-1103	3µm	50mm	4.6mm		SIL-99-1108	5µm	50mm	4.6mm
SIL-99-1104	3µm	100mm	4.6mm		SIL-99-1109	5µm	100mm	4.6mm

Custom column dimensions also available



Corporate Headquarters:  
Transgenomic, Inc.  
12325 Emmet Street, Omaha, NE 68164  
phone: +1(402) 452.5400 Fax: (402) 452 5401  
email: [info@transgenomic.com](mailto:info@transgenomic.com)  
On the web at: [www.transgenomic.com](http://www.transgenomic.com)

# NEW

# Powerful