

PART OF THE **DURAN**® BOTTLE SYSTEM

ISOLATED


DURAN® DOUBLE WALLED BOTTLE



**DURAN
WHEATON
KIMBLE**


Excellence in your hands

DURAN® GLS 80® DOUBLE WALLED BOTTLE WIDE MOUTH



The DURAN® double walled, wide mouth bottles GLS 80® incorporate an integral jacket that isolates the contents from the external environment. Heated or cooled liquids can be circulated through the jacket to control the temperature within the screw topped DURAN® bottle. The DURAN® double walled bottles offer a sealable and more flexible alternative to open topped jacketed beakers.

Cat. No.	Capacity (mL)	Thread	d (OD) (mm)	h (mm)	Pack Unit
with screw cap and pouring ring from PP (blue)					
24 256 44 51	500	80	110	175	1
24 256 54 56	1 000	80	110	275	1

**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz
 Tel +49 6131 880 96-0, Fax +49 6131 880 96-20
 e-mail: info@mz-at.de, www.mz-at.de



> Find your nearest **distributor** on our global network:
www.DWK-LifeSciences.com/DURAN/distributors

