


PART OF THE **DURAN**® BOTTLE SYSTEM

PURE

DURAN® PURE BOTTLE



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de




**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® PURE BOTTLE



with DIN thread



High-quality Type I borosilicate glass bottles GL 25 / GL 45 from 25 ml – 20 l, developed for the needs of the pharmaceutical industry. DURAN® PURE Bottles are closed with a protective cover to prevent contamination. DWK Life Sciences provides CoAs with each delivery, assures 100 % traceability and change control management. DURAN® PURE Bottles are available with the documentation required by regulatory authorities and meet the requirements of the international pharmacopeias.


Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs.

Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
81 801 14 01	25	25	36	70	10
81 801 24 06	100	45	56	100	10
81 801 36 08	250	45	70	138	10
81 801 44 07	500	45	86	176	10
81 801 54 03	1 000	45	101	225	10
81 801 63 05	2 000	45	136	260	10
81 801 73 01	5 000	45	182	330	1
81 801 86 06	10 000	45	227	410	1
81 801 91 05	20 000	45	288	505	1

DURAN® PURE BOTTLE AMBER



with DIN thread



High-quality Type I borosilicate glass bottles GL 25 / GL 45, amber coloured from 25 ml – 20 l, developed for the needs of the pharmaceutical industry. DURAN® PURE Bottles are closed with a protective cover to prevent contamination. Very uniform, durable and chemically resistant amber colouring due to the use of innovative technology offers protection for light sensitive compounds. DWK Life Sciences provides CoAs with each delivery, assures 100 % traceability and change control management. DURAN® PURE Bottles are available with the documentation required by regulatory authorities and meet the requirements of the international pharmacopeias.

Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of API's.

Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
81 806 14 06	25	25	36	70	10
81 806 24 02	100	45	56	100	10
81 806 36 04	250	45	70	138	10
81 806 44 03	500	45	86	176	10
81 806 54 08	1 000	45	101	225	10
81 806 63 01	2 000	45	136	260	10
81 806 73 06	5 000	45	182	330	1
81 806 86 02	10 000	45	227	410	1
81 806 91 01	20 000	45	288	505	1




> Find your nearest **distributor** on our global network:
www.DWK-LifeSciences.com/DURAN/distributors



DURAN® PURE GLS 80® BOTTLE



with GLS 80® thread



High-quality Type I borosilicate glass bottles GLS 80® from 500 ml – 20 l, developed for the needs of the pharmaceutical industry. DURAN® PURE Bottles are closed with a protective cover to prevent contamination. The 80 mm wide outer diameter of the bottle mouth allows easy filling and emptying of bottles. DWK Life Sciences provides CoAs with each delivery, assures 100% traceability and change control management. DURAN® PURE Bottles are available with the documentation required by regulatory authorities and meet the requirements of the international pharmacopeias.


Typical applications: manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs.

Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
81 860 44 03	500	GLS 80®	101	148	10
81 860 54 08	1 000	GLS 80®	101	218	10
81 860 63 01	2 000	GLS 80®	136	248	10
81 860 73 06	5 000	GLS 80®	182	310	1
81 860 86 02	10 000	GLS 80®	227	385	1
81 860 91 01	20 000	GLS 80®	288	480	1

DURAN® PURE GLS 80® BOTTLE AMBER



with GLS 80® thread



High-quality Type I borosilicate glass bottles GLS 80®, amber coloured from 500 ml – 20 l, developed for the needs of the pharmaceutical industry. DURAN® PURE Bottles are closed with a protective cover to prevent contamination. The 80 mm wide outer diameter of the bottle mouth allows easy filling and emptying of bottles. Very uniform, durable and chemically resistant amber colouring due to the use of innovative technology offers protection for light sensitive compounds. DWK Life Sciences provides CoAs with each delivery, assures 100% traceability and change control management. DURAN® PURE Bottles are available with the documentation required by regulatory authorities and meet the requirements of the international pharmacopeias.

Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs


Cat. No.	Capacity (mL)	DIN Thread (GL)	d (OD) (mm)	h (mm)	Pack Unit
81 866 44 09	500	GLS 80®	101	148	10
81 866 54 05	1 000	GLS 80®	101	218	10
81 866 63 07	2 000	GLS 80®	136	248	10
81 866 73 03	5 000	GLS 80®	182	310	1
81 866 86 08	10 000	GLS 80®	227	385	1
81 866 91 07	20 000	GLS 80®	288	480	1



DURAN® PURE PREMIUM LIP SEAL CAP



PFA



High-quality plastic screw cap with GL 45 thread size, developed for the needs of the pharmaceutical industry. The DURAN® PURE Premium Lip Seal Cap is manufactured in ISO Class 7 clean rooms and made from PFA with lip seal design. Designed for use with DURAN® PURE GL 45 Bottles. The screw cap is compatible with autoclaving and other common sterilization methods. DWK Life Sciences provides CoAs with each delivery, assures 100% traceability and change control management. The DURAN® PURE Premium Lip Seal Cap is available with the documentation required by regulatory authorities and meets the requirements of the international pharmacopeias.


Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs.

Cat. No.	DIN Thread (GL)	d (OD) (mm)	h (mm)	Colour	Pack Unit
89 886 28 02	45	53	27	colourless	5

DURAN® PURE PREMIUM SCREW CAP



PFA, with PTFE coated silicone seal



High-quality plastic screw caps with GL 25 / GL 45 thread sizes, developed for the needs of the pharmaceutical industry. DURAN® PURE Premium Screw Caps are manufactured in ISO Class 7 clean rooms and made from PFA with highly inert PTFE coated silicone seal. Designed for use with DURAN® PURE GL 25 / 45 Bottles. The screw caps are compatible with autoclaving and other common sterilization methods. DWK Life Sciences provides CoAs with each delivery, assures 100% traceability and change control management. DURAN® PURE Premium Screw Caps are available with the documentation required by regulatory authorities and meet the requirements of the international pharmacopeias.


Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs

Cat. No.	DIN Thread (GL)	d (OD) (mm)	h (mm)	Colour	Pack Unit
89 108 13 06	25	32	22	colourless	5
89 108 86 79	45	51	26	colourless	5

DURAN® PURE GLS 80® SCREW CAP



PSU, with PTFE coated silicone seal



High-quality plastic screw cap with GLS 80® thread size, developed for the needs of the pharmaceutical industry. The DURAN® PURE GLS 80® Screw Cap is made from PSU with highly inert PTFE coated silicone seal. Designed for use with DURAN® PURE GLS 80® Bottles. The screw cap is compatible with autoclaving and other common sterilization methods. DWK Life Sciences provides CoAs with each delivery, assures 100% traceability and change control management. The DURAN® PURE GLS 80® Screw Cap is available with the documentation required by regulatory authorities and meets the requirements of the international pharmacopeias.

Typical applications: Manufacturing of pharmaceutical products, storage of bulk drug intermediates and ingredients, transportation of APIs.

Cat. No.	DIN Thread (GL)	d (OD) (mm)	h (mm)	Colour	Pack Unit
89 116 91 09	GLS 80®	88	40	white	5

