


PART OF THE DURAN® BOTTLE SYSTEM

SMART

DURAN® YOUTILITY® BOTTLE




**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® YOUTILITY LABORATORY BOTTLE



GL 45




The specially shaped gripping zones on both sides of the bottle enable easier and safer handling. With the new DURAN® YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories. The slimmer DURAN® YOUTILITY bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators. A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive YOUTILITY bottle labels. Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes. Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (PP) and a GL 45 cap (PP).

| Cat. No. | Capacity (mL) | DIN Thread (GL) | d (OD) (mm) | h (mm) | Pack Unit |
|--|---------------|-----------------|-------------|--------|-----------|
| with screw-cap and pouring ring from PP | | | | | |
| 21 881 28 54 | 125 | 45 | 55 | 124 | 4 |
| 21 881 36 53 | 250 | 45 | 66 | 158 | 4 |
| 21 881 44 52 | 500 | 45 | 78 | 193 | 4 |
| 21 881 54 57 | 1 000 | 45 | 93 | 253 | 4 |

DURAN® YOUTILITY LABORATORY BOTTLE AMBER



GL 45, USP <660> and USP <671> (Spectral Transmission) compliant



The specially shaped gripping zones on both sides of the bottle enable easier and safer handling. With the new DURAN® YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories. The slimmer DURAN® YOUTILITY bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators. A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive YOUTILITY bottle labels. Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes. The glass bottle body is moulded from the tried and tested DURAN® borosilicate 3.3 pharmacopoeial Type 1 neutral glass. DURAN® glass offers very good chemical resistance and high temperature resistance. Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (PP) and a GL 45 cap (PP).

| Cat. No. | Capacity (mL) | DIN Thread (GL) | d (OD) (mm) | h (mm) | Pack Unit |
|--|---------------|-----------------|-------------|--------|-----------|
| with screw-cap and pouring ring from PP | | | | | |
| 21 886 28 59 | 125 | 45 | 55 | 124 | 4 |
| 21 886 36 58 | 250 | 45 | 66 | 158 | 4 |
| 21 886 44 57 | 500 | 45 | 78 | 193 | 4 |
| 21 886 54 53 | 1 000 | 45 | 93 | 253 | 4 |




> Find your nearest **distributor** on our global network:
www.DWK-LifeSciences.com/DURAN/distributors



DURAN® YOUTILITY SCREW CAP FROM PP



GL 45




The DURAN® YOUTILITY Screw Cap GL 45 is manufactured from a food-grade polypropylene (PP). Ergonomically shaped screw cap with optimised grooves and ridges for a more efficient and easier tightening or removal. The faster opening and closing thread of the YOUTILITY screw cap is fully compatible with DIN GL 45 bottle threads. The optimised cap sealing system ensures a liquid tight seal. A pre-defined labelling area on the cap is compatible with the dedicated DURAN® self-adhesive YOUTILITY labels.

| Cat. No. | DIN Thread (GL) | d (OD) (mm) | h (mm) | Colour | Pack Unit |
|---------------------|-----------------|-------------|--------|--------|-----------|
| screw cap | | | | | |
| 29 229 28 02 | 45 | 53 | 25 | cyan | 10 |
| pouring ring | | | | | |
| 29 241 28 08 | 45 | | 4 | cyan | 16 |

DURAN® YOUTILITY BOTTLE TAG



GL 45, from silicone





The silicone GL 45 Bottle Tags can be securely attached around the neck of the YOUTILITY bottle for easy customisation and unambiguous bottle identification. The elastic Bottle Tags will fit around any DIN GL 45 bottle neck.

| Cat. No. | Thread | Colour | Pack Unit |
|--------------|--------|--------------------------------|-----------|
| 29 243 29 04 | 45 | eight colours, two pieces each | 16 |
| 29 243 28 18 | 45 | red | 20 |
| 29 243 28 26 | 45 | orange | 20 |
| 29 243 28 34 | 45 | yellow | 20 |
| 29 243 28 42 | 45 | green | 20 |
| 29 243 28 59 | 45 | blue | 20 |
| 29 243 28 67 | 45 | purple | 20 |
| 29 243 28 75 | 45 | black | 20 |
| 29 243 28 83 | 45 | white | 20 |




DURAN® YOUTILITY LABELS

Printable self-adhesive



The robust DURAN® YOUTILITY identification labels are suitable for use with the YOUTILITY bottle, plus many other types of laboratory glassware, such as GL 45 bottles, beakers and conical flasks. The tear resistant, white polyester labels use a high performance, peelable adhesive. An easy-to-use label creator web app. is available (www.duranlabels.com) to design, and print your own label designs. Ideal for use in fridges, freezers, autoclaves, incubators, and water baths, without the risk of the labels falling off. Can be printed with office printers, copiers, or marked with technical lab pens. Chemically resistant to many typical laboratory chemicals, and solvents. Wide -40 to +150 °C thermal performance range.

| Cat. No. | Description | Pack Unit |
|--------------|------------------|----------------|
| 29 401 02 03 | white, polyester | 1 x 100 labels |



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de



> Find your nearest **distributor** on our global network:
www.DWK-LifeSciences.com/DURAN/distributors