

# GraceSmart™ Columns



## High Quality HPLC Phases at Exceptional Value

- High purity phases
- Efficient and reproducible separations
- General use selectivity
- Exceptional value


GraceSmart™ HPLC phases use high-purity silica and are monomerically bonded with uniform coverage. This translates into symmetrical peaks for acids/bases, and predictable reversed-phase selectivity. Whether routine analysis or new method development, use GraceSmart™ columns to get premium performance at exceptional value.

GraceSmart™ Phase Specifications				
Phase	Encapping	Surface Area	Pore Size	Particle Size
C18	TMS, non-polar	220m <sup>2</sup> /g	120Å	3, 5µm

hplc columns | small molecule

## Reproducible Methods Start with Reproducible Columns

Our tightly controlled silica synthesis and bonding keep capacity factor and selectivity variations to a minimum. The advanced packing methods deliver consistently high column-to-column efficiency. Every column is individually tested and shipped with a quality assurance chromatogram.




Consistent batch-to-batch capacity for acidic, basic, and neutral compounds.

### GraceSmart™ RP18 Columns


Packing	Format	i.d. x Length	Part No.
C18, 5µm	Analytical	4.6 x 150mm	5138812
	Analytical	4.6 x 250mm	5138810
	Analytical	2.1 x 150mm	5138811
C18, 3µm	Analytical	2.1 x 250mm	5138813
	Analytical	4.6 x 150mm	5141752
	Analytical	4.6 x 100mm	5141753
	Analytical	2.1 x 150mm	5141754
	Analytical	2.1 x 100mm	5141755
	Analytical	2.1 x 50mm	5141756
Analytical	4.6 x 50mm	5141811	

## Expect Competitive Performance

Six sigma techniques and lean manufacturing are key to our cost competitive position, without compromising quality or performance. GraceSmart™ columns show similar or better efficiency and asymmetries for challenging base and acid components in comparison to industry leading columns.



GraceSmart™ columns compare favorably to leading competitor columns.








AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH  
 Barcelona-Allee 17 • D-55129 Mainz  
 Tel +49 6131 880 96-0  
 Fax +49 6131 880 96-20  
 e-mail: info@mz-at.de  
 www.mz-at.de



# GraceSmart™ Columns

hplc columns | small molecule



## more applications

To view our complete searchable chromatogram database visit [www.discoverysciences.com/chromdb/](http://www.discoverysciences.com/chromdb/)

