


# Periodic Table of The Analytical Fittings



**AUTHORIZED DISTRIBUTOR**  
 MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz  
 Tel +49 6131 880 96-0, Fax +49 6131 880 96-20  
 e-mail: info@mz-at.de, www.mz-at.de



LEGEND		FITTINGS		TOOLS		ICONS	
1 Pressure Rating	2 Tubing Size	3 Fitting / Female Part Number	4 Torque Required	5 Head Style	6 Fitting / Ferrule Material	7 Tool Category	8 Tool Part Number
9 Torque Required	10 Head Style						

<b>3750 PSI 1/16</b> M-644-03 / M-650 2.3 in-lbs Headless Micro Knurl PEEK	<b>15000 PSI 1/32</b> UH-904 3.4 in-lbs Headless Knurl PEEK Blend	<b>3750 PSI 1/16</b> M-660 / M-650 2.3 in-lbs Headless Micro Knurl PEEK	<b>4000 PSI .025</b> F-125 3.4 in-lbs Standard Micro Knurl PEEK	<b>2500 PSI 1/8</b> P-331 / P-350 3.4 in-lbs Standard Knurl PEEK	<b>TOOL N-290</b> 3.4 in-lbs Micro Headless Knurl	<b>5000 PSI 360 µm</b> F-124S 3.4 in-lbs Standard Micro Knurl PEEK	<b>5000 PSI 1/32</b> F-126S 3.4 in-lbs Standard Micro Knurl PEEK	<b>5000 PSI 1/32</b> F-126H 3.4 in-lbs Micro Knurl PEEK	<b>5000 PSI 360 µm</b> F-124H 3.4 in-lbs Headless Micro Knurl PEEK	<b>15000 PSI 1/32</b> PK-126 3.4 in-lbs Standard Micro Knurl PEEK Blend	<b>15000 PSI 1/32</b> PK-126H 3.4 in-lbs Headless Micro Knurl PEEK Blend	<b>15000 PSI 360 µm</b> PK-124 3.4 in-lbs Standard Micro Knurl PEEK Blend	<b>30000 PSI 1/32</b> VHP-700 20 in-lbs 4 mm Hex SS	<b>25000 PSI 1/32</b> VHP-920 3 in-lbs 4 mm Hex SS / PEEK Blend	<b>30000 PSI 1/32</b> VHP-900 / VHP-900-01 20 in-lbs 4 mm Hex SS
--	---	---	---	--	--	--	--	---	--	---	--	---	---	---	--

<b>4000 PSI 1/16</b> F-100 3.4 in-lbs Diamond Knurl PCTFE	<b>5000 PSI 1/16</b> F-130 3.4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/32</b> F-354 / LT-132 4 in-lbs Flush Nut SS / PEEK	<b>5000 PSI 1/32</b> LT-110 / LT-132 4 in-lbs Standard Knurl PEEK	<b>6000 PSI 190 µm</b> F-200 / F-148 3.4 in-lbs Single Wing Delrin / PCTFE	<b>6000 PSI 360 µm</b> F-200 / F-151 3.4 in-lbs Single Wing Delrin / PCTFE	<b>6000 PSI 1/32</b> F-200 / F-113 3.4 in-lbs Single Wing Delrin / PEEK	<b>6000 PSI 1/16</b> F-200 / F-142 3.4 in-lbs Single Wing Delrin / PEEK	<b>9000 PSI 1/16</b> F-195 / F-192 4 in-lbs Standard Knurl PEEK	<b>9000 PSI 1/16</b> F-198 / F-192 4 in-lbs Headless Knurl PEEK	<b>9000 PSI 1/16</b> F-287 / F-192 4 in-lbs Standard Knurl PEEK / PEEK Blend	<b>11000 PSI 1/16</b> VHP-3200 3.4 in-lbs Diamond Knurl SS / PEEK Blend	<b>12000 PSI 1/16</b> PK-120BLK 8 in-lbs Standard Knurl PEEK Blend	<b>20000 PSI 1/16</b> U-400 / U-401 20 in-lbs 1/4" Hex SS	<b>20000 PSI 1/16</b> U-420 / U-401 20 in-lbs 1/4" Hex SS	<b>25000 PSI 1/16</b> VHP-320 10 in-lbs 1/4" Hex SS / PEEK Blend	<b>30000 PSI 1/16</b> VHP-200 / VHP-200-01 20 in-lbs 1/4" Hex SS
---	---	--	---	--	--	---	---	---	---	--	---	--	---	---	--	--

<b>5000 PSI 1/16</b> F-120 3.4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/16</b> F-127 3.4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/16</b> F-354 / LT-100 4 in-lbs Flush Nut SS / PEEK	<b>5000 PSI 1/16</b> LT-110 / LT-100 4 in-lbs Standard Knurl PEEK	<b>6000 PSI 190 µm</b> F-300 / F-148 3.4 in-lbs Double Wing PEEK / PCTFE	<b>6000 PSI 360 µm</b> F-300 / F-151 3.4 in-lbs Double Wing PEEK / PCTFE	<b>6000 PSI 1/32</b> F-300 / F-113 3.4 in-lbs Double Wing PEEK	<b>6000 PSI 1/16</b> F-300 / F-142 3.4 in-lbs Double Wing PEEK	<b>9000 PSI 1/16</b> F-196 / F-192 4 in-lbs Standard Knurl PEEK	<b>9000 PSI 1/16</b> F-350 / F-192 4 in-lbs Flush Nut SS / PEEK Blend	<b>9000 PSI 1/16</b> F-284 / F-192 4 in-lbs Headless Knurl PEEK / PEEK Blend	<b>10000 PSI 1/16</b> LT-110 / LT-135 7 in-lbs Standard Knurl PEEK / PEEK Blend	<b>16500 PSI 1/32</b> PK-110 / PK-132 8 in-lbs Standard Knurl PEEK Blend	<b>16500 PSI 1/16</b> PK-110 / PK-100 8 in-lbs Standard Knurl PEEK Blend	<b>20000 PSI 1/16</b> U-320 / U-321 20 in-lbs 1/4" Hex SS	<b>25000 PSI 1/16</b> VHP-325 10 in-lbs 1/4" Hex SS / PEEK Blend	<b>30000 PSI 1/16</b> VHP-205 / VHP-200-01 20 in-lbs 1/4" Hex SS
---	---	--	---	--	--	--	--	---	---	--	---	--	--	---	--	--

<b>1000 PSI 1/16</b> P-249 2.3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> P-349 2.3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> P-331 / P-359 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> F-356 / P-259 3.4 in-lbs Flush Nut SS / ETFE	<b>1350 PSI 1/16</b> P-255 / P-259 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>1500 PSI 1/8</b> F-156 / P-360 3.4 in-lbs Female Knurl PEEK	<b>1500 PSI 1/8</b> P-331 / P-360 3.4 in-lbs Standard Knurl PEEK	<b>1850 PSI 1/16</b> F-356 / P-260 3.4 in-lbs Flush Nut SS / PEEK	<b>1850 PSI 1/16</b> P-255 / P-260 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/32</b> F-356 / P-248 3.4 in-lbs Flush Nut SS / ETFE	<b>2500 PSI 1/32</b> P-255 / P-248 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>2500 PSI 1/16</b> F-356 / P-250 3.4 in-lbs Flush Nut SS / PEEK	<b>2500 PSI 1/16</b> P-255 / P-250 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/8</b> F-156 / P-350 3.4 in-lbs Female Knurl PEEK	<b>2500 PSI 1/8</b> P-331 / P-350 3.4 in-lbs Standard Knurl PEEK	<b>4500 PSI 1/8</b> F-364 / LT-200 4 in-lbs Flush Nut SS / PEEK	<b>20000 PSI 1/8</b> C-235 / C-236 40 in-lbs Female Hex SS
---	--	---	---	--	--	--	---	---	---	--	---	---	--	--	---	--

<b>1000 PSI 1/8</b> P-100 3.4 in-lbs Diamond Knurl PCTFE	<b>1000 PSI 1/8</b> F-364 / P-359 3.4 in-lbs Flush Nut SS / ETFE	<b>1000 PSI 1/8</b> P-336 / P-359 3.4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> LT-115 / P-259 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> P-287 / P-259 3.4 in-lbs Headless Knurl PPS / ETFE	<b>1500 PSI 1/8</b> P-381 / P-360 3.4 in-lbs Standard Knurl PPS / PEEK	<b>1500 PSI 1/8</b> F-364 / P-360 3.4 in-lbs Flush Nut SS / PEEK	<b>1500 PSI 1/8</b> P-336 / P-360 3.4 in-lbs Headless Knurl PEEK	<b>1850 PSI 1/16</b> LT-115 / P-260 3.4 in-lbs Standard Knurl PEEK	<b>1850 PSI 1/16</b> P-287 / P-260 3.4 in-lbs Headless Knurl PPS / PEEK	<b>2500 PSI 1/32</b> LT-115 / P-248 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>2500 PSI 1/32</b> P-287 / P-248 3.4 in-lbs Headless Knurl PPS / ETFE	<b>2500 PSI 1/16</b> LT-115 / P-250 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/16</b> P-287 / P-250 3.4 in-lbs Headless Knurl PPS / PEEK	<b>2500 PSI 1/8</b> F-364 / P-350 3.4 in-lbs Flush Nut SS / PEEK	<b>2500 PSI 1/8</b> P-336 / P-350 3.4 in-lbs Headless Knurl PEEK	<b>4500 PSI 1/8</b> C-235 / LT-200 40 in-lbs 5/16" Hex SS / PEEK
--	--	---	---	---	--	--	--	--	---	---	---	--	---	--	--	--

<b>1000 PSI 1/8</b> F-156 / P-359 3.4 in-lbs Female Knurl PEEK / ETFE	<b>1000 PSI 1/8</b> LT-215 / P-359 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>1000 PSI 1/8</b> P-387 / P-359 3.4 in-lbs Headless Knurl PPS / ETFE	<b>1350 PSI 1/16</b> P-232 / P-259 3.4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> P-420 / P-259 3.4 in-lbs Female Knurl PEEK / ETFE	<b>1500 PSI 1/8</b> LT-215 / P-360 3.4 in-lbs Standard Knurl PEEK	<b>1500 PSI 1/8</b> P-387 / P-360 3.4 in-lbs Headless Knurl PPS / PEEK	<b>1850 PSI 1/16</b> P-232 / P-260 3.4 in-lbs Standard Knurl PEEK	<b>1850 PSI 1/16</b> P-420 / P-260 3.4 in-lbs Female Knurl PEEK	<b>2500 PSI 1/32</b> P-232 / P-248 3.4 in-lbs Headless Knurl PEEK / ETFE	<b>2500 PSI 1/32</b> P-420 / P-248 3.4 in-lbs Female Knurl PEEK / ETFE	<b>2500 PSI 1/16</b> P-232 / P-250 3.4 in-lbs Headless Knurl PEEK	<b>2500 PSI 1/16</b> P-420 / P-250 3.4 in-lbs Female Knurl PEEK	<b>2500 PSI 1/8</b> LT-215 / P-350 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/8</b> P-387 / P-350 3.4 in-lbs Standard Knurl PPS / PEEK	<b>4500 PSI 1/8</b> LT-215 / LT-200 4 in-lbs Standard Knurl PEEK	<b>4500 PSI 1/8</b> LT-210 / LT-200 4 in-lbs Double Wing PEEK
---	--	--	--	--	---	--	---	---	--	--	---	---	---	--	--	---

<b>1000 PSI 1/16</b> P-229 2.3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> P-329 2.3 in-lbs Standard Knurl PEEK	<b>1000 PSI 2.5 mm</b> P-317 / P-366 3.4 in-lbs Standard Knurl PEEK	<b>1350 PSI 1/16</b> P-213 / P-259 3.4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> P-219 / P-259 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> P-217 / P-259 3.4 in-lbs Standard Knurl PEEK / ETFE	<b>2500 PSI 1/32</b> P-213 / P-248 3.4 in-lbs Headless Knurl PEEK / ETFE	<b>2500 PSI 1/32</b> P-219 / P-248 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/16</b> P-213 / P-250 3.4 in-lbs Headless Knurl PEEK	<b>2500 PSI 1/8</b> P-219 / P-250 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1/8</b> P-317 / P-352 3.4 in-lbs Standard Knurl PPS / PEEK	<b>2500 PSI 1/8</b> P-319 / P-352 3.4 in-lbs Standard Knurl PEEK	<b>2500 PSI 1.8 mm</b> P-317 / P-355 3.4 in-lbs Standard Knurl PPS / PCTFE	<b>2500 PSI 1.8 mm</b> P-319 / P-355 3.4 in-lbs Standard Knurl PEEK / PCTFE	<b>5000 PSI 2.0 mm</b> P-357 / P-357-2 3.4 in-lbs Standard Knurl PEEK	<b>9000 PSI 1/16</b> F-197 / F-192 4 in-lbs Standard Knurl PEEK	<b>20000 PSI 1/16</b> U-450 / U-401 20 in-lbs 5/16" Hex SS
---	--	---	--	--	--	--	---	---	--	--	--	--	---	---	---	--

For more information, please visit [www.kinesis-group.com](http://www.kinesis-group.com) or contact your local Kinesis office:

- UK**  
Kinesis Ltd  
+44 (0)1480 212122  
9 Orion Court  
Ambuscade Road  
Colmworth Business Park  
St Neots, Cambridgeshire  
PE19 8YX, UK
- USA**  
Kinesis Inc  
+1 (856) 719-8787  
701 South Route 73  
Suite E  
West Berlin  
NJ 08091  
USA
- Germany**  
Kinesis GmbH (Abimed)  
+49 (0)2173 89 05 0  
Raiffeisenstr. 3  
40764 Langenfeld  
Deutschland
- Australia**  
Kinesis Pty  
+61 (0)7 3829 3996  
Unit 4, Jadsaw Court  
5-11 Daintree Drive  
Redlands Business Park  
Redland Bay, Qld 4165  
Australia
- India**  
Kinesis India Pvt Ltd  
+91 22 661 41 668  
ANM House, Plot No. A-141,  
Road No. 23  
Wagle Industrial Area  
Thane (W) 400604  
India

**MZ ANALYSENTECHNIK**  
 MZ-Analysentechnik GmbH  
 Wöhlerstraße 2-6 • D-55120 Mainz  
 Tel +49 6131 880 66 19  
 Fax +49 6131 880 66 20  
 e-mail: info@mz-at.de  
 www.mz-at.de

HAND TIGHTENING/EXTENDER TOOLS								PRE-SET TORQUE TOOLS				
<b>F-345</b> 3.4 in-lbs Flush Nut 10-32	<b>F-346</b> 3.4 in-lbs Flush Nut 1/4-28	<b>N-290</b> 3.4 in-lbs Micro Headless Knurl	<b>P-277</b> 3.4 in-lbs Standard Micro Knurl	<b>P-297</b> 3.4 in-lbs Headless Knurl	<b>P-299</b> 3.4 in-lbs Standard Knurl	<b>P-399</b> 3.4 in-lbs Standard Micro Knurl	<b>VHP-1001</b> 10 in-lbs 1/4" Hex	<b>VHP-9000</b> 2-10 in-lbs 4 mm Hex	<b>VHP-1000</b> 10 in-lbs 1/4" Hex	<b>VHP-2000</b> 14 in-lbs 1/4" Hex	<b>VHP-3000</b> 10 in-lbs 8 mm Hex	<b>P-1000</b> 4 in-lbs Standard Knurl

ADJUSTABLE TORQUE DRIVER AND TORQUE DRIVER EXTENDER TOOLS							
<b>VHP-4000</b> 2-10 in-lbs 1/4" Drive	<b>N-291</b> 2-10 in-lbs Micro Headless Knurl	<b>P-268</b> 2-10 in-lbs 1/4" Hex	<b>P-278</b> 2-10 in-lbs Female Nut Knurl	<b>P-279</b> 2-10 in-lbs Micro Knurl	<b>P-291</b> 2-10 in-lbs Standard Knurl	<b>P-292</b> 2-10 in-lbs Headless Knurl	<b>F-347</b> 2-10 in-lbs Flush Nut 1/4-28



**Fluidics for Life**  
 Fluidic Connections | Valves | Pumps | Degassers | Column Hardware | Manifolds