

## PuriFlash Dry Load Columns

EMPTY OR PARTIALLY FILLED WITH SILICA OR CELITE

- Self-contained, no adapters or special fittings
- Hand tight to 15 bar / 215 PSI
- Compatible with all Flash systems

Based on a proven design, the PuriFlash Dry Load columns are simple to load and the luer fittings provide a universal fit for all Flash systems.

The empty PuriFlash Dry Load columns allow you to add your own silica or celite. The HC Series is supplied with high capacity silica for maximum adsorption. The Celite Series is supplied with cellite and recommended for non-adsorptive loading.



Part number	Description	Price	Unit
PF-DLE/4G	PuriFlash Dry Load Column, Empty, 7ml, 12 x 65mm (FO8880)	\$80	20/PK
PF-DLE/12G	PuriFlash Dry Load Column, Empty, 30ml, 21 x 80mm (FO8890)	\$105	20/PK
PF-DLE/25G	PuriFlash Dry Load Column, Empty, 45ml, 21 x 130mm (FO8900)	\$135	20/PK
PF-DLE/40G	PuriFlash Dry Load Column, Empty, 75ml, 27 x 130mm (FL9521)	\$165	20/PK
PF-DLE/60G	PuriFlash Dry Load Column, Empty 115 ml, 37 x 110 mm	\$165	10/PK
PF-DLE/80G	PuriFlash Dry Load Column, Empty, 145ml, 31 x 193mm	\$95	5/PK
PF-DLE/100G	PuriFlash Dry Load Column, Empty 235 ml, 62 x 78 mm	\$235	5/PK
PF-DLE/120G	PuriFlash Dry Load Column, Empty, 215ml, 36 x 212mm	\$100	5/PK
PF-DLE/200G	PuriFlash Dry Load Column, Empty, 375ml, 60 x 133mm	\$125	5/PK
PF-DLE/300G	PuriFlash Dry Load Column, Empty, 575ml, 60 x 205mm	\$185	5/PK
PF-DLE/KA	PuriFlash Dry Load Column, Empty, 2 Each: 4g, 12g, 25g, 40g	\$35	KIT
PF-DLSIHC08/4G	PuriFlash Dry Load Column, HC Silica, 3g, 7ml, 12 x 65mm	\$95	20/PK
PF-DLSIHC08/12G	PuriFlash Dry Load Column, HC Silica, 11g, 30ml, 21 x 80mm	\$140	20/PK
PF-DLSIHC08/25G	PuriFlash Dry Load Column, HC Silica, 20g, 45ml, 21 x 130mm	\$157	20/PK
PF-DLSIHC08/40G	PuriFlash Dry Load Column, HC Silica, 32g, 75ml, 27 x 130mm	\$189	20/PK
PF-DLSIHC08/K	PuriFlash Dry Load Column Kit, HC Silica, 4 each: 3g, 11g, 20g, 32g	\$125	KIT
PF-DLCET08/12G	PuriFlash Dry Load Column, Celite, 5g, 30ml, 21 x 80mm	\$140	20/PK
PF-DLCET08/25G	PuriFlash Dry Load Column, Celite, 10g, 45ml, 21 x 130mm	\$180	20/PK
PF-DLCET08/40G	PuriFlash Dry Load Column, Celite, 15g, 75ml, 27 x 130mm	\$210	20/PK
PF-CLDET08/4G	PuriFlash Dry Load Column, Celite, 2g, 7ml, 12 x 65mm	\$95	20/PK
PF-DLP/4G	Top Frit Plunger for 12mm ID columns	\$10	EA
PF-DLP/25G	Top Frit Plunger for 21mm ID columns (12 & 25g)	\$10	EA
PF-DLP/40G	Top Frit Plunger for 27mm ID columns	\$10	EA
PF-DLP/80G	Top Frit Plunger for 31mm ID columns	\$20	EA
PF-DLP/120G	Top Frit Plunger for 36mm ID columns	\$25	EA
PF-DLP/200G	Top frit Plunger for 60mm ID columns (200 & 300g)	\$30	EA
FQ5191	Frits for PuriFlash Dry Load Column, 4g	\$44	20/PK
FQ5461	Frits for PuriFlash Dry Load Column, 12g	\$84	20/PK
FQ5471	Frits for PuriFlash Dry Load Column, 25g	\$84	20/PK
FQ5921	Frits for PuriFlash Dry Load Column, 40g	\$106	20/PK
JO0261	Frits for PuriFlash Dry Load Column, 80g	\$71	10/PK
JO0271	Frits for PuriFlash Dry Load Column, 120g	\$106	10/PK
JO0281	Frits for PuriFlash Dry Load Column, 200g	\$100	5/PK
JO0291	Frits for PuriFlash Dry Load Column, 300g	\$100	5/PK