

CONVENIENCE

Direct Adaptive HPLC Column Connectors



Direct Adaptive (Di-Ad™) HPLC Column Connectors were designed for Speed and Convenience when changing between different column brands on your instrument. They are very easy to use; one can quickly remove and replace columns without tools and the need to cut ferrules and adjust tubing that can ultimately and undesirably affect system volume. The spring-loaded, Di-Ad™ pilot length automatically adjusts to become compatible with any column 10-32 threaded, end-fitting. With a pressure rating up to 10,000 psi.

Using the Direct Adaptive Fittings can be considered a "Best Practice" in HPLC as it is an error-free process that can reduce Peak Broadening from poor connections. The time savings in easily switching columns is a great advantage and the peace of mind that your column and instrument connection will not affect your chromatography.

Di-Ad™ Fittings

Direct Adaptive HPLC
Column Connectors

- » Reduce Peak Tailing from Poor Connections
- » Pressure Rating Stainless Steel Tubing: 10,000 psi
- » Connection with all Brands of Columns is Easy
- » Universal, Self Adapting Pilot Length



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

Different internal depths can cause band broadening if tubing is not adjusted properly, Di-Ad™ eliminates this.


MANUFACTURED BY:

MICROSOLV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

customers@mtc-usa.com
www.mtc-usa.com



Convenient, Direct Adaptive HPLC Column Connectors

MicroSolv Brand HPLC Column Connectors made with Stainless Steel, laser cut tubing

Product Number	Description	Length	Single or Double	Internal Diameter Size
49905-10-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Double End Fittings	0.005"
49905-10-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Single End Fitting	0.005"
49905-20-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	20cm	Double End Fittings	0.005"
49905-20-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	20cm	Single End Fitting	0.005"
49905-30-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Double End Fittings	0.005"
49905-30-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Single End Fitting	0.005"
49910-10-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Double End Fittings	0.010"
49910-10-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Single End Fitting	0.010"
49910-20-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	20cm	Double End Fittings	0.010"
49910-20-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	20cm	Single End Fitting	0.010"
49910-30-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Double End Fittings	0.010"
49910-30-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Single End Fitting	0.010"
49930-10-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Double End Fittings	0.030"
49930-10-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	10cm	Single End Fitting	0.030"
49930-30-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Double End Fittings	0.030"
49930-30-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	30cm	Single End Fitting	0.030"
49930-50-DD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	50cm	Double End Fittings	0.030"
49930-50-SD	Di-Ad™ Direct Adaptive HPLC Column Connector. 1 EA.	50cm	Single End Fitting	0.030"

Other Sizes and Formats are Available Upon Request.

MANUFACTURED BY: