



1	5'-CMP	2	5'-AMP	3	5'-UMP
	100 ppm		100 ppm		100 ppm
4	5'-IMP	5	5'-GMP		
	100 ppm		100 ppm		

HPLC Conditions

Column	CAPCELL PAK NH ₂ UG80 S5 : 4.6 mm i.d.×250 mm
Mobile phase	50mmol/L (NH ₄) ₂ HPO ₄ (pH3.0)
Flow rate	1.0 mL/min
Temperature	40 °C
Detect	UV 254 nm
Inj. vol.	10 µL