



# TECHNICAL SPECIFICATIONS

**SKU:** ABALONASE™

**Form:** Clear aqueous solution.

**Volumes:** 10mL, 25mL, 50mL and 100mL.

**β-Glucuronidase activity:** > 50,000 units/mL (Phenolphthalein-β-D-glucuronide as substrate), > 20,000,000 units/g protein.

**Sulfatase activity:** none.

**Storage Temperature:** +4 to +8°C.

---

**SKU:** ABALONASE™ +

**Form:** Clear aqueous solution.

**Volumes:** 10mL, 25mL, 50mL and 100mL.

**β-Glucuronidase activity:** > 50,000 units/mL (Phenolphthalein-β-D-glucuronide as substrate), > 20,000,000 units/g protein.

**Sulfatase activity:** > 400 U/mL.

**Storage Temperature:** +4 to +8°C